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H-Cr Series

Reagent Strips For Urinalysis User's Guide

Rev:12/2018

Product Name: Reagent Strips For Urinalysis
Specification: 100strips/bottle, 50 strips/bottle,25 strips/bottle


Products Type	Test Item
H13-Cr	Urobilinogen, Bilirubin, Ketone (acetoacetic acid), Creatinine,Blood, Protein, Microalbumin, Nitrite, Leukocytes, Glucose, Specific Gravity, pH and Ascorbic Acid
H12-Cr	Urobilinogen, Bilirubin, Ketone (acetoacetic acid), Creatinine,Blood, Protein, Microalbumin, Nitrite, Leukocytes, Glucose, Specific Gravity, pH
H2-Cr	Microalbumin and Creatinine

The detection of albumin possesses a significance for a variety of diseases in early detection.
(1) Value for hypertension. Hypertensive patients' urinary albumin excretion rate is significantly higher than normal urine, elevated urinary albumin is an important indicator of cardiovascular disease progression.
(2) Micro-albumin can predict the development and early diagnosis of diabetic nephropathy. The earliest clinical signal of diabetic nephropathy is the presence of micro-proteinuria,therefore, measurement of miroalbumin will help patients with diabetes to take measures in advance to protect renal function.
(3) Microalbumin detection is sensitive indicator of diabetic microvascular disease.

Products Type	Applicable Instruments
H13-Cr	DIRUI H - 50/100/300/500 Urine Analyzer
H12-Cr	DIRUI H - 50/100/300/500 Urine Analyzer
H2-Cr	DIRUI H - 50/100/300/500 Urine Analyzer

Sample Requirement
Collect fresh urine in a clean and dry container. Don't centrifuge the urine. Mix the sample well before taking the test. The urine test must be taken within an hour. All sample must always be taken and kept under sanitary conditions.
Note: The preservatives will not prevent the deterioration of ketones, bilirubin or urobilinogen. The growth of bacteria in the long-term storage sample may affect the test results of glucose, pH, nitrite and blood.

[Test Method]
Room temperature for test:(25±5)℃
Visual Reading
1.Immerse the reagent area of the strip in the urine specimen and remove quickly.
2.Run the edge of the strip against the rim of the container to remove the excess urine.
3.Hold the strip horizontally and compare the result on the strip with the colour chart on the bottle label . Make note of the result. For a semi-quantitative result, read the result according to the time specified on the colour chart. For a qualitative results, the strip should be read between 1-2 minutes after dipping. If a positive result is obtained, repeat the test and compare with the colour chart at the time specified. Colour changes beyond 2 minutes are of no diagnostic value.

Instrumental Testing
Follow the directions given in appropriate instrument-operating manual.

[Storage Condition and Validity]
Storage Condition: The strips must be stored in dry place at the temperature between 2℃-30℃. Avoid direct sunlight. As to protect the activity of the reagent, must eliminate the humidity, sunlight and heat from the ambience environment.
Validity: Stored in dry place and avoid sunlight, with the temperature between 2℃-30℃, the sealed validity is 2 years; After opened seal, cover the tap closely, stored in dry, avoid sunlight with the temperature between 2℃-30℃, the validity is 1 month.

Reaction Principles
Glucose: The glucose oxidized by glucose oxidase catalyzes the formation of glucuronic acid and peroxide hydrogen. Peroxide hydrogen releases neo-ecotypes oxide [O] under the function of peroxidase. [O] oxidizes iodide potassium , which makes the colour changes.
Bilirubin:The direct bilirubin and dichlorobenzene diazonium produce azo dyes in a strongly acid medium.
Ketone: The acetoacetate and sodium nitroprusside cause reaction in alkaline medium, which produces violet colour.
Specific Gravity: Electrolyte (M^+X^-) in the form of salt in urine reacts with poly methyl vinyl ether and maleic acid(-COOH),which are weak acid ionic exchanger. The reaction produces hydrogenous ionogen, which reacts with pH indicator that causes the colour change.
Blood: Hemoglobin acts as peroxidase. It can cause peroxidase release neo-ecotypes oxide(O). (O) oxidizes the indicator and make the colour change subsequently.
pH: The method of pH indicator is applied.
Protein: This is based on the protein-error-of-indicator principle. Anion in the specific pH indicator attracted by cation on protein molecule makes the indicator further ionized, which changes its colour.
Urobilinogen: Urobilinogen and diazonium produce pink azo dyes under the function of strong acid medium.
Nitrite: Nitrite in the urine and aromatic amino sulphanilamide are diazotized to form a diazonium compound. The diazonium compound reacting with tetrahydro benzo(h) quinolin 3-phenol causes the colour change.
Leukocytes: Granulocyte leukocytes in urine contains esterases that catalyze the hydrolysis of the pyrrole amino acid ester to liberate 3-hydroxy 5-pheny pyrrole. This pyrrole reacting with diazonium forms a purple colour.
Microalbumin: sulfone phthalein dye have high sensitivity to microalbumin by the method of protein error.
Ascorbic Acid: Ascorbic acid, with 1,2-dihydroxy alkenes, under the alkaline condition, deoxidizes the blue 2,6-dichloroindophenolate into colourless N-(p-pheno)-2,6- dichloro-P-amine phenol.
Creatinine:ceratinine and 3,5-dinitro-benzoic acid react in strong alkaline condition,and generate colored complex.

Result Explanation
Glucose
The test is for specificity of glucose. There is no false positive result occurred in reagent strip, caused by any substance in urine. When the ascorbic acid concentration≥2.8mmol/L or acetoacetic acid concentration≥1.0mmol/L, the sample of glucose concentration is 3～7mmol/L may occur false negative result.
Bilirubin
Normally, even the most sensitive method can't detect bilirubin in urine. It is abnormal to have little bilirubin in urine, which requires further inspection. Medicines that dyes urine red and anything that shows red itself in an acid medium e.g., phenazopyridine may affect the test result. High concentration of the ascorbic acid may cause false negative result.
Ketone
The reagent strip reacts with acetoacetic acid and acetone in urine. It doesn't do with β-hydro butyric acid. Normal urine specimens usually conduct negative results in the test. False positive results may occur in highly pigmented urine or those containing a large amount of levodopa metabolites.
Specific Gravity
The reagent strip for Specific Gravity allows the urine specimens specific gravity between 1.000 and 1.030.In general, the mean error between the results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.005 may be added to readings from urines with pH equal or greater than 6.5.Urine reading instrument can automatically make these adjustments in strip-readings. The urine nonionic constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated specific gravity readings may occur in the presence of moderate quantities of protein(1.75g/L).
Blood
'Trace' reaction may vary among the patients. Clinical judgments are required for individual cases. The presence of green spots (intact erythrocytes) or green color (haemoglobin/myoglobin) on the reagent area within 60 seconds indicates for further diagnostic check. Blood is often found in the urine of the menstruating females. Haemoglobin 150 μ g/L-620 μ g/L is approximately equivalent to 5-15 cells/ μ L intact erythrocytes.The reagent strip is highly sensitive to haemolobin and thus can be used as a supplementary to the microscopic examination. The sensitivity of the strip might be reduced in urine with a large amount of specific gravity. The strips are equally sensitive to myglobin as to haemoglobin. Certain oxidizing contaminants, such as hypochlorite, may lead to false positive results. Microbial peroxidase associated with urinary tract infection may also produce a false positive result.Ascorbic acid less than 5.0mmol/L in urine may not influence the result of the test.
pH
The strip tests for pH values are generally in the range of 5.0-8.5 visually and 5.0-9.0 instrumentally.
Protein and microalbumin
Adult normal urinary albumin concentration is less than 20mg / L. The results within 20 ~ 200mg / L indicate clinical micro-albuminuria; the results greater than 200mg / L indicate clinical proteinuria. When albumin excretion rate within 30 ~ 300mg/24 hours, it is shown as the precursor of micro-albuminuria. Urinary albumin excretion may be increased due to a temporary exercise, urinary tract infection, acute febrile illness.
The Microalbumin test can accurately and specifically determine albumin. A strip result of 15 mg/dL is indicative of clinical albuminuria. The test is not affected by other proteins at concentrations at least nine times greater than the excretion rate considered to be abnormal.
Albumin reagent strip is used for testing urinary volume albumin, and reagent strip has low sensitivity to mucinogen and globulin that it can be tested when more than 60mg/dL. Visible haematuria (≥5mg/dL) may also occur false positive result.

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The following possibilities occurs when micro albumin emerges in urine:
(1) Occasional urinary micro-albumin may symbolize physiological proteinuria, such as functional proteinuria, or orthostatic proteinuria caused by diet, exercise, mental stress
(2) The clinical significance is very important if the consecutive urine microalbumin, or micro-albumin and glucose, or micro-albumin and blood emerge in urine at the same time.
The Protein test pad is not specific for a particular protein，and proteins other than albumin can cause a positive response. The test is less sensitive to mucoproteins and globulins, which are generally detected at levels of 60 mg/dL or higher.
A visibly bloody urine (≥5 mg/dL) may cause falsely elevated results.

Urobilinogen

The reagent strips can detect urobilinogen in low amount as 3 μmol/L (approximately 0.2 Ehrlich unit/dL) in urine. A result of 33 μ mol/L in urine indicate the critical value, representing the transition from normal to abnormal, which requires further check on patients and specimens. The negative results are not final to determine the absence of urobilinogen.

Nitrite

Gram-negative bacteria in urine converts nitrate(derived from foods)into nitrite.The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink color development should be taken as positive result. The degrees of colour development the numbers of bacteria are not in direct proportion. The negative result doesn't mean the existence of bacteria in a large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite. (2)when urine has not remained in the bladder long enough(four hours up)to let the nitrate covert into nitrite.(3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 2.8mmol/L ascorbic acid or less won't interfere the test result.

Leukocytes

Test area react with esterase in leucocytes (granulocytic leukocytes). Normal urine specimens generally yield negative result; positive results (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance; however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random specimens from females due to contamination of the specimen by vaginal discharge. Elevated glucose concentrations (160mmol/L) or high specific gravity may cause decreased test results.

Ascorbic Acid

The test area can detect the ascorbic acid in urine. Through the ascorbic acid detection, we will know the level of ascorbic acid in the body and the effect degree that the ascorbic acid bring to the test on glucose, bilirubin, blood and nitrite. It will reduce the sensitivity when the oxidant (such as potassium permanganate, hypochlorite) in the urine.

Creatinine

Normal adult urinary creatinine concentration: 0.6 ~ 2.0g/24 hours (strip test result: 4.4 ~ 17.7mmol / L). Random urine creatinine test results vary widely, which may be in the range of from 0.9mmol / L to 26.5mmol / L. High level happens to concentrated urine and morning urine (strip test results may be higher than 17.7mmol / L). Due to xcessive drinking and other polyuria to dilute urine, which will lead to a typical low concentration of urine (test results ≤ 50mg / dL) Visible haematuria (≥ 50mg / dL) or cimetidine may cause creatinine false test results elevated.

Aulbumin—creatinine Ratio

Adult normal urine albumin level is less than 30mg albumin / g creatinine (3.4mg albumin / mmol creatinine), sample rate result is in the range of 30 ~ 300mg / g (3.4 ~ 33.9mg/mmol) can be defined as microalbuminuria, and proteinuria when > 300mg / g (> 33.9mg/mmol)

Result Table

Conventional uints and international units available

Test Item	Abb.	Print or Display Result			
		Conventional Unit		International Unit	
Microalbumin	MA	10 mg/L	80 mg/L	10 mg/L	80 mg/L
		30 mg/L	150 mg/L	30 mg/L	150 mg/L
Creatinine	Cr	10 mg/dL	200 mg/dL	0.9 mmol/L	17.7mmol/L
		50 mg/dL	300 mg/dL	4.4 mmol/L	26.5 mmol/L
		100 mg/dL		8.8 mmol/L	
Aulbumin-creatinine Ratio	A:C	<30 mg/g (Normal)		<3.4 mg/mmol(Normal)	
		30-300 mg/g (Abnormal)		3.4-33.9 mg/mmol(Abnormal)	
		>300 mg/g (Highly abnormal)		>33.9 mg/mmol(Highly abnormal)	

Remark: The shaded parts are abnormal results.

Test Method Limits

Like all the other laboratory tests, definitive diagnostic or therapeutic decisions should not be made or based on any single result or method.
Specific performance characteristics are based on clinical and analytical studies. In clinical sample , the sensitivity depends upon several factors; the variability of colour perception, the presence or absence of inhibitory factors typically found in urine, specific gravity, pH, and the lighting conditions when the product is read visually. Each colour block or instrumental display value represents a range of values. Because of sample and reading variability, sample with analyte concentrations that fall between nominal levels may give results at either level. Results at levels greater than the second positive level for the Protein, Glucose, Ketone, and Urobilinogen tests will usually be within one level of the true concentration. Exact agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical system of the instruments.

Points For Attention

1.In vitro diagnosis for professional only.
2. Reagent strips must be stored in original cylinder; do not take out strips except use them immediately; removing the desiccant is forbidden; tighten the cap immediately after take out strips;
3. The expired should not be used; Please make sure whether the strip is over the validity, and test it by quality control if the deterioration of strips occurs, which leads color shallower or deeper, or any question occurs with the unexpected results.
4. Water can not be used as negative control.
5. One time use only for each strip.
6. Please read the instruction carefully before the use.
7. Please refer to the lab biological treatment of harmful materials to dispose the waste strips.
8. Do not place them in refrigerator, however they should be stored in dry circumstance; avoid direct sunlight; do not touch the reactive area of strip; ambient moisture, light and heat must be prevented in order to protect the reactivity of the strips.

Test Range and Sensitivity of H-Cr Series Strips

Item	Sensitivity	Instrumental Test Range	Visual Test Range
Glucose(mmol/L)	2.8-5.6	Neg.-56	
Protein(g/L)	0.15-0.3	Neg.-3.0	Neg. -20.0
Microalbumin(mg/L)	30-80	10-150	10-150
Ketone(acetoacetic acid) (mmol/L)	0.5-1.0	Neg.-7.8	Neg.-16
Blood (Ery/ μ L)	5-15	Neg.-200	
Bilirubin (μ mol/L)	8.6-17	Neg.-103	
Nitrite(μ mol/L)	13-22	Neg.- Pos.	
Leukocytes (leuko/ μ L)	5-15	Neg.-500	
Urobilinogen (μ mol/L)	3.4-17	3.4-135	
Ascorbic acid(mmol/L)	0.6-1.4	0-5.7	
Creatinine(mmol/L)	——	0.9-26.5	
PH	——	5.0-9.0	5.0-8.5
Specific Gravity	——	1.005-1.030	1.000-1.030

REACTIVE INGREDIENTS

Microalbumin	sulfone phthalein buffer nonreactive ingredients	2.2%w/w 96.0%w/w 1.8%w/w
Protein	tetrabromphenol blue buffer nonreactive ingredients	0.1% w/w 97.4% w/w 2.5% w/w
Blood	diisopropylbenzene dihydro peroxide tetramethylbenzidine buffer nonreactive ingredients	26.0% w/w 1.5% w/w 35.3% w/w 37.2 %w/w
Glucose	glucose oxidase (microbial.123U) peroxidase(horseradish. 203U) potassium iodide buffer nonreactive ingredients	1.7% w/w 0.2 % w/w 0.1% w/w 71.8% w/w 26.2% w/w
Ketone	sodium nitroprusside nonreactive ingredients buffer	5.7% w/w 29.9% w/w 64.4% w/w
Leukocytes	pyrrole amino acid ester diazonium salt buffer nonreactive ingredients	4.3% w/w 0.4% w/w 92.6% w/w 2.7% w/w
Nitrite	p-arsanilicacid-N-(1-Naphthol) -ethylenediamine tetrahydroquinoline buffer nonreactive ingredients	1.3% w/w 0.9%w/w 89.6% w/w 8.2% w/w
Specific Gravity	bromthymol blue poly(methyl vinyl ether co maleic anhydride) sodium hydroxide	4.8%w/w 90.2%w/w 5.0% w/w
pH	bromcresol green bromthymol blue nonreactive ingredients	3.3% w/w 55.0% w/w 41.7% w/w
Bilirubin	2,4-dichloroaniline diazonium salt buffer nonreactive ingredients	0.6%w/w 57.3%w/w 42.1%w/w
Urobilinogen	fast B blue buffer nonreactive ingredients	0.2%w/w 98.0%w/w 1.8%w/w
Ascorbic acid	2,6-dichloroindophenolate hydrate buffer nonreactive ingredients	0.8%w/w 40.7%w/w 58.5%w/w
Creatinine	3,5-dinitro-benzoicaid buffer nonreactive ingredients	4.8%w/w 85.2%w/w 10 %w/w

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Rev: 12/2018

[Product Name]
Common Name: H13-Cr Reagent Strips for Urinalysis
Trade Name: Reagent Strips for Urinalysis
[Specification]
100 strips/bottle, 50 strips/bottle, 25 strips/bottle

[Intended Use]
H13-Cr Reagent Strips for Urinalysis is used for qualitative and semi-quantitative test of urobilinogen, bilirubin, ketone (acetacetic acid), blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid, microalbumin, and creatinine in urine.
This product is used for in vitro diagnosis by professionals.

[Testing Principle]
Glucose: The glucose oxidized by glucose oxidase catalyzes the formation of glucuronic acid and peroxide hydrogen. Peroxide hydrogen releases neo-ecotypes oxide [O] under the function of peroxidase. [O] oxidizes iodide potassium, which makes the color changes.
Bilirubin: The direct bilirubin and dichlorobenzene diazonium produce azo dyes in a strongly acid medium.
Ketone: The acetacetic acid and sodium nitroprusside cause reaction in alkaline medium, which produces violet compound.
Specific Gravity: Electrolyte (M+X-) in the form of salt in urine reacts with poly methyl vinyl ether and maleic acid (-COOH),which are weak acid ionic exchanger. The reaction produces hydrogenous ionogen, which reacts with pH indicator that causes the color change.
Blood: Hemoglobin acts as peroxidase. It can cause peroxidase release neo-ecotypes oxide (O). (O) oxidizes the indicator and make the color change subsequently.
pH: The method of pH indicator is applied.
Protein: This is based on the protein-error-of-indicator principle. Anion in the specific pH indicator attracted by cation on protein molecule makes the indicator further ionized, which changes its color.
Urobilinogen: Urobilinogen and diazonium produce red violet dyes under the function of strong acid medium.
Nitrite: Nitrite in the urine and aromatic amino sulphanilamide are diazotized to form a diazonium compound. The diazonium compound reacting with tetrahydro benzo(h) quinolin 3-phenol causes the red azo dye.
Leukocyte: Granulocyte leukocytes in urine contains esterases that catalyze the hydrolysis of the pyrrole amino acid ester to liberate 3-hydroxy 5-pheny pyrrole. This pyrrole reacting with diazonium forms a purple azo dye.
Microalbumin: Sulfone phthalein dye have high sensitivity to microalbumin by the method of protein error.
Ascorbic Acid: Ascorbic acid, with 1,2-dihydroxy alkenes, under the alkaline condition, deoxidizes the blue 2,6-dichloroindophenolate into colorless N-(p-phenol)-2,6- dichloro-P-amine phenol.
Creatinine: Creatinine and 3,5-dinitro-benzoic acid react in strong alkaline condition, and generate colored compound.


[Main Ingredients]
Leukocyte: 4.3% w/w pyrrole amino acid ester; 0.4% w/w diazonium salt; 92.6% w/w buffer; 2.7% w/w nonreactive ingredients.
Nitrite:1.3% w/w p-arsanilicacid-N-(1-Naphthol)-ethylenediamine; 0.9%w/w tetrahydroquinoline; 89.6% w/w buffer; 8.2% w/w nonreactive ingredients.
Protein: 0.1% w/w tetrabromphenol blue; 97.4% w/w buffer; 2.5% w/w nonreactive ingredients.
Blood: 26.0% w/w diisopropylbenzene dihydro peroxide; 1.5% w/w tetramethylbenzidine; 35.3% w/w buffer; 37.2 %w/w nonreactive ingredients.
Urobilinogen: 0.2%w/w fast B blue; 98.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Bilirubin: 0.6%w/w 2,4-dichloroaniline diazonium salt; 57.3%w/w buffer; 42.1%w/w nonreactive ingredients.
pH: 3.3% w/w methyl red; 55.0% w/w bromthymol blue; 41.7% w/w nonreactive ingredients.
Specific Gravity: 4.8%w/w bromthymol blue; 90.2%w/w poly(methyl vinyl ether co. maleic anhydride).
Glucose: 1.7% w/w glucose oxidase (microbial.123U); 0.2 % w/w peroxidase(horseradish. 203U); 0.1% w/w tetramethyl benzidine; 71.8% w/w buffer; 26.2% w/w nonreactive ingredients.
Ketone: 5.7% w/w sodium nitroprusside; 29.9% w/w buffer; 64.4% w/w nonreactive ingredients.
Microalbumin: 2.2%w/w sulfone phthalein dye; 96.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Ascorbic Acid: 0.8%w/w 2,6-dichloroindophenolate hydrate; 40.7%w/w buffer; 58.5%w/w nonreactive ingredients.
Creatinine: 4.8%w/w 3,5-dinitro-benzoicaid; 85.2%w/w buffer; 10 %w/w nonreactive ingredients.

[Storage Conditions and Shelf Life]
Storage Condition: The strips must be stored in dry place at the temperature between 2℃~30℃, and do not keep in the refrigerator; Avoid direct sunlight and do not touch the reaction area of the strip; As to protect the activity of the reagent, must eliminate the humidity, sunlight and heat from the ambience environment.
Shelf Life: Stored in dry place and avoid sunlight, with the temperature between 2℃~30℃, the sealed shelf Life is 2 years; After opened, cover the tap closely, stored in dry, avoid sunlight with the temperature between 2℃~30℃, the shelf Life is 1 month.

[Applicable Instrument]
DIRUI H-800 Automatic Urine Analyzer

[Sample Requirements]
Collect fresh urine in a clean and dry container. Don't centrifuge the urine. Mix the sample well before the test. The urine test must be taken within two hours All samples must always be taken and kept under sanitary conditions.
Note: The preservatives will not prevent the deterioration of ketones, bilirubin or urobilinogen. The growth of bacteria in the long-term storage sample may affect the test results of glucose, pH, nitrite and blood.

[Test Method]
Visual Reading:
Room temperature for test:(25±5)℃
1. Immerse the reagent area of the strip in urine sample and take it out quickly and immediately.
2. Run the edge of the strip against the rim of the container to remove the excess urine.
3. Hold the strip up horizontally and compare the result on the strip with the color chart on the bottle label closely. Make a note of the result. For a semi-quantitative result, take the result according to the time specified on the color chart. For a qualitative result, the strip should be read within 1-2 minutes after dipping. If a positive result is obtained, repeat the test and compare with the color chart at the time specified. Color changes beyond 2 minutes are of no diagnostic value.

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[Result Explanation]
Glucose
The test is for specificity of glucose. There is no false positive result occurred in reagent strip, caused by any other substance in urine. When the ascorbic acid concentration≥2.8mmol/L or acetoacetic acid concentration≥1.0mmol/L, the sample of glucose concentration is 3～7mmol/L may occur false negative result.
Bilirubin
Normally, even the most sensitive method can't detect bilirubin in urine. It is abnormal to have little bilirubin in urine, which requires further inspection. Medicines that dyes urine red and anything that shows red itself in an acid medium e.g., phenazopyridine may affect the test result. High concentration of the ascorbic acid may cause false negative result.
Ketone
The reagent strip acetone reacts with acetoacetic acid of urine. It doesn't do with β-hydro butyric acid. Normal urine sample usually shows negative results in the test. False positive results may occur in highly pigmented urine or those containing a large amount of levodopa metabolites.
Specific Gravity
The reagent strip for Specific Gravity allows the urine sample specific gravity between 1.000 and 1.030. In general, the mean error between the results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.005 may be added to readings from urines with pH equal or greater than 6.5. Urine reading instrument can automatically make these adjustments in strip-readings. The urine nonionic constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated

DIRUI® DIRUI INDUSTRIAL CO.,LTD.
IVD CE EC REP EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands

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specific gravity readings may occur in the presence of moderate quantities of protein(1g/L-1.75g/L).

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Nitrite

Gram-negative bacteria in urine converts nitrate(derived from foods)into nitrite. The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink color development should be taken as
itive result. The degrees of color development the numbers of bacteria are not in direct proportion. The negative result doesn't mean the existence of bacteria in a large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite. (2)when urine has not remained in the bladder long enough(four hours up)to let the nitrate covert into nitrite. (3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 1.4mmol/L ascorbic acid or less won't interfere the result.

Leukocyte

Test area reacts with esterase in leucocytes (granulocytic leukocytes). Normal urine sample generally yield negative result; positive results (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance; however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random sample from females due to contamination of the sample by vaginal discharge. Elevated glucose concentrations (160mmol/L) or high specific gravity may cause decreased test results.

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Albumin-Creatinine Ratio

Urine albumin content of adult is less than 30 mg albumin/g creatinine (3.4 albumin/mmol creatinine). A ratio of 30-300mg/g (3.4-33.9mg/mmol) indicates microalbuminuria. A ratio of more than 300mg/g (>33.9mg/mmol) indicates albuminuria.

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[Performance Index]

Test Range and Sensitivity of H13-Cr Series Strips

Item	Sensitivity	Instrument Test Range	Visual Test Range
Glucose(mmol/L)	2.8-5.6	Neg.-56	
Protein(g/L)	0.15-0.3	Neg.-3.0	Neg.-20.0
Microalbumin(mg/L)	30-80	10-150	10-150
Ketone (Acetacetic Acid)(mmol/L)	0.5-1.0	Neg.-7.8	Neg.-16
Blood(Ery/uL)	5-15	Neg.-200	
Bilirubin(umol/L)	8.6-17	Neg.-103	
Nitrite(umol/L)	13-22	Neg.-Pos.	
Leukocyte(Leuko/uL)	5-15	Neg.-500	
Urobilinogen(umol/L)	3.4-17	3.4-135	
Ascorbic Acid(mmol/L)	0.6-1.4	0-5.7	
Creatinine(mmol/L)	—	0.9-26.5	
pH	—	5.0-9.0	5.0-8.5
Specific Gravity	—	1.005-1.030	1.000-1.030

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H14-Ca Urinalysis Reagent Strip User's Guide

H-800 Automatic Urine Analyzer Suited Strips

Rev:11/2019

[Product Name]
Common Name: H14-Ca Urinalysis Reagent Strips
Trade Name: Urinalysis Reagent Strips
[Specification]
100 strips/bottle, 50 strips/bottle, 25 strips/bottle

[Intended Use]
H14-Ca Urinalysis Reagent Strips is used for qualitative and semi-quantitative test of urobilinogen, bilirubin, ketone (acetacetic acid), blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid, microalbumin, creatinine and Ca in urine.
This product is used for in vitro diagnosis by professionals.

[Testing Principle]
Glucose: The glucose oxidized by glucose oxidase catalyzes the formation of glucuronic acid and peroxide hydrogen. Peroxide hydrogen releases neo-ecotypes oxide [O] under the function of peroxidase. [O] oxidizes iodide potassium, which makes the color changes.
Bilirubin: The direct bilirubin and dichlorobenzene diazonium produce azo dyes in a strongly acid medium.
Ketone: The acetacetic acid and sodium nitroprusside cause reaction in alkaline medium, which produces violet compound.
Specific Gravity: Electrolyte (M+X-) in the form of salt in urine reacts with poly methyl vinyl ether and maleic acid (-COOH),which are weak acid ionic exchanger. The reaction produces hydrogenous ionogen, which reacts with pH indicator that causes the color change.
Blood: Hemoglobin acts as peroxidase. It can cause peroxidase release neo-ecotypes oxide (O). (O) oxidizes the indicator and make the color change subsequently.
pH: The method of pH indicator is applied.
Protein: This is based on the protein-error-of-indicator principle. Anion in the specific pH indicator attracted by cation on protein molecule makes the indicator further ionized, which changes its color.
Urobilinogen: Urobilinogen and diazonium produce red violet dyes under the function of strong acid medium.
Nitrite: Nitrite in the urine and aromatic amino sulphanilamide are diazotized to form a diazonium compound. The diazonium compound reacting with tetrahydro benzo(h) quinolin 3-phenol causes the red azo dye.
Leukocyte: Granulocyte leukocytes in urine contains esterases that catalyze the hydrolysis of the pyrrole amino acid ester to liberate 3-hydroxy 5-pheny pyrrole. This pyrrole reacting with diazonium forms a purple azo dye.
Microalbumin: Sulfone phthalein dye have high sensitivity to microalbumin by the method of protein error.
Ascorbic Acid: Ascorbic acid, with 1,2-dihydroxy alkenes, under the alkaline condition, deoxidizes the blue 2,6-dichloroindophenolate into colorless N-(p-phenol)-2,6- dichloro-P-amine phenol.
Creatinine: Creatinine and 3,5-dinitro-benzoic acid react in strong alkaline condition, and generate colored compound.
Ca: Ca reacts with o-cresolphthalein complexone produce purple red, the color shade is proportional to Ca concentration.


[Main Ingredient]
Leukocyte: 4.3% w/w pyrrole amino acid ester; 0.4% w/w diazonium salt; 92.6% w/w buffer; 2.7% w/w nonreactive ingredients.
Nitrite:1.3% w/w p-arsanilicacid-N-(1-Naphthol)-ethylenediamine; 0.9%w/w tetrahydroquinoline; 89.6% w/w buffer; 8.2% w/w nonreactive ingredients.
Protein: 0.1% w/w tetrabromphenol blue; 97.4% w/w buffer; 2.5% w/w nonreactive ingredients.
Blood: 26.0% w/w diisopropylbenzene dihydro peroxide; 1.5% w/w tetramethylbenzidine; 35.3% w/w buffer; 37.2 %w/w nonreactive ingredients.
Urobilinogen: 0.2%w/w fast B blue; 98.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Bilirubin: 0.6%w/w 2,4-dichloroaniline diazonium salt; 57.3%w/w buffer; 42.1%w/w nonreactive ingredients.
pH: 3.3% w/w methyl red ; 55.0% w/w bromthymol blue; 41.7% w/w nonreactive ingredients.
Specific Gravity: 4.8%w/w bromthymol blue; 95.2%w/w poly(methyl vinyl ether co. maleic anhydride).
Glucose: 1.7% w/w glucose oxidase (microbial.123U); 0.2 % w/w peroxidase(horseradish. 203U); 0.1% w/w tetramethyl benzidine; 71.8% w/w buffer; 26.2% w/w nonreactive ingredients.
Ketone: 5.7% w/w sodium nitroprusside; 29.9% w/w buffer; 64.4% w/w nonreactive ingredients.
Microalbumin: 2.2%w/w sulfone phthalein dye; 96.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Ascorbic Acid: 0.8%w/w 2,6-dichloroindophenolate hydrate; 40.7%w/w buffer; 58.5%w/w nonreactive ingredients.
Creatinine: 4.8%w/w 3,5-dinitro-benzoicaid; 85.2%w/w buffer; 10 %w/w nonreactive ingredients.
Ca: 2.5%w/w o-cresolphthalein complexon; 87.5%w/w buffer; 10 %w/w nonreactive

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[Applicable Instrument]
DIRUI H-800 Automatic Urine Analyzer

[Sample Requirement]
Collect fresh urine in a clean and dry container. Don't centrifuge the urine. Mix the sample well before the test. The urine test must be taken within two hour. All samples must always be taken and kept under sanitary conditions.
Note: The preservatives will not prevent the deterioration of ketones, bilirubin or urobilinogen. The growth of bacteria in the long-term storage sample may affect the test results of glucose, pH, nitrite and blood.

[Test Method]
Room temperature for test:(25±5)℃
Visual Reading:
1. Immerse the reagent area of the strip in urine sample and take it out quickly and immediately.
2. Run the edge of the strip against the rim of the container to remove the excess urine.
3. Hold the strip up horizontally and compare the result on the strip with the color chart on the bottle label closely. Make note of the result. For a semi-quantitative result, take the result according to the time specified on the color chart. For a qualitative result, the strip should be read within 1-2 minute after dipping. If a positive result is obtained, repeat the test and compare with the color chart at the time specified. Color changes beyond 2 minutes are of no diagnostic value.

Instrumental Testing:
Follow the directions given in appropriate instrument-operating manual.

[Result Explanation]
Glucose
The test is for specificity of glucose. There is no false positive result occurred in reagent strip, caused by any other substance in urine. When the ascorbic acid concentration≥2.8mmol/L or acetoacetic acid concentration≥1.0mmol/L, the sample of glucose concentration is 3～7mmol/L may occur false negative result.
Bilirubin
Normally, even the most sensitive method can't detect bilirubin in urine. It is abnormal to have little bilirubin in urine, which requires further inspection. Medicines that dyes urine red and anything that shows red itself in an acid medium e.g., phenazopyridine may affect the test result. High concentration of the ascorbic acid may cause false negative result.
Ketone
The reagent strip acetone reacts with acetoacetic acid of urine. It doesn't do with β-hydro butyric acid. Normal urine sample usually shows negative results in the test. False positive results may occur in highly pigmented urine or those containing a large amount of levodopa metabolites.
Specific Gravity
The reagent strip for Specific Gravity allows the urine sample specific gravity between 1.000 and 1.030. In general, the mean error between the results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.005 may be added to readings from urines with pH equal or greater than 6.5. Urine reading instrument can automatically make these adjustments in strip-readings. The urine nonionic constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated specific gravity readings may occur in the presence of moderate quantities of protein(1g/L-7.5g/L).

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Blood
'Trace' reaction may vary among the patients. Clinical judgments are required for individual cases. The presence of green spots (intact erythrocytes) or green color (haemoglobin/myoglobin) on the reagent area within 60 seconds indicates for further diagnostic check. Blood is often found in the urine of the menstruating females. Haemoglobin 150μg/L-620μg/L is approximately equivalent to 5-15 cells/μL intact erythrocytes. The reagent strip is highly sensitive to haemolobin and thus can be used as a supplementary to the microscopic examination. The sensitivity of the strip might be reduced in urine with a large amount of specific gravity. The strips are equally sensitive to myglobin as to haemoglobin. Certain oxidizing contaminants, such as hypochlorite, may lead to false positive results. Microbial peroxidase associated with urinary tract infection may also produce a false positive result. Vitamine in urine may not influence the result of the test.

pH
The strip tests for pH values are generally in the range of 5.0-8.5 visually and 5.0-9.0 instrumentally.

Urobilinogen
The reagent strips can detect urobilinogen in low amount as 3μmol/L (approximately 0.2 Ehrlich unit/dL) in urine. A result of 33μmol/L in urine indicate the critical value, representing the transition from normal to abnormal, which requires further check on patients and samples. The negative results are not final to determine the absence of urobilinogen.

Nitrite
Gram-negative bacteria in urine converts nitrate(derived from foods)into nitrite. The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink color development should be taken as positive result. The degrees of color development the numbers of bacteria are not in direct proportion. The negative result doesn't mean the existence of bacteria in a large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite. (2)when urine has not remained in the bladder long enough(four hours up)to let the nitrate covert into nitrite. (3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 1.4mmol/L ascorbic acid or less won't interfere the result.

Leukocyte
Test area react with esterase in leucocytes (granulocytic leukocytes). Normal urine sample generally yield negative result; positive results (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance; however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random sample from females due to contamination of the sample by vaginal discharge. Elevated glucose concentrations (160mmol/L) or high specific gravity may cause decreased test results.

Ascorbic Acid
The test area can detect the ascorbic acid in urine. Through the ascorbic acid detection, we will know the level of ascorbic acid in the body and the effect degree that the ascorbic acid bring to the test on glucose, bilirubin, blood and nitrite. It will reduce the sensitivity when the oxidant (such as potassium permanganate, hypochlorite)present in the urine.

Protein and Microalbumin
Albumin concentration in normal adult urine is less than 20 mg/L. A result of 20-200 mg/L indicates clinical microalbumin. A result of higher than 200 mg/L indicates clinical albuminuria. A albumin discharge rate of 30-300 mg/24 hour indicates a warning of urine microalbumin symptom. A albumin discharge rate of higher than 300 mg/24 hour indicates a warning of urine albumin symptom. Urine albumin discharge may be elevated by temporary sport, urinary tract infection or acute fever disease.
Microalbumin strip is intended for detecting urine microalbumin, and is sensitive to microalbumin. Its sensitivity for other protein is 9 times lower than microalbumin. The test is less sensitive to mucoproteins and globulins, which are generally detected at levels of 60 mg/dL or higher. A visibly bloody urine (≥5 mg/dL) may cause falsely elevated results.
Microalbumin in urine may lead to:
(1)Occasionally presented urine microalbumin may be functional albuminuria or postural proteinuria caused by physiological albuminuria, such as diet, sport, mental stress, etc.
(2)Continuously presented urine microalbumin, microalbumin and glucose, or microalbumin and blood are with great clinical significance.

Creatinine
Normal urine creatinine concentration of adult is 0.6-2.0g/24h (about 4.4-17.7 mmol/L by test strip). High deviation exists on the test results of random urine sample creatinine, with the range from 0.9mmol/L to 26.5mmol/L. concentrated and morning urine have higher level of creatinine (maybe more than 17.7mmol/L tested by test strip). The concentration of the component required to test can be diluted due to over-ingestion of water or other substance, will result in typical low concentration urine (the tested result ≤50mg/dL). Visible hematuria(≥5mg/dL) or Cimitidine may result in creatinine testing result spurious rising.

Albumin-Creatinine Ratio
Urine albumin content of adult is less than 30 mg albumin/g creatinine (3.4 albumin/mmol creatinine). A ratio of 30-300mg/g (3.4-33.9mg/mmol) indicates microalbuminuria. A ratio of more than 300mg/g (＞33.9mg/mmol) indicates albuminuria.

Result Table
Routine system of units and international system of units can be selected in urine analyzer test.

Test Item	Ab.	Print or Display Result			
		Routine System of Units		International System of Units	
Microalbumin	MA	10 mg/L	80 mg/L	10 mg/L	80 mg/L
		30 mg/L	150 mg/L	30 mg/L	150 mg/L
Creatinine	Cr	10 mg/dL	200 mg/dL	0.9 mmol/L	17.7 mmol/L
		50 mg/dL	300 mg/dL	4.4 mmol/L	26.5 mmol/L
		100 mg/dL		8.8 mmol/L	
Albumin-Creatinine Ratio	A: C	＜30 mg/g(normal)		＜3.4 mg/mmol(normal)	
		30-300 mg/g(abnormal)		3.4-33.9 mg/mmol(abnormal)	
		＞300 mg/g(highly abnormal)		＞33.9 mg/mmol(highly abnormal)	

Remark: shaded area indicates abnormal result.
Ca:A great deal of Mg(more than 10mmol/L) will seduce urine Ca rising.

[Limitation of Test Method]
Like all lab tests, diagnosis result and therapeutic schedule can not be made according to any single diagnosis method.
Reagent strips' application is based on clinical analysis. In clinical sample, the sensitivity depends upon several factors: the variability of color perception, specific gravity, pH, and the lighting conditions change when the product is read visually. Each colour block or instrumental display value represents a range of values. Because of sample and reading variability, sample with analyte concentrations that fall between nominal levels may give results at either level. Results at levels greater than the second positive level for the Protein, Glucose, Ketone, and Urobilinogen tests will usually be within one level of the true concentration. Exact agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical system of the instruments.

[Performance Index]
Test Range and Sensitivity of H14-Ca Series Strips

Item	Sensitivity	Instrument Test Range	Visual Test Range
Glucose(mmol/L)	2.8-5.6	Neg.-56	
Protein(g/L)	0.15-0.3	Neg.-3.0	Neg.-20.0
Microalbumin(mg/L)	30-80	10-150	10-150
Ketone (Acetacetic Acid)(mmol/L)	0.5-1.0	Neg.-7.8	Neg.-16
Blood(Ery/uL)	5-15	Neg.-200	
Bilirubin(umol/L)	8.6-17	Neg.-103	
Nitrite(umol/L)	13-22	Neg.-Pos.	
Leukocyte(Leuko/uL)	5-15	Neg.-500	
Urobilinogen(umol/L)	3.4-17	3.4-135	
Ascorbic Acid(mmol/L)	0.6-1.4	0-5.7	
Creatinine(mmol/L)	—	0.9-26.5	
Ca(mmol/L)	2.0-2.5	1.0-10	
pH	—	5.0-9.0	5.0-8.5
Specific Gravity	—	1.005-1.030	1.000-1.030

[Points for Attention]
1.In vitro diagnosis for professional only.
2.Reagent strips must be stored in original bottle, do not take out strips except use them immediately, removing the desiccant is forbidden, tighten the cap immediately after taking out strips.
3.The expired should not be used, please make sure whether the strip is over the validity, and test it by quality control if the deterioration of strips occurs, which leads color shallower or deeper, or any question occurs with the unexpected results.
4.Water can not be used as negative control.
5.One time use only for each strip.
6.Please read the Guide carefully before use.
7.Please refer to the lab biological treatment of harmful materials to dispose the waste strips.
8.Do not place them in refrigerator, however they should be stored in dry circumstance, avoid direct sunlight, do not touch the reactive area of strip, ambient moisture, light and heat must be prevented in order to protect the reactivity of the strips.

Notes on symbols and marks

LOT	Batch code		Expiry Date
	Single use	IVD	In Vitro Diagnostic Use
	Manufactured by		Store at
	Please read package insert	CE	These test strips conform to the directive 98/79/EC(IVD-directive)
EC REP	Authorised Representative	REF	Catalogue number

DIRUI®
DIRUI Industrial Co.,Ltd.
95 Yunhe Street, New & High Tech.
Development Zone
Changchun, Jilin 130012 P.R. China
Tel :+86(431)85100409
Fax:+86(431)85172581
E-mail:dirui@dirui.com.cn
Http://www.dirui.com.cn

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H-Ca Series Urinalysis Reagent Strip User's Guide

Rev: 04/2019

[Product Name]
Common Name: H-Ca Series Urinalysis Reagent Strips
Trade Name: Urinalysis Reagent Strips
[Specification]
100 strips/bottle, 50 strips/bottle, 25 strips/bottle

[Intended Use]
H-Ca Series Urinalysis Reagent Strips(including H12-Ca,H13-Ca,H14-Ca Urinalysis Reagent Strips) are used for qualitative and semi-quantitative test of urobilinogen, bilirubin, ketone (acetacetic acid), blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid, microalbumin, creatinine and Ca in urine.
This product is used for in vitro diagnosis by professionals.

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
[Main Ingredient]
Leukocyte: 4.3% w/w pyrrole amino acid ester; 0.4% w/w diazonium salt; 92.6% w/w buffer; 2.7% w/w nonreactive ingredients.
Nitrite:1.3% w/w p-arsanilicacid-N-(1-Naphthol)-ethylenediamine; 0.9%w/w tetrahydroquinoline; 89.6% w/w buffer; 8.2% w/w nonreactive ingredients.
Protein: 0.1% w/w tetrabromphenol blue; 97.4% w/w buffer; 2.5% w/w nonreactive ingredients.
Blood: 26.0% w/w diisopropylbenzene dihydro peroxide; 1.5% w/w tetramethylbenzidine; 35.3% w/w buffer; 37.2 %w/w nonreactive ingredients.
Urobilinogen: 0.2%w/w fast B blue; 98.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Bilirubin: 0.6%w/w 2,4-dichloroaniline diazonium salt; 57.3%w/w buffer; 42.1%w/w nonreactive ingredients.
pH: 3.3% w/w bromcresol green; 55.0% w/w bromthymol blue; 41.7% w/w nonreactive ingredients.
Specific Gravity: 4.8%w/w bromthymol blue; 90.2%w/w poly(methyl vinyl ether-alt-maleic anhydride), 5%w/w sodium hydroxide.
Glucose: 1.7% w/w glucose oxidase (microbial.123U); 0.2 % w/w peroxidase(horseradish. 203U); 0.1% w/w potassium iodide; 71.8% w/w buffer; 26.2% w/w nonreactive ingredients.
Ketone: 5.7% w/w sodium nitroprusside; 29.9% w/w buffer; 64.4% w/w nonreactive ingredients.
Microalbumin: 2.2%w/w sulfone phthalein dye; 96.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Ascorbic Acid: 0.8%w/w 2,6-dichloroindophenolate hydrate; 40.7%w/w buffer; 58.5%w/w nonreactive ingredients.
Creatinine: 4.8%w/w 3,5-dinitro-benzoicaid; 85.2%w/w buffer; 10 %w/w nonreactive ingredients.
Ca: 2.5%w/w o-cresolphthalein complexon; 87.5%w/w buffer; 10 %w/w nonreactive

[Storage Condition and Validity]
Storage Condition: The strips must be stored at the temperature between 2℃-30℃. As to protect the activity of the reagent, must eliminate the humidity, sunlight and chemical reagent.
Validity: Stored in dry place and avoid sunlight, with the temperature between 2℃-30℃, the sealed validity is 2 years; After opened seal, cover the tap closely, stored in dry, avoid sunlight with the temperature between 2℃-30℃, the validity is 1 month.

[Applicable Instrument]
DIRUI H-100/300/500 Urine Analyzer

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[Test Method]
Room temperature for test:(25±5)℃
Visual Reading:
1. Immerse the reagent area of the strip in urine sample and take it out quickly and immediately.
2. Run the edge of the strip against the rim of the container to remove the excess urine.
3. Hold the strip up horizontally and compare the result on the strip with the color chart on the bottle label closely. Make note of the result. For a semi-quantitative result, take the result according to the time specified on the color chart. The pH and Protein can be read at any time within 60 seconds after dipping. For a qualitative result, the strip should be read within 1-2 minute after dipping. If a positive result is obtained, repeat the test and compare with the color chart at the time specified. Color changes beyond 2 minutes are of no diagnostic value.

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Urobilinogen

The reagent strips can detect urobilinogen in low amount as 3μmol/L (approximately 0.2 Ehrlich unit/dL) in urine. A result of 33μmol/L in urine indicate the critical value, representing the transition from normal to abnormal, which requires further check on patients and samples. The negative results are not final to determine the absence of urobilinogen.

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Creatinine

Normal urine creatinine concentration of adult is 0.6-2.0g/24h (about 4.4-17.6 mmol/L by test strip). High deviation exists on the test results of random urine sample creatinine, with the range from 0.9mmol/L to 26.5mmol/L. concentrated and morning urine have higher level of creatinine (maybe more than 17.7mmol/L tested by test strip). The concentration of the component required to test can be diluted due to over-ingestion of water or other substance, will result in typical low concentration urine (the tested result ≤50mg/dL). Visible hematuria(≥5mg/dL) or Cimitidine may result in creatinine testing result spurious rising.

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Result Table

Routine system of units and international system of units can be selected in urine analyzer test.

Test Item	Ab.	Print or Display Result			
		Routine System of Units		International System of Units	
Microalbumin	MA	10 mg/L	80 mg/L	10 mg/L	80 mg/L
		30 mg/L	150 mg/L	30 mg/L	150 mg/L
Creatinine	Cr	10 mg/dL	200 mg/dL	0.9 mmol/L	17.7 mmol/L
		50 mg/dL	300 mg/dL	4.4 mmol/L	26.5 mmol/L
		100 mg/dL		8.8 mmol/L	
Albumin-Creatinine Ratio	A: C	＜30 mg/g(normal)		＜3.4 mg/mmol(normal)	
		30-300 mg/g(abnormal)		3.4-33.9 mg/mmol(abnormal)	
		＞300 mg/g(highly abnormal)		＞33.9 mg/mmol(highly abnormal)	

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[Limitation of Test Method]

Like all lab tests, diagnosis result and therapeutic schedule can not be made according to any single diagnosis method.

Reagent strips' application is based on clinical analysis. In clinical sample, the sensitivity depends upon several factors: the variability of color perception, specific gravity, pH, and the lighting conditions change when the product is read visually. Each colour block or instrumental display value represents a range of values. Because of sample and reading variability, sample with analyte concentrations that fall between nominal levels may give results at either level. Results at levels greater than the second positive level for the Protein, Glucose, Ketone, and Urobilinogen tests will usually be within one level of the true concentration. Exact agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical system of the instruments.

[Performance Index]

Test Range and Sensitivity of H-Ca Series Strips

Item	Sensitivity	Instrument Test Range	Visual Test Range
Glucose(mmol/L)	2.8-5.6	Neg.-56	
Protein(g/L)	0.15-0.3	Neg.-3.0	Neg.-20.0
Microalbumin(mg/L)	30-80	10-150	10-150
Ketone (Acetacetic Acid)(mmol/L)	0.5-1.0	Neg.-7.8	Neg.-16
Blood(Ery/uL)	5-15	Neg.-200	
Bilirubin(μmol/L)	8.6-17	Neg.-103	
Nitrite(μmol/L)	13-22	Neg.-Pos.	
Leukocyte(Leuko/uL)	5-15	Neg.-500	
Urobilinogen(μmol/L)	3.4-17	3.4-135	
Ascorbic Acid(mmol/L)	0.6-1.4	0-5.7	
Creatinine(mmol/L)	0.9-4.4	0.9-26.5	
Ca(mmol/L)	2.0-2.5	1.0-10	
pH		5.0-9.0	5.0-8.5
Specific Gravity	—	1.005-1.030	1.000-1.030

[Points for Attention]

1.In vitro diagnosis for professional only.

2.Reagent strips must be stored in original bottle, do not take out strips except use them immediately, removing the desiccant is forbidden, tighten the cap immediately after taking out strips.

3.The expired should not be used, please make sure whether the strip is over the validity, and test it by quality control if the deterioration of strips occurs, which leads color shallower or deeper, or any question occurs with the unexpected results.

4.Water can not be used as negative control.

5.One time use only for each strip.

6.Please read the Guide carefully before use.

7.Please refer to the lab biological treatment of harmful materials to dispose the waste strips.

8.Do not place them in refrigerator, however they should be stored in dry circumstance, avoid direct sunlight, do not touch the reactive area of strip, ambient moisture, light and heat must be prevented in order to protect the reactivity of the strips.

Notes on symbols and marks

LOT	Batch code		Expiry date
	Single use	IVD	In Vitro Diagnostic Use
	Manufactured by		Store at
	Please read package insert	CE	These test strips conform to the directive 98/79/EC(IVD-directive)
EC REP	Authorised Representative	REF	Catalogue number

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FUS series Urinalysis Strip User's Guide

FUS series Automatic Urinalysis System Suited Strips

Rev:06/2015

In order to provide reliable results, please read the user's guide carefully before use.

[Product Name]
FUS series Urinalysis Strips

[Specification]
125 strips/bottle, 100 strips/bottle, 50 strips/bottle, 25 strips/bottle

[Intended Use]
FUS series Urinalysis Strips is used for qualitative or semi-quantitative test of urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid, microalbumin, creatinine, Ca.
Test items for different types of products are shown below:

Products Type	Test Item
FUS-10 Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH
FUS-11 Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid
FUS-11MA Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, pH, nitrite, leukocyte, glucose, specific gravity, microalbumin
FUS-12MA Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, microalbumin, ascorbic acid
FUS-13Cr Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, pH, specific gravity, nitrite, glucose, leukocyte, ascorbic acid, microalbumin, creatinine
FUS-14Ca Urinalysis Strips	urobilinogen, bilirubin, ketone, blood, protein, pH, specific gravity, nitrite, glucose, leukocyte, ascorbic acid, microalbumin, creatinine, Ca

[Testing Principle]
Glucose: The glucose oxidized by glucose oxidase catalyzes the formation of glucuronic acid and peroxide hydrogen. Peroxide hydrogen releases neo-ecotypes oxide [O] under the function of peroxidase. [O] oxidizes iodide potassium, which makes the color changes.
Bilirubin: The direct bilirubin and dichlorobenzene diazonium produce azo dyes in a strongly acid medium.
Ketone: The acetacetic acid and sodium nitroprusside cause reaction in alkaline medium, which produces violet compound.
Specific Gravity: Electrolyte (M+X-) in the form of salt in urine reacts with poly methyl vinyl ether and maleic acid (-COOH),which are weak acid ionic exchanger. The reaction produces hydrogenous ionogen, which reacts with pH indicator that causes the color change.
Blood: Hemoglobin acts as peroxidase. It can cause peroxidase release neo-ecotypes oxide (O). (O) oxidizes the indicator and make the color change subsequently.
pH: The method of pH indicator is applied.
Protein: This is based on the protein-error-of-indicator principle. Anion in the specific pH indicator attracted by cation on protein molecule makes the indicator further ionized, which changes its color.
Urobilinogen: Urobilinogen and diazonium produce red purplish dyes under the function of strong acid medium.
Nitrite: Nitrite in the urine and aromatic amino sulphanilamide are diazotized to form a diazonium compound. The diazonium compound reacting with tetrahydro benzo(h) quinolin 3-phenol causes the red azo dye.
Leukocyte: Granulocyte leukocytes in urine contains esterases that catalyze the hydrolysis of the pyrrole amino acid ester to liberate 3-hydroxy 5-pheny pyrrole. This pyrrole reacting with diazonium forms a purple azo dye.
Microalbumin: Sulfone phthalein dye have high sensitivity to microalbumin by the method of protein error.
Ascorbic Acid: Ascorbic acid, with 1,2-dihydroxy alkenes, under the alkaline condition, deoxidizes the blue 2,6-dichloroindophenolate into colorless N-(p-phenol)-2,6- dichloro-P-amine phenol.
Creatinine: Creatinine and 3,5-dinitro-benzoic acid react in strong alkaline condition, and generate colored compound.
Ca: Ca reacts with o-cresolphthalein complexone produce purple red, the color shade is proportional to Ca concentration.


[Main Ingredient]
Leukocyte: 4.3% w/w pyrrole amino acid ester; 0.4% w/w diazonium salt; 92.6% w/w buffer; 2.7% w/w nonreactive ingredients.
Nitrite:1.3% w/w p-arsanilicacid-N-(1-Naphthol)-ethylenediamine; 0.9%w/w tetrahydroquinoline; 89.6% w/w buffer; 8.2% w/w nonreactive ingredients.
Protein: 0.1% w/w tetrabromphenol blue; 97.4% w/w buffer; 2.5% w/w nonreactive ingredients.
Blood: 26.0% w/w diisopropylbenzene dihydro peroxide; 1.5% w/w tetramethylbenzidine; 35.3% w/w buffer; 37.2 %w/w nonreactive ingredients.
Urobilinogen: 0.2%w/w fast B blue; 98.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Bilirubin: 0.6%w/w 2,4-dichloroaniline diazonium salt; 57.3%w/w buffer; 32.1%w/w caffeine; 10.0%w/w nonreactive ingredients.
pH: 3.3% w/w methyl red; 55.0% w/w bromthymol blue; 41.7% w/w nonreactive ingredients.
Specific Gravity: 4.8%w/w bromthymol blue; 90.2%w/w poly(methyl vinyl ether co. maleic anhydride); 5.0%w/w sodium hydroxide.
Glucose: 1.7% w/w glucose oxidase (microbial.123U); 0.2 % w/w peroxidase(horseradish. 203U); 0.1% w/w potassium iodide; 71.8% w/w buffer; 26.2% w/w nonreactive ingredients.
Ketone: 5.7% w/w sodium nitroprusside; 29.9% w/w buffer; 64.4% w/w nonreactive ingredients.
Microalbumin: 2.2%w/w sulfone phthalein dye; 96.0%w/w buffer; 1.8%w/w nonreactive ingredients.
Ascorbic Acid: 0.8%w/w 2,6-dichloroindophenolate hydrate; 40.7%w/w buffer; 58.5%w/w nonreactive ingredients.
Creatinine: 4.8%w/w 3,5-dinitro-benzoicaid; 85.2%w/w buffer; 10 %w/w nonreactive ingredients.
Ca: 2.5%w/w o-cresolphthalein complexon; 87.5%w/w buffer; 10 %w/w nonreactive ingredients.

[Storage Conditions and Shelf Life]
Storage Conditions: The strips must be stored in dry place at the temperature between 2℃-30℃, and protected from direct sunlight and chemical reagents.
Shelf Life: Stored in dry place and protected from sunlight, with the temperature between 2℃-30℃, the sealed shelf life is 2 years; After opened, cover the tap closely, stored in dry, protected from sunlight with the temperature between 2℃-30℃, the shelf life is 1 month.

[Applicable Instrument]
DIRUI H-800/900/1000/1100/1200、HC-900 Urine analyzer
DIRUI FUS-2000 Urinalysis Hybrid

[Sample Requirements]
Collect fresh urine in a clean and dry container. Don't centrifuge the urine. Mix the sample well before the test. The urine test must be taken within two hour. All samples must always be taken and kept under sanitary conditions.
Note: The preservatives will not prevent the deterioration of ketones, bilirubin or urobilinogen. The growth of bacteria in the long-term storage sample may affect the test results of glucose, pH, nitrite and blood.

[Test Methods]
Room temperature for test:(25±5)℃
Visual Reading:
1. Immerse the reagent area of the strip in urine sample and take it out quickly and immediately.
2. Run the edge of the strip against the rim of the container to remove the excess urine.
3. Hold the strip up horizontally and compare the result on the strip with the color chart on the bottle label closely. Make note of the result. For a semi-quantitative result, take the result according to the specified time on the color chart. For a qualitative result, the strip should be read within 1-2 minute after dipping. If a positive result is obtained, repeat the test and compare with the color chart at the specified time. Color changes beyond 2 minutes are of no diagnostic value.

Instrumental Reading:
Follow the directions given in appropriate instrument user manual.

[Result Explanation]
Glucose
The test is for specificity of glucose. There is no false positive result occurred in reagent strip, caused by any other substance in urine. When the ascorbic acid concentration≥2.8mmol/L or acetoacetic acid concentration≥1.0mmol/L, the sample of glucose concentration is 3～7mmol/L may occur false negative result.
Bilirubin
Normally, even the most sensitive method can't detect bilirubin in urine. It is abnormal to have little bilirubin in urine, which requires further inspection. Medicines that dyes urine red and anything that shows red itself in an acid medium e.g., phenazopyridine may affect the test result. High concentration of the ascorbic acid may cause false negative result.
Ketone
The reagent strip acetone reacts with acetoacetic acid of urine. It doesn't do with acetone or β-hydro butyric acid. Normal urine sample usually shows negative results in the test. False positive results may occur in highly pigmented urine or those containing a large amount of levodopa metabolites.
Specific Gravity
The reagent strip for Specific Gravity allows the urine sample specific gravity between 1.000 and 1.030. In general, the mean error between the results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.005 may be added to readings from urines with pH equal or greater than 6.5. Urine reading instrument can automatically make these adjustments in strip-readings. The urine nonionic constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated

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specific gravity readings may occur in the presence of moderate quantities of protein(1g/L-7.5g/L).

Blood

'Trace' reaction may vary among the patients. Clinical judgments are required for individual cases. The presence of green spots (intact erythrocytes) or green color (haemoglobin/myoglobin) on the reagent area within 60 seconds indicates for further diagnostic check. Blood is often found in the urine of the menstruating females. Haemoglobin 150μg/L-620μg/L is approximately equivalent to 5-15 cells/μL intact erythrocytes. The reagent strip is highly sensitive to haemolobin and thus can be used as a supplementary to the microscopic examination. The sensitivity of the strip might be reduced in urine with a large amount of specific gravity. The strips are equally sensitive to myglobin as to haemoglobin. Certain oxidizing contaminants, such as hypochlorite, may lead to false positive results. Microbial peroxidase associated with urinary tract infection may also produce a false positive result. Vitamine in urine may not influence the result of the test.

pH

The strip tests for pH values are generally in the range of 5.0-8.5 visually and 5.0-9.0 instrumentally.

Urobilinogen

The reagent strips can detect urobilinogen in low amount as 3μmol/L (approximately 0.2 Ehrlich unit/dL) in urine. The normal content of Urobilinogen is 3μmol/L-16μmol/L. A result of 33μmol/L in urine indicate the critical value, representing the transition from normal to abnormal, which requires further check on patients and samples. The negative results are not final to determine the absence of urobilinogen.

Nitrite

Gram-negative bacteria in urine converts nitrate(derived from foods)into nitrite. The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink color development should be taken as positive result. The degrees of color development the numbers of bacteria are not in direct proportion. The negative result doesn't mean the existence of bacteria in a large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite. (2)when urine has not remained in the bladder long enough(four hours up)to let the nitrate covert into nitrite. (3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 1.4mmol/L ascorbic acid or less won't interfere the result.

Leukocyte

Test area react with esterase in leucocytes (granulocytic leukocytes). Normal urine sample generally yield negative result; positive results (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance; however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random sample from females due to contamination of the sample by vaginal discharge. Elevated glucose concentrations (160mmol/L) or high specific gravity may cause decreased test results.

Ascorbic Acid

The test area can detect the ascorbic acid in urine. Through the ascorbic acid detection, we will know the level of ascorbic acid in the body and the effect degree that the ascorbic acid bring to the test on glucose, bilirubin, blood and nitrite. It will reduce the sensitivity when the oxidant (such as potassium permanganate, hypochlorite)present in the urine.

Protein and Microalbumin

Albumin concentration in normal adult urine is less than 20 mg/L. A result of 20-200 mg/L indicates clinical microalbumin. A result of higher than 200 mg/L indicates clinical albuminuria. A albumin discharge rate of 30-300 mg/24 hour indicates a warning of urine microalbumin symptom. A albumin discharge rate of higher than 300 mg/24 hour indicates a warning of urine albumin symptom. Urine albumin discharge may be elevated by temporary sport, urinary tract infection or acute fever disease.

Microalbumin strip is intended for detecting urine microalbumin, and is sensitive to microalbumin. Its sensitivity for other protein is 9 times lower than microalbumin. The test is less sensitive to mucoproteins and globulins, which are generally detected at levels of 60 mg/dL or higher. A visibly bloody urine (≥5 mg/dL) may cause falsely elevated results.

Microalbumin in urine may lead to:

(1)Occasionally presented urine microalbumin may be functional albuminuria or postural proteinuria caused by physiological albuminuria, such as diet, sport, mental stress, etc.

(2)Continuously presented urine microalbumin, microalbumin and glucose, or microalbumin and blood are with great clinical significance.

Creatinine

Normal urine creatinine concentration of adult is 0.6-2.0g/24h (about 4.4-17.7 mmol/L by test strip). High deviation exists on the test results of random urine sample creatinine, with the range from 0.9mmol/L to 26.5mmol/L. concentrated and morning urine have higher level of creatinine (maybe more than 17.7mmol/L tested by test strip). The concentration of the component required to test can be diluted due to over-ingestion of water or other substance, will result in typical low concentration urine (the tested result ≤50mg/dL). Visible hematuria(≥5mg/dL) or Cimitidine may result in creatinine testing result spurious rising.

Albumin-Creatinine Ratio

Urine albumin content of adult is less than 30 mg albumin/g creatinine (3.4 albumin/mmol creatinine). A ratio of 30-300mg/g (3.4-33.9mg/mmol) indicates microalbuminuria. A ratio of more than 300mg/g (＞33.9mg/mmol) indicates albuminuria.

Ca:A great deal of Mg(more than 10mmol/L) will seduce urine Ca rising.

Result Table

Routine system of units and international system of units can be selected in urine analyzer test

Test Item	Ab.	Print or Display Result			
		Routine System of Units		International System of Units	
Microalbumi n	MA	10 mg/L	80 mg/L	10 mg/L	80 mg/L
		30 mg/L	150 mg/L	30 mg/L	150 mg/L
Creatinine	Cr	10 mg/dL	200 mg/dL	0.9 mmol/L	17.7 mmol/L
		50 mg/dL	300 mg/dL	4.4 mmol/L	26.5 mmol/L
		100 mg/dL		8.8 mmol/L	
Albumin-Creatinine Ratio	A: C	＜30 mg/g(normal)		＜3.4 mg/mmol(normal)	
		30-300 mg/g(abnormal)		3.4-33.9 mg/mmol(abnormal)	
		＞300 mg/g(highly abnormal)		＞33.9 mg/mmol(highly abnormal)	

Remark: shaded area indicates abnormal result.

[Limitation of Test Method]

Like all the other laboratory tests, definitive diagnosis or therapeutic decisions should not be made or based on any single result or method.

The application of the reagent strips is based on the study of clinical analysis. For the clinical urine samples, the sensitivity depends on the following factors: color variability, SG, pH, lighting condition, etc. Visual reading and instrumental reading represent a range of values. Because of urine sample and reanding variability, there is a deviation between the value of tested analyte and actual value. For the tests of protein, glucose, ketone and urobilinogen, the deviation of positive value above the second positive level should be within a "+". Because of the inherent difference of sensitization between eyes and optical system, visual reading and instrumental reading might not coincide completely.

[Product Performance Index]

Test Range and Sensitivity of FUS series Urinalysis Strips

Item	LOD	Instrumental Test Range	Visual Test Range
Glucose(mmol/L)	5.6	Neg.-56	
Protein(g/L)	0.2	Neg.-3.0	Neg.-20.0
Microalbumin1(g/L)	0.15	Neg.-0.15	Neg.-0.15
Ketone(mmol/L)	0.5	Neg.-7.8	Neg.-16
Blood(ca CELLS/μL)	10	Neg.-200	
Bilirubin(μmol/L)	17	Neg.-103	
Nitrite(μmol/L)	18	Neg.- 36	
Leukocyte(ca CELLS/μL)	15	Neg.-500	
Urobilinogen(μmol/L)	3.4	3.4-135	
Ascorbic Acid(mmol/L)	0.6	0-5.7	
Creatinine(mmol/L)	0.9	0.9-26.5	
Ca(mmol/L)	1.0	1.0-10	
pH	—	5.0-9.0	5.0-8.5
Specific Gravity	—	1.005-1.030	1.000-1.030
Microalbumin2(mg/L)	10	10-150	10-150

Notes:

Microalbumin 1 is used for FUS-11MA and FUS-12MA Urinalysis Strips.
Microalbumin 2 is used for FUS-13Cr and FUS-14Ca Urinalysis Strips.

Specificity

Detect the Anti-VC Standard Solution presented in the following table using the strips, the test results should not be negative.

Item	Anti-VC Standard Solution Concentration	VC Concentration
NIT	36μmol/L	1.4mmol/L
GLU	14mmol/L	2.8mmol/L

[Points for Attention]

1.In vitro diagnosis for professional only.
2.Reagent strips must be stored in original bottle, do not take out strips except use them immediately. Do not remove the desiccants, cap the bottle immediately and tightly after taking out strips.
3.The expired should not be used, please make sure whether the strip is over the validity, and test it by quality control if the deterioration of strips occurs, which leads color shallower or deeper, or any question occurs with the unexpected results.
4.Water can not be used as negative control.
5.Each strip can be used only once.
6.Please read the User's Guide carefully before use.
7.Please dispose the used strips according to Treatment Regulations of Lab Biohazard Materials.
8.Do not place them in refrigerator, however they should be stored in dry place and protected from direct sunlight, do not touch the reactive area of strip, ambient moisture, light and heat must be prevented in order to protect the reactivity of the strips.

Notes on symbols and marks

LOT	Batch code		Use by
	Single use	IVD	In Vitro Diagnostic Use
	Manufactured by		Store at
	Please read package insert	CE	These test strips conform to the directive 98/79/EC(IVD-directive)
EC REP	Authorised Representative	REF	Catalogue number

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DIRUI®

Urinalysis Control (Negative)

【Name】Urinalysis control
【Volume】8ml/bottle、4ml/bottle
【Intended Use】
Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin,Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.
【Test Principle】
Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.
【Main Composition】
Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinagen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.
Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.
【Storage and Validity】
Its validity period is 12 months when stored at 2℃～ 8℃, sealed and protected from sunlight. Its validity period is 1 month when stored at 2℃～ 8℃, sealed and protected from sunlight after the first use.
【Applied Analyzer】Urine Analyzer
【Test Method】
1. Place the Control at room temperature(18℃-25℃) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2.Using method on semi-automatic urine analyzer: Adjust the analyzer to QC status; drip the Control on the strip, the strip side with control should be placed upward; removethe extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer.
3.Using method on automatic urine analyzer: Adjust the analyzer to QC status; pour the Control into a clean tube; place the tube into corresponding position of the analyzer.
【Result Explanation】
All QC results within the target range indicates qualified, otherwise, the result is unqualified.
【Limits】
Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.
【Matters Need Attention】
1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3.Tighten the cap immediately after use, and store it at 2℃-8℃.
【Reference】
Refer to the supplied control target list.
The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might becaused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.
【Documents】ZL01 1 28009.3 Urinalysis Control and Preparation
【Instruction Approved and Modified Date】11/2014
【Symbols】

Store at | Please read package insert | Use by | LOT Batch code | IVD In Vitro Diagnostic Use | Manufactured by
EC REP Authorised Representative | CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD) | REF Catalogue number

DIRUI®

Urinalysis Control (Negative)

【Name】Urinalysis control
【Volume】8ml/bottle、4ml/bottle
【Intended Use】
Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin,Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.
【Test Principle】
Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.
【Main Composition】
Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinagen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.
Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.
【Storage and Validity】
Its validity period is 12 months when stored at 2℃～ 8℃, sealed and protected from sunlight. Its validity period is 1 month when stored at 2℃～ 8℃, sealed and protected from sunlight after the first use.
【Applied Analyzer】Urine Analyzer
【Test Method】
1. Place the Control at room temperature(18℃-25℃) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2.Using method on semi-automatic urine analyzer: Adjust the analyzer to QC status; drip the Control on the strip, the strip side with control should be placed upward; removethe extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer.
3.Using method on automatic urine analyzer: Adjust the analyzer to QC status; pour the Control into a clean tube; place the tube into corresponding position of the analyzer.
【Result Explanation】
All QC results within the target range indicates qualified, otherwise, the result is unqualified.
【Limits】
Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.
【Matters Need Attention】
1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3.Tighten the cap immediately after use, and store it at 2℃-8℃.
【Reference】
Refer to the supplied control target list.
The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might becaused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.
【Documents】ZL01 1 28009.3 Urinalysis Control and Preparation
【Instruction Approved and Modified Date】11/2014
【Symbols】

Store at | Please read package insert | Use by | LOT Batch code | IVD In Vitro Diagnostic Use | Manufactured by
EC REP Authorised Representative | CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD) | REF Catalogue number

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DIRUI®

Urinalysis Control (Positive)

【Name】Urinalysis control
【Volume】8ml/bottle、4ml/bottle
【Intended Use】
Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin,Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.
【Test Principle】
Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.
【Main Composition】
Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinagen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.
Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.
【Storage and Validity】
Its validity period is 12 months when stored at 2℃～ 8℃, sealed and protected from sunlight. Its validity period is 1 month when stored at 2℃～ 8℃, sealed and protected from sunlight after the first use.
【Applied Analyzer】Urine Analyzer
【Test Method】
1. Place the Control at room temperature(18℃-25℃) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2.Using method on semi-automatic urine analyzer: Adjust the analyzer to QC status; drip the Control on the strip, the strip side with control should be placed upward; removethe extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer.
3.Using method on automatic urine analyzer: Adjust the analyzer to QC status; pour the Control into a clean tube; place the tube into corresponding position of the analyzer.
【Result Explanation】
All QC results within the target range indicates qualified, otherwise, the result is unqualified.
【Limits】
Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.
【Matters Need Attention】
1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3.Tighten the cap immediately after use, and store it at 2℃-8℃.
【Reference】
Refer to the supplied control target list.
The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might becaused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.
【Documents】ZL01 1 28009.3 Urinalysis Control and Preparation
【Instruction Approved and Modified Date】07/2013
【Symbols】

Store at | Please read package insert | Use by | LOT Batch code | IVD In Vitro Diagnostic Use | Manufactured by
EC REP Authorised Representative | CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD) | REF Catalogue number

DIRUI®

Urinalysis Control (Positive)

【Name】Urinalysis control
【Volume】8ml/bottle、4ml/bottle
【Intended Use】
Urinalysis control(referred to as control) is used for quality control of urinalysis strips and analyzers. It can conduct quality control for 13 items of strips, such as Glucose, Bilirubin,Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leucocytes, Microalbumin, Creatinine and Ca.
【Test Principle】
Chemistry reagent contained in the Control reacts with the ingredient within the urinalysis strips, thus, strips color changed.
【Main Composition】
Positive: phosphate buffer 0.2%w/w, glucose 1.0%w/w, sodium chloride 0.5%w/w, hemoglobin 0.1%w/w, albumin 0.7%w/w, ethyl acetoacetate 0.9%w/w, sodium nitrite 0.3%w/w, esterase 0.5%w/w, urea 2%w/w, creatinine 0.2%w/w, Ca 0.1%w/w, bilirubin substitute 0.1%w/w, urobilinagen substitute 0.1%w/w, other non-reactive substances and stabilizers 93.3%w/w.
Negative: urea 2%w/w, sodium chloride 0.5%w/w, phosphate buffer 0.2%w/w, other non-reactive substances and stabilizers 97.3%w/w.
【Storage and Validity】
Its validity period is 12 months when stored at 2℃～ 8℃, sealed and protected from sunlight. Its validity period is 1 month when stored at 2℃～ 8℃, sealed and protected from sunlight after the first use.
【Applied Analyzer】Urine Analyzer
【Test Method】
1. Place the Control at room temperature(18℃-25℃) for certain time until it returns to normal temperature. And mix the Control evenly through bottom up and down slightly.
2.Using method on semi-automatic urine analyzer: Adjust the analyzer to QC status; drip the Control on the strip, the strip side with control should be placed upward; removethe extra control with absorbent paper, the absorbent paper can not touch the strip pads(avoid cross-contamination); place the strip on the waiting place of analyzer.
3.Using method on automatic urine analyzer: Adjust the analyzer to QC status; pour the Control into a clean tube; place the tube into corresponding position of the analyzer.
【Result Explanation】
All QC results within the target range indicates qualified, otherwise, the result is unqualified.
【Limits】
Bilirubin and Urobilinogen used for the Control are replaced by chemical substances. Therefore, there is a trace difference in color between the result of control and the result of urine.
【Matters Need Attention】
1. The Control follows outside of the control bottle can not be collected back to avoid cross-contamination.
2. Control-skin contact should be avoided. Rinse the contact part with plenty of water if accidental skin contact occur.
3.Tighten the cap immediately after use, and store it at 2℃-8℃.
【Reference】
Refer to the supplied control target list.
The supplied reference is obtained through repeat tests. For the same lab, the results might be different when testing on different dates. The difference might becaused by environment, instrument, reagent or test method difference. Even so, the test results should be within the supplied reference range.
【Documents】ZL01 1 28009.3 Urinalysis Control and Preparation
【Instruction Approved and Modified Date】07/2013
【Symbols】

Store at | Please read package insert | Use by | LOT Batch code | IVD In Vitro Diagnostic Use | Manufactured by
EC REP Authorised Representative | CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD) | REF Catalogue number

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(Concentrated type)

Instruction

Product Name
Cleaning liquid
Volume
200ml
Intended Use
It is used for cleaning and rinsing of tubing system. It will effectively clean the remained urine and other test-interfering particles to avoid cross-contamination among urine sample.
Main Ingredient
Lowfoaming surfactant
For Instrument
Automatic Urine Analyzer
Usage
1.Pour the detergent into a clean barrel with volume of 10L.
2.Rinse the detergent bottle with distilled water, and pour the water into the 10L barrel.
3.Dilute the detergent to 10L level with distilled water and mix it.
4.Connect the 10L barrel with Automatic Urine Analyzer, please refer to the Automatic Urine Analyzer user manual for connection method.
Important Notice
1.This detergent is used for Automatic Urine Analyzer only, not applicable to other instruments.
2.It should be used together with Automatic Urine Analyzer under room temperature.
3.The detergent contains chemical medicine, like preservative, etc. which can not be swallowed. Rinse with water immediately if the detergent contacts your skin.
Storage Condition and Period of Validity
Stored at 2℃~30℃, sealed, the validity is 1 year.Its validity after diluted is 1 month at 2℃~30℃.
Manual Approved and Modified Date
08/2013
Symbol Explanation

Please read instruction	Storage Temperature	Expiry	REF Catalogue number
Manufacturer	EC REP Authorised Representative	IVD In Vitor Diagnosis	
CE European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD)		LOT Lot No	

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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DIRUI®

H Series Urine Analyzer
Cleaning Liquid for Refractometer and Turbidimeter

Instruction

【Product Name】 H Series Urine Analyzer Cleaning Liquid for Refractometer and Turbidimeter

【Volume】 50ml

【Intended Use】 To clean the refractometer and turbidimeter of H series Urine Analyzer

【Main Composition】 Surfactant 、5%W/W sodium hypochlorite

【Storage Conditions and Guarantee Period】 Kept airtight, stored in a dark and dry place at 2℃-30℃, the shelf life is 12 months. After unsealed, stored at 2℃-30℃, the shelf life is 60 days.

【Applicable Instrument】 H series Automatic Urine Analyzer

【Usage】 Pour at least 8mL cleaning liquid into a tube, and insert the tube into the STAT gate. Press "Rinse Refractometer and Turbidimeter" for rinsing.

【Matters Needing Attention】

1. Avoid contact with skin and eyes. Wash it away with plenty of water if careless contact occurs. Seek medical advice immediately when necessary.
2. It should be used when returned to room temperature. Tighten it immediately after use.
3. Avoid direct sunlight, and avoid freezing during storage.
4. Sample and waste liquid are potentially biologically infectious. The operator should obey the regulations of safe lab operation, and dispose the waste according to the local rules of medical waste, infectious waste, and industrial waste, etc.

【Date of Instruction Approval and Modification】 07/2013

【Explanation for Symbols】

Please Read the Instruction	Storage Temperature	Manufacturer
LOT Lot Number	Expiration Date	IVD In vitro Diagnostic Reagent
CE Meet the Requirements of IVDD	EC REP	Authorised Representative
R31-36/38 S2-25 Danger Level	REF	Catalogue Number

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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DIRUI®

H series urine analyzer

control liquid for Specific gravity

Instructions

【Product name】H series urine analyzer control liquid for Specific gravity
【Dimension】8ml/vial
【Intended purpose】H series urine analyzer control liquid for Specific gravity is used for DIRUI's H series urine analyzer's refractometer's quality control.
【Principle】Using refraction method, refractometer tests the quality control liquid.
【Composition】carbamide（4.5%w/w）、Sunset yellow.（0.1%w/w）、Hydrazine yellow（0.1%w/w）、reatinine（0.125%w/w）、Sodium Chloride（4.5%w/w）、potassium chloride（1.5%w/w）. others（89.18%w/w）.
【Storage and expiration】Stored in dry place and protected from sunlight, with the temperature between 2℃-8℃, the sealed shelf life is 12 months. After opened, stored at the temperature between 2℃-8℃, the shelf life is 30 days.
【Applicable analyzer】H series urine analyzer.
【Reference（refrence range）】

refrence range	Vertical
Level1	1.003～1.007
Level2	1.026～1.034
Level3	1.046～1.054

【Procedure】According to H series urine analyzer user manual.
【Result explanation】If the test result of the control liquid falls in the range of reference, it means the instrument is normal; if not in the range, it means the status of instrument is not normal. Please rerun the test or calibrate refractometer.
【Limitation of testing method】The turbidity control is only applicable to DIRUI's H series urine analyzer.
【Warning and precautions】
1、Please resume the product to the room temperature （18℃～25℃）and blend it well before use.
2、keep it sealed in cooling condition under temperature 2℃～8℃.
3、Avoid contact with skin and eyes. In case of skin contact, rinse with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, induce vomiting and seek medical advice.
【Product performance】

Level	Accuracy
Level1	SG＝1.005±0.002
Level2	SG＝1.030±0.004
Level3	SG＝1.050±0.004

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Please read package insert	Store at	Manufactured by	LOT	Batch code	
IVD	In Vitro Diagnostic Use	EC REP	Authorised Representative	REF	Catalogue Number
CE	The calibration liquid for Turbidity conform to the directive 98/79/EC (IVD-directive)			Use by	

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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control liquid for Turbidity

Instructions

【Product name】H series urine analyzer control liquid for Turbidity
【Dimension】 8ml/vial
【Intended purpose】H series urine analyzer control liquid for Turbidity is used for DIRUI's H series urine analyzer turbidity quality control.
【Principle】Using scattering theory, turbidimeter tests control liquid.
【Composition】Fu MA hydrazine turbidity quality control liquid （800NTU）.
【Storage and expiration】 Stored in dry place and protected from sunlight, with the temperature between 2℃-8℃, the sealed shelf life is 12 months. After opened, stored at the temperature between 2℃-8℃, the shelf life is 30 days.
【Applicable analyzer】H series urine analyzer.
【Procedure】According to H series urine analyzer user manual.
【Reference（Reference range）】

Level	Reference range
Level1	Clear
Level2	Cloudy

【Result explanation】If the test result of the turbidity control falls in the range of reference, it means the instrument is normal; if not in the range, it means the status of instrument is not normal.Please rerun the test or calibrate turbidimeter.
【Limitation of testing method】The turbidity control is only applicable to DIRUI's H series urine analyzer
【Warning and precautions】Please blend well and avoid contact with skins. In case of skin contact, rinse with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, induce vomiting and seek medical advice.
【Product performance】

Level	Accuracy
Level1	Turbidity=200±30NTU
Level2	Turbidity=700±30NTU

【Date of Instruction Approval and Modification】01/2014
【Explanation for Symbols】

Please read package insert	Store at	Manufactured by	LOT Batch code
Use by	IVD In Vitro Diagnostic Use	EC REP Authorised Representative	
CE The control liquid for Turbidity conform to the directive 98/79/EC (IVD-directive)		REF Catalogue number	

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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Color control

Instructions

【Product name】H series urine analyzer color control
【Dimension】8ml/vial
【Intended purpose】H series urine analyzer color control is used for DIRUI's H series urine analyzer color quality control.
【Principle】Urine analyzer's color sensor tests the color control.
【Composition】Red：Amaranth；Blue： Light blue；Green：Amaranth、Light blue.
【Storage and expiration】Stored in dry place and protected from sunlight, with the temperature between 2℃-8℃, the sealed shelf life is 12 months. After opened, stored at the temperature between 2℃-8℃, the shelf life is 30 days.
【Applicable analyzer】H series urine analyzer
【Procedure】According to H series urine analyzer user manual.
【Reference（reference range）】

Level	Reference range
Red	Red
Green	Green
Blue	Blue

【Result explanation 】If the test result of the color control falls in the range of reference, it means the instrument is normal; if not in the range, it means the status of instrument is not normal.
【Limitation of testing method】The color control is only applicable to DIRUI's H series urine analyzer.
【Warning and precautions】
1、Please resume the product to the room temperature （18℃～25℃）and blend it well before use.
2、keep it sealed in cooling condition under temperature 2℃～8℃.
3、Avoid contact with skin and eyes. In case of skin contact, rinse with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, induce vomiting and seek medical advice.
【Date of Instruction Approval and Modification】01/2014
【Explanation for Symbols】

Please read package insert	Store at	Manufactured by	LOT Batch code
Use by	IVD In Vitro Diagnostic Use	EC REP Authorised Representative	
CE The Color control conform to the directive 98/79/EC (IVD-directive)	Keep in Dark Places	REF Catalogue Number	

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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DIRUI®

H series urine analyzer

calibration liquid for Specific gravity

Instructions

【Product name】H series urine analyzer calibration liquid for Specific gravity
【Dimension】8ml/vial
【Intended purpose】H series urine analyzer calibration liquid for Specific gravity is used for DIRUI's H series urine analyzer's refractometer's calibration.
【Principle】Using refraction method, refractometer tests the calibration liquid.
【Composition】carbamide（3.6%w/w）、Sodium Chloride（3.6%w/w）、potassium chloride（1.2%w/w）、creatinine（0.001%w/w）、Sunset yellow（0.08%w/w）、Hydrazine yellow（0.08%w/w）、 others（91.44%w/w）.
【Storage Conditions and Guarantee Period】 Stored in dry place and protected from sunlight, with the temperature between 2℃-8℃, the sealed shelf life is 12 months. After opened, stored at the temperature between 2℃-8℃, the shelf life is 30 days.
【Applicable analyzer】H series urine analyzer.
【Procedure】According to H series urine analyzer user manual.
【Result explanation】After refractometer calibration, the screen will show “calibration Ok” or “calibration not OK”. “Calibration OK” means the instrument is in normal status; “Calibration not OK” means the calibration failed, please rerun calibration or clean the refractometer.
【Limitation of testing method】The calibration liquid is only applicable to DIRUI's H series urine analyzer.
【Warning and precautions】
1、Please resume the product to the room temperature 18℃～25℃ and blend it well before use.
2、keep it sealed in cooling condition under temperature 2℃～8℃.
3、Avoid contact with skin and eyes. In case of skin contact, rinse with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, induce vomiting and seek medical advice.
【Product performance】Accuracy：SG＝1.040±0.005
【Date of Instruction Approval and Modification】01/2014
【Explanation for Symbols】

Please read package insert	Store at	Manufactured by	LOT Batch code
IVD In Vitro Diagnostic Use	EC REP Authorised Representative		REF Catalogue Number
CE The calibration liquid for Turbidity conform to the directive 98/79/EC (IVD-directive)			Use by

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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DIRUI Turbilitor Calibration Liquid

Instructions

【Product name】H series urine analyzer calibration liquid for Turbidity
【Dimension】8ml/vial
【Intended purpose】H series urine analyzer calibration liquid for Turbidity is used for DIRUI's H series urine analyzer turbidimeter calibration.
【Principle】Using scattering theory, turbidimeter tests calibration liquid.
【Composition】Fu MA hydrazine turbidity quality control liquid （800NTU）.
【Storage Conditions and Guarantee Period】 Stored in dry place and protected from sunlight, with the temperature between 2℃-8℃, the sealed shelf life is 12 months. After opened, stored at the temperature between 2℃-8℃, the shelf life is 30 days.
【Applicable analyzer】H series urine analyzer.
【Procedure】According to H series urine analyzer user manual.
【Result explanation】After turbidimeter calibration, the screen will show “calibration Ok or “calibration not OK”. “Calibration OK” means the instrument is in normal status; “Calibration not OK” means the calibration failed, please rerun calibration or clean the turbidimeter.
【Limitation of testing method】The turbidity calibration liquid is only applicable to DIRUI's H series urine analyzer.
【Warning and precautions】Please blend well and avoid contact with skins. In case of skin contact, rinse with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, induce vomiting and seek medical advice.
【Product performance】accuracy：Turbidity=400±30NTU
【Date of Instruction Approval and Modification】01/2014
【Explanation for Symbols】

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Instruction

【Product name】Standard Solution

【Package specification】30mL, 60mL, 125mL, 500mL, 2L, 5L, 2×125mL, 3×125mL, 4×125mL;

【Intended use】This product is used to calibrate Urine Sediment Analyzer to guarantee to get accurate test results.

【Test principle】Flat flow cell imaging principle

【Main components】QC blood: 0.02%~0.1%, disodium hydrogen phosphate: 6.7mmol/L.
Note: fixed value is different in different lot number (see the label for values of each batch).

【Storage conditions and shelf life】Stored at 2℃~8℃, dry, sealed and protected from sunlight. The shelf life refers to the label. It shall be sealed and kept in a dark place at the temperature of 2℃~8℃, and the open-vial shelf life is 7 days.

【Test method】Refer to the user manual of the applicable instrument.

【Sample requirement】Fresh urine sample without centrifugation.

【Product performance indices】

1. Accuracy
 Standard solution particles :900particles/μL~1,300particles/μL; relative bias shall be within the range of ±8.0%.
2. Homogeneity
 CV≤10.0%

【Precautions】

1. Before use, gently shake it and avoid over exertion;
2. Avoid contact with skin and eyes;
3. Rinse eyes with the flow of warm water for a few minutes if the product splashes into eyes(ensuring eyelids open). Seek medical care immediately if pain or swelling happens. Rinse with soap water and water once skin contact occurs.

【Bibliography】

1. Ferris JA urine microscopy comparison & standard: Lab Med14: 659-662,1983.
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

【Instruction Approved and Modified Date】

Approval date: 02/2017, Modified date: 01/2018

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Please read the Specification	Storage temperature	Manufacturer	In Vitro Diagnostic Reagent
Batch No.	Upward	Expiration date	Kept in a dark place
Comply with In-Vitro Diagnostic Devices Directive (98/79/EC)			Authorised Representative

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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【Product name】Positive Control

【Package specification】

Positive control: 40mL, 60mL, 125mL, 500mL, 2L, 5L, 2×125mL, 3×125mL, 4×125mL;

【Intended use】This product is used for QC during test.

【Test principle】Flat flow cell imaging principle

【Main components】QC blood: 0.02%~0.1%, disodium hydrogen phosphate: 6.7mmol/L.

Note: fixed value is different in different lot number (see the label for values of each batch).

【Storage conditions and shelf life】Stored at 2℃~8℃, dry, sealed and protected from sunlight. The shelf life refers to the label. It shall be sealed and kept in a dark place at the temperature of 2℃~8℃, and the open-vial shelf life is 30 days.

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【Product performance indices】

1. Accuracy: Positive control particles: 900particles/μL~1,300particles/μL; relative bias shall be within the range of ±8.0%.
2. Homogeneity: CV≤10.0%

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【Instruction Approved and Modified Date】

Approval date: 02/2017, Modified date: 01/2018

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Please read the Specification	Storage temperature	Manufacturer	In Vitro Diagnostic Reagent
Batch No.	Upward	Expiration date	Kept in a dark place
Comply with In-Vitro Diagnostic Devices Directive (98/79/EC)		Authorised Representative	

DIRUI INDUSTRIAL CO.,LTD.
95 Yunhe Street, New & High Tech.
Development Zone Changchun, Jilin 130012 P.R.China
Tel: +86 431 85100409 Fax: +86 431 85172581
E-mail:dirui@dirui.com.cn http://www.dirui.com.cn

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
Urine Sediment Analyzer Detergent

User Manual

【Product Name】Urine Sediment Analyzer Detergent
【Packaging Specification】300mL, 500mL
【Intended Use】Rinse the pipeline and flow cell
【Test Principle】Flow cell technique and high-speed imaging technique
【Main Ingredient】Sodium hypochlorite
【Storage Condition and Shelf life】Stored at 2℃~30℃ in a dry place, protected from sunlight, and sealed. The shelf life is 12 months. After opened, the shelf life is 30 days.
【For Instrument】FUS Series Urine Sediment Analyzer and Urinalysis Hybrid
【Usage】See User Manual of FUS Series Urine Sediment Analyzer and Urinalysis Hybrid.
【Sample Requirement】Fresh urine sample without centrifugation.
【Performance Indicator】pH is 12.10±0.5 at（25±1）℃
【Matters Needing Attention】

1.Avoid contact with skin and eyes;
2.Flush eyes with the flow of warm water for a few minutes if the product splashed into eyes (ensuring eyelids opened). Seek medical advice immediately if pain or swelling happened. Flush with soap water and water once skin contact occurs.

R31: Exposed to acid will release toxic gas　R36/38: Irritating to eyes and skin
S2: Out of reach for children　S25: Keep away from eyes
Paste hazard label, danger and safety instructions on DIRUI product comply with European Community regulations.
【Bibliography】
1.Ferris JA Comparison and standards of urine microscopy Lab Med 14:659-662,1983
2.Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.
【Symbol Explanation】

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2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition.

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1.Avoid contact with skin and eyes;
2.Rinse eyes with the flow of warm water for a few minutes if the product splashed into eyes (ensuring eyelids opened). Seek medical care immediately if pain or swelling happened. Rinse with soap water and water once skin contact occurs.

【Bibliography】

1.Ferris JA Comparison and standards of urine microscopy Lab Med 14:659-662,1983
2.Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

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Please read the User Manual	Stored at	Manufacturer
IVD In vitro diagnostic reagent	LOT Lot Number	Valid until
CE Conform to EU in vitro diagnostic medical equipment directives98/79/EC		
EC REP Authorised Reprsentative	REF	Catalogue number

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1. Ferris JA Comparison and standards of urine microscopy Lab Med 14:659-662, 1983
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition.

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DIRUI®

Urine Sediment Analyzer Reagent -- **Sheath** **User Manual**

【Name】Sheath for Urine Sediment Analyzer
【Specification】5L/10L/15L/20L
【Intended Use】Used to wrap the sample to form a sheath flow to conduct urine physical component counting
【Test Principle】Flat flow cell imaging principle
【Ingredients】PBS buffer, sodium chloride, EDTA, non-ionic surfactant
【Storage Conditions and Shelf Life】Stored in the shade at 2℃～30℃ and kept dry and airtight, the shelf life is 18 months. The open vial shelf life is 60 days.
【For Instrument】Urine Sediment Analyzer
【Usage】See user manual of urine sediment analyzer.
【Sample Requirement】Fresh urine sample without centrifugation.
【Performance Indicator】pH is about 7.55 ±0.20 at（25±1）℃
【Important Notice】

1. Avoid contact with skin and eyes;
2. Rinse eyes with the flow of warm water for a few minutes if the product splashes into eyes (ensuring eyelids opened). Seek medical care immediately if pain or swelling happens. Rinse with soap water and water once skin contact occurs.

【Bibliography】

1. Ferris JA Comparison and Standards of Urine Microscopy Lab Med 14:659-662,1983
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

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Please read the User Manual	Stored at	Manufacturer
IVD In vitro diagnostic reagent	LOT Lot Number	Valid until
CE Conform to EU in vitro diagnostic medical equipment directives98/79/EC		
EC REP Authorised Representative	REF Catalogue number	

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