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everall7

VILLACRYL STC HOT

Heat-curing acrylic resin
for crowns and bridges veneering

EN | PL | RU | UA

Everall7 Sp. z o.o.
ul. Augustowska 14 | 02-981 Warszawa | Polska
T: +48 22 858 82 72 | biuro@everall7.pl
www.everall7.pl

CE 1639

VILLACRYL STC HOT

INSTRUCTIONS FOR USE - EN

Heat-curing acrylic resin for crowns and bridges veneering

Product characteristics: heat-cured acrylic resin for crowns and bridges facing

Classification accordance EN ISO 10477: Type 2, Class 1

Indications for use: acrylic faced crowns and bridges

Clinical benefits:

Restoration of biting and chewing functions lost as a result of partial or complete loss of teeth.

Target group or target groups of patients and intended users of the devices

Users: Professionals in the dental technology industry - certified dental technicians producing crowns and bridges.

Patients: People of different ages and sexes, with partial or complete lack of teeth, eligible for use of acrylic crowns or bridges

Contraindications: Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Side effects: Possible irritation and inflammation the oral cavity caused by fungi and bacteria adhering to the crown or bridge in the absence of appropriate hygiene.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

Dental technicians

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

Patients

Irritation, inflammation and allergic reactions caused by residual monomer washed out of the crown.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

Colour shades

Dentine: A1; A2; A3; A3,5; A4; B1; B2; C2; C4; D2

Enamel: E1; E2; E3; E4

Assortment

REF	Set
V210A1Z01	VILLACRYL STC HOT color A1 (80 g powder + 40 ml liquid)
V210A2Z02	VILLACRYL STC HOT color A2 (80 g powder + 40 ml liquid)
V210A3Z03	VILLACRYL STC HOT color A3 (80 g powder + 40 ml liquid)
V210A35Z04	VILLACRYL STC HOT color A3,5 (80 g powder + 40 ml liquid)
V210A4Z05	VILLACRYL STC HOT color A4 (80 g powder + 40 ml liquid)
V210B1Z06	VILLACRYL STC HOT color B1 (80 g powder + 40 ml liquid)
V210B2Z07	VILLACRYL STC HOT color B2 (80 g powder + 40 ml liquid)
V210C2Z08	VILLACRYL STC HOT color C2 (80 g powder + 40 ml liquid)
V210C4Z09	VILLACRYL STC HOT color C4 (80 g powder + 40 ml liquid)
V210D2Z12	VILLACRYL STC HOT color D2 (80 g powder + 40 ml liquid)
V210Z13	VILLACRYL STC HOT dentine powder 40g color A2, A3, A3,5, A4, B2, C2 + dentine powder 20g color A1, B1, C4, D2 + enamel powder 20g color E1, E2, E3, E4 + liquid 200ml
	Liquid
V210L02	VILLACRYL STC HOT liquid 200ml

Chemical composition

Powder: acrylic polymers, benzoyl peroxide, pigments

Liquid: methyl methacrylate, cross-linking agent

Compatible products

For use with metals and dental alloys.

Method of use

Attention: In the purpose of the proper selection of tooth shade color it is recommended using the table attached to the instruction. The table is only a proposal of application color blends.
Color should be selected individually and in dubious cases the own sample should be make.

Preliminary works

The metallic structure of crown or bridge shall be made by traditional method. The chewing surfaces shall be left metallic. The part designed for facing shall have retention elements (e.g. balls).

When structure is polished, places designed for facing shall be sand blasted with the jet of sand (aluminum oxide), grain size 250 µm, under pressure 4-6 bar and masked using product Villacryl **Opaker constituting simultaneously the first layer of color reconstructed tissue. It is possible to mix powders of opaker creating individually first layer of color. The prepared structure is ready for working** model and model with wax the shape of the future facing. Wax shall be removed for acrylate by traditional method. The metal part of structure shall be embedded in cl[illegible]m (Stodent II is recommended) in the lower part of polymerization box so that the face[illegible]mains exposed. The gypsum surface shall be isolated with gypsum-gypsum insulation after[illegible]ng. Then place **the second part of the box, the so called "backing", and it shall be also fi[illegible]ss II gypsum.** After 20 minutes, when gypsum is bounded to the proper hardness, pol[illegible] box shall be opened and wax shall be scalded with boiling water. The form shall be [illegible] gypsum-acrylate isolator (e.g. Izo-Sol). When the polymerization box is cooled to [illegible]mperature, ca **23 °C, the acrylic dough can be applied.**

Recommended following mixing proportions: 2,4 g of power and 1 g o[illegible]

Dosing and mixing

It is recommended to use gloves from polyethylene HDPE for protection o[illegible]

Prepare two chemically resistant containers (glass, polyethylene, etc.) [illegible] powder with **liquid. Pour proper amount of monomer to the first container, add powder** [illegible] color till saturation is obtained, and thoroughly mix. It is possible to mix [illegible] dentine creating individually the appropriate shade of color. After 2 minutes pour t[illegible] container the proper amount of monomer, add powder (enamel) and thoroughly mix [illegible]. Containers shall be accurately covered. The acrylic paste binding level shall be controlled after short period of time.

Application

Paste can be put to the mould after ca 8-10 minutes from mixing powder with liquid. Thoroughly mix the dentine part, squeezing with hands protected by gloves from HDPE, place in mould and cover with polyethylene foil. Close the polymerization box and place under the press. Carefully press till the box edges are completely folded. Take the box from the press, remove the foil, cut **off obliquely with sharp knife dentine part toward the secant edge. Apply the prepared previously** enamel layer in place of removed acrylate, cover with foil and press the polymerization box till 3000 kg pressure. Open the box and check correctness of dentine transition into enamel. If everything is OK, it is recommended to remove the excess of the acrylate pressed aside, again assemble the box, this time without foil separator, and again press till 3000 kg pressure. Take the box from the press after 10 minutes and place in polymerization frame.

Polymerization

The polymerization box fastened in frame place in water with temperature 80-90 °C and heat for 30 minutes.

Processing

The polymerization box after cooling to ambient temperature - shall be opened and the prepared gypsum element shall be released. Processing and polishing by standard method.

Safety instructions Warnings

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Precautions

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: The area of product contact with skin wash thoroughly with soap with water.

P501 Dispose of contents/container to an approved waste disposal plant.

Inflammable liquid (contains methyl methacrylate) - flash point +10 °C. Keep away from fire and radiation sources. While working with the material do not: consume meals or smoke tobacco. Avoid breathing vapour of liquid and direct contact of product components or the acrylic dough on the skin and eyes. Work should be carried out in good ventilated premises. Use standard ventilation. In case of contact of the powder particles with eyes the mechanical irritation may occur. In case of contact acrylic powder or liquid with eyes wash it with water and contact physician. The area of contact the product components or wet acrylic dough with skin wash thoroughly with soap water. In case of accidental swallowing contact the physician. Liquid dangerous to the environment. Do not pour into drainage system.

Hygiene, sterilization, disinfection

Product

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

The crowns made of materials

There is a risk of transmission of infection. The crown should be cleaned and disinfected by the dentist before being placed in the patient's mouth for the first time.

For the disinfection of crown made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water. Provide the appropriate information to the dentist.

Information to be provided to patients / prosthesis users:

If any symptoms of an allergic reaction or irritation occur, discontinue using the crown immediately and consult a doctor.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the material after polymerization.

For the hygiene of crowns and bridges, popular and generally available toothpastes, threads, interdental brushes can be used, it is recommended to use a pressure device of the WATERPIK type.

The material is not resistant to long-term effects of alcoholic beverages.

The material may change color or lose its original physical properties when exposed to chemicals or if it is misused.

Report any damage of the material to the dentist. Repairs may be only performed by a dental technician.

If a piece of damaged crown is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures.

Storage and handling

Liquid is sensitive to influence of temperature and UV radiation under the influence of which the uncontrolled reaction of polymerization may occur. Product components shall be stored in original **containers, in dark and good ventilated compartment with temp. 5-25 °C, place not accessible by children. Do not use the product after expiry date specified on the sales packaging and labels of** the product ingredients.

Disposal considerations

Dispose of in accordance with applicable regulations. Carton packaging, label instructions - the possibility of recycling, dispose of as communal not dangerous waste. Powder component and a packaging are not dangerous. Dispose of as communal not dangerous waste. In case of plastic container the possibility of recycling. Liquid and liquid packaging - hazardous waste. Pass to an authorized waste contractor.

Problem prevention (2)

See website: **www.everall7.pl**

Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Instructions in written form, oral or information relayed during Producers presentation result from actual level of dental technique and technology. It should be treated as obligatory information.

Information/s obtained from third persons (due to their different level of authorization) do not release the user from personal check whether the preparation is intended to the designated appli. cation. Due to the fact that the application takes place without the possibility of control from the producer's side; therefore it limits the producers liability to the value of provided material.

We shall be very grateful for all comments related to products, please mail it to our address.

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VILLACRYL OPAKER — INSTRUCTIONS FOR USE - EN

Universal metal masking material
The product characteristics: Acrylic resin for masking metal elements in dentures, crowns, and bridges.
Classification in accordance with EN ISO 10477: Type 2, Class 1
Indications for use: skeletal dentures, crowns, and bridges faced with acrylic
Clinical benefits:
Restoration of biting and chewing functions lost as a result of partial or complete loss of teeth.
Target group or target groups of patients and intended users of the devices
Users: Professionals in the dental technology industry - certified dental technicians producing dentures, crowns and bridges.
Patients: People of different ages and sexes, with partial or complete lack of teeth, eligible for use of denatures, acrylic or bridges
Contraindications: Do not use this product on persons allergic to any of thecomponents of the product. Do not use the product by pregnant women and nursing mothers in the group of dental technicians.
Side effects:
Dental technicians
Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.
Patients
No side effect are known at the moment.
Data necessary to identify the product and the contents of the package
Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.
Package content - see table below
Assortment

REF	Set
V210Z10	Villacryl Opaker 4x7 g powder (the grey, yellow, brown, pink colours) + 2x12 ml liquid + 12 ml primer
V210Z12	Villacryl Opaker 2x7 g powder (the pink colour) + 12 ml liquid + 12 ml primer

Chemical composition
Powder: acrylic polymer, benzoyl peroxide, pigments
Liquid: acrylic monomers, N,N dimethyl-p-toluidine
Primer: silane, isopropyl alcohol
Product compatibility
For use with dental metals and alloys as well as acrylic materials for the preparing of dentures, crowns and dental bridges.
Method of use
Preparation: Cleanse the metal surface of any leftovers of the protective mass or metal oxides using a prosthetic sandblaster. Prepare the deansed surface by sandblasting it with pure aluminium oxide with the granulation of 250 µm and under the pressure of 4 bar. Clean the prepared surface with acetone, apply **Opaker Primer**, and spread it evenly with compressed air (do not use oil compressors). **Opaker Primer** is used for the chemical joining of acryl-metal, and it prevents the creation of the so-called "randkluft."
Mixing: Pour a few drops of the **Opaker Liquid** monomer to a separate container, and add some **Opaker Powder** to saturate it. Stir for about 30 seconds till you get a homogenous consistency. The **Opaker Powder** powders can be mixed with one another in order to get the required colour.
Applying: Spread on the surface of the metal with a thin, opaque layer of the prepared product with the help of a small brush, and wait for about 5 minutes.

Curing: Put the painted metal in a pressure polymerizer with the water at the temperature of 95 °C for 10 minutes under the pressure of 2 bar. Curing the opaque product can be also performed in a stream of water over a flame of a blowpipe. One ought be careful not to burn the product layer. Otherwise it will become chalky, and it will not adhere to the metal surface. If the product iscured in a stream of water, it keeps a permanent in colour; it adheres to metal; and it covers it completely.
Treatment and finishing: After curing remove the Opaker surface with alcohol in order to remove the unpolymerized layer. Material cured in water has a brighter colour than cured in a stream of hot air.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H319 Causes eye irritation.
H335 May cause irritation of the respiratory tract.
H336 May cause drowsiness or dizziness.
Precautions
P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.
P261 Avoid breathing in dust/smoke/gas/mist/vapours/spray.
P271 Use only outdoors or in a well-ventilated area.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 If on skin: Wash thoroughly with soap with water.
P304+P340 IF INHALED: Take the victim to fresh air and ensure conditions for rest in the position allowing for free breathing.
P305+P351+P338 IF IN EYES: In case of contact with the acrylic powder or liquid with eyes, wash them with water and contact a doctor.
Remove contact lenses if present and easy to do - continue rinsing.
P501 Dispose of the contents/container to an approved waste disposal plant.
An Inflammable primer (contains Isopropyl alcohol) - flash point +15 °C.
An Inflammable liquid (contains) - flash point +10 °C.
Keep away from fire and radiation sources.
While working with the material, do not consume meals or smoke tobacco.
Avoid breathing the vapour of the liquid, and the direct contact of the product components or the wet acrylic dough with the skin and eyes.
Work should be carried out in well ventilated premises.
Use workplace ventilation.
Wash thoroughly the areas of contact of the product components with skin with soap and water. In case of accidental consumption, contact the physician. Do not pour into the drainage system.
Hygiene, sterilization, disinfection
Product
The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.
Denture framework / metal element of the crown masked by material
It does not require sterilization or disinfection before applying the correct acrylic layer of the restoration.
Information to be provided to patients / prosthesis users:
Not applicable
Storage and handling
Product components are to be stored in original containers, in a dark and well ventilated room with temp. 5-25 °C, in a place not accessible to children. Do not use the product after the expiry date. Unused material and packaging are dangerous for the environment. Remove according to local regulations.
Disposal considerations
Dispose of in accordance with applicable regulations. Carton packaging, labels, instructions - the possibility of recycling; dispose of as communal waste not constituting a danger. The powder component and powder packaging are not dangerous waste.

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Dispose of as communal, not dangerous waste. In case of packaging- the possibility of recycling. The liquid and liquid packaging - hazardous waste. Pass on to an authorized waste contractor. The primer and primer packaging - hazardous waste. Pass on to an authorized waste contractor.
Problem prevention (2)
See website: **www.everall7.pl**
Additional information
Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Instructions in the written and oral form, and information relayed during a presentation by the Producer result from the current level of dental technique and technology. They ought to be treated as applicable information. The information obtained from third persons (due to their iidifferent level of authorisation) does not relase the user from checking personally whether the product is intended for the designated application. Due to the fact that the application takes place without the possibility of control on the producer's part, it limits the producer's liability to the value of provided material.

We kindly request all the comments connected with our products to be sent to our address.

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