MEDRAD® Twist & Go

Sterile High Pressure Connector Tubing Products
Produits de tubulure de raccord haute pression stérile
Tubo conector de alta presión estéril

Instructions for Use / Mode d'emploi / Instrucciones de Uso

Catalog Number Référence catalogue Número de Catálogo	Product Name Nom du Produit Nombre del Producto
TAG C50 HPCT	Sterile High Pressure Connector Tubing, 50 cm (20 in) Tubulure de raccord haute pression stérile, 50 cm Tubo conector de alta presión estéril de 50 cm (20")
TAG C75 HPCT	Sterile High Pressure Connector Tubing, 75 cm (30 in) Tubulure de raccord haute pression stérile, 75 cm Tubo conector de alta presión estéril de 75 cm (30")
TAG C120 HPCT	Sterile High Pressure Connector Tubing, 120 cm (48 in) Tubulure de raccord haute pression stérile, 120 cm Tubo conector de alta presión estéril de 120 cm (48")
TAG C150 HPCT	Sterile High Pressure Connector Tubing, 150 cm (60 in) Tubulure de raccord haute pression stérile, 150 cm Tubo conector de alta presión estéril de 150 cm (60")

©2013-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, MEDRAD Twist & Go, and Twist & Go are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the MEDRAD® Twist & Go sterile high pressure connector tubing products can be found on page 2 of these Instructions for Use.

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Introduction: Read the information contained in these instructions for use. Understanding the information will assist you in using the ***MEDRAD® Twist & Go*** (Twist & Go) High Pressure Connector Tubing properly.

Important Safety Notice: The ***Twist & Go*** High Pressure Connector Tubing is intended to be used by qualified medical professionals with proper training and experience in angiographic procedures and the use of injectors from Bayer.
As with any injector, operator vigilance is required when using an Injection System. Visual inspections and good practice during set-up and use and compliance with the operation manual remain the primary methods to limit the risk of air being present during fluid delivery.

Indications for Use: The ***Twist & Go*** High Pressure Connector Tubing is intended to be used in the delivery of contrast media and common flushing solutions. The device is indicated for single use only with injectors from Bayer with ***Twist & Go*** Syringes.

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⚠ WARNING

Air Embolism Hazard. Serious patient injury or death may result.

- Ensure that one operator is designated as being responsible for purging the High Pressure Connector Tubing. Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.
- Ensure that patient is not connected while purging air from High Pressure Connector Tubing.
- Remove all air from the High Pressure Connector Tubing.
- Use caution when removing air from the syringe or tubing connections. Component damage may result from the use of tools during air removal.
- Securely tighten the high pressure connector tubing to syringe tip connection. Component damage may result from the use of tools during tightening.
- Do not use tubing to manipulate the injector position.
- Avoid contact of and damage to syringe tip, tubing and connections when manipulating injector position.

⚠ WARNING

Bloodborne Contamination Hazard. Serious patient and/or worker injury or death may result.

- Use caution when removing air from the syringe or tubing connections. Component damage may result from the use of tools during air removal.
- Do not use tools to securely tighten the high pressure connector tubing to syringe tip connection.
- Avoid contact of and damage to syringe tip, tubing and connections when manipulating injector position.
- For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess, or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.
- Do not exceed 1200 psi (8273 kPa).
- Do not use tubing to manipulate injector position.

Biological Contamination Hazard. Serious patient and/or worker injury or death may result.

- Properly discard disposables after use if contamination may have occurred during setup or use.

⚠ CAUTION

Environmental Contamination Hazard. Patient and/or worker injury may result.

- Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.
- Follow sterile technique principles, specifically, maintain sterility of the syringe tip and high pressure connector tubing.

INSTALLING, FILLING AND PURGING A HIGH PRESSURE CONNECTOR TUBING

1. Remove the dust cap from the syringe tip if attached.
2. Insert the ***Twist & Go*** High Pressure Connector Tubing onto the ***Twist & Go*** syringe tip.
3. Turn the ***Twist & Go*** (A) nut clockwise to securely tighten the high pressure connector tubing to the syringe tip (B).
4. Purge all air from the high pressure connector tubing. Use the injector to push contrast and flushing solutions out until all of the air has been removed from the High Pressure Connector Tubing. Gentle tapping at any connection point may be needed to facilitate air removal.

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A glossary of the symbols used on the MEDRAD® Avanta sterile disposable products can be found on page 2 of these Instructions for Use.

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Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

Vincent Medical (Dong Guan) Mfg. Co., Ltd.
45-46 ShaBu Industrial Zone
Qiao Long District
Tang Xia Town, Dong Guan City
Guangdong, China 523730

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht
The Netherlands

LATEX

CE 2797

DN-258026 Rev. A February 25, 2019

English	5	Français	6
Deutsch	7	Italiano	8
Español (LatAm)	9	Русский	10
Svenska	11	日本語	12
Português (Brasileiro)	13	Nederlands	14
Hrvatski	15	Türkçe	16
Dansk	17	Български	18
Suomi	20	Ελληνικά	21
Magyar	22	Srpski	23
Česky	24	Norsk	25
Polski	26	Română	27
Slovenský	28	Español (España)	29
한국어	30	Slovenščina	31
Português (Europeu)	32		

© 2010-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, MEDRAD Avanta, Avanta, MEDRAD FluiDots, and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the sterile Disposable Sets properly.

Important Safety Notice: The MEDRAD® Avanta Sterile Disposable Sets are intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD® Avanta Fluid Management Injection System and angiographic procedures and techniques.

Indications for Use: The MEDRAD® Avanta Sterile Disposable Sets are specifically intended for use in the x-ray angiography environment. They are designed to administer intravascular radio-opaque contrast compounds, and common flushing agents at various volumes and flow rates into humans for use in diagnostic and interventional angiographic procedures performed in cardiology, radiology, and vascular surgery.

Contraindications: For use only with the MEDRAD® Avanta Fluid Management Injection System. The MEDRAD® Avanta Sterile Disposable Sets are not intended to be used for chemotherapy and are not intended to administer fluids other than intravascular contrast agents and common flushing solutions.

Restricted Sale: Rx Only.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

⚠ **For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only.** Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

⚠ **Air embolization can cause death or serious injury to the patient.** Do not connect a patient to the injector, or attempt an injection until all air has been removed from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators (where applicable) to verify that fluid is present and to reduce the chance of air embolization during an injection.

⚠ **Do not attempt to aspirate fluid from the waste port stopcock when open to the PIV.** Doing so may increase the risk of inducing air into the fluid path.

⚠ **Do not attempt to aspirate fluid from the hemodynamic port.** Doing so may increase the risk of inducing air into the fluid path.

⚠ **Before connecting the Single-Patient Sterile Disposable Set (SPAT) to the patient, ensure all stopcocks and open ports are closed to air and all air has been removed from the fluid path before injection.** Improper manipulation of the waste port stopcock as well as any fluid port left open to air may increase the risk of introducing air into the fluid path.

⚠ **Biological contamination can result from reusing disposable items or failure to follow sterile technique during setup or use.** Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

⚠ **Do not use if sterile package is opened or damaged.** Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

⚠ **Patient or operator injury may result from mishandling of the sharp spike.** Use care in handling and inserting the spike into the contrast bottle and the saline bag/bottle.

⚠ **Storage of filled syringes can promote bacterial growth.** Syringes from Bayer are designed to be filled just prior to the procedure. Discard unused filled syringes.

⚠ **Patient injury could result if syringe is not properly engaged.** Improper engagement may cause air embolization or under-volume delivery.

⚠ **Patient or operator injury may result from contrast media leaks or ruptures to the syringe or connector tubing.** Use of greater pressure or occlusions in the fluid path may result in leaks or ruptures to the syringe or connector tubing. Ensure that the fluid path is open; do not exceed appropriate pressure.

⚠ **Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe.** Do not remove the plunger to fill the syringe.

⚠ **Inadvertent aspiration can occur if the Multi-Patient disposable contrast tubing set is kinked or blocked during a syringe fill process.** Verify multi-patient contrast tubing is free of kinks before attempting to fill the syringe.

⚠ **Patient injury can result if the patient is connected during any contrast or saline purge.** Ensure the patient is disconnected before any purge.

⚠ **Blood can contaminate the multi-patient portion of the disposable set.** Ensure that blood does not appear past the pressure isolation valve in the single patient disposable set.

⚠ **Ensure vent cap on spike(s) are open if using bottles.** If closed, a vacuum may be created and air may be introduced.

⚠ **Do not exceed the five patients per day usage requirement for the Multi-patient Disposable Set (MPAT).** This is the usage approved by appropriate regulatory bodies.

⚠ **Use of non-Bayer supplied disposables, including administration sets and additions to administration sets, such as but not limited to bleed back control devices and pressure transducers may cause patient injury if not properly connected or flushed. These devices must be compatible with your system.** Refer to manufacturer's instructions for proper use of these devices.

⚠ **Use of extension tubing connected to the Single Patient Disposable Set may cause patient injury and will degrade system performance.** Do not connect extension tubing to the Single Patient Disposable Set.

Cautions

⚠ **Component damage may occur if disposables are not installed properly.** Assure all connections are secure. This will help minimize leaks, disconnections, air introduction, and component damage. Do not use tools to over tighten connections or to assist in removal of disposables.


⚠ **Do not use methyl alcohol on disposable sets.** Use of methyl alcohol will degrade the integrity of the disposable sets.

⚠ **All contrast and saline media should be used in accordance with the manufacturer's indication for use.**

NOTE: To perform the following procedures, the operator must select the Case Setup screens on the DCU. Refer to the injector operation manual for additional information.

Multi-Patient Sterile Disposable Set (MPAT) Installation (usable for up to five patients). The Multi-Patient Disposables Set consists of the following items. (Also see MPAT Setup illustration):

Vented Contrast Spike	2nd HPCT w/ male luer
Contrast Drip Chamber	Saline Spike
Low Pressure Contrast Tube (LPCT)	Saline Drip Chamber
High Pressure Contrast Tube (HPCT)	Low Pressure Saline Tube (LPST)
Contrast Flow Control Valve (CFCV)	Post Drip Chamber Saline Tube
150 mL Multi-Patient Syringe	

General Setup Instructions

Be sure that the patient is disconnected and the injector is on.

1. Install syringe. (See Install Multi-Patient Syringe below).
2. Install saline tubing (See Install Saline Tubing below).
3. Install contrast tubing. (See Install Contrast Tubing below).
4. Rotate head to upright position, fill syringe and remove dust caps.
5. Purge multi-patient disposables. Catch the discharged fluid in a suitable container and tap disposable connection points and pressure jacket to remove air.
6. Rotate head downward.
7. Attach MEDRAD® Avanta Sterile Hand Controller, and maintain sterility, if applicable.

NOTE: Refer to Operation Manual for MEDRAD® Avanta Fluid Management Injection System for further information.

After step 7 in the General Setup Instructions is complete, go to SPAT Installation.

Install Multi-Patient Syringe

NOTE: Ensure patient is not connected.

Select Multi-Patient Syringe from Case Setup.

1. Retract piston.
2. Insert syringe into pressure jacket and raise drop front arms.
3. Engage the piston to the plunger.

Install Saline Tubing

NOTE: Ensure patient is not connected.

Select Multi-Patient Disposables from Case Setup, then select Install Saline Tubing.

1. Install drip chamber.
2. Route tubing through pump and close door.
3. Route tubing around peg.
4. Insert tubing into gross air detector and close door.
5. Spike saline bag.
6. Fill saline drip chamber completely and snap into position before use.

NOTE: To differentiate the saline and contrast lines of the MPAT, the contrast line will be identified as the line connected to the stopcock that is placed in the CFCV.

Install Contrast Tubing

Select Install Contrast Tubing from Case Setup.


1. Install drip chamber.
2. Connect tubing to syringe.
3. Install stopcock.
4. Insert tubing into gross air detector and close door.
5. Spike contrast bottle.
6. Fill contrast drip chamber completely and snap into position before use.

Go to step 4 in General Setup Instructions.

Single-Patient Sterile Disposable Set (SPAT) Installation (must be replaced for each patient use). The Single-Patient Disposable Set consists of the following items (Also see SPAT Setup illustration):

Rotating luer **(A)**	High Pressure Saline Tube w/ male luer connector **(D)**
Stopcock for waste port **(B)**	High Pressure Contrast Tube w/ female luer connector **(E)**
Pressure Isolation Valve (PIV) Assembly / Hemodynamic pressure monitoring port **(C)**	

Single-Patient Sterile Disposable Set (AVA 500 SPAT Angio) is indicated for use when a user attachable transducer is unnecessary and does not contain item (C) pictured below.

Select Single-Patient Disposables from Case Setup.

Be sure that the patient is disconnected.

1. Remove all dustcaps.
2. Connect the SPAT saline tube to MPAT saline tube.
3. Connect the SPAT contrast tube to MPAT contrast tube.

To purge air out of SPAT:

NOTE: Ensure patient is not connected.

When using AVA 500 SPAT L continue with step 4. When using AVA 500 SPAT ANGIO proceed to step 5.

4. Attach transducer to hemodynamic port, verify the transducer stopcock is open, and the waste port stopcock is off to waste port.
 a. Attach 20cc saline syringe to the transducer port.
 b. Hold SPAT parallel to the patient with the PIV and 20cc syringe connection facing upward at a 45° angle.
 c. Forcefully flush through the transducer and PIV assembly until all air is expelled from both components. Flick or tap PIV to facilitate air removal (Repeat as necessary).
 d. Close transducer port stopcock. Remove 20cc syringe.
5. Initiate contrast and saline purge, tap SPAT to MPAT saline connection point and wasteport stopcock to remove air while purging is in progress.
6. Repeat steps above as needed, until all air is purged.
7. Carefully inspect the entire SPAT and MPAT to verify all air is removed before connecting to patient

NOTE: Disconnect the SPAT from the MPAT immediately following each procedure.

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