MEDRAD® VFlow

Disposable Hand Controller
Contrôleur manuel jetable
Controlador manual desechable

Instructions for Use / Mode d'emploi / Instrucciones de Uso

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©2013-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Mark 7 Arterion, VFlow, MEDRAD Mark 7 Arterion, and MEDRAD VFlow are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the MEDRAD® VFlow sterile disposable products can be found on page 2 of these Instructions for Use.

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

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Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the **MEDRAD VFlow** Sterile Hand Controller properly.

Important Safety Notice: The MEDRAD VFlow Sterile Hand Controller is intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD® Mark 7 Arterion Injection System and angiographic procedures and techniques.

Indications for Use: The MEDRAD VFlow Sterile Hand Controller is specifically intended for use with the MEDRAD Arterion Injection System for angiographic studies.

Contraindications: None known.

Restricted Sale: Rx Only.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

⚠ WARNING

Environmental Contamination Hazard - Serious patient or worker injury or death may result.

- Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Cross Contamination Hazard - Serious patient and/or worker injury or death may result.

- For devices labelled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

⚠ CAUTION

Environmental Contamination Hazard - Minor or moderate patient and/ or worker injury may result.

- Properly discard disposables after use, in accordance with hospital hazard waste disposal procedures. Do not reuse disposables.

NOTICE


Electro-Mechanical Hazard - Equipment Damage may result.

- Do not use tools to over tighten connections or to assist in the removal of disposables.

NOTE: Refer to the MEDRAD Mark 7 Arterion injector operation manual for a full set of instructions, warnings and precautions.

The MEDRAD VFlow Hand Controller is a sterile device intended for single patient use. The Hand Controller works for Fixed Flow Rate and Variable Flow Rate injections. When in the Variable Flow Rate injection mode, the flow rate increases incrementally as the Hand Controller plunger (A) is depressed, and decreases as the Hand Controller plunger is released. In the Fixed Flow Rate injection mode, the Hand Controller acts as a start switch, and release of the device stops the injection.


The Hand Controller button (B) is non-functional and will beep from the Injector head and Display Control Unit when depressed.

Hand Controller

Hand Controller Installation

1. Ensure that the injector is in the disarmed state.
2. Open the hand controller pouch using sterile technique.
3. Using sterile technique, remove the hand controller from the pouch.
4. Plug the hand controller into the underside of the Injector Head (see drawing below). Listen for an audible “click” to confirm proper connection.

NOTE: If the hand controller is damaged or does not work properly, discontinue use and discard.

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©2010-2011, 2013- 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Avanta, and MEDRAD Avanta are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the MEDRAD® Avanta sterile disposable products can be found on page 2 of these Instructions for Use.

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
USA
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

English	7	Deutsch	8
Français	9	Italiano	10
Svenska	11	Nederlands	12
Dansk	13	Suomi	14
Português (Brasileiro)	15	Español (España)	16
Norsk	17	日本語	18
Srpski	19	Български	20
Hrvatski	21	Česky	22
Ελληνικά	23	Magyar	24
Polski	25	Română	26
Русский	27	Slovenský	28
Türkçe	29	简体中文	30
Português (Europeu)	31	Español (LatAm)	32
한국어	33	Lietuvių	34
Tiếng Việt	35	Slovenščina	36

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Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the MEDRAD Avanta Sterile Hand Controller properly.

Important Safety Notice: The MEDRAD Avanta Sterile Hand Controller is intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD Avanta Fluid Management Injection System and angiographic procedures and techniques.

Indications for Use: The MEDRAD Avanta Sterile Hand Controller is specifically intended for use with the MEDRAD Avanta Fluid Management Injection System in the x-ray angiography environment.

Contraindications: The MEDRAD Avanta Sterile Hand Controller is only intended to be used with MEDRAD Avanta Fluid Management Injection System.

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⚠ Warnings

- **Air embolization can cause death or serious injury to the patient.** Do not connect a patient to the injector, or attempt an injection until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions of the use of FluiDots indicators (where applicable) to reduce the chance of air embolism.
- **Biological contamination can result from reusing disposable items or failure to follow sterile technique.** Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

NOTE: If the MEDRAD Avanta Sterile Hand Controller is used without the optional MEDRAD Avanta Hand Controller Sterile Sheath, discard the hand controller after each patient.

- **Do not use if sterile package is opened or damaged.** Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.
- **Patient or operator injury may result from contrast media leaks or ruptures to the syringe or connector tubing.** Use of greater pressure or occlusions in the fluid path may result in leaks or ruptures to the syringe or connector tubing. Ensure that the fluid path is open; do not exceed appropriate pressure.

⚠ Caution

- **Component damage may occur if not installed properly.** Assure all connections are secure; do not overtighten. This will help minimize leaks, disconnections, and component damage.


NOTE: Refer to the injector operation manual and the MEDRAD Avanta Hand Controller Sterile Sheath IFU for additional information.

Hand Controller Installation

The MEDRAD Avanta Sterile Hand Controller can be used for up to five patients. Sterility can be maintained between patients by using MEDRAD Avanta Sterile Hand Controller Sheaths (available for purchase).

1. Ensure that the injector is in the disarmed state.
2. Open the hand controller bag using sterile technique.
3. Using sterile technique, remove the hand controller from the bag.
4. Plug the hand controller into the Fluid Control Module – Hand Controller Connector (see drawing below). Listen for an audible “click” to confirm proper connection.

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5. Follow instructions provided in the MEDRAD Avanta Fluid Management Injection System Operations manual for hand controller use.

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