Qwik-Fit Syringe for MR Injectors
Seringue Qwik-Fit Syringe pour Systèmes d'injection IRM
Jeringa Qwik-Fit para inyectores MR

Instructions for Use / Mode d'emploi / Instrucciones de Uso

Injector Injecteur Inyector	Catalog Number Référence catalogue Número de Catálogo	Product Name Nom du Produit Nombre del Producto
MEDRAD® Spectris	**SQK 65VS** **SQK 65VA**	**Disposable MRI Kit with Check Valve /** **Kit IRM jetable avec valve anti-reflux /** **Kit de desechables para RM con válvula antirretorno**
MEDRAD® Spectris Solaris	**SSQK 65/115VS**	**Disposable MRI Kit for 65/115 MR Injector System /** **Kit IRM jetable pour système d'injection IRM 65/115 /** **Kit de desechables para sistema inyector de RM para 65/115**
	SSQK 115	**Saline Syringe & Tubing for MR Injector System /** **Seringue de sérum physiologique et tubulure pour système d'injection IRM /** **Jeringa y tubo de solución salina para sistema inyector de RM**
MEDRAD® Spectris Solaris EP	**SSQK 65/115VS**	**Disposable MRI Kit for 65/115 MR Injector System /** **Kit IRM jetable pour système d'injection IRM 65/115 /** **Kit de desechables para sistema inyector de RM para 65/115**
	SSQK 115	**Saline Syringe & Tubing for MR Injector System /** **Seringue de sérum physiologique et tubulure pour système d'injection IRM /** **Jeringa y tubo de solución salina para sistema inyector de RM**

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

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Bayer, the Bayer Cross, MEDRAD, Spectris Solaris, MEDRAD Spectris Solaris, Spectris Solaris EP, MEDRAD Spectris Solaris EP, FluiDots, MEDRAD FluiDots, and Qwik-Fit Syringe are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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⚠ Warnings

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators (where applicable) to reduce the chance of air embolism.

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Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used Visually inspect contents and package before each use

Patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged. Ensure the alignment marks on the syringe and injector head are properly aligned and the piston and plunger are interlocked. Improper engagement may cause air embolization or under-volume delivery.

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Operator injury could result from using a sharp spike. Use care in handling and inserting the spike into the bottle.

⚠ Warnings

Serious injury may result from syringe failure. Do not retract pistons with connector tubing installed. Retracting the pistons with the connector tubing installed on syringes will create a vacuum in the syringe due to the check valve in the connector tubing. This vacuum may accelerate the plunger rapidly toward the tip of the syringe when it is removed from the injector causing the syringe to break.

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Syringe Installation:

1. Ensure that the piston(s) on the injector head are fully retracted.
2. Remove the new syringe(s) from the package and insert into the injector head Syringe Interface ensuring that the Keyed Wings at the syringe base and openings in the Syringe Interface are aligned.


 Note that the syringe is keyed to properly fit one way.

3. Twist the syringe 1/4 turn clockwise until the syringe clicks, locks, and the alignment marks on the syringe and injector head are aligned.

Syringe Filling:


NOTE: One person should fill and arm the injector. If you must change the injector operator, be certain that the second operator verifies that the syringe is properly filled and that any air is eliminated.

1. Tilt the injector head up.
2. Advance the plunger to full forward position.
3. Remove the protective cover(s) from the syringe luer tip and attach a sterile filling device (spike or Female/Female Adapter (FFA)) onto the tip of the syringe.
4. Transfer the contrast media into the syringe by one of following methods:
 a. **Spike(s)** - The SQK 65S, SQK 65VS, and SSQK 65/ 115VS kits contain two spikes. Use the smaller spike to load contrast media syringe A and the larger spike to load saline syringe B. The syringe and spike included in the SSQK 115 package are intended for use with saline.
 b. **FFA** - The FFA is included in the SQK 65VA kit. Usethe FFA to load contrast from a prefilled syringe. Prepare the prefilled syringe by removing the protective closure. Attach the prefilled syringe to the FFA. The FFA may also be purchased separately. The catalog number for this item is FFA 50.

NOTE: Use a filling device to reduce the volume and size of air bubbles drawn into the syringe during filling. It is more difficult to remove the air bubbles if you use smaller diameter tubes, such as a catheter-over-needle, or a tube longer than 10 in. (25 cm.)

5. Aspirate contrast media into the syringe.
6. The syringe has FluiDots indicators to help you detect the presence of contrast media in the syringe. Verify that the Qwik-Fit Syringe disposable FluiDots indicators are round in the filled portion of the syringe. The rounded shape of the FluiDots indicators varies according to the type of contrast media, but an oblong shape indicates the presence of air. Rounded FluiDots indicators do not indicate the total absence of air bubbles in the syringe tip.

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
1. Remove the filling device and expel air from the syringe.

 Do not hit the syringe to remove air bubbles. Reverse the plunger 3 to 5 ml then expel the remaining air.

2. The MEDRAD Spectris MR Injector uses two syringes and low pressure connector tubing with a 'T' connector. Attach the straight portion to the contrast media syringe A and the extension to the saline syringe B. Do not reverse the piston plunger on syringe A after you have attached the tube.

3. Verify that the tubing is not kinked or obstructed.
4. Use the manual knob to prime the connector tubing and make certain that the air is expelled.


 An absence of flow is an obvious indication of air in the fluid path.

5. Verify that all the air has been expelled from the syringe and fluid path.
6. Attach the connector tube to the vascular entry device.
7. Tilt the injector head down.

Reinstalling a Syringe


If you remove a syringe from the injector and then wish to reinstall it, perform the following steps:

1. Insert the end of the syringe in the horizontal cutouts in the injector head.
2. Advance the piston until it is past the plunger feet and the piston and plunger interlock.
3. Rotate the syringe 1/4 turn clockwise until the syringe locks and alignment marks are positioned.
4. Proceed as normal by aspirating and dislodging any air bubbles.

Removing a Syringe

To detach the syringe when the procedure has been completed, rotate the syringe 1/4 turn counterclockwise and remove it from the injector head. Discard all used components.

NOTE: Do not use unless the plunger feet are properly engaged with the piston.

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©2012-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Twist & Go, MEDRAD Twist & Go, Mark 7 Arterion, MEDRAD Mark 7 Arterion, FluiDots, and MEDRAD FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

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A glossary of the symbols used on the MEDRAD® Twist & Go sterile disposable products can be found on page 2 of these Instructions for Use.

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Introduction: Read the information contained in these instructions for use. Understanding the information will assist you in using the MEDRAD *Twist & Go* syringe properly.

Important Safety Notice: The MEDRAD *Twist & Go* syringe is intended to be used by qualified medical professionals with proper training and experience in angiographic procedures and the use of the *MEDRAD® Mark 7 Arterion Injection System*. The operations manual is intended as instructions on the proper use of the *MEDRAD Mark 7 Arterion* injector and syringe.

As with any injector, operator vigilance is required when using the MEDRAD Mark 7 Arterion injector/syringe. Visual inspections and good practice during set-up and use and compliance with the operator's manual remain the primary methods to limit the risk of air being present during fluid delivery.

Indications for Use: The MEDRAD *Twist & Go* Syringe, Quick Fill Tube and other disposables supplied from Bayer are specifically intended for single use only with the MEDRAD *Mark 7 Arterion Injection System* for angiographic studies.

Contraindications: This device is not intended to be used for chemotherapy and is not intended to administer fluids other than intravascular contrast agents and common flushing solutions.

Restricted Sale: Rx Only

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred)..

⚠ WARNING

Air Embolism Hazard. Serious patient injury or death may result.

- Ensure that one operator is designated as being responsible for filling and refilling the syringe. Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.
- Ensure that patient is not connected while purging air from syringe, or engaging or advancing plunger.
- Orient the Injector Head to the Purge (upright) position during filling of syringe and purging of air.
- Purge all air out of syringe and any and all disposables after filling.
- Tap syringe after filling to facilitate air removal.
- Verify that the *MEDRAD®* FluiDots indicators are rounded to ensure that fluid is present in the syringe.

⚠ WARNING

Bloodborne Contamination Hazard. Serious patient and/or worker injury or death may result.

- Use caution when removing air from the syringe. Component damage may result from the use of tools during air removal.
- Press the **End Case** button on the Display Control Unit. Select **Yes** to acknowledge that you want to end the case, and that the patient has been disconnected from the system.
 - Alternatively, turn the Manual Knob counter-clockwise to retract the syringe plunger 2ml.
- Do not exceed 1200 psi (8273 kPa)

Biological Contamination Hazard. Serious patient and/or worker injury or death may result.

- Properly discard disposables after use or if contamination may have occurred during setup or use.

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

⚠ CAUTION

Environmental Contamination Hazard. Patient and/or worker injury may result.

- **Do not use if sterile package is opened or damaged.** Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.
- Follow sterile technique principles, specifically, maintain sterility of the syringe tip and plunger and syringe barrel internal surface, and Quick Fill Tube.
- Do not scrape off dried, potentially contaminated contrast into recesses during syringe installation.
- Do not reuse disposables.


Bacterial Contamination Hazard. Patient and/or worker injury may result.

- Use filled syringes immediately.
- Do not store filled syringes for later use.
- Discard previously filled unused syringes.

INSTALLING A SYRINGE

Prior to installing a syringe, ensure that the system is on and the pressure jacket is installed.

1. Ensure the syringe piston is fully retracted. To retract the piston, press the **Enable** button (A), and then press reverse arrows (B) on the **Fill Strip** (C).

NOTE: Your finger placement on the Fill Strip determines the speed at which the syringe retracts or advances. Move your finger away from the **Enable** button to increase the speed.

2. Open the syringe package and remove the syringe.
3. Insert the syringe in the Pressure Jacket. Install the syringe with the raised syringe alignment key (D) aligned with the triangle (E) on the Pressure Jacket.
4. Maintain syringe tip sterility, and raise and completely close the drop front (F).
5. On the Injector Head, press the **Enable** button, and then press the forward arrows (G) on the **Fill Strip** to fully advance the plunger in the syringe.


See Next Page for Filling, Purging, and Removing Syringe

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1. Orient the Injector Head to the Purge (upright) position during filling of syringe and purging of air.
2. Remove the Quick Fill Tube from the syringe package.
3. Remove the dustcap from the syringe tip and set aside maintaining sterility.
4. Attach the short end of the Quick Fill Tube to the syringe tip.
5. Insert the long end of the Quick Fill Tube into the fluid source (usually contrast media). Raise the contrast bottle until the Quick Fill Tube is fully inserted into the contrast.
 NOTE: Use a Quick Fill Tube supplied by Bayer or equivalent device to reduce the volume and size of air bubbles drawn into the syringe during filling. It is more difficult to remove the air bubbles if you use smaller diameter tubes or a tube longer than 10 in. (25 cm.)
6. On the Injector Head, press the **Enable** button (A), and then press the reverse arrows (B) on the **Fill Strip** (C) until the system fills the syringe with the desired contrast volume.
 - Alternatively, on the Injector Head, press the **Enable** button and then press and release the **Auto-Fill** button (I). The MEDRAD *Mark 7 Arterion* fills the syringe with the preconfigured contrast volume at the preconfigured speed. The volume and speed are configured from the Display Control Unit Options tab.
7. If required, use your free hand to gently tap the base of the pressure jacket to facilitate migration of air bubbles remaining within the syringe and on the syringe plunger to the syringe tip.
8. Carefully observe the FluiDots indicators to ensure that fluid is present in the syringe. Verify that the FluiDots indicators are round in the filled portion of the syringe. The rounded shape of the FluiDots indicators varies according to the type of contrast media, but an oblong shape indicates the presence of air. Rounded FluiDots indicators do not indicate the total absence of air bubbles in the syringe tip.

Empty Syringe Filled Syringe

9. Purge all air out of syringe after filling. Turn the Manual Knob (J) clockwise to remove air from the syringe.
10. Remove the Quick Fill Tube from the syringe tip.
11. Visually confirm that all air bubbles have been removed from the syringe. Tap syringe after filling to facilitate air removal.
12. The system is ready to accept installation of the MEDRAD *Twist & Go* High Pressure Connector Tubing (TAGHPCT). If you are not installing the TAGHPCT at this time, replace the sterile dustcap. See operations manual for TAGHPCT specifications, connecting and purging.

REMOVING A SYRINGE

1. Disconnect the disposable tubing set from the vascular entry device, such as a catheter or sheath. The disposable tubing set does not need to be disconnected from the syringe.
2. Press the **End Case** button on the Display Control Unit. Select **Yes** to acknowledge that you want to end the case, and that the patient has been disconnected from the system.
 - Alternatively, turn the Manual Knob counter-clockwise to retract the syringe plunger 2ml.
3. Open the drop front.
4. Rotate the syringe 1/4 turn clockwise and gently pull the syringe out of the Pressure Jacket. Discard the syringe with disposable tubing set into a bio-hazard container.
 NOTE: Once the syringe is removed from the injector and the Injector Head is rotated into the Purge position, the injector sounds an audible beep three times, and the piston automatically retracts to the home position. Auto Retract must be enabled for this feature to function.

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Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
USA
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

© 2011-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Mark V ProVis, MEDRAD Mark V ProVis, FluiDots, MEDRAD FluiDots, and Quik-Fit Syringe are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

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A glossary of the symbols used on the listed sterile disposable products can be found on page 2 of these Instructions for Use.

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ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for use: The contents of this package are intended to be used in the delivery of contrast media. These devices are indicated for single-use only with injectors from Bayer.

Contraindications: These devices are not intended for multiple-use, drug infusion, chemotherapy, or any other use for which the device is not indicated.

Restricted Sale: Federal (USA) law restricts these devices for sale by or on the order of a physician.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

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Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators to reduce the chance of air embolism.

Biological contamination can result from reusing disposable items or failure to follow aseptic technique. Properly discard disposable items after use.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Bacterial contamination can occur if syringes are used to store contrast media. Use loaded syringes immediately. Do not store loaded syringes for later use. Discard unused loaded syringes.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or Operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures to the connector tube.

⚠ Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.

Refer to the Injector Operation Manual for further instructions.

MEDRAD Mark IV, MEDRAD Mark V, MEDRAD Mark V Plus, MEDRAD Mark V Provis, or MEDRAD CT-100/200/300 Series Injectors:

1. Ensure that the piston on the head is fully retracted. To retract the piston(s), press the Reverse Load Switch on the injector head.
2. Set the Mechanical (safety) Stop to zero in the full forward position (applies to MEDRAD Mark IV and MEDRAD CT-100 Series Injectors only). Position the empty Pressure Jacket in the top of the turret. If there is a syringe in the top pressure jacket, remove and discard.
3. Remove the syringe from its package and slide it in the top Pressure Jacket (the flat side of the flange must face upward). This positions the nut with its cutout facing up for easier hub attachment.
4. Fully rotate the turret to place the syringe plunger in front of the piston.

Note: The turret must be fully rotated to the closed position to ensure piston/syringe plunger alignment to prevent contrast leakage.

The following are general guidelines only. For additional instructions, refer to the Operation Manual for your Liebel Flarsheim injector.

Liebel Flarsheim 6000 and 6000 CT Injectors:

1. Ensure that the piston on the head is fully retracted.
2. Unlock pressure jacket and rotate downward.
3. Remove syringe from its package and insert syringe into pressure jacket, making sure to position alignment key with pressure jacket slot.
4. Rotate pressure jacket plate upward and latch by turning knob fully clockwise.
5. Advance piston to engage syringe plunger.

Liebel Flarsheim Angiomat Injector:


1. Press the LOAD button on the console. Press the LOAD button on the head to fully retract the plunger.
2. Set the mechanical stop at zero.
3. Open the turret door so the pressure jacket is ready to accept a syringe. If there is a syringe in the pressure jacket, remove and discard.
4. Remove syringe from its package and insert syringe into pressure jacket.
5. Close and secure the turret.
6. Press the UNLOAD button to advance the piston forward.

Syringe Filling:

Note: One person should fill and arm the injector. If you must change the injector operator, be certain that the second operator verifies that the syringe is properly filled and that any air is eliminated.

1. Tilt the injector head up.
2. Press and hold the forward motion control until the piston is fully advanced.
3. Remove the protective covers from the syringe luer tip and attach a sterile filling device onto the tip of the syringe (Quick Fill Tube or other filling device).
4. Transfer the contrast media into the syringe using a Quick Fill Tube. Raise the contrast bottle until the Quick Fill Tube is fully inserted. The end of the Quick Fill Tube should be near the bottom of the contrast bottle.

 Note: Use a Quick Fill Tube or equivalent device to reduce the volume and size of air bubbles drawn into the syringe during filling. It is more difficult to remove the air bubbles if you use smaller diameter tubes, such as a catheter-over-needle, or a tube longer than 10 in. (25cm.)
5. Aspirate contrast media into the syringe.
6. Press and hold the REVERSE piston motion key until the syringe is filled with the volume of contrast media desired.
7. Use the FluiDots indicators to help you detect the presence of contrast media in the syringe. Verify that the FluiDots indicators are round in the filled portion of the syringe. The rounded shape of the FluiDots indicators varies according to the type of contrast media, but an oblong shape indicates the presence of air. Rounded FluiDots indicators do not indicate the total absence of air bubbles in the syringe tip

Empty Syringe *Filled Syringe*

Expelling Air and Attaching the Connector Tube:

After you fill the syringe, you must expel all air from the syringe and then attach the connector tubing. To do so:

1. Remove the filling device and expel air from the syringe by manually advancing the piston so that the contrast media is in the tip of the syringe.

 Gently tap the base of the pressure jacket to dislodge the air bubbles. Expel remaining air.
2. Attach the connector tube to the syringe:
 a. **FasTurn** - Place the connector tube hub in the syringe nut's cutout. Turn the syringe nut clockwise approximately 2 1/4 turns until secure while holding the hub stationary.
 b. **Push-N-Turn** - Place the connector tube hub in the syringe nut's cutout. Push the syringe nut towards the injector head until it snaps; turn the nut clockwise approximately 1 turn until secure.

 NOTE: Position the hubs so that when the syringe nut is secure, the wings are across the cutout, not with one wing in the cutout.
3. Verify that the tubing is not kinked or obstructed.
4. Use the manual knob to advance the syringe plunger to provide a slow flow of contrast media through the connector tubing.

 An absence of flow is an obvious indication of air in the fluid path
5. Prime the connector tubing and make certain that the air is expelled.
6. Verify that all of the air has been expelled from the syringe and fluid path.
7. Attach the connector tube to the vascular entry device.
8. Tilt the injector head down.

Removing a Syringe:

To remove the syringe when the procedure has been completed, fully retract the piston and rotate the turret counter-clockwise. Remove the syringe from the pressure jacket. Discard all used components.

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