GYNEMED

Bromelain in Dulbecco´s PBS

Material included in the package:
Bromelain in Dulbecco´s PBS
Product code: 4 GM 501BROM10
- 1 x10 ml Bromelain in Dulbecco´s PBS

Material not included in the package:
- Incubator (no CO_2)
- Sperm Cup
- LAF bench (ISO 5 environment)
- Microscope
- Test tubes

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for Bromelain in Dulbecco´s PBS is available upon request and can also be downloaded from our website.
- Bromelain in Dulbecco´s PBS is manufactured and tested according to the following specifications:

pH (at 37°C)	7.20-7.50
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.5
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %)	≥80
Sperm Survival test survival after 4 hours exposure of density selected spermatozoa to the test medium in %	≥80

EN: Intended use/Intended users:
- Bromelain in Dulbecco's PBS is designed for liquefaction of viscous semen samples prior to semen analysis and preparation for further IVF treatment.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- Bromelain solved in Dulbecco's PBS
- Bromelain in Dulbecco's PBS is an enzyme preparation and has a specific activity of 10 IU/ml.

Instructions for use:
- Warm the Bromelain in Dulbecco's PBS to 37°C.
- Dilute the semen sample with the same volume of Bromelain (for highly viscous, tenacious ejaculate we recommend to previously grind the ejaculate roughly!)
- Swing the semen solution carefully
- Incubate the ejaculate for approx. 10 minutes at 37°C.
- Use the liquefied semen sample for evaluation. Attention: To calculate the sperm concentration the dilution 1+1 (1:2) of semen with Bromelaine must be accounted.

Precautions and warnings:
- Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination.
- Only for the intended use.
- The long term safety of ICSI, including the use of Bromelain, on children born following this procedure is unknown.
- Since the device in its intended use does not come into direct contact with the patient, toxicity as well as sensitization is highly unlikely. However, caution may be advised with ejaculates from patients with known latex and associated food allergies, since Bromelain shares possibly antigen epitopes (Levy and Leynadier, 1997).
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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Storage instructions and stability:
- The shelf life is 9 month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be resterilized after opening.
- Do not use after expiry date.

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	Consult instructions of use
	Temperature limitation
LOT	Lot number
	Use by
REF	Catalogue number
	Do not use if package is damaged
STERILE A	Sterilized using aseptic processing techniques (filtration)
	Keep away from sunlight
	Do not resterilize

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GYNEMED

GM501 Gradient 100%

Material included in the package:
Product code: 4 GM 501G-100-50
- 1 x 50 ml GM501 Gradient 100 %

Product code: 4 GM 501G-100-100
- 1 x 100 ml GM501 Gradient 100 %

Product code: 4 GM 501G-100-250
- 1 x 250 ml GM501 Gradient 100 %

Product code: 4 GM 501G-100-500
- 1 x 500 ml GM501 Gradient 100 %

Material not included in the package:
- Incubator or water bath at 37°C (optional)
- LAF bench (ISO 5 environment)
- Test tubes
- 3 cc syringes with 21 G needle
- Centrifuge (must be able to operate for up to 30 minutes at 400 g)
- Sperm washing medium (e.g. GM501 SpermActive)

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Gradient 100% is available upon request and can also be downloaded from our website.
- GM501 Gradient 100% is manufactured according to the following specifications:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Osmolality (mOsm/kg)	300-330
Density (g/ml)	1.1150-1.1250
Viscosity (cP)	< 1.75
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.5
Sperm Survival test survival after 4 hours exposure of density selected spermatozoa to the test medium in % survival after 24 hours exposure of density selected spermatozoa to the test medium in %	≥ 80 ≥ 75
Mouse Embryo Assay (MEA)	not MEA tested

Calculation of g-forces:
- The g-force of your centrifuge can be calculated using this formula:

g = 1.118 x r x rpm2 or rpm = Square root {g/ (1.118 x r)}
r = radius of centrifuge in mm
rpm = rotations per minute/1000
Example 1
r = 150 mm
rpm = 1200 rotations per minute
g = 1.118 x 150 x 1.44 = 242 g
Example 2
r = 150 mm
g = 300 g
rpm = SQR {300/ (1.118 x 150)} = 1.33
rpm = 1330 rotations per minute

EN: Intended use/Intended users:
- GM501 Gradient 100 % is an isotonic solution for semen preparation with a density of approximately 1.12 g/ml. GM501 Gradient 100 % can be used in combination with IUI, IVF and ICSI.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- GM501 Gradient 100 % consists of silane-coated colloidal silica particles suspended in HEPES-buffered medium.

Instructions for preparation of gradients:
- Mix the density gradient bottles by 5 bottle inversions before use.
- We advice to produce a 2-phasic gradient system (45% and 90%). You may prefer a different mixing ratio (e.g. 40% and 80%) or a multi layer gradient (45%-70%-90%).
- Mix 9 parts of GM501 Gradient 100 % with 1 parts of washing medium to produce the 90% gradient.
- Mix 4.5 parts of GM501 Gradient 100 % with 5.5 parts of washing medium to produce the 45 % gradient.

Note: Gradients should be prepared and repacked under sterile conditions (e.g. LAF-bench, ISO Class 5). For optimal results prepare the gradient media maximum 24 hours prior to use. Mix well after diluting the GM501 Gradient 100 %.

Instruction for use with fresh semen samples:
1. Before use warm all components of the system and the samples to 37°C or to room temperature.
2. Mix the density gradient bottles by 5 bottles inversions before use.
3. Pipet 2.5 ml of the lower density gradient (45 %) into a sterile disposable centrifuge tube.
4. Using a 3 cc syringe with a 21 G needle, layer 2.5 ml of the higher density gradient (90 %) under the lower density gradient (45 %) free of air bubbles. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle at the bottom of the centrifuge tube and slowly dispensing the higher density gradient. This two layer gradient is stable for about two hours.
5. Gently place 2.5 ml of liquefied semen onto the upper layer using a transfer pipette or syringe.
6. Centrifuge at 350-400 g for 15-18 minutes. In case, no pellet is visible after this step, centrifuge for another 3 minutes.
7. Aspirate the supernatant.
8. Using a syringe, resuspend the pellet with 2-3 ml of fresh washing medium.
9. Centrifuge at 300 g for 8-10 minutes. In case you want to gain higher sperm concentrations it is advisable to centrifuge for the whole 10 minutes.
10. Aspirate the supernatant and repeat the last two steps.
11. Finally remove the supernatant and resususpend the peltet again in appropriat medium for the use with the subsequent procedure of assisted reproductive medicine (e.g. IVF, ICSI, IUI).

Instruction for use with frozen semen samples:
1. Before use warm all components of the system and the samples to 37°C or to room temperature.
2. Mix the density gradient bottles by 5 bottle inversions before use.
3. Pipet 1 ml of the lower density gradient (45 %) into a sterile disposable centrifuge tube.
4. Using a 3 cc syringe with a 21 G needle, layer 1 ml of the higher density gradient (90 %) under the lower density gradient (45 %) free of air bubbles. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle at the bottom of the centrifuge tube and slowly dispensing the higher density gradient. These two layers of density are stable for about two hours.
5. Gently place maximum 0.5 ml of thawed semen onto the upper layer using a transfer pipette or syringe.
6. Centrifuge for 15-20 minutes at 350 g. In case, no pellet is visible after this step, centrifuge for another 3 minutes.
7. Aspirate the supernatant down to no less than the 0.5 ml mark above the pellet.
8. Using a syringe, resuspend the pellet with 2-3 ml of fresh washing medium.
9. Centrifuge at 300 g for 8-10 minutes. In case you want to gain higher sperm concentrations it is advisable to centrifuge for the whole 10 minutes.
10. Aspirate the supernatant and repeat the last two steps.
11. Finally remove the remaining liquid to have the pellet resuspended again in the desired amount for the use with the subsequent procedure of assisted reproductive medicine (e.g. IVF, ICSI, IUI).

In order to achieve better seperation of the spermatozoa the sample should not be diluted, but rather the centrifugal force should be increased (not higher than 500 g).

Precautions and warnings:
- All human organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (e.g. LAF-bench, ISO Class 5) to avoid possible contamination.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:
- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.

Storage instructions and stability:
- The shelf life is 18 months from time of manufacture.
- Store products between 2-8°C.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Do not freeze before use.
- Do not use after expiry date.
- Keep away from (sun)light.
- Open and close bottles under aseptic conditions (e.g. LAF bench, ISO Class 5)
- Content cannot be re-sterilized after opening.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

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GYNEMED

GM501 Gradient 45% • 90%

Material included in the package:
GM501 Gradient 45 %
Product code: 4 GM 501G-45-10
- 1 x 10 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-50
- 1 x 50 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-100
- 1 x 100 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-250
- 1 x 250 ml GM501 Gradient 45 %

Material included in the package:
GM501 Gradient 90 %
Product code: 4 GM 501G-90-10
- 1 x 10 ml GM501 Gradient 90 %

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Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Gradient 45 % and 90 % is available upon request and can also be downloaded from our website.
- GM501 Gradient 45 % and 90 % is manufactured according to the following specifications:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Osmolality (mOsm/kg)	310-340 (GM501 Gradient 45 %) 320-350 (GM501 Gradient 90 %)
Density (g/ml)	1.1050-1.1150 (GM501 Gradient 90 %)
Viscosity (cP)	< 1.65 (GM501 Gradient 90 %)
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.5
Sperm Survival test survival after 4 hours exposure of density selected spermatozoa to the test medium in %	≥ 80
survival after 24 hours exposure of density selected spermatozoa to the test medium in %	≥ 75
Mouse Embryo Assay (MEA)	not MEA tested

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EN: Intended use/Intended users:
- GM501 Gradient 45 % and 90 % is a gradient system for semen preparation. GM501 Gradient 45 % and 90 % can be used in combination with IUI, IVF and ICSI.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- GM501 Gradient 45 % and 90 % consists of silane-coated colloidal silica particles suspended in HEPES-buffered medium.

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Instruction for use with frozen semen samples:
1. Before use warm all components of the system and the samples to 37°C or to room temperature.
2. Mix the density gradient bottles by 5 bottle inversions before use.
3. Pipet 1 ml of the lower density gradient (45 %) into a sterile disposable centrifuge tube.
4. Using a 3 cc syringe with a 21 G needle, layer 1 ml of the higher density gradient (90 %) under the lower density gradient (45 %) free of air bubles. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle at the bottom of the centrifuge tube and slowly dispensing the higher density gradient. These two layers of density are stable for about two hours.
5. Gently place maximum 0.5 ml of thawed semen onto the upper layer using a transfer pipette or syringe.
6. Centrifuge for 15-20 minutes at 350 g. In case, no pellet is visible after this step, centrifuge for another 3 minutes.
7. Aspirate the supernatant down to no less than the 0.5 ml mark above the pellet.
8. Using a syringe, resuspend the pellet with 2-3 ml of fresh washing medium.
9. Centrifuge at 300 g for 8-10 minutes. In case you want to gain higher sperm concentrations it is advisable to centrifuge for the whole 10 minutes.
10. Aspirate the supernatant and repeat the last two steps.
11. Finally remove the remaining liquid to have the pellet resuspended again in the desired amount for the use with the subsequent procedure of assisted reproductive medicine (e.g. IVF, ICSI, IUI).

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GYNEMED

GM501 Collagenase

Material included in the package:
Product code: 4 GM 501COLL
- 1 x 3 ml GM501 Collagenase

Material not included in the package:
- Incubator (no CO_2)
- LAF-bench (ISO 5 environment)
- Pipettes
- Syringe
- Cannula
- Microscope
- Test tubes
- HEPES-buffered sperm processing medium (e. g. GM501 SpermAir)

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Collagenase is available upon request and can also be downloaded from our website.
- GM501 Collagenase is manufactured and tested according to the following specifications:

pH (at 37°C)	7.20-7.70
Osmolality (mOsm/kg)	270-310
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Sperm Survival test survival after exposure of density selected spermatozoa to the test medium in %	≥ 80

NOTE:
Performing LAL-endotoxin- and MEA-tests is not possible with these media as the activity of the Collagenase inactivates the enzymes (LAL) and damages the mouse embryos (MEA), which are used for these assays, respectively.

EN: Intended use/Intended users:
- GM501 Collagenase is designed for digestion of human testicular tissue obtained by testicular biopsy (TESE). Using Collagenase, it is possible to degrade testicular tissue in single cells to facilitate the isolation of sperm cells, that can be used for the microinjection into oocytes (ICSI).
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alanyl-Glutamine, water, non-essential and essential Amino Acids, HEPES, Phenolred.
- It has an activity of 1000 CDU/ml (Collagen Digestion Units). The Collagenase is obtained from culture filtrates of *Clostridium histolyticum*.

Instructions for use:
1. Extract 1.5 ml GM501 Collagenase out of original bottle and fill into a 5 ml round-bottom tube.
2. Warm the GM501 Collagenase at 37°C. GM501 Collagenase is HEPES-buffered. Incubation in a CO_2-Incubator will lower the pH.
3. For digestion of testicular tissue carefully pick up the chosen tissue pieces with a fine syringe cannula. For easier handling, if necessary, fill the tissue suspension into a 60 mm petri dish. Let adhesive transport or cryo medium drop off as much as possible and transfer into Collagenase tubes.
4. Close the tube completely. Now put the tube for 60 minutes in the incubator or ideally in a heat cabinet for digestion of the tissue. Slight agitation every 20-30 min will facilitate the formation of a single cell suspension.
5. Suspend the digested tissue by carefully pipetting up and down. Under ideal conditions a suspension of single testicular tissue cells and free semen cells has been formed. If coarse tissue pieces are still visible, repeat step 4 for a further 20 to 30 minutes.
6. Now centrifuge the tissue cell suspension and wash twice with 1–2 ml HEPES-buffered sperm processing medium (e. g. GM501 SpermAir). Discard the obtained supernatant. Alternatively, the cell suspension can be processed using a density gradient system (e. g. GM501 Gradient).
7. After the last centrifugation re-suspend the pellet in a small volume of 30–80 µl HEPES-buffered sperm processing medium. Add a few µl of this suspension into a dish.
8. Control by microscope whether viable semen cells are present.
9. Continue IVF treatment according to internal standard procedures.

Precautions and warnings:
- GM501 Collagenase does not contain any antibiotics.
- The long term safety of ICSI including the use of Collagenase, on children following this procedure is still under investigation.
- Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:
- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.
- The GM501 Collagenase may appear cloudy, if warmed up to 37°C. This may occur occasionally without loss of efficiency and is not evidence of microbial contamination. It will disappear, if GM501 Collagenase is cooled down to room temperature.

Storage instructions and stability:
- The shelf life is 9 month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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GYNEMED GmbH & Co. KG