GYNEMED

GM501 SpermStore

Material included in the package:
Product code: 4 SCP-20
- 1 x 20 ml GM501 SpermStore

Material not included in the package:
- Sperm freezing straws
- Sterile pipettes
- Freezing tank with liquid nitrogen
- LAF bench (ISO 5 environment)
- Sperm wash medium

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 SpermStore is available upon request and can also be downloaded from our website.
- GM501 SpermStore is manufactured and tested according to the following specifications:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Sperm Survival test survival after 4 hours exposure of density selected spermatozoa to the test medium in %	≥80
Mouse Embryo Assay (MEA)	not MEA tested

EN: Intended use/Intended users:
- GM501 SpermStore is a medium for freezing human sperm including epididymal or testicular sperm.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- GM501 SpermStore is a ready-to-use HEPES buffered cryopreservation medium which also contains physiologic salts, glycine, dextrose monohydrate, lactate, glycerol, sucrose and human serum albumin (4.00 g/Liter) to protect sperm from damage due to the freezing procedure.

GM501 SpermStore and Spermpreparation:
- GM501 SpermStore can be used in combination with GM501 Gradient.

Instructions for use:
- Ensure all media are well mixed before use.

Before freezing:
- In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted using a suitable sperm washing medium.
- GM501 SpermStore needs to be warmed to room temperature before use to avoid cold-shock.

Freezing:
1. Allow the semen to liquefy at room temperature for 30 minutes.
2. Mix 1 ml of sperm with 0.7 ml of GM501 SpermStore. Add the medium in drops while gently swirling.
3. Leave the mixture for 10 minutes at room temperature for equilibration.
4. Aspirate the sample/medium mixture into the freezing straws, leaving approximately 1.5 cm of air at the end of the straws.
5. Seal the straws.
6. Dry off individually with a lint-free wipe.
7. Shake to move the air-bubble to the centre of the straw.
8. Place the straws horizontally on a styrofoam board in a liquid nitrogen bath to allow for freezing in vapour phase. Leave for (at least) 15 minutes.
9. Transfer straws quickly into liquid nitrogen and store at -196°C.

Thawing:
1. Remove as many straws as required from the liquid nitrogen.
2. Place the samples in tap water for 5 minutes (room temperature or 37°C).
3. Cut off the end of the straw, place the open end inside a container (e.g. a test tube) and tap the straw lightly against the side of the container to allow complete evacuation of the mixture.
4. Dilute the thawed sperm in a suitable insemination medium (at least 3 ml per 0.5 ml semen) and mix thoroughly.
5. Centrifuge for 15 minutes at 300-350 x g.
6. Resuspend pellet in a suitable insemination medium.

After thawing:
- If necessary, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium.

Precautions and warnings:
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Product does not contain any antibiotics. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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Storage instructions and stability:
- The shelf life is 18 months from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Do not use after expiry date.

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The media should be used in the order displayed above.

Material not included in the package:
- 4-well dishes
- Freezing tank with liquid nitrogen
- Water bath (able to hold 37°C)
- Pipettes
- Forceps
- Vitrification device
- LAF bench (ISO 5 environment)
- Microscope
- Lab timer

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for the media are available upon request and can also be downloaded from our website.
- GM501 VitriStore Freeze/VitriStore Thaw are manufactured and tested according to the following specifications:

pH	7.20-7.40
Osmolality (mOsm/kg)	270-290 (VitriStore Freeze PV) 805-850 (VitriStore Thaw 1) 535-565 (VitriStore Thaw 2) 405-435 (VitriStore Thaw 3) 270-290 (VitriStore Thaw 4)
Sterility	sterile - SAL 10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %)	≥80

EN: Intended use/Intended users:
- GM501 VitriStore Freeze/VitriStore Thaw are a set of ready-to-use media for vitrification and thawing of human embryos.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- GM501 VitriStore Freeze/VitriStore Thaw are DMSO/ethylene glycol based vitrification media that also contain PBS, sucrose, Ficoll and human serum albumin (10-20 g/Liter).
- GM501 VitriStore Freeze/VitriStore Thaw do not contain antibiotics.

Instructions for use:
- Ensure all media are well mixed before use. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps
- In a 4-well dish fill the first well with 300 µl of VitriStore Pre-vitrification Medium, the second with 250 µl VitriStore Freeze Medium 1 and the third with 250 µl VitriStore Freeze Medium 2. Next open as many packs of vitrification devices as will be required for the vitrification step. Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Up to 5 vitrification cycles (of the same patient!) can be performed with one media set-up. Do not use the same media for different patients!

Freezing preparation
Transfer the embryos from the blastocyst cell culture medium into each of the VitriStore Freeze solutions using the following scheme:

Stage	VitriStore Pre-vitrification	VitriStore Freeze 1	VitriStore Freeze 2	Temperature
Early Blastocyst, Morulae	2 min.	2 min.	30 sec.	Room temperature
Blastocyst-expanded Bl.	2 min.	3 min.	30 sec.	37°C
Blastocyst-expanded Bl + artificial shrinkage*	2 min.	2 min.	30 sec.	Room temperature

**Before starting the vitrification procedure, in order to reduce the negative effect of the blastocoel, expanded blastocysts should be collapsed by reducing artificially with a glass pipette the volume of the blastocoel. (Vanderzwalmen et. al., 2002; Son et al., 2003; Hiraoka, 2004)*
Vitrification
1. Using an attenuated pipette or an equally suitable device, deposit maximum 2 blastocysts in a volume of approximately 0.3 µl of VitriStore Freeze 2 in the gutter of the tip of the vitrification straw.
2. Place the vitrification straw in the outer sheath and seal it as indicated in the instructions for use of the vitrification device.
3. Plunge the sealed device into the liquid nitrogen.

Thawing
1. In a 4-well dish, fill the first well with 1 ml VitriStore Thaw Medium 1, the second with VitriStore Thaw Medium 2 (250-300 µl), The third with VitriStore Thaw Medium 3 (250-300 µl), the fourth with VitriStore Thaw Medium 4 (250-300 µl). Preheat the media to 37°C.
2. Remove the vitrification straw from the outer sheath as indicated in the instructions for use of the vitrification device.
3. Immediately plunge the vitrification straw into pre-heated VitriStore Thaw Medium 1 (37°C) and leave in Thawing 1 Medium for 3 minutes.
4. Transfer into VitriStore Thaw Medium 2 (37°C) and leave in this medium for 2 minutes.
5. Transfer into VitriStore Thaw Medium 3 (37°C) and leave in this medium for 2 minutes.
6. Finally transfer into VitriStore Thaw Medium 4 (37°C) and wash for at least 1 minute.
7. Transfer into blastocyst culture medium for continued cell culture.

GM501 VitriStore Freeze/VitriStore Thaw and Embryoculture:
- GM501 VitriStore Freeze/VitriStore Thaw can be used in combination with GM501 Media (Washing- and IVF-Media) before freezing and after thawing.

Precautions and warnings:
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (LAF-bench ISO Class 5) to avoid possible contamination.
- Only for the intended use.
- The long term safety of embryo vitrification on children born following this procedure is unknown.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:
- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.

Storage instructions and stability:
- The shelf life is 12 months from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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GYNEMED

GM501 GentleVit

GentleVit Freeze Kit
Material included in the package:

Product code: 4 GVF_KIT1

- 1 x 5 ml GentleVit Pre-vitrification Medium (GPI) (4 GPI1005)
- 1 x 1 ml GentleVit Freeze Medium 1 (GF1) (4 GF1001)
- 1 x 1 ml GentleVit Freeze Medium 2 (GF2) (4 GF2001)
- 1 x 1 ml GentleVit Freeze Medium 3 (GF3) (4 GF3001)
- 1 x 1 ml GentleVit Freeze Medium 4 (GF4) (4 GF4001)
- 1 x 1 ml GentleVit Freeze Medium 5 (GF5) (4 GF5001)

GentleVit Thaw Kit
Material included in the package:

Product code: 4 GVT_KIT1

- 1 x 5 ml GentleVit Thaw Medium 1 (GT1) (4 GT1005)
- 1 x 1 ml GentleVit Thaw Medium 2 (GT2) (4 GT2001)
- 1 x 1 ml GentleVit Thaw Medium 3 (GT3) (4 GT3001)
- 1 x 1 ml GentleVit Thaw Medium 4 (GT4) (4 GT4001)
- 1 x 1 ml GentleVit Thaw Medium 5 (GT5) (4 GT5001)
- 1 x 1 ml GentleVit Thaw Medium 6 (GT6) (4 GT6001)

The media should be used in the order displayed above (the bottles may be in a different order in the box). One kit will provide sufficient medium for approximately 3-4 procedures.

Up to 5 vitrification cycles (of the same patient) can be performed with one media set-up. Do not use the same media for different patients!

Material not included in the package:

- Well dishes
- Freezing tank with liquid nitrogen
- Water bath (able to hold 37°C)
- Pipettes
- Forceps
- Vitrification device
- LAF bench (ISO 5 environment)
- Microscope
- Labtimer

Product specifications and quality control:

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for the media are available upon request and can also be downloaded from our website.
- GM501 GentleVit Freeze/Thaw are manufactured and tested according to the following specifications (according to Ph. Eur./ USP):

pH	7.20-7.50 (release criteria: 7.20-7.40)
Osmolality (mOsm/kg)	270-295 (GPI/GT6: release criteria: 270-290) 805-865 (GT3: release criteria: 805-850) 535-565 (GT4) 405-435 (GT5)
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
MEA (1-cell assay, blastocysts after 96 h in %)	≥80

EN: Intended use/Intended users:

- GM501 GentleVit Freeze/GentleVit Thaw are a set of ready-to-use media for vitrification and warming of human oocytes and embryos.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:

- GM501 GentleVit Freeze/GentleVit Thaw are HTF-based and contain HEPES, sucrose, human serum albumin (12-20 g/liter). Cooling media also contain DMSO, Ethylene Glycol (EG) and Ficoll.
- GM501 GentleVit Freeze/GentleVitThaw does not contain antibiotics.

Vitrification of Oocytes:

Preparation for Vitrification

- Ensure that all media bottles of the kit are well mixed before use and warmed to room temperature (20-25°C). Optionally, GT1 can be warmed to 37°C. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps:

- In a 6-well dish, fill the first well with 250-300 µl of GentleVit Pre-vitrification medium, the second with Freeze Medium **1**, the third with Freeze Medium **2** and continue doing so until Freeze Medium **5**.
- Open the necessary number of vitrification devices, taking into account that 1 device can hold 2-3 oocytes, in a maximum volume load of 1 µl (check the instructions for the device you are using). Conveniently place the separate parts of the device on the workbench for easy access later in the procedure.

Vitrification protocol:

- Oocytes are sequentially exposed to the following media:

	GPI	GF1	GF2	GF3	GF4	GF5
	2 min.	3 min.	3 min.	3 min.	5-6 min.	60 sec.*
DMSO/EG (%)	0	1.25	2.5	5	10	20

***Note:** The complete process of placing the oocyte in "Freeze Medium 5", loading the oocyte on the vitrification device in maximum 1µl GF5, inserting the device in the outer straw and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen.*

Thawing of Oocytes: (Zygote to Blastocyste):

Preparation for Thawing

- Ensure that all media bottles of the kit are well mixed before use and warmed to room temperature (20-25°C). One exception is Warming 1, which may be warmed to 37°C before use. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps:

- For thawing of 4-cell embryos till blastocyst, the following media are **NOT** required:

GM501 GentleVit Thaw Medium 1
GentleVit Freeze Medium 2

- In a 6-well dish, fill the first well with 500-800µl of Thaw Medium **1**, the second with Thaw Medium **2**, the third with Thaw Medium **3**, the fourth with Thaw Medium **4**, the fifth with Thaw Medium **5** (only for zygotes) and the last with Thaw Medium **6**.

Thawing protocol:

- In the first warming step, ensure homogeneous temperature by swirling the straw gently in GT1.
- Oocytes are sequentially exposed to the following media:

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***Note:** Wash for 1-2 min before transfer to culture medium.*

Vitrification of Embryos (Zygote to Blastocyst):

Preparation for Vitrification

- Ensure that all media bottles of the kit are well mixed before use and warmed to room temperature (20-25°C). We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps:

- For vitrification of embryos, the following media are **NOT** required:

GentleVit Freeze Medium 1
GentleVit Freeze Medium 2

- In a 6-well dish, fill the first well with 250-300µl of GentleVit Pre-vitrification medium, the second with Freeze Medium **3**, the third with Freeze Medium **4** and the last with Freeze Medium **5**.
- Open the number required of vitrification devices, taking into account that 1 device can hold 1-2 embryos in a maximum volume load of 1µl (check the instructions for the device you are using). Conveniently place the separate parts of the device on the workbench for easy access later in the procedure.

Vitrification protocol:

- Embryos are sequentially exposed to the following media:

	GPI	GF3	GF4	GF5
Zygotes	2 min.	5 min.	5 min. 30 sec.	40-60 sec.*
4-cell to blastocyst	2 min.	5 min.	4. min.	40-60 sec.*
DMSO/EG (%)	0	5	10	20

***Note:** The complete process of placing the oocyte in "Freeze Medium 5", loading the oocyte on the vitrification device in maximum 1µl GF5, inserting the device in the outer straw and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen.*

Thawing of Embryos (Zygote to Blastocyst):

Preparation for Thawing

- Ensure that all media bottles of the kit are well mixed before use and warmed to room temperature (20-25°C). One exception is Warming 1, which may be warmed to 37°C before use. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps:

- For thawing of 4-cell embryos till blastocyst, the following media are **NOT** required:

GM501 GentleVit Thaw Medium 1
GentleVit Freeze Medium 2

- In a 6-well dish, fill the first well with 500-800µl of Thaw Medium **1**, the second with Thaw Medium **2**, the third with Thaw Medium **3**, the fourth with Thaw Medium **4**, the fifth with Thaw Medium **5** (only for zygotes) and the last with Thaw Medium **6**.

Thawing protocol:

- In the first warming step, ensure homogeneous temperature by swirling the straw gently in GT1.
- Embryos are sequentially exposed to the following media:

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***Note:** Wash for 1-2 min before transfer to culture medium.*

Precautions and warnings:

- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (LAF-bench ISO Class 5) to avoid possible contamination.
- Only for the intended use.
- The long term safety of oocyte and embryo vitrification on children born following this procedure is unknown.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.

Storage instructions and stability:

- The shelf life is 12 month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C). Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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GentleVit Freeze Medium 1
GentleVit Freeze Medium 2

- In a 6-well dish, fill the first well with 250-300µl of GentleVit Pre-vitrification medium, the second with Freeze Medium **3**, the third with Freeze Medium **4** and the last with Freeze Medium **5**.
- Open the necessary number of vitrification devices, taking into account that 1 device can hold 1-2 embryos in a maximum volume load of 1µl (check the instructions for the device you are using). Conveniently place the separate parts of the device on the workbench for easy access later in the procedure.

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GYNEMED

GM501 Cult

Material included in the package:
GM501 Cult
Product code: 4 GM 501H-20
- 1 x 20 ml GM501 Cult

Product code: 4 GM 501H-50
- 1 x 50 ml GM501 Cult

Material included in the package:
GM501 Cult with Gentamicin
Product code: 4 GM 501H+G-20
- 1 x 20 ml GM501 Cult with Gentamicin

Product code: 4 GM 501H+G-50
- 1 x 50 ml GM501 Cult with Gentamicin

Material included in the package:
GM501 Cult with Gentamicin and Phenolred
Product code: 4 GM 501H+PR+G-20
- 1 x 20 ml GM501 Cult with Gentamicin and Phenolred

Product code: 4 GM 501H+PR+G-50
- 1 x 50 ml GM501 Cult with Gentamicin and Phenolred

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pH (at 37°C, 6% CO_2)	7.20-7.45
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %, 96 h exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- The media are ready-to-use bicarbonate-buffered culture media, designed for fertilisation and for human embryo culture from day 1 to blastocyst stage. It can also be used for embryo transfer.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
GM501 Cult
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50% solution), Na-Pyruvate, EDTA, Alanyl-Glutamine, water, non-essential and essential Amino Acids, human serum albumin (10.00 g/liter).

GM501 Cult with Gentamicin
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50% solution), Na-Pyruvate, EDTA, Alanyl-Glutamin, Wasser, nicht-essentielle und essentielle Aminosäuren, Humanalbuminlösung (10,00 g/liter).

GM501 Cult with Gentamicin and Phenolred
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50% solution), Na-Pyruvate, EDTA, Alanyl-Glutamin, water, non-essential and essential Amino Acids, human serum albumin (10.00 g/liter), Gentamicin (10.00 mg/liter), Phenolred.

Instructions for use:
- The medium must be equilibrated over night in a humidified CO_2-incubator (at 5-7 % CO_2, 37°C).
- The media are designed for embryo culture under oil as well as in open systems.

Precautions and warnings:
- GM501 Cult with Gentamicin and GM501 Cult with Gentamicin and Phenolred contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Parmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination, even when GM501 Cult medium contains Gentamicin.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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	Consult instructions of use
8 °C / 2 °C	Temperature limitation
LOT	Lot number
	Use by
REF	Catalogue number
	Do not use if package is damaged
STERILE A	Sterilized using aseptic processing techniques (filtration)
	Keep away from sunlight
	Do not resterilize

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GYNEMED

GM501 Flush

Material included in the package:
Product code: 4 GM 501F-20
- 1 x 20 ml GM501 Flush

Product code: 4 GM 50F-50
- 1 x 50 ml GM501 Flush

Product code: 4 GM 501F-100
- 1 x 100 ml GM501 Flush

Product code: 4 GM 50F-250
- 1 x 250 ml GM501 Flush

Product code: 4 GM 50F-500
- 1 x 500 ml GM501 Flush

Material not included in the package:
- Incubator (no CO_2)
- Petri-dishes
- LAF bench (ISO 5 environment)
- Microscope
- Test tubes
- Rubber free syringe

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Flush is available upon request and can also be downloaded from our website.
- GM501 Flush is manufactured and tested according to the following specifications:

pH (at 37°C)	7.20-7.60 (Release criteria: 7.20-7.40)
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %, 60 min exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- Cell culture medium for human oocyte pick-up.
- GM501 Flush is a ready-to-use medium for flushing the ovarian follicles during aspiration and/or oocyte pick-up intended for extra corporeal fertilization procedures.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alanyl-Glutamin, water, non-essential and essential, Amino Acids, HEPES, 2.5 IU/ml Heparin.
- Gentamicin sulfate can be added by the user (10 mg/liter).

Instructions for use:
- GM501 Flush needs to be warmed at 37°C over night before use (no CO_2, with closed lid). GM501 Flush is HEPES-buffered. Incubation in a CO_2-incubator will lower the pH.

Oocyte collection:
- Keep GM501 Flush medium at 37°C or as close to this as possible.
- Take about 20-50 ml medium at a time, using a rubber free syringe. Assemble needle and tubing set according to the manufacturers instructions.
- Attach syringe containing GM501 Flush medium to the appropriate Luer, allowing connection to the lumen of the needle.
- The person holding the syringe should try to apply warmth to it by holding it in the palm of the hand covering as much of the syringes surface as possible.
- During aspiration of follicles, if an ovum is not extracted after a follicle has been emptied, flushing is indicated. The collapsed follicle may be refilled with GM501 Flush, after which it may be aspirated once again. This may be done repeatedly, after refilling the syringe with fresh aspiration medium each time, until an ovum is recovered or the surgeon decides to abandon the follicle (after three attempts at flushing it may be sensible to move on, but the ultimate decision remains with the surgeon).

Precautions and warnings:
- All human organic material should be considered potentially infectious. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination.
- GM501 Flush contains heparin which is derived from porcine intestinal mucosa.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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	Consult instructions of use
	Temperature limitation
LOT	Lot number
	Use by
REF	Catalogue number
	Do not use if package is damaged
STERILE A	Sterilized using aseptic processing techniques (filtration)
	Keep away from sunlight
	Do not resterilize

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GYNEMED

GM501 SpermAir

Material included in the package:
Product code: 4 GM 501AIR+PR+G-02
- 1 x 2 ml GM501 SpermAir

Product code: 4 GM 501AIR+PR+G-20
- 1 x 20 ml GM501 SpermActive

Product code: 4 GM 501AIR+PR+G-50
- 1 x 50 ml GM501 SpermActive

Material not included in the package:
- Incubator (no CO_2)
- Petri-dishes
- LAF bench (ISO 5 environment)
- Microscope
- Test tubes

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 SpermAir is available upon request and can also be downloaded from our website.
- GM501 SpermAir is manufactured and tested according to the following specifications:

pH (at 37°C)	7.20-7.50
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %, 60 min exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- GM501 SpermAir is a ready-to-use medium designed for all human sperm-preparation (e.g. density gradient centrifugation), sperm-washing and swim-up techniques as well as for testicular tissue.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50% solution), Na-Pyruvate, EDTA, Alanyl-Glutamine, 21 mM HEPES, water, non-essential and essential Amino Acids, human serum albumin (5.00 g/liter), Gentamicin (10.00 mg/liter), Phenolred.

Instructions for use:
- Do not equilibrate GM501 SpermAir in a CO_2-incubator, just warm it up to 37°C. GM501 SpermAir is HEPES-buffered. Incubation in a CO_2-incubator will lower the pH.

Precautions and warnings:
- GM501 SpermAir contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination, even when GM501 SpermAir Medium contains Gentamicin.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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GYNEMED

GM501 PVP

Material included in the package:
Product code: 4 PVP 0001
- 1 x 1 ml GM501PVP

Product code: 4 PVP 0001-5
- 5 x 0.2 ml GM501 PVP

Material not included in the package:
- Incubator at 37°C (no CO_2)
- Petri-dishes
- LAF bench (ISO 5 environment)
- Microscope
- Pipettes
- Mineral Oil (e.g. GM501 Mineral Oil)
- HEPES-buffered sperm processing

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- GM501 PVP has been evaluated by an independent test laboratory using the ISO 10993-5/USP30-NF25 Agar Diffusion test for cytotoxicity. The product was considered non-cytotoxic.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 PVP is available upon request and can also be downloaded from our website.
- GM501 PVP is manufactured and tested according to the following specifications:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Osmolality (mOsm/kg)	280-320
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 1.00
Viscosity at 25°C (cP)	> 170
MEA (1-cell assay, blastocysts after 96 h in %)	≥80

EN: Intended use/Intended users:
- GM501 PVP is a ready-to-use medium to reduce the motility of sperm, making it easier to catch them with an ICSI pipette.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- GM501 PVP is 10 % dialysed PhEur grade polyvinylpyrrolidone solved in an aqueous solution containing physiologic salts, HEPES, lactate, pyruvate, glucose and human albumin solution (4.00 g/liter).

Used Abbreviations:
- ICSI - Intra cytoplasmatic Sperm Injection
- PVP - Polyvinyl pyrrolidone

Instructions for use:
- Evaporation of the PVP solution on the slide or dish and the presence of human albumin may occasionally cause a small amount of precipitation in the medium.

Standard procedure:
1. Warm the PVP solution to 37°C.
2. Place a small drop of PVP solution (5 µl - 10 µl) in a dish and cover with GM501 MineralOil.
3. Add a small volume (1 µl - 2 µl) of washed sperm into the centre of the PVP droplet.
4. Wait for a few minutes to allow sperm to migrate to the periphery of the droplet.
5. Select and recover the spermatozoa for injection.

Alternative procedure (with extra washing step):
1. Warm the PVP solution and HEPES-buffered sperm processing medium to 37°C.
2. Place a small drop of PVP solution (5 µl - 10 µl) and 1 or more small drops HEPES-buffered sperm processing medium in a dish and cover with GM501 MineralOil.
3. Add a small volume (1 µl - 2 µl) of washed sperm into the centre of the PVP droplet.
4. Wait for a few minutes to allow sperm to migrate to the periphery of the droplet.
5. Select the spermatozoa for injection and nick (break) the tail of the spermatozoon with the tip of the pipette.
6. Transfer the spermatozoon into one of the HEPES-buffered sperm processing medium droplets and wash by transferring the sperm cell in and out of the sperm processing medium several times.
7. Aspirate the sperm cell into the pipette and use for ICSI procedure.

Precautions and warnings:
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- GM501 PVP medium does not contain any antibiotics. Always work under strict hygienic conditions (e.g. LAF-bench, ISO Class 5) to avoid possible contamination.
- The long term safety of ICSI, including the use of PVP, on children born following this procedure is still under investigation.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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Storage instructions and stability:
- The shelf life is 9month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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	Consult instructions of use
8 °C / 2 °C	Temperature limitation
LOT	Lot number
	Use by
REF	Catalogue number
	Do not use if package is damaged
STERILE A	Sterilized using aseptic processing techniques (filtration)
	Keep away from sunlight
2 STERILIZE	Do not resterilize

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GYNEMED

GM501 Wash

Material included in the package:
GM501 Wash
Product code: 4 GM 501W-20
- 1 x 20 ml GM501 Wash

Product code: 4 GM 501W-50
- 1 x 50 ml GM501 Wash

Product code: 4 GM 501W-500
- 1 x 500 ml GM501 Wash

Material included in the package:
GM501 Wash with Phenolred and Gentamicin
Product code: 4 GM 501W+G+PR-20
- 1 x 20 ml GM501 Wash with Phenolred and Gentamicin

Product code: 4 GM 501W+G+PR-50
- 1 x 50 ml GM501 Wash with Phenolred and Gentamicin

Material not included in the package:
- CO_2-incubator at 37°C
- Petri-dishes
- LAF bench (ISO 5 environment)
- Microscope
- Test tubes

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for the media are available upon request and can also be downloaded from our website.
- The media are manufactured and tested according to the following specifications:

pH (at 37°C, 6% CO_2)	7.20-7.50
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %, 60min exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- The media are ready-to-use media designed for washing procedures of human oocytes and embryos and any short-term handling procedures outside the incubator like washing after Hyaluronidase treatment (denudation), ICSI, polar body or blastomere biopsy and others.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
GM501 Wash
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alanyl-Glutamine, water, non-essential and essential Amino Acids, 15 mM HEPES, human serum albumin (5.00 g/liter).

GM501 Wash with Phenolred and Gentamicin
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrat, Na-Lactat (50 % Lösung), Na-Pyruvat, EDTA, Alanyl-Glutamin, Wasser, nicht-essentielle und essentielle Aminosäuren, 15 mM HEPES, Humanalbuminlösung (5.00 g/Liter).

GM501 Wash with Phenolred and Gentamicin
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alanyl-Glutamine, water, non-essential and essential Amino Acids, 15 mM HEPES, human serum albumin (5.00 g/liter), Gentamicin (10 mg/liter), Phenolred.

Instructions for use:
- The medium must be equilibrated over night in a humidified CO_2-incubator (at 5-7 % CO_2, 37°C).

Precautions and warnings:
- GM501 Wash with Phenolred and Gentamicin contains the antibiotic: Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 enviroment, e.g. LAF-bench) to avoid possible contamination, even when GM501 Wash medium contains Gentamicin.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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Storage instructions and stability:
- The shelf life is 6 month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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GYNEMED

GM501 Hyaluronidase

Material included in the package:
Product code: 4 HY 0001-5
- 5 x 1 ml GM501 Hyaluronidase

Product code: 4 HY 0010
- 1 x 10 ml GM501 Hyaluronidase

Material not included in the package:
- Incubator at 37°C (no CO_2)
- Petri dishes
- LAF bench (ISO 5 environment)
- Pipette suitable for denudation of oocytes
- Mineral Oil (e.g. GM501 Mineral Oil)
- Washing medium (e.g. GM501 Wash)
- Microscope

Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Hyaluronidase is available upon request and can also be downloaded from our website.
- GM501 Hyaluronidase is manufactured and tested according to the following specifications:

pH	7.30-7.60
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 1.00
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96h in %, 30 min exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- GM501 Hyaluronidase is a ready-to-use solution designed to facilitate the mechanical removal of cumulus cells. Hyaluronidase digests the extracellular matrix in the cumulus-oocyte complex consisting of hyaluronic acid.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- 80 IU/ml pharmaceutical grade hyaluronidase from bovine origin solved in HEPES-buffered medium, containing 4.00 g/liter human serum albumin.

Instructions for use:
- GM501 Hyaluronidase contains HEPES; no CO_2 incubation is required, just warm it up to 37°C (incubation in a CO_2 incubator will lower the pH below 7).
- Depending on the dishes, 100 µl (micodrops) or 400-600 µl (center-/4-well) of GM501 Hyaluronidase should be used for denudation.
- Additionally prepare 3-5 drops (100 µl) of washing medium (e.g. GM501 Wash) for oocyte washing.
- The dishes should be covered with Mineral Oil (e.g. GM501 Mineral Oil).
- Place the oocyte in the Hyaluronidase (up to max. 5 oocytes) for about 30 seconds.
- Using a fine pipette, transfer the partially denuded oocytes to the first washing drop.
- Remove the corona by repeatedly pipetting the oocytes.
- Use the other drops to further wash the denuded oocyte.

Precautions and warnings:
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Product does not contain any antibiotics. Always work under strict hygienic conditions (e.g. LAF-bench, ISO Class 5) to avoid possible contamination.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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Bovine sourced Hyaluronidase:
- The Pharmaceutical grade hyaluronidase used in this product is sourced from bovine testis. It is certified not to contain any neuronal tissue (Certificate of Suitability/ CEP).
- Hyaluronidase is sourced from animals determined „fit for human consumption“ and originating from countries with „negligible BSE risk“ (Resolution No. 18 „Recognition of the Bovine Spongiform Encephalopathy Risk Status of Member Countries (2014)“, adopted by the OIE).
- According to the WHO guidelines on Tissue Infectivity Distribution in Transmissible Spongiform Encephalopathies (2010), testes from bovine source are classified as Category IC material (i.e. “Tissues with no detected infectivity”).

Storage instructions and stability:
- The shelf life is 12 months from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- Do not use after expiry date.

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	Consult instructions of use
8°C / 2°C	Temperature limitation
LOT	Lot number
	Use by
REF	Catalogue number
	Do not use if package is damaged
STERILE A	Sterilized using aseptic processing techniques (filtration)
	Keep away from sunlight
	Do not resterilize

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GYNEMED

GM501 SpermActive

Material included in the package:
GM501 SpermActive
Product code: 4 GM 501SA-020
- 1 x 20 ml GM501 SpermActive

Product code: 4 GM 501SA-050
- 1 x 50 ml GM501 SpermActive

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Product specifications and quality control:
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 SpermActive is available upon request and can also be downloaded from our website.
- GM501 SpermActive is manufactured and tested according to the following specifications:

pH (at 37°C, 6 % CO_2)	7.20-7.50
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL10^{-3} (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25
Mouse Embryo Assay (MEA) (1-cell assay, blastocysts after 96 h in %, 60 min exposure at zygote stage)	≥80

EN: Intended use/Intended users:
- GM501 SpermActive is a ready-to- use medium designed for all human sperm preparation, sperm washing and swim up techniques like density-gradient centrifugation as well as for testicular tissue.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:
- NaCl, KCl, KH_2PO_4, $MgSO_4 7H_2O$, $NaHCO_3$, $CaCl_2 2H_2O$, D(+)-Glucose anhydrous, Na-Lactate, Na-Pyruvate, EDTA, Alanyl-Glutamine, water, non-essential and essential Amino Acids, 15 mM HEPES, human serum albumin (5.00 g/liter), Gentamicin (10.00 mg/liter), Phenolred

Instructions for use:
- GM501 SpermActive must be equilibrate over night in a humidified CO_2-incubator (at 5-7 % CO_2, 37°C).

Precautions and warnings:
- GM501 SpermActive contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination, even when GM501 SpermActive Medium contains Gentamicin.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

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GYNEMED GmbH & Co. KG