AUSTRIA
Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

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en english

These instructions for use can also be downloaded from: www.bauerfeind.com/downloads

Dear Customer,

thank you for placing your trust in Bauerfeind. The VenoTrain compression stockings you have purchased represent a high quality medical device. Excellent workmanship, perfect fit and use of yarns that are uniquely soft and gentle to the skin, contribute to ensuring a pleasant yet therapeutic wearing experience. For a successful therapy, we recommend wearing your VenoTrain compression stockings daily. Consult your physician on a regular basis. This will help to ensure long term therapeutic success.

If there is a CE mark on the sewn-in label on your stockings, they are a standard product. If your name is printed on the sewn-in label and there is no CE mark, the stockings are a custom-made product in the sense of § 3 No. 8 of the German Medical Devices Act.

Please read these **user instructions** carefully as they contain important information about use, indications and contraindications, as well as notes regarding fitting and care. Keep the instructions for use so that you can consult them later if any questions arise.

Intended purpose

- Long-term treatment of phlebological diseases in the case of stockings
- Long-term treatment of lymphatic diseases in the case of arm sleeves and stockings

Indications

- Prophylaxis (prevention) of vein problems
- Primary (hereditary vein weakness) and secondary varicosis (as collateral circulation in the condition after deep leg vein thrombosis has occurred)
- Guide vein insufficiency
- Varicosis during pregnancy
- Thrombosis prophylaxis (prevention of a thrombosis)
- Thrombophlebitis (superficial inflammation of the vein) and condition after phlebitis has healed
- Deep leg vein thrombosis and condition after thrombosis (blood clot formation in the veins)
- Post-thrombotic syndrome (chronic venous insufficiency as a result of a deep leg vein thrombosis)
- Conditions resulting from chronic venous insufficiency (e.g. discoloration and hardening of the skin)
- After sclerotherapy (vein obliteration) and venous surgery
- Edemas (tissue swelling due to increased fluid retention) during pregnancy
- Post-traumatic (as a result of injury) and post-operative edemas
- Lymphedemas (tissue swelling due to obstruction of lymph drainage) and lipedemas (swelling of the fatty tissue)
- Angiodysplasias (congenital vascular malformations)

The compression class needs to be selected by the consultant physician. Regular medical checkups are an important part of therapeutic success.

Permissible conditions and places of use

In accordance with the indications

Application risks

⚠ Caution*

There are no known application risks if the product is used properly, provided your legs have been measured by a trained specialists and you have an understanding of the correct way of donning the garment. In the case of stockings with a silicone top band, skin irritation can occur in patients who are intolerant to this material.

Contraindications

- Advanced peripheral arterial occlusion (circulatory disturbance in the leg arteries)
- Decompensated cardiac insufficiency (myocardial weakness)
- Septic phlebitis (bacterial inflammation of the vein)
- Phlegmasia coerulea dolens (acute and extreme venous occlusion with blue discoloration and severe swelling of the leg)
- Weeping dermatosis (skin rash)
- Material intolerances
- Paresthesia of the extremity (numbness and reduced pain perception with advanced peripheral neuropathy)

The manufacturer accepts no liability – in particular product liability – in the event of failure to comply with the above contraindications.

Application instructions

- The initial fitting of the product must only be performed by medically trained specialists.
- It is best to put on your stockings first thing in the morning, as soon as you get up.
- Wear rubber gloves (e.g. Bauerfeind gloves) to put on and take off your VenoTrain compression stockings as this makes it easier to distribute the material over the leg.
- Your legs and hands must be dry. **Please note: compression stockings are fitted and not pulled on!**
- Avoid grasping the fabric with sharp fingernails. Make sure that you take off all jewelry before starting.
- Please avoid over-stretching the stockings, and do not pull up the product by the top band.
- Protect your VenoTrain compression stockings against contact moisturizers as this may have a damaging effect on the fabric.
- We recommend applying moisturizers in the evening after you have removed your compression stockings.
- Avoid damaging the knitted fabric (e.g. as a result of wearing unsuitable shoes) and do not cut off any loose threads or labels, as this may affect the durability of the garment.
- If you experience any unexpected problems with the compression stockings, take them off as a precaution and consult a physician.
- Standard products can be stored for up to 36 months, after which medical compression stockings can be worn for up to 6 months. The product's storage life under proper storage conditions is indicated by an "hourglass" symbol printed on the label on the packet. Custom-made medical compression stockings are intended for immediate use and are therefore not given a storage life date.

Cleaning instructions

- Avoid dry cleaning and ironing.
- Prolong the life of your stockings by daily handwashing or machine washing on a delicate wash cycle (please note the care instructions on the label of your stockings).
- The product can be machine dried provided a delicate cycle is chosen (refer to the information on the label). However, direct heat should be avoided (radiators, sunshine). When handwashing the stockings, we recommend drying them between two towels (do not wring).
- Only use a mild washing powder / liquid and avoid the use of washing additives like fabric softener, optical brighteners or stain removers.
- Please wipe the inside of the silicone top band with a damp cloth regularly (e.g. using Bauerfeind silicone top band cleaning cloths).
- Store your VenoTrain compression stockings in a dry place away from direct light.

Please refer to the care instructions on the fabric label of your compression stocking.

Information on the material used

Bauerfeind is authorized according to the Oeko-Tex Standard 100 to use the "Confidence in Textiles" trade mark, as we only use materials which have been tested for harmful substances. Further information on material content can be found on the label sewn into every compression stocking.

Notes on reuse

VenoTrain compression stockings are intended for your treatment only. For this reason, do not pass the VenoTrain on to other people.

Warranty

Warranty is provided in accordance with the statutory provisions The legal provisions of the country in which the product was purchased apply. Country-specific warranty regulations that govern the relationship between the retailer and the buyer also apply where relevant. If you believe that a claim may be made under the warranty, please first contact the retailer from which you purchased the product directly.

Please do not make any changes to the product yourself. This particularly applies to individual adjustments made by specialists. Please follow our instructions for use and care. These are based on our experience gained over many years and ensure the functionality of our medically effective products for a long time. Only aids that function optimally can help you. Failure to follow the instructions may also restrict the warranty.

Disposal

The compression stockings can be disposed of in household waste. If deviating national statutory provisions exist, these must be applied.

Version: 2017-12

For more details, please consult your physician or medical retailer.

Donning tips

Product with closed toe option

1) Turn the stocking inside-out up to the heel (fold it over to the left). Pull the foot part over your toes and then over the upper part of your foot.
2) Fold the leg part back over and pull it over your heel at the same time. Make sure that the heel is in the correct position.
3) Now pull the stocking up gradually and distribute it uniformly until it is free of creases. Avoid pulling on the top band when donning the garment.

Product with open toe option

1) Use the blue VenoTrain donning aid supplied.
2) Gather the leg part together and pull the foot part over the donning aid. Make sure that the heel is in the correct position.
3) Now pull the stocking up gradually and distribute it uniformly until it is free of creases. Avoid pulling on the top band when donning the garment. You can now remove the donning aid by withdrawing it from the open toe.

Please follow the instructions on the foldout sheet of this brochure!

Putting on your VenoTrain compression stockings is made much easier by the VenoTrain glider donning aid.

*Note on the risk of personal injury (risk of injury, harm to health, accidents and property damage)

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