dimeda®
SURGICAL INSTRUMENTS

To whom it may concern

STATEMENT

Herewith, «DIMEDA Instrumente GmbH», having its principal business office at Gänsäcker 54+58, 78532 Tuttllingen, Germany, represented by its General Manager, Mr. Bastian Milkau certifies the following information:
«Instruction Manual for reusable surgical instruments» covers the following range of our products:
scalpel knives, knife handles, scalpel blades, knives, scissors, punches, dissectors, nasal files, drilling instruments, raspatories, skinny raspatories, extractors with accessories, rib-cutting forceps, bone forceps, surgical snares, surgical spoons, curettes, bone files, bone saws with accessories, drills, bone pliers, bone chisels, bone punches, tractors,bone forceps, screwdrivers with accessories, plaster knives, plaster dilator, hand plaster saws, screw taps, sleeves with accessories, plates and plate sets with accessories, needles, suture needles, needle holders, awl, staples, clamping applicators, ligature instrument with accessories, pulmonary spatulas, spatulas, dilators with accessories, spreaders, mouth dilators with accessories, operating retractors with accessories, speculum with accessories, hooks with accessories, depressors, elevators with accessories, scalers, artery clamps with accessories, bulldog clamps, clamps for umbilical cord, clamps with accessories, surgical tissue clamps, tweezers, biopsy forceps, forceps others, ear mirror with accessories, goniometers, dynamometers, cephalometer, esthesiometer, tonometer with accessories, caliper, percussion hummers with accessories, camertone with accessories, pin wheel, caliper gauge, stethoscopes with accessories, head lamp, head light, probes with accessories, forehead mirrors with accessories, skin markers, anchors with accessories, hummers with accessories, microsurgical instruments with accessories, obstetrical forceps, urethrotomes with accessories, uterine dilators with accessories, vaginal speculum with accessories, metallic catheters, surgical clamps for penis, additional accessories for gynecology and urology, dental forceps, chisel, tweezers, curettes, periodontometer, excavators, amalgam plugger, polish instrument, filling instrument, additional instruments for filling, separator, mould forceps, matrix, mouth spoon, modeling instrument, modeling carver, hummers, wire dispenser, additional dental accessories,

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«Instructions for Use and Sterilization Instructions for Orthopaedic Implants», «Instruction Manual for reusable surgical instruments» cover cutting instrument for needles, wires with accessories.

«Instruction Manual for reusable surgical instruments», «Instruction for Use Suction – Irrigation Instruments» cover tubes for aspiration/ irrigation with accessories.

«Instruction Manual for Sterile Containers with Accessories» cover sterilization container with accessories.

«Instruction Manual for Sterile Containers with Accessories», «Instruction Manual for reusable surgical instruments» spread for medical containers with accessories.

«Instruction Manual for reusable surgical instruments» Instructions for Use and Sterilization Instructions for Orthopaedic Implants» cover additional accessories for surgical instruments, orthopedic and for osteosynthesis.

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Products

This instruction manual is valid for the following products:

REF

Suctioning and rinsing instruments
art. No.
04.126.00 - 04.184.03
22.210.00 - 22.375.12
40.268.11 - 40.333.20
44.099.15 - 44.105.24
46.327.05 - 46.348.19
70.590.04 - 70.810.14
76.085.01 - 76.093.11

Important Information

Read this instruction manual carefully before every application and keep it easily accessible for all users or the respective specialist staff. Carefully read the warnings marked with this symbol. Improper use of the products may result in serious injuries to the patient, the users or third parties.

1 Scope

The instruments must be used according to their intended use in the medical fields and by respectively trained and qualified staff only. The treating physician and/or user is responsible for choosing the equipment for specific applications and/or operative use, for the appropriate training and information, and for the sufficient experience regarding the handling of the equipment.

2 Intended use / Indication

Suction-irrigation instruments can be connected to a suction-irrigation pump for the suction of body fluids or for rinsing and cleaning body cavities during a surgical intervention or during treatment.

Contraindication

The instruments must be used according to their intended use in the medical fields and by respectively trained and qualified staff only. The instruments of Dimeda GmbH may not be used in the areas of cardiothoracic surgery and neurosurgery.

Materials

Material 1.4301 (AISI 304)
Material 1.4305 (AISI 303)

3 Precautions and Warnings

Attention!
The suctioning- and rinsing instruments are only designed for surgical use and must not be used for any other purpose. Improper handling and maintenance as well as misappropriated use may result in the premature wear and tear of the instruments.

Application period
The suctioning and rinsing instrument intended application duration is 60 minutes. The product must be replaced after 60 minutes.

Material intolerance
Under no circumstances must the instruments be used if the user or specialist staff become aware of the patient being intolerant to the material.

Functional Impairment
Surgical instruments corrode and become impaired in their functionality if they come into contact with aggressive substances. It is therefore necessary to observe the storage and sterilisation instructions.

Operating Conditions
The aforementioned products require correct maintenance and care in order to guarantee that the products operate safely. In addition to this, functionality testing and a visual check should be performed prior to each application. For this reason, please pay attention to the respective chapters in this inst

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Storage
There are no specific storage requirements concerning the products. Nevertheless, we recommend storing medical products in a clean and dry environment.

Creutzfeldt Jakob Disease
With regard to the reprocessing of medical devices that have been used on patients or suspected patients suffering from or suspected of suffering from Creutzfeldt-Jacob disease (CJD) or its variant (vCJD), the requirements specified in the corresponding appendix of the guidelines for hospital hygiene and infection prevention and the requirements specified by publications in the Federal Health Gazette must be adhered to. The medical devices that were used on this group of patients must be disposed of by incineration (European Waste Catalogue EAK 18 01 03) without risk. Dry heat, ethanol, formaldehyde and glutaraldehyde have a fixing but no inactivating effect on TSE pathogens. Of the sterilization methods available, only steam sterilization (especially 134°C, 18 minutes) has been shown to have a limited effect.

4 Liability and Warranty

As a manufacturer, Dimeda Instrumente GmbH is not liable for consequential damage resulting from improper use or handling. This particularly applies to use which is not compliant with the defined intended use, or non-compliance with the instructions on preparation and sterilisation. This also applies to repairs or changes to the product which are not carried out by authorized staff of the manufacturer. These disclaimers also apply to warranty services.

5 Sterility

State upon Delivery
Medical products are delivered in a non-sterile condition and need to be prepared and sterilised by the user prior to the first application and any subsequent application according to the following instructions.

6 Reprocessing

Warnings
- Frequent reprocessing impairs the quality of the products.
- Any tap water to be used must comply with COUNCIL DIRECTIVE 98/83/EC of 3 November 1998 on the quality of water intended for human consumption.
- The detergents and disinfectants used in order to maintain validation are listed in this instruction for reprocessing. Should an alternative detergent or disinfectant (listed at the RKI or VAH) be used, the responsibility lies with the reprocessor.
- Reassemble disassembled products prior to sterilisation.

Use Site
The first steps of a proper reprocessing take place in the operating theatre. Coarse contaminations must be removed prior to storing the instruments if possible. For this purpose, the instruments should be rinsed under cold tap water (<40°C). If this procedure is not sufficient to remove the obvious soiling, a soft plastic brush can be used to remove soiling.
Whenever possible, dry removal (moistened, closed system) should be the method of choice. A drying of any residues should be avoided! Wherever possible, dry disposal is to be preferred, since with wet disposal the prolonged lying of the medical devices in solutions can lead to material damage (e.g. corrosion). Long periods of waiting until the reprocessing, for instance overnight or over the weekend, must be avoided with both types of removal (<60 minutes).

Transport
The products must be disposed of in a dry state immediately (<60 min) after use, if possible. This means that the products have to be transported in a closed container from the place of application to the purification, so that the products do not dry up.

Preparing the Decontamination
The products must be disassembled prior to the following reprocessing steps and/or must be exposed to the following reprocessing steps in an open condition, where possible. Rinse residue must be avoided. The products must be reprocessed in appropriate screen baskets or rinsing shields (choose size according to product). The products must be positioned in the cleaning basket at a minimum clearance from one another. Avoid overlapping so that the damaging of the products during the cleaning process can be excluded.

Pre-cleaning
Rinse the products under cold tap water suitable for use as drinking water (<40°C) until all visible contaminations are removed. Remove persistent dirt with a soft brush. Movable parts on the instrument need to be moved. Cavities, lumen, gaps and slits need to be rinsed intensively (>60 sec) with cold tap water suitable for use as drinking water (<40°C) using a water pistol (or similar). Insertion of the products into an ultrasonic bath (<40°C) with an alkaline cleaner (0.5% Neodisher® MediClean forte, Dr. Weigert), time of treatment 5 min. at a frequency of about 35 kHz. For this please observe the instructions of the cleaning agent manufacturer. Rinse instruments briefly (<15 sec) under cold water. Movable parts must be moved. Respectively rinse cavities, lumen, gaps and slits again (>30 sec.) with cold city water (<40°C) using a water pistol (or similar).

Cleaning/disinfection
Automated cleaning and/or disinfection process
(*Miele Disinfector G7835 CD as per ISO 15883*):
- Pre-clean for 1 minute
- Drain water
- Pre-clean for 4 minutes
- Drain water
- Clean for 6 minutes at 58°C +/- 1°C using 0.5 % alkaline detergent (*0,5 % Neodisher Mediclean forte*)
- Drain water
- Neutralize for 3 minutes (0.1% *0,1 % NeodisherZ*) with cold tap water suitable for use as drinking water <40°C
- Drain water
- Clean for 2 minutes with FD water <40°C

Please observe the special instructions of the manufacturer of the automatic cleaning machine.

Automated Disinfection
Automated thermal disinfection in a cleaning and disinfection device taking into consideration the national requirements for the A_0 value; for instance, A_0 value 3000: >5 minutes at >95°C.

Automated Drying
Automated drying in accordance with the drying operation of the cleaning and disinfection device for at least 30 minutes at 92°C +/- 2°C.

Sterilisation
(Typ B Autoklav by Tuttmauer as per DIN EN 13060)

Sterilisation of products with a fractionated pre-vacuum method (according to DIN EN ISO 17665-1) taking into consideration the respective national requirements. The sterilisation of the products must be conducted in suitable sterilisation packaging according to DIN EN ISO 11607-1 and EN 868
The sterilisation must be completed using a fractionated pre-vacuum method with the following parameters:
134°C,
5 minutes hold time,
3 pre-vacuum cycles
Drying in vacuum for at least 20 minutes
The instruction manual of the manufacturer of the autoclave and the recommended directions for maximum loading with goods to be sterilised must be observed. The autoclave must be installed, maintained, validated and calibrated in accordance with requirements.

Additional Information
The reprocessor is responsible for ensuring that the actual reprocessing, including the used equipment, materials and the staff involved in the reprocessing facility, achieves the desired results. This typically requires the validation and routine monitoring of the method and the equipment used.

7 Lifespan of the Products

Frequent reprocessing has little effect on the surgical instruments. Product life is limited by wear due to proper use and damage to the instrument. The product may no longer be used under the following conditions, among others:
Corrosion, damage, fractures, cracks, deformation, porosity, functional limitations, products with unrecognizable or missing marking. For this reason, the conditioner must comply with the relevant instructions for functional testing (section 8).

8 Functional Inspection

Check products after reprocessing and before sterilization with regard to the following aspects:

- Cleanliness
- Damage, including but not limited to signs of corrosion (rust, pitting), discoloration, deep scratches, flaking, cracking and wear.
- Proper function
- Missing or removed (ground off) part numbers.
- Do not use improperly functioning or defective and excessively worn products, products with

unrecognizable markings, missing or removed (ground) part numbers.
- Check products for perfect surfaces, correct assembly and functionality. Do not use severely damaged products.
- Reassemble dismantled products before sterilization.

9 Service and Repair

⚠ Service and Repair

Do not carry out any repairs or changes to the product yourself. Authorized staff of the manufacturer are solely responsible for such work. Should you wish to make complaints or queries, or offer us any advice regarding our products, please feel free to contact us.

⚠ Returns

Defective or non-compliant products must go through the entire reprocessing process before being sent back for repairs/service.

10 Packaging, Storage and Disposal

Standard packaging of the products for sterilisation according to ISO 11607 and EN 868.
Store sterile products in a dry, clean, and dust-free environment, secured against damage, at moderate temperatures.
The medical products of the manufacturer should be stored and kept in single packaging, boxes or protective containers. Please handle the instruments with care during transportation, storage and reprocessing. The user and/or specialist staff intended for this is responsible for ensuring that the sterile state is maintained after the sterilisation.
The disposal of the products, packaging as well as the accessories must be performed in accordance with current national rules and laws. No specific instruction regarding this matter is provided by the manufacturer.

11 Combination products

Suction and rinsing pumps:

Please refer to the manufacturer's instructions to ensure that these are equipped with the compatible connections and that the suction and rinsing performance corresponds to the area of application.

Connections:

04.150.27/04.150.31/04.150.80 to be used with:

04.150.98, 04.150.99

04.150.80 to be used with:

04.150.90, 04.150.91, 04.150.92, 04.150.93,

04.150.97

04.165.00 to be used with:

04.165.01, 04.165.02, 04.165.03, 04.165.04

To be used with a hose diameter of 6mm:

04.165.11, 04.165.12, 04.165.13, 04.165.14, 04.171.03, 04.171.04, 04.171.05, 04.171.06, 04.171.15, 04.171.25, 04.172.02 ,04.172.03, 04.172.04, 04.172.05, 04.172.06, 04.172.15, 04.172.25

To be used with a hose diameter of 6-7mm:

04.160.21, 04.161.21

To be used with a hose diameter of 6-9mm:

04.126.00, 04.130.10, 04.130.30, 04.130.38, 04.140.06, 04.140.08, 04.141.35, 04.143.31, 04.145.30, 04.147.33, 04.147.90, 04.148.27, 04.149.27, 04.150.22, 04.150.27, 04.150.27AT, 04.150.27VOM, 04.150.31, 04.150.80, 04.150.99VOM, 04.152.33, 04.152.99, 04.153.03, 04.153.04, 04.153.24, 04.154.16, 04.154.26, 04.154.27, 04.156.04, 04.156.05, 04.165.00, 04.170.01, 04.170.02, 04.170.03, 04.182.01, 04.182.02, 04.184.00, 04.184.01, 04.184.02, 04.184.03, 40.268.11, 40.269.07, 40.269.09, 40.269.11, 40.269.12, 40.269.15, 40.270.15, 40.270.17, 40.270.20, 40.270.25, 40.270.30, 40.270.40, 40.270.50, 40.271.15, 40.271.20, 40.271.25, 40.271.30, 40.271.40, 40.271.50, 40.272.15, 40.272.20, 40.272.25, 40.272.30, 40.272.40, 40.272.50, 40.273.15, 40.273.20, 40.273.25, 40.273.30, 40.273.40, 40.273.50, 40.275.05, 40.275.06, 40.275.07, 40.275.08, 40.275.09, 40.275.10, 40.275.12, 40.275.15, 40.282.30, 40.282.35, 40.282.45, 40.283.02, 40.283.03, 40.283.70, 40.283.71, 40.283.72, 40.283.73, 40.283.74, 40.283.75, 40.284.70, 40.286.13, 40.286.15, 40.286.18, 40.286.22, 40.287.13, 40.287.15, 40.287.18, 40.287.22, 40.288.13, 40.288.22, 40.289.13, 40.289.15, 40.289.18, 40.289.22, 40.289.45, 40.290.03, 40.290.04, 40.290.05, 40.290.06, 40.290.07, 40.290.08, 40.290.09, 40.290.10, 40.290.11, 40.290.12, 40.290.14, 40.290.16, 40.290.21, 40.292.03, 40.292.04, 40.292.05, 40.292.06, 40.292.07, 40.292.08, 40.292.08TU, 40.292.09, 40.292.10, 40.292.10TU, 40.292.11, 40.292.12, 40.292.12TU, 40.292.15, 40.293.06, 40.293.08, 40.293.09, 40.293.10, 40.293.12, 40.295.06, 40.295.07, 40.295.08, 40.295.09, 40.295.10, 40.295.11, 40.295.12, 40.295.15, 40.295.46, 40.295.50, 40.295.51, 40.295.99, 40.296.06, 40.296.16, 40.296.18, 40.296.30, 40.296.32, 40.297.06, 40.297.07, 40.297.08, 40.297.09, 40.297.10, 40.297.11, 40.297.12, 40.297.15, 40.306.01, 40.306.02, 40.306.03, 40.306.04, 40.306.05, 40.306.10, 40.306.11, 40.306.12, 40.307.45, 40.307.55, 40.308.00, 40.308.01, 40.308.04, 40.308.06, 40.308.08, 40.309.00, 40.309.01, 40.309.11, 40.310.00, 40.310.01, 40.310.02, 40.310.03, 40.310.04, 40.313.15, 40.313.20, 40.313.25, 40.313.30, 40.320.09, 40.320.11, 40.321.09, 40.321.15, 40.323.03, 40.323.04, 40.323.05, 40.323.06, 40.323.07, 40.323.08, 40.323.09, 40.323.12, 40.324.02, 40.324.03, 40.324.04, 40.324.05, 40.324.06, 40.324.07, 40.324.08, 40.324.09, 40.324.10, 40.324.12, 40.325.03, 40.325.04, 40.325.05, 40.325.06, 40.325.07, 40.325.08, 40.325.09, 40.325.12, 40.330.09, 40.330.17, 40.331.09, 40.331.17, 40.332.09, 40.332.17, 40.333.09, 40.333.20, 44.101.03, 44.101.05, 44.101.07, 44.101.13, 44.101.15, 44.101.17, 44.102.02, 44.102.03, 44.102.20, 44.102.24, 44.103.01, 44.103.02, 44.103.03, 44.103.04, 44.103.11, 44.103.12, 44.103.13, 46.327.05, 46.327.07, 46.327.09, 46.341.02, 46.341.03, 46.341.05, 46.348.19, 70.671.01, 70.671.02, 70.671.03, 70.810.50

To be used with a hose diameter of 8-9mm:

04.155.30, 04.155.35, 04.155.99, 40.276.04

To be used with a hose diameter of 9-13mm:

04.158.00

Pipes with Luer Lock connection:

22.210.00, 22.220.03, 22.220.04, 22.220.05, 22.220.06, 22.370.02, 22.370.03, 22.370.04, 22.370.05, 22.370.12, 22.370.15, 22.370.23, 22.370.37, 22.372.06, 22.373.04, 22.375.08, 22.375.12, 40.294.05, 40.305.06, 40.305.07, 40.305.08, 40.305.09, 40.305.10, 40.305.12, 40.305.15, 40.305.16, 40.305.18, 40.305.20, 40.305.25, 40.305.30, 40.306.06, 40.306.07, 40.309.02, 40.309.03, 40.309.04, 40.309.05, 40.332.20, 44.099.15, 44.099.17, 44.099.20, 44.105.16, 44.105.18, 44.105.20, 44.105.22, 44.105.24, 46.330.25, 70.590.04, 70.592.01, 70.592.02, 70.592.03, 70.592.04, 76.087.47, 76.087.48, 76.087.49, 76.087.57, 76.087.58, 76.087.59, 76.087.60, 76.087.61, 76.087.62, 76.087.63, 76.087.64, 76.087.65, 76.087.66, 76.087.67, 76.087.68, 76.087.69, 76.087.81, 76.087.85, 76.087.89

Pipes with Luer Lock connection / Suction tube with hose attachement Ø 6-9mm

46.330.30, 46.330.40, 46.331.25, 46.331.30, 46.331.40, 46.332.25, 46.332.30, 46.332.40, 46.333.25, 46.333.30, 46.333.40, 46.334.25, 46.334.30, 46.334.40

Pipes/grips with Luer Lock connection:

76.085.01, 76.085.02, 76.085.03, 76.085.04, 76.085.05, 76.085.06, 76.085.07, 76.085.08, 76.085.09, 76.085.10, 76.085.11, 76.085.12, 76.086.03, 76.086.31, 76.086.32, 76.086.58

Pipes/grips with Luer Lock connection / Hose attachement with Ø 6-9mm

76.085.90, 76.085.91

To be used with a hose diameter of 13mm:

76.087.45, 76.087.46, 76.093.01, 76.093.02, 76.093.03, 76.093.04, 76.093.05, 76.093.06, 76.093.07, 76.093.08, 76.093.09, 76.093.10, 76.093.11

76.087.05/76.087.43/76.087.45/ 76.087.46 To be used with:

76.086.70, 76.086.72, 76.086.74, 76.086.87, 76.086.88, 76.087.00, 76.087.01, 76.087.02, 76.087.03, 76.087.04, 76.087.07, 76.087.08, 76.087.09, 76.087.10, 76.087.11, 76.087.12, 76.087.14, 76.087.15, 76.087.16, 76.087.17, 76.087.19, 76.087.20, 76.087.21, 76.087.22, 76.087.23, 76.087.24, 76.087.25, 76.087.26, 76.087.27, 76.087.28, 76.087.29, 76.087.31, 76.087.32, 76.087.33, 76.087.34, 76.087.37, 76.087.38, 76.087.40, 76.087.41, 76.087.50, 76.087.51, 76.087.52, 76.087.53, 76.087.54, 76.087.55, 76.087.70, 76.087.71, 76.087.72, 76.087.73, 76.087.74, 76.087.75, 76.087.76, 76.087.77, 76.087.78, 76.087.79, 76.087.87, 76.087.88, 76.088.07, 76.088.08, 76.088.17, 76.088.18, 76.088.30, 76.088.31, 76.088.32, 76.088.33, 76.088.34

70.810.50 to be used with:

70.810.03, 70.810.04, 70.810.06, 70.810.08, 70.810.10, 70.810.12, 70.810.14

12 Description of Symbols Used

⚠	Attention!
	Observe the instruction manual
REF	Item number
LOT	Lot designation
CE	CE labelling
	Indication of a non-sterile product
	Name and address of the manufacturer

dimeda®
SURGICAL INSTRUMENTS
DIMEDA Instrumente GmbH
Gänsäcker 54 - 58
D-78532 Tuttlingen/Germany
i.A.

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Instruction manual

For reusable surgical instruments CE

Valid from:	10.12.2020	Version:	5

dimeda® SURGICAL INSTRUMENTS

Manufacturer

Dimeda Instrumente GmbH
Gänsäcker 54+58
78532 Tuttlingen, Germany
Phone:+49 (0) 7462 / 9461-0
Fax: +49 (0) 7462 / 9461-33
http://www.dimeda.de
info@dimeda.de

Products

This instruction manual is valid for all reusable surgical instruments, which

- are one-piece instruments
- optionally have simple joints or
- include simple moving parts
- are optionally composed of multiple replaceable parts (e.g. handle and multiple attachments)

Important Information

Read this instruction manual carefully before every application and keep it easily accessible for all users or the respective specialist staff.

Carefully read the warnings marked with this symbol. Improper use of the products may result in serious injuries to the patient, the users or third parties.

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The instruments intended application duration is 60 minutes. The product must be replaced after 60 minutes.

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3 Precautions and Warnings

Attention!

The instruments are only designed for surgical us
other purpo

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maintenance as well as misappropriated use may result in the premature wear and tear of the instruments.

Material intolerance

Under no circumstances must the instruments be used if the user or specialist staff become aware of the patient being intolerant to the material.

Functional Impairment

Surgical instruments corrode and become impaired in their functionality if they come into contact with aggressive substances. It is therefore necessary to observe the storage and sterilisation instructions.

Operating Conditions

The aforementioned products require correct maintenance and care in order to guarantee that the products operate safely. In addition to this, functionality testing and a visual check should be performed prior to each application. For this reason, please pay attention to the respective chapters in this instruction manual.

Combination with other products

Should the products be reassembled after disassembly, individual parts must not be replaced with parts from other manufacturers! If the intended purpose of the product entails certain parts being exchanged (e.g. different attachments), no parts from different manufacturers must be used! We recommend to also purchase other accessories (e.g. detergents) at Dimeda Instrumente GmbH.

Storage

There are no specific storage requirements concerning the products. Nevertheless, we recommend storing medical products in a clean and dry environment.

Creutzfeldt Jakob Disease

With regard to the reprocessing of medical devices that have been used on patients or suspected patients suffering from or suspected of suffering from Creutzfeldt-Jacob disease (CJD) or its variant (vCJD), the requirements specified in the corresponding appendix of the guidelines for hospital hygiene and infection prevention and the requirements specified by publications in the Federal Health Gazette must be adhered to. The medical devices that were used on this group of patients must be disposed of by incineration (European Waste Catalogue EAK 18 01 03) without risk. Dry heat, ethanol, formaldehyde and glutaraldehyde have a fixing but no inactivating effect on TSE pathogens. Of the sterilization methods available, only steam sterilization (especially 134°C, 18 minutes) has been shown to have a limited effect.

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5 Sterility

State upon Delivery

Medical products are delivered in a non-sterile condition and need to be prepared and sterilised by the user prior to the first application and any subsequent application according to the following instructions.

6 Reprocessing

Warnings

- Frequent reprocessing impairs the quality of the products.
- Any tap water to be used must comply with COUNCIL DIRECTIVE 98/83/EC of 3 November 1998 on the quality of water intended for human consumption.
- The detergents and disinfectants used in order to maintain validation are listed in this instruction for reprocessing. Should an alternative detergent or disinfectant (listed at the RKI or VAH) be used, the responsibility lies with the reprocessor.
- Reassemble disassembled products prior to sterilisation.
- The preparation may only be carried out by qualified medical personnel. Machine reprocessing must be qualified and validated by the user. The washer-disinfectors must fully comply with the requirements of DIN EN ISO 15883.

Use Site

The first steps of a proper reprocessing take place in the operating theatre. Coarse contaminations must be removed prior to storing the instruments if possible. For this purpose, the instruments should be rinsed under cold tap water (<40°C). If this procedure is not sufficient to remove the obvious soiling, a soft plastic brush can be used to remove soiling.
Whenever possible, dry removal (moistened, closed system) should be the method of choice. A drying of any residues should be avoided! Wherever possible, dry disposal is to be preferred, since with wet disposal the prolonged lying of the medical devices in solutions can lead to material damage (e.g. corrosion). Long periods of waiting until the reprocessing, for instance overnight or over the weekend, must be avoided with both types of removal (<60 minutes).

Transport

The products must be disposed of in a dry state immediately (<60 min) after use, if possible. This means that the products have to be transported in a closed container from

the place of application to the purification, so that the products do not dry up.

Preparing the Decontamination

The products must be disassembled prior to the following reprocessing steps and/or must be exposed to the following reprocessing steps in an open condition, where possible. Rinse residue must be avoided. The products must be reprocessed in appropriate screen baskets or rinsing shields (choose size according to product). The products must be positioned in the cleaning basket at a minimum clearance from one another. Avoid overlapping so that the damaging of the products during the cleaning process can be excluded.

Pre-cleaning

Rinse the products under cold tap water suitable for use as drinking water (<40°C) until all visible contaminations are removed. Remove persistent dirt with a soft brush. Movable parts on the instrument need to be moved. Cavities, lumen, gaps and slits need to be rinsed intensively (>60 sec) with cold tap water suitable for use as drinking water (<40°C) using a water pistol (or similar). Insertion of the products into an ultrasonic bath (<40°C) with an alkaline cleaner (0.5% Neodisher® MediClean forte, Dr. Weigert), time of treatment 5 min. at a frequency of about 35 kHz. For this please observe the instructions of the cleaning agent manufacturer. Rinse instruments briefly (<15 sec) under cold water. Movable parts must be moved. Respectively rinse cavities, lumen, gaps and slits again (>30 sec.) with cold city water (<40°C) using a water pistol (or similar).

Cleaning/disinfection

Automated cleaning and/or disinfection process

(Miele Disinfector G7835 CD as per ISO 15883):

- Pre-clean for 1 minute
- Drain water
- Pre-clean for 4 minutes
- Drain water
- Clean for 6 minutes at 58°C +/- 1°C using 0.5 % alkaline detergent (0,5 % Neodisher Mediclean forte)
- Drain water
- Neutralize for 3 minutes (0.1% 0,1 % NeodisherZ) with cold tap water suitable for use as drinking water <40°C
- Drain water
- Clean for 2 minutes with FD water <40°C

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Automated Disinfection

Automated thermal disinfection in a cleaning and disinfe[illegible] considerat[illegible] the A_0 valu[illegible]

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Automated Drying

Automated drying in accordance with the drying operation of the cleaning and disinfection device for at least 30 minutes at 92°C +/- 2°C.

Sterilisation

(Typ B Autoklav by Tuttmauer as per DIN EN 13060)

Sterilisation of products with a fractionated pre-vacuum method (according to DIN EN ISO 17665-1) taking into consideration the respective national requirements. The sterilisation of the products must be conducted in suitable sterilisation packaging according to DIN EN ISO 11607-1 and EN 868

The sterilisation must be completed using a fractionated pre-vacuum method with the following parameters:

134°C,

5 minutes hold time,

3 pre-vacuum cycles

Drying in vacuum for at least 20 minutes

The instruction manual of the manufacturer of the autoclave and the recommended directions for maximum loading with goods to be sterilised must be observed. The autoclave must be installed, maintained, validated and calibrated in accordance with requirements.

⚠ **Additional Information**

The reprocessor is responsible for ensuring that the actual reprocessing, including the used equipment, materials and the staff involved in the reprocessing facility, achieves the desired results. This typically requires the validation and routine monitoring of the method and the equipment used.

7 Lifespan of the Products

Frequent reprocessing has little effect on the surgical instruments. Product life is limited by wear due to proper use and damage to the instrument. The product may no longer be used under the following conditions, among others:

Corrosion, damage, fractures, cracks, deformation, porosity, functional limitations, products with unrecognizable or missing marking. For this reason, the conditioner must comply with the relevant instructions for functional testing (section 8).

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- Cleanliness
- Damage, including but not limited to signs of corrosion (rust, pitting), discoloration, deep scratches, flaking, cracking and wear.
- Proper function
- Missing or removed (ground off) part numbers.
- Do not use improperly functioning or defective and excessively worn products, products with unrecognizable markings, m or removed (ground) part numbers.
- Check products for perfect surfaces, correct assembly a functionality. Do not use sev damaged products.
- Reassemble dismantled proc before sterilization.

9 Service and Repair

⚠ **Service and Repair**

Do not carry out any repairs or changes the product yourself. Authorized staff of manufacturer are solely responsible for work. Should you wish to make compla queries, or offer us any advice regardin products, please feel free to contact us.

⚠ **Returns**

Defective or non-compliant products mu through the entire reprocessing process before being sent back for repairs/servic

10 Packaging, Storage and Disposal

Standard packaging of the products for sterilisation according to ISO 11607 and 868.

Store sterile products in a dry, clean, an dust-free environment, secured against damage, at moderate temperatures.

The medical products of the manufactur should be stored and kept in single packaging, boxes or protective containe Please handle the instruments with care during transportation, storage and reprocessing. The user and/or specialist intended for this is responsible for ensur that the sterile state is maintained after t sterilisation.

The disposal of the products, packaging well as the accessories must be perform accordance with current rules and laws. specific instruction regarding this matter provided by the manufacturer.

11 Description of Symbols Us

⚠	Attention!
	Observe the instruction manual
REF	Item number
LOT	Lot designation
CE	CE labelling
	Indication of a non-steril product
	Name and address of the manufacturer

dimeda®
SURGICAL INSTRUMENTS
DIMEDA Instruments GmbH
Gänsäcker 61 [illegible]
D-78532 Tuttlingen / Germany

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Instructions For Use

1.	Intended use – important information for the user	2
2.	Product description	3
3.	Intended purpose	3
4.	Implant selection	3
5.	Indications	4
6.	Contraindications	4
7.	Potential adverse effects	4
8.	Correct handling	4
9.	Postoperative care:	5
10.	Compatibility:	5
11.	Single-use products	6
12.	Device reprocessing	6
13.	Validated machine reprocessing methods	6
a)	Machine cleaning, disinfection and drying in the washer-disinfector	6
b)	Visual inspection	7
c)	Packaging	7
d)	Validated sterilisation method	7
14.	Symbols for labelling	7
15.	Storage	7
16.	Guarantee policy	8
17.	Bibliography	8

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Please read these instructions in full and observe them.

1. Intended use – important information for the user

These instructions for use contain general instructions on the application and use osteosynthesis implants from Dimeda Instrumente GmbH.
Purchasing this implant provides you with a high quality product whose appropriate and cor handling are described below.
To minimise risks and unnecessary strain on the patient, please carefully read and save th instructions for use.
⚠ The products are non-sterile disposable medical devices (single-use products ⊗). Be use, the implants must be reprocessed in accordance with the specifications defined in these II

Item name:

Orthopaedic implants incl. accessories:
- Drill wires, Steinmann pins,
- Intramedullary pins
- K-wires
- Bone wires (cerclage)

Dimensions/specifications: See table below

Item number	Item name	UMDNS code	Dimensions
33.xxx.xx	K-wires / drill wires made of implant steel	16-104	L=50-600mm D=0.6-3.0 mm
33.xxx.xx	Steinmann pins made of implant steel	16-078	L=50-600mm D=3.2-6.0mm
33.xxx.xx	Bone soft wire (cerclage) made of implant steel	16-104	L=10m D=0.2-2.0mm

Material:

Material	Material specification	Standard
Implant steel 1.4441	X2CrNiMo 18-15-3	DIN ISO 5832-1

Before using an Dimeda implant, the user must carefully read the following instructions use and comply with its recommendations, warnings and instructions.
Dimeda cannot be held liable for complications arising due to the use of implants/instrume outside of the control of Dimeda, including but not limited to product selection and deviations fr use/handling and surgical technique.

⚠ The implants are **not** intended for use on the central nervous/circulatory system and m not be used for this purpose.

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2. Product description

Orthopaedic implants offer orthopaedic and trauma surgeons a means for fixing bones precisely. They play a supportive role in the treatment and healing of bone fractures (osteosynthesis and correction of degenerative diseases). However, implants are not suitable for replacing normal body structures or bearing the whole body weight.

3. Intended purpose

- K-wires (drill wires)
 intended for use in the closed reduction and fixation of fractures by means of a rotating drill wire (K-wire). This includes surgical fracture treatment by means of percutaneous intramedullary splints (e.g., at the metacarpals) or percutaneous "pinning", fracture fixation by K-wire insertion, if possible with fixation of the wire in the contralateral cortical bone.
- Steinmann pins
 intended for use in the extension treatment of fractures. Extension treatment is based on applying continuous longitudinal traction to the injured limb. Depending on the fracture to be treated, a drill wire must be applied transversely through the bone, and longitudinal traction is applied with the aid of a metal clamp and a variable weight.
- Bone wire (cerclage wire)
 intended for fracture care by means of simple cerclage wiring as stand-alone treatment. The bone is encircled by the soft wire, and tension is applied by twisting.
- Intramedullary pins
 intended for treating shaft fractures and some metaphyseal and epiphyseal fractures.

4. Implant selection

The following factors must be considered in fracture treatment:

Selection of the correct implant:

Choosing the proper implant is extremely important for successful treatment. The potential for success is increased by selecting the proper implant size and shape. The characteristics of human bone and soft tissue pose restrictions on the size and strength of implants. If a strong bone union is to be achieved, the patient requires adequate external assistance. Physical strain and load on the fracture site must be limited to prevent delayed healing and/or late sequelae.

Patient-related factors

a) Weight:
Patient overweight or obesity can adversely affect the implant and its stability. The load limits must be considered individually depending on the implant site and purpose.

b) Occupation or activities:
Occupational activities pose a risk to the healing process if they are associated with the application of considerable force and hence physical loads on the body. In that case, immobilisation is necessary to ensure healing.

c) Senility, mental illness or alcoholism:
There is a risk of affected patients ignoring certain necessary limitations and precautions, leading to the failure of the product or other complications.

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2. Product description

Orthopaedic implants offer orthopaedic and trauma surgeons a means for fixing bones precisely. They play a supportive role in the treatment and healing of bone fractures (osteosynthesis and correction of degenerative diseases). However, implants are not suitable for replacing normal body structures or bearing the whole body weight.

3. Intended purpose

- K-wires (drill wires)
 intended for use in the closed reduction and fixation of fractures by means of a rotating drill wire (K-wire). This includes surgical fracture treatment by means of percutaneous intramedullary splints (e.g., at the metacarpals) or percutaneous "pinning", fracture fixation by K-wire insertion, if possible with fixation of the wire in the contralateral cortical bone.
- Steinmann pins
 intended for use in the extension treatment of fractures. Extension treatment is based on applying continuous longitudinal traction to the injured limb. Depending on the fracture to be treated, a drill wire must be applied transversely through the bone, and longitudinal traction is applied with the aid of a metal clamp and a variable weight.
- Bone wire (cerclage wire)
 intended for fracture care by means of simple cerclage wiring as stand-alone treatment. The bone is encircled by the soft wire, and tension is applied by twisting.
- Intramedullary pins
 intended for treating shaft fractures and some metaphyseal and epiphyseal fractures.

4. Implant selection

The following factors must be considered in fracture treatment:

Selection of the correct implant:

Choosing the proper implant is extremely important for successful treatment. The potential for success is increased by selecting the proper implant size and shape. The characteristics of human bone and soft tissue pose restrictions on the size and strength of implants. If a strong bone union is to be achieved, the patient requires adequate external assistance. Physical strain and load on the fracture site must be limited to prevent delayed healing and/or late sequelae.

Patient-related factors

a) Weight:
Patient overweight or obesity can adversely affect the implant and its stability. The load limits must be considered individually depending on the implant site and purpose.

b) Occupation or activities:
Occupational activities pose a risk to the healing process if they are associated with the application of considerable force and hence physical loads on the body. In that case, immobilisation is necessary to ensure healing.

c) Senility, mental illness or alcoholism:
There is a risk of affected patients ignoring certain necessary limitations and precautions leading to the failure of the product or other complications.

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d) Degenerative diseases and nicotine use:
In some cases, degenerative disease may be so advanced at the tim implantation that it may substantially decrease the expected useful life o implant. In that case, the products serve only as a means of delaying or tempo alleviating the disease.

e) Sensitivity to foreign bodies:
Where sensitivity or allergy to the materials used in the implant is suspe appropriate tests must be conducted prior to implant selection and placement.

5. *Indications*

Fixation of bones and bone fragments following successful reduction.

6. *Contraindications*

- All concomitant diseases that may pose a risk to fixation or the success of surge negatively affect them, such as obesity, impairment of the circulation, etc.
- Poor bone quality or quantity that poses a risk for or adversely affects the secure fix of implants
- Severe muscular, neurological or vascular disease that poses a risk for or adve affects the success of the procedure/surgery
- Allergy patients who are susceptible to developing reactions when coming into co with materials used in the implant
- Acute or chronic, local or systematic infections
- Nicotine use that may pose a risk for the success of the procedure/surgery due to del bone/wound healing
- Mental conditions that make it impossible to understand and follow physician instruc and/or to participate in the rehabilitation program (e.g., alcohol or drug use, Parkins disease, Alzheimer's disease, etc.)

7. *Potential adverse effects*

The adverse effects listed below have been described in the literature:

- Loosening of implants
- Wound infection (skin and deep wound infection)
- Vascular complications
- Nonunions
- Nerve damage
- Inflammations
- Metal allergies

8. *Correct handling*

Implants must be treated with the care required when handling medical devices. If the shap the implant has to be altered, it should not be bent excessively, bent backward, notche scratc tions, in combination with improper handling and application, lead to surface defects and/or structural changes in the material and to product damage a failure.

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9. Postoperative care:

Doctors must inform their patients about the load restrictions of the implant and give instruction on postoperative behaviour and on gradually increasing physical loads. Failure to do so ca generate malalignment, delayed bone healing, implant failure, infection, thrombophlebitis and/c wound haematoma.
The doctor takes the final decision on when to remove the implant. If possible and appropriat for the individual patient, we recommend removing fixation products after healing is complete This particularly applies to young and active patients. The risk of negative effects, such a secondary infection, allergies, material fatigue fractures, implant failure and/or impaired bloo circulation, increases with the length of time the implant remains in the body.

10. Compatibility:

Dimeda does not recommend the combination of Dimeda products with those of othe manufacturers since the designs, materials, mechanics and construction are not mutuall harmonised. Dimeda assumes no liability for any complications arising from combinin components or from using third-party medical devices in combination.
Unless described otherwise, mixing different implant metals is not recommended. Mixing o metals may lead to galvanic corrosion and the release of ions. This may cause inflammatory response, metal sensitivity reactions and/or long-term adverse systemic effects. In addition, the corrosion process can reduce the mechanical strength of the implant.

⚠ ***Warning:***

- The implants are intended for single use. They must not be reused after being used once Dimeda assumes no liability in case of non-compliance.
- Implants that have come into contact with blood, soft tissue, bone or body fluids must no be reused and must be disposed of following the disposal guidelines for contaminate products. Contamination residues on the implants may lead to patient or user injuries o infections.
- Before use, the user must check the implants for visible damage such as cracks, fracture or damaged tips. Damaged implants must not be placed.
- The products must be handled exclusively by trained medical staff. They must be use and implanted exclusively by appropriately trained doctors.
- Incorrectly selected, positioned or sized implants or incorrect fixation may lead to unusua stress that may negatively affect the life of the implants.
- The implants must be used exclusively for the intended indication. These implants must not be used for other indications (off-label use).

⚠ ***Precautions:***

- The patient should be regularly monitored and tested for infection for as long as the implant remains in the body.
- These implants have been developed for temporary use and should be removed afte complete healing of the fracture.
- The implants must not be machined or modified.
- Unless the implant is labelled as “MR Safe” or “MR Conditional”, there are hazard associated with the use of Dimeda products in an MRI environment. This includes but i not limited to:
 - Heating and/or migration of the implant

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- Artefacts created by the implant

Dimeda implants have not been tested for MRI compatibility.

11. Single-use products

Products intended for single use must not be reused.
The reuse or clinical reprocessing of single-use products may compromise the structural int of the device and lead to device failure, causing patient injury, illness or death. Furthermor reuse or clinical reprocessing of single-use devices increases the risk of contamination, e.g to the transmission of pathogens from one patient to another. This can also lead to patient i illness or death.
Dimeda advises against the clinical reprocessing of contaminated implants. Dimeda de contaminated by blood, tissue and/or body fluids or substances must not be reused and mu disposed of in accordance with hospital guidelines and protocol. Even if components a externally intact after use, minor defects and invisible material damage may cause ma fatigue.

12. Device reprocessing

Dimeda products are supplied in non-sterile condition and must be cleaned, disinfecte steam sterilised prior to surgical use.
Before reprocessing, the implants must be removed from the original packaging; the pro are not intended to be sterilised in their original packaging.

Dimeda recommends machine reprocessing with a standard cleaning program using a wa disinfector in accordance with ISO 15883-2.
Manual reprocessing methods are unsuitable due to the implant design (threads, holes, dri etc.) and have therefore not been validated.

13. Validated machine reprocessing methods

Validated machine reprocessing methods include:

a) Machine cleaning, disinfection and drying in the washer-disinfector
b) Visual check
c) Packaging
d) Validated sterilisation method

a) Machine cleaning, disinfection and drying in the washer-disinfector

The validated cleaning and disinfection method refers to the Miele standard program "DES-V TD" in the Miele G7835 CD washer-disinfector.

The specifications of the washer-disinfector manufacturer regarding appropriate and professi operation and loading as well as maintenance of the washer-disinfector must be strictly obse

Procedure:

- Pr[illegible]ter for 1 minute
- Cl[illegible]nd the alkaline cleaner "Neodisher Mediclean forte, 0.5%" for 5 minutes at 55°C ± 5°C

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- Neutralisation with the "Neodisher Z" neutraliser for 2 minutes
- Rinsing with deionised water for at least 1 minute
- Thermal disinfection at 55°C for 5 minutes
- Drying at 60°C ± 5°C for 30 minutes

b) Visual inspection

After cleaning, the implants must be visually inspected for cleanliness and damage. Implants that are not clean must be cleaned again, and damaged implants must be separated out and discarded.

c) Packaging

The instruments must be packaged for sterilisation in accordance with 11607-1. The validated sterilisation method applies to double foil bags.

d) Validated sterilisation method

The validated sterilisation method applies to the autoclave Tuttnauer Type B 3870 EHS.

Sterilisation:

- 2 fractionated pre-vacuum phases
- Holding time: min. 5 minutes, max. 7 minutes at 134°C
- Drying for at least 10 minutes

The unit manufacturer's operating instructions as well as instructions for use and maintenance must be strictly followed.

14. Symbols for labelling

Symbol	Explanation	Symbol	Explanation
	Observe instructions for use		Manufacturer
	For single use only		Batch name
	Attention!		Product is supplied non-sterile
	Manufacturing date		Protect from heat and direct sunlight

15. Storage

Must be stored at room temperature in a clean and dry location and protected from humidity and direct sunlight.

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16. Guarantee policy

Dimeda Instrumente GmbH is responsible for each individual product being p inspected and packaged with the utmost care. Since ***Dimeda*** has no influence or cor correct and professional application, ***Dimeda*** cannot be held responsible for complicatio failure of an application. ***Dimeda*** sets and individual products are mutually compatible. U responsible for verifying the mutual compatibility of the products before use.
Dimeda employees are not authorised to change the above conditions or extend liability into additional product-related commitments.

17. Bibliography

"Hygiene requirements for the reprocessing of medical devices", RKI (Robert Koch Instit

For additional information, please contact:

Dimeda Instrumente GmbH
Gänsäcker 54 + 58
78532 Tuttlingen
Germany
info@dimeda.de
Tel.: +49 7462 94 61 – 13
Fax: +49 7462 94 61 – 33
www.dimeda.de

Notified body:
CE 0123
TÜV-Süd Product Service GmbH
Ridlerstraße 65
80339 München, Germany
meineanfrage@tuev-sued.de
+49 89 50084747
www.tuev-sued.de

dimeda®
SURGICAL INSTRUMENTS
DIMEDA Instrumente GmbH
Gänsäcker 54 - 58
D-78532 Tuttlingen / Germany
i.A.

Last revised: April 2019
Instructions for use for orthopaedic implants (Revision 12)

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Instruction manual

For sterile containers with accessories CE

Valid from:	23.04.2018	Version:	1

dimeda®
SURGICAL INSTRUMENTS

Manufacturer:

Dimeda Instrumente GmbH
Gänsäcker 54 + 58
D-78532 Tuttlingen
Tel: +49 7462 94610
www.dimeda.de
info@dimeda.de

Products

These operating instructions are valid for all sterile container systems and accessories of Dimeda Instrumente GmbH with the article number range:

REF 88.XXX.XX

Important note

Read these instructions for use carefully before each use and keep them easily accessible for the user or the respective qualified personnel.

Read the warnings marked with this symbol carefully. Improper use of the products may result in serious injuries to the patient, users or third parties.

1 Scope of application

The products listed above may only be used by suitably trained and qualified personnel. The products may only be used in a sterile environment.

2 Precautions and Warnings

- Do not use steel wool or abrasive detergents.
- Do not use cleaning solutions with iodine or high chlorine content.
- Do not place contaminated or used Dimeda Instrumente medical devices in a case for cleaning in the washer-disinfector. Contaminated/used Dimeda Instrumente products need to be cleaned separately before putting them into the sterilizing container. Cases are designed as organization containers for steam sterilization, as storage containers for medical devices and as organization containers during surgery.
- Mechanical cleaning is preferable as it leads to a more effective result. Machine cleaning and disinfection provide greater safety in the process.
- Alkaline cleaning agents (pH>10) are not suitable for all materials. The Robert Koch Institute points out that potential problems can be caused such as increased wear of the following materials: aluminium, silicone elastomers, adhesive joints, soldered joints of silver and tin, sealing materials, plastic coatings, fibre-optic cables and optical surfaces with anti-reflection coating.
- Defective products must have undergone the entire reprocessing treatments before being returned for repair or complaint. Proof of decontamination must be enclosed with the return shipment. - The sterilization parameters only apply to adequately pre-cleaned components.
- The parameters performed only apply to properly installed, maintained and calibrated conditioning systems that meet the requirements of ISO 15883 and ISO 17665 standards.
- Patients, which are to be considered as "at high risk" regarding Creutzfeldt-Jakob Disease (CJD) and related infections, shall be operated with disposable instruments. Instruments which were used for operations on patients with suspected or evidenced CJD disease have to be disposed of and/or to be handled according to current valid national recommendations.

3 Limits of clinical preparation

- Repeated/frequent preparation in accordance with this instruction has only a minor effect on the service life of the containers.
- The service life of a sterilization container is essentially determined by wear and tear and damage caused by the application.
- In case of proper use 4 times a week on average, a service life of at least 10 years can be reached.

4 Application

- Dimeda Instrumente sterile container systems combine approved filter technologies, tested materials and design features to create a reliable container system. They are reusable container systems that offer a wide range of dimensions and equipment to provide effective packaging, storage and transport of the instruments to be sterilized. The container systems are perfectly suitable for fractionated vacuum processes.

Intended purpose

The Dimeda Instrumente sterile container systems are designed for the sterilization of medical instruments. The containers enable the sterilization and storage of the products inside until they are finally used. Depending on the model, the containers are available with perforated and non-perforated tub bottoms and perforated lids. The 1/1, ¾ and ½ containers are also available with a safety lid.

Combination products

Dimeda Instrumente sterile container systems consist of sterile containers, sieve baskets and filters. Accessories for the container systems are also available. A sieve basket of the appropriate size should be used for the respective container size. Possible combinations of the various container designs are described below. A detailed overview of combinable products can be found in section COMBINATION PRODUCTS

- **Standard container**
In the lid and, if applicable, in the tub, there are filter holders below/above the perforations. Disposable paper filters or permanent filters must be inserted into these filter holders before sterilization.
The lid of the standard container sizes 1/1, ½ and ¾ can additionally be covered with a safety lid if required. This safety lid enables additional protection against contamination during storage or transport of the sterile container.

- **Sieve baskets**
For each container size matching sieve baskets in different heights, corresponding lids and matching feet are available.

- **Security seals**
Security seals are attached to the outside of the closures by passing the seal through the opening of the spring closure system and then locking the seal. The seal breaks when the closure system is opened/raised.

- **Silicone mats**
The baskets sieve are to be put into the containers and can be equipped with a silicone mat if required.

- **Indicator labels**
The indicator inside the container changes colour during steam sterilization at 134°C. Please observe the shelf life of the labels according to the manufacturer's instructions. The indicator labels may only be used for the intended purpose. Non-compliance with the specifications may falsify the result.

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Notes on the use of paper filters

- Paper filters are intended for single use only.
- Paper filters are manufactured according to EN ISO 11307-1.
- Paper filters must not be labelled with adhesive material (e.g. for documentation of the runs), as the adhesive material may contain harmful substances. In addition, the germ barrier is destroyed by glueing.

The paper filters used must have a size that makes sure that the perforation in the container lid is completely covered.

Notes on the use of permanent filters

- PTFE filters are designed for multiple use.
- Permanent filters must not be labelled with adhesive material (e.g. for documentation of the runs), as the adhesive may contain harmful substances. In addition, the germ barrier is destroyed by glueing.
- In case of coarse soiling, the filter must be removed and then cleaned.

The permanent filter must have a size that makes sure that the perforation in the container lid is completely covered.

5 Handling and preparation

General information

The Dimeda Instrumente sterilization containers consist of an aluminium alloy whose surface is anodized for corrosion protection. Aggressive cleaning agents, metal brushes or scouring cloths can permanently damage this surface and must therefore not be used. If these instructions are not followed, warranty shall be excluded.

The sterilization containers may only be handled by instructed or trained personnel in order to prevent damage to the containers, closures, seals and sterile filters/cassettes.

The sterilisation containers are also available with coloured lids to facilitate allocation to the individual disciplines and specialist departments. Sterilisation indicators and coloured identification labels provide information on content, location and condition.

In accordance with the normative specifications and recommendations, suitable measures (e.g. sealing, process indicators) must be taken to ensure that sterilised and non-sterilised sterilisation containers cannot be confused. Only intact seals ensure that the sterilization container has not been opened without permission.

Preparation for cleaning

- Separation of container tray and lid
- Remove the contents of the container (sieve basket, instruments, etc.).
- Remove the filter holder/cassette from the inside of the lid and, if applicable, from the tub section (for containers with bottom perforation).
- For disposable paper filters: Dispose of disposable filter
- Removal of disposable seals and indicator labels

Commissioning of a new container

- The container must be thoroughly cleaned before it is used for the first time.
- The container shall be processed in a validated mechanical cleaning and disinfection process.
- For this purpose, a neutral cleaning agent should be used in the machine.
- After reprocessing in a mechanical cleaning and disinfection process, the products must be steam sterilised in a fractionated steam sterilisation process at 134°C.
- In addition, all moving parts on the container must be regularly maintained with an approved instrument care oil.
- After cleaning, suitable new filters must be used (see filter change).

6 Preparation

Additional Information

- The urban water that is used must comply with Council Directive 98/83/EC of 3 November 1998 on the quality of water intended for human consumption.
- The reprocessing instructions shall specify the detergents and disinfectants used for validation. If an alternative detergent and disinfectant (RKI or VAH listed) is used, the responsibility lies with the person in charge.
- Reassemble dismantled products before sterilisation.

Manual pre-cleaning

- For aluminium containers and lids, mild, neutral cleaning agents shall be used or alternatively chemical products which have been expressly approved by manufacturers for the treatment of aluminium products. If necessary, the products must be tested for suitability by using the appropriate procedure. A soft, suitable sponge should be used for manual cleaning. Scrubbing sponges should not be used as they destroy the surfaces and thus the passivation and will lead to the loss of warranty.
- After cleaning, careful rinsing with low-salt water (e.g. deionised water) and sufficient drying is required.
- Do not use metal brushes or abrasive cleaners.
- Finally, disinfection must be carried out in accordance with the relevant hygiene requirements.
- Ultrasound
 If the pre-cleaning by sponge and rinsing with a water jet gun has left behind any visually visible impurities, a pre-cleaning by ultrasound must be carried out.
- Dip container and sieves into ultrasonic bath filled with water <40°C, 0.5% alkaline detergent (Neodisher FA, Dr. Weigert) for 5 minutes and clean for 10 minutes. Rinse containers and sieves with a water jet gun (4 bar) > 10 sec.

Cleaning / Disinfection

Automatic cleaning/disinfection process:

(Washing machine: Washer - Disinfector G 7735 CD (Miele):

Step 1: 1 minute Pre-cleaning with cold tap water Drinking water quality <40°C

Step 2: Water drain

Step 3: 3 minutes pre-cleaning with cold tap water Drinking water quality <40°C

Step 4: Water drain

Step 5: 5 minutes of cleaning at 55°C ± 5°C with 0.5% alkaline cleaner (Neodisher FA, Dr. Weigert).

Step 6: Water drain

Step 7: 3 minutes neutralisation with cold tap water Drinking water quality <40°C

Step 8: Water drain

Step 9: 2 minutes rinse with cold tap water Drinking water quality <40°C

The special instructions of the automatic cleaning machine manufacturer must be observed.

Automatic disinfection

Automatic thermal disinfection in washer-disinfectors, taking into account the national requirements for the A0 value 3000:

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Automatic drying

Automatic drying according to the automatic drying process of the washer-disinfector for at least 30 minutes (at 60°C ± 5°C in the rinsing room). If necessary, subsequent manual drying with lint-free cloth and blowing out the corners with sterile, oil-free compressed air.

Filter change

After changing the filter, the filter holder must be pressed into the correct position until it snaps audibly into place. Dimeda Instrumente lids may only be used with Dimeda Instrumente filter holders.

- Disposable paper sterile filters must be replaced before each new sterilization.
- The suitability and accuracy of the filters can only be guaranteed if the Dimeda Instrumente filters are used.
- Warranty services can only be provided if the original Dimeda Instrumente filters are always used.
- PTFE filters have been tested for a service life of 1200 cycles and must then be replaced.

Attention

Combine only original Dimeda Instrumente component such as lids, tubs, filters, seals, cassettes and filter holders with each other, so that the tightness and the germ barrier are not put at risk. Otherwise, Dimeda Instruments does not assume any warranty.

Sterilization

Sterilization of container systems using fractionated pre-vacuum processes (according to DIN EN ISO 17665-14), taking into account the respective national requirements. Products must be sterilised in suitable sterilisation packaging in accordance with EN ISO 11607 and DIN EN 868.

EU standard

The sterilization shall be performed with a fractionated pre-vacuum method with the following parameters:
134°C
≥ 5 minutes holding time,
3 pre-vacuum cycles
Vacuum drying for at least 20 minutes.
The instructions for use provided by the manufacturer and recommended guidelines for maximum loading with sterilization material must be observed. The autoclave must be properly installed, maintained, validated and calibrated.

Container loading

The total weight of the container load should not exceed the following quantities, otherwise satisfactory sterilization cannot be guaranteed.

Model	Max. Loading in kg
1/1 (Full-) Size Container	9,0 kg
¾ Size-Container	7,0 kg
½ Size-Container	5,0 kg
Flat -Container	1,5 kg
Mini-Container	1,0 kg
Dental-Container	1,8 kg

Placement in sterilizer

The containers are designed in such a way that they can be used in any commercial large sterilizer for sterilization with moist heat. Please note that heavy containers are to be positioned at the bottom of the sterilization chamber. Due to their design, the containers can be stacked on top of each other safely and easily during sterilization without slipping. Stacking is only recommended for sterilization cycles using fractionated vacuum processes. The stacking height should not exceed 46 cm to achieve effective air removal and steam penetration. Follow the instructions of the sterilizer manufacturer.

CAUTION!

Note the following when sterilizing: Never pack the container in another outer packaging. Never cover the perforation fields in the lid and base with any foil packaging or the like, as this would obstruct the air and steam flow into the container. This would result in a vacuum-induced deformation of the container due to insufficient pressure equalization and the sterility of the contents of the container is not guaranteed. When loading and unloading the steriliser and for any transports please always use the carrying handles and never the lid.

Sequence control

- Operate the loaded sterilizer according to the sterilizer manufacturer's specifications for the selected sterilization cycle (with reference to temperature and sterilization time) taking into account the validation results.
- To avoid condensate accumulating in the container, the container should cool down completely on the sterilization trolley.
- After each sterilisation, the sterile material must be assessed and released in accordance with the internal instructions and the validation results. This is consistently done by qualified employees with.

Additional information

It is the responsibility of the person in charge for conditioning to ensure that the actual treatment carried out with the equipment, materials and personnel, used in the treatment facility, achieves the desired results. This usually requires validation and routine monitoring of the process and the equipment used.

7 Examination

The sterilization containers must be checked for proper functioning before each use. Damage to the closures, seals, filter holders, filters, cassettes and bent and dented parts mean that the sterilization containers must be repaired and must not be used. Do not use damaged sterilization containers.

- All moving parts on the container must be treated with an approved instrument care oil.
- If damage to the seals is found, they must be replaced immediately.
- The seals should not be treated with spray, oil or solvent. For cleaning and care, an occasional wipe with a damp cloth is sufficient.
- If damage is found to the sterilization containers, they must be inspected, repaired or replaced.
- Spare parts can be obtained from Dimeda Instrumente

Service and repair

Do not carry out any repairs or modifications to the product yourself. Only authorized personnel of the manufacturer is responsible and intended for this. Should you have any complaints or remarks regarding our products, please contact us.

Return

Defective or non-compliant products must have undergone the entire reprocessing treatments before being returned for repair / service.

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8 Storage, transport and disposal

Storage

Please refer to DIN 58953-9 (Application Technology of Sterilization Containers) for the storage period of medical devices in sterilization containers. The storage time usually depends on the storage conditions and must be determined by the responsible hygiene personnel. In case of a particularly high requirement for asepsis or in case of deviations from the specified storage conditions, shorter storage periods or additional packaging should be used..

Storage conditions:

- Temperature: 15 – 26°C
- Air humidity. 30 – 50%
- Air pressure: normal atmospheric pressure

Container loading in different ways, storage times and storage conditions are to be determined by the responsible hygiene personnel. The Dimeda Instrumente sterile containers were tested for a storage period of 6 months by applying Bacillus subtilis spore suspension. Therefore, a storage period of 6 months can be guaranteed. The containers must be stored under protected conditions (e.g. in closed cabinets) dust-protected, clean, dry and free from vermin.

Transport

The sterile containers should only be transported using the carrying handles provided for this purpose.

Disposal

The products, packaging material and accessories must be disposed of in accordance with national regulations and laws. The manufacturer does not give any specific instructions for this.

9 Liability and Warranty

Dimeda Instrumente GmbH, as manufacturer, is not liable for consequential damages resulting from improper use or handling. This also applies to repairs or modifications to the product carried out by unauthorized personnel of the manufacturer. These exclusions of liability also apply to warranty services.

10 Symbol description

Attention!

Follow the instructions for use

REF Item number

LOT Batch indication

CE CE mark

Non sterile

Manufacturer

Store in a dry place

Temperature limitation

Protect from sunlight

11 Applied standards

In order to ensure the safety of the sterilization containers during production and handling, the following standards have been observed:

DIN EN 868-2
Packaging materials and systems for medical devices to be sterilized - Part 8:
Reusable sterilisation containers for steam sterilisers to EN 285; requirements and test methods

EN ISO 11607-1
Packaging for medical devices to be sterilized in the final packaging - Part 1:
Requirements for materials, sterile barrier systems and packaging systems

DIN 58952-2
Sterilization; packaging materials for sterile goods, metal sterilization trays

DIN 58952-3
Sterilization; packaging material for sterile goods,
Sterilizing sieve trays made of metal

DIN 58953-9
Sterilization; Sterile Supply - Part 9:
Application technology for sterilization containers

EN ISO 14937
Sterilisation of health care products - General requirements for characterisation of a sterilising agent and for development, validation and routine control of a sterilisation process for medical devices

EN ISO 17665-1
Sterilization of health care products - Moist heat - Part 1:
Requirements for the development, validation and control of the application of a sterilisation process for medical devices.

DIN EN 285
Sterilization - Steam Sterilizers - Large Sterilizers

To ensure sterility safety, tests were carried out by an independent and accredited testing laboratory. The purpose of these tests was to validate a moist heat sterilization process for the reusable Dimeda Instrumente sterile container systems.

Based on the results, we therefore stipulate the sterilization procedures indicated on page 3 of this instruction manual.

12 Materials

Our sterilization containers are made of an anodized aluminium alloy and the accessories are made of stainless instrument steel.

13 Delivery condition

NON STERILE The Dimeda Instrumente sterile container systems are delivered in a non-sterile condition and must be prepared and sterilized by the user according to the following instructions before the first and each further application.

DIMEDA INSTRUMENTE GMBH ACCEPTS NO LIABILITY IF IT CAN BE PROVEN THAT THIS CUSTOMER INFORMATION HAS BEEN VIOLATED.

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