dimeda®
SURGICAL INSTRUMENTS

To whom it may concern

STATEMENT

Herewith, «DIMEDA Instrumente GmbH», having its principal business office at Gänsäcker 54+58, 78532 Tuttllingen, Germany, represented by its General Manager, Mr. Bastian Milkau certifies the following information:
«Instruction Manual for reusable surgical instruments» covers the following range of our products:
scalpel knives, knife handles, scalpel blades, knives, scissors, punches, dissectors, nasal files, drilling instruments, raspatories, skinny raspatories, extractors with accessories, rib-cutting forceps, bone forceps, surgical snares, surgical spoons, curettes, bone files, bone saws with accessories, drills, bone pliers, bone chisels, bone punches, tractors,bone forceps, screwdrivers with accessories, plaster knives, plaster dilator, hand plaster saws, screw taps, sleeves with accessories, plates and plate sets with accessories, needles, suture needles, needle holders, awl, staples, clamping applicators, ligature instrument with accessories, pulmonary spatulas, spatulas, dilators with accessories, spreaders, mouth dilators with accessories, operating retractors with accessories, speculum with accessories, hooks with accessories, depressors, elevators with accessories, scalers, artery clamps with accessories, bulldog clamps, clamps for umbilical cord, clamps with accessories, surgical tissue clamps, tweezers, biopsy forceps, forceps others, ear mirror with accessories, goniometers, dynamometers, cephalometer, esthesiometer, tonometer with accessories, caliper, percussion hummers with accessories, camertone with accessories, pin wheel, caliper gauge, stethoscopes with accessories, head lamp, head light, probes with accessories, forehead mirrors with accessories, skin markers, anchors with accessories, hummers with accessories, microsurgical instruments with accessories, obstetrical forceps, urethrotomes with accessories, uterine dilators with accessories, vaginal speculum with accessories, metallic catheters, surgical clamps for penis, additional accessories for gynecology and urology, dental forceps, chisel, tweezers, curettes, periodontometer, excavators, amalgam plugger, polish instrument, filling instrument, additional instruments for filling, separator, mould forceps, matrix, mouth spoon, modeling instrument, modeling carver, hummers, wire dispenser, additional dental accessories,

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«Instructions for Use and Sterilization Instructions for Orthopaedic Implants», «Instruction Manual for reusable surgical instruments» cover cutting instrument for needles, wires with accessories.

«Instruction Manual for reusable surgical instruments», «Instruction for Use Suction – Irrigation Instruments» cover tubes for aspiration/ irrigation with accessories.

«Instruction Manual for Sterile Containers with Accessories» cover sterilization container with accessories.

«Instruction Manual for Sterile Containers with Accessories», «Instruction Manual for reusable surgical instruments» spread for medical containers with accessories.

«Instruction Manual for reusable surgical instruments» Instructions for Use and Sterilization Instructions for Orthopaedic Implants» cover additional accessories for surgical instruments, orthopedic and for osteosynthesis.

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Instruction manual

For reusable surgical instruments CE

Valid from:	10.12.2020	Version:	5

dimeda® SURGICAL INSTRUMENTS

Manufacturer

Dimeda Instrumente GmbH
Gänsäcker 54+58
78532 Tuttlingen, Germany
Phone:+49 (0) 7462 / 9461-0
Fax: +49 (0) 7462 / 9461-33
http://www.dimeda.de
info@dimeda.de

Products

This instruction manual is valid for all reusable surgical instruments, which

- are one-piece instruments
- optionally have simple joints or
- include simple moving parts
- are optionally composed of multiple replaceable parts (e.g. handle and multiple attachments)

Important Information

Read this instruction manual carefully before every application and keep it easily accessible for all users or the respective specialist staff.

Carefully read the warnings marked with this symbol. Improper use of the products may result in serious injuries to the patient, the users or third parties.

1 Scope

The instruments must be used according to their intended use in the medical fields and by respectively trained and qualified staff only. The treating physician and/or user is responsible for choosing the equipment for specific applications and/or operative use, for the appropriate training and information, and for the sufficient experience regarding the handling of the equipment.

The instruments intended application duration is 60 minutes. The product must be replaced after 60 minutes.

2 Contraindication

The instruments must be used according to their intended use in the medical fields and by respectively trained and qualified staff only. The instruments of Dimeda GmbH may not be used in the areas of cardiothoracic surgery and neurosurgery.

3 Precautions and Warnings

Attention!

The instruments are only designed for surgical us
other purpo

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maintenance as well as misappropriated use may result in the premature wear and tear of the instruments.

Material intolerance

Under no circumstances must the instruments be used if the user or specialist staff become aware of the patient being intolerant to the material.

Functional Impairment

Surgical instruments corrode and become impaired in their functionality if they come into contact with aggressive substances. It is therefore necessary to observe the storage and sterilisation instructions.

Operating Conditions

The aforementioned products require correct maintenance and care in order to guarantee that the products operate safely. In addition to this, functionality testing and a visual check should be performed prior to each application. For this reason, please pay attention to the respective chapters in this instruction manual.

Combination with other products

Should the products be reassembled after disassembly, individual parts must not be replaced with parts from other manufacturers! If the intended purpose of the product entails certain parts being exchanged (e.g. different attachments), no parts from different manufacturers must be used! We recommend to also purchase other accessories (e.g. detergents) at Dimeda Instrumente GmbH.

Storage

There are no specific storage requirements concerning the products. Nevertheless, we recommend storing medical products in a clean and dry environment.

Creutzfeldt Jakob Disease

With regard to the reprocessing of medical devices that have been used on patients or suspected patients suffering from or suspected of suffering from Creutzfeldt-Jacob disease (CJD) or its variant (vCJD), the requirements specified in the corresponding appendix of the guidelines for hospital hygiene and infection prevention and the requirements specified by publications in the Federal Health Gazette must be adhered to. The medical devices that were used on this group of patients must be disposed of by incineration (European Waste Catalogue EAK 18 01 03) without risk. Dry heat, ethanol, formaldehyde and glutaraldehyde have a fixing but no inactivating effect on TSE pathogens. Of the sterilization methods available, only steam sterilization (especially 134°C, 18 minutes) has been shown to have a limited effect.

4 Liability and Warranty

As a manufacturer, Dimeda Instrumente GmbH is not liable for consequential damage resulting from improper use or handling. This particularly applies to use which is not compliant with the defined intended use, or non-compliance with the instructions on preparation and sterilisation. This also applies to repairs or changes to the product which are not carried out by authorized staff of the manufacturer. These disclaimers also apply to warranty services.

5 Sterility

State upon Delivery

Medical products are delivered in a non-sterile condition and need to be prepared and sterilised by the user prior to the first application and any subsequent application according to the following instructions.

6 Reprocessing

Warnings

- Frequent reprocessing impairs the quality of the products.
- Any tap water to be used must comply with COUNCIL DIRECTIVE 98/83/EC of 3 November 1998 on the quality of water intended for human consumption.
- The detergents and disinfectants used in order to maintain validation are listed in this instruction for reprocessing. Should an alternative detergent or disinfectant (listed at the RKI or VAH) be used, the responsibility lies with the reprocessor.
- Reassemble disassembled products prior to sterilisation.
- The preparation may only be carried out by qualified medical personnel. Machine reprocessing must be qualified and validated by the user. The washer-disinfectors must fully comply with the requirements of DIN EN ISO 15883.

Use Site

The first steps of a proper reprocessing take place in the operating theatre. Coarse contaminations must be removed prior to storing the instruments if possible. For this purpose, the instruments should be rinsed under cold tap water (<40°C). If this procedure is not sufficient to remove the obvious soiling, a soft plastic brush can be used to remove soiling.
Whenever possible, dry removal (moistened, closed system) should be the method of choice. A drying of any residues should be avoided! Wherever possible, dry disposal is to be preferred, since with wet disposal the prolonged lying of the medical devices in solutions can lead to material damage (e.g. corrosion). Long periods of waiting until the reprocessing, for instance overnight or over the weekend, must be avoided with both types of removal (<60 minutes).

Transport

The products must be disposed of in a dry state immediately (<60 min) after use, if possible. This means that the products have to be transported in a closed container from

the place of application to the purification, so that the products do not dry up.

Preparing the Decontamination

The products must be disassembled prior to the following reprocessing steps and/or must be exposed to the following reprocessing steps in an open condition, where possible. Rinse residue must be avoided. The products must be reprocessed in appropriate screen baskets or rinsing shields (choose size according to product). The products must be positioned in the cleaning basket at a minimum clearance from one another. Avoid overlapping so that the damaging of the products during the cleaning process can be excluded.

Pre-cleaning

Rinse the products under cold tap water suitable for use as drinking water (<40°C) until all visible contaminations are removed. Remove persistent dirt with a soft brush. Movable parts on the instrument need to be moved. Cavities, lumen, gaps and slits need to be rinsed intensively (>60 sec) with cold tap water suitable for use as drinking water (<40°C) using a water pistol (or similar). Insertion of the products into an ultrasonic bath (<40°C) with an alkaline cleaner (0.5% Neodisher® MediClean forte, Dr. Weigert), time of treatment 5 min. at a frequency of about 35 kHz. For this please observe the instructions of the cleaning agent manufacturer. Rinse instruments briefly (<15 sec) under cold water. Movable parts must be moved. Respectively rinse cavities, lumen, gaps and slits again (>30 sec.) with cold city water (<40°C) using a water pistol (or similar).

Cleaning/disinfection

Automated cleaning and/or disinfection process

(Miele Disinfector G7835 CD as per ISO 15883):

- Pre-clean for 1 minute
- Drain water
- Pre-clean for 4 minutes
- Drain water
- Clean for 6 minutes at 58°C +/- 1°C using 0.5 % alkaline detergent (0,5 % Neodisher Mediclean forte)
- Drain water
- Neutralize for 3 minutes (0.1% 0,1 % NeodisherZ) with cold tap water suitable for use as drinking water <40°C
- Drain water
- Clean for 2 minutes with FD water <40°C

Please observe the special instructions of the manufacturer of the automatic cleaning machine.

Automated Disinfection

Automated thermal disinfection in a cleaning and disinfe[illegible] considerat[illegible] the A_0 valu[illegible]

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Automated Drying

Automated drying in accordance with the drying operation of the cleaning and disinfection device for at least 30 minutes at 92°C +/- 2°C.

Sterilisation

(Typ B Autoklav by Tuttmauer as per DIN EN 13060)

Sterilisation of products with a fractionated pre-vacuum method (according to DIN EN ISO 17665-1) taking into consideration the respective national requirements. The sterilisation of the products must be conducted in suitable sterilisation packaging according to DIN EN ISO 11607-1 and EN 868

The sterilisation must be completed using a fractionated pre-vacuum method with the following parameters:

134°C,

5 minutes hold time,

3 pre-vacuum cycles

Drying in vacuum for at least 20 minutes

The instruction manual of the manufacturer of the autoclave and the recommended directions for maximum loading with goods to be sterilised must be observed. The autoclave must be installed, maintained, validated and calibrated in accordance with requirements.

⚠ **Additional Information**

The reprocessor is responsible for ensuring that the actual reprocessing, including the used equipment, materials and the staff involved in the reprocessing facility, achieves the desired results. This typically requires the validation and routine monitoring of the method and the equipment used.

7 Lifespan of the Products

Frequent reprocessing has little effect on the surgical instruments. Product life is limited by wear due to proper use and damage to the instrument. The product may no longer be used under the following conditions, among others:

Corrosion, damage, fractures, cracks, deformation, porosity, functional limitations, products with unrecognizable or missing marking. For this reason, the conditioner must comply with the relevant instructions for functional testing (section 8).

8 Functional testing

Check products after reprocessing and before sterilization with regard to the following aspects:

- Cleanliness
- Damage, including but not limited to signs of corrosion (rust, pitting), discoloration, deep scratches, flaking, cracking and wear.
- Proper function
- Missing or removed (ground off) part numbers.
- Do not use improperly functioning or defective and excessively worn products, products with unrecognizable markings, m or removed (ground) part numbers.
- Check products for perfect surfaces, correct assembly a functionality. Do not use sev damaged products.
- Reassemble dismantled proc before sterilization.

9 Service and Repair

⚠ **Service and Repair**

Do not carry out any repairs or changes the product yourself. Authorized staff of manufacturer are solely responsible for work. Should you wish to make compla queries, or offer us any advice regardin products, please feel free to contact us.

⚠ **Returns**

Defective or non-compliant products mu through the entire reprocessing process before being sent back for repairs/servic

10 Packaging, Storage and Disposal

Standard packaging of the products for sterilisation according to ISO 11607 and 868.

Store sterile products in a dry, clean, an dust-free environment, secured against damage, at moderate temperatures.

The medical products of the manufactur should be stored and kept in single packaging, boxes or protective containe Please handle the instruments with care during transportation, storage and reprocessing. The user and/or specialist intended for this is responsible for ensur that the sterile state is maintained after t sterilisation.

The disposal of the products, packaging well as the accessories must be perform accordance with current rules and laws. specific instruction regarding this matter provided by the manufacturer.

11 Description of Symbols Us

⚠	Attention!
	Observe the instruction manual
REF	Item number
LOT	Lot designation
CE	CE labelling
	Indication of a non-steril product
	Name and address of the manufacturer

dimeda®
SURGICAL INSTRUMENTS
DIMEDA Instruments GmbH
Gänsäcker 61 [illegible]
D-78532 Tuttlingen / Germany

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To whom it may concern

STATEMENT

Herewith, «DIMEDA Instrumente GmbH», having its principal business office at Gänsäcker 54+58, 78532 Tuttllingen, Germany, represented by its General Manager, Mr. Bastian Milkau certifies that our Instruction for use «DIMEDA Cables for Electrosurgery», «Bipolar and Monopolar Electrodes» covers the following products of Class I: electrosurgery instruments with accessories.

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ENGLISH Bipolar and Monopolar Electrodes

⚠ Please read all information contained in this insert attentively.
Incorrect handling and care, as well as misuse, can lead to premature wear of surgical instruments or risks to patients and users.

Intended Use
Bipolar electrodes for open surgery
The bipolar electrodes (REF 90.502.17, 90.502.18) serve cutting and coagulation of biological tissue. The fully assembled instrument (if assembly is needed) has to be – with the appropriate cable - to monopolar or bipolar output of an HF generator.
Only the defined parameters has to be used.

Maximum output voltage of the generator, U_{max}:
500 Vp , 250 Vp (depending on the type, see catalogue)

Appropriate connecting cables for bipolar electrodes for open surgery:
Dimeda bipolar cable REF 90.190.30 – 90.190.50; 90.193.30 – 90.193.50; 90.415.30 – 90.418.50; 90.425.30 – 90.428.50; 92.706.01 – 92.810.06

Bipolar and monopolar electrodes for MIS
The bipolar electrodes for MIS (REF 92.305.03 – 92.305.06 und 92.303.01 - 92.303.03) are intended to be used for coagulation of biological tissue.
The electrodes are intended for use in minimally invasive and especially laparoscopic surgical procedures. The instrument has to be inserted through a trocar sleeve with the appropriate diameter.
The fully assembled instrument (if assembly is needed) has to be connected- with the appropriate cable - to monopolar or bipolar output of an HF generator. When indicated, monopolar or accordingly bipolar coagulation or cutting current can be selectively applied.

Maximum output voltage of the generator, U_{max}:
for bipolar electrodes for MIS:
500 Vp , 300 Vp , 250 Vp (depending on the type, see catalogue)

for monopolar electrodes for MIS:
2kV

Suitable cables for bipolar electrodes for MIS:
Dimeda bipolar cable REF 90.190.30 – 90.190.50; 90.193.30 – 90.193.50; 90.415.30 – 90.418.50; 90.425.30 – 90.428.50; 92.706.01 – 92.810.06.

Suitable cables for monopolar electrodes for MIS:
Dimeda Monopolar Cable REF 90.190.30 – 90.190.50; 90.193.30 – 90.193.50; 90.415.30 – 90.418.50; 90.425.30 – 90.428.50; 92.706.01 – 92.810.06.

⚠ *Instruments for electrosurgery must only be used by persons who have been specially trained or instructed in this.*

Contraindications
- Do not use the instrument if, in the opinion of the attending physician, the risks to the patient outweigh the benefits.

Incidents that have been reported in connection with the use of electrosurgical systems
- Unintended activation with resulting tissue injury in the wrong location and/or damage to the equipment.
- Fire in connection with surgical drapes and other inflammable materials.
- Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
- Explosions caused by sparks in the proximity of inflammable gases.
- Perforation of organs. Sudden severe bleedings.

Use and safety instructions
Non-observance of these use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
- When using electrosurgery in patients with pacemakers or other active implants, special requirements apply (e.g. low HF-current, patient monitoring). In any case, a cardiologist or appropriate medical specialist must be consulted.
- Before initial use and any further use, all instruments must be completely cleaned, disinfected and sterilised and their function must be checked.
- It is very important to check every surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas such as blades, tips, notches, locking and blocking devices, as well as all movable parts, insulations and ceramic elements must be checked carefully.
- Never use any damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- When temporarily not in use, the instrument must be placed electrically insulated from the patient.
- Activate electrosurgical current only if the contact areas are in full view and have good contact with the tissue that needs to be treated. Do not touch any other metallic instruments, trocar sleeves, optics or similar objects during use.
- Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

⚠ Applies for monopolar mode of operation:
Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.

Assembly
For assembly and disassembly of the instrument follow the pictogram, which is available upon request, or can be downloaded on www.dimeda.de.

Special instructions for the hysterectomy loop electrode

Insert the instrument through the cannula with the loop pulled inside the shaft.
Thoroughly check correct placement of the loop around the uterus prior to activation of the electrical current; the entire tissue surrounded by the loop will be cut.
After correct placement, maintain tension of the loop to avoid uncontrolled slip of the loop. Maintain some safety distance to the ligature.
During electrosurgical cutting, pull the loop continuously inside the shaft. Cutting will start at the uninsulated part of the loop wire.
The loop wire is a replaceable, sterilisable product for single use that should be replaced following each procedure.

eprocessing
ue to the product design, the materials used and the tended purpose, it is not possible to define a limit with egard to the maximum possible number of reprocessing ycles. The serviceable life of the instruments is etermined by their function as well as by a careful andling.
struments for electrosurgery are by their nature subject increased wear depending on the type and time of use.

reparation and transport
mmediately after each use, clean the instruments with a oft brush under cold tap water until all visible ontamination is removed. Do not use fixation agents or ot water (>40°C). Storage and transport of the struments to the reprocessing location must take place a sealed container.
omplex instruments must be taken apart for cleaning nd disinfection in accordance with pictogram.

lachine reprocessing
leaning
lace the instruments in a basket on the insert module or n the inserts of the MIS module and start the cleaning rocess.
Prerinse with cold water for 1 min
Discharge
Prerinse with cold water for 3 min.
Discharge
Wash at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent for 5 min.
Discharge
Neutralise with warm tap water (>40°C) and a neutralising agent for 3 min.
Discharge
Rinse with warm tap water (>40°C) for 2 min.
).Discharge

isinfection
achine-operated thermal disinfection must be carried ut under observation of the national requirements garding the A0 value (see ISO 15883).

rying
ry the outside of the instruments by carrying out a drying cle of the cleaning/disinfection machine.
necessary, manual drying may additionally be carried ut using a lint-free cloth. Dry cavities by blowing with erile compressed air.

anual reprocessing
ltrasonic pre-cleaning
The instruments are placed in an ultrasonic bath with 0.5% enzymatic cleaning detergent and treated with ultrasound for 15 minutes at 40°C/104°F.
Remove the instrument and rinse them completely with cold water to remove the cleaning detergent.

eaning
repare a cleaning bath according to the manufacturer's structions.
Rinse products with cold tap water (<40°C) until all visible contamination has been removed. Remove adhering dirt by using a soft brush.
Place products in the prepared cleaning bath so that they are completely submersed. Observe residence time according to the manufacturer's instructions.
Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
The following step only applies to channels and the insides of tubes: Push the brush into and out of the tubes at least six times. Rinse the tubes with DI water. Repeat the procedure.
Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Drying
Manual drying is carried out using a lint-free cloth and sterile compressed air, in particular for drying cavities and channels.

Functional test and packaging
Perform visual inspection for cleanliness and integrity; if required, perform an assembly and functional test.
If necessary, repeat reprocessing until the instrument is visually clean. Packaging must comply with the ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation
Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) under observation of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of at least 132°C and at most 137°C
- Exposure time: at least 3 min.
- Drying time: at least 10 min.

⚠ If contamination with prions (CJD) is suspected, differing national guidelines are to be followed and longer holding times (i.e. 15 min.) may apply.

Storage
Sterilised instruments must be stored in a dry, clean and dust-free environment. The applicable national guidelines must be followed.

Repairs
Never attempt to perform repairs yourself. Service and repair work must only be performed by persons trained and qualified accordingly. If you have any question regarding these matters, contact either the manufacturer or your medico-technical department.
⚠ *Defective products must complete the entire reprocessing process before being returned for repair.*

Information on the validation of the reconditioning
The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):
Neodisher FA by Dr. Weigert (alkaline)
Endozime by Ruhof (enzymatic)
Cleaning agents (manual cleaning):
Cidezyme, Enzol Enzym detergent, Johnson&Johnson
Disinfectants (manual disinfection):
Cidex OPA, Johnson&Johnson
Neutralising agent:
Neodisher Z by Dr. Weigert
Cleaning and disinfection device:
Miele Desinfector G 7735 CD
Miele insert module E 327-06
Miele MIS module E 450

For details, see report.
SMP GmbH # 01707011901 (machine cleaning)
MDS GmbH # 135196-10 (manual cleaning, sterilisation)
Nelson Labs # 200432706-02 (sterilisation)
MDS GmbH Testbericht 084183-10 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling
During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.
This applies particularly to blades, fine tips and other sensitive areas.

Disposal
Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

Warranty
Dimeda Instrumente GmbH exclusively supplies tested and faultless products to its customers.
All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which have been r... compared to the original product, misused or ... used improperly.

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Attention

Refer to instructions for use

CE-Mark and registration number of the Notified Body
TÜV Süd
Product Service GmbH
Ridler-Straße 65

Manufacturer
Production date

Attention:According to US-laws, this device must only be sold by a doctor or on the instruction of a doctor.

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REF

90.190.xx
90.192.xx
90.193.xx
90.290.xx
90.291.xx
90.415.xx
90.416.xx
90.417.xx
90.418.xx
90.425.xx
90.426.xx
90.427.xx
90.428.xx
92.706.xx
92.712.xx
92.810.xx

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ENGLISH **Dimeda Cables for electrosurgery**

Caution

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Indications for Use

Cables for electrosurgery are designed either to conduct electrical power from the output of a high-frequency generator to the instrument or to connect a neutral electrode with the generator.

Attention: *Instruments and cables for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.*

Use and safety instructions

- All cables have to be completely cleaned, disinfected and sterilised before initial use and any other use.
- The cables have to be submitted to a visual inspection and a functional test before each use.
- Never use damaged cables.
- Store the cables in a clean, cool and dry place. Protect them from damage.
- Handle the cables with care, and do not pinch them.
- Make sure that the cable is correctly inserted in the socket or connected with the plug.

Reprocessing

Preparation and transport

Remove coarse dirt from the cables immediately after each use. Do not use fixation agents or hot water (>40°C).

Machine reprocessing

Cleaning

Place the cables in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Prerinse for 1 min. with cold water
2. Discharging
3. Prerinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C)
10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the cables by carrying out a drying cycle of the cleaning/disinfection machine.
If necessary, manual drying may additionally be carried out using a lint free cloth. Dry cavities by blowing with sterile, compressed air.

Manual reprocessing

Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submersed. Observe residence time according to the manufacturer's instructions.
3. Clean the cable in the bath manually using a soft brush. Brush all surfaces several times.
4. Rinse the products thoroughly with running tap water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the cables in the disinfectant bath and observe the specified residence time. Rinse the cables thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Dry manually by using a lint free cloth and, in particular, for drying cavities, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness.
If necessary, repeat the reprocessing process until the instrument is optically clean.
Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 11140 / ISO 17665) in consideration of the respective national requirements.

- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):
Neodisher FA by Dr. Weigert (alkaline)
Endozime by Ruhof (enzymatic)
Cleaning agents (manual cleaning):
Enzol Enzym, detergent by Johnson&Johnson
Neutralising agent:
Neodisher Z by Dr. Weigert
Cleaning and disinfection device:
Miele G 7736 CD
Miele insert module E 327-06
Miele MIS module E 450

For details, see report.
SMP GmbH # 01707011901-2 (machine cleaning)
Northview Laboratories #P8H066 (manual cleaning, sterilisation)
Nelson Labs # 200432706-02 (sterilisation)
MDS GmbH Testbericht 084183-10 (Sterilisation)

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Handling

During transport, cleaning, care, sterilisation and storage, all cables should be handled with maximum care.

Warranty

Dimeda Instruments GmbH exclusively supplies tested and faultless products to the customers.
All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexpert way.

Explanation of symbols

LOT — Batch code

REF — Reference Number

Attention!

Refer to Instructions for Use

Manufacturer
Production date Year-Month

LATEX — Free of natural rubber latex

NON STERILE — Unsterile

℞ ONLY — Attention: According to US-laws, this device must only be sold by a doctor or on the instruction of a doctor.

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