MEDRAD® Stellant

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For use with injector software versions 105.6 or higher and 901.4 or higher and with Certegra® Workstation software version 101.05 or higher

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©2017-2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer

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Bayer, the Bayer Cross, MEDRAD, Stellant, MEDRAD Stellant, Certegra, MEDRAD FluiDots, and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

EC REP **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

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ENGLISH

Important Safety Notice
This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Please follow your institution's policy for media administration and aseptic technique.

NOTE: If cleaning of the multi-patient kit is desired, use only ethanol-based disinfectants. Exposure to other cleaning agents or disinfectants may lead to component failure.

NOTE: Refer to the injector system operation manual for additional operating instructions.

Indications for Use
The contents of this package are intended to be used in the delivery of contrast media or saline. These devices are indicated for use only with MEDRAD® Stellant CT Injection System.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

⚠ Warning

Follow aseptic technique when setting up and using disposables. Do not touch the end of the spike or connectors. Do not remove protective covers until ready to make connections. Properly discard the patient tube after each patient injection and the syringes and injector tubing within 12 hours. Contamination can result from failure to follow aseptic technique.

Do not use the multi-patient kit (syringes and Injector tubing) for longer than 12 hours. These disposables have been designed and validated for only 12 hours. Extended use of the multi-patient kit poses risks to the patient including contamination and injury due to device malfunction.

Do not exceed the media manufacturer's recommended use time. Refer to the media manufacturer's "Instructions for Use" for specific indications and use time.

For devices labeled for single use, please note: This product is intended for single use only. Do not re-sterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

For devices labeled for multi-use, please note: Do not re-sterilize or reprocess. Potential device failure includes significant component deterioration and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be reprocessed or re-sterilized.

Do not fill or inject unless the syringe is properly engaged. Patient injury can result if syringe is not properly engaged.


Do not re-install syringes onto an injector. Syringe sterility will be compromised, and patient infection can result. Use loaded syringes immediately. Discard unused loaded syringes. Contamination can occur if syringes are used to store contrast media.

⚠ Warning

Visually inspect contents and package before each use. Do not use contents if sterile package is opened or damaged. Patient or operator injury can result.


Do not exceed pressures identified on the front of the package; ensure that the fluid path is open. Use of greater pressures or occlusions in the fluid path can result in fluid leaks or tubing ruptures and patient or operator injury.

1 /


Configure Settings:
1. Press Setup/Fluid Delivery Setup button.
2. Set LPCT type to P8100 LPCT.
3. Set Auto Load Purge Volume to OFF to reduce fluid waste during syringe filling and refilling.
4. Set Fluid Control to OFF to decrease resistance felt when turning manual knobs and reduce leaking of fluid.
5. Press OK to save settings.

1 / 12h

Install Syringes:
1. Insert and push syringes until they snap into place.

MEDRAD® FluiDots
Round FluiDots: Syringe filled with Fluid


Ellipsoidal FluiDots: Syringe filled with Air

⚠ Caution

Ensure all connections are secure; do not over-tighten. This will help minimize leaks, disconnection, and component damage.

Install Fluid Transfer Set:
1. Remove white caps from both syringes.
2. Remove orange caps from fluid transfer sets.
3. Install fluid transfer sets onto syringe tips.

Install Connector Tube:
1. Remove white caps from fluid transfer set.
2. Remove small orange caps from connector tube.
3. Install the end of the connector tube that looks like a "T" onto Syringe A (Contrast).
4. Install opposite tube onto Syringe B (Saline).

Install Use-Time Label:
1. Write the syringe and injector tubing installation and expiration time on the label.

NOTE: Maximum use time is 12-hours.


2. Slide the label onto the transfer set or connector tube, so it is visible.

1 /

Install Single Patient Tube:
1. Remove a new patient tube from packaging.
2. Remove orange caps from connector tubing and patient tube.
3. Install patient tube onto the injector tubing.

NOTE: Leave white cap on patient tube until ready to connect to patient catheter.

Dust cap on: new patient tube


Dust cap off: used patient tube

⚠ Warning

Use care in handling and inserting spike into contrast bottle and saline bag. The spike is sharp and can cause personal injury.

Spike Fluids (AS NEEDED):
1. Remove cap from spike and fluid source.
2. Without touching the spike end, insert spike into fluid source.
3. Hang the spiked fluid source.
4. Repeat steps 1-3 for other fluid source.

Fill Syringes (AS NEEDED):
1. Fill syringes either automatically or manually (See Below).
2. Install heat maintainers.
3. Ensure air is expelled from syringes.

AUTOMATIC FILL:
Program a new or recall an existing protocol and press Auto Load on injector head.
Press Fill A or Fill B as needed.

NOTE: For first-time auto-fill, air remains in the syringe.

MANUAL FILL:
Press Move Piston on injector head and use forward and reverse piston controls to fill syringes to desired amount.

NOTE: Ensure final piston movement is in the forward direction. If the Syringe Undocked error appears, manually advance piston to re-dock syringe.

Prime Tubing:

NOTE: A whistling sound during priming is normal.

NOTE: While priming, tap the check valves on the patient tube to remove air bubbles.


1. Prime tubing.
2. Confirm that air is expelled from tubing.

AUTOMATIC PRIME:
Press Prime to prime injector and patient tubing.

MANUAL PRIME:
Either press Move Piston and use forward controls or turn manual knobs on injector head.

⚠ Warning

Clear all air from the syringes and tubing before connecting a patient to the injector. Air embolization can cause death or serious injury to the patient.


Connect to the Patient:
1. Rotate injector head down.
2. Confirm there is no air in syringes or tubing.
3. Remove cap from patient tube and connect to the patient.
4. Press I Checked for Air button on the injector head.

Prepare & Start Injection:
1. Ensure protocol is correct for current patient.
2. Lock the protocol at the control room monitor.
3. Press Arm.
4. Press yellow Start button.

AFTER EACH PATIENT
Discard Patient Tube:
1. Disconnect patient tube from the patient catheter.
2. Disconnect patient tube from injector tubing.
3. Discard patient tube in accordance with facility policy.
4. Immediately install new patient tube.
5. Repeat 1-per-patient steps as needed within the 12-hour use life.

WITHIN 12-HOUR USE TIME
Discard multi-patient kit:
1. Remove syringe heat maintainers.
2. Remove syringes. Rotate syringes 1/4 turn counter clockwise and pull out.
3. Remove used fluid sources.
4. Discard syringes, tubing, and fluid sources in accordance with facility policy.

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