
Apyx® MEDICAL

Energy Elevating Lives®

EN - Instructions for Use Disposable Split Adult Return Electrode Pad

REF APYX-ESREC

Languages:

English (USA) / EN-US – English (USA)	Español / ES - Spanish	한국어 / KO – Korean	Русский / RU – Russian
العربية / AR - Arabic	Eesti keel / ET – Estonian	Lietuvis / LT – Lithuanian	Slovensky / SK - Slovak
Български / BG - Bulgarian	Français / FR - French	Latvijas / LV – Latvian	Svenska / SV - Swedish
中国人 / CN - Chinese	Hrvatski / HR – Croatian	Nederlands / NL – Dutch	ไทย / TH – Thai
Česky / CS - Czech	Magyar / HU - Hungarian	Norsk / NO - Norwegian	Türkçe / TR - Turkish
Dansk / DA - Danish	Italiano / IT - Italian	Português / PT - Portuguese	Українська / UA Ukrainian
Deutsch / DE - German	ქართული / KA – Georgian	Română / RO – Romanian	Uzbek / UZ – Uzbek
Ελληνικά / EL - Greek	Қазақ / KK – Kazakh	Српски / RS - Serbian	Tiếng Việt / VI – Vietnamese

SYMBOLS GLOSSARY: Symbol Glossary is provided at www.apyxmedical.com/symbolglossary

Federal Law (USA) restricts this device to sale by or on the order of a physician. For a listing of indications, contraindications, precautions, and warnings, please refer to this document.

The device is not made with any natural rubber latex.

WARNING: This product can expose you to chemicals including Di (Iso-Nonyl) phthalate (DINP), which is known to the State of California to cause cancer. For more information, go to: www.P65Warnings.ca.gov.

INDICATIONS FOR USE/INTENDED USE:

The Disposable Split Adult Return Electrode Pad is a disposable, single-use, neutral electrode, which provides a return path for high frequency electrical current to the ESU device.

DEVICE DESCRIPTION AND OPERATING PRINCIPLE:

Disposable Split Return Electrodes (also known as grounding pads and dispersive electrodes) are single use, non-sterile, disposable devices. They are used in conjunction with an electrosurgical generator and accessories, such as electrosurgical handpieces, pencils, and electrodes. Grounding pads are used in monopolar surgery to disperse electrosurgical current delivered by a monopolar electrode to low density (reducing the risk of patient burns) and return the electrosurgical current back to the generator.

RF energy, which is generated from the generator travels through the interconnecting cable to an electrosurgical accessory (e.g.: handpiece), where the energy is delivered to an electrode to coagulate tissue. The electrical current passing through the patient's tissue seeking a return ground is directed to the grounding pad that is attached to the patient. The current will flow through the ground pad back to the generator, completing the circuit. This provides safety to the patient by dispersing the RF energy over the large surface area of the grounding pad to avoid burning the patient due to high current density.

The electrode is placed on or under the patient, and the cable is attached to the appropriate port on the electrosurgical generator. Grounding pads are attached directly to the patient's skin by a gel adhesive. This grounding pad is a split return electrode, and it is used with compatible generators to provide additional safety in that the contact quality is monitored.

CONTRAINDICATIONS: None

WARNINGS:

This medical device cannot be effectively cleaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection, or other risks of device failure to the patient.

TIONS:

- When using a split return electrode, the NEM (neutral electrode monitor or contact quality monitor) confirms the total resistance is within the preset safety range. A split pad must be used to utilize this safety feature.
- If the patient is repositioned for the surgical procedure, verify that the pad remains in proper contact with the skin and the cable remains connected to the pad. Failure to do so may result in patient burns.
- During the surgical procedure, ensure that no object obstructs the pad site.
- If higher than normal power settings are required, a problem may exist. Before increasing the power setting, verify that the pad has full contact with patient's skin. Check the cable and connectors and inspect all active accessories.
- Do not relocate the pad after initial application.
- Do not use damaged pads.
- Do not cut or modify the pad in any way.
- Review and follow any additional instructions, warnings, and cautions in the Instruction Manual for the generator in use.

INSTRUCTIONS FOR USE:

A. Inspection and Handling

Upon removal of the Return Electrode Pad from pouch visually inspect the pad for any damage that may include a missing clear protective backing; rips or tears; holes in the pad; damaged cord or cord connection (as applicable); missing or damaged connection tab; gel sticking to the clear backing when removed from the pad; adhesive missing from the perimeter on the clothed portion of the patient contact side of the pad. The pad should be handled to prevent inducing damage. After the removal of the clear protective backing, handle the pad to ensure that it remains flat as it is placed onto the patient's skin.

B. Selection and Preparation of Placement Site

1. Select a well-vascularized, muscular, and convex site in proximity to the surgical site. Do not place over scar, inflamed skin, fatty tissue, bony prominences, metal prostheses, ECG electrodes and cables, pacemakers, or where fluids may pool.
2. Shave all hair from the selected site. Remove any excess moisture.

C. Placement of Grounding Pad

1. Plug pad into generator.
2. Remove the backing from pad.
3. Apply to skin starting from one end of the pad and smooth out the opposite end.
4. Apply adequate pressure to the entire pad area to ensure adequate contact with the skin.

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D. Removal and Disposal of Grounding Pad

1. At the completion of the surgical procedure, remove the pad slowly and with care to avoid skin trauma.
2. Disconnect the cord from the generator.
3. Dispose of pad in accordance with local procedures.

Technical Description

General Description:


- ✓ Reduces electrosurgical current delivered by an electrode to low density and returns current back to electrosurgical generator
- ✓ Energy source: Compatible, electrosurgical generator is required

Specifications

Surface Area & Size:

- ✓ Pad size: 10.5cm (105mm) x 21.5cm (215mm)
- ✓ Cable length: 9.85 ft (3 meters)

Environmental Information

Patient Population

- ✓ Minimum weight: 25 lbs (11.4 kg)

Site of Use and Site Conditions

- ✓ Place on a well-vascularized, muscular, and convex site in proximity to the surgical site
- ✓ Do not place over scar, inflamed skin, fatty tissue, bony prominences, metal prostheses, ECG electrodes and cables, pacemakers, or where fluids may pool.
- ✓ Shave all hair from the selected site. Dry skin thoroughly.

Intended User Profile

Education:

- ✓ Trained Licensed Physician

Knowledge:

- ✓ Read and understand the information in the supplied IFU
- ✓ Understands electrosurgery and electrosurgical techniques

rmissible impairments:

- ✓ Mild reading vision impairment or corrected vision to 20/20

Intended Conditions for Use

Environment:

- ✓ Doctor office, Surgery center, Hospital, or similar environment
- ✓ Intended for professional use only
- ✓ Ambient Temperature +10° C to +40° C
- ✓ Indoor or contained area use

Disposal Information:

- ✓ Dispose of used grounding pad as biohazard medical waste.
- ✓ Do not reprocess the grounding pad.

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