Becton Dickinson
Infusion Therapy AB
Medication and Procedural Solution
Heather.Hagvik@bd.com

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Sterile single use Peripheral IV Catheter with Injection Port and Needle Guard BD Venflon™ Pro Safety is designed for infusion / injection and to provide access to the patient's peripheral blood vessels for infusion of fluids, blood transfusion and blood pressure control. These catheters can be used for patients of all categories, taking into account the size of the veins and the individual characteristics of the patient, the suitability of the catheter for the injected solutions and the duration of therapy.

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Instructions for use.

1. Carefully remove the needle cap without bending the needle and inspect the entire catheter device;
2. Align the needle and catheter tip;
3. Make sure that the bevel of the needle is directed upwards during insertion. Perform venipuncture. (18-22G Needle with BD Instaflash™ Technology: Reverse blood flow is visible in the catheter).
4. Verify that blood is flowing into the return flow control chamber.
5. After making sure that the blood is flowing back, tilt the entire catheter with the needle down slightly and make sure that the tip of the catheter is in the vein.
6. While holding the return flow control chamber, advance the catheter over the needle and into the vein. When removing the needle, push the vein proximal to the tip of the cannula and withdraw the needle while holding the needle retainer.
7. Do not bend the needle while removing it.
8. As the catheter slides off the needle, the plastic sleeve of the needle will begin to stretch from the needle retainer.
9. Once the needle is fully withdrawn, the protective cap will cover the tip of the needle.
10. Immediately place the needle in an appropriate sharps container.
11. Check if the port cover is closed.
12. Secure the catheter and apply a sterile dressing per protocol.

Indications for use. This device is intended for use with an auto injector set to a maximum pressure of 2240 kPa.
If injection is done through a port system, the port chamber must be compatible to work with the auto injector.

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Becton Dickinson Infusion Therapy AB, Sweden /
Becton Dickinson Infusion Therapy AB