PPL2122143 (2019-02-23)

Armada 14

PTA Catheter

CE
2797

- (GB) **INSTRUCTIONS FOR USE**
- (BG) **ИНСТРУКЦИИ ЗА УПОТРЕБА**
- (CZ) **NÁVOD NA POUŽITÍ**
- (DK) **BRUGSANVISNING**
- (NL) **GEBRUIKSAANWIJZINGEN**
- (EST) **KASUTUSJUHISED**
- (FIN) **KÄYTTÖOHJEET**
- (F) **MODE D'EMPLOI**
- (D) **GEBRAUCHSANWEISUNG**
- (GR) **ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ**
- (H) **HASZNÁLATI ÚTMUTATÓ**
- (I) **ISTRUZIONI PER L'USO**
- (LV) **LIETOŠANAS INSTRUKCIJAS**
- (LT) **NAUDOJIMO INSTRUKCIJOS**
- (PL) **INSTRUKCJA UŻYTKOWANIA**
- (P) **INSTRUÇÕES DE UTILIZAÇÃO**
- (RO) **INSTRUCȚIUNI DE UTILIZARE**
- (RUS) **ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ**
- (SK) **NÁVOD NA POUŽITIE**
- (SLO) **NAVODILA ZA UPORABO**
- (E) **INSTRUCCIONES DE USO**
- (S) **BRUKSANVISNING**
- (TR) **KULLANIM TALİMATLARI**
- (N) **BRUKSANVISNING**

Abbott

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Armada 14
PTA Catheter

I. DESCRIPTION

The Armada 14 Percutaneous Transluminal Angioplasty (PTA) Catheter is a dual lumen catheter with a balloon located near the distal atraumatic tip. One lumen is used for inflating the balloon and is accessed via the side leg port. The second lumen, starting at the straight entry port, allows access to the distal tip of the catheter for guide wire insertion (max 0.014"). The balloon has two radiopaque markers for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating section of the balloon and help in balloon placement.

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How the system is supplied

This device is EtO sterilized. Non-pyrogenic. Do not use if the package is open or damaged. Packaged with every Armada 14 PTA Catheter is a protective tubing that is positioned over the balloon, a packaging stylet and a balloon compliance card.

Storage

The Armada 14 PTA Catheter should be stored in a dark, dry and cool place. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that the device is used prior to the use by date on the package label.

II. INDICATIONS

The device is indicated to dilate stenoses in femoral, popliteal, infra popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The 2.0 to 4.0 mm balloon diameters are also indicated for post-dilatation of balloon-expandable stents up to 40 mm and self-expanding stents up to 80 mm in the vessels listed above.

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IV. WARNINGS/PRECAUTIONS

- This device should only be used by physicians who are experienced and have a thorough understanding of the clinical and technical aspects of PTA.
- One-time use only – do not resterilize! This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and/or chemical characteristics introduced under conditions of repeated use, cleaning, and/or resterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or spaces and diminished safety and/or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to patient and/or user.
- Do not use if inner package is damaged or opened.
- Employ aseptic techniques during removal from the package and during use.
- Any use for procedures other than those indicated in these instructions is not recommended.
- Use prior to the use by date.
- Carefully inspect the catheter prior to use to verify that it has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used.
- Precautions to prevent or reduce blood clotting should be taken when any catheter is used.
- Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.
- When the system is introduced into the vascular system, it should be manipulated only under high quality fluoroscopy.
- The Armada 14 PTA Catheter must always be introduced, moved and or withdrawn over a guide wire (max. 0.014").
- Never attempt to move the guide wire when the balloon is inflated.
- Never use air or any gaseous medium to inflate the balloon.
- Do not advance the Armada 14 PTA Catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the Armada 14 PTA Catheter through a smaller sized sheath introducer than indicated on the label.
- The size of the inflated balloon should be selected not to exceed the diameter of the artery immediately distal, or proximal, to the stenosis.
- Inflation in excess of the rated burst pressure may cause the balloon to rupture. Use of a pressure monitoring device is recommended.
- If a distal protection device is used, follow the manufacturer's instruction for use. Allow and maintain adequate distance between the Armada 14 PTA Catheter and the distal protection device to avoid engagement.
- Rated burst pressure and balloon fatigue testing of the Armada 14 PTA balloons within deployed stents has demonstrated the following:
 - The 2.0 to 4.0 mm balloon diameters can safely post-dilate balloon expandable stents up to 40 mm in length.
 - The 2.0 to 4.0 mm balloon diameters can safely post-dilate self expanding stents up to 80 mm in length.

The safety of using additional balloon diameters and/or lengths to post-dilate stents has not been established.

- When post-dilating stents, use a balloon length that is appropriate for the deployed stent length.

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- Abrupt closure
- Access site hematoma
- Aneurysm
- Angina
- Arrythmias
- Arteriovenous fistula
- Bleeding complications which may require transfusion
- Cerebral ischemia/transient ischemic attack (TIA)
- Death
- Embolism (air, tissue, thrombotic, systemic or device component)
- Fever/pyrogenic reaction
- Hypersensitivity or allergic reaction to contrast agents and drug reactions
- Hypertension/hypotension
- Infection
- Ischemia, including tissue ischemia, steal syndrome and necrosis
- Leg edema
- Myocardial ischemia or infarction
- Nausea and vomiting
- Neuropathies or nerve injury
- Occlusion
- Organ failure (single, multiple)
- Pain
- Palpitations
- Pseudoaneurysm
- Renal failure/insufficiency
- Restenosis
- Stroke/cerebrovascular accident (CVA)
- Vascular complications, including entry site, which may require vessel repair
- Vascular thrombosis
- Vessel injury, e.g. dissection, perforation
- Vessel spasm

VI. MATERIALS REQUIRED

- 1000 µ/500 cc Heparinized Normal Saline (HepNS)
- 4F sized Introducer sheath
- One inflation device (with pressure monitoring – recommended)
- One 3 way stopcock
- 0.014" guide wire
- 50% contrast diluted 1:1 with normal saline
- Two to three 10–20 cc syringes

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The diameter of the expanded balloon should not exceed that of the artery just distal, or proximal, to the stenosis. Verify that the selected accessories accommodate the balloon catheter as labeled.

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1) Fill an inflation device/syringe with a mixture of contrast medium and normal saline.
2) Remove the stylet and protective tubing from the balloon.
3) Attach the inflation device/syringe to the connector of the balloon lumen. Hold the catheter with the distal tip pointing downwards.

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4) Open the stopcock to the catheter; pull negative for 30 seconds; release to neutral.
5) Close the stopcock to the catheter; purge the inflation device/syringe of all air.
6) Repeat steps 4 through 6 until all air is expelled.
7) Slowly release device pressure to neutral.

VIII. INTRODUCTION AND DILATATION

Introduction of the system

1) The PTA catheter is designed to be introduced percutaneously using the Seldinger technique.
2) When the catheter is ready for introducing into the vascular system, the balloon protection tubing should be removed completely from the catheter.
3) Place the prepared catheter over a pre-positioned guide wire and advance the tip to the introduction site. It is advantageous to use the balloon catheter with an introducer to facilitate entry.

 Note: Perform all further catheter manipulations under fluoroscopy.
4) Position the catheter with the center of the balloon in the middle of the stenosis. The radiopaque marker bands indicate the length of the balloon.

Inflation of the balloon

5) When an acceptable position has been reached, inflate the balloon to achieve the desired dilatation. Caution: Do not exceed the rated burst pressure. Higher pressure can damage the balloon or catheter or overdistend the selected artery.
6) Deflate the balloon by aspirating the inflation syringe or inflation device.
7) Maintaining a vacuum in the balloon, withdraw the catheter.

 Note: Gentle counterclockwise twisting motion of the balloon may ease withdrawal through the sheath or from the percutaneous entry site. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath as one unit.
8) Please note that if multiple balloon-inflations/deflations have taken place, some resistance can occur upon device withdrawal.
9) The results of PTA should be checked by angiography.

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Abbott Vascular has exercised reasonable care in the manufacture of this device. Abbott Vascular excludes all warranties, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device as well as factors relating to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond Abbott Vasculars' control directly affect this device and the results obtained from its use. Abbott Vascular shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

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Percutaneous Transluminal Angioplasty (PTA) Catheter

I. DESCRIPTION

The Armada 18 Percutaneous Transluminal Angioplasty (PTA) Catheter is an over-the-wire (OTW) catheter with a balloon located near the distal atraumatic tip. One lumen is used for inflating the balloon and is accessed via the side leg port. The second lumen, starting at the straight entry port, allows access to the distal tip of the catheter for guide wire insertion (max 0.018"). The balloon has two radiopaque markers for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating section of the balloon and help in balloon placement.

The balloon is inflated using the side leg port, at which point the balloon expands to a known diameter at specific pressure. The working pressure range for the balloon is between the nominal pressure and the rated burst pressure. All balloons distend to sizes above the nominal diameter at pressures greater than the nominal pressure.

The balloon diameter and length are printed on the luer. Refer also to the package label for information about catheter length, guide wire compatibility, nominal and rated burst pressure, and sheath compatibility.

How the system is supplied

This device is EtO sterilized. Non-pyrogenic. Do not use if the package is open or damaged. Packaged with every Armada 18 PTA Catheter is a protective tubing that is positioned over the balloon, a packaging stylet and a balloon compliance card.

Storage

The Armada 18 PTA Catheter should be stored in a dark, dry and cool place. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that the device is used prior to the use by date on the package label.

II. INDICATIONS

The Armada 18 is indicated to dilate stenosis in femoral, popliteal, infra-popliteal, tibial, peroneal, and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. In addition, the device is also indicated for post-dilatation of balloon expandable and self-expanding stents.

III. CONTRAINDICATIONS

- Inability to cross lesion with a guide wire
- Use in the coronary arteries.

IV. WARNINGS / PRECAUTIONS

- This device should only be used by physicians who are experienced and have a thorough understanding of the clinical and technical aspects of PTA.
- One-time use only – do not resterilize! This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to patient and / or user.
- Do not use if inner package is damaged or opened.
- Employ aseptic techniques during removal from the package and during use.
- Any use for procedures other than those indicated in these instructions is not recommended.
- Use prior to the use by date.
- Carefully inspect the catheter prior to use to verify that it has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used.
- Precautions to prevent or reduce blood clotting should be taken when any catheter is used.
- Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.
- When the system is introduced into the vascular system, it should be manipulated only under high quality fluoroscopy.
- The Armada 18 PTA Catheter must always be introduced, moved and or withdrawn over a guide wire (max. 0.018").
- Never attempt to move the guide wire when the balloon is inflated.
- Never use air or any gaseous medium to inflate the balloon.
- Do not advance the Armada 18 PTA Catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the Armada 18 PTA Catheter through a smaller sized sheath introducer than indicated on the label.
- The size of the inflated balloon should be selected not to exceed the diameter of the artery immediately distal, or proximal, to the stenosis.
- Inflation in excess of the rated burst pressure may cause the balloon to rupture. Use of a pressure monitoring device is recommended.
- If a distal protection device is used, follow the manufacturer's instruction for use. Allow and maintain adequate distance between the Armada 18 PTA Catheter and the distal protection device to avoid engagement.
- Rated burst pressure and balloon fatigue testing of the Armada 18 PTA balloons within deployed stents has demonstrated that the Armada 18 can safely post-dilate balloon expandable and self-expanding stents.
- When post-dilating stents, use a balloon length that is appropriate for the deployed stent length.

V. POTENTIAL COMPLICATIONS

The following complications may occur as a result of PTA, but may not be limited to:

- Abrupt closure
- Allergic reaction (contrast medium, drug, or stent material)
- Aneurysm, pseudoaneurysm, or arteriovenous fistula
- Angina or coronary ischemia
- Arrhythmias (including premature beats, bradycardia, atrial or ventricular tachycardia, atrial or ventricular fibrillation)
- Bleeding complications requiring transfusion or surgical intervention
- Death
- Detachment and / or implantation of a component of the system
- Embolization, arterial or other (air, tissue, plaque, thrombotic material, device)
- Emergent or urgent surgery
- Fever
- Hematoma or hemorrhagic event with or without surgical repair
- Hyperperfusion syndrome
- Hypotension or hypertension
- Infection
- Ischemia or infarction of tissue or organ not covered under other adverse events
- Myocardial infarction
- Pain (limb or catheter site)
- Peripheral nerve injury
- Pulmonary embolism
- Renal failure or insufficiency
- Restenosis of vessel
- Shock
- Stroke, cerebrovascular accident (CVA), or transient ischemic attack (TIA)
- Target limb loss (amputation of toe, foot, and / or leg)
- Vascular thrombosis or occlusion at puncture site, treatment site, or remote site
- Venous thromboembolism
- Vessel dissection, perforation, or rupture
- Vessel spasm or recoil
- Worsening claudication or rest pain

VI. MATERIALS REQUIRED

- 1000 µ/500 cc Heparinized Normal Saline (HepNS)
- Introducer sheath size according to label
- One inflation device (with pressure monitoring–recommended)
- One 3 way stopcock
- 0.018" guide wire maximum
- Contrast diluted 1:1 with normal saline
- Two to three 10–20 cc syringes

VII. SELECTION AND PREPARATION OF DEVICE AND COMPATIBILITY WITH ACCESSORIES

Selection of balloon size and compatibility of system to accessories

The diameter and length of the expanded balloon should approximate the size of the artery just distal or proximal to the lesion and length of the stenosis. Verify that the selected accessories accommodate the balloon catheter as labeled.

Preparation of PTA catheter

It is very important to check before use that the packaging has not been damaged in a way that might have rendered the catheter unsterile. It is also important at the same time to verify that the chosen catheter is the correct one for the planned procedure. Do not inflate the PTA catheter during preparation of the device.

Perform the following steps to remove all air and verify the integrity of the PTA catheter:

1) Fill an inflation device / syringe with a mixture of contrast medium and normal saline.
2) Remove the stylet and protective tubing from the balloon.
3) Attach the inflation device / syringe to the connector of the balloon lumen. Hold the catheter with the distal tip pointing downwards.
4) Open the stopcock to the catheter; pull negative for 30 seconds; release to neutral.
5) Close the stopcock to the catheter; purge the inflation device / syringe of all air.
6) Repeat steps 4 through 6 until all air is expelled.
7) Open the stopcock to the catheter; pull negative and maintain negative.

VIII. INTRODUCTION AND DILATATION

Introduction of the system

1) The PTA catheter is designed to be introduced percutaneously using the Seldinger technique.
2) Prior to introducing the catheter into the vascular system, confirm that the balloon protection tubing has been removed completely and the catheter is under negative pressure.
3) Place the prepared catheter over a pre-positioned guide wire and advance the tip to the introduction site. It is advantageous to use the balloon catheter with an introducer to facilitate entry.

 Note: Perform all further catheter manipulations under fluoroscopy.

 Note: If unusual resistance is felt at any point during delivery, do not force passage. Resistance may indicate a problem and the use of excessive force may result in device damage and vessel injury. Maintain guide wire placement and remove the catheter and sheath as a **single unit**.
4) Position the catheter with the center of the balloon in the middle of the stenosis. The radiopaque marker bands indicate the length of the balloon.

Inflation of the balloon

5) When an acceptable position has been reached, inflate the balloon to achieve the desired dilatation.

 Caution: Do not exceed the rated burst pressure. Higher pressure can damage the balloon or catheter or over distend the selected artery.
6) Deflate the balloon by aspirating the inflation syringe or inflation device.
7) Maintain vacuum during withdrawal.
8) Gentle counterclockwise twisting motion of the balloon may ease withdrawal through the sheath or from the percutaneous entry site. Do not twist the catheter more than three (3) full turns. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath as one unit.
9) The results of PTA should be checked by angiography.

IX. PRODUCT INFORMATION DISCLOSURE

Abbott Vascular has exercised reasonable care in the manufacture of this device. Abbott Vascular excludes all warranties, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device as well as factors relating to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond Abbott Vascular's control directly affect this device and the results obtained from its use. Abbott Vascular shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Armada is a trademark of the Abbott Group of Companies.

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©2013 Abbott