# Perclose ProGlide

**Suture-Mediated Closure (SMC) System**

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# Perclose ProGlide Suture-Mediated Closure (SMC) System

**INSTRUCTIONS FOR USE**

**Table of Contents**

1.0 CAUTION
2.0 DEVICE DESCRIPTION
Figure 1: Perclose ProGlide Suture-Mediated Closure System
3.0 HOW SUPPLIED
4.0 INDICATIONS
5.0 CONTRAINDICATIONS
6.0 WARNINGS
7.0 PRECAUTIONS
8.0 SPECIAL PATIENT POPULATIONS
9.0 POTENTIAL ADVERSE EVENTS
10.0 THE PERCLOSE PROGLIDE SMC SYSTEM CLINICAL PROCEDURE
10.1 Examination and Selection of Products
10.2 Access Site and Puncture Considerations
10.3 Single SMC Device Placement, Including Optional Maintaining Wire Access and Pre-Close Techniques
10.3.1 Maintaining Wire Access During Knot Advancement (Closing over the wire)
10.3.2 Pre-Close Technique
10.3.3 Knot Advancement Using the Suture Trimmer
10.4 Multiple SMC Device Placement, Utilizing Pre-Close and Maintaining Wire Access Techniques
10.5 Suture Breakage
10.6 Post-Procedure Patient Management
10.7 Recommendation for Patient Ambulation and Discharge
11.0 PRODUCT INFORMATION DISCLOSURE
12.0 PATENTS AND TRADEMARKS

**TO ENSURE PROPER DEPLOYMENT AND USE OF THIS DEVICE AND TO PREVENT INJURY TO PATIENTS, READ ALL INFORMATION CONTAINED IN THESE INSTRUCTIONS FOR USE.**

## 1.0 CAUTION

This device should only be used by physicians (or allied healthcare professionals, authorized by, or under the direction of, such physicians) who are trained in diagnostic and / or interventional catheterization procedures and who have been trained by an authorized representative of Abbott Vascular. Prior to use, the operator must review the Instructions for Use and be familiar with the deployment techniques associated with the use of this device.


## 2.0 DEVICE DESCRIPTION

The Perclose ProGlide Suture-Mediated Closure (SMC) System is designed to deliver single monofilament polypropylene suture to close femoral vessel puncture sites following diagnostic or interventional catheterization procedures. This Perclose ProGlide SMC device is composed of a plunger, handle, guide, and sheath. The Perclose ProGlide tracks over a standard 0.038" (0.97 mm) (or smaller) guide wire. A hemostasis valve restricts the blood flow through the sheath with or without the guide wire in place. The guide houses the needles, and the foot, and precisely controls the placement of these needles around the puncture site. The handle is used to stabilize the device during use. The plunger advances the needles and is used to retrieve the suture. A marker lumen is contained within the guide, with the intraluminal port of the lumen positioned at the distal end of the guide. Proximally, the marker lumen exits from the body of the device. The marker lumen allows a pathway for back-bleeding (obtaining mark) from the femoral vessel to ensure proper device positioning.

The Perclose Suture Trimmer is included, and is designed to position the tied suture knot to the top of the access site and to trim the trailing limbs of suture.

The Perclose ProGlide SMC System is designed for use in closing common femoral artery and vein access sites. **The Perclose ProGlide SMC System is depicted in Figure 1.**

**Figure 1: Perclose ProGlide Suture-Mediated Closure System**

## 3.0 HOW SUPPLIED

The Perclose ProGlide SMC device and accessories are provided sterile and non-pyrogenic in unopened, undamaged packages. Products are sterilized with ethylene oxide and intended for single use only. This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. Store in a cool, dry place.

**Perclose ProGlide SMC System Includes:**

One (1) Perclose ProGlide SMC device
One (1) Perclose Suture Trimmer

## 4.0 INDICATIONS

The Perclose ProGlide SMC System is indicated for the percutaneous delivery of suture for closing the common femoral artery and vein access site of patients who have undergone diagnostic or interventional catheterization procedures:

- For access sites in the common femoral artery using 5F to 21F sheaths. For sheath sizes greater than 8F, at least two devices and the pre-close technique are required.
- For access sites in the common femoral vein using 5F to 24F sheaths. For sheath sizes greater than 8F, at least one device and the pre-close technique is required.

## 5.0 CONTRAINDICATIONS

There are no known contraindications to the use of this device. Attention is drawn to Sections 6.0 WARNINGS and 7.0 PRECAUTIONS.

## 6.0 WARNINGS

Do not use the Perclose ProGlide SMC device or accessory if the packaging or sterile barrier has been previously opened or damaged or if the components appear to be damaged or defective.

DO NOT RESTERILIZE OR REUSE. The Perclose ProGlide SMC device and accessory are intended for single use only.

Do not use the Perclose ProGlide SMC System if the sterile field has been broken where bacterial contamination of the sheath or surrounding tissues may have occurred, since such a broken sterile field may result in infection.

Do not use the Perclose ProGlide SMC System if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma. Perform a femoral angiogram to verify the location of the puncture site. **NOTE:** This may require both a Right Anterior Oblique (RAO) and Left Anterior Oblique (LAO) angiogram to adequately visualize where the sheath enters the femoral artery or vein.

Do not use the Perclose ProGlide SMC System if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a hematoma or retroperitoneal bleed.

Do not use the Perclose ProGlide SMC System if the puncture site is located in the superficial femoral artery or the profunda femoris artery, or the bifurcation of these vessels, since such puncture sites may result in a pseudoaneurysm, intimal dissection, or an acute vessel closure (thrombosis of small artery lumen). Perform a femoral angiogram to verify the location of the puncture site. **NOTE:** This may require both a Right Anterior Oblique (RAO) and Left Anterior Oblique (LAO) angiogram to adequately visualize where the sheath enters the femoral artery or vein.

## 7.0 PRECAUTIONS

- Prior to use, inspect the Perclose ProGlide SMC System to ensure that the sterile packaging has not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
- As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the Perclose ProGlide SMC System. Employ appropriate groin management, as per hospital protocol, post procedure and post hospital discharge to prevent infection.
- Use a single wall puncture technique. Do not puncture the posterior wall of the vessel.
- Do not deploy the Perclose ProGlide SMC device at an angle greater than 45 degrees, as this may cause a cuff miss.
- There are no reaccess restrictions if previous access site repairs were achieved with Abbott Vascular SMC devices.
- If significant blood flow is present around the Perclose ProGlide SMC device, do not deploy needles. Remove the Perclose ProGlide SMC device over a 0.038" (0.97 mm) (or smaller) guide wire and insert an appropriately sized introducer sheath.
- When pushing the plunger assembly to advance the needles, stabilize the device to ensure the device does not twist or move forward during deployment. Twisting the device could lead to needle deflection resulting in a cuff miss. Do not use excessive force or repeatedly push the plunger assembly. Excessive force on the plunger during deployment could potentially cause breakage of the device, which may necessitate intervention and / or surgical removal of the device and vessel repair.
- Do not apply excessive force to the lever when returning the foot to its original position **(marked #4)** down to the body of the device. Do not attempt to remove the device without closing the lever. Excessive force on the lever of the device or attempting to remove the device without closing the lever could cause breakage of the device and / or lead to vessel trauma, which may necessitate intervention and / or surgical removal of the device and vessel repair.
- **Do not advance or withdraw the Perclose ProGlide SMC device against resistance until the cause of that resistance has been determined** (see Section 10.3 SINGLE SMC DEVICE PLACEMENT). **Excessive force used to advance or torque the Perclose ProGlide SMC device should be avoided, as this may lead to significant vessel damage and / or breakage of the device, which may necessitate intervention and / or surgical removal of the device and vessel repair.**

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- If excessive resistance in advancing the Perclose ProGlide SMC device is encountered, withdraw the device over a 0.038" (0.97 mm) (or smaller) guide wire and reinsert the introducer sheath or use manual compression.
- Remove the Perclose ProGlide sheath before tightening the suture. Failure to remove the sheath prior to tightening the suture may result in detachment of the tip of the sheath.
- In using this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing damage due to application of surgical instruments such as clamps, forceps or needle holders.
- For closure of access sites using a 5F–8F procedural sheath, use manual compression in the event bleeding from the femoral access site persists after the use of the Perclose ProGlide SMC device.
- For closure of access sites using a procedural sheath greater than 8F, depending on the type of interventional catheterization procedure, use manual compression methods, compression assisted devices, and / or other appropriate treatment methods in the event bleeding from the femoral access site persists after the use of Perclose ProGlide SMC devices.

**8.0 SPECIAL PATIENT POPULATIONS**

The safety and effectiveness of the Perclose ProGlide SMC devices have not been established in the following special patient populations:

- Patients with introducer sheaths < 5F or > 21F in the femoral artery during the catheterization procedure.
- Patients with introducer sheaths < 5F or > 24F in the femoral vein during the catheterization procedure.
- Patients with small femoral arteries or veins (< 5 mm in diameter).
- Patients with access sites above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks.
- Patients having access in vessels other than the common femoral artery or vein.
- Patients having a hematoma, pseudoaneurysm or arteriovenous fistula present prior to sheath removal.
- Patients with femoral artery calcium which is fluoroscopically visible at access site.
- Patients with severe claudication, iliac or a femoral artery diameter stenosis greater than 50% or previous bypass surgery or stent placement in the vicinity of access site.
- Patients with access sites in vascular grafts.
- Patients with prior intra-aortic balloon pump at access site at any time prior.
- Patients with ipsilateral arterial access sites punctured and compressed within 48 hours of closure. **NOTE:** The previous / initial puncture site may have the potential to re-bleed due to an unstable clot and / or anticoagulants, even if the new puncture site is successfully closed with a Perclose ProGlide SMC device.
- Patients where there is difficulty inserting the introducer sheath or greater than one ipsilateral arterial puncture at the start of the catheterization procedure.
- Patients with antegrade punctures.
- Patients receiving glycoprotein IIb/IIIa inhibitors before, during, or after the catheterization procedure.
- Patients who are pregnant or lactating.
- Patients with bleeding diathesis or coagulopathy.
- Patients younger than 18 years of age.
- Patients who are morbidly obese (Body Mass Index $\geq$ 40 kg/m$^2$).
- Patients with active systemic or cutaneous infection or inflammation.

Before considering early discharge, assess the patient for the following clinical conditions:

- Conscious sedation
- Anticoagulation, thrombolytic, or antiplatelet therapy
- Unstable cardiac status
- Hematoma at the closure site
- Hypotension
- Pain while walking
- Bleeding at the closure site
- Any comorbid condition requiring observation

The presence of any of the above factors has generally led to the deferral of early discharge recommendations.

**9.0 POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with the use of suture-mediated closure devices may include, but are not limited to, the following:

- Allergic reaction or hypersensitivity to device components
- Anemia
- Arterial stenosis / occlusion
- Arteriovenous fistula
- Bleeding / hemorrhage
- Bruising / hematoma
- Death
- Deep vein thrombosis
- Device entrapment
- Device failure / malfunction / misplacement
- Diminished pulses distal to closure site
- Embolism
- Hypotension / hypertension
- Infection / sepsis
- Inflammation
- Intimal tear / dissection
- Ischemia distal to closure site
- Nerve injury
- Numbness
- Pain
- Perforation
- Pseudoaneurysm
- Pulmonary embolism
- Retroperitoneal hematoma / bleeding
- Thrombus formation
- Vascular injury
- Vasoconstriction / vasospasm
- Vasovagal episode
- Wound dehiscence

**10.0 THE PERCLOSE PROGLIDE SMC SYSTEM CLINICAL PROCEDURE**

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the Perclose ProGlide SMC System. The techniques and procedures described below are not intended as a substitute for the operator's experience and judgment in treating any specific patients.

**10.1 Examination and Selection of Products**

1. Select the Perclose ProGlide SMC device(s) for closure of 5F to 21F introducer sheath access sites in the common femoral artery and 5F to 24F introducer access sites in the common femoral vein.
2. After carefully inspecting the packaging of the Perclose ProGlide SMC device and accessories for damage to the sterile barrier, remove the device from the package.
3. Exercise care when using additional instruments, such as clamps, forceps or needle holders during device handling, to reduce the possibility of accidental device breakage or damage to the suture.
4. Verify marker lumen patency by flushing the marker lumen with saline until the saline exits the marker port. **Do not use the Perclose ProGlide SMC device if the marker lumen is not patent.**

**10.2 Access Site and Puncture Considerations**

1. An extremely deep tissue tract can influence needle trajectory, which prevents the Perclose ProGlide SMC device needles from engaging the cuffs, or secure knot tying, as the Suture Trimmer may not be able to advance the knot to the vessel wall for complete apposition before locking the knot. An extremely deep tissue tract may require a long access needle and / or upon inserting the Perclose ProGlide device, requires compression of the subcutaneous tissue to be able to obtain pulsatile flow.
2. Before inserting the access needle, use of ultrasound guidance to visualize the access site or fluoroscopy to visualize the femoral head is recommended. When using the femoral head as a reference point, target the middle of the femoral head as the puncture site. Performing a femoral angiogram through the introducer sheath (or procedural sheath) to verify that the access site is in the common femoral artery or vein is recommended before anticoagulants are given.
3. Puncture the anterior wall of the common femoral artery or vein at an angle of approximately 45 degrees. Avoid side wall or posterior wall punctures.
4. Prior to deployment of the Perclose ProGlide SMC device, perform a femoral angiogram to evaluate the access site for vessel size, calcium deposits, tortuosity, and for disease or dissections of the vessel wall to avoid device cuff misses (device needles not engaging with the cuffs) and / or posterior wall suture placement and possible ligation of the anterior and posterior walls of the vessel. The puncture should be proximal to the bifurcation of the superficial femoral artery and the profunda femoris branch and distal to the inferior margin of the inferior epigastric artery.
5. There are **no** reaccess restrictions after using Abbott Vascular closure devices.

**10.3 Single SMC Device Placement, Including Optional Maintaining Wire Access and Pre-Close Techniques**

**The following instructions detail the deployment sequence to close the access site of a catheterization procedure performed with a single device. For venous sheath sizes greater than 8F, at least one device and the pre-close technique is required. Follow Section 10.3.2 to close with one device using the pre-close technique.**

1. Place a 0.038" (0.97 mm) (or smaller) guide wire through the procedural (or introducer) sheath. Remove the procedural sheath while applying pressure on the groin to maintain hemostasis.
2. Backload the device over the guide wire until the guide wire exit port of the device sheath is just above the skin line. Remove the guide wire before the exit port crosses the skin line.
3. Continue to advance the device in the vessel until brisk pulsatile flow of blood is evident from the marker lumen. Position the device at a 45-degree angle. Deploy the foot by lifting the lever **(marked #1)** on the top of the handle. **Do not deploy the foot in the artery unless brisk pulsatile flow of blood ("mark") is evident from the marker lumen. In the vein, the flow of blood may not be pulsatile or blood may only fill the marker lumen.**
4. Gently pull the device back to position the foot against the vessel wall. If proper position of the foot has been achieved, tactile sensation will be felt. In the artery, blood marking will cease or be significantly reduced to a slight drip. In the vein, there may be no change in blood marking. If marking does not stop or significantly change, evaluate the angiogram for device position in the vessel, vessel size, calcium deposits, tortuosity, disease and for location of the puncture to ensure footplate is not in bifurcation or side branch. Reposition the device to stop blood marking or reinsert the wire, remove the device to hold manual compression or insert a new sheath.
5. While maintaining device position, stabilize the device with your free hand (the one not used to deploy the device) to maintain the gentle retraction and to ensure the device does not twist or move forward during deployment. Use your other hand to deploy needles by pushing on the plunger assembly **(in the direction marked #2)** until you visually confirm that the collar of the plunger makes contact with the proximal end of the body. Do not use excessive force or repeatedly push the plunger assembly. After visually confirming contact with the body of the device only **one time**, this step is complete.

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6. Using your thumb as a fulcrum on the handle, gently disengage the needles by pulling the plunger assembly back **(in the direction marked #3)** and completely remove the plunger and needles from the body of the device. One suture limb will be attached to the anterior needle. The posterior needle will be free of suture. Pull back on the plunger until the suture is taut, which confirms that all the sutures have been fully retracted from the body of the device.
7. Do not attempt to redeploy the needles if the suture limb is not attached to the anterior needle. Reinsert the guide wire, and remove the device over the wire. Insert a new Perclose ProGlide SMC device to complete the procedure.
8. Use the QuickCut suture-trimming mechanism located on the handle to cut the suture from the anterior needle distal of the link. Use of a new, sterile scalpel or scissors is optional.
9. Relax the device and then return the foot to its original position by pushing the lever **(marked #4)** down to the body of the device. Do not attempt to remove the device without closing the lever.
10. Withdraw the Perclose ProGlide SMC device until the guide wire port exits the skin line.
11. Grasp the suture adjacent to the device sheath and pull the suture ends through the distal end of the proximal guide. The rail suture limb is blue and is the longer of the two suture limbs. This rail suture limb will be used to advance the knot. The shorter, non-rail limb is white tipped and will be used to lock the knot.
12. Continue with Knot Advancement (Step 10.3.3), if closing over the wire, continue with Section 10.3.1.

#### 10.3.1 Maintaining Wire Access During Knot Advancement (Closing over the wire)

**If the operator chooses to maintain wire access, reinsert the guide wire after Section 10.3 Step 11. Follow Steps A–H when closing over the wire.**

A. There should be sufficient guide wire exposed out of the guide wire exit port before removing the device.
B. Wrap the rail limb (long blue limb) of suture around your left index finger, low, close to skin level.
C. Remove the Perclose ProGlide SMC device with the right hand, while maintaining an adequate length of guide wire inside the vessel. This allows placement of another Perclose ProGlide SMC device or a sheath in the event that hemostasis is not obtained.
D. While removing the device with the right hand, simultaneously advance the knot to the access site by applying slow, consistent increasing tension to the rail suture limb, keeping the suture coaxial to the tissue tract. (Do not advance the knot with the Suture Trimmer until the wire has been completely removed from the patient.)
E. Assess the site for adequate hemostasis. If bleeding is controlled, the operator should then remove the guide wire. Next, use the Suture Trimmer, as described in Step 13, to advance and tighten the knot until complete closure is achieved. Confirm the security of the knot by having the patient cough and / or bend his / her leg. Active testing for hemostasis is only for 5F-8F sheaths. For sheaths greater than 8F, active confirmation should not be performed; only visual confirmation of hemostasis should be deployed. Additionally, patients may be able to move freely in bed without head of bed or leg restrictions if the close is successful.
F. If the wire is still in place and suture breakage occurs during knot advancement, or if hemostasis is not achieved, another Perclose ProGlide SMC device can be used to complete the procedure or a sheath can be inserted. Either remove the broken suture limbs or cut the suture limbs close to the knot using the Suture Trimmer (as in Section 10.3.3 Step J) or a new, sterile scalpel or scissor.
G. Care should be taken to avoid excessive force if the insertion of another Perclose ProGlide SMC device or introducer sheath is required. To avoid resistance, use a sheath small enough to avoid undue force.
H. If hemostasis cannot be achieved after the wire is removed, apply manual compression.

#### 10.3.2 Pre-Close Technique

**The Perclose ProGlide suture can be placed around the access site at the beginning of the procedure and knot advancement can be placed on hold until the procedure is complete. Follow Steps A–F when using a pre-close technique.**

A. Complete Section 10.3 Steps 1–11 to remove suture from the distal end of the proximal guide.
B. Immediately place a shodded hemostat or clamp to hold the two suture limbs together at the distal end of the non-rail suture limb (shortest blue limb with white tip). To prevent knot advancement or locking of the knot, care must be taken not to pull on the individual suture limbs until the clamp is securely holding the two limbs together.
C. After securing the suture limbs and before inserting the procedural sheath, gently pull on the clamp until the suture is taut to remove any suture slack from the tissue tract.
D. Place the clamped suture under a sterile towel during the procedure. **NOTE:** The monofilament suture can be damaged by opening and closing the clamp. If you want to attach the suture to the drape, it is recommended that you use a second clamp with the tip placed through the handle of the first clamp and attach the second clamp to the drape.
E. At the end of the catheterization, reinsert the guide wire into the procedural sheath.
F. Follow Section 10.3.1 to close the access site over the guide wire. Note: In Step C, the procedural sheath will be removed instead of a Perclose ProGlide device.

#### 10.3.3 Knot Advancement Using the Suture Trimmer

A. Securely wrap the rail (longer, blue) limb of the suture around your left forefinger, close to the skin.
B. Do not tighten the suture around the sheath. Completely remove the device or the sheath (if the device was deployed at the beginning of the catheterization procedure) from the vessel, while simultaneously pulling gently on the rail limb. Always keep the suture coaxial to the tissue tract.
C. With the rail (longer, blue) suture limb securely wrapped around your left forefinger, place the rail limb into the Suture Trimmer utilizing the following steps: Retract the Thumb Knob on the Suture Trimmer with the right hand. Place the Suture Trimmer under the suture limb making an "x" or a "cross" between the suture limb and the mid-point of the Suture Trimmer. Slide the Suture Trimmer back to load the suture into the Suture Gate located at the distal end of the Suture Trimmer. Keeping the Thumb Knob retracted, turn the Suture Trimmer coaxial to the suture and then release the Thumb Knob to capture the suture in the Suture Gate. Releasing the Thumb Knob before the suture is coaxial to the Suture Trimmer can cause the suture to be caught within the sliding mechanism at the distal tip and damage the suture. Once the suture is loaded correctly, the Suture Trimmer should slide easily.
D. With the suture limb and Suture Trimmer coaxial to the tissue tract, move the knot to the vessel surface by advancing the Suture Trimmer with the right hand while placing slow, consistent increasing tension on the rail suture with the left forefinger. Avoid quick or jerking type movements with the suture limbs. The Suture Trimmer and suture limbs should always remain coaxial to the tissue tract. The Thumb Knob should be at 12 o'clock (facing the ceiling) and the Suture Trimmer should not be rotated.
E. With the rail (longer, blue) suture limb securely wrapped around the left forefinger, place the Suture Trimmer under the left thumb to assume a single-handed position and complete knot advancement with slow, consistent increasing tension until the suture is taut (guitar string tightness).
F. With the Suture Trimmer in place and suture taut, tighten the knot by gently pulling the non-rail (shorter, white tipped) suture limb, keeping it coaxial to the tissue tract.
G. Hemostasis of the access site is achieved when the knot is fully advanced to the vessel surface, the slack is gently pulled from the knot with the non-rail limb while the Suture Trimmer holds tension on the rail limb of the suture, and the tissue is in complete apposition.
H. Remove the Suture Trimmer from the tissue tract, relax tension on the suture, and test for hemostasis by having the patient cough or bend his / her leg. If hemostasis has not been achieved, assume the single-handed position for 20 seconds, or until hemostasis is achieved. Secure the knot again by gently pulling on the white tipped non-rail suture limb. Do not apply excessive pressure to the suture.
I. Confirm the security of the knot by having the patient cough and / or bend his / her leg. Active testing for hemostasis is only for 5F-8F sheaths. For sheath greater than 8F, active confirmation should not be performed; only visual confirmation of hemostasis should be deployed. Additionally, patients may be able to move freely in bed without head of bed or leg restrictions if the close is successful.
J. Once hemostasis is achieved, use the Suture Trimmer to trim the sutures below the skin. While holding both suture limbs together and pulled taut, load both suture limbs into the Suture Trimmer (as described above in Step C) and advance the Suture Trimmer to the vessel surface. Trim the sutures by pulling back on the red Trimming Lever. Keep the Trimming Lever pulled back while removing the Suture Trimmer and trimmed suture limbs from the tissue tract. If only one suture limb has been loaded and trimmed, repeat the same technique on the other suture limb.
K. If hemostasis cannot be achieved, apply manual compression until hemostasis is achieved.

### 10.4 Multiple SMC Device Placement, Utilizing Pre-Close and Maintaining Wire Access Techniques

**The following introductions detail the deployment sequence for closing the access site of a catheterization procedure with multiple devices using the pre-close technique.**

1. Place a 0.038" (0.97 mm) (or smaller) guide wire through the introducer sheath. Remove the introducer sheath while applying pressure on the groin to maintain hemostasis.
2. Backload the SMC device over the guide wire until the guide wire exit port of the device sheath is just above the skin line. Remove the guide wire before the exit port crosses the skin line.
3. Place the first Perclose ProGlide SMC devices over the guide wire. Continue to advance the device until brisk pulsatile flow of blood is evident from the marker lumen. In the vein, the flow of blood may not be pulsatile or blood may only fill the marker lumen. The device lever **(marked #1)** and logo should be facing the ceiling (12 o'clock).
4. Rotate the device approximately 30 degrees towards the patient's right side. Position the device at a 45-degree angle. Deploy the foot by lifting the lever **(marked #1)** on top of the handle. **Do not deploy the foot in the artery unless brisk pulsatile flow of blood ("mark") is evident from the marker lumen. In the vein, flow of blood may not be pulsatile or blood may only fill the marker lumen.**
5. Gently pull the device back to position the foot against the vessel wall. If proper position of the foot has been achieved, tactile sensation will be felt. In the artery, blood marking will cease or be significantly reduced to a slight drip. In the vein, there may be no change in blood marking. If marking does

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not stop or significantly change, evaluate the angiogram for device position in the vessel, vessel size, calcium deposits, tortuosity, disease and for location of the puncture to ensure footplate is not in bifurcation or any side branch. Reposition the device to stop blood marking (maintain the 30 degree rotation) or reinsert the wire and evaluate the situation before continuing the procedure.

6. While maintaining the device position, stabilize the device with your free hand (the one not used to deploy the device) to maintain the gentle retraction and to ensure the device doesn't twist or move forward during deployment. Use your other hand to deploy needles by pushing on the plunger assembly **(in the direction marked #2)** until you visually confirm that the collar of the plunger makes contact with the proximal end of the body. Do not use excessive force or repeatedly push the plunger assembly. After visually confirming contact with the body of the device only **one time**, this step is complete.
7. Using your thumb as a fulcrum on the handle, gently disengage the needles by pulling the plunger assembly back **(in the direction marked #3)** and completely remove the plunger and needles from the body of the device. One suture limb will be attached to the anterior needle. The posterior needle will be free of suture. Pull back on the plunger until the suture is taut, which confirms that the suture has been fully retracted from the body of the device.
8. Do not attempt to redeploy the needles if the suture limb is not attached to the anterior needle. Reinsert the guide wire, and remove the device over the wire. Insert a new Perclose ProGlide SMC device to complete the procedure.
9. Use the QuickCut suture-trimming mechanism located on the handle to cut the suture from the anterior needle distal of the link. Use of a new, sterile scalpel or scissors is optional.
10. Relax the device and then return the foot to its original position by pushing the lever **(marked #4)** on top of the device, down to its original position. Do not attempt to remove the device without closing the lever.
11. Withdraw the Perclose ProGlide SMC device until the guide wire port exits the skin line.
12. Slightly rotate the device until you can see the two suture limbs in the bend of the distal guide. Grasp the sutures adjacent to the sheath. While holding the two suture limbs **together**, gently pull both suture ends through the distal end of the proximal guide.
13. **Immediately** place a shodded hemostat or clamp to hold the two suture limbs together at the distal end of the non-rail suture limb (shortest limb with white tip). To prevent knot advancement or locking of the knot, care must be taken not to pull on the individual suture until the clamp is securely holding the two suture limbs together.
14. Gently pull on the clamp until the suture is taut to **remove any suture slack from the tissue tract.** Place the clamped suture on the **right side** of the patient under a sterile towel. It is important to identify which suture is deployed first as this is the knot that needs to be advanced first at the end of the procedure. **NOTE:** The monofilament suture can be damaged by opening and closing the clamp. If you want to attach the suture to the drape, it is recommended that you use a second clamp with the tip placed through the handle of the first clamp and attach the second clamp to the drape. It is important to remember this suture was placed first when suture tying occurs at the end of the procedure.
15. Reinsert the guide wire. There should be sufficient guide wire inside the vessel and exposed out of the guide wire exit port for device exchange.
16. Remove the Perclose ProGlide SMC device, while holding compression above the puncture site and maintaining adequate length of guide wire inside the vessel. This allows placement of the second Perclose ProGlide SMC device.
17. Repeat Steps 2–13 with the second Perclose ProGlide device. **NOTE:** In Step 4, the second device should be rotated approximately 30 degrees towards the patient's left side.
18. After removing any excess suture slack from the tissue track, place the clamped suture for the second device on the left side of the patient under a sterile towel. It is important to identify which suture was placed first and which suture was placed second. At the completion of the procedure, the suture knots will be advanced in the order they were placed. The knot from the first device placed on the right side of the patient would be advanced, followed by the knot from the second device placed on the left side of the patient.
19. **Knot advancement will be placed on hold at this point and the sutures set aside under sterile towels until the operator proceeds with the procedure.**
20. When utilizing a pre-close technique, the Perclose ProGlide SMC device is exchanged for an appropriately sized introducer sheath.
21. After completing the procedure, a hydrophilic or general purpose guide wire should be advanced into the vessel. Maintain adequate length of guide wire in both the vessel and exposed out the guide wire exit port to ensure guide wire access is maintained until hemostasis is achieved.
22. Heavily irrigate the secured Perclose ProGlide sutures with heparinized saline to remove any dry blood.
23. Remove the clamp from the first suture (patient's right side). The rail suture limb is blue and is the longer of the two suture limbs. This rail suture limb will be used to advance the knot. The shorter, non-rail suture limb is white tipped and will be used to lock the knot.
24. Wrap the rail suture limb securely around your left forefinger, close to skin level, and hold the suture coaxial to the tissue tract. **While maintaining guide wire access,** carefully remove the entire sheath system and simultaneously pull the rail limb with slow, consistent increasing tension. Avoid quick or jerking type movements with the suture limbs. Manual pressure should be applied proximal to the puncture site for hemostasis, while the sheath is removed and during initial suture advancement.
25. **DO NOT lock the knot at this point.** Due to the size of the access site, use of Suture Trimmer (Section 10.3, Step 13) may be needed to approximate the tissue edges. However, **DO NOT lock or excessively tighten the knot while the guide wire is still in the vessel.** Again place the suture limbs on the right side of the patient for easy identification as the first suture deployed.
26. Remove the clamp from the second suture (patient's left side) and advance the knot using the same technique and **maintaining guide wire access. Again, DO NOT lock the knot.** Place the suture limbs on the left side of the patient for easy identification as the second suture placed.
27. Assess for hemostasis. If brisk bleeding is observed, advance the first (patient's right side) suture again and then advance the second (patient's left side) suture again. Multiple knot advancements are common when closing larger sheath sizes. However, **DO NOT lock or excessively tighten the knot while the guide wire is still in the vessel. Until the wire is removed, some bleeding will be visible, but it should not be brisk.**
28. If acceptable hemostasis is not observed, additional Perclose ProGlide SMC devices may be deployed at this point. Repeat Steps 2–11 with a third Perclose ProGlide device. **NOTE:** In Step 4, the third device should not be rotated. The device will be deployed with the logo facing the ceiling / 12 o'clock. After deployment of this device, advance the knot in the same fashion. **DO NOT lock the knot or excessively tighten the knot while the guide wire remains in the vessel.**
29. Assess the site for adequate hemostasis. **If bleeding is controlled, the operator should then remove the guide wire.** With the rail suture limb (longer, blue limb) securely wrapped around the left forefinger, again advance the first suture (patient's right side) and then place the Suture Trimmer under the left thumb to assume a single-handed position and complete knot advancement with slow, consistent increasing tension until the suture is taut (guitar string tightness). With the Suture Trimmer in place and the suture taut, tighten the knot by gently pulling the non-rail (shorter, white tipped) suture limb keeping it coaxial to the tissue tract.
30. **DO NOT cut the suture.** Follow the same steps to advance the second suture (patient's left side), coaxial to the tissue tract and lock the knot but **DO NOT cut the suture.** If applicable, advance the third sutures and lock the knots in the order that they were placed (first, second, third). Assess bleeding. If hemostasis is deemed adequate, cut the suture tails below the surface of the skin using the Suture Trimmer (Section 10.3, Step 13 J) or a new, sterile scalpel or scissors.

### 10.5 Suture Breakage

1. If suture breakage occurs **before the knot** is tightened, if a wire is still in place, discard the suture material, and use another Perclose ProGlide SMC device to complete the procedure.
2. If suture breakage occurs **after a knot** has been advanced and / or tightened, and a wire is still in place, another Perclose ProGlide SMC device can be used to complete the procedure or a sheath can be inserted. Care should be taken to avoid excessive force if the reintroduction of another device or introducer sheath is required. To avoid resistance, use an introducer sheath small enough to be introduced without undue force.
3. In all cases, if another Perclose ProGlide SMC device or introducer sheath cannot be inserted, use manual compression to obtain hemostasis.
4. To prevent suture breakage, always pull on the suture limbs with slow, consistent increasing tension. Avoid quick or jerking type movements with the suture limbs.
5. To prevent damage to the suture and subsequent suture breaks, the Suture Trimmer and suture limbs should always remain coaxial to the tissue tract. The Thumb Knob should be at 12 o'clock (facing the ceiling) and the Suture Trimmer should not be rotated. When loading the suture into the Suture Trimmer, keep the Thumb Knob retracted until the suture and Suture Trimmer are coaxial, then release the Thumb Knob to capture the suture in the Suture Gate.

### 10.6 Post-Procedure Patient Management

1. Apply an appropriate dressing to the access site.
2. Assess the access site as per hospital standard of care.

### 10.7 Recommendation for Patient Ambulation and Discharge

Patients who have undergone a diagnostic or interventional procedure using 5F–8F sheaths may be ambulated two hours after the Perclose ProGlide SMC device procedures.

For patients who have undergone an interventional catheterization procedure using sheaths greater than 8F, post-procedure ambulation and patient care is at the discretion of the physician.

In determining whether to ambulate or discharge an individual patient, it is important to consider all clinical factors including, but not limited to, anticoagulation regimen, antiplatelet and thrombolytic agents administered, oozing or bleeding from the access site, venous access site hemostasis, the general cardiovascular condition of the patient, anesthetic levels, and the overall clinical condition of the patient.

## 11.0 PRODUCT INFORMATION DISCLOSURE

Abbott Vascular Inc. has exercised reasonable care in the manufacturing of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability of fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc. directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

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**12.0 PATENTS AND TRADEMARKS**

This product and / or its use may be covered by one or more of the following United States Patents: 6,136,010; 6,190,396; 6,358,258; 6,558,399; 6,730,102; 6,746,457; 6,964,668; 7,001,400; 7,029,487; 7,048,747; 7,094,246; 7,147,646 ; 7,201,762; 7,235,087; EP 721313; EP 1093346. Other U.S. patents pending. Foreign patents issued and pending.

QuickCut, Perclose, and ProGlide are trademarks of the Abbott Group of Companies.

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
### CAUTION

This device is restricted to sale by or on the order of physicians (or other healthcare professionals, authorized by, or under the direction of, such physicians) who are trained in diagnostic and therapeutic catheterization procedures and who have been trained by an authorized representative of Abbott Vascular. Prior to use, the operators must review the instructions for use and be familiar with the deployment techniques associated with the use of this device.

### DEVICE DESCRIPTION


The StarClose SE Vascular Closure System is designed to deliver a nitinol (nickel and titanium) clip to close femoral artery puncture sites following percutaneous catheterization procedures. The StarClose SE Vascular Closure System consists of the StarClose SE Clip Applier (see **Figure 1**) and a 6F Exchange System (see **Figure 2**). An implantable clip is mounted on the Clip Applier, which delivers the clip through the Exchange Sheath for extravascular closure of puncture sites.

**Figure 1 - StarClose SE Clip Applier**

1. Locator Wings
2. Flex-guide
3. Clip Delivery Tube
4. Sheath Splitter
5. Side-arm Slot
6. Stabilizer
7. Deployment Button
8. Thumb Advancer
9. Safety Release
10. Plunger
11. Access Ports

**Figure 2 - 6F Exchange System**

A. Exchange Sheath
   1. Exchange Hub with Hemostasis Valve
   2. Sheath Tube
B. Dilator
   3. Dilator Hub
C. J Tip Guide Wire

### INDICATIONS

The StarClose SE Vascular Closure System is indicated for the percutaneous delivery of an extravascular clip for closure of femoral artery access sites following catheter-based procedures.

### CONTRAINDICATIONS

The StarClose SE Vascular Closure System is contraindicated for use in patients with known hypersensitivity to nickel-titanium.

### WARNINGS

Do not use the StarClose SE Vascular Closure System if the packaging or sterile barrier has been previously opened or damaged or if the components appear to be damaged or defective.

DO NOT RESTERILIZE OR REUSE the StarClose SE Vascular Closure System if the sterile field has been broken where bacterial contamination of the sheath or surrounding tissues may have occurred, since such a broken sterile field may result in infection.

Do not use the StarClose SE Vascular Closure System if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma. Perform a femoral angiogram to verify the location of the puncture site.

Do not use the StarClose SE Vascular Closure System if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal hematoma.

Do not use the StarClose SE Vascular Closure System if the puncture site is located in the superficial femoral artery or the profunda femoris artery, since such puncture sites may result in a pseudoaneurysm, intimal dissection, or an acute vessel closure (thrombosis of small artery lumen). Perform a femoral angiogram to verify the location of the puncture site.

### PRECAUTIONS

1. The safety and effectiveness of the StarClose SE Vascular Closure System have not been established in the following patient populations:
   - Patients with introducer sheaths < 5F or > 6F during the catheterization procedure
   - Patients with ipsilateral arterial access sites punctured and compressed within 3 months before the catheterization procedure
   - Patients having arterial access in vessels other than the common femoral artery
   - Patients having a hematoma, pseudoaneurysm, or arteriovenous fistula present prior to sheath removal
   - Patients with femoral artery calcium, which is fluoroscopically visible at access site
   - Patients with small femoral arteries (< 5 mm in diameter)
   - Patients with severe claudication, iliac or femoral artery diameter stenosis greater than 50%, or previous bypass surgery or stent placement in the vicinity of access site
   - Patients with access sites in vascular grafts
   - Patients with prior intra-aortic balloon pump at access site at any time prior
   - Patients with ipsilateral femoral venous sheath during catheterization procedure
   - Patients with which there is difficulty inserting the introducer sheath or greater than one ipsilateral arterial puncture at the start of the catheterization procedure
   - Patients with intra-procedural bleeding around access site.
   - Patients younger than 18 years of age
   - Patients who are pregnant or lactating
   - Patients with bleeding diathesis or coagulopathy
   - Patients who are morbidly obese (Body Mass Index > 35 kg/m$^2$)
   - Patients with hypertension (systolic BP > 180 mmHg or diastolic BP > 110 mmHg) unresponsive to medical therapy
   - Patients with active systemic or cutaneous infection or inflammation
   - Patients with access sites above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks
2. The StarClose SE Vascular Closure System is provided sterile and non-pyrogenic in unopened undamaged packaging. Products are sterilized with ethylene oxide and intended for single use only. Do not resterilize. Store in a cool, dry place.

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3. Prior to use, ensure that the StarClose SE Vascular Closure System and sterile packaging have not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
4. As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the StarClose SE Vascular Closure System. Employ appropriate groin management, as per hospital protocol, post-procedure and post-hospital discharge to prevent infection.
5. Use a single wall puncture technique. Do not puncture the posterior wall of the artery.
6. Do not use the StarClose SE Vascular Closure System to close vessels with diameters less than 5 mm.
7. Do not deploy the Clip in areas of calcified plaque.
8. The StarClose SE Vascular Closure System can **ONLY** be used with the StarClose Exchange System (included in the StarClose SE Vascular Closure System packaging).
9. **Do not advance or withdraw the StarClose SE Vascular Closure Device against resistance until the cause of that resistance has been determined.** Excessive force used to advance or torque the StarClose SE device should be avoided, as this may lead to significant vessel damage and / or breakage of the device, which may necessitate interventional and / or surgical removal of the device and vessel repair.

## ADDITIONAL CONSIDERATIONS

The StarClose Clip has been shown to be MR Conditional immediately following implantation. A patient with this implant can be scanned safely immediately after clip placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.


In non-clinical testing, the StarClose Clip produced a temperature rise of 0.5°C at maximum MR system-reported whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 Tesla MR system using a transmit / receive body coil.

The MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the StarClose Clip. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary.

## CLOSURE PROCEDURE

1. It is necessary to create a 5 – 7 mm skin incision at the sheath site to accommodate the insertion of the Clip Delivery Tube into the tissue tract. This should be done at the beginning of the procedure prior to the administration of anticoagulants and antiplatelet agents, if possible. Blunt dissection by single spread with surgical instrument in the skin incision is recommended, especially in those patients with scar tissue from previous catheterization procedures.
2. Perform a femoral angiogram through the side port of the procedural sheath to determine the location of the arteriotomy site, the vessel size, and the presence of disease (calcified plaque, stenosis, chronic or acute occlusion), tortuosity or presence of arterial wall dissection.
3. Re-prep the access site, placing clean towels around the access site, and wear new sterile gloves prior to handling the Clip Applier and proceed with the closure procedure.
4. Prepare the access site for closure:
   a. Insert a guide wire and exchange for a StarClose Exchange Sheath.
5. **CLICK 1:** Connect the StarClose Exchange Sheath to the Clip Applier **(Figure 3)**.

**Figure 3**


   a. Approximate the Clip Applier to the Hub of the StarClose Exchange Sheath. Hold the bottom of the Sheath Hub with the left hand and insert the distal tip of the Flex-guide into the Hub using the right hand.
   b. Slowly advance the Flex-guide down through the Sheath. Make sure not to severely bend the Flex-guide during insertion and advancement.

**Note:** When connecting the Sheath Hub to the Clip Applier, hold the Hub up at the same angle as the tissue tract above the skin surface. This action allows a firm connection to be made without pushing against the skin.

   c. An audible "click" should be heard. Check to make sure the Hub-to-Clip Applier engagement is secure by gently pulling on the Sheath Hub.
6. **CLICK 2:** Deploy Locator Wings and initiate splitting of Sheath **(Figure 4)**.

**Figure 4**

   a. Place the left hand on the Stabilizer to stabilize the device at the angle of the tissue tract. Orient the Clip Applier so the front side faces upward (the side with the product logo).
   b. Assume a syringe grip on the proximal end of the device with the right hand by placing the index finger and the middle finger on the proximal posts, and the thumb on the Plunger.
   c. Retract the device out of the tissue 3 to 4 cm with the right hand. Firmly depress the Plunger with the right thumb until the number "2" appears in the number window. This opens and locks the Locator Wings in the expanded position, and advances the Thumb Advancer and the Clip Delivery Tube approximately 3 cm to initiate splitting of the Sheath. The initial splitting of the Sheath can be seen above the skin level.

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
d. An audible "click" should be heard. Check to make sure the number "2" is completely visible in the number window. This indicates that the Locator Wings are opened and in the locked position.

7. **CLICK 3:** Advance the Thumb Advancer and deliver the Clip to the access site (**Figure 5**).

**Figure 5**

a. Maintain the left hand on the Stabilizer to stabilize the device at the angle of the tissue tract, gently retract the device with the right hand until slight resistance is felt. The goal is to place the Locator Wings against the anterior surface of the arterial wall to locate the access site without applying tension on the artery.

b. Assume a syringe grip on the distal end of the device with the right hand by placing the index finger on the distal post, the middle two fingers in the loop of the Stabilizer, and the thumb on the Thumb Advancer.


c. While maintaining apposition of the Locator Wings on the anterior surface of the arterial wall and keeping the Flex-guide straight, advance the Thumb Advancer using the pad of the right thumb until the number "3" appears in the number window. This advances the Clip Delivery Tube over the Flex-guide while splitting the Sheath from the Hub to the distal tip.

d. An audible "click" should be heard. Check to make sure the number "3" is completely visible in the number window. This indicates that the Thumb Advancer is fully advanced and in the locked position.

**Note:** Do not advance the Thumb Advancer against excessive resistance.

If excessive resistance is experienced during advancement of the Thumb Advancer, manually collapse the Locator Wings and remove the device following the steps below (**Figure 6**):

**Figure 6**


- Retract the Thumb Advancer as far proximal as possible until resistance is felt. This is to decrease any interaction between the Clip Delivery Tube and the Flex-guide.
- Slide the Safety Release in the direction of the arrow as seen in **Figure 6** and on the device to collapse the Locator Wings. The Locator Wings are fully collapsed when the number "2" on the Plunger exits the number window completely.
- Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers.
- Provide counter-traction with the left hand on the patient's body, and assertively pull the device out with the right hand.
- It is important not to rock or twist the device as this may cause arterial damage.

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8. **CLICK 4:** Deploy Clip to close the arteriotomy **(Figure 7)**.

**Figure 7**

a. While stabilizing the body of the device with the right hand and supporting the Clip Delivery Tube with the left hand, raise the body of the device to a 60 to 75 degree angle.


b. Support the Clip Delivery Tube with the left hand and gently push the device down on top of the artery to seat the Clip Delivery Tube on top of the access site. Pulsation of the artery may be felt.

c. While maintaining downward pressure and device stability, depress the Deployment Button with the right thumb. An audible "click" should be heard. This releases the Clip to close the arteriotomy and collapses the Locator Wings simultaneously.

d. Maintain the downward pressure for 2 – 3 seconds. Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers. Provide counter-traction with the left hand on the patient's body, and remove the device with the right hand.

**Note:** If resistance is met upon removal of the device, follow the steps below to remove the device **(Figure 8)**:

**Figure 8**

- Locate the two Access Ports on the lateral and medial sides of the device approximate to the location of the product logo.
- Insert the distal end of the Dilator (or an 18-Gauge Thin-wall Needle) into the Access Ports, one at a time, to release the locking mechanism on the Thumb Advancer.
- An audible "click" should be heard. Check that the number "3" exits the number window completely and the Thumb Advancer is freed from the locked position. Repeat the above steps if necessary.
- Retract the Thumb Advancer as far proximal as possible until resistance is felt. This is to decrease any interaction between the Clip Delivery Tube and the Locator Wings.
- Slide the Safety Release to collapse the Locator Wings and reset the Plunger. The Locator Wings are fully collapsed when the number "2" on the Plunger exits the number window completely.
- Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers.
- Provide counter-traction with the left hand on the patient's body, and assertively pull the device out with the right hand.
- It is important not to rock or twist the device as this may cause arterial damage.

9. Check for hemostasis.

a. Use sterile gauze to firmly depress the area around the puncture site to express any blood in the tissue tract.

b. While supporting the patient's leg, have him / her bend his / her knee toward the chest, and have the patient lift his / her head off the table and cough.

c. Adjunctive compression can be applied for up to 5 minutes, as needed.

d. If hemostasis is not achieved, apply manual compression until hemostasis is achieved.

**POST-PROCEDURE PATIENT MANAGEMENT**

- Apply an appropriate dressing to the puncture site.
- Assess the insertion site as per hospital protocol.

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In determining whether to ambulate or discharge an individual patient, it is important to consider all clinical factors including, but not limited to, anticoagulation regimen, antiplatelet and thrombolytic agents administered, oozing or bleeding from the access site, venous access site hemostasis, the general cardiovascular condition of the patient, anesthetic and medication levels, and the overall clinical condition of the patient.

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## PRODUCT INFORMATION DISCLOSURE

Abbott Vascular Inc. has exercised reasonable care in the manufacture of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc., directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

## HOW SUPPLIED

The StarClose SE Vascular Closure System and accessories are provided sterile and non-pyrogenic in unopened, undamaged packaging. Products are sterilized with ethylene oxide and intended for single use only. This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. The device and the primary packaging do not contain latex. Store in a cool, dry place.

**StarClose SE Vascular Closure System** contains:

One (1) StarClose SE Clip Applier
One (1) 6F Exchange Sheath (11 cm)
One (1) 6F Dilator
One (1) 0.038" (0.97 mm) J Tip Guide Wire

## PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 5,674,231; 5,810,846; 6,197,042; 6,277,140; 6,391,048; 6,461,364; 6,623,510; 6,626,918; 6,632,238; 6,695,867; 6,719,777; 6,749,621; 6,780,197; 6,942,674; 7,001,398; 7,144,411. Other patents pending.

StarClose SE is a trademark of the Abbott Group of Companies.

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# Prostar XL

## Percutaneous Vascular Surgical System

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## ENGLISH

**PROSTAR XL**
**PERCUTANEOUS VASCULAR SURGICAL (PVS) SYSTEM**
**INSTRUCTIONS FOR USE**

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### CAUTION

This device should only be used by physicians (or allied healthcare professionals, authorized by, or under the direction of, such physicians) who are trained in diagnostic and therapeutic catheterization procedures and who have been trained by an authorized representative of Abbott Vascular.

Prior to use, the operator must review the Instructions for Use and be familiar with the deployment techniques associated with the use of this device.

### DEVICE DESCRIPTION

The Prostar XL Percutaneous Vascular Surgical (Prostar XL PVS) device is designed to deliver polyester suture(s) to close femoral artery puncture sites following catheterization procedures. The Prostar XL device has two sutures and four needles.

The Prostar XL device consists of a sheath, which contains two pairs of sutured needles, a needle guide that precisely controls the placement of needles around the puncture site and a rotating barrel that receives deployed needles. The Prostar XL sheaths have J-shaped tips at the distal end.


A marker lumen is contained within the barrel of the device with the intra-arterial port of the lumen positioned in the needle guide.

Proximally, the marker lumen exits from the hub of the device. The marker lumen allows a pathway for back-bleeding from the femoral artery and ensures proper device positioning. The barrel rotates independently from the central core and is designed to prepare the subcutaneous track. Barrel rotation is accomplished by depressing the interlocks exiting from the hub. The Prostar XL devices track over a standard 0.038" (0.97 mm)(or smaller) guide wire.

The Prostar XL PVS System consists of a 10F Prostar XL device and a Perclose Knot Pusher. The Perclose Knot Pusher is designed to advance the tied suture knot to the arteriotomy.

The Prostar XL device and the Perclose Knot Pusher are depicted in **Figure 1.**

**Figure 1**

Prostar XL Percutaneous Vascular Surgical System

A. Prostar device
B. Perclose Knot Pusher

### INDICATIONS FOR USE

The Prostar XL Percutaneous Vascular Surgical Systems are intended for percutaneous delivery of sutures for closing the common femoral artery access site, reducing the time to hemostasis and ambulation (patient walks ten feet/three meters) for patients who have undergone catheterization procedures.

The Prostar XL PVS System reduces the time to discharge in patients having undergone diagnostic catheterization procedures without complicating clinical conditions (refer to PRECAUTIONS and SPECIAL PATIENT POPULATIONS).

The Prostar XL PVS System is designed for use in conjunction with 8.5 to 24F sheaths.

### CONTRAINDICATIONS

There are no known contraindications to the use of this device. Attention is drawn to the Warnings, Precautions, and Special Patient Populations.

### WARNINGS

Do not use the Prostar XL PVS device or accessories if the packaging or sterile barrier has been previously opened or damaged, or if the components appear to be damaged or defective.

DO NOT RESTERILIZE OR REUSE. The Prostar XL PVS device and accessories are intended for single use only.

Do not use the Prostar XL PVS System if the sterile field has been broken where bacterial contamination of the sheath or surrounding tissues may have occurred, since such a broken sterile field may result in infection.

Do not use the Prostar XL PVS System if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma. Perform a femoral angiogram to verify the location of the puncture site.

Do not use the Prostar XL PVS System if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal hematoma.

Do not use the Prostar XL PVS System if the puncture site is located in the superficial femoral artery or the profunda femoris artery, since such puncture sites may result in a pseudoaneurysm, intimal dissection, or an acute vessel closure (thrombosis of small artery lumen). Perform a femoral angiogram to verify the location of the puncture site.

### PRECAUTIONS

1. The Prostar XL PVS System is provided sterile and non-pyrogenic in unopened undamaged packages. Products are sterilized with ethylene oxide and intended for single use only. Do not resterilize. Store in a cool, dry place.
2. Prior to use, inspect the Prostar XL PVS System to ensure that the sterile packaging has not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
3. As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the Prostar XL PVS System. Employ appropriate groin management, as per hospital protocol, post procedure and post hospital discharge to prevent infection.
4. Use a single wall puncture technique. Do not puncture the posterior wall of the artery. Avoid posterior wall suture placement.
5. Do not insert the Prostar XL PVS device into the femoral artery at an angle greater than 45 degrees to the longitudinal plane of the artery.
6. An experienced operator should use accepted surgical technique as warranted by surgical circumstances to provide adequate knot security.
7. There are no reaccess restrictions if previous arteriotomy repairs were achieved with an Abbott Vascular Suture Mediated Device.
8. **Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined** (see Prostar XL Device Placement section). **Excessive force used to advance or torque the Prostar XL PVS device should be avoided, as this may lead to significant vessel**

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**damage and / or breakage of the device, which may necessitate intervention and / or surgical removal of the device and vessel repair.**

9. If excessive resistance in advancing the Prostar XL PVS device is encountered, withdraw the Prostar XL device over a 0.038" (0.97 mm) (or smaller) guide wire and reinsert the introducer sheath, or use conventional compression therapy.
10. In the event suture breakage occurs after an initial knot has been tied, care should be taken to avoid excessive force if the reintroduction of another Prostar XL device or introducer sheath is required. Any resistance to introduction should result in advancement of an introducer sheath small enough to be introduced without undue force.
11. If significant blood flow is present around the barrel of the Prostar XL device, do not deploy the needles. Remove the Prostar XL device over a 0.038" (0.97 mm) (or smaller) guide wire and insert an appropriately sized introducer sheath.
12. Remove the Prostar XL sheath before tightening the suture. Failure to remove the sheath prior to tightening the suture may result in detachment of the tip of the sheath.
13. In using this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing damage due to application of surgical instruments such as forceps or needle holders.
14. Do not attempt to re-deploy Prostar XL needles after the needles have been "backed-down" into the sheath (refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section).
15. If using the "pre-close" technique (placement of the Prostar XL device sutures prior to dilating the access site beyond 10F), sutures should be gradually tightened as the introducer sheath is removed to maintain hemostasis.
16. Use conventional compression methods in the event bleeding from the femoral access site persists, after the use of the Prostar XL PVS device.
17. The Prostar XL device should be advanced using fluoroscopic guidance to prevent possible vessel damage.

## SPECIAL PATIENT POPULATIONS

The safety and effectiveness of the Prostar XL PVS System have not been established in the following patient populations:

- Patients with introducer sheaths < 8.5F or > 24F during the catheterization procedure.
- Patients with ipsilateral arterial access sites punctured and compressed within 48 hours of closure.
- Patients having arterial access in vessels other than the common femoral artery.
- Patients having a hematoma, pseudoaneurysm or arteriovenous fistula present prior to sheath removal.
- Patients with femoral artery calcium, which is fluoroscopically visible at access site.
- Patients with small femoral arteries (< 5 mm in diameter).
- Patients with severe claudication, iliac or femoral artery diameter stenosis greater than 50% or previous bypass surgery or stent placement in the vicinity of access site.
- Patients with access sites in vascular grafts.
- Patients with prior intra-aortic balloon pump at access site at any time prior.
- Patients with ipsilateral femoral venous sheath during the catheterization procedure.
- Patients with whom there is difficulty inserting the introducer sheath or greater than one ipsilateral arterial puncture at the start of the catheterization procedure.
- Patients with intra-procedural bleeding around access site.
- Patients receiving glycoprotein IIb/IIIa inhibitors before, during, or after the catheterization procedure.
- Patients younger than 18 years of age.
- Patients who are pregnant or lactating.
- Patients with bleeding diathesis or coagulopathy.
- Patients who are morbidly obese where less than one third of the access needle is above the skin line.
- Patients with active systemic or cutaneous infection or inflammation.
- Patients with access sites above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks.
- Patients with antegrade punctures.

## CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of Prostar XL PVS System. The techniques and procedures described below are not intended as a substitute for the operator's experience and judgment in treating any specific patients.

### Examination and Selection of Products

1. The outer pouch packaging of the Prostar XL PVS System provides the sterile barrier. After careful inspection of the packaging to ensure there is no damage to the sterile barrier, remove the device from the package.
2. Exercise care during device handling to reduce the possibility of accidental device breakage.
3. Exercise care when using additional surgical instruments such as forceps, hemostats or needle holders during device handling, to reduce the possibility of accidental device breakage.
4. Verify marker port patency by flushing the lumen with saline until the saline exits from the marker port. **Do not use the Prostar XL PVS device if the marker lumen is not patent.**

### Arterial Puncture Considerations

1. Puncture the anterior wall of the common femoral artery at an angle of approximately 45 degrees.
2. Avoid sidewall or posterior wall femoral artery punctures.
3. Puncture locations are ideally located in the common femoral artery below the level of the inguinal ligament and above the common femoral artery bifurcation.

### Prostar XL Device Placement

**The following instructions detail the deployment sequence to close the access site of a catheterization procedure performed through an applicable sized sheath.**

When utilizing a "pre-close" technique (placement of the Prostar XL device sutures prior to dilating the access site beyond 10F) refer to step 7 in the Suture Management section for more information.

1. Prior to Prostar XL device placement, perform a femoral angiogram through the introducer sheath to verify that the access site is in the common femoral artery.
2. Evaluate the femoral artery site for size, calcium deposits, and tortuosity, to avoid posterior wall suture placement and possible ligation of the anterior and posterior walls of the femoral artery.
3. Re-prep the access site, placing clean towels around the access site, and wearing new sterile gloves prior to handling the device and proceeding with the closure procedure.
4. While the introducer sheath remains in place, use a scalpel to slightly extend the incision and forceps to dilate subcutaneous tissue.

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
5. Place a 0.038" (0.97 mm) (or smaller) guide wire through the introducer sheath. Remove the introducer sheath while applying pressure on the groin to maintain hemostasis.
6. Carefully back-load the Prostar XL device over the guide wire until the guide wire exit port is just above the skin line. Remove the guide wire. Continue to advance the Prostar XL device until the barrel is at skin level.
7. Unlock the hub by depressing the interlocks with thumb and forefinger. Once the hub is unlocked, rotate the hub while gently advancing the barrel at a **45-degree angle**, or less.
8. A steady, continuous drip of blood from the dedicated marker lumen occurs when the Prostar XL device is properly positioned (see **Figure 2**). Marking from the lumen(s) containing suture(s) may occur but should not be used as the indicator for proper positioning and needle deployment.

   **Do not clamp the suture lumen with a hemostat or other instrument. Doing so will prevent suture deployment.**

**Figure 2**


9. Lock the hub back in place.
   - 9.1. If a continuous drip of blood (luminal marking) from the dedicated marker lumen is not apparent, withdraw the device to expose the marker port. Flush the marker lumen to verify patency, and then continue to gently advance the Prostar XL device while rotating the barrel.
   - 9.2. If continuous marking is still not achieved, remove the Prostar XL device and follow conventional compression protocol, or replace the Prostar XL device with an appropriately sized introducer sheath.
   - 9.3. Do not deploy needles until a continuous drip of blood is evident from the dedicated marker lumen.

### Needle Deployment


1. Confirm that the interlocks are re-engaged and correctly aligned with the locking indents in the hub (locked).
2. With the left hand, hold the hub of the device in position at an **angle of 45 degrees** (or less).
3. With the right hand, rotate the handle counterclockwise, to unlock the handle.
4. Ensure blood marking is maintained.
5. Pull the handle away from the hub to deploy the needles (see **Figure 3**).

**Figure 3**

6. If resistance to rotating the handle is encountered, **do not attempt to deploy the needles.** Significant resistance is an indication that the hub is not properly positioned. Refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section to ensure proper hub positioning.
7. **Do not clamp the suture lumen with a hemostat or other instrument. Doing so will prevent suture deployment.**
8. Continue to pull the handle until the needle tips emerge at the top of the barrel.
9. While steadily holding the device in position, confirm that all four needles are visible in the hub.
10. If significant resistance is encountered prior to needle tips emerging at the top of the barrel, or if all four of the needles are not deployed, terminate deployment. Refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section to perform needle back-down procedure.
11. Using a hemostat, remove the posterior needle(s) followed by the anterior needle(s) using the funnel-shaped hub as a fulcrum to facilitate needle removal (see **Figure 4**).


**Figure 4**

### Suture Management

1. Once needles have been removed from the hub, remove slack in the suture(s) by pulling the suture ends to evenly matched lengths and tensioning until resistance is felt. Cut the suture ends close to the needles. Dispose of needles in accordance with hospital policy.
2. Withdraw the Prostar XL device while maintaining access to the sutures exiting the hub during this step. Create a "bow string" effect with the exposed suture(s) by bending the Prostar XL device sheath away from the operator and applying tension to the suture ends exiting the hub (see **Figure 5**). This ensures that the sutures are not tied around the sheath.

**Figure 5**

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3. Grab the anterior suture ends (one white and one green) adjacent to the topside of the sheath and pull the suture ends through the distal end of the barrel. Place this pair of suture ends toward the patient's head.
4. Grab the posterior suture ends (one white and one green) adjacent to the underside of the sheath and pull the suture ends through the distal end of the barrel. Ensure that this pair of suture ends comes under the device and towards you. Place this pair of suture ends toward the patient's feet.
5. Identify the two ends of a single suture by color, and then tension the ends using a gentle seesawing motion.
6. Continue to withdraw the Prostar XL device until the guide wire exit port exits the skin line. Reinsert a 0.038" (0.97 mm) (or smaller) guide wire into the visible guide wire exit port so that guide wire access can be maintained until hemostasis is verified. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

Access sites larger than 10F require the use of the "pre-close" technique.

7. When utilizing a "pre-close" technique, the Prostar XL device is exchanged for an appropriately sized introducer sheath or an additional Prostar XL device if additional sutures are to be placed around the same access site.
   - Sliding, self-locking surgical knots should be tied prior to introducer sheath removal according to the Knot Advancement section below.

Access sites larger than 10F may require the use of an additional Prostar XL device.

8. Use of an additional Prostar XL device:
   - 8.1 If additional sutures are to be placed, position an additional Prostar XL device rotated 45 degrees relative to the first Prostar XL device to establish alternating suture positions and repeat all the steps above beginning at **Prostar XL Device Placement**.
   - 8.2 Identify and secure the sutures from the first device separately from the sutures of the additional device.

**Knot advancement**


1. Identify the green suture ends. Place one green suture end in each hand.
2. Gently tension the suture to create a longer suture end. (The short end will be your rail limb of suture, for the green knot).
3. Tie a sliding, self-locking surgical knot using the green suture.
4. Lay the green suture aside.
5. Identify the white suture ends. Place one white suture end in each hand.
6. Gently tension the suture to create a longer suture end. The short end will be your rail limb of suture, for the white knot.
7. Tie a sliding, self-locking surgical knot using the white suture.
8. Securely wrap the white rail suture around your left index finger. Saturate the suture with saline. Gently pull on the white rail suture first, keeping the suture coaxial to the tissue tract. Completely remove the Prostar XL device from the artery, leaving the guide wire in the artery.

   **Note:** If using the "pre-close" technique, sutures should be gradually tightened as the introducer sheath is removed to maintain hemostasis.

   - 8.1 Do not compress the femoral access site while withdrawing the Prostar XL device from the tissue tract. Applying excessive force to the Prostar XL sheath, due to site compression, during removal of the device may result in a needle guide break. Please refer to **NEEDLE GUIDE BREAK PROCEDURE** section for further instruction.
   - 8.2 Do not tighten the suture around the Prostar XL sheath. Tightening the suture around the Prostar XL sheath during device removal may result in a needle guide break. Please refer to **NEEDLE GUIDE BREAK PROCEDURE** section for further instructions.

9. Saturate the sutures with saline. Load the rail suture into the Knot Pusher and continue to advance the knot forward to the arteriotomy (see **Figure 6**).

**Figure 6**

10. With the rail suture limb securely wrapped around the left forefinger, place the Knot Pusher under the left thumb to assume a single-handed position to complete knot advancement.
11. Do not tension the non-rail suture to tighten the knot.
12. Remove the Knot Pusher from the tissue tract and from the suture without using the thumb knob.
13. Lay the white suture down, being sure to separate the rail and non-rail ends.
14. Repeat steps 8–13 with the green suture.
15. Repeat steps 1–14 for second set of sutures if additional device is used.
16. If hemostasis has been achieved, remove the guide wire from the patient. Continue with step 18.
17. If hemostasis has not been achieved:
    - 17.1. Repeat steps 9–14.
    - 17.2. Exchange the Prostar XL device with an appropriately sized introducer sheath. Care should be taken to avoid excessive force if the reintroduction of another introducer sheath is required. To avoid resistance, use an introducer sheath small enough to be introduced without undue force, but large enough to maintain hemostasis.
    - 17.3. Apply conventional compression therapy.
18. Once again, securely wrap the white rail suture around your left index finger.
19. Load the white rail suture into the Knot Pusher and advance to the level of the arteriotomy.
20. With the rail suture limb securely wrapped around the left forefinger, place the Knot Pusher under the left thumb to assume a single-handed position and to complete knot advancement.
21. With the Knot Pusher in place, tighten the knot by gently pulling on the white non-rail suture.
22. Remove the Knot Pusher from the tissue tract and from the suture without using the thumb knob.
23. Repeat these same steps 18–22 with the green suture.
24. If hemostasis is not complete, repeat the single-handed position for 20 seconds. After 20 seconds gently pull on the non-rail suture to tighten the knot. (Repeat first with the white suture and next with the green suture). Do not apply excessive pressure to the Knot Pusher or suture.
25. Once hemostasis is achieved, trim the sutures below the skin.

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
## TECHNIQUE FOR NEEDLE BACK-DOWN

The following describes a safety feature ("needle back-down") that permits the physician to return the needles into the sheath. This feature provides the option of exchanging the Prostar XL device with another Prostar XL device or an introducer sheath so that the patient may be treated with conventional compression therapy.

1. If resistance to counterclockwise rotation of the handle is felt, **do not attempt to deploy the needles.** Check that the interlocks are correctly aligned with the locking indents in the hub. The star on the edge of the hub should be centered directly between the interlocks.
2. If the needles are not easily deployed, back the needles down into the sheath prior to removal of the device.
3. Manually remove the coiled sections of the suture lumens to uncover the suture loops.
4. Use a hemostat to grasp the pull rod close to the hub. Advance the pull rod 1 cm into the core (see **Figure 7**) and then gently pull the slack out of the uncovered suture loops.

**Figure 7**

5. Repeat step 4 until the handle snaps into position on the proximal end of the device.
6. Pull back on the suture loops exiting the suture lumens to ensure that all slack in the sutures has been removed.
7. Prior to device removal, use fluoroscopy to verify the needle has been returned to the needle guide. **The needle tips should be as close as possible to the proximal edge of the radiopaque sheath ring before removal of the Prostar XL device.**
8. Perform the same "needle back-down" procedure in the event that all needles are not deployed. **DO NOT REMOVE ANY DEPLOYED NEEDLES.**
9. Do not attempt to re-deploy the Prostar XL device after the needles have been "backed-down". Replace the Prostar XL device with another Prostar XL device, an introducer sheath or use conventional compression therapy.

## SUTURE BREAKAGE

1. If suture breakage occurs prior to completing the initial knot, discard the suture material and remove the device over the guide wire. Use another Prostar XL PVS device to complete the procedure or replace the introducer sheath.
2. If the sutures are inadvertently tangled or removed prior to knot tying, discard the suture material and remove the Prostar XL device over the guide wire. Use another Prostar XL device to complete the procedure.
3. If suture breakage occurs after an initial knot has been tied, another Prostar XL device may be introduced or the introducer sheath may be replaced. If replacement of the introducer sheath is attempted, care should be taken to avoid excessive force during introduction. Any resistance to introduction should result in advancement of an introducer sheath small enough to be introduced without undue force, but large enough to maintain hemostasis.

## NEEDLE GUIDE BREAK PROCEDURE

The following describes a safety procedure that allows the physician to remove the Prostar XL device from the patient without surgical intervention, in the event that the Needle Guide of the Prostar XL device breaks.

**Procedure for Needle Guide Break pre needle deployment:**

1. Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined. **Advancing or withdrawing the Prostar XL device against resistance may cause the device to break.**
2. **Note:** though the guide will appear to be broken, the sheath remains attached to the device via the holder stop. The Prostar XL device has a redundant system that prevents the sheath from becoming completely detached from the device.
3. To remove a Prostar XL device with a Needle Guide Break that occurred during device insertion, evaluate if the Guide Wire Exit Port can be visualized above the skin.
4. If the Guide Wire Exit Port cannot be visualized:
   - 4.1. Manually remove the coiled sections of the suture lumens to uncover the suture loop.
   - 4.2. Pull back on the suture loops exiting the suture lumens to ensure that all slack in the sutures has been removed.
   - 4.3. Clamp the suture lumens with hemostats.
   - 4.4. Remove the device until the Guide Wire Exit Port is visualized.
5. Insert a 0.038" (0.97 mm) (or smaller) guide wire into the Guide Wire Exit Port and replace the Prostar XL device with another Prostar XL device, an introducer sheath, or use conventional compression therapy. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

**Procedure for Needle Guide Break post needle deployment:**

1. Do not compress the femoral access site while withdrawing the Prostar XL device from the tissue tract. **Applying excessive force to the Prostar XL sheath, due to site compression, during removal of the device may result in a needle guide break.**
2. Do not tighten the suture around the sheath. **Tightening the suture around the sheath during Prostar XL device removal may result in a needle guide break.**
3. Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined. **Advancing or withdrawing the Prostar XL device against resistance may cause the device to break.**
4. **Note:** though the guide will appear to be broken, the sheath remains attached to the device via the holder stop. The Prostar XL device has a redundant system that prevents the sheath from becoming completely detached from the device.
5. To remove the Prostar XL device when a Needle Guide Break has occurred during device removal, evaluate if the Guide Wire Exit Port can be visualized above the skin.
   - 5.1. If the Guide Wire Exit Port cannot be visualized, apply similar and simultaneous backward tension to both the circular handle and the device hub.
   - 5.2. Remove the device until the Guide Wire Exit Port is visualized.

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6. Insert a 0.038" (0.97 mm) (or smaller) guide wire into the Guide Wire Exit Port and replace the Prostar XL device with another Prostar XL device, an introducer sheath, or utilize conventional compression therapy. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

**POST PROCEDURE PATIENT MANAGEMENT**

1. Apply an appropriate dressing to the puncture site.
2. Assess the insertion site as per hospital protocol.

**COMPLICATIONS**

Potential complications resulting from procedures associated with the use of the Prostar XL PVS System include, but are not limited to:

- Localized vessel wall trauma which may lead to surgical repair
- Arterial thrombus
- Deep vein thrombosis
- Pseudoaneurysm
- Hematoma
- Local infection
- Nerve injury
- Local pulse deficits or ischemia
- Blood loss that may lead to blood transfusion
- Local discomfort
- Wound dehiscence.

**RECOMMENDATION FOR PATIENT AMBULATION AND DISCHARGE**

1. Patients may be ambulated after Prostar XL device procedures with normalized (150–180 seconds) Activated Clotting Times (ACTs).
2. Patients who have undergone Prostar XL device procedures may be eligible for earlier discharge from the hospital relative to conventional compression (i.e., manual or mechanical methods).
3. Before considering early discharge, assess for the following clinical conditions:
   - Conscious sedation
   - Anticoagulation, thrombolytic or antiplatelet therapy
   - Unstable cardiac status
   - Hematoma at the closure site
   - Hypotension
   - Pain while walking
   - Bleeding at the closure site
   - Any comorbid condition requiring observation.

   The presence of any of the above factors has generally led to the deferral of early discharge recommendations.

**PRODUCT INFORMATION DISCLOSURE**

Abbott Vascular Inc. has exercised reasonable care in the manufacture of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc. directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

**HOW SUPPLIED**

**Prostar XL Percutaneous Vascular Surgical System**

Each System Includes:
One (1) Prostar XL Percutaneous Vascular Surgical device
One (1) Perclose Knot Pusher

The Prostar XL PVS device and accessories are provided sterile and non-pyrogenic in unopened undamaged packages. Products are sterilized with ethylene oxide and intended for single use only. This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. Store in cool, dry place.

**PATENTS AND TRADEMARKS**

This product and / or its use may be covered by one or more of the following United States Patents: 5,417,699; 5,527,322; 5,613,974; 5,746,755; 5,779,719; 5,792,152; 5,797,929; 5,860,991; 5,902,311; 5,921,994; 6,036,699; 6,117,145; 6,206,893; 6,355,050; 6,517,553; 7,029,480; 7,445,626; EP 673228; EP 727965; EP 910288; EP 1236437. Other U.S. patents pending. Foreign patents issued and pending.

Prostar and Perclose are trademarks of the Abbott Group of Companies.

**Please refer to the end of this booklet for Graphical Symbols for Medical Device Labeling.**

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©2018 Abbott