Becton Dickinson
Infusion Therapy AB
Medication and Procedural Solution
Heather.Hagvik@bd.com

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Instructions for use. Purpose: Sterile intravenous peripheral single use catheter BD Venflon™ I is intended for access to peripheral veins of the patient's circulatory system for rehydration, parenteral nutrition, drug administration, blood transfusion, blood sampling and monitoring. Sterile in undamaged packaging. Precautions: Reuse may lead to infection or other disease / injury. Do not try to reinstall a partially or completely removed needle. Be careful to maintain the integrity of the catheter. Don't use scissors. Dispose of the device according to your hospital's protocol. Use only ISO compatible Luer connections. Non-ISO Luer lock connections may cause fluid / blood leakage during catheter handling. Possible complications associated with peripheral catheter insertion may include minor bleeding, hematoma, and failed insertion attempts. Directions for use: Medical personnel should be trained in the practice of venipuncture and catheter insertion and be aware of the associated complications. PTFE. Non-pyrogenic. DEHP or natural rubber latex are not part of the material.

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Becton Dickinson Infusion Therapy AB, Sweden
Becton Dickinson Infusion Therapy AB