Medtronic

Melody™

Melody™: Transcatheter Pulmonary Valve • **Melody™:** Транскатетърна пулмонална клапа • **Melody™:** Katetrizačně implantovaná pulmonální chlopeň • **Melody™:** Transkateter lungeklap • **Melody™:** Transkatheter-Pulmonalklappe • **Melody™:** Διακαθετηριακή πνευμονική βαλβίδα • **Melody™:** Válvula pulmonar transcatéter • **Melody™:** transkateteraalne pulmonaalklapp • **Melody™:** katetrilla implantoitava keuhkovaltimoläppä • **Melody™ :** Valve pulmonaire transcathéter • **Melody™:** Transkateterski plućni zalistak • **Melody™:** Transzkatéteres pulmonalis billentyű • **Melody™:** Valvola polmonare transcatetere • **„Melody™“:** transkateterinis plaučių vožtuvas • **Melody™:** transkatetra plaušu artērijas vārstulis • **Melody™:** Транскатетерски пулмонарен залистоκ • **Melody™:** Transkatheterpulmonalisklep • **Melody™:** Transkateter pulmonalklaff • **Melody™:** Przezcewnikowa zastawka pnia płucnego • **Melody™:** válvula pulmonar transcateter • **Melody™:** Válvula pulmonar transcateter • **Melody™:** Valvă pulmonară transcateter • **Melody™:** клапан легочной артерии транскатетерный • **Melody™:** Pulmonálna chlopňa pre transkatétrové zavedenie • **Melody™:** Transkatetrska pljučna zaklopka • **Melody™:** transkateterski pulmonalni zalistak • **Melody™:** Kateterburen pulmonalisklaff • **Melody™:** Transkateter Pulmoner Kapak • **Melody™:** транскатетерний клапан легеневої артерії

Ensemble™: Transcatheter Valve Delivery System • **Ensemble™:** Носеща система за транскатетърна клапа • **Ensemble™:** Systém pro katetrizační zavedení chlopně • **Ensemble™:** Transkateter klapindføringssystem • **Ensemble™:** Transkatheter-Klappeneinführsystem • **Ensemble™:** Διακαθετηριακό σύστημα προώθησης και τοποθέτησης βαλβίδας • **Ensemble™:** Sistema dispensador de la válvula transcatéter • **Ensemble™:** transkateteraalse klapi paigaldussüsteem • **Ensemble™:** katetrilla implantoitavan läpän asetin • **Ensemble™ :** Système d'implantation pour valve transcathéter • **Ensemble™:** Sustav za uvođenje transkateterskog zaliska • **Ensemble™:** Transzkatéteres billentyűbevezető rendszer • **Ensemble™:** Sistema di posizionamento per valvola transcatetere • **„Ensemble™“:** transkateterinio vožtuvo įvedimo sistema • **Ensemble™:** transkatetra vārstuļa ievadīšanas sistēma • **Ensemble™:** Транскатетерски систем за воведување залистоκ • **Ensemble™:** Plaatsingssysteem voor transkatheterkleppen • **Ensemble™:** Transkateter klaffinnføringssystem • **Ensemble™:** Przezcewnikowy system wprowadzający zastawki • **Ensemble™:** sistema de colocação da válvula transcateter • **Ensemble™:** Sistema de colocação da válvula transcateter • **Ensemble™:** Sistem de plasare a valvei transcateter • **Ensemble™:** устройство доставки клапана транскатетерного • **Ensemble™:** Systém pre transkatétrové zavedenie chlopne • **Ensemble™:** Sistem za vstavljanje transkatetrske zaklopke • **Ensemble™:** sistem za dopremanje transkateterskog zaliska • **Ensemble™:** Kateterburet klaffinföringssystem • **Ensemble™:** Transkateter Kapak Taşıma Sistemi • **Ensemble™:** система доставки транскатетерного клапана

Ensemble™ II: Transcatheter Valve Delivery System • **Ensemble™ II:** Носеща система за транскатетърна клапа • **Ensemble™ II:** Systém pro katetrizační zavedení chlopně • **Ensemble™ II:** Transkateter klapindføringssystem • **Ensemble™ II:** Transkatheter-Klappeneinführsystem • **Ensemble™ II:** Διακαθετηριακό σύστημα προώθησης και τοποθέτησης βαλβίδας • **Ensemble™ II:** Sistema dispensador de la válvula transcatéter • **Ensemble™ II:** transkateteraalse klapi paigaldussüsteem • **Ensemble™ II:** katetrilla implantoitavan läpän asetin • **Ensemble™ II :** Système d'implantation pour valve transcathéter • **Ensemble™ II:** Sustav za uvođenje transkateterskog zaliska • **Ensemble™ II:** Transzkatéteres billentyűbevezető rendszer • **Ensemble™ II:** Sistema di posizionamento per valvola transcatetere • **„Ensemble™ II“:** transkateterinio vožtuvo įvedimo sistema • **Ensemble™ II:** transkatetra vārstuļa ievadīšanas sistēma • **Ensemble™ II:** Транскатетерски систем за воведување залистоκ • **Ensemble™ II**: Plaatsingssysteem voor transkatheterkleppen • **Ensemble™ II:** Transkateter klaffinnføringssystem • **Ensemble™ II:** Przezcewnikowy system wprowadzający zastawki • **Ensemble™ II:** sistema de colocação da válvula transcateter • **Ensemble™ II:** Sistema de colocação da válvula transcateter • **Ensemble™ II:** Sistem de plasare a valvei transcateter • **Ensemble™ II:** устройство доставки клапана транскатетерного • **Ensemble™ II:** Systém pre transkatétrové zavedenie chlopne • **Ensemble™ II:** Sistem za vstavljanje transkatetrske zaklopke • **Ensemble™ II:** sistem za dopremanje transkateterskog zaliska • **Ensemble™ II:** Kateterburet klaffinföringssystem • **Ensemble™ II:** Transkateter Kapak Taşıma Sistemi • **Ensemble™ II:** система доставки транскатетерного клапана

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання

CE
0344

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1 Device description

The implant system consists of the Melody transcatheter pulmonary valve, model PB10 (stented bovine jugular vein valve) and either the Ensemble transcatheter valve delivery system, model NU10 or the Ensemble II transcatheter valve delivery system, model ENS10.

1.1 Melody transcatheter pulmonary valve (TPV)

The TPV consists of a heterologous (bovine) jugular vein valve sutured within a laser-welded, platinum-iridium stent with gold brazing of the welds. The TPV is available in a 16 mm bovine jugular vein (nominal length of 30 mm) and an 18 mm bovine jugular vein (nominal length of 28 mm). A final sterilization step is performed using a sterilant that contains 1% glutaraldehyde and 20% isopropyl alcohol, and in which the TPV is preserved and packaged until used. Adequate rinsing with isotonic saline solution must be performed before implantation to reduce the glutaraldehyde concentration.

See Table 1 for TPV sizing recommendations.

Table 1. TPV sizing

TPV catalog number	Diameter of intended implant location	Maximum expanded inner diameter of the TPV
PB1016	18 mm to 20 mm	20 mm
PB1018	18 mm to 22 mm	22 mm

1.2 Ensemble and Ensemble II transcatheter valve delivery system (delivery system)

Each delivery system consists of a balloon-in-balloon catheter with a retractable polytetrafluoroethylene (PTFE) sheath large enough to cover the TPV after crimping. The delivery system has a 22 Fr crossing profile. At inflation, the inner balloon is half the diameter of the outer balloon. Both balloons are made of nylon. The delivery system's sheath has a side-port used to flush the system and a hemostatic sleeve over the sheath to minimize bleeding at the insertion site. The delivery system has a polyether block amide (Pebax, Arkema Corporation) distal obturator that is conical in shape. The delivery system is compatible with a 0.889 mm (0.035 in) guidewire. The Ensemble II delivery system also includes radiopaque marker bands on the delivery system beneath the balloons. The marker bands are intended to aid in visibility of balloon location under fluoroscopy.

See Table 2 for delivery system sizes.

Table 2. Delivery system sizes

Ensemble delivery system catalog number	Ensemble II delivery system catalog number	Outer balloon outer diameter
NU1018	ENS1018	18 mm
NU1020	ENS1020	20 mm
NU1022	ENS1022	22 mm

2 Indications

The Melody TPV is indicated for use in patients with the following clinical conditions:

- Patients with regurgitant prosthetic right ventricular outflow tract (RVOT) conduits or bioprostheses with a clinical indication for invasive or surgical intervention, OR
- Patients with stenotic prosthetic RVOT conduits or bioprostheses where the risk of worsening regurgitation is a relative contraindication to balloon dilatation or stenting

3 Contraindications

- Venous anatomy unable to accommodate a 22 Fr size introducer sheath
- Implantation of the TPV in the left heart
- RVOT unfavorable for good stent anchorage
- Severe RVOT obstruction, which cannot be dilated by balloon
- Obstruction of the central veins
- Clinical or biological signs of infection
- Active endocarditis
- Known allergy to aspirin or heparin
- Pregnancy

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4 Warnings and precautions

Physician training: Physicians are required to have successfully completed the Physician/Center training program provided by Medtronic before using the Melody TPV and either the Ensemble or Ensemble II delivery system.

Device lifetime: The Melody TPV is used in the palliative treatment of RVOT conduit or bioprosthetic valve dysfunction with a goal of delaying the next surgical intervention. Some patients will require valve reintervention or replacement, the timing of which is dependent on several factors.

4.1 Warnings

4.1.1 TPV

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Do not resterilize the TPV by any method.

Exposure of the device and container to irradiation, steam, ethylene oxide, or other chemical sterilants will render the device unfit for use.

Do not use the device if:

- The device has been dropped, damaged, or mishandled in any way
- The Use By date has elapsed
- Each tamper-evident seal is broken
- The serial number tag on the TPV does not match the serial number on the container label
- The shipping temperature indicator window inside the shelf carton is black. If the shipping temperature indicator is black, the TPV is not suitable for clinical use.
- The storage solution does not completely cover the device or there is evidence of leakage

Do not expose the device to solutions other than the storage and rinsing solutions.

Do not add antibiotics to either the storage or the rinse solution. Do not apply antibiotics to the device.

Do not allow the device to dry. Maintain tissue moisture with irrigation or immersion.

Do not attempt to repair a damaged device.

Do not handle or use forceps to manipulate the valve leaflet tissue.

Do not use forceps to manipulate the stent.

Do not overexpand the device beyond the maximum recommended size as shown in Table 1, as this may result in a regurgitant TPV.

4.1.2 Delivery system

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Do not dilate the balloon appreciably beyond the diameter of the obstructive vessel and the intended implant site.

Do not use air or any gaseous substance as a balloon-inflation medium.

Do not advance the guidewire, balloon-dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.

Do not remove the guidewire from the delivery system at any time during the procedure.

4.2 Precautions

- Rinsing procedures of the TPV must be strictly followed to avoid necrosis of the tissue at the site of implantation.
- Exposure to glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to or breathing of the chemical vapor. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water for a minimum of 15 minutes. In the event of eye contact, flush with water for a minimum of 15 minutes and seek medical attention immediately.
- The sealed delivery system packaging should be inspected prior to opening. If the seal is broken or the packaging has been damaged, sterility cannot be assured.
- Proper functioning of the delivery system depends on its integrity. Use caution when handling the delivery system. Damage may result from kinking, stretching, or forceful wiping of the delivery system.
- This delivery system is not recommended for pressure measurement or delivery of fluids.

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- Maintain tight catheter connections and use aspiration before proceeding to avoid air introduction into the system.
- The delivery system must be carefully flushed to avoid the introduction of air bubbles.
- Before crimping (reducing) the size of the TPV on the balloon, the orientation and function of the valve must be tested.
- Do not take off the labels. These determine the orientation of the TPV before the crimping procedure on the balloon is completed.
- The crimping procedure must be carried out carefully. While crimping, the orientation of the valve must be known at all times. No change of orientation should occur as the TPV is mounted on the balloon.
- The crimping procedure must be carried out carefully. Do not place excessive pressure on the device during crimping as this can overstress the stent and lead to stent fracture.
- The inflation diameter of the balloon used during TPV delivery should approximate the diameter of the obstructive vessel and the intended implant site.
- Use of 2 appropriate size inflation devices (1 for each balloon) with pressure gauges is highly recommended during this procedure when inflating the balloon to deliver the TPV.
- The TPV is rigid and may make navigation through vessels difficult.
- Do not, under any circumstances, advance any portion of the delivery system against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- Balloon deployment should be conducted under fluoroscopic guidance with appropriate X-ray equipment.
- It is very important to ensure the balloons are completely deflated before pulling the delivery system back into the sheath.
- If resistance is felt upon removal, the delivery system and guidewire should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the guidewire and delivery system as a unit and withdrawing both together, using a gentle twisting motion combined with traction.

5 Potential complications/adverse events

The following list includes potential procedural complications that may result from implantation of the TPV:

- Rupture of the RVOT conduit
- Compression of a coronary artery
- Perforation of a major blood vessel
- Embolization or migration of the TPV
- Perforation of a heart chamber
- Arrhythmias
- Allergic reaction to contrast media
- Cerebrovascular events (transient ischemic attack [TIA], cerebrovascular accident [CVA])
- Infection/sepsis
- Fever
- Hematoma
- Radiation-induced erythema, blistering, or peeling of the skin
- Pain, swelling, or bruising at the catheterization site

The following list includes potential device-related adverse events that may occur following TPV implantation:

- Stent fracture[1]
- Stent fracture resulting in recurrent obstruction
- Endocarditis
- Embolization or migration of the TPV
- Valvular dysfunction (stenosis or regurgitation)
- Paravalvular leak
- Valvular thrombosis
- Pulmonary thromboembolism
- Hemolysis

6 Patient information

6.1 Anticoagulation/antiplatelet information

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on each patient's condition.

[1] The term "stent fracture" refers to the fracture of the Melody TPV. However, in subjects with multiple stents in the RVOT it is difficult to definitively attribute stent fractures to the Melody frame versus another stent.

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6.2 Identification and management of stent fractures

The potential for stent fracture should be considered in all patients who undergo TPV placement, regardless of the previously discussed or subsequently characterized risk factors.

6.2.1 Identification

For patients found to have a substantial increase in the degree of RVOT obstruction, the possibility of an associated stent fracture should be considered and evaluated. Radiographic assessment of the stent with chest radiography or fluoroscopy should be included in the routine postprocedural evaluation of patients who receive a TPV.

6.2.2 Management

If a stent fracture is detected, continued monitoring of the stent should be performed in conjunction with clinically appropriate hemodynamic assessment. In patients with stent fracture and significant associated RVOT obstruction or regurgitation, reintervention should be considered in accordance with usual clinical practice.

Reintervention may include implantation of an additional Melody TPV or surgical conduit replacement. Note that limited data are available in the Melody clinical studies on reimplantation of another Melody TPV within the original Melody TPV.

6.3 Endocarditis

Endocarditis is a potential adverse event associated with all bioprosthetic valves (Chapter 5). Patients should make their health care providers aware that they have a bioprosthetic valve before any procedure.

A low incidence of suspected endocarditis has been reported in patients implanted with the Melody TPV. Unexplained, prolonged fever may be an indication of infection, and patients with these conditions should be advised to seek medical attention.

Prophylactic antibiotic therapy is recommended for patients implanted with a Melody TPV undergoing dental procedures.

6.4 Traceability labels

Product traceability labels are provided for hospital use only and are not to be returned to Medtronic.

6.5 Registration information

A patient registration form is included in TPV packages for countries that require patient registration.

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU.

After implantation, please complete all requested information. The serial number is located on both the package and the identification tag attached to the TPV. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

Medtronic will provide an Implanted Device Identification Card to the patient. The card contains the name and telephone number of the patient's physician as well as information that medical personnel would require in the event of an emergency. Patients should be encouraged to carry this card with them at all times.

6.6 MRI safety information

Nonclinical testing and modeling has demonstrated that the Melody TPV is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial gradient magnetic field of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the Melody TPV is expected to produce a maximum in vivo temperature rise of less than 2.1°C after 15 minutes of continuous scanning.

MR image quality may be compromised if the area of interest is in the same area, or relatively close to the position of the device. In nonclinical testing, the image artifact caused by the device extends approximately 3 mm from the Melody TPV when imaged with a spin echo pulse sequence and 6 mm when imaged with a gradient echo pulse sequence and a 3.0 T MRI System. The lumen of the device was obscured.

For deployment of a Melody TPV inside a bioprosthesis or stent, consult the MRI labeling pertaining to the bioprosthesis or stent for additional image artifact information. It may be necessary to optimize MR imaging parameters for the presence of this implant.

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Scanning under the conditions defined above may be performed immediately after implantation.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

7 How supplied

7.1 Packaging

The TPV is chemically sterilized and provided **sterile** and **nonpyrogenic** in a sealed glass container with a screw cap. Sterility is compromised if each tamper-evident seal is broken, the container is damaged, or leakage is evident. The outside of the container is **nonsterile** and should not be placed in the sterile field.

The delivery system is **nonpyrogenic** and is sterilized with ethylene oxide gas and packaged in a double pouch. The delivery system is sterile if the pouches are undamaged and unopened. Do not use the delivery system if the outer pouch is damaged. The delivery system should never be stored in only the inner pouch. The inner pouch does not provide a sterile barrier. The outer surfaces of the outer pouch are **nonsterile** and must not be placed in the sterile field.

7.2 Storage

Store the TPV at 15°C to 25°C (59°F to 77°F). Store the delivery system at room temperature and away from direct sunlight. The delivery system should never be stored in only the inner pouch. The inner pouch does not provide a sterile barrier.

Appropriate inventory control should be maintained so that devices with the earliest Use By dates are preferentially used to avoid expiration. Devices must be used by the Use By date shown on the product labels.

8 Instructions for use

The following is a sequential outline of the catheterization/implant procedure. The type of diagnostic catheters, guidewires, dilation balloons, sizing balloons, or other tools needed are at the discretion of the operator.

8.1 Access site preparation and preimplant diagnostics

1. Perform sterile preparation and draping of the access site.
2. Gain arterial and venous access.

 Note: If femoral venous access is not available (for example, due to stenosis or trauma), an alternative venous access site, such as the internal jugular vein or subclavian vein, can be used for access.
3. Administer heparin to achieve a target ACT of >250 seconds.
4. Introduce a catheter into the arterial sheath and advance it into the ascending aorta. Perform an aortogram to demonstrate the coronary arteries are not adjacent to the RVOT, and that there is no risk of coronary compression when a stent or a TPV is implanted.

 If the coronaries appear to be in close proximity to the implant site, and coronary compression appears to be possible, further investigations should be done before moving forward with TPV implantation.
5. Advance an angiographic catheter to the right ventricle (RV) or proximal part of the RVOT for angiography. The angiographic projections obtained will be based on the relative position of the RVOT.

 If assessing a bioprosthesis, adjust the fluoroscopy angles so the bioprosthesis annulus and/or commissures are in one plane or in profile. Consider the features of the bioprosthesis when determining the optimal placement for the TPV.
6. To assess for suitability of TPV implantation, follow either Step 6a (for conduits) or Step 6b (for bioprosthetic valves).
 a. **For conduits:** Obtain angiographic measurements of the intended implantation site to assess suitability of the conduit for TPV implantation. If the angiographic measurements are unclear or if there is a question of the conduit being compliant, a low-pressure sizing balloon (<811 kPa [<8 atm]) may be used to further assess the current condition of the conduit. If the narrowest conduit dimension (usually in the lateral projection) is ≤18 mm and less than the original conduit diameter, predilate the site of conduit obstruction (the implantation site) to facilitate optimal relief of conduit obstruction. The size of the predilation balloon should be at least 2 mm greater than the narrowest diameter of the conduit in any projection, no more than 110% of the nominal conduit diameter. Perform another conduit angiogram to ensure the conduit is intact. If there is no conduit injury, repeat the predilation steps, using the same guidelines with increasingly larger balloons until appropriate implant diameter has been reached.
 b. **For bioprosthetic valves:** Obtain angiographic measurements of the intended implantation site to assess suitability for TPV implantation. If angiographic measurements are unclear or if there is a question of the inner diameter (ID) of the bioprosthesis being too large for implant (>22 mm), a low-pressure sizing balloon (<811 kPa [<8 atm]) may be used to further assess the ID. If the current ID of the bioprosthesis is <22 mm and is ≥2 mm less than the manufacturer's published ID, use a balloon to predilate the bioprosthesis. The balloon should not exceed the manufacturer's published ID of the bioprosthesis. If the bioprosthesis ID measures >22 mm, the bioprosthesis should be considered too large for Melody TPV implantation.

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Note: If obstruction is noted, further preparation of the implantation site (for example, predilation or pre-stenting) may be required to facilitate optimal relief of obstruction before implanting the TPV.

7. Perform final sizing using a balloon diameter that clearly demonstrates the narrowest portion of the conduit or bioprosthesis (at a pressure of <811 kPa [<8 atm]). If the coronary arteries are an acceptable distance from the implant site, the anatomy is considered suitable for Melody TPV implantation. See Table 1 for sizing recommendations.

 If questions about the coronary anatomy remain, repeat an aortic root angiogram or perform a selective coronary angiogram of the coronary system close to the conduit while simultaneously inflating the sizing balloon or largest predilation balloon to full expansion. If a coronary artery appears to be compressed by the balloon, the patient has unsuitable anatomy for the TPV, and the procedure should be abandoned.

8. Place a guidewire across the RVOT with its tip located as far distal in the pulmonary arterial bed as possible. Remove the diagnostic catheter. Leave the dilator in place and prepare the TPV and the delivery system.

8.2 Preparation of the TPV

1. Before opening, carefully examine the jar and lid for damage, leakage, or broken seals. **The jar should contain enough sterilant to cover the TPV. Rinse the TPV for a minimum of 2 minutes to reduce the glutaraldehyde concentration from the TPV as directed in the following steps.**
 a. Using aseptic technique, prepare 3 sterile bowls; 1 remaining empty and 2 containing isotonic saline solution (500 mL) for rinsing.
 b. Using aseptic technique, remove the TPV by grasping the serial number tag with atraumatic forceps and lifting it from the jar. The outside of the jar is **nonsterile**. Do not allow the TPV to come into contact with the outside of the jar.
 c. A serial number tag is sutured to the outflow aspect of the TPV. Verify that the serial number on the tag matches the jar label serial number and the serial number on the patient registration form. If any differences in serial number are noted, do not use the TPV. Do not detach the tag from the TPV until implantation is imminent.
 d. Drain the residual storage solution from the valve into the empty discard bowl (bowl 1) by holding the TPV with the serial tag (outflow end) downward.
 e. Transfer the empty TPV to the rinse bowl (bowl 2). Fill the TPV with rinse solution and alternately empty and fill by inverting and swirling, emptying and filling the valve for 1 minute, then empty the solution from the valve into the bowl.
 f. Transfer the empty TPV to the rinse bowl (bowl 3) and repeat step e for a minimum of 1 minute. Leave the TPV in the rinse bowl until implantation is imminent to prevent the tissue from drying.
 g. Empty the rinse solution from the TPV before loading the TPV onto the delivery system.
2. Remove the serial number tag by cutting the suture, which attaches the tag to the TPV.

8.3 Delivery system loading and placement of the TPV

1. Carefully examine the delivery system to confirm it was not damaged in shipment, and the balloon size is suitable for the intended procedure. The size of the delivery system to be used is based on the prepared implant site, as described in Table 1. The delivery system used should contain a TPV delivery balloon that is 0 to 2 mm larger than the waist at the implant site, but no more than 110% of the nominal conduit or bioprosthetic valve diameter.
2. Flush the delivery system guidewire lumen and the side port. Prepare the balloons by completely deflating with a syringe. Connect inflation syringes to the inner and outer balloon lumens. Inflation devices with pressure monitoring capabilities are recommended. Use a solution ratio of one-third to two-thirds contrast to saline in the inflation syringes.
3. Reduce the size of the TPV while crimping it using mandrels of decreasing size. It is recommended to use a 2.5 mL sterile syringe for crimping to an intermediate size prior to final crimping onto the balloon catheter. When reduced to the intermediate size, slide the TPV over the tip of the delivery system and center over the balloons.
4. Verify that the outflow end of the TPV is oriented toward the distal end of the delivery system. The outflow end of the TPV can be detected by the blue suture used to attach the valve at this location.
5. Gently crimp the TPV onto the balloon using finger pressure and a rolling action to exert equal pressure on all sides of the TPV to elongate the stent. Crimp only until no movement is felt on the delivery system.

 Avoid bending or twisting the TPV.

 Ensure the blue suture is adjacent to the blue tip of the delivery system to ensure proper orientation of the TPV for delivery.

6. Carefully slide the sheath over the TPV and balloons, ensuring that the crowns at the inflow end of the TPV do not get caught on the sheath as it is advanced over the TPV. Flush the sheath using the sidearm to remove air from the delivery system. Continue advancing the sheath and flushing until the sheath fits snugly over the proximal end of the blue tip.

 Note: Check to ensure that the blue suture is adjacent to the blue tip of the delivery system.

7. Remove the venous dilator from the access site, and advance the delivery system over the guidewire. Carefully introduce the delivery system through the skin and advance it toward the implantation site. Should bleeding be noted at the venous access site, advance the sleeve on the proximal end of the shaft into the vein to stop the bleeding. Advance the delivery system to the site

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where the TPV is intended to be deployed. This requires manipulations of the delivery system and the guidewire. **Maintain adequate guidewire position at all times.**

8. Uncover the TPV once it has reached the level of the implantation site. Hold the shaft of the delivery system in place while pulling back on the outer sheath. The TPV is fully uncovered when the flushing port of the outer sheath is aligned with the proximal marker on the delivery system shaft. While uncovering, ensure that there is no inadvertent movement of the TPV, which could possibly lead to the loss of the appropriate TPV position.
9. Once the TPV is fully uncovered, connect a syringe with pure contrast to the flush port of the delivery system. Inject a small amount of contrast to confirm good position of the TPV and full retraction of the sheath. Further minor adjustments of the position of the TPV are still possible at this point.
10. Inflate the inner balloon followed by the outer balloon. Once the TPV is properly deployed against the implantation site, deflate the balloons and carefully remove the delivery system.
11. Measure the RV-PA pressure. If a significant gradient exists, perform balloon dilation with a high-pressure balloon in order to minimize the gradient.
12. Compare the RV pressure with the systemic pressure measured through the arterial approach.
13. Perform a contrast injection into the proximal main PA (or distal RVOT) to demonstrate valve function and position. Ensure that the valve is not held open by the guidewire and the delivery system as this would give the false impression of pulmonary regurgitation.
14. Remove the catheters, guidewires, and sheaths while maintaining hemostasis.

 Note: The delivery system is for single use and should only be inflated once. Do not use the delivery system to expand the TPV after initial implantation. The TPV may be expanded after placement by using a different balloon catheter. To minimize the risk of creating regurgitation, do not exceed the maximum recommended size as shown in Table 1.

 Note: Circumferential tearing of the delivery balloon catheter prior to complete expansion of the TPV may cause the balloon to become tethered to the TPV, requiring surgical removal. In case of rupture of an adequately sized balloon after TPV expansion, the balloon can be withdrawn and a new balloon catheter can be exchanged over a guidewire to complete expansion of the TPV.

Table 3. Melody TPV sizing chart

Delivery system size – inner balloon / outer balloon	**Inner balloon maximum applied pressure** RBP		**Outer balloon applied pressures**		**Corresponding TPV outside diameter (mm) (balloon inflated)**
	atm	**kPa**	**atm**	**kPa**	
Size 18 mm Inner: 9 mm × 3.5 cm / Outer: 18 mm × 4 cm	**5**	**506**	1	101	17.9
			2	203	18.6
			3	304	19.4
			4 RBP	**405**	**20.1**
Size 20 mm Inner: 10 mm × 3.5 cm / Outer: 20 mm × 4 cm	**5**	**506**	1	101	19.7
			2	203	20.7
			3	304	21.7
			4 RBP	**405**	**22.4**
Size 22 mm Inner: 11 mm × 3.5 cm / Outer: 22 mm × 4 cm	**4.5**	**456**	1	101	21.8
			2	203	22.8
			3 RBP	**304**	**24.1**

Notes:

- Do not exceed **bolded** pressure values for either the inner or outer balloon of the delivery system size.
- RBP = rated burst pressure = maximum applied pressure
- atm = atmosphere
- kPa = kilopascal

9 Return of explanted TPVs

Medtronic is interested in obtaining recovered TPVs. Specific pathological studies of the explanted TPV will be conducted under the direction of a consulting pathologist. A written summary of the findings will be returned to the physician. To obtain a product return kit,

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contact a Medtronic distribution center or a Medtronic Sales Representative. If a kit is not available, place the explanted TPV in a container of glutaraldehyde or 10% buffered formalin immediately after excision. For further instructions on the return of an explanted device, contact a Medtronic Sales Representative.

9.1 Product incidents

Report any product incidents to a Medtronic Sales Representative or to the manufacturer, Medtronic B.V. The manufacturer's contact information is located on the last page of the IFU.

10 Disclaimer of warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

ALTHOUGH THE MELODY TRANSCATHETER PULMONARY VALVE MODEL PB10, ENSEMBLE TRANSCATHETER VALVE DELIVERY SYSTEM MODEL NU10, AND ENSEMBLE II TRANSCATHETER VALVE DELIVERY SYSTEM MODEL ENS10, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© 2019 Medtronic
The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.