Medtronic

Riptide™

Aspiration System

fr	Système d'aspiration
de	Aspirationssystem
it	Sistema per aspirazione
es	Sistema de aspiración
sv	aspirationssystem
nl	Aspiratiesysteem
pt-pt	Sistema de aspiração
fi	Aspiraatiojärjestelmä
da	Aspirationssystem
el	Σύστημα αναρρόφησης
cs	Aspirační systém
hu	Aspirációs rendszer
ru	Аспирационная система
pl	System do aspiracji
tr	Aspirasyon Sistemi
no	Aspirasjonssystem
sk	Aspiračný systém
ro	Sistem de aspirație
bg	Система за аспирация
et	Aspiratsioonisüsteem
hr	Sustav za aspiraciju
lt	aspiracinė sistema
lv	aspirācijas sistēma
sl	Aspiracijski sistem
pt-br	Sistema de aspiração
ko	흡인 시스템
ar	نظام شفط

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English en
Instructions for Use

Riptide™ Aspiration System

DEVICE DESCRIPTION
The Riptide™ Aspiration System is comprised of the following devices:

• React™ 68 Catheter	REACT-68
• React™ 71 Catheter	REACT-71
• Riptide™ Large Bore Aspiration Tubing	MAT-110-110
• Riptide™ Aspiration Pump	MAP-1000EU
• Riptide™ Collection Canister with Intermediate Tubing	MAC-1200

The Riptide™ Aspiration System is designed to remove occlusive thrombus from the cerebral vasculature using continuous aspiration. The catheter shaft has a hydrophilic coating that spans the distal 40 cm to reduce friction during use. The catheter acts as a conduit for providing direct aspiration to the occlusion site generated by the Riptide™ Aspiration Pump. The catheter is introduced through a guide catheter or sheath and into the intracranial vasculature and guided over a neurovascular guidewire and/or microcatheter to the site of the primary occlusion. The catheter is connected to the Riptide™ Aspiration Pump through the aspiration tubing. The Intermediate Tubing connects the Riptide™ Collection Canister to the Riptide™ Aspiration Pump. The Riptide™ Aspiration System in this configuration may be used to aspirate thrombus from the occluded vessel.

Note: The catheter and aspiration tubing are Applied Parts.

INTENDED PURPOSE / INDICATIONS FOR USE
The Riptide™ Aspiration System is intended for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

CONTRAINDICATION
There are no known contraindications.

WARNINGS
- The Riptide™ Aspiration System should only be used by physicians who have received appropriate training in interventional neuro-endovascular techniques for treatment of acute ischemic stroke. The Riptide™ Aspiration System is designed for setup by clinical support staff under the supervision of a trained physician.
- Do not advance, retrieve, or use any component of the Riptide™ Aspiration System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter against resistance may result in damage to the device or vessel.
- For vessel safety, do not perform more than three attempts at aspiration in the same vessel using the Riptide™ Aspiration System.
- Do not use devices in conjunction with or in place of those devices specified above when using the Riptide™ Aspiration System.

PRECAUTIONS
- The catheter and the aspiration tubing are supplied sterile and are intended for single use only. Do not resterilize or re-use. Resterilization and/or re-use may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location.
- Use the catheter, aspiration tubing, and Riptide™ Collection Canister with Intermediate Tubing prior to the "Use-by date".
- The Riptide™ Aspiration Pump and Riptide™ Collection Canister with Intermediate Tubing are supplied non-sterile and are intended for use outside of the sterile field only.
- The Riptide™ Collection Canister with Intermediate Tubing is intended for single use only. Dispose in accordance with standard biologic waste disposal procedures.
- Do not use kinked or damaged devices. Do not use open or damaged packages.
- Inspect the catheter before use to verify that its size and condition are suitable for the specific procedure.
- Use the catheter under fluoroscopic guidance to confirm placement of the catheter proximal to the thrombus prior to applying aspiration.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the aspiration tubing flow switch is in the open position for only the minimum time needed to remove thrombus.
- The recommended continuous duration of aspiration is a maximum of 90 seconds.
- If repositioning of the catheter is necessary during the revascularization procedure, such reposition should be performed over an appropriate neurovascular guidewire and/or microcatheter using standard techniques.
- Do not use automated high-pressure contrast injection equipment with the catheter because it may damage the device.
- Medical management and post stroke care should follow AHA/ASA guidelines and institutional protocols.
- As in all surgical interventions, monitoring of intra-procedural blood loss is recommended so that appropriate management may be instituted.
- Possible malfunctions of the Riptide™ Aspiration System include, but are not limited to, issues associated with vacuum for the aspiration tubing, Riptide™ Aspiration Pump, and Riptide™ Collection Canister with Intermediate Tubing that may lead to decreased therapeutic response during mechanical thrombectomy.
- Operators should take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

POTENTIAL ADVERSE EVENTS
- acute occlusion
- allergic reaction and anaphylaxis from contrast media
- arrhythmia
- complication at puncture site
- complication of radiation exposure (e.g., alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia)
- death
- decreased therapeutic response
- embolism
- fistula
- foreign body in patient
- hypertension/hypotension
- infection
- inflammation
- intracranial hemorrhage
- ischemia
- nerve damage
- neurological deficits including stroke
- organ failure/shock
- pain
- stenosis
- vessel spasm, thrombosis, dissection, or perforation

PROCEDURE
1. Refer to Warnings, Precautions, and Potential Adverse Events prior to use.
2. As each component of the Riptide™ Aspiration System is used, remove the component from the packaging and inspect for damage or for kinks.
3. Prepare a guide catheter or sheath according to the device manufacturer's *Instructions for Use*.
4. Place the guide catheter or sheath into the appropriate artery for access to the occlusion site.

CATHETER PREPARATION AND USE
1. Confirm vessel diameter and ensure the catheter is appropriate for treatment of the vessel.

Recommended Sizing Guidelines	
Model	Recommended Vessel Diameter (mm)
	Minimum
REACT-68	≥2.7
REACT-71	≥2.7

2. Prepare the catheter for use by flushing the packaging hoop and device with heparinized saline.
3. Remove the catheter from the packaging hoop.
4. Inspect the catheter for kinks or damage. If any damage is observed, replace with a new device.

 NOTE: When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).
5. Connect a rotating hemostasis valve to the proximal hub of a guide catheter or sheath.
6. Insert the catheter through the rotating hemostasis valve. If required, the provided split y-introducer may be used to facilitate the introduction of the catheter through the rotating hemostasis valve, inserting the tip through the introducer starting at the split end and penetrating the valve with the tubular introducer end. After the catheter has passed through the rotating hemostasis valve, the split y-introducer may be removed by backing it out of the valve, pulling gently on the split ends to separate the remaining portion. Discard the separated introducer.
7. Using conventional catheterization techniques under fluoroscopic guidance advance the catheter into the target vessel over an appropriate neurovascular guidewire and/or microcatheter. Position the catheter proximal to the thrombus. Remove the guidewire and/or microcatheter from the catheter.

 NOTE: Access to the site of occlusion may be aided by the triaxial use of a compatible inner microcatheter.


Figure 1: Riptide™ Aspiration System Setup Diagram

Riptide™ Aspiration System PREPARATION AND USE
1. Refer to **Figure 1** (above) for the correct arrangement of devices for proper setup of the Riptide™ Aspiration System.
2. Refer to the Riptide™ Aspiration Pump *User's Manual* for pump setup, operating instructions and troubleshooting if required.
3. Refer to the Riptide™ Collection Canister *Quick Reference Guide* for assembly setup.
4. Place the Riptide™ Aspiration Pump on a stable, flat horizontal surface outside of the sterile field within reach of the aspiration tubing and providing sufficient slack to the patient table. Locate the Riptide™ Aspiration Pump such that a minimum of 6 inches (approximately 15 cm) clearance exists in all directions to adjacent objects or other common cath lab equipment, powered or un-powered.
5. Attach the Riptide™ Aspiration Pump's power cord to the designated mains power specified in the Riptide™ Aspiration Pump *User's Manual*. Confirm that the power cord is securely plugged into the back of the pump enclosure. If required, attach the potential equalization connector to the potential equalization plug.
6. The Riptide™ Collection Canister lid should be pre-assembled to the cup as it is removed from its packaging. If it is not or becomes separated, secure the lid to the cup by snapping it tightly on the lip.
7. The large drain port on the Riptide™ Collection Canister lid should be pre-assembled with a closure cap as it is removed from its packaging. If it is not or becomes separated, seal the port with the large tethered cap on the lid by pressing it in place firmly. The small caps attached to the lid may remain tethered and unused until needed.
8. Place the Riptide™ Collection Canister into the receptacle on the Riptide™ Aspiration Pump and press firmly into place. The horizontal port marked with the embossed label as labeled "PATIENT" should point away from the pump enclosure to enable unobstructed attachment of the aspiration tubing.
9. The Intermediate Tubing should be pre-assembled to the Riptide™ Collection Canister lid as it is removed from its packaging. If it is not or becomes separated, secure the connector at one end of the supplied Intermediate Tubing over the vertical central port of the Riptide™ Collection Canister marked with the embossed label as "VACUUM" by pressing down firmly. Either connector end of the Intermediate Tubing may be used.
10. Attach other end of the Intermediate Tubing to the vacuum inlet of the Riptide™ Aspiration Pump by pressing down firmly.
11. Open the aspiration tubing pouch with the contents accessible toward the sterile field. Remove the tubing assembly from the pouch into the sterile field direction. Place or hold the distal end with the flow control switch on the patient table and pass the proximal end of aspiration tubing from the sterile field to the Riptide™ Collection Canister located on the Riptide™ Aspiration Pump.
12. Attach the connector at the proximal end of the aspiration tubing to the horizontal Riptide™ Collection Canister port marked with the embossed label as labeled "PATIENT". Press the connector firmly onto the port.
13. Slide the actuator of the flow control switch on the aspiration tubing to OFF to close the tubing to flow.
14. Press the start button on the Riptide™ Aspiration Pump and allow it to run for at least one minute. Confirm normal operation and the presence of vacuum as indicated on the gauge.
15. Using the vacuum control knob on the Riptide™ Aspiration Pump, adjust the vacuum to a reading of minimum 20 inHg (68 kPa) on the gauge but not exceeding 25 inHg (85 kPa). Adjust the vacuum by rotating the vacuum control knob until the gauge reads the desired vacuum level. To increase the vacuum, turn the knob clockwise. To decrease the vacuum, turn the knob counter-clockwise.

 Note: It is normal for the indicator needle of the gauge to oscillate or vibrate in response to pump pulsation. Allow the range of the needle to stabilize and adjust the vacuum control so that the midpoint of the motion aligns with the desired setting.
16. The Riptide™ Aspiration System is now ready to aspirate.
17. Connect luer fitting at the distal end of the aspiration tubing to the side port of the rotating hemostasis valve attached to the catheter.
18. Advance the distal tip of the catheter to engage with the proximal end of the thrombus or in accordance with the physicians preferred technique.
19. To begin aspiration, slide the flow control switch toward the ON position. To stop the aspiration flow, slide the actuator of the flow control switch toward the OFF position.
20. If necessary, additional passes with catheter may be used to further remove thrombus at the discretion of the physician.
21. After withdrawal from the patient, disconnect the aspiration tubing from the catheter; replace it with 5cc or 10cc syringe, and aspirate approximately 5cc of blood from the catheter to remove any thrombus that may remain.
22. Obtain a post-treatment angiogram in accordance with the physicians preferred technique.
23. At the conclusion of the aspiration procedure, turn OFF the Riptide™ Aspiration Pump. The catheter, the aspiration tubing, and the Riptide™ Collection Canister with Intermediate Tubing are single use. Remove and discard in accordance with standard biological waste disposal procedures.
24. Clean and prepare the Riptide™ Aspiration Pump for storage as described in the *User's Manual*.

WARRANTY DISCLAIMER
Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2019 Medtronic.
Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.
Micro Therapeutics, Inc.
Medtronic B.V.