en

STERRAD™ Instrument Tray Mats

INSTRUCTIONS FOR USE

INDICATIONS FOR USE
The STERRAD Instrument Mats are designed to protect surgical instruments in containment devices for sterilization in prevacuum steam, ethylene oxide and STERRAD Sterilization Systems.

GENERAL INFORMATION
1. Cleaning of the mats after the sterilization process is recommended after each use. A neutral pH or mild detergent is recommended.
2. Only STERRAD Instrument Mats have been validated by Advanced Sterilization Products for use with the APTIMAX™ Trays. (See product offering below.)

PRECAUTIONS
1. To enable maximum penetration/diffusion, do not stack the APTIMAX Trays in any sterilizer. In addition, do not overload trays and do not exceed a total weight of 4.4 kg (9.7 lbs.) per tray.
2. Discontinue use of the STERRAD Instrument Mats when they become cracked or damaged in any way.
3. To avoid cycle cancellations in the STERRAD System, ASP recommends not exceeding the following amount of mats per load:

STERRAD 100S System	2774 cm^2 (430 sq. in.)
STERRAD 50 System	1387 cm^2 (215 sq. in.)
STERRAD 200 System	5548 cm^2 (860 sq. in.)
STERRAD NX™ System	1387 cm^2 (215 sq. in.)
STERRAD 100NX™ System	2250 cm^2 (348 sq. in.)

4. If combinations of mats are used resulting in a cycle cancellation in the STERRAD System, please reduce the number of mats and reprocess the load.

Mat REF	Fits Tray REF
99205	13832
99206	13833
99207	13834
99208	13835
99209	13829
99210	13830
99211	13831
99212	13836
99213	13837

For a glossary of symbols used, go to eifu.asp.com.

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© ASP 2021