ASP

APTIMAX™

Instrument Tray Holders

Porte-Instruments
Instrumenten-Tray Halterungen
Sistema Di Bloccaggio Per Strumenti
Bandejas Soporte para Instrumental
Instrumenthouders
Fixadores Para Bandejas De Instrumentais
Holdere til Instrumentbakker
Instrumentinpitimet
Holdere for Instrumentbrett
Instrumenthållare
Βάσεις Απόθεσης Δίσκων Εργαλείων
滅菌トレーホルダー
Uchwyty do Tac na Instrumenty
Tálcák Műszertartói
Držiaky kontajnerov na nástroje
Instrumendi Korvide Raamid
Držáky Kontejnerů na Nástroje
Alet Tepsisi Tutaçları
Držači instrumenata
Держатели/разделители для лотков
Suporturi pentru tăvile de instrumentar

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APTIMAX Instrument Tray Holders

INSTRUCTIONS FOR USE

- REF: 99400, 99401, 99402, 99403, 99404, 99405, 99406, 99407
- Align holder as desired for the instrument set up and press holder into bottom of the tray.
- After insertion, check to ensure that the holder is secure before further processing.
- Holders are validated for use in pre-vacuum steam, ethylene oxide, and STERRAD™ Sterilization processes. For information on specific sterilization cycles that the holders can be processed in, please refer to the Instructions for Use for the APTIMAX Instrument Trays.
- Choose an appropriate size holder for the instrument. The instrument should be secure in the holder but not tight (e.g., the instrument should move freely within the holder but not fall out of the holder when the tray is inverted).
- Cleaning of the holders after the sterilization process is recommended after each use.
- Instrument Holding Bar, REF 99407, is designed to also fit the lid of APTIMAX Instrument Trays (REF 13826, 13827, 13828, 13829, 13830, 13831, 13832, 13833, 13834, 13835, 13836, 13837, 13838) to prevent instruments from bumping/touching the lid.
- Holders are not designed to fit in the 135XX, 136XX, and 137XX series of STERRAD Instrument Tray bases or lids.

For a glossary of symbols used, go to eifu.asp.com.

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