Medtronic

SELECTSITE™ C304-HIS

Deflectable catheter system
Système de cathéter orientable
Steuerbares Führungskathetersystem
Sistema de catéter flexible
Buigbare-kathetersysteem
Sistema a catetere orientabile
Styrbart katetersystem
Böjbart katetersystem
Σύστημα καμπτόμενου καθετήρα
Bøjeligt katetersystem
Sistema de cateter moldável
Hajlékony, irányítható katéterrendszer
System cewnika z odchylaną końcówką
Systém ovladatelného katétru
Гибкая катетерная система
Systém ohýbateľného katétra
Bükülebilir kateter sistemi
Система на огъваем катетър
Sistema de cateter defletível

Technical manual • Manuel technique • Technisches Handbuch • Manual técnico • Technische handleiding • Manuale tecnico • Teknisk håndbok • Teknisk manual • Τεχνικό εγχειρίδιο • Teknisk håndbog • Manual técnico • Műszaki leírás • Instrukcja techniczna • Technická příručka • Техническое руководство • Technická príručka • Teknik el kitabı • Техническо ръководство • Manual técnico

! USA **Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

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

English

1. Device description

The device features a guide wire to access the vein, a valve to reduce blood loss during the implant procedure, a deflectable catheter to introduce a transvenous device, a catheter dilator to facilitate deflectable catheter passage, and a guide catheter slitter to remove the deflectable catheter.

The deflectable catheter is designed for placement of transvenous devices in or near the bundle of His. It features a deflecting distal section, controlled by the deflectable catheter handle, and an out-of-plane curve on the distal tip. The body of the deflectable catheter is radiopaque for visibility on fluoroscopy.

1.1. Contents of package

Dispose of all single-use accessories according to local environmental requirements. The deflectable catheter delivery system is supplied sterile.

Each package contains the following items:

- 1 C304-HIS deflectable catheter
- 1 catheter dilator
- 1 valve
- 1 guide catheter slitter
- 1 guide wire with tip straightener
- product documentation

2. Indications for use

The device is indicated to provide a pathway through which diagnostic and therapeutic transvenous devices are introduced within the chambers and coronary vasculature of the heart and for introducing balloon catheters into the coronary sinus or leads into vessels of the left heart via the coronary sinus.

3. Contraindications

Use of the device is contraindicated for patients with obstructed or inadequate vasculature and for right ventricular use in patients with tricuspid valvular disease or a mechanical tricuspid heart valve.

4. Warnings and precautions

Transvenous device compatibility – Use the deflectable catheter only with compatible transvenous devices. No test data is available to demonstrate compatibility of the deflectable catheter with any device not manufactured by Medtronic.

- For the C304-HIS catheters, no test data is available to demonstrate compatibility of the transvenous devices with an outer diameter larger than 1.5 mm (4.6 Fr), with the exception of the 1.85 mm (5.6 Fr) dilator included with this package.

Consequences of using the deflectable catheter with incompatible devices may include the inability to deliver the transvenous device or damage to the transvenous device during delivery.

Single use – This product is intended for single use only.

Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized.

Vessel and tissue damage – Use care when passing a catheter through vessels and tissue.

Thrombus formation – Thrombogenicity evaluations were conducted using a heparinized model. If your patient cannot be adequately anticoagulated, it is unknown whether thrombus formation may occur with this product.

Inspecting the sterile package – Inspect the sterile package with care before opening it.

- Contact a Medtronic representative if the seal or package is damaged.
- Do not store this product above 40°C (104°F).
- Do not use the product after its expiration date.

Handling the deflectable catheter – Handle the deflectable catheter with care at all times.

- Do not kink, stretch, or severely bend the deflectable catheter.
- Do not use surgical instruments to grasp the deflectable catheter.
- Avoid contact with liquids other than isopropyl alcohol, blood, saline, or contrast solution.

Using the valve – Use the valve to impede the back flow of venous blood during the implant procedure.

- Ensure that the flush port stopcock is closed before attaching the valve to the deflectable catheter hub.

External defibrillation equipment – Keep external defibrillation equipment nearby for immediate use during acute lead system testing, the implant procedure, or whenever arrhythmias are possible or intentionally induced during post-implant testing.

5. Potential complications

The following are known potential adverse events associated with the use of this product.

Note: Usage of this product may result in adverse events, which may lead to injury, death, or other serious adverse reactions.

- Air embolism
- Allergic reaction to contrast media
- Arteriovenous fistula formation
- Bleeding at the insertion site
- Brachial plexus injury
- Cardiac tamponade
- Dislodgement
- Dissection
- Endocarditis
- Heart block
- Hematoma formation
- Hemothorax
- Infection
- Irregular heart beat
- Mediastinal widening
- Perforation
- Pneumothorax
- Subclavian artery puncture
- Thrombophlebitis
- Thrombosis
- Valve damage
- Vascular occlusion
- Vessel damage

6. Directions for use

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The following procedures are provided for information only. Each physician must apply the information in these instructions according to professional medical training and experience.

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6.1. Preparing the deflectable catheter for insertion

Warning: The deflectable catheter must be thoroughly flushed with saline and free of air prior to use to avoid an air embolism in the patient.

Warning: Do not use a power injection syringe to inject contrast solution.

1. Remove the components of the sterile package within the sterile field.
2. Verify that the deflectable catheter deflects properly by rotating the handle clockwise (Figure 1). Each time the handle is rotated, the distal section is further deflected. After the rotating has stopped, the distal section will hold the amount of deflection for the number of times the handle was rotated.


Figure 1. Deflecting and straightening the deflectable catheter

3. Verify that the deflectable catheter straightens properly by rotating the handle counterclockwise.
 Note: After the deflectable catheter has been deflected and straightened, the deflectable catheter may not be completely straight. If necessary, insert a catheter dilator and carefully manipulate the distal section of the deflectable catheter to straighten it as much as possible. While the curve on the distal end is fixed, deflecting and inserting the dilator will affect the curvature. The catheter dilator prevents the deflectable catheter shape from becoming deformed.
4. Attach the valve to the hub of the deflectable catheter (Figure 2).


Figure 2. Attaching the valve to the deflectable catheter

1. Valve
2. Stopcock
3. Flush port

5. Flush the deflectable catheter through the valve flush port with saline or heparinized saline prior to use.
6. Flush the catheter dilator.
7. Insert the catheter dilator through the valve and into the deflectable catheter. The distal end of the catheter dilator should extend past the distal end of the deflectable catheter.
8. Set the deflectable catheter assembly aside within the sterile field.

6.2. Inserting the deflectable catheter assembly

Warning: After the deflectable catheter assembly has been advanced into the right atrium, slowly remove or retract the catheter dilator to minimize the risk of perforation.

Caution: To avoid damaging tissue and the tip of the deflectable catheter assembly, keep the catheter dilator fully inserted into the deflectable catheter during passage to the right atrium.

Caution: If the subclavian approach is required for deflectable catheter assembly insertion, avoid techniques that may damage the deflectable catheter assembly.

- The insertion site should be as far lateral as possible to avoid clamping the deflectable catheter assembly between the clavicle and the first rib.
- Do not force the deflectable catheter assembly if significant resistance is encountered during insertion or passage.
- Do not use techniques such as adjusting the patient's posture to facilitate passage. If resistance is encountered, it is recommended to use an alternate venous entry site.

The deflectable catheter assembly may be inserted through several different venous routes, including the right or left cephalic vein or the external or internal jugular vein. Use the cephalic vein whenever possible to avoid damage to the deflectable catheter assembly in the first rib/clavicular (thoracic inlet) space.

Insert the deflectable catheter assembly:

1. Gain venous access.
2. Straighten the J-shaped tip of the guide wire by advancing the tip straightener over the end of the guide wire.
3. Insert the guide wire into the venous system, using the tip straightener.
4. Remove the tip straightener and advance the guide wire until the distal tip reaches the right atrium.
5. Place the deflectable catheter assembly over the proximal end of the guide wire.
6. Advance the deflectable catheter assembly through the insertion vein and then to the right atrium.
7. Remove the catheter dilator, leaving the guide wire and deflectable catheter in place.

6.3. Positioning the deflectable catheter in the right atrium or right ventricle

1. If positioning the deflectable catheter in the right ventricle, pass the guide wire through the tricuspid valve. Then use the guide wire to pass the deflectable catheter through the tricuspid valve.
2. Position the deflectable catheter toward the desired location in the right atrium or right ventricle.
3. Use fluoroscopy to verify that the deflectable catheter is in the desired position.

6.4. Positioning the deflectable catheter in the coronary sinus

1. Rotate the deflectable catheter posteriorly and to the patient's left to enter the coronary sinus ostium. Use fluoroscopy and inject small amounts

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of contrast solution to verify that the deflectable catheter is in the desired position.

2. Advance the deflectable catheter into the coronary sinus.

6.5. Obtaining a venogram

When placing a lead in the coronary sinus, venograms should be obtained to indicate the anatomy of the coronary vasculature.

- When using the SelectSite C304-HIS catheter, refer to the product documentation packaged with a 1.85 mm (5.6 Fr) or smaller venogram balloon for further instructions on obtaining a venogram.

6.6. Inserting the transvenous device into the deflectable catheter and positioning the transvenous device

Insert the transvenous device into the deflectable catheter. After obtaining the desired amount of deflection for placement of the transvenous device, refer to the product documentation packaged with the transvenous device for further positioning instructions.

If necessary, use the deflectable catheter to reposition the transvenous device.

6.7. Removing the deflectable catheter from the transvenous device

Warning: Before attempting to remove the deflectable catheter, undeflect the catheter as much as possible to avoid vessel damage during removal.

Warning: Keep the cutting surface of the slitter in line with the hub and catheter shaft axis. Do not move the catheter or slitter in a back and forth motion. This ensures that the hub and guide catheter wall are cleanly slit by the slitter blade, preventing transvenous device damage or dislodgements.

Caution: During and after slitting the deflectable catheter, handle the slit edges of the deflectable catheter carefully to prevent transvenous device damage.

Note: Remove the deflectable catheter with care. Removing the deflectable catheter may dislodge the transvenous device. If the transvenous device dislodges, refer to the product documentation supplied with the transvenous device for repositioning instructions.

When the transvenous device is in final position, the deflectable catheter must be removed from the transvenous device before surgical closure. The deflectable catheter is removed from the transvenous device using a slitter.

Remove the deflectable catheter from the transvenous device:

1. Remove the valve from the deflectable catheter hub (Figure 3).



Figure 3. Removing the valve

2. Follow the appropriate directions as provided in Appendix A for information on slitting.
3. Discard the deflectable catheter and slitter after the deflectable catheter has been removed from the body of the transvenous device.


7. Detailed device description

7.1. Specifications

Parameter		Model C304-HIS
Guide wire	Material:	Stainless steel
	Length:	120 cm (47.2 in)
	Diameter:	0.09 cm (0.035 in)
Deflectable catheter	Material:	Polyether block amide
	Usable length:	43 cm (16.9 in)
	Inner diameter:	1.9 mm (5.7 Fr)
	Outer diameter:	2.8 mm (8.4 Fr)

Deflectable catheter reach		44 mm (1.73 in)
Curve diameter		79 mm (3.11 in)

Catheter dilator	Material:	Polyethylene
	Usable length:	57.9 cm (22.8 in)
	Outer diameter:	1.85 mm (5.6 Fr)
Valve	Inner diameter:	5.1 mm (15.2 Fr) maximum
Deflectable catheter slitter		1.3 mm (4.0 Fr) - 2.0 mm (6.0 Fr)
Compatible leads		1.5 mm (4.6 Fr) maximum

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8. Medtronic disclaimer of warranty

For complete disclaimer of warranty information, see the accompanying disclaimer of warranty document.

9. Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate address or telephone number listed on the back cover.

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© 2020 Medtronic
The following are trademarks of Medtronic:
Medtronic USA, Inc.
Medtronic B.V.