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INSTRUCTIONS FOR USE

MitraClip™ and TriClip™ ACCESSORIES

Stabilizer, Ref No. SZR01ST
Support Plate, Ref No. PLT01ST
Lift, Ref No. LFT01ST

WARNING: Read all instructions carefully. Failure to follow these instructions, warnings and precautions may lead to device damage, user injury or patient injury. Failure to properly clean and / or re-sterilize prior to use may result in infections, malfunction of the MitraClip and TriClip Accessories, user injury or patient injury. Use of the MitraClip and TriClip Accessories should be restricted to those physicians trained in the proper use of either the MitraClip or TriClip System and procedure.

TABLE OF CONTENTS

INTENDED USE 4
PATIENT AND REUSABLE ACCESSORY PREPARATION PRIOR TO PROCEDURAL USE 4
CONTRAINDICATIONS 4
WARNINGS 4
PRECAUTIONS 4
POTENTIAL COMPLICATIONS AND ADVERSE EVENTS 4
HOW SUPPLIED 4
STORAGE AND HANDLING 4
MITRACLIP AND TRICLIP REUSABLE ACCESSORIES 4
DESCRIPTION 4
DIMENSIONS 4
MATERIAL LIST 5
REQUIRED CLEANING EQUIPMENT NOT INCLUDED 5
STABILIZER CLEANING AND STERILIZATION 5
MANUAL CLEANING 5
AUTOMATIC CLEANING 5
STERILIZATION 5
SUPPORT PLATE CLEANING 5
LIFT CLEANING 5
FUNCTIONAL INSPECTION 6
PATENTS 6

INTENDED USE
Reusable Accessories are intended to be used only with the MitraClip System or TriClip System product family during the MitraClip or TriClip procedure.

PATIENT AND REUSABLE ACCESSORY PREPARATION PRIOR TO PROCEDURAL USE
Follow the instructions regarding patient and / or reusable accessory preparation provided with the applicable MitraClip or TriClip System IFU.

NOTE: The Stabilizer, Support Plate and Lift are provided non-sterile. The Support Plate and Lift must be cleaned prior to use. The Stabilizer must be cleaned and sterilized prior to use. Follow the cleaning and / or sterilization instructions provided in this Instructions for Use. Ensure Lift and Support Plate are covered completely by sterile drape during the MitraClip or TriClip procedure. Use towels as necessary to minimize direct contact between the patient and all surfaces of both the Lift and Support Plate.

CONTRAINDICATIONS
None

WARNINGS
Read all instructions carefully. Failure to follow these instructions, warnings and precautions may lead to MitraClip and TriClip Accessories damage, MitraClip or TriClip System damage, user injury, or patient injury.

Use universal precautions for biohazards and sharps while handling the MitraClip and TriClip Accessories to avoid user injury.

Use of the MitraClip and TriClip System Accessories should be restricted to those physicians trained in the proper use of either the MitraClip or TriClip System and procedure.

Failure to properly clean and / or re-sterilize prior to use may result in infections, malfunction of the MitraClip or TriClip Accessories, user injury or patient injury.

Use of automatic washers to clean the Support Plate and Lift may result in damage or malfunction.

PRECAUTIONS
Inspect all product prior to use. Do not use if the package is open or damaged, or if product is damaged.

POTENTIAL COMPLICATIONS AND ADVERSE EVENTS
Failure to properly clean, sterilize (Stabilizer only) and inspect the MitraClip and TriClip Accessories per the instructions provided in this Instructions for Use may result in the following:

Allergic reaction	Fever
Failure to deliver Implant to the intended site	Valve injury
Infection	Vessel perforation or laceration

HOW SUPPLIED

Contents:
SZR01ST: One (1) Stabilizer

PLT01ST: One (1) Support Plate

LFT01ST: One (1) Lift

NOTE: MitraClip and TriClip Accessories are packaged separately.

Non-sterile:
The Stabilizer, Support Plate and Lift are provided non-sterile. The Support Plate and Lift must be cleaned prior to use. The Stabilizer must be cleaned and sterilized prior to use. Follow the cleaning and / or sterilization instructions provided in this Instructions for Use.

STORAGE AND HANDLING

NOTE: Proper care and handling of Stabilizer, Lift and Support Plate is required to ensure that each reusable accessory maintains its appearance, corrosion resistance, and proper function over time.

New Reusable Accessory Device
Handle with care. Store in original packaging prior to initial use. Keep dry. Keep away from sunlight.

NOTE: The Stabilizer, Support Plate and Lift are provided non-sterile. The Support Plate and Lift must be cleaned prior to use. The Stabilizer must be cleaned and sterilized prior to use.

Reusable Accessory Device after Cleaning and / or Sterilization
Stabilizer: Handle with care. Store in double wrapped standard sterile drape. Keep dry. Keep away from sunlight.

Support Plate and Lift: Handle with care. After cleaning, cover device to avoid contamination. Keep dry. Keep away from sunlight.

MITRACLIP AND TRICLIP ACCESSORIES

DESCRIPTION
The reusable accessories used in conjunction with either the MitraClip or TriClip System include: 1) a Stabilizer, 2) a Support Plate and 3) Lift. The Stabilizer is provided separately as a non-sterile reusable device and must be cleaned and sterilized prior to each use. During the procedure the Stabilizer is used in the sterile field to support and position the Steerable Guide Catheter and Clip Delivery System. The Support Plate and Lift are provided separately as non-sterile reusable devices and must be cleaned prior to each use. The Lift and Support Plate are used outside the sterile field to provide a stable platform for the Stabilizer and MitraClip or TriClip System during the procedure.

DIMENSIONS

Component	Approximate Dimension	Packaged Weight
Stabilizer	*520 mm (L) x 110 mm (W) x 210 mm (H)*	*~ 7 lbs (~ 3 kg)*
Support Plate	*(Unfolded) 660 mm (L) x 508 mm (W) x 9 mm (H)*	*~ 11 lbs (~ 5 kg)*
Lift	*(Fully extended) 427 mm (L) x 353 mm (W) x 290 mm (H)*	*~ 7 lbs (~ 3 kg)*

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
Figure 2: Support Plate



Figure 3: Lift

Number	Name
Stabilizer	
1	*Threaded Hole (7)*
Support Plate	
2	*Hinge*
3	*Plexiglas Plate (2)*
Lift	
4	*Black Mat*
5	*Orange Mat*
6	*Leg (4)*
7	*Leg Extension (4)*
8	*Rubber Tip (4)*

MATERIAL LIST

Component	Material	Quantity
Stabilizer		
Stabilizer	*304 Stainless Steel*	*1*
Support Plate		
Base	*Clear Polycarbonate*	*2*
Hinge	*Steel*	*1*
#8-32x1/4" Bolt	*Stainless Steel*	*9*
Lift		
Rubber Leg Tip	*Rubber*	*4*
Leg Extension, Long	*Hard Black Anodized Aluminum Tubing*	*2*
Leg Extension, Short	*Hard Black Anodized Aluminum Tubing*	*2*
Snap Button	*Zinc-Plated Steel*	*4*
Stool	*Stainless Steel*	*1*
Adhesive	*Cyanoacrylate*	*N/A*
Black Mat	*Butyl Rubber*	*1*
Adhesive	*Silicone Adhesive*	*N/A*
Orange Mat	*Silicone*	*1*

REQUIRED CLEANING EQUIPMENT NOT INCLUDED

Manual Cleaning

- Enzol (or equivalent enzymatic detergent)
- LF2100 (or equivalent low-foam non-enzymatic detergent)
- Appropriate sized non-abrasive nylon bristle brushes (or equivalent)
- Appropriate sized non-abrasive nylon bristle channel brushes (or equivalent), (Stabilizer only)
- Standard polyester wipers
- Clean, dry towels

Automated Cleaning (Stabilizer Only)

- Automatic Washer
- Enzol (or equivalent enzymatic detergent)
- LF2100 (or equivalent low-foam non-enzymatic detergent)

STABILIZER CLEANING AND STERILIZATION

CAUTION: Failure to properly clean and dry the device may lead to a reduction in device life. DO NOT use steel wool, wire brushes, pipe cleaners, or abrasive detergents and DO NOT use saline to clean or soak the device. DO NOT use chlorine bleach.

WARNING: Failure to properly clean and dry the device may lead to inadequate sterilization.

- All of the following steps should be followed.
- Follow detergent manufacturer instructions for proper handling and detergent solution preparation.
- Cleaning may become difficult if debris is allowed to dry on the device. Cleaning is recommended immediately after use.

MANUAL CLEANING

Pre-cleaning:

1. If attached, unthread Fasteners from Stabilizer and discard.

CAUTION: Failure to remove Fasteners may result in device damage and / or inadequate cleaning.

2. Prepare LF2100 (non-enzymatic) or Enzol (enzymatic) (or equivalent to either detergent) detergent solution.
3. Rinse the Stabilizer under running tap water at room temperature (~22° C) for a minimum of 1 minute.

Cleaning:

1. Completely submerge the Stabilizer in the detergent solution for a minimum of 3 minutes.
2. Brush all surfaces of the Stabilizer with a non-abrasive nylon bristled brush for a minimum of three (3) minutes, re-wetting the brush as necessary.
3. Clean threaded holes (7) with appropriately sized nylon bristle channel brush.
 a. Thoroughly brush each threaded hole for a minimum of thirty (30) seconds each.
4. Immediately following cleaning with detergent, thoroughly rinse Stabilizer under running tap water at room temperature (~22° C) for a minimum of 2 minutes.
5. Wipe all surfaces with a clean dry towel immediately.

Drying:

1. Prior to sterilization thoroughly air dry the Stabilizer until no visible moisture is observed.
2. After Stabilizer is completely dry, prepare device for steam sterilization. Double wrap the device with standard sterilization wrap.

AUTOMATED CLEANING

1. If attached, unthread Fasteners from Stabilizer and discard.

CAUTION: Failure to remove Fasteners may result in device damage and / or inadequate cleaning.

Stabilizer cleaning may also be performed with an appropriate automatic washer (see equipment list above).

2. For cleaning with an automatic washer, the parameters in Table 1 should be followed at a minimum.

Table 1: Stabilizer Automatic Washer Parameters

Detergent	Cycle Parameters									
	Pre-wash		Detergent Soak		Wash		Rinse		Dry	
	Time (min.)	Temperature	Time (min.)	Temperature	Time (min.)	Temperature	Time (min.)	Temperature	Time (min.)	Temperature
LF2100, Enzol, or equivalent to either detergent	2:00	Cold	1:00	Hot	2:00	150° F (66° C)	1:15	180° F (82° C)	8:00	240° F (116° C)

3. After Stabilizer is completely dry, prepare device for steam sterilization. Double wrap the device with standard sterilization wrap.

STERILIZATION

- The sterilization method for the Stabilizer is Prevacuum Steam Sterilization.

CAUTION: DO NOT use other sterilization methods, e.g., flash, ‡Steris™, or ‡Sterrad™.

- The minimum recommended exposure temperature and time to achieve a sterility assurance level (SAL) of 10^{-6} using the Prevacuum method is as follows:

Table 2: Stabilizer Steam Sterilization Parameters

Method	Cycle Parameters		
	Temperature	Time	Drying Time
Prevacuum Steam Sterilization	270° F to 275° F (132° C to 135° C)	4 Minutes	45 Minutes

Storage Prior to Use:

1. After device is sterilized, follow instructions for storage and handling as outlined in the STORAGE section of this Instructions for Use.

SUPPORT PLATE CLEANING

CAUTION: Failure to properly clean, dry and store the device may lead to a reduction in device life. DO NOT use steel wool, wire brushes, pipe cleaners, or abrasive detergents and DO NOT use saline to clean or soak the device. DO NOT use chlorine bleach.

- All of the following steps should be followed.
- Follow detergent manufacturer instructions for proper handling and detergent solution preparation.
- Cleaning may become difficult if debris is allowed to dry on the device. Cleaning is recommended immediately after use.

Pre-cleaning:

1. Rinse the Support Plate under running tap water at room temperature (~22° C) for a minimum of 30 seconds.

Cleaning:

2. Prepare LF2100 (non-enzymatic) or Enzol (enzymatic) (or equivalent to either detergent) detergent solution.
3. Fully saturate standard polyester wiper with detergent solution.
4. Wipe all surfaces of the Support Plate with saturated wiper for a minimum of 2 minutes, re-wet wiper as necessary to maintain saturation.
5. Wipe the entire hinge for minimum of 30 seconds re-wetting the wiper as necessary to maintain saturation.
 a. Open and close the plate as necessary to access all surfaces of the hinge.
6. Immediately following cleaning with detergent, thoroughly rinse Support Plate under running tap water at room temperature (~22° C) for a minimum of 2 minutes.
7. Wipe all surfaces with a clean dry towel immediately.

Drying:

8. Allow the Support Plate to air dry until no visible moisture is observed.

Storage Prior to Use:

9. After device is completely dry, follow instructions for storage and handling as outlined in the STORAGE section of this Instructions for Use.

LIFT CLEANING

CAUTION: Failure to properly clean, dry and store the device may lead to a reduction in device life. DO NOT use steel wool, wire brushes, pipe cleaners, or abrasive detergents and DO NOT use saline to clean or soak the device. DO NOT use chlorine bleach.

- All of the following steps should be followed.
- Follow detergent manufacturer instructions for proper handling and detergent solution preparation.
- Cleaning may become difficult if debris is allowed to dry on the device. Cleaning is recommended immediately after use.

Pre-cleaning:

1. Remove telescoping legs (black) from Lift.

CAUTION: Failure to remove telescoping legs may result in inadequate cleaning.

2. Rinse the Lift and each telescoping leg under running tap water at room temperature (~22° C) for a minimum of 30 seconds each.

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Cleaning:

3. Prepare LF2100 (non-enzymatic) or Enzol (enzymatic) (or equivalent to either detergent) detergent solution.
4. Fully saturate standard polyester wiper with detergent solution.
5. Wipe all surfaces of the Lift with saturated wiper for a minimum of 2 minutes, re-wet wiper as necessary to maintain saturation.
 a. During cleaning brush the exposed surface of the black mat on each side for a minimum of 30 seconds using a non-abrasive nylon bristle brush.
6. Wipe all surfaces of each telescoping leg with saturated wiper for a minimum of 30 seconds.
 a. If necessary, use an appropriate sized nylon bristle brush to access all surfaces.
7. Immediately following cleaning with detergent, thoroughly rinse Lift under running tap water at room temperature (~22° C) for a minimum of 2 minutes.
8. Immediately following cleaning with detergent, thoroughly rinse each telescoping leg under running tap water at room temperature (~22° C) for a minimum of 30 seconds each.
9. Wipe all surfaces of the Lift and telescoping legs with a clean dry towel immediately.

Drying:

10. Allow the Lift and telescoping legs to air dry until no visible moisture is observed.

Re-assembly:

11. After the Lift and telescoping legs are completely dry, re-assemble the telescoping legs to the Lift.

Storage Prior to Use:

12. After device is completely dry, follow instructions for storage and handling as outlined in the STORAGE section of this Instructions for Use.

FUNCTIONAL INSPECTION

STABILIZER

Inspection should be performed within the sterile field prior to each use.

i. Ensure that the Stabilizer is free of corrosion (i.e. rust).
ii. Assemble the sterilized Stabilizer by threading the Fasteners into the Stabilizer.
iii. Ensure that the Fasteners can be fully threaded into the Stabilizer holes.

SUPPORT PLATE

Inspection should be performed outside the sterile field prior to each use.

i. Ensure that the hinge is free of corrosion and securely attached to each polycarbonate piece (2) of the Support Plate.
ii. Ensure the Support Plate is free of cracks and / or gouges.
iii. Open and close the plate to ensure the hinge functions properly.

LIFT

Inspection should be performed outside the sterile field prior to each use.

i. Ensure metal surfaces are free of corrosion.
ii. Ensure rubber leg tips are free of cracks and / or slits.
iii. Ensure that each leg extension moves freely up and down along the leg and can maintain desired position.
iv. Ensure both mats are securely adhered on the Lift.

CAUTION: DO NOT use any of the above devices if damage is detected.

PATENTS

This product and / or its use are covered by one or more of the following United States Patents: 8,057,493; 7,736,388; 7,682,369; 7,666,204; 7,655,015; 7,608,091; 7,604,646; 7,563,267; 7,288,097; 7,226,467; 7,048,754; 6,770,083; 6,752,813; 6,629,534; 6,461,366. Other U.S. patents pending. Foreign patents issued and pending.

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

THIS PAGE INTENTIONALLY LEFT BLANK

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© 2019 Abbott. All Rights Reserved.