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EN F7830S BIPOLAR CABLE, LENGTH 200CM – STERILE DIRECTIONS FOR USE

DESCRIPTION
Threshold cable to use during the application of permanent endocavitary electrodes for permanent pacemaker implantation.
The cable, model **F7830S,** is designed for connecting any type of electrode plug for temporary or permanent epimiocardic or endocavitary stimulation to any type of device designed to detect the Sensing and Pacing thresholds.
The length of cable, supplied sterile, is 200 cm which allows for keeping the measuring device at a distance from the sterile field.
The model **F7830S** cable consists of one bipolar flexible copper cables sheathed in grey polyurethane. The cable can be connected by a 2mm red/black Touch-proof plug (E-DIN 428802-2) at one end and a red/black protected crocodile clip at the other.
The cable model **F7830S** is also supplied in the variants **F7830/AF**, **F7830DIN**, **F7831S**, **F7832S** and **F7837S**.
Model **F7830/AF** is made up by protected red and black crocodiles and by 4mm female red and black sockets.
Model **F7830DIN** is made up by protected red and black crocodiles and DIN 1.5 mm connection.
In model **F7831S** the ends are the same as model F7830S, whereas the pair of flexible copper cables sheathed in grey polyurethane is in blue colour.
The cable model **F7832S** consists of a pair of flexible copper cables sheathed in grey polyurethane 400 cm length.
The model **F7837S** it is the same as F7830S but with black PVC DEHP-FREE insulation.

CONTENTS OF THE BOX
The cardboard box on sale contains ten primary packets.
Each primary packet is sterile and contains one bipolar cable with red/black terminals in double shrink-wrapped bags in paper and polyethylene, except for the F7837S cable which is packaged in one shrink-wrapped bags in paper and polyethylene.

STERILIZATION
The product is disposable and sterilized with EO gas - do not reuse. Reusing the product may lead to:
- Possible presence of biological residues that may cause cross-infections
- Alteration of materials
- Loss of initial functional features of product

INDICATIONS
For detecting the Sensing and Pacing thresholds during the installation of temporary or permanent epimiocardic or endocavitary electrodes.

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WARNINGS
ATTENTION: in models that have cable longer than 3 meters, the cable may be subject to electromagnetic disturbances that, when particularly intense, could affect the operation of the device to which it is connected. If necessary remove noise sources from the cable.
All the patient-supporting devices must be compliant with the standards in force:
- Check the packaging before use. If wet on the inside, no longer intact, or damaged in any way, then the sterility can no longer be guaranteed
- After having checked the integrity of the packaging, check that the connecting wire between the terminals is integral, with no interruptions; in the presence of damage or visible defects do not use the product and return it to FIAB
- External cardiac stimulators, and the equipment for recording endocavitary potentials, must all be complying with the regulations in force.

NB: before making connections consult the instruction manual pertaining to the instrument to be used.
- All the electrical instruments must be opportunely earthed.
- Make sure that there is a defibrillator on hand.
- Medical devices must be connected and activated by qualified personnel only.

PREPARATION AND USE
PERMANENT CARDIAC PACING
- Prepare the subcutaneous pocket in which the permanent pacemaker will be inserted.

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- Position the catheter in the cardiac chamber which is to be stimulated.
- Extract the cable from the package after having controlled its integrity.
- Insert the plugs in the analyzer sockets of a PSA (Pacing Sensing Analyzer).
- Connect the red alligator clip to the subcutaneous pocket and the black alligator clip to the pin of the plug of the lead in position.

ATTENTION: In the case of a bipolar lead, the black alligator clip must be connected to the terminal of the catheter plug, and the red alligator clip must be applied to the ring pole of the plug itself.

- Activate the "SENSING" circuit of the PSA and measure the maximum value expressed in mV of the spontaneous electrical activity of the chamber under examination. If this value is too low, it will be opportune to change the position of the catheter.
- Activate the "PACING" circuit of the PSA and set a stimulation frequency of approx. 10 beats more than the spontaneous frequency of the chamber under examination. Set the duration of the stimulus at 0,5 m/sec. Begin the stimulation with a stimulus amplitude value of at least 5V or 5mA, and while continuing to stimulate, search for the minimum amplitude capable of picking up the cardiac chamber in a stable manner.

TEMPORARY ELECTRO-STIMULATION

- Position the electrode for temporary cardiac pacing (intravenous lead or epimocardial heart wire/s).
- Insert the red plug in the "+ " socket (usually red) and the black plug in the "-" socket (usually black) of the external cardiac stimulator.
- Connect the red alligator clip to the distal ring electrode, and the black alligator clip to the electrode of the tip of the intravenous lead. In the case of hearth wire the references "+" and "-" refer to the cranio-caudal position of one electrode with respect to the other. Consequently the red alligator clip will be connected to the electrode in the cranial position while the black alligator clip will be connected to the electrode in the caudal position.
- Set the "sensing" of the external cardiac stimulator on the minimum value.
- Set the duration of the stimulus at 0,5 m/sec. and the amplitude at 10 V or 10 mA.
- Begin the stimulation and gradually reduce the amplitude in order to ascertain the minimum threshold capable of creating stable chamber stimulation.
- Set the amplitude at a value equal to least three times the minimum threshold.
- Gradually increase the "sensing" value in order to ascertain the maximum value above which the spontaneous activity is no longer detected by the stimulator, and set an intermediate value

POTENTIAL COMPLICATIONS

The product must never be used if visible damage is present along the cable and/or its terminals.

The lead or heart wire, whichever is inserted, and to which the cable is connected, constitutes a direct connection with low resistance from the outside towards the heart. For this reason its is indispensable for all the equipment used to be compliant with the regulations in force. Particular attention must be paid to the earthing of all the equipment positioned near the patient.

STORAGE

The product must be stored in a suitable environment. The temperature and humidity limits are shown on its label.

GUARANTEE - NOTICE

FIAB SpA guarantees that the product complies with Directive 93/42/EEC and has been manufactured according to the procedures of FIAB Quality System certified ISO 13485. No responsibility may be ascribed to the producer who shall not be held liable for medical costs, direct or indirect damage due to lacking function or malfunction of the above products, when used differently from the instructions for use. We recommend to report opportunely any malfunction or defect of the product to FIAB Quality Assurance Service.

WASTE DISPOSAL

Waste coming from hospitals must be disposed of in accordance with regulations in force.

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EN

SERIES F7833 – F7834 – F7838
THRESHOLD CABLE FOR CARDIAC STIMULATION – STERILE
DIRECTIONS FOR USE

DESCRIPTION
Threshold cable to use during the application of permanent endocavitary electrodes for permanent pacemaker implantation.
The cable is designed for connecting any type of electrode plug for temporary or permanent epimiocardic or endocavitary stimulation to any type of device designed to detect the Sensing and Pacing thresholds.
The length of cable, supplied sterile, allows for keeping the measuring device at a distance from the sterile field.
The cable consists of one bipolar flexible copper cables sheathed in grey polyurethane. The cable can be connected by a 2mm red/black plug with protective cap at one end and a red/black protected crocodile clip at the other.
F7833 cable can be also produced with a blue insulating covering, changing its reference to **F7834** and with a black PVC DEHP-FREE insulation covering, changing its reference to **F7838**.

CONTENTS OF THE BOX
The cardboard box on sale contains ten primary packets.
Each primary packet is sterile and contains one bipolar cable with red/black terminals in double shrink-wrapped bags in paper and polyethylene, except for the F7838 cable which is packaged in one shrink-wrapped bags in paper and polyethylene.

STERILIZATION
The product is disposable and sterilized with EO gas – do not reuse. Reusing the product may lead to:
- Possible presence of biological residues that may cause cross-infections;
- Alteration of materials;
- Loss of initial functional features of product;

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
INDICATIONS
For detecting the Sensing and pacing thresholds during the installation of temporary or permanent epimiocardic or endocavitary electrodes.

CONTRAINDICATIONS
No contraindications exist if used as indicated by the manufacturer.
Do not use as an extension cable for temporary cardiac pacing or monitoring of the patient's heartbeat.
Do not use beyond the duration of the implantation procedure.

WARNINGS
ATTENTION: in models that have cable longer than 3 meters the cable may be subject to electromagnetic disturbances that, when particularly intense, could affect the operation of the device to which it is connected. If necessary remove noise sources from the cable.

All the patient-supporting devices must be compliant with the standards in force.
- Check the packaging before use. If wet on the inside, no longer intact, or damaged in any way, then the sterility can no longer be guaranteed
- After having checked the integrity of the packaging, check that the connecting wire between the terminals is integral, with no interruptions; in the presence of damage or visible defects do not use the product and return it to FIAB
- External cardiac stimulators, and the equipment for recording endocavitary potentials, must all be complying with the regulations in force.

NB: before making connections consult the instruction manual pertaining to the instrument to be used.
- All the electrical instruments must be opportunely earthed.
- Make sure that there is a defibrillator on hand.
- Medical devices must be connected and activated by qualified personnel only.
- After removing the protective caps, the accidental contact with the metallic part of the pin (free terminals far from the patient) could cause dangerous arrhythmias due to induced currents transmitted to the cardiac tissue through the bare pin; so the cable-lead connection must be done after connecting the cable to the device, on the contrary cable-lead disconnection must be done before disconnecting the cable to the device (see figure).

PREPARATION AND USE
PERMANENT CARDIAC PACING
- Prepare the subcutaneous pocket in which the permanent pacemaker will be inserted.
- Position the catheter in the cardiac chamber which is to be stimulated.
- Extract the cable from the package after having controlled its integrity.
- Insert the plugs in the analyzer sockets of a PSA (Pacing Sensing Analyzer).
- Connect the red alligator clip to the subcutaneous pocket and the black alligator clip to the pin of the plug of the lead in position.

ATTENTION: In the case of a bipolar lead, the black alligator clip must be connected to the terminal of the catheter plug, and the red alligator clip must be applied to the ring pole of the plug itself.
- Activate the “SENSING” circuit of the PSA and measure the maximum value expressed in mV of the spontaneous electrical activity of the chamber under examination. If this value is too low, it will be opportune to change the position of the catheter.
- Activate the “PACING” circuit of the PSA and set a stimulation frequency of approx. 10 beats more than the spontaneous frequency of the chamber under examination. Set the duration of the stimulus at 0,5 m/sec. Begin the stimulation with a stimulus amplitude value of at least 5V or 5mA, and while continuing to stimulate, search for the minimum amplitude capable of picking up the cardiac chamber in a stable manner.

TEMPORARY ELECTRO-STIMULATION
- Position the electrode for temporary cardiac pacing (intravenous lead or epimocardial heart wire/s).
- Insert the red plug in the “+ “ socket (usually red) and the black plug in the “-“ socket (usually black) of the external cardiac stimulator.
- Connect the red alligator clip to the distal ring electrode, and the black alligator clip to the electrode of the tip of the intravenous lead. In the case of hearth wire the references “+” and “-“ refer to the cranio-caudal position of one electrode with respect to the other. Consequently the red alligator clip will be connected to the electrode in the cranial position while the black alligator clip will be connected to the electrode in the caudal position.
- Set the “sensing” of the external cardiac stimulator on the minimum value.
- Set the duration of the stimulus at 0,5 m/sec. and the amplitude at 10 V or 10 mA.
- Begin the stimulation and gradually reduce the amplitude in order to ascertain the minimum threshold capable of creating stable chamber stimulation.
- Set the amplitude at a value equal to least three times the minimum threshold.
- Gradually increase the “sensing” value in order to ascertain the maximum value above which the spontaneous activity is no longer detected by the stimulator, and set an intermediate value

POTENTIAL COMPLICATIONS
The product must never be used if visible damage is present along the cable and/or its terminals.
The lead or heart wire, whichever is inserted, and to which the cable is connected, constitutes a direct connection with low resistance from the outside towards the heart. For this reason its is indispensable for all the equipment used to be compliant with the regulations in force. Particular attention must be paid to the earthing of all the equipment positioned near the patient.

STORAGE
The product must be stored in a suitable environment. The temperature and humidity limits are shown on its label.

GUARANTEE - NOTICE
FIAB SpA guarantees that the product complies with Directive 93/42/EEC and has been manufactured according to the procedures of FIAB Quality System certified ISO 13485. No responsibility may be ascribed to the producer who shall not be held liable for medical costs, direct or indirect damage due to lacking function or malfunction of the above products, when used differently from the instructions for use. We recommend to report opportunely any malfunction or defect of the product to FIAB Quality Assurance Service.

WASTE DISPOSAL
Waste coming from hospitals must be disposed of in accordance with regulations in force.

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