CorVocet™

BIOPSY SYSTEM

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUKSANVISNING

BRUGSANVISNING

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INSTRUCTIONS FOR USE

DESCRIPTION OF THE PRODUCT

The Merit Corvocet™ Biopsy System is a sterile single patient use device comprised of a core needle biopsy device and optional coaxial introducer. It is available in several needle gauge sizes and lengths. The top and rear actuator buttons on the device are color coded according to the various gauge sizes, Yellow=20 gauge, Pink=18 gauge, Purple=16 gauge and Green=14 gauge.

INDICATIONS FOR USE

The disposable Corvocet Biopsy System is intended for use in obtaining core biopsy samples from soft tissues such as liver, kidney, prostate, spleen, breast, lung, lymph nodes and various soft tissue tumors. It is not intended for use in bone.

CONTRAINDICATIONS

Physician judgment is required when considering biopsy on patients with bleeding disorders, or receiving anticoagulant medications.

℞ Only **Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physcian trained and/or experienced in the use of this device.

WARNINGS

- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.
- Unusual force applied to the needle or unusual resistance against the needle may cause the needle to deform or bend. A bent or deformed needle may not function properly.
- Post-biopsy care may vary with the biopsy technique utilized and the individual patient's physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.
- The collection of multiple needle cores may help to ensure the detection of any cancer tissue. A "negative" biopsy in the presence of suspicious radiographic findings does not preclude the presence of carcinoma.
- The Corvocet Biopsy System is not intended for use in bone.

PRECAUTIONS

- This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific physiology being biopsied.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Inspect the product for damage. If damaged, appropriately discard the entire device and prepare a new instrument.
- These instructions for Corvocet are NOT meant to define or suggest any medical or surgical technique. The individual physician is responsible for the proper procedure and techniques to be used with this product.
- Never test the product by firing into the air. Damage may occur to the device and could result in patient and/or user injury.
- The introduction of the needle into the body should be carried out under imaging control (ultrasound, X-Ray, CT). This product has not been tested for MR imaging compatibility.

NOTE: If collecting multiple samples, inspect the needle for a damaged point, bent shaft or other imperfections after each sample is collected. Do not use the needle if any imperfection is noted.

POTENTIAL COMPLICATIONS

Potential complications associated with core biopsy procedures are site specific and include, but are not limited to:


- Air Embolism
- Bleeding
- Unintended Organ Injury
- Track Seeding
- Infection
- Perforation
- Hematoma
- Hemorrhage
- Adjacent tissue injury
- Pain
- Hemoptysis
- Hemothorax
- Non-target tissue, organ or vessel perforation
- Pneumothorax

Air embolism is a rare but serious potential complication of lung biopsy procedures. Rapid deterioration of neurological status and/or cardiac arrhythmia may be indicative of air embolism. Prompt diagnosis and treatment must be considered if the patient exhibits signs or symptoms of air embolism.


INSTRUCTIONS FOR USE

BIOPSY DEVICE PREPARATION


1. Using aseptic technique remove the biopsy device from its package. Remove protective needle sheath that is secured to device and pull back slightly on the top trigger to remove packaging clip.

2. Adjust device throw length to desired level (usable range is 10 mm – 25 mm). Each "click" is 0.5 mm. The colored line on the adjustment dial indicates the penetration depth.


3. Energize (prime) device by pulling back on the top or rear trigger until you feel a click. When the device is primed, a yellow indicator will show in the indicator window (labeled "ready").

4. Advance the rear or top trigger until you feel a slight increase in resistance.

NOTE: Advancing the trigger past this point of increased resistance will cause the device to fire.


5. Device is ready to fire (primed).

BIOPSY PROCEDURE (WITHOUT OPTIONAL COAXIAL INTRODUCER)

1. Using aseptic technique, prepare site as required. Adequate anesthesia should be considered prior to procedure.
2. Verify instrument is ready to fire (primed).
3. Advance the needle tip of the biopsy device to a position proximal to the desired biopsy site.


NOTE: Once fired, the biopsy needle will travel to the preset penetration depth beyond the biopsy needle tip.

4. While maintaining the device's position and needle orientation, depress the rear trigger or push the top trigger distally to cause the needle to automatically advance.

NOTE: For devices that include a safety, it must be disengaged before the device will fire. The safety will toggle to red when ready to fire.

5. Remove the device from the patient.
6. Pull proximal on the top or rear trigger until you feel a click. See Figure below. This motion expels the specimen.

NOTE: This will prime the biopsy device.

7. Remove the specimen.
8. If additional biopsies are required, advance the rear or top trigger until you feel a slight increase in resistance.

NOTE: Advancing the triggers past this point of increased resistance will cause the device to fire.

9. Device is ready to fire (primed). Repeat steps 2-8 for all additional biopsies.

NOTE: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. Do not use needle if any imperfection is noted.

10. After your last sample is collected it is recommended that you disengage the device by dry firing it before you discard it.

BIOPSY PROCEDURE (WITH OPTIONAL CORVOCET COAXIAL INTRODUCER)


Select the appropriate Corvocet Coaxial Introducer to correspond with the gauge and length of the Corvocet Biopsy System in use. Refer to catalog for product numbers and descriptions.

NOTE: An optional blunt tip stylet may be included with 18 and 20 gauge Corvocet Biopsy Systems. The blunt tip stylet may be used to manipulate through soft tissue and around vasculature or other organs to minimize the risk of unintentional damage to these areas.

1. Using aseptic technique, prepare site as required. Adequate anesthesia should be considered prior to procedure.
2. Using aseptic technique, remove the Corvocet Coaxial Introducer from its package.
3. If preferred, advance the supplied depth stop to the predetermined placement depth.


NOTE: The depth stop should be adjusted so that the coaxial introducer is in proper position when the depth stop is in contact with the skin. This will help stabilize the Corvocet Coaxial Introducer.

4. Using image guidance, insert the tip of the Corvocet Coaxial Introducer proximal to the lesion to be biopsied using the depth stop as an aid for proper placement and adjust as necessary.
5. Hold the guiding cannula needle hub. Squeeze the stylet hub and pull proximally to remove the stylet from the outer cannula. Leave the cannula in place as a guide for the placement of the Corvocet biopsy device.

6. Verify instrument is ready to fire (primed).
7. Place the Corvocet biopsy device all the way through the guiding cannula to the lesion to be biopsied.

NOTE: There should not be a gap between the hub of the Corvocet Coaxial Introducer and the proximal end of the biopsy device.

8. Ensure device needle tip is at the correct location.

NOTE: Once fired the biopsy needle will travel to the preset penetration depth beyond the biopsy needle tip.

9. While maintaining the device's position and needle orientation, depress the rear trigger or push the top trigger distally to cause the needle to automatically advance.

NOTE: For devices that include a safety it must be disengaged before the device will fire.

10. Remove the device from the patient.
11. Pull proximal on the top or rear trigger until you feel a click. This motion expels the sample.

NOTE: This will prime the biopsy device.

12. Remove the specimen.
13. If additional biopsies are required, advance the rear or top trigger until you feel a slight increase in resistance.

NOTE: Advancing the triggers past this point of increased resistance will cause the device to fire.

14. Device is ready to fire (primed). Repeat steps 6-12 for all additional biopsies.

NOTE: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. Do not use needle if any imperfection is noted.

15. After your last sample is collected it is recommended that you disengage the device by dry firing it before you discard it.

Recommendation: When collecting multiple samples, wipe the core tissue biopsy needle with sterile moist gauze prior to reinsertion into the guiding cannula. This will aid in proper movement of the Corvocet biopsy needle within the cannula.

⚠	Caution: Consult accompanying documents
2	Single use only
	Do not use if package is damaged

REF	Catalog Number
LOT	Lot Number
	Use By
STERILE EO	Sterilized Using Ethylene Oxide

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MERITMEDICAL®

CE 2797

Manufacturer:
Merit Medical Systems, Inc. 1600 West Merit Parkway, South Jordan, Utah 84095
U.S.A. 1-801-253-1600 U.S.A Customer Service 1-800-356-3748

EC REP

Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222

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INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

The Achieve Programmable Automatic Biopsy System is a spring loaded, core needle biopsy system with delay or automatic firing options. The automatic firing mode releases the cannula and stylet in rapid sequence and captures a tissue sample with the push of a button. The delay firing method is used to verify the position of the sample notch in the target area prior to depressing the A button, which activates the cutting cannula, capturing the tissue sample.


Figure 1: Achieve Programmable Automatic Biopsy System

INTENDED USE

Intended for use in obtaining core biopsy samples from soft tissue such as kidney, liver, prostate, spleen, lymph nodes, and various soft tissue masses. Not intended for use in bone.

The Achieve Programmable Automatic Biopsy System is also indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. It is designed to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality (e.g., malignancy). When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

WARNING: Test firing the needle may damage the cutting cannula and/or stylet tip and is not recommended. If test firing is performed, the cutting cannula may be protected by holding the cannula near the clear taper with moderate pressure between the thumb and fingers during firing. Keep fingers or other objects from in front of the needle as it will advance 2.5 cm very rapidly when activated. Never test fire the instrument in an unsafe direction.

In the automatic mode, the Achieve automatic biopsy device is designed to fire the cannula and stylet in rapid succession. Please exercise caution before proceeding as the device will advance 2.5 cm when activated.

Ensure the device is fully discharged before disposing of the device.

Only clinicians, who have been properly trained in this procedure, preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

Good medical judgment should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders.

Post-biopsy care may vary with the biopsy technique utilized and the individual patient's physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.

CAUTION

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Do not use if package is damaged.

Do not resterilize. After resterilization, the sterility of the device is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprossessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to the potential adverse effects on components that are influenced by thermal and/or mechanical changes.

NOTE: If collecting multiple samples, inspect the needle for a damaged point, bent shaft or other imperfections after each sample is taken.

NOTE: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.

If biasing the device prior to obtaining a sample, take care to avoid bending the stylet notch.

STEP 1: COAXIAL INTRODUCER PLACEMENT
During the pre-procedure set up, insert the biopsy device through the coaxial introducer needle to verify compatibility.

Adjust the sliding skin stop on the coaxial introducer needle to the desired location and lock in place. Insert the coaxial introducer needle until the sliding skin stop rests on the skin surface. Use imaging, if desired, to verify the tip is at or near the edge of the region to be sampled. Remove the inner stylet.


Figure 2: Coaxial Introducer Needle

AUTOMATIC & DELAY MODE

STEP 2: PREPARATION OF ACHIEVE FOR USE
To charge the Achieve biopsy device, simply pull back on the charging handle and release twice, first to charge the cutting cannula, second to charge the inner stylet. **See Figures 3 through 5.**

FIGURE 3:
Non-charged and after firing.

FIGURE 4:
Pull back the charging handle once to retract the cannula which exposes the notched stylet.

FIGURE 5:
Pull back the charging handle a second time to retract the notched stylet.

STEP 3: POSITION THE NEEDLE
Position the needle at or near the edge of the region to be sampled.

STEP 4: FIRE THE INSTRUMENT
Two firing sequence methods are available once the Achieve biopsy device is charged.
A = Automatic mode
D = Delay mode (Stylet Only)

1. Delay Mode – Depression of the **D** button releases the inner stylet into the sample area. The operator uses this delay firing method to depressing the **A** button, which activates the cutting cannula, capturing the tissue sample.
 ition of the sample notch in the target area prior
 NOTE: If desired, a slight bias may be imparted to the stylet prior to firing by forcing the stylet notch against the lesion. Take care to avoid bending the stylet notch.
2. Automatic Mode – TThe operator can choose the fully sequential mode, omitting the delay feature, simply by depressing the **A** button first. This releases the stylet and cutting cannula in sequence, capturing the tissue sample.

STEP 5: BIOPSY SPECIMEN
To retrieve the tissue sample, remove the device from the patient, hold it with the **A** and **D** buttons facing up, and pull back and release the charging handle once. This action will retract the outer cutting cannula, exposing the stylet notch with the tissue sample. **See Figure 4**.

STEP 6: OBTAINING MULTIPLE SAMPLES
After the tissue sample is removed, additional samples from the same patient may be obtained. To obtain additional samples, pull the charging handle a second time, which will retract the stylet within the cutting cannula, placing the Achieve device in a fully charged position. **See Figure 5**.

MANUAL STYLET INSERTION MODE
Reasons for performing a biopsy in this manner include: a desire to push through a dense lesion slowly with the combined strength of the stylet and cannula, and a desire to visualize the placement of the needle within the lesion prior to firing.

STEP 1: POSITION THE NEEDLE
Position the needle such that the tip of the needle extends through the target area.

STEP 2: REVEAL THE STYLET NOTCH
While holding the handle steady, pull back on the charging handle one time to charge the cannula and reveal the stylet notch. Take care to maintain the position of the stylet tip while charging the cannula.

STEP 3: FIRE CANNULA
NOTE: If desired, a slight bias may be imparted to the stylet prior to firing by forcing the stylet notch against the lesion. Take care to avoid bending the stylet notch. Press the A button to fire the cannula through the target site.

STEP 4: BIOPSY SPECIMEN
To retrieve the tissue sample, remove the device from the patient, hold it with the **A** and **D** buttons facing up, and pull back and release the charging handle once. This action will retract the outer cutting cannula, exposing the stylet notch with the tissue sample.

STEP 5: OBTAINING MULTIPLE SAMPLES
After the tissue sample is removed, additional samples from the same patient may be obtained. To obtain additional samples, press the **A** button while holding the outer cutting cannula adjacent to the handle. This action will place the Achieve device in the uncharged position, and prevent damage to the stylet or cannula tips.

SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
LOT	Lot Number
REF	Catalog Number
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Use if Package is Damaged
	Single Use
℞ ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.
MD	Medical Device
Sterile Package	Sterile Package

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INSTRUCTIONS FOR USE

DEVICE DESCRIPTION


The TEMNO™ A.C.T. Biopsy System is a single patient use core biopsy device. It is available in several needle gauge sizes and lengths. The plunger is color coded according to the various gauge sizes (e.g., gray=22G, green=20G, orange=18G, white=14G).

INDICATIONS FOR USE

Intended for use in obtaining core biopsy samples from soft tissue such as breast, kidney, liver, lung, thyroid, and various soft tissue masses. Not intended for use in bone.


STEP 1: COAXIAL TEMNO PLACEMENT

During the pre-procedure set-up, insert the biopsy device through the coaxial introducer needle to verify compatibility. Adjust sliding stopper on coaxial needle to desired calibration level. Lock in place. Insert coaxial and stylet until sliding stopper rests on skin surface. Image coaxial needle to make sure the tip is proximal to the lesion. Remove inner stylet.


STEP 2: SAMPLE NOTCH SIZE

Select sample notch size by turning the white coaxial bolt on the Temno Biopsy Needle from a maximum of 19 mm to a minimum of 9 mm as shown by the indicator window. The amount of sampling notch exposure and needle advancement are measured from the distal end of the coaxial introducer.

STEP 3: TEMNO NEEDLE SPRING CHARGE


Charge Temno needle by pulling back hard on the plunger until a firm click is felt indicating that the needle spring is charged.

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
STEP 4: TEMNO NEEDLE PLACEMENT

With the needle stylet fully retracted, insert the Temno needle into the coaxial introducer and gently connect the coaxial's hub onto the Temno's white coaxial bolt until snug.

STEP 5: CAPTURE BIOPSY SPECIMEN
Advance the stylet forward, with thumb on plunger, to deliver the specimen sampling notch into the area to be biopsied. Release the cutting cannula by fully depressing the plunger with thumb. This will cut a tissue core.

STEP 6: RETRIEVE BIOPSY SPECIMEN
Withdraw Temno needle from coaxial introducer. Pull back on plunger as in Step 3 above until firm click is felt. Carefully advance the stylet forward to expose specimen notch and biopsy core. Remove tissue sample.

STEP 7: OBTAINING MULTIPLE SAMPLES
Additional tissue samples may be obtained by repeating Step 4 through Step 6 above. Note that the distal tip of the coaxial introducer remains proximal to target area.

CAUTION: Test firing the needle may cause damage to the cutting cannula and is not recommended. If test firing is performed, the cutting cannula may be protected by holding the cannula with moderate pressure between the thumb and fingers during firing.

NOTE: The needle stylet will be completely covered by the cutting cannula after charging the device.

CAUTION: For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

WARNING:
- Not intended for use in bone.
- Do not use if package is damaged.

Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
LOT	Lot Number
REF	Catalog Number
STERILE EO	Sterilized using Ethylene Oxide
	Do Not Use if Package is Damaged
2	Single Use
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
MD	Medical Device
Sterile Package	Sterile Package

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
TEMNO®
Coaxial Introducer Needle

INSTRUCTIONS FOR USE

INDICATIONS FOR USE:
Temno® Coaxial Introduer Needles are intended to be used as needle guides in soft tissue biopsy procedures.

NOTE: Not intended for use in bone. Not compatible with MRI.

CAUTION: Use with imaging equipment such as ultrasound, X-ray, CT, etc.

Do not use if package is damaged.

DIRECTIONS FOR USE:
1. During the pre-procedure set-up, insert the biopsy device through the coaxial introducer needle to verify compatibility.
2. Inspect the needle for damage before use.
3. Adjust sliding skin stopper on coaxial introducer needle to desired calibration and lock in place.
4. Insert introducer needle until sliding stopper rests on skin surface. Image coaxial needle, if desired, to make sure the tip is at or near the edge of the region to be sampled.
5. Pull the stylet out of the cannula.
6. Insert a procedural biopsy needle into the introducer needle until the biopsy device touches the hub on the introducer needle. **NOTE:** Automated devices must be prepared to fire prior to insertion.
7. Perform the biopsy.
8. Withdraw the biopsy needle from the introducer needle and remove the specimen.
9. If another biopsy of the same tissue is required, repeat steps 6 and 7. **NOTE:** If changing the angle of the introducer needle is required, reinsert the stylet prior to repositioning the needle. Remove the stylet prior to taking the biopsy.
10. Remove the introducer needle from biopsy site.

CAUTION: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
LOT	Lot Number
REF	Catalog Number
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Use if Package is Damaged
2	Single Use
℞ ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.

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
TEMNO Evolution®

INSTRUCTIONS FOR USE


DEVICE DESCRIPTION
The TEMNO Evolution® biopsy system is a single patient use core biopsy device. It is available in several needle gauge sizes and lengths. The plunger is color coded according to the various gauge sizes (e.g., green=20G, orange=18G, blue=16G, white=14G)."

INDICATIONS FOR USE
The TEMNO Evolution® biopsy device is used to remove, by cutting, a specimen of tissue for microscopic evaluation.


STEP 1. COAXIAL TEMNO PLACEMENT
During the pre-procedure set-up, insert the biopsy device through the coaxial introducer needle to verify compatibility. Adjust sliding stopper on coaxial needle to desired calibration level. Lock in place. Insert coaxial and stylet until sliding stopper rests on skin surface. Image coaxial needle to make sure the tip is on the surface of the lesion. Remove inner stylet.


STEP 2. TEMNO NEEDLE SPRING CHARGE
If using a 10 mm Notch Exposure, apply the spacer onto the TEMNO needle and press fit firmly onto luer taper of TEMNO. Charge TEMNO needle by pulling back on plunger until firm click is felt indicating that the needle spring is charged. It is possible to load the needle in two different positions.
- First Click = 10 mm Notch Exposure (Use Spacer)
- Second Click = 20 mm Notch Exposure (Do not use Spacer)


STEP 3. TEMNO NEEDLE PLACEMENT
With stylet fully retracted, (i.e. specimen notch completely covered by the cannula as in STEP 2) insert TEMNO needle into coaxial needle. Press fit coaxial needle hub onto luer taper of Temno or luer taper of spacer for 10 mm notch.

STEP 4. TEMNO NEEDLE EXECUTION
Press plunger with thumb to advance specimen notch into area to be sampled. Image again to insure that notch is inside the lesion. Press plunger past stop with thumb, automatically triggering the cannula to close over the notch trapping the specimen.

STEP 5. TEMNO NEEDLE EXECUTION
Draw needle out. Pull back on plunger to recharge the spring. Push plunger to first stop exposing the specimen notch. Remove tissue specimen from notch in stylet. Be careful not to push plunger past first stop during this step.

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CAUTION: For single use only. Re-use may result in a non-functional product or contribute to cross contamination. Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

WARNING
- Not intended for use in bone.
- Do not use if package is damaged.

SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
LOT	Lot Number
REF	Catalog Number
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Use if Package is Damaged
	Single Use
℞ ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Caution: Consult accompanying documents

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South Jordan, Utah 84095 U.S.A.