DiaClon ABO/D + Reverse Grouping

English

B001228 v.01 2020-10-07

A, B, DVI-, ctl/A_1, B

Determination of the ABO/Rh blood groups combined with reverse grouping

Product-Identification: 50092

ID-Card ″DiaClon ABO/D + Reverse Grouping″

4 x 12	REF 001234
24 x 12	REF 001237
60 x 12	REF 001236
112 x 12	REF 001235

INTRODUCTION

ABO blood group typing, using anti-A and anti-B test sera, is known as direct or forward grouping test.

Reverse grouping uses red cell reagents of known ABO antigen specificity to indicate the presence or absence of anti-A and anti-B isoagglutinins, the results of which determine the reverse group. Discrepancies between forward and reverse grouping require further investigation [1].


Classification of blood groups must be based on both forward and reverse grouping.

The ID-Card ″DiaClon ABO/D + Reverse Grouping″ allows combined testing of forward and reverse grouping as well as RhD determination.

REAGENTS IVD

ID-Card ″DiaClon ABO/D + Reverse Grouping″ contains monoclonal anti-A [cell line A5], anti-B [cell line G½] and anti-D [cell lines LHM 59 / 20 (LDM3) + 175-2] within the gel matrix. The microtube ctl is the negative control. Two microtubes with ″neutral″ gel serve for reverse grouping with A_1 and B cells.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspensions.
- Test cell Reagents: ID-DiaCell A_1 and B in a 0.8% suspension, in 10 ml vials, ready-to-use.

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When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1,500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

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PREPARATION OF BLOOD SAMPLE

a) Red cell suspension (for ABO/D determination)
Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

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b) Plasma or Serum for reverse grouping
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C, or in accordance with local / national policies / guidelines.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow the test cell reagent to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of ˝ID-DiaCell A_1˝ to microtube 5 (A_1).
4. Pipette 50 µl of ˝ID-DiaCell B˝ to microtube 6 (B).
5. Pipette 50 µl of the patient serum or plasma to both microtubes 5 and 6. An incubation for 10 minutes at room temperature is recommended (see REMARKS, point 4).
6. Pipette 10 or 12.5 µl of the patient's red cell suspension to the microtubes 1–4 (A, B, D, ctl).
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the reactions.

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Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (anti-A, anti-B, anti-AB).In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under ˝REMARKS 1˝.

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C) Reactions for reverse grouping

A_1	B	Blood group
+ to ++++	negative	**B**
negative	+ to ++++	**A**
+ to ++++	+ to ++++	**O**
negative	negative	**AB**

If questionable reactions are obtained, repeat reverse grouping with 4 red cell reagents (A_1, A_2, B and O).

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The anti-D used has been selected so as not to react with DVI variants.

Weak D may give a negative reaction. If all weak / partial D's are required to be detected, all D negative results must be retested. Note that most guidelines do not recommend further testing for weak or partial-D in patients.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under ˝REMARKS 1˝.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/D + Reverse Grouping˝ and the reverse grouping application have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns (only for forward grouping) for which ABO/RH1 phenotypes have been previously determined by a reference method:

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ABO Reverse Grouping

Application	Total number of samples	Sensitivity	Specificity
Reverse A_1	3735	100%	100%
Reverse B	3735	100%	100%

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REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Then proceed as under ˝PREPARATION OF BLOOD SAMPLE˝ and ˝TEST PROCEDURE˝.
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections B, C and D.
 - If the negative control remains positive, the results of the ABO/Rh determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
2. Where a discrepancy occurs between the results of ABO typing and reverse grouping, consult the ˝DiaMed ABO discrepancy chart˝ for appropriate information.
3. Full forward and reverse grouping requires the use of anti-A, -B, -AB, and A_1, A_2, B and O cells. The ID-Card ˝DiaClon ABO/D + Reverse Grouping˝ does not contain anti-AB and allows the use of A_1 and B cells only. It should be used only for confirmation of previously fully tested samples with established blood group status or in accordance with local or national guidelines / recommendations for blood grouping.
4. For reverse grouping, an incubation of at least 10 minutes at 18–25 °C prior to centrifugation will enhance the reactions and minimise repeat testing due to weak isoagglutinins.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the serum or the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein

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1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for in vitro diagnostic medical devices. (2009/886/EC).

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
Reverse grouping uses red cell reagents of known ABO antigen specificity to indicate the presence or absence of anti-A and anti-B isoagglutinins, the results of which determine the reverse group. Discrepancies between forward and reverse grouping require further investigation [1]. Classification of blood groups must be based on both forward and reverse grouping.
The ID-Card ″DiaClon ABO/D + Reverse Grouping for Patients″ allows combined testing of forward and reverse grouping as well as RhD determination.

REAGENTS IVD

ID-Card ″DiaClon ABO/D + Reverse Grouping for Patients″ contains monoclonal anti-A [cell line A5], anti-B [cell line LM306/686 (LB-2)] and anti-D (cell lines LHM 59/20 (LDM3)/175-2] within the gel matrix. The microtube ctl is the negative control. Two microtubes with ″neutral″ gel serve for reverse grouping with A_1 and B cells.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
- Test cell reagents: ID-DiaCell A_1 and B in a 0.8% suspension, in 10 ml vials, ready-to-use.

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1. Dispense 0.5 ml of ID-Diluent 2 into a clean glass tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

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b) Plasma or Serum for reverse grouping
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of ˝ID-DiaCell A1˝ to microtube 5 (A_1).
4. Pipette 50 µl of ˝ID-DiaCell B˝ to microtube 6 (B).
5. Pipette 50 µl of the patient serum or plasma to both microtubes 5 and 6. An incubation for 10 minutes at room temperature is recommended (see ˝REMARKS 4˝).
6. Pipette 10 or 12.5 µl of the patient's red cell suspension to the microtubes 1–4 (A, B, DVI-, ctl).
7. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
8. Read and record the reactions.

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Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (anti-A, anti-B, anti-AB). In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

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C) Reactions for reverse grouping

A1	B	Blood group
+ to ++++	negative	**B**
negative	+ to ++++	**A**
+ to ++++	+ to ++++	**O**
negative	negative	**AB**

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* ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial-D types as appropriate for the category of sample being tested.

The anti-D used has been selected so as **not to react** with DVI variants.

Weak D may give a negative reaction. If all weak D's are required to be detected, all D negative results must be retested. Note that most guidelines do not recommend further testing for weak or partial-D in patients.

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REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Then proceed as under ˝PREPARATION OF BLOOD SAMPLE˝ and ˝TEST PROCEDURE˝.
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections B, C and D.
 - If the negative control remains positive, the results of the ABO/Rh determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
2. Where a discrepancy occurs between the results of ABO typing and reverse grouping, consult the ˝Bio-Rad Quick Reference Chart, Resolving ABO Discrepancies˝ for appropriate information.
3. Full forward and reverse grouping requires the use of anti-A, -B, -AB, and A_1, A_2, B and O cells. The ID-Card ˝DiaClon ABO/D + Reverse Grouping for Patients˝ does not contain anti-AB and allows the use of A_1 and B cells only. It should be used only for confirmation of previously fully tested samples with established blood group status or in accordance with local or national guidelines/recommendations for blood grouping.
4. For reverse grouping, an incubation of at least 10 minutes at 18–25 °C prior to centrifugation will enhance the reactions and minimise repeat testing due to weak isoagglutinins.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/D + Reverse Grouping for Patients˝ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RhD phenotypes have been previously determined by a reference method:

ABO and RhD forward grouping

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: A5	4385	100%	100%
Anti-B, Cell line: LB-2	4559	100%	100%
Anti-D (DVI-), Cell lines: LHM59/20 (LDM3), 175-2	4679*	99.72%**	100%

* These samples include 81 samples with weak/variant expressions of the D antigen.
** 62 weak/variant samples have been detected directly and the 8 DVI variants tested were not detected.

Remark: The Anti-D well is not conceived to detect every weak/variant expressions of the D antigen.

Reverse Grouping

Applications	Total number of samples	Sensitivity	Specificity
A_1	3735	100%	100%
B	3735	100%	100%

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BIBLIOGRAPHY

1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al,: The gel test: A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al,: The gel test: A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/ EC on common technical specifications for *in vitro* diagnostic medical devices (2009/886/EC).

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English B001028 v.01 2020-11-05

A, B, AB, DVI+, ctl, DAT

Determination of the ABO/Rh blood groups with direct antiglobulin test (DAT) for newborns

Product-Identification: 50961

ID-Card ″DiaClon ABO/Rh for Newborns DVI+″

4 x 12	REF 001047
24 x 12	REF 001048
60 x 12	REF 001049
112 x 12	REF 001050

INTRODUCTION

Anti-A, anti-B and anti-AB test sera are necessary to detect the presence or absence of A/B antigens on human red cells. Since A and B antigens are not fully developed at birth, weaker reactions may occur with red cells of newborns than of adults and subgroups often cannot be identified.

The serum of adults contains antibodies directed against the A and B antigens absent from their own red cells. Both antibodies appear after the first 4 to 6 months of life. As a result, reverse grouping is not usually undertaken on newborn blood samples. Confirmation of the newborn's blood group is indicated when the A and B antigen expression is fully developed (2–4 years).

The D antigen as well as weak D is fully developed at birth. The determination of the RhD status of the newborn's blood group is important if the mother is RhD negative.

A direct antiglobulin test (DAT) on newborn blood samples has become a standard procedure, since it is of importance to know if the newborn's red cells have been coated with maternal antibodies in-utero.


REAGENTS IVD

ID-Card ″DiaClon ABO/Rh for Newborns DVI+″ contains monoclonal anti-A [cell line LM 297/628 (LA-2)], anti-B [cell line LM 306/686 (LB-2)], anti-AB [cell lines ES131 (ES-15) + Birma-1 + ES-4] and anti-D (cell lines ESD-1M + 175-2), within the gel matrix.
The microtube (ctl) is the negative control.
The anti-human globulin serum is a blend of rabbit anti-IgG and monoclonal anti-C3d (cell line C139-9).

Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″DiaClon ABO/Rh for Newborns DVI+″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspensions.
(see related package insert)

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SAMPLE MATERIAL

Draw blood samples using acceptable phlebotomy techniques. Cord or heel prick samples may be used. It is normally not necessary to wash the cells before use (see ˝REMARKS˝). Where cord samples are used, care must be taken to avoid contamination with Wharton's jelly.

Preferably, blood samples should be drawn into Citrate, EDTA or CPD-A. For reliable results, use of freshly collected blood is indicated.

PREPARATION OF BLOOD SAMPLE

Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:

Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells or 20 µl of whole blood, mix gently.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to all 6 microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for blood groups ABO

Anti-A	Anti-B	Anti-AB	Blood group
+ to ++++	negative	+ to ++++	**A**
negative	+ to ++++	+ to ++++	**B**
+ to ++++	+ to ++++	+ to ++++	**AB**
negative	negative	negative	**O**

In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

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The anti-D test sera were selected so as **to react** with DVI variants. Some DVI variants can give very weak reaction.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under ″REMARKS 1″.

D) Reactions for the direct antiglobulin test (DAT)

A negative reaction indicates the absence of detectable antibodies on the newborn's red cells.
A positive reaction (± to ++++) indicates that the newborn's red cells are sensitized (red cells coated with antibodies).

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ″DiaClon ABO/Rh for Newborns DVI+″, as well as the DAT application, have been evaluated according to the Common Technical Specification (CTS) [4] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes and DAT have been previously determined by a reference method:

ABO/RH1 grouping

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: LM297/628 (LA-2)	4566	100%	100%
Anti-B, Cell line: LM306/686 (LB-2)	4566	100%	100%
Anti-AB, Cell lines: ES131 (ES-15), Birma-1, ES-4	4737	100%	100%
Anti-D(DVI+), Cell lines: ESD-1M, 175-2	3724*	99.90%**	100%

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DAT

Applications	Total number of samples	Sensitivity	Specificity
Anti-IgG + C3d	700	95.74%	100%

REMARKS

1. The negative control (ctl) must always show a negative reaction. If the negative ctl is positive, proceed as follows:
 - Wash the red cells first with isotonic saline solution (or ID-Diluent 2) before preparing the 0.8% cell suspension.
 - Then proceed as under ″Preparation of blood sample″ and ″Test procedure″.

2. If the negative control remains positive, safe interpretation of the ABO and RhD groups cannot be made and further investigation is required.

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BIBLIOGRAPHY

1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. American Association of Blood Banks: Technical Manual; 11th ed. 1993.
3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. American Association of Blood Banks: Technical Manual; 11th ed. 1993.
3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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INTRODUCTION

According to Mollison [1], the frequencies of the different ABO blood groups in the Caucasian population are approximately as follows:

O 46 %
A 41 %
B 9 %
AB 4 %

To detect the presence or absence of the A/B antigens on red cells, antibodies against the corresponding antigens, anti-A and anti-B are used which can be of human or monoclonal origin. ABO forward typing should not be considered complete without reverse grouping whereby the patient's serum is tested against known A_1, A_2, B and O red cells.

Approximately 85% of the Caucasian population are RhD positive [1]. The expression ″Rh positive″ or ″Rh negative″ is based on the presence or absence of the RhD antigen on the red cells. This may be determined using anti-D test serum which can be of human or monoclonal origin. The sensitivity of the ID-System allows a direct detection of most weak D's.

The ID-Card ″DiaClon ABO/D″ provides the complete profile for ABO/RhD in one single procedure step, including the confirmation of RhD. The first anti-D detects the presence of the D variant DVI, the second anti-D is negative for the DVI variant.
Some DVI variants can give very weak reactions.

REAGENTS IVD

ID-Card ″DiaClon ABO/D″ contains monoclonal anti-A [cell line A5], anti-B [cell line G½], anti-AB [cell lines: ES-131 (ES-15) + Birma-1 + ES-4], anti-D [cell lines ESD-1M + 175-2], anti-D [cell lines LHM59/20 (LDM3) + 175-2] within a gel matrix. The microtube ctl is the negative control.
Preservative: < 0.1% NaN_3.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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* Weaker reactions than +++ may indicate A or B subgroups and further investigations should be performed. See ″Bio-Rad Quick Reference Chart, Resolving ABO Discrepancies″.

In the presence of weak or very weakly expressed antigens the reaction can be negative.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under ″REMARKS 2″.

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* For safe interpretation, it is advised to further investigate the sample using ˝ID-Partial D Typing˝ (REF 001451).
** ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial D types as appropriate for the category of sample being tested.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under ˝REMARKS 2˝.

The reaction strength of both anti-D may vary with weak D's. If the detection of all D weak phenotypes are required, proceed to further tests of negative reactions such as with the ID-Card ˝Anti-D˝ with human antibodies and red cell suspension in ID-Diluent 1.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/D˝ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: A5	4385	100%	100%
Anti-B, Cell line: G1/2	4385	100%	100%
Anti-AB, Cell lines: ES131 (ES-15), Birma-1, ES-4	4737	100%	100%
Anti-D(DVI+), Cell lines: ESD-1M, 175-2	3724*	99.90%	100%
Anti-D(DVI-), Cell lines: LHM59/20(LDM3), 175-2	4679**	99.72%	100%

* These samples include 57 samples with weak (Weak/variant) expression of the RH1 antigen. 54 weak/variant samples have been detected directly. All the DVI variants were detected.
** These samples include 81 samples with weak (Weak/variant) expression of the RH1 antigen. 62 weak/variant samples have been detected directly. The DVI variants were not detected.

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REMARKS

1. The monoclonal anti-B in the ID-Card ˝DiaClon ABO/D˝ has been shown not to react with red cells possessing an acquired B-antigen.
2. If the negative control ctl is positive, wash the red cells first with isotonic saline solution (or ID-Diluent 2) before preparing the red cell suspension. Where the ctl microtube subsequently shows a negative result, the reactions can be interpreted as under B and C.
3. If the ctl remains positive, the results of the ABO/RhD determination should be considered invalid and further investigation is required.
4. Some examples of DVI red cells have weak expression in serological tests, possessing lower numbers of DVI antigen sites than other DVI cells. Reaction strengths for DVI red cells can therefore vary. For safe interpretation of a positive result, it is advised to further investigate the sample using the ID Partial-D Typing kit (REF 001451; see related package insert).

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1. Mollison P. L., Engelfriet C. P. and Contreras M.: Blood Transfusion in Clinical Medicine. 9th ed.1993; Blackwell Scientific Publications, Oxford.
2. Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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English B001334 v.01 2020-09-25

A, B, AB, DVI+, DVI-, ctl

Blood group determination ABO and RhD with human antibodies

Product-Identification: 50031
ID-Card ″ABO/RhD″

4 x 12 REF 001334
24 x 12 REF 001337
60 x 12 REF 001336
112 x 12 REF 001335

INTRODUCTION

According to Mollison [1], the frequencies of the different ABO blood groups in the Caucasian population are as follows:

O..... 46%
A 41%
B 9%
AB..... 4%

To detect the presence or absence of the A/B antigens on red cells, antibodies against the corresponding antigens, anti-A, anti-B and anti-AB are used which can be of human or monoclonal origin. ABO forward typing should not be considered complete without reverse grouping, whereby the patient's serum is tested against known A_1, A_2, B, and O red cells.

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
The expression ″RhD positive″ or ″RhD negative″ is based on the presence or absence of the RhD antigen on the red cells. To detect the presence or absence of the RhD antigen on red cells, anti-D test serum is used which can be of human or monoclonal origin. The sensitivity of the ID-System allows a direct detection of weak D, previously known as D^u. In most countries, it is prescribed that RhD determination is confirmed by a second test, using another reagent.

The ID-Card ″ABO/RhD″ provides the best profile for ABO/RhD in one single procedure step: including the confirmation of RhD. The first anti-D is of human origin and may detect the presence of the D variant DVI. The second anti-D is monoclonal and is negative for the DVI variant.

REAGENTS IVD

ID-Card ″ABO/RhD″ contains a mixture of polyclonal antibodies from human serum and monoclonal anti-A (cell line: A5), polyclonal anti-B, anti-AB and anti-D of human origin and one monoclonal anti-D (cell-line 175-2) within the gel matrix. The microtube ctl is the negative control.
Preservative: < 0.1 % NaN_3.

Caution: the source materials from which these products were manufactured, were found non reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: modified bromelin solution for red cell suspensions.

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1. Dispense 0.5 ml of ID-Diluent 1 into a clean glass tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.
3. Incubate for 10 minutes at room temperature (18–25 °C).

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 µl or 12.5 µl of the red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

INTERPRETATION OF THE RESULTS

A) Principle [2]
Positive: agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: compact button of cells on the bottom of the microtube.

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Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (anti-A, anti-B, anti-AB). In the presence of weak or very weakly expressed antigens the reaction can be negative.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under ˝REMARKS 1˝.

C) Reactions for RhD

Polyclonal anti-D	Monoclonal anti-D	Blood group
+++ to ++++	+++ to ++++	RhD positive
± to ++*	± to ++**	RhD weak
± to ++*	negative	Variant DVI
negative	negative	RhD negative

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The reaction strength of both anti-D may vary with weak D's, see ˝REMARKS 2˝.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under ˝REMARKS 1˝.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the polyclonal antibodies present in the ID-Cards ˝ABO/RhD˝ have been evaluated according to the Common Technical Specification (CTS) [4] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RhD phenotypes have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A	2992	100%	100%
Anti-B	2992	100%	100%
Anti-AB	2992	100%	100%
Anti-D(DVI+)	3049*	99.96%	100%
Anti-D(DVI-), Cell line: 175-2	1065**	99.53%	100%

* These samples include 57 samples with variant (weak/partial) expression of the RH1 antigen. 56 variant samples have been detected directly. All the DVI variants were detected.
** These samples include 57 samples with variant (weak/partial) expression of the RH1 antigen. 50 variant samples have been detected directly. The DVI variants were not detected.

Remark: No Anti-D/RH1 reagent detects all variant (weak/partial) expressions of the RH1 antigen.

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REMARKS

1. The negative control (ctl) must always show a negative reaction. If it is positive, proceed as follows:
 - Wash the red blood cells first with isotonic saline solution (or ID-Diluent 2) before preparing the cell suspension.
 - Then proceed as under ″PREPARATION OF BLOOD SAMPLE″ and ″TEST PROCEDURE″.

 Where the ctl microtube subsequently shows a negative result, the reactions can be interpreted as described in sections B and C. If the ctl remains positive, the results of the ABO/Rh Determination should be considered invalid, and further investigation is required.
2. The D antigen consists of at least 9 epitopes and probably more [3]. Weak D (Dweak) cells possess all known epitopes, but have fewer antigen sites per cell than RhD positive cells. Partial D phenotypes (D variants, D categories), lack one or more epitopes, to which an allo-anti-D may be produced following immunisation by pregnancy or transfusion with RhD positive cells. The polyclonal anti-D used in the ID-Card will react with both weak D and partial D phenotypes, except for R_0^{Har} (RH33) which may not always be detected. The monoclonal anti-D will give a positive result with weak D cells, but will be negative with cells of the DVI category.
3. If used for testing of patients samples, note that most guidelines do not recommend further testing for weak or partial D phenotypes.

LIMITATIONS

a) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either in vivo or *in vitro*, may react with all human antisera. Further investigation of such reactions is required.
b) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions for red cells other than ID-Diluent 1 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant results.

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These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therei

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B001344 v.01 2020-09-25

A, B, DVI-, DVI-/ctl, DAT

Blood grouping ABO and double determination of RhD, combined with direct anti-human globulin test (DAT)

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With the ID-Card ″DiaClon ABO/D + DAT″ a complete profile for ABO forward typing and RhD status can be obtained in one single test procedure, including a direct anti-human globulin test (DAT). The two anti-D antibodies in the ID-card are negative for the variant DVI. The test procedure is thus suitable for patients and not recommended for donors.


REAGENTS IVD

ID-Card ″DiaClon ABO/D + DAT″ contains monoclonal anti-A [cell line A5], anti-B [cell line G½], anti-D [cell lines LHM 50/3 (LDM1) + TH-28 + RUM-1] anti-D [cell lines LHM59/20 (LDM3) + 175-2], within a gel matrix. The microtube ctl is the negative control. The anti-human globulin serum is a blend of rabbit anti-IgG and monoclonal anti-C3d (cell line C139-9).
Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″DiaClon ABO/D + DAT″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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The raw materials of the anti-human globulin serum are of particular lots used only in this product with its indicated lot number.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspensions

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Weaker reactions than +++ may indicate A or B subgroups and further investigations should be performed. See Bio-Rad ABO-Identification chart. In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under ˝REMARKS 2˝.

C) Reactions for RhD

+++ to ++++	± to ++*	negative
RhD positive	RhD weak	negative

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The anti-D test sera were selected so as **not to react** with DVI variants.

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The reaction strength of both anti-D may vary with weak D's. If the detection of all D weak phenotypes are required, proceed to further tests of negative reactions such as with the ID-card ˝Anti-D˝ with human antibodies and red cell suspension in ID-Diluent 1. If used for testing of patients samples, note that most guidelines do not recommend further testing for weak or partial D phenotypes.

D) Reactions for the direct anti-human globulin test (DAT)
A negative reaction indicates absence of detectable IgG and/or C3d antibodies on the red cells.
A positive reaction (± to ++++) indicates that the patient's red cells are sensitised (red cells coated with antibodies and/or C3d).

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/D + DAT˝ and the DAT application have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes and DAT have been previously determined by a reference method:

ABO/RH1 grouping

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: A5	4385	100%	100%
Anti-B, Cell line: G1/2	4385	100%	100%
Anti-D(DVI-), Cell lines: LHM50/3 (LDM1), TH28, RUM-1	3965*	99.69%	100%
Anti-D(DVI-), Cell lines: LHM59/20 (LDM3), 175-2	4679**	99.72%	100%

* These samples include 81 samples with weak (Weak/variant) expression of the RH1 antigen. 63 weak/variant samples have been detected directly. The DVI variants were not detected.

** These samples include 81 samples with weak (Weak/variant) expression of the RH1 antigen. 62 weak/variant samples have been detected directly. The DVI variants were not detected.

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DAT

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-IgG + C3d	700	95.74%	100%

REMARKS

1. The monoclonal anti-B in the ID-Card ˝DiaClon ABO/D + DAT˝ has been shown not to react with red cells possessing an acquired B-antigen.
2. If the negative control ctl is positive, wash the red cells first with isotonic saline solution (or ID-Diluent 2) before preparing the red cell suspension. Where the ctl microtube subsequently shows a negative result, the reactions can be interpreted as under B, C and D.
3. If the ctl remains positive, the results of the ABO/Rh Determination should be considered invalid and further investigation is required.

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BIBLIOGRAPHY

1. Mollison P. L., Engelfriet C. P. and Contreras M.: Blood Transfusion in Clinical Medicine. 9th ed.1993; Blackwell Scientific Publications, Oxford.
2. Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109 –113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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B001037 v.01 2020-11-05

A, B, AB, DVI+, CDE, ctl

Determination of the ABO/Rh Blood groups

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O..... 46%
A..... 41%
B..... 9%
AB..... 4%

To detect the presence or absence of the A / B antigens (ABO1 / ABO2) on red cells, antibodies against the corresponding antigens, anti-A and anti-B are used which can be of human or monoclonal origin. ABO forward typing should not be considered complete without reverse grouping whereby the patient's serum is tested against known A_1, A_2, B and O red cells.

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Current recommendations for RhD typing suggest that for transfusion recipients and antenatal patients, anti-D reagents should not detect the DVI phenotype. The D antigen consists of many epitopes, and in the 9 epitope model, DVI red cells lack all but three. This means that individuals possessing the DVI phenotype may produce an anti-D to the missing epitopes after immunisation by fetal or transfused RhD positive red cells. To ensure that appropriate therapeutic measures are instigated, DVI patient's red cells should be assigned ″Rh negative″ status. Conversely, donor bloods should be tested with anti-D that does detect DVI and assigned ″Rh positive″ status, to avoid the unit being transfused to an RhD negative or partial-D patient.

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
The ID-Card ″DiaClon ABO/Rh for Donors″ offers a complete test profile for forward typing of ABO plus RhD and CDE status in one single test procedure.

REAGENTS IVD

ID-Card ″DiaClon ABO/Rh for Donors″ contains monoclonal anti-A [cell line: A5], anti-B [cell line: G½], anti-AB [cell lines: ES 131 (ES-15), Birma-1, ES-4], anti-D [cell lines: ESD-1M und 175-2] and anti-CDE (cell lines: MS-24, MS-201, MS-26, MS-80) within the gel matrix. The microtube (ctl) is the negative control.
Preservative: < 0.1% NaN_3.

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ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the donor's red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for blood groups ABO

Antisera			Blood group
Anti-A	**Anti-B**	**Anti-AB**	
+++ to ++++	negative	+++ to ++++	**A**
negative	+++ to ++++	+++ to ++++	**B**
+++ to ++++	+++ to ++++	+++ to ++++	**AB**
negative	negative	negative	**O**

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The anti-D test sera were selected so as **to react** with DVI variants. Some DVI variants can give very weak reactions.

If the detection of all weak D phenotypes is required, samples with negative reactions should be retested using ID-Cards ˝Anti-D˝ human (6 x anti-D), see related package insert.

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Important: the microtube ctl must show a negative reaction. If the ctl is positive, the CDE determination is not valid. Repeat the test as described under ˝Remarks˝.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/Rh for Donors˝ have been evaluated according to the Common Technical Specification (CTS) [4] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: A5	4385	100%	100%
Anti-B, Cell line: G1/2	4385	100%	100%
Anti-AB, Cell lines: ES131 (ES-15), Birma-1, ES-4	4737	100%	100%
Anti-D(DVI+), Cell lines: ESD-1M, 175-2	3724*	99.90%**	100%
Anti-CDE, Cell lines: MS-24, MS-26, MS-201, MS-80	2034	100%	100%

* These samples include 57 samples with weak (Weak/variant) expression of the RH1 antigen.
** 54 weak/variant samples have been detected directly. All the DVI variants were detected.

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REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Then proceed as under ″PREPARATION OF BLOOD SAMPLE″ and ″TEST PROCEDURE″.
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections B, C and D.
 - If the negative control remains positive, the results of the ABO/Rh determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
2. Mutations in the blood group gene may give rise to weak or variant forms of an antigen, which may result in unexpectedly weak or negative results. DiaClon anti-C is specific for the C (RH2) antigen and does not react with the C^w (RH8) antigen. The C antigen may be very weak in Cw+ individuals as well as other phenotypes exhibiting depressed C expression. Very weak or negative reactions may therefore be seen in these rare circumstances.

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

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DiaClon ABD-Confirmation for Donors

English B001134 v.01 2020-10-07

A, B, DVI+/A, B, DVI+

Confirmation of blood groups ABO and RhD for donors

Product-Identification: 51051

ID-Card ″DiaClon ABD-Confirmation for Donors″

4 x 12	REF 001134
24 x 12	REF 001133
60 x 12	REF 001136
112 x 12	REF 001135

INTRODUCTION

Compatibility testing must include the examination of previous transfusion records for the recipient's serological history [1]. If previous testing was done, results of current testing must be compared with the interpretation of the previously determined blood group. Concurrence between both test results gives some assurance that there have been no identification errors and that the tests have been performed and interpreted correctly.

Current recommendations for RhD typing suggest that for transfusion recipients and antenatal patients, anti-D reagents should not detect the DVI phenotype. The D antigen consists of many epitopes, and in the 9 epitope model, DVI red cells lack all but three. This means that individuals possessing the DVI phenotype may produce an anti-D to the missing epitopes after immunisation by fetal or transfused RhD positive red cells. To ensure that appropriate therapeutic measures are instigated, DVI patients' red cells should be assigned ″Rh negative″ status. Conversely, donor bloods should be tested with anti-D that does detect DVI and assigned ″Rh positive″ status, to avoid the unit being transfused to an RhD negative or partial-D patient.


The ABD confirmation card can be used for the ABO/RhD blood group control of donors.

REAGENTS IVD

ID-Card ″DiaClon ABD-Confirmation for Donors″ contains monoclonal anti-A (cell line LM297/628 [LA-2]), anti-B (cell line LM306/686 [LB-2]) and anti-D (cell lines ESD-1M, 175-2) within the gel matrix.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspensions.

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1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 10 or 12.5 µl of the donor's red cell suspension of the first sample to the first 3 microtubes of the ID-Card. A second donor's sample can be pipetted into the last 3 microtubes.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

Note: If only one donor is to be tested, remove the aluminium foil from the first 3 microtubes. The remaining microtubes may be used at a later date for a subsequent sample.

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Weaker reactions than +++ may indicate A or B subgroups and futher investigations should be performed. For correct interpretation, a complete grouping test should be performed. In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

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If all weak D's are required to be detected, all D negative results must be retested. For the complete procedure see box insert ID-Cards ˝ABO/Rh˝ or ˝Anti-D˝ with polyclonal reagents of human origin.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABD-Confirmation for Donors˝ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: LM297/628 (LA-2)	4566	100%	100%
Anti-B, Cell line: LM306/686 (LB-2)	4566	100%	100%
Anti-D(DVI+), Cell lines: ESD-1M,175-2	3724*	99.90%**	100%

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LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) ABD-confirmation does not replace the complete ABO/RhD determination. It should only be used as a confirmation of previously determined blood groups where, according to general rules, ABO-typing is confirmed with reverse grouping and ABO/Rh-typing is validated by the presence of a negative control test.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Holland P.V. ed.; Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington,V.A: American Association of Blood Banks, 1984:
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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1. Holland P.V. ed.; Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington,V.A: American Association of Blood Banks, 1984:
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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English

B001025 v.01 2021-02-15

A, B, AB, DVI-, CDE, ctl
Determination of the ABO/Rh blood groups

Product-Identification: 50012
ID-Card ″DiaClon ABO/Rh for Patients″

4 x 12. REF 001044
24 x 12. REF 001043
60 x 12. REF 001046
112 x 12. REF 001045

INTRODUCTION

According to Mollison [1], the frequencies of the different ABO blood groups (ISBT symbol ABO, number 001) in the Caucasian population are as follows:

O..... 46%
A 41%
B 9%
AB..... 4%

Approximately 85% of the Caucasian population are RhD positive [1].

To detect the presence or absence of the A/B antigens (ABO1/ABO2) on red cells, antibodies against the corresponding antigens, anti-A and anti-B are used which can be of human or monoclonal origin. ABO forward typing should not be considered complete without reverse grouping whereby the patient's serum is tested against known A_1, A_2, B and O red cells.

The expression ″Rh positive″ or ″Rh negative″ is based on the presence or absence of the RhD antigen (RH1) on the red cells. This may be determined using anti-D test serum which can be of human or monoclonal origin. The sensitivity of the ID-System allows a direct detection of most examples of the weak D phenotype.

Current recommendations for RhD typing suggest that for transfusion recipients and antenatal patients, anti-D reagents should not detect the DVI phenotype. The D antigen consists of many epitopes, and in the 9 epitope model, DVI red cells lack all but three. This means that individuals possessing the DVI phenotype may produce an anti-D to the missing epitopes after immunisation by fetal or transfused RhD positive red cells. To ensure that appropriate therapeutic measures are instigated, DVI patient's red cells should be assigned ″Rh negative″ status.
Conversely, donor bloods should be tested with anti-D that does detect DVI and assigned ″Rh positive″ status, to avoid the unit being transfused to an RhD negative or partial-D patient.

The Rh antigens C (RH2) and/or E (RH3) may be present on red cells with or without the presence of D (RH1) antigens. The presence of either C and/or E antigens can be indicated by the use of an anti-CDE reagent.


The ID-Card ″DiaClon ABO/Rh for Patients″ offers a complete test profile for forward typing of ABO plus RhD and CDE status in one single test procedure.

REAGENTS IVD

ID-Card ″DiaClon ABO/Rh for Patients″ contains monoclonal anti-A [cell line A5], anti-B [cell line G½], anti-AB [cell lines ES 131 (ES-15), Birma-1, ES-4], anti-D [cell lines LHM 59/20 (LDM 3) and 175-2] and anti-CDE (cell lines MS-24, MS-201, MS-26, MS-80) within the gel matrix. The microtube (ctl) is the negative control.
Preservative: < 0.1% NaN_3.

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STABILITY

ID-Card ″DiaClon ABO/Rh for Patients″ is stable up to the expiry date indicated on the label, with an intact seal of the aluminium foil and stored at 18–25 °C in the position indicated on the outer packaging. Do not store near any heat, air conditioning sources of ventilation outlets. Do not freeze.

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ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspensions. *(See related package insert)*

FURTHER MATERIALS REQUIRED

- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Tubes for suspensions
- ID-Working table
- ID-Centrifuge 6, 12 or 24

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample. Samples that cannot be analysed rapidly should be stored between +2 °C and +8 °C and tested within 48 hours. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used. For automated testing please refer to the user manual of the corresponding instrument.

The reliability of the results depends on correct compliance to the Good Laboratory Practices for reagents and samples.

PREPARATION OF BLOOD SAMPLE

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 25 µl of packed cells, mix gently. Where whole blood from a segment of the blood bag is used directly, add 50 µl of blood cells to the 0.5 ml of ID-Diluent 2.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the patients' red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

INTERPRETATION OF THE RESULTS

A) Principle:
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

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B) Reactions for blood groups ABO

Anti-A	Anti-B	Anti-AB	Blood group
+++ to ++++	negative	+++ to ++++	**A**
negative	+++ to ++++	+++ to ++++	**B**
+++ to ++++	+++ to ++++	+++ to ++++	**AB**
negative	negative	negative	**O**

Weaker reactions than +++ may indicate A or B subgroups and further investigations should be performed. For correct interpretation, a complete grouping test should be performed (forward and reverse grouping). In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the ABO determination is not valid. Repeat the test as described under ˝REMARKS˝.

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Weak D may give a negative reaction. If all weak D's are required to be detected, all D negative results must be retested. The monoclonal anti-D used in the ID-Card ˝DiaClon ABO/Rh for Patients˝ does not react with red cells of the partial DVI phenotype. For the complete procedure see box insert ID-Cards ˝ABO/Rh˝ or ˝Anti-D˝ with polyclonal reagents of human origin. Note that most guidelines do not recommend further testing for weak or partial-D in patients.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the RhD determination is not valid. Repeat the test as described under ˝REMARKS˝.

D) Reactions for CDE
Presence of the antigens RhD, C or E is indicated by positive reactions of +++ to ++++.

Important: the microtube ctl must show a negative reaction. If the ctl is positive, the CDE determination is not valid. Repeat the test as described under ˝REMARKS˝.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABO/Rh for Patients˝ have been evaluated according to the Common Technical Specification (CTS) on reagents used for determining ABO/RH1 blood group systems. Evaluation was performed with samples coming from donors, patients and newborns for whom ABO/RH1 group/phenotypes has been previously determined by a reference method. Total number of tested samples exceeded the CTS requirements.

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A/Anti-B/Anti-AB	3363	100%	100%
Anti-DVI-	3416*	99.7%**	100%
Anti-CDE	1633	100%	100%

* These samples include 54 samples with weak expression of the RH1 antigen.
** 85% of RHW1 (weak D) tested antigens were detected.

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ˝DiaClon ABO/Rh for Patients˝ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

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REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Then proceed as under ″Preparation of blood sample″ and ″Test procedure″.
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections B, C and D.
 - If the negative control remains positive, the results of the ABO/Rh Determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
2. The anti-D test sera were selected so as **not to react** with DVI variants.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

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DiaClon ABD-Confirmation for Patients

English B001254 v.01 2020-10-07

A, B, DVI-/A, B, DVI-
monoclonal
Confirmation of blood groups ABO and RhD

Product-Identification: 50053

ID-Card ″DiaClon ABD-Confirmation for Patients″

4 x 12	REF 001254
24 x 12	REF 001257
60 x 12	REF 001256
112 x 12	REF 001255

INTRODUCTION

Compatibility testing must include the examination of previous transfusion records for the recipient's serological history [1].If previous testing was done, results of current testing must be compared with the interpretation of the previously determined blood group. Concurrence between both test results gives some assurance that there have been no identification errors and that the tests have been performed and interpreted correctly.

Current recommendations for RhD typing suggest that for transfusion recipients and antenatal patients, Anti-D reagents should not detect the DVI phenotype. The D antigen consists of many epitopes, and in the 9 epitope model, DVI red cells lack all but three. This means that individuals possessing the DVI phenotype may produce an Anti-D to the missing epitopes after immunisation by fetal or transfused RhD positive red cells. To ensure that appropriate therapeutic measures are instigated, DVI patients' red cells should be assigned ″Rh negative″ status. Conversely, donor bloods should be tested with Anti-D that does detect DVI and assigned ″Rh positive″ status, to avoid the unit being transfused to an RhD negative or partial-D patient.

The ABD confirmation card can be used for the ABO/RhD blood group control of patients.

REAGENTS IVD

ID-Card ″DiaClon ABD-Confirmation for Patients″ contains monoclonal Anti-A [cell line LM 297 / 628 (LA-2)], Anti-B [cell line LM 306 / 686 (LB-2)] and Anti-D (cell lines TH-28, RUM-1, LDM1] within the gel matrix.

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1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension of the first sample to the first 3 microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

Note: If only one patient is to be tested, remove the aluminium foil from the first 3 microtubes. The remaining microtubes may be used at a later date for a subsequent sample.

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Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (Anti-A, Anti-B, Anti-AB). In the presence of weak or very weakly expressed antigens the reaction can be negative. The Anti-B of monoclonal origin does not react with the acquired B antigen.

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The Anti-D test sera were selected so as **not to react** with DVI variants.

Note that most guidelines do not recommend further testing for weak or partial-D in patients.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon ABD Confirmation for Patients˝ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes have been previously determined by a reference method:

ABO/RH1 grouping

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: LM297/628 (LA-2)	4566	100%	100%
Anti-B, Cell line: LM306/686 (LB-2)	4566	100%	100%
Anti-D(DVI-), Cell lines: TH28,RUM-1,LDM1	3965*	99.69%	100%

* These samples include 81 samples with weak (Weak/variant) expression of the RH1 antigen.
** 63 weak/variant samples have been detected directly. The DVI variants were not detected.

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LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
b) ABD-confirmation does not replace the complete ABO/Rh determination. It should only be used as a confirmation of previously determined blood groups where, according to general rules, ABO-typing is confirmed with reverse grouping and ABO/Rh-typing is validated by the presence of a negative control test.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant results.

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Product-Identification: 50110
ID-Card ″DiaClon Rh-Subgroups + K″

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According to Issitt [1] their frequencies in the Caucasian population are as follows:

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
The ID-Card ″DiaClon Rh-Subgroups + K″ offers a complete profiling of the Rh phenotype and includes K typing.

REAGENTS IVD

ID-Card ″DiaClon Rh-Subgroups + K″ contains monoclonal antibodies anti-C (cell line MS-24), anti-c (cell line MS-33), anti-E (cell line MS-260), anti-e (cell lines MS-16, MS-21, MS-63) and anti-K (cell line MS-56) within each gel matrix. The microtube (ctl) is the negative control.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspensions.

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PREPARATION OF BLOOD SAMPLE

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:

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B) Reactions for Rh phenotypes and K
- A positive reaction (+ to ++++) indicates presence of the corresponding antigen. (see ˝REMARKS 2. + 3.˝)
- Reactions of ≤ 2+ may indicate the presence of weak or variant forms of the antigen.
- A negative reaction indicates absence of the corresponding antigen.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the antigen determination is not valid. Repeat the test as described under ˝REMARKS 1.˝.

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3. Kp^a (K3) antigen in *cis* weakens the expression of K antigens.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
c) Use of suspension solutions for red cells other than ID-Diluent 2 may modify the reactions.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regulary according to GLP procedures.
f) Too heavy or too weak red cell suspensions can cause aberrant results.
g) A pinkish color in the gel may be observed in the anti-e microtube when tested with an e-/RH:-5 quality control sample. This slight coloration may impact the automatic reading of this negative reaction by an instrument. Therefore, as long as the reaction fits with a negative reaction (refer to ″ID-Card – Interpretation Guide″), it can be edited and corrected as negative.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ″DiaClon Rh-Subgroups + K″ have been evaluated according to the Common Technical Specification (CTS) [4] on reagents used for determining RH/KEL1 phenotype.
Evaluation was performed with samples coming from donors, patients and newborns for which RH/KEL1 phenotypes have been previously determined by a reference method. Total number of tested samples exceeded the CTS requirements.

Antibodies	Total number of samples	Sensitivity	Specificity
RH2/RH3/RH4/RH5	1871	100%	100%
KEL1	1871	100%	100%

Reproducibility

Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ″DiaClon Rh-Subgroups + K″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

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English

B002224 v.01 2021-02-15

C, c, E, e, K, ctl
monoclonal
Determination of Rh phenotypes and K

Product-Identification: 50710
ID-Card ˝DiaClon Rh + K Pheno II˝

4 x 12. REF 002224
24 x 12. REF 002227
60 x 12. REF 002226
112 x 12. REF 002225

INTENDED USE

This ID-Card is used for the determination of the C (RH2), c (RH4), E (RH3), e (RH5), K (KEL1) antigens on human red blood cells. For the use on semi or fully automated systems, please refer to the according user manual.

INTRODUCTION

Besides the RhD (RH1) antigen, other important antigens of the Rh system (ISBT, number 004) are: C (RH2), E (RH3), c (RH4) and e (RH5).

According to Reid, M. [5] their frequencies in the Caucasian population are as follows:

C 68%
c..... 80%
E 29%
e..... 98%

Appropriate antigen-positive red cells may stimulate antibody production in antigen C (RH2), c (RH4), E (RH3) and e (RH5) negative individuals. The determination of the Rh phenotypes can, therefore be important during pregnancy, for previously transfused patients and for patients with known irregular antibodies [2].


Approximately 9% of the Caucasian population are K (KEL1, K, K1) positive. The K antigen is strongly immunogenic. Anti-K has been reported as the cause of hemolytic transfusion reactions, both immediate and delayed, and hemolytic disease of the newborn.

The ID-Card ˝DiaClon Rh + K Pheno II˝ offers a complete profiling of the Rh phenotype and includes K typing.

REAGENTS IVD

ID-Card ˝DiaClon Rh + K Pheno II˝ contains monoclonal antibodies anti-C (cell lines MS-273, P3x25513G8), anti-c (cell line 951), anti-E (cell lines MS-80, MS-258), anti-e (cell lines MS-62, P3GD512) and anti-K (cell line AEK4) within each gel matrix. The microtube (ctl) is the negative control.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspensions. (see related package insert)

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1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 25 µl of packed cells or 50 µl of whole blood, mix gently. (see ″LIMITATIONS″, g)

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil. (see ″LIMITATION″, a)

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the patients' red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for Rh phenotypes and K
- A positive reaction (+ to ++++) indicates presence of the corresponding antigen. (see ″REMARKS 2. + 3.″)
- Reactions of ≤ 2+ may indicate the presence of weak or variant forms of the antigen.
- A negative reaction indicates absence of the corresponding antigen.

Important: The microtube ctl must show a negative reaction. If the ctl is positive, the antigen determination is not valid. Repeat the test as described under ″REMARKS 1.″.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ″DiaClon Rh + K Pheno II″ have been evaluated according to the Common Technical Specifications on reagents used for determining RH/KEL1 phenotype.
Evaluation was performed with samples coming from donors, patients and newborns for which RH/KEL1 phenotypes have been previously determined by a reference method. Total number of tested samples exceeded the CTS requirements.

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-RH2 (anti-C), anti-RH4 (anti-c), anti-RH3 (anti-E), anti-RH5 (anti-e)	1210	100%	100%
Anti-KEL1 (anti-K)	1210	100%	100%

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ″DiaClon Rh + K Pheno II″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

REMARKS

1. The negative control must always show a negative reaction.
 - If the negative control is positive, wash the red cells in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension.
 - Proceed as under ″PREPARATION OF BLOOD SAMPLE″ and ″TEST PROCEDURE″.
 - If the negative control subsequently shows a negative result, the reactions can be interpreted as described in sections A and B.
 - If the negative control remains positive, the results of the phenotype determination should be considered invalid and further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.

2. Mutations in the blood group gene may give rise to weak or variant forms of an antigen, which may result in unexpectedly weak or negative results. DiaClon anti-C is specific for the C (RH2) antigen and does not react with the C^w (RH8) antigen. The C antigen may be very weak in Cw+ individuals as well as other phenotypes exhibiting depressed C expression. Very weak or negative reactions may therefore be seen in these rare circumstances.

3. Kp^a (K3) antigen in cis weakens the expression of K antigens.

LIMITATIONS

a) Do not use ID-Cards which show signs of drying, have bubbles, damaged seals or drops of gel in the reaction chamber or on the under side of the aluminium foil. Cards that have drops of antiserum or gel in the upper part of the tube can be used. Cards that shows drops of antiserum in the reaction chamber with no other signs of deterioration can be used, but the cards should be centrifuged before use.
b) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
c) Use of suspension solutions for red cells other than ID-Diluent 2 may modify the reactions.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regulary according to GLP procedures.
f) Too heavy or too weak red cell suspensions can cause aberrant results.
g) The use of whole blood rather than packed cells may give rise to a cell suspension that is too weak and could lead to false negative reactions.

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4 x 12. REF 002134
24 x 12. REF 002137
60 x 12. REF 002136
112 x 12. REF 002135

INTENDED USE

This ID-Card is used for the determination of the C (RH2), C^w (RH8), c (RH4), E (RH3), e (RH5) and K (KEL1) antigens on human red blood cells, for manual use as well as on instruments. For the use on semi or fully automated systems, please refer to the appropriate user manual.

SUMMARY

- Besides the RhD (RH1) antigen, other important antigens of the Rh system (ISBT number 004) are: C (RH2), C^w (RH8), E (RH3), c (RH4) and e (RH5).
- In the Kell System (ISBT number 006), the K (KEL1) antigen is important.

According to Reid, M. [1] their frequencies are as follows:

	Antigen frequency (%)		
	Caucasians	**Blacks**	**Others**
C	68	27	Asians: 93
C^w	2	1	Finns: 4, Latvians: 9
c	80	98	Asians: 47
E	29	22	Asians: 39
e	98	98	Asians: 96
K	9	2	Asians: rare, Iranian Jews: 12, Arabs: 25

- Antigen-positive red cells may stimulate antibody production in antigen C (RH2), C^w (RH8), c (RH4), E (RH3), e (RH5) and K (KEL1) negative individuals.
- The determination of the Rh and K (KEL1) phenotypes can therefore be important during pregnancy, for previously transfused patients and for patients with known irregular antibodies [2].
- The K (KEL1) antigen is strongly immunogenic. Anti-K has been reported as the cause of hemolytic transfusion reactions, both immediate and delayed, and hemolytic diseases of the newborn.
- Mutations in the blood group gene may give rise to weak or variant forms of an antigen, which may result in unexpectedly weak or negative results with anti-C, -c, -E, -e or -K.
- DiaClon anti-C is specific for the C (RH2) antigen and does not react with the C^w (RH8) antigen.
- The C (RH2) antigen may be very weak in C^w+ individuals as well as in other phenotypes exhibiting depressed C (RH2) expression. Very weak or negative reactions may therefore be seen in these rare circumstances.
- Kp^a (K3) antigen in cis weakens the expression of Kell antigens.

PRINCIPLE OF THE TEST

The ID-Card ″DiaClon Rh-Subgroups + C^w + K″ is designed for detection of the presence or absence of the C (RH2), C^w (RH8), c (RH4), E (RH3), e (RH5) and K (KEL1) antigens on red blood cells. The test combines the principles of agglutination and gel filtration. The red cell suspension is distributed into the microtubes and centrifuged. Non agglutinated red cells are collected at the bottom of the microtubes while the agglutinates are dispersed throughout the length of the gel, depending upon their size. Their position in the gel determines the intensity of the reaction.

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REAGENTS

ID-Card ″DiaClon Rh-Subgroups + C^w + K″ contains monoclonal antibodies in the gel matrix:
- Anti-C (cell lines MS-273, P3x25513G8)
- Anti-C^w (cell line MS-110),
- Anti-c (cell line 951)
- Anti-E (cell lines MS-80, MS-258)
- Anti-e (cell lines MS-62, P3GD512)
- Anti-K (cell line AEK4)

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MATERIAL PROVIDED

- ID-Card ″DiaClon Rh-Subgroups + C^w + K″

MATERIAL REQUIRED BUT NOT PROVIDED

- ID-Diluent 2, Id-nº 05761
- ID-Card ″Control Card A″, Id-nº 52020; ctl x 6
- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Centrifuge

STORAGE AND HANDLING

- Store at 18–25 °C.
- Do not store near any heat, air conditioning sources or ventilation outlets.
- Store in an upright position.
- Shelf life: see expiry date on the label.

WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Strict adherence to the procedures and recommended equipment are essential. The equipment should be checked regularly according to the GLP procedures.
- Do not use reagents after the expiration date.
- Do not freeze or expose reagents to excessive heat.
- Used test material must be discarded as hazardous materials. Manage waste according to local, state and/or national regulations.
- All ID-Reagents and test samples must be brought to room temperature (18-25 °C) prior to use.
- Use reagents as provided.
- The bovine albumin used for the production of this reagent is from BSE-free sources.
- Contains sodium azide, which may react with lead or copper plumbing to form explosive azides.
- Do not use ID-Cards which show signs of drying, have bubbles, damaged seals or drops of gel in the reaction chamber or under of the aluminium foil. Cards that have drops of antiserum or gel in the upper part of the tube can be used. Cards that show drops of antiserum in the reaction chamber with no other signs of deterioration can be used, but the cards should be centrifuged before use.

SAMPLES

For optimal results, the determination should be performed using a freshly drawn sample or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into CPD, EDTA or citrate anticoagulant. Samples drawn into plain tubes (no anticoagulant) may be used for manual use only.

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RED BLOOD CELL SUSPENSIONS

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 25 µl of packed cells, mix gently.

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
TEST PROCEDURE

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from all microtubes by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the 5% red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

CONTROL PROCEDURE

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

INTERPRETATION OF THE RESULTS

Reaction Grading Guide

-	A compact pellet of RBCs** with a smooth surface at the bottom of the well with no visible agglutination.
+	A pellet of RBCs at the bottom of the well with agglutinated RBCs visible in the lower half of the gel column.
++	Agglutinated RBCs distributed throughout the entire length of the gel column, with no line of RBCs on the top of the well.
+++	Most agglutinated RBCs concentrated at the top of the gel or upper half of the gel column.
++++	Agglutinated RBCs concentrated as a line on the top of the gel column with a few agglutinated RBCs just underneath the gel surface.
***dp: Double Population**	Agglutinated RBCs as line at the top of the gel or dispersed in upper part of the gel and non-agglutinated RBCs forming a pellet at the bottom of the well.

**RBCs = Red Blood Cells

A) Principle [3]
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells at the bottom of the microtube.

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B) Reactions for ″Rh phenotypes and K″

- A positive reaction (+++ to ++++) indicates presence of the corresponding antigen.
- Positive reactions ≤ ++ with anti-C, anti-c, anti-E, anti-e or anti-K may indicate the presence of weak or variant forms of the antigen.
- In case positive reactions ≤ ++ are obtained with anti-C^w, further investigations must be carried out.
- A negative reaction indicates absence of the corresponding antigen.
- It is strongly recommended to include the ″Control Card A″ (Id-nº 52020) when using the ID-Card ″DiaClon Rh-Subgroups + C^w + K ″ (Id-nº 52000). In any case, if only positive reactions are obtained with all the 6 antisera, the test must be repeated and the ID-Card ″Control Card A″ (Id-nº 52020) processed in parallel according to the ″TEST PROCEDURE″ described above.
 - A negative reaction obtained with the ″Control Card A″ validates the result, obtained by using the ″DiaClon Rh-Subgroups + C^w + K″.
 - In case a positive reaction is obtained with the ″Control Card A″ (Id-nº 52020), the antigen determination is not valid. Repeat the test by washing the cells in warm (37 °C) isotonic saline solution or ID-Diluent 2 before preparing the cell suspension for the ID-Cards ″DiaClon Rh-Subgroups + C^w + K″ and ″Control Card A″. If the control remains positive, the results of the antigen determinations must be considered invalid and further investigations following recommended techniques must be undertaken to ascertain the reason before valid antigen typing can be assured. If the negative control subsequently shows a negative result, the reactions can be interpreted.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Card ″DiaClon Rh-Subgroups + C^w + K″ have been evaluated according to the Common Technical Specifications (CTS) [4] on reagents used for determining RH/KEL1 phenotype.

Evaluation was performed with samples coming from donors, patients and newborns for which RH/KEL1 phenotypes have been previously determined by a reference method. Total number of tested samples exceeded the CTS requirements.

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-RH2 (anti-C), anti-RH4 (anti-c), anti-RH3 (anti-E), anti-RH5 (anti-e)	1290	100%	100%
Anti-KEL1 (anti-K)	1290	100%	100%

Evaluation for Anti-RH8 (Anti-C^w) was performed with samples coming from donors and patients for which phenotypes have been previously determined by a reference method.

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-RH8 (anti-C^w)	340	100%	100%

Reproducibility

Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ″DiaClon Rh-Subgroups + C^w + K″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

LIMITATIONS

Erroneous and abnormal results may be caused by:

a) Fibrin residues in the red cell suspension which may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
b) Use of suspension solutions for red cells other than ID-Diluent 2.
c) Bacterial or other contamination of materials used.
d) Too heavy or too weak red cell suspensions.
e) Patient medication or disease yielding a cross-reaction.
f) Incomplete resuspension of the red blood cells suspension.
g) Contamination between microtubes due to pipetting errors.
h) Use of procedure other than the one described above.
i) Use of cards showing signs of drying, air bubbles or damages foil stripes.
j) Grossly haemolysed, lipemic or icteric blood samples, blood samples with abnormally high concentrations of protein or blood samples from patients who received plasma expanders of high molecular weight.

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BIBLIOGRAPHY

1. Reid M E, Lomas-Francis C, Olsson M.L. THE BLOOD GROUP ANTIGEN FactsBook. 3rd ed. Oxford: Elsevier; 2012
2. Daniels G. Human Blood Groups. publ.: Willey-Blackwell. Third edition, 2013.
3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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REF	Catalog reference
LOT	Batch number
IVD	*In vitro* diagnostic
	Consult instructions for use
	Expiry date (YYYY-MM-DD)
	Storage temperature
	Legal manufacturer
	Consult **downloads.bio-rad.com** to download the latest version of these instructions for use

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*These products may only be used with the following analyzers authorized by Bio-Rad: IH-1000, IH-500, Saxo ID-Reader and Swing TwinSampler.**

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1. Reid M E, Lomas-Francis C, Olsson M.L. THE BLOOD GROUP ANTIGEN FactsBook. 3rd ed. Oxford: Elsevier; 2012.
2. Daniels G. Human Blood Groups. publ.: Willey-Blackwell. Third edition, 2013.
3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/ EC on common technical specifications for *in vitro* diagnostic medical devices (2009/886/EC).

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INTRODUCTION


Transfusions are made with ABO, RhD and compatible blood, in certain cases Rh phenotypes and Kell status are also taken into account. Other blood group antigens are generally not considered. However, when a clinically significant antibody is present, appropriate antigen negative blood should be used. For such cases, the rapid availability of fully typed donor blood is of great advantage. The ID-System facilitates complete antigen profiling.

The ID-Card ″ID-Antigen Profile I″ is suitable for the determination of the antigens P1 (P1), Le(a) (LE1), Le(b) (LE2), Lu(a) (LU1) and Lu(b) (LU2); the ″ID-Antigen Profile II″ for the antigens k (KEL2), Kp(a) (KEL3), Kp(b) (KEL4), Jk(a) (JK1) and Jk(b) (JK2) and the ″ID-Antigen Profile III″ for the antigens M (MNS1), N (MNS2), S (MNS3), s (MNS4), Fy(a) (FY1) and Fy(b) (FY2).

In the ID-Cards ″ID-Antigen Profile I + II″, the gel suspensions contain the corresponding antibodies, requiring only the addition of the red cell suspension (in ID-Diluent 1). In the ID-Card ″ID-Antigen Profile III″, the first 2 microtubes contain neutral gel and the last 4 gel with polyspecific anti-human globulin (AHG). The corresponding antibodies (ID-test sera, specially adapted for the ID-System) are added after the red cell suspension (in ID-Diluent 2).

REAGENTS IVD

- ID-Card ″ID-Antigen Profile III″, 2 microtubes with neutral gel and 4 microtubes with polyspecific anti-human globulin (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix (to be used with the specific ID-test sera). Preservative: < 0.1% NaN_3.

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
Stability: see expiry date on label.

- ID-Card ″ID-Antigen Profile I″, containing monoclonal anti-P1 antibodies (cell line 650), monoclonal anti-Le^a (cell line P3N20V3 + LM112/161), monoclonal anti-Le^b (cell line LM129/181), anti-Lu^a and anti-Lu^b, polyclonal antibodies from human serum, within the gel matrix. The microtube ctl is the negative control. Preservative: < 0.1% NaN_3.
- ID-Card ″ID-Antigen Profile II″, containing anti-k, anti-Kp^a and anti-Kp^b, polyclonal antibodies from human serum, monoclonal anti-Jk^a (cell line MS-15) and monoclonal anti-Jk^b (cell line MS-8), within the gel matrix. The microtube ctl is the negative control. Preservative: < 0.1% NaN_3.
- ID-test sera (for ″ID-Antigen Profile III″), ready-to-use, in 1.4 ml or 5.0 ml vials: anti-M monoclonal [cell line LM110/140 (LM-1)] and anti-N monoclonal (cell line 1422 C 7); anti-S, anti-s, anti-Fy^a and anti-Fy^b, polyclonal antibodies from human serum. Preservatives: the antibiotics trimetoprime and sulfamethoxazole.

Caution: The source materials from which these products were manufactured, were found non reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

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ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: modified bromelin solution for red cell suspension.
- ID-Diluent 2: modified LISS for red cell suspension.

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PREPARATION OF BLOOD SAMPLE

I. ID-Antigen Profile I and II

Prepare a 5% red cell suspension in ID-Diluent 1 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 1 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed red cells, mix gently.
3. Incubate the red cell suspension for 10 minutes at room temperature (18–25 °C).

Use within 15 minutes after incubation.

II. ID-Antigen Profile III

Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed cells or 20 µl of whole blood, mix gently.

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CONTROLS

Known positive and negative samples should be included in accordance to the relevant guidelines of quality assurances.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow ID-Cards to reach room temperature before use.
Allow the reagents to reach room temperature (18–25 °C) before use.

I. ID-Antigen Profile I and II

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 10 or 12.5 µl of the red cell suspension in ID-Diluent 1 to all microtubes.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

II. ID-Antigen Profile III

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension in ID-Diluent 2 to all microtubes.
4. Add 50 µl of the ID-test sera to the appropriate microtubes.
5. Incubate the ID-Card for 10 minutes at room temperature (18–25 °C).
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for Antigen Profiles

- Positive reactions of + to ++++ indicate presence of the corresponding antigen.
- A double population must also be considered as positive. However, a double population may also indicate the presence of cells positive and negative for the corresponding antigen (eg Posttransfusion, if the antigen configuration of the transfused blood was different from the patient's).
- Negative reactions indicate absence of the corresponding antigen.
- Strong ++++ reactions are very rare.
- Anti-Jk^{a} may react more strongly with Jk (a+b-) than with Jk (a+b+) red blood cells (dosage effect).

REMARKS

I. ID-Antigen Profile I and II

- The negative control must always show a negative reaction.
- If the negative control is positive, wash the red cells first in warm isotonic saline solution or ID-Diluent 2, before preparing the red cell suspension. Then repeat the negative control test.
- If the negative control is negative, proceed to the required antigen test.
- If the negative control remains positive, further investigations following recommended techniques should be undertaken to ascertain the reason, before valid antigen typing can be assured.
- ID-Antigen Profile I: Serologically, the P1 phenotype exhibits variation in strength on the red cells, which appears to be under genetic control. The rare dominant inhibitor of Lutheran and other red cell antigens also inhibits the expression of P1.

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II. ID-Antigen Profile III

Prior to testing for the presence of an antigen, it should be assured that the red cells are free from in vivo or *in vitro* coating by autoantibodies and/or complement components, which may react with the AHG, producing falsely positive results. Proceed to the direct antiglobulin test (DAT) as follows:

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Take off the aluminium foil by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension (in ID-Diluent 2) to a microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.
 - If the DAT is negative, proceed to the antigen test.
 - If the DAT is positive, wash the cells with warm isotonic saline solution, ID-Diluent 2 or other recommended techniques, before preparing the cell suspension in ID-Diluent 2. Then repeat the DAT.
 - If the DAT is negative, proceed to the antigen test.
 - If it remains positive, the autoantibody should be eluted following recommended techniques before proceeding to the antigen test.
 - Glycophorin B (GPB) carries the S (MNS3), s (MNS4) antigenic determinants and also expresses the ˝N˝ antigen. Some anti-N sera may therefore cross react with red cells carrying double dose of the M antigen. Unexpected results may also occur when variants of the M and N antigens are present, for example, the rare Mg (MNS11) antigen which is encoded by a gene which produces virtually no M or N antigen. Red cells can therefore be erroneously classified NN or NM when they are genotypically NMg or MMg. Mc (MNS8) is another rare allele of M and N which reacts with the majority of anti-M sera and some anti-N sera. Me (MNS13) is a determinant shared by M and He antigens and some anti-M sera also contain an anti-Me component. These potential sources of discrepant results of MN typing can be of particular importance in parentage testing.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either in vivo or *in vitro* may react with all human antisera.Further investigation of such reactions is required.
c) Certain drugs are known to cause positive antiglobulin reactions.
d) Some pathological conditions are also reported as causing positive antiglobulin tests.
e) Bacterial or other contamination of materials used can cause false positive or false negative results.
f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
h) Use of suspension solutions for red cells other than ID-Diluent 1 and 2 may modify the reactions.
i) Use of other test sera than ID-test sera may modify the reactions.
j) Too heavy or too weak red cell suspensions can cause aberrant results.
k) Lewis antigens may be lost as result of disease. Lewis antigens are often greatly reduced during pregnancy.
l) The very rare Fy^x gene is responsible for the production of extremly weak Fy^b antigens that may give a negative reaction with Fy^b antisera.
m) Some cells which are k (K2) and Kp^a (K3) positive may react weakly to negative with anti-k test sera.

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Product-Identification: 50150
ID-Card ˝DiaClon Anti-H˝ 1 x 12. REF 001911

INTRODUCTION


The H antigen is present on all human red cells, except those with the phenotype O_h (Bombay blood group). H is the precursor substance for A and B antigens. The order of reactivity of anti-H with different ABO-blood groups tends to be:
$O > A_x > A_3 > A_2 > B > A_2B > A_1 > A_1B$.

The ID-Card ˝DiaClon Anti-H˝ consists of 6 microtubes with anti-H, monoclonal, within the gel matrix, providing an easy and standardised test procedure for clear results and safe interpretation.

For the classification of A subgroup status anti-H should be used in conjunction with anti-A, anti-AB and anti-A_1 (available in different configurations on ID-Cards, please see the catalog for ID-System).

REAGENTS IVD

ID-Card ˝DiaClon Anti-H˝ contains 6 microtubes with anti-H monoclonal (cell lines: H-86/50 + H-89/8), within the gel matrix. Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: modified bromelin solution.

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B) Reactions for A subgroups

	Anti-A	Anti-AB	Anti-A_1	Anti-H
A_1	++++	++++	++ to ++++	neg to +
A_{int}	++++	++++	++ to +++	++ to +++
A_2	+++ to ++++	+++ to ++++	neg to +	++ to ++++
A_3	+ to +++*	++ to ++++	neg	++ to ++++
A_x	neg to +	+ to ++	neg	++ to ++++
A_1 B	++++	++++	++ to ++++	neg to +
A_2 B	+++ to ++++	+++ to ++++	neg to +	neg to ++
A_3 B	++ to ++++*	++++	neg	neg to ++
A_x B	neg to +	++++	neg	+ to +++

* Mixed field reactions with anti-A may occasionally be observed, to be interpreted as positive.

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LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the serum or the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 1 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

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Product-Identification: 50140
ID-Card ″Anti-A_1 Absorbed″ 1 x 12. REF 001811

INTRODUCTION

The principle subgroups of A are A_1 and A_2. They are genetically determined via transferases which produce both quantitative and qualitative differences between red cells of the A_1 and A_2 phenotypes. Red cells of both phenotypes react strongly with anti-A.

The serological distinction between the two may be made with anti-A_1 prepared from human group B serum.


The reagent ″Anti-A_1 Absorbed″ agglutinates A_1 cells but reacts negatively or only weakly with A_2 cells. Approximately 80% of group A or AB individuals have red cells that are agglutinated by anti-A_1 and are thus classified as A_1 or A_1B.

For the classification of A subgroup status, anti-A_1 should be used in conjunction with anti-A, anti-AB and anti-H (available in different configurations on ID-Cards, please see the catalog for ID-System).

REAGENTS IVD

ID-Card ″Anti-A_1 Absorbed″ with 6 microtubes containing polyclonal antibodies Anti-A_1 from human serum within the gel matrix. Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: Modified Bromelin solution for red cell suspension.

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1. Dispense 0.5 ml of ID-Diluent 1 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

Incubate the red cell suspension for 10 minutes at room temperature.

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B) Reactions for A subgroups antigens

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* Double populations may occasionally be observed, and should be interpreted as positive.

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English

B001711 v.01 2021-02-15

Product-Identification: 52020
ID-Card ″Control Card A″

1 x 12 REF 001711
4 x 12 REF 001714

INTENDED USE

This ID-Card is a negative control to be used with the below mentioned ID-Cards for manual use as well as on instruments; see corresponding box inserts:

- ″DiaClon Anti-C^w″, Id-n° 52010, B001321.
- ″DiaClon Rh-Subgroups + C^w + K″, Id-n° 52000, B002134.

For the use on semi or fully automated systems, please refer to the appropriate user manual.

SUMMARY

The ID-Card ″Control Card A″ is designed to validate the reactions obtained with the above mentioned ID-Cards on red blood cells.

PRINCIPLE OF THE TEST

The test combines the principles of agglutination and gel filtration. The red cell suspension is distributed into the microtubes and centrifuged. Non agglutinated red cells are collected at the bottom of the microtube while the agglutinates are dispersed throughout the length of the gel, depending upon their size. Their position in the gel determines the intensity of the reaction.

REAGENTS

ID-Card ″Control Card A″ contains 6 microtubes with antibody free gel matrix.

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WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Strict adherence to the procedures and recommended equipment essential. The equipment should be checked regularly according to the GLP procedures.
- Do not use reagents after the expiration date.
- Do not freeze or expose reagents to excessive heat.
- Used test material must be discarded as hazardous materials. Manage waste according to local, state and/or national regulations.
- All ID-Reagents and test samples must be brought to room temperature (18–25°C) prior to use.
- Use reagents as provided.
- The bovine albumin used for the production of this reagent is from BSE-free sources.
- Contains sodium azide, which may react with lead or copper plumbing to form explosive azides.
- Do not use ID-Cards which show signs of drying, have bubbles, damaged seals or drops of gel in the reaction chamber or under of the aluminium foil. Cards that have drops of antiserum or gel in the upper part of the tube can be used. Cards that show drops of antiserum in the reaction chamber with no other signs of deterioration can be used, but the cards should be centrifuged before use.

SAMPLES

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into CPD, EDTA or citrate anticoagulant. Samples drawn into plain tubes (no anticoagulant) may be used for manual use only.

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1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 25 µl of packed cells, mix gently.

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
TEST PROCEDURE

1. Identify the ID-Card with the unique patient or donor numbers / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the 5% red cell suspension to the appropriate microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

Unused wells can be utilized up to six centrifugations if foil seal is intact.

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INTERPRETATION OF THE RESULTS

Reaction Grading Guide

-	A compact pellet of RBCs** with a smooth surface at the bottom of the well with no visible agglutination.
+	A pellet of RBCs at the bottom of the well with agglutinated RBCs visible in the lower half of the gel column.
++	Agglutinated RBCs distributed throughout the entire length of the gel column, with no line of RBCs on the top of the well.
+++	Most agglutinated RBCs concentrated at the top of the gel or upper half of the gel column.
++++	Agglutinated RBCs concentrated as a line on the top of the gel column with a few agglutinated RBCs just underneath the gel surface.
*dp: Double Population	Agglutinated RBCs as line at the top of the gel or dispersed in upper part of the gel and non-agglutinated RBCs forming a pellet at the bottom of the well.

**RBCs = Red Blood Cells

A) Principle [1]
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells at the bottom of the microtube.

B) Reactions for ˝Control Card A˝
- A negative reaction validates the results obtained with the same samples using the ID-Card ˝DiaClon Rh-Subgroups + + + K˝, Id-n° 52000, and/or ˝DiaClon Anti-C^w˝, Id-n° 52010.
- A positive reaction with the ˝Control Card A˝ invalidates the results obtained on the same samples using the ID-Card ˝DiaClon Rh-Subgroups + C^w + K˝, Id-n° 52000, and/or ˝DiaClon Anti-C^w˝, Id-n° 52010.
- In case a positive control result is obtained, the cells must be washed in warm (37 °C) isotonic saline solution or ID-Diluent 2 before preparing the cell suspensions for the ID-Cards tested, ˝DiaClon Rh-Subgroups + C^w + K˝ or ˝DiaClon Anti-C^w˝ in addition to the ˝Control Card A˝ and then the tests must be repeated. If the control remains positive, the results of the antigen determinations must be considered invalid and further investigations following recommended techniques must be undertaken to ascertain the reason before valid antigen typing can be assured. If the negative control subsequently shows a negative result, the reactions can be interpreted.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the ID-Card ˝Control Card A˝ have been evaluated internally. Evaluation was performed with donor and patient samples which had reacted negatively with a ctl reagent in a reference method.

Test	Sensitivity		Specificity	
Control Card A	N=0	N/A	N=340	100%

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-card ˝Control Card A˝ have been evaluated internally. Neither false positive nor false negative results were observed.

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LIMITATIONS

Erroneous and abnormal results may be caused by:
a) Fibrin residues in the red cell suspension which may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are at the bottom of the microtube after centrifugation.
b) Use of suspension solutions for red cells other than ID-Diluent 2.
c) Bacterial or other contamination of materials used.
d) Too heavy or too weak red cell suspensions.
e) Patient medication or disease yielding a cross-reaction.
f) Incomplete resuspension of the red blood cells suspension.
g) Contamination between microtubes due to pipetting errors.
h) Use of procedure other than the one described above.
i) Use of cards showing signs of drying, air bubbles or damaged foil stripes.
j) Grossly haemolysed, lipemic or icteric blood samples, blood samples with abnormally high concentrations of protein or blood samples from patients who received plasma expanders of high molecular weight.

BIBLIOGRAPHY

1. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
2. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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INTENDED USE

This ID-Card is used for the determination of the C^{w} antigen (RH8) on human red blood cells, for manual use as well as on instruments. For the use on semi or fully automated systems, please refer to the appropriate user manual.

SUMMARY

C^{w} (RH8) is a relatively low frequency antigen within the Rh system (ISBT, number 004).
According to Reid, M. [1] the frequency is as follows:

Antigen frequency (%)			
	Caucasians	**Blacks**	**Others**
C^{w} (RH8)	2	1	Finns: 4, Latvians: 9

Antigen-positive red cells may stimulate antibody production in C^{w} (RH8) negative individuals, and once formed, the antibody may cause destruction of corresponding red cells.

The C^{w} antigen (RH8) occurs most often with a *DCe* gene complex [2].

PRINCIPLE OF THE TEST

The ID-Card ″DiaClon Anti-C^{w}″ is designed for the detection of the presence or absence of the C^{w} antigen (RH8) on red blood cells. The test combines the principles of agglutination and gel filtration. The red cell suspension is distributed into the microtubes and centrifuged. Non agglutinated red cells are collected at the bottom of the microtube while the agglutinates are dispersed throughout the length of the gel, depending upon their size. Their position in the gel determines the intensity of the reaction.

REAGENTS

ID-Card ″DiaClon Anti-C^{w}″ contains 6 microtubes with monoclonal anti-C^{w} (cell line MS-110) in the gel matrix.

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WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Strict adherence to the procedures and recommended equipment are essential. The equipment should be checked regularly according to the GLP procedures.
- Do not use reagents after the expiration date.
- Do not freeze or expose reagents to excessive heat.
- Used test material must be discarded as hazardous materials. Manage waste according to local, state and/or national regulations.
- All ID-Reagents and test samples must be brought to room temperature (18–25 °C) prior to use.
- Use reagents as provided.
- The bovine albumin used for the production of this reagent is from BSE-free sources.
- Contains sodium azide, which may react with lead or copper plumbing to form explosive azides.
- Do not use ID-Cards which show signs of drying, have bubbles, damaged seals or drops of gel in the reaction chamber or under of the aluminium foil. Cards that have drops of antiserum or gel in the upper part of the tube can be used. Cards that show drops of antiserum in the reaction chamber with no other signs of deterioration can be used, but the cards should be centrifuged before use.

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TEST PROCEDURE

1. Identify the ID-Card with the unique patient or donor numbers / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the 5% red cell suspension to the appropriate microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 Minutes in the ID-Centrifuge.
5. Read and record the results.


Unused wells can be utilized up to 6 centrifugations if foil seal is intact.

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B) Reactions for ˝Anti-C^{w}˝
- A positive reaction (+++ to ++++) indicates presence of the corresponding antigen.
- In case positive reactions ≤ ++ are obtained, further investigations must be carried out.
- A negative reaction indicates absence of the corresponding antigen.
- In case a positive reaction is obtained, the test should be repeated and the ID-Card ˝Control Card A˝ (Id-n° 52020) processed in parallel according to the TEST PROCEDURE described above.
- A negative reaction obtained with the ˝Control Card A˝ validates the result obtained by using the ˝DiaClon Anti-C^{w}˝.
- In case a positive reaction is obtained with the ˝Control Card A˝, the antigen determination is not valid. Repeat the test by washing the cells in warm (37 °C) isotonic saline solution or ID-Diluent 2 before preparing the cell suspension for the ID-Cards ˝DiaClon Anti-C^{w}˝ and ˝Control Card A˝. If the control remains positive, the results of the C^{w} (RH8) determination must be considered invalid and further investigations following recommended techniques must be undertaken to ascertain the reason before valid antigen typing can be assured. If the negative control subsequently shows a negative result, the reactions can be interpreted.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performance of the monoclonal antibody present in the ID-Card ˝DiaClon Anti-C^{w}˝ has been evaluated in conformity with the application field and comparable to other approved products.
Evaluation for Anti-RH8 (Anti-C^{w}) was performed with samples coming from donors and patients for which phenotypes have been previously determined by a reference method.

Antibody	Total number of samples	Sensitivity	Specificity
Anti-RH8 (Anti-C^{w})	340	100%	100%

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ˝DiaClon Anti-C^{w}˝ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

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English

B001211 v.01 2020-09-25

human
Determination of the RhD and weak D antigen

Product-Identification: 50160
ID-Card ″Anti-D″

1 x 12 REF 001211
24 x 12 REF 001213

INTRODUCTION

According to Mollison [1], approximately 85% of the Caucasian population are RhD positive.


The expression ″RhD positive″ or ″RhD negative″ is based on the presence or absence of the RhD antigen on the red cells. To detect the presence or absence of the RhD antigen on red cells, anti-D test serum is used which can be of human or monoclonal origin.

The ID-Card ″Anti-D″ allows the testing of 6 samples. It is especially useful for confirmation of RhD determination, verification of negative reactions and when doubtful results with other methods (such as microplates, automated procedures) have been observed.

REAGENTS IVD

ID-Card ″Anti-D″ contains polyclonal anti-D of human origin within the gel matrix.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 1: modified bromelin solution for red cell suspensions.

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1. Dispense 0.5 ml of ID-Diluent 1 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.
3. Incubate the red cell suspension for 10 minutes at room temperature (18–25 °C).

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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Most D variant phenotypes will be detected with the polyclonal anti-D, except for R_0^{Har} (RH33) which may not always be detected. D variant phenotypes include weak expression of full D antigen and partial D types, which lack epitopes associated with the D antigen. Some DVI variants can give very weak reactions. Weak D may give a negative reaction. If all weak D's are required to be detected, all D negative results must be retested (using ˝ID-DiaClon Anti-D˝ for verification of weak D by IAT).

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the polyclonal antibody present in the ID-Cards ″Anti-D″ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which RHD phenotype have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-D	3049*	99.96%	100%

* These samples include 57 samples with variant (weak/partial) expression of the RH1 antigen. 56 variant samples have been detected directly.

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REMARKS

Prior to testing for the presence of antigens, it must be assured that the red cells are free of enzyme reactive autoantibodies.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 1 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

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ЛИТЕРАТУРА
1. Mollison, P.L., Engelfriet, C.P. and Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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B001251 v.01 2020-09-25

Determination of the RhD

Product-Identification: 50171
ID-Card ″DiaClon Anti-DVI neg.″

1 x 12 REF 001251
24 x 12 REF 001252

INTRODUCTION

The expression ″RhD positive″ or ″RhD negative″ is based on the presence or absence of the RhD antigen on the red cells. According to Mollison [1], approximately 85% of the Caucasian population are RhD positive. The ″DiaClon Anti-DVI neg.″ is a blend of monoclonal antibodies.

The ID-Card ″DiaClon Anti-DVI neg.″ allows the testing of 6 samples. It is especially useful for confirmation of RhD determination, verification of negative reactions and when doubtful results with other methods (such as microplates, automated procedures) have been observed.

REAGENTS IVD

ID-Card ″DiaClon Anti-DVI neg.″ with 6 microtubes containing a blend of monoclonal anti-D (cell lines TH-28, MS-26, 175-2) within the gel matrix.
Preservative: < 0.1% NaN_3.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Card ″DiaClon Anti-DVI neg.″ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which RhD phenotypes have been previously determined by a reference method:

RhD grouping

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-D (DVI-) Cell lines: TH-28, MS-26, 175-2	1131*	99.27% (953/960)	100% (171/171)

* These samples include 57 samples with variant (weak/partial) expression of the RH1 antigen. 47 variant samples have been detected directly. The DVI variants were not detected.

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REMARKS

1. Most examples of weak D will be detected directly with the ID-Card ″DiaClon Anti-DVI neg.″.
2. When the detection of very weak D is required, samples with negative reactions should be retested using the ID-Card ″Anti-D″ (human) or ″DiaClon Anti-D″ (monoclonal) for confirmation of weak D by IAT (see related package insert). For safe interpretation of a positive result, it is advised to further investigate the sample using the ID Partial-D Typing kit (REF 001451; see related package insert).
3. Negative control:
 Where the ID-Card ″DiaClon Anti-DVI neg.″ is used as a confirmation following the use of other monoclonal ID-Cards containing a negative control, and where this latter was found to be negative, no further control is required. In all other cases a negative control has to be performed.

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Mollison P. L., Engelfried C. P., and Contreras, M.: Blood Transfusion in Clinical Medicine. 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109 –113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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1. Mollison P. L., Engelfried C. P., and Contreras, M.: Blood Transfusion in Clinical Medicine. 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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Product-Identification: 51180
ID-Card ″DiaClon Anti-P_1″ 1 x 12 REF 007212

INTRODUCTION

The P blood group system (ISBT, number 003) has the symbol P1 and contains only one antigen, P1 [1]. P1 positive cells are designated P_1, and P1 negative cells, P_2. P and Pk are within the globoside collection; P (ISBT, number 209) is GLOB1 and Pk is GLOB2. However there are serological relationships between all these antigens. The symbol p is used to denote the ″null″ phenotype.

The phenotypes and frequencies in the caucasian population are:


Phenotype	Antigens	Frequency
P_1	PP1Pk	79%
P_2	PPk	21%
p	-	very rare
P_1^k	P1Pk	very rare
P_2^k	Pk	very rare

Serologically, the P_1 phenotype exhibits variation in strength on the red cells, which appears to be under genetic control. The rare dominant inhibitor of Lutheran and other red cell antigens also inhibits the expression of P1.

The ID-Card ″DiaClon Anti-P_1″ allow testing of 6 samples simultaneously.

REAGENTS IVD

ID-Card ″DiaClon Anti-P_1″ with 6 microtubes containing monoclonal anti-P_1 antibodies (cell line 650) within the gel matrix. Preservative: < 0.1% NaN_3.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension (mix gently before use) to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for the determination of the antigen P1
- Positive reactions of + to ++++ indicate presence of the corresponding antigen. Strong ++++ reactions are very rare.
- Negative reactions indicate absence of the corresponding antigen.

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LIMITATIONS

a) ID cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions for red cells other than ID-Diluent 1 may modify the reactions.
f) Too heavy red cell suspensions can cause aberrant results.

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ID-Card Fy^{a}, ID-Card Fy^{b}

English B007270 v.01 2021-02-15

Determination of the Fy^{a} (FY1) and Fy^{b} (FY2) antigens with the ID-Card Fy^{a}, ID-Card Fy^{b}, test serum ID-Anti-Fy^{a} and test serum ID-Anti-Fy^{b}

Product-Identification:

ID-Card ″ID-Card Fy^{a}″	Id-n°: **50350**	6 x Fy^{a}	1 x 12	REF 007270
Test serum ″ID-Anti-Fy^{a}″	Id-n°: **09210**	vial	1 x 3.6 ml	REF 007271
Test serum ″ID-Anti-Fy^{a}″	Id-n°: **09210**	vial	1 x 5.0 ml *	REF 007272
ID-Card ″ID-Card Fy^{b}″	Id-n°: **50360**	6 x Fy^{b}	1 x 12	REF 007280
Test serum ″ID-Anti-Fy^{b}″	Id-n°: **09310**	vial	1 x 3.6 ml	REF 007281
Test serum ″ID-Anti-Fy^{b}″	Id-n°: **09310**	vial	1 x 5.0 ml *	REF 007282

* Ideal for automation

INTRODUCTION

According to Mollison [1], the phenotype frequency of the Duffy blood groups in the Caucasian population is as follows:

Fy a+ b-± 18%
Fy a+ b+.....± 49%
Fy a- b+± 33%

Both Duffy antigens are fully developed at birth. Fy^{a} frequently causes immunisation while Fy^{b} is rarely the cause of immunisation.

Six genes are described in the Duffy system, (ISBT symbol FY, number 008) responsible for the production of the antigens Fy^{a}, Fy^{b}, Fy3, Fy5, Fy6.

Fy (a-b-) is the result of homozygosity for an apparently silent gene, Fy, and may determine the Fy4 antigen. Fy^{x} is responsible for the production of very weak Fy^{b} antigens that are only detected if strong anti-Fy^{b} antisera are used [2].


Both Fy^{a} and Fy^{b} are destroyed by enzyme treatment of the red cells and disappear when stored at low pH.

Anti-Fy^{a} and anti-Fy^{b} are practically always of immune origin. They both cause hemolytic transfusion reactions and can, although not frequently, cause haemolytic disease of the newborn (HDN).

The ID-Cards ″ID-Card Fy^{a}″ and ″ID-Card Fy^{b}″ allow testing of 6 samples simultaneously.

REAGENTS IVD

- ID-Card ″ID-Card Fy^{a}″ with 6 microtubes containing polyspecific anti-human globulin serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within the gel matrix.
 Preservative: < 0.1% NaN_3.
- ID-Card ″ID-Card Fy^{b}″ with 6 microtubes containing polyspecific anti-human globulin serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within the gel matrix.
 Preservative: < 0.1% NaN_3.

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
- Test serum "ID-Anti-Fy^{a}", vial containing 3.6 ml and 5 ml anti-Fy^{a}, polyclonal antibodies from human serum. Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

- Test serum "ID-Anti-Fy^{b}", vial containing 3.6 ml and 5 ml anti-Fy^{b}, polyclonal antibodies from human serum. Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow the reagent to reach room temperature (18–25 °C) before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension (mix gently before use) to the appropriate microtube.
4. Add 50 µl of the corresponding ID-test serum.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for Fy^{a} and Fy^{b}
Positive reactions of + to +++ indicate the presence of the corresponding antigen.
Strong ++++ positive reactions are rare.
Negative reactions indicate the absence of the corresponding antigen.

REMARKS

1. Cells that are direct antiglobulin test (DAT) positive due to *in vivo* coating with IgG will be agglutinated by the AHG within the gel matrix. A DAT should be performed on the sample prior to antigen typing to eliminate this possibility.
2. Repeated centrifugation of the ID-Cards does not affect the performance of the product.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either *in vivo* or *in vitro* may react with all human antisera. Further investigation of such reactions is required.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
g) Use of other test sera than ID-test sera may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.
i) The very rare Fy^{x} gene is responsible for the production of extremly weak Fy^{b} antigens that may give a negative reaction with Fy^{b} antisera.

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Determination of the S (MNS3) and s (MNS4) antigens with the ID-Card S, ID-Card s, test serum ID-Anti-S and test serum ID-Anti-s

Product-Identification:

ID-Card ″S″	Id-n°: **50330**	6 x S	1 x 12 REF 007130
Test serum ″ID-Anti-S″	Id-n°: **09010**	vial	1 x 3.6 ml REF 007131
Test serum ″ID-Anti-S″	Id-n°: **09010**	vial	1 x 5.0 ml * REF 007132
ID-Card ″s″	Id-n°: **50340**	6 x S	1 x 12 REF 007140
Test serum ″ID-Anti-s″	Id-n°: **09110**	vial	1 x 3.6 ml REF 007141
Test serum ″ID-Anti-s″	Id-n°: **09110**	vial	1 x 5.0 ml * REF 007142

* Ideal for automation

INTRODUCTION

The antigens S (MNS3) and s (MNS4) are carried on the human glycophrin B (GPB) molecule and are part of the MNS system (ISBT number 002).

According to Daniels [1], the phenotype frequency of the antithetical antigens S and s in the Caucasian population is as follows:


SS..... 11%
Ss 44%
ss 45%

Ss antibodies are not frequently encountered, however they have both been implicated in severe hemolytic disease of the newborn and hemolytic transfusion reactions, and are therefore clinically important. Antigen typing of paternal or fetal neonatal or donor bloods may also therefore be appropriate, according to the clinical situation.

The ID-Cards ″S″ and ″s″ allow testing of 6 samples simultaneously to determine the S and s phenotype of donor blood prior to crossmatching, or for confirmation of S and s phenotype of patient blood after antibody identification.

REAGENTS IVD

- ID-Card ″S″ with 6 microtubes containing polyspecific anti-human globulin serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ″s″ with 6 microtubes containing polyspecific anti-human globulin serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

- Test serum ″ID-Anti-S″. Vial containing 3.6 ml and 5 ml anti-S, polyclonal antibodies from human serum. Preservatives: the antibiotics trimethoprim and sulfamethoxazole.
- Test serum ″ID-Anti-s″. Vial containing 3.6 ml and 5 ml anti-s, polyclonal antibodies from human serum. Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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Allow the reagents to reach room temperature (18–25 °C) before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension (mix gently before use) in the appropriate microtube.
4. Add 50 µl of the corresponding ID-test serum.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for S and s antigens
- Positive reactions of + to ++++ indicate presence of the corresponding antigen.
- Negative reactions indicate absence of the corresponding antigen.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the polyclonal antibodies present in ID-Cards ″S″ and ″s″, Test sera ″ID-Anti-S″ and ″ID-Anti-s″ are verified by tests with positive and negative samples to ensure conformity with product specifications. The results obtained with the ID-Cards ″S″ and ″s″, Test sera ″ID-Anti-S″ and ″ID-Anti-s″ were compared directly with an established device already CE marked:

Test	Total number of samples	Sensitivity	Specificity
Anti-S Polyclonal, human	207	100%	100%
Anti-s Polyclonal, human	207	100%	100%

Results obtained showed that performance characteristics of each test of the tests ″Anti-S″ and ″Anti-s″ are in conformity with the application field and comparable to other approved products.

REMARKS

1. Cells that are direct antiglobulin test (DAT) positive due to in vivo coating with IgG will be agglutinated by the AHG within the gel matrix. A DAT should be performed on the sample prior to antigen typing to eliminate this possibility.
2. Repeated centrifugation of the ID-Card does not affect the performance of the product.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either in vivo or *in vitro* may react with all human antisera. Further investigation of such reactions is required.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
g) Use of other test sera than ID-test sera may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

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English B002121 v.01 2020-11-05

monoclonal
Determination of the K Antigen (KEL1)

Product-Identification: 50200
ID-Card ″DiaClon Anti-K″ 1 x 12 REF 002121

INTRODUCTION

The Kell system (ISBT number 006) is currently known to contain 22 antigens, numbered K1 to K25 (K8, K9 and K15 are now obsolete).

According to Mollison [1], the frequencies of the antithetical antigens K (KEL1) and k (KEL2) in the Caucasian population are as follows:

KK..... 0.2 %
Kk 8.8 %
kk..... 91.0 %

The K antigen is very immunogenic and therefore antibodies against the antigen are regularly found. They have been associated both with transfusion reactions and hemolytic disease of the newborn (HDN).


The clinical importance of anti-K has resulted in the systematic determination of the K antigen in both donors and patients. The ID-Card ″DiaClon Anti-K″ can be used to determine the K antigen status of donor blood prior to compatibility testing, for confirmation of the K antigen of patient blood after antibody identification, or for paternal/fetal antigen typing in allo-immunised pregnancies.

The ID-Card ″DiaClon Anti-K″ allows testing of six samples simultaneously.

REAGENTS IVD

ID-Card ″DiaClon Anti-K″ with 6 microtubes containing anti-K monoclonal (cell line MS-56) within the gel matrix.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspension.
(see related package insert)

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Take off the aluminium foil from as many microtubes as needed by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to all microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for Kell

- A positive reaction of + to ++++ indicates the presence of the corresponding antigen.
- A negative reaction indicates the absence of the corresponding antigen.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the monoclonal antibodies present in the ID-Cards ″DiaClon Anti-K″ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for which KEL1 phenotype have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-K (Anti-KEL1)	N=1498	100%	100%

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REMARKS

1. IgG coating of the red cells (positive DAT), in the absence of an autoagglutinin reacting in saline will normally not give rise to a false positive reaction with monoclonal reagents.
2. Certain disease states, e.g. multiple myeloma and cold agglutinin disease (CAD) are known to give rise to spontaneous aggregation of red cells. In most cases, washing cells with warm isotonic saline solution prior to testing resolves the problem.
3. Kp^a (K3) antigen in cis weakens the expression of K antigens.

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BIBLIOGRAPHY

1. Mollison P. L., Engelfriet C. P. and Contreras M.: Blood Transfusion in Clinical Medicine, 9th ed. 1993; Blackwell Scientific Publications, Oxford.
2. Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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Anti-K and Anti-k

English B002111 v.01 2020-09-25

human
Determination of the K- (KEL1) and k- (KEL2) Antigens

Product-Identification:

ID-Card ″Anti-K″	Id-n°: **50190**	CE 0123	1 x 12 REF 002111
ID-Card ″Anti-k″	Id-n°: **50260**	CE	1 x 12 REF 007251

INTRODUCTION

The KEL (Kell) system (ISBT number 006) is currently known to contain 25 antigens, numbered K1 to K28 (K8, K9 and K15 are now obsolete).

According to Mollison [1], the phenotype frequencies of the antithetical antigens K (KEL1) and k (KEL2) in the Caucasian population are as follows:

KK..... 0.2%
Kk 8.8%
kk 91.0%

The antibodies of the Kell system have been associated with transfusion reactions and hemolytic disease of the newborn (HDN).


The clinical importance of anti-K has resulted in the systematic determination of the K and Cellano antigens in both donors and patients. The ID-Cards ″Anti-K″ and ″Anti-k″ can be used to determine the antigen status of donor blood prior to crossmatching, for confirmation of the antigen after antibody identification, or for paternal/fetal antigen typing in allo-immunised pregnancies.

The ID-Cards ″Anti-K″ and ″Anti-k″ allow testing of six samples simultaneously.

REAGENTS IVD

- ID-Card ″Anti-K″ with 6 microtubes containing polyclonal anti-K antibodies from human serum within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ″Anti-k″ with 6 microtubes containing polyclonal anti-k antibodies from human serum within the gel matrix. Preservative: < 0.1% NaN_3.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow ID-Cards to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Take off the aluminium foil from as many microtubes as needed by holding the ID-Card in the upright position.
3. Add 10 µl or 12.5 µl of the red cell suspension to the appropriate microtubes.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the reactions.

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B) Reactions for K and k
Positive reactions of + to ++++ indicate presence of the corresponding antigen.
Negative reactions indicate absence of the corresponding antigen.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the human antibodies present in the ID-Cards ″Anti-K″ and ″Anti-k″ have been evaluated according to the Common Technical Specification (CTS) [5] for *in vitro* diagnostic medical devices. Evaluation was performed with samples coming from donors, patients and newborns for ″Anti-K″ and donors for ″Anti-k″ for which KEL phenotype have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-K (Anti-KEL1)	N=1460	100%	100%
Anti-k (Anti-KEL2)	N=332	100%	N/A*

* According to Mollison [1], the KK phenotype frequency in the Caucasian population is 0.2%.

REMARKS

1. Prior to testing for the presence of antigens, it must be assured that the red cells are free of enzyme reactive autoantibodies.
2. Repeated centrifugation of the ID-Cards does not affect the performance of the product.
3. Kp^a (K3) antigen in cis weakens the expression of K-antigens.

LIMITATIONS

a) Some cells which are k (K2) and Kp^a (K3) positive may react weakly to negative with anti-k test sera.
b) Cells that have become polyagglutinable, due to cryptoantigen exposure e.g. T antigen, either in vivo or *in vitro*, may react with all human antisera. Further investigation of such reactions is required.
c) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
f) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
g) Use of suspension solutions for red cells other than ID-Diluent 1 may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

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DiaClon Anti-Jk^{a} and DiaClon Anti-Jk^{b}

English B007321 v.01 2021-02-15

monoclonal
Determination of the Jk^{a} (JK1) and Jk^{b} (JK2) Antigens

Product-Identification:

ID-Card ″DiaClon Anti-Jk^{a}″	Id-n°: **50270**	1 x 12. REF 007321
ID-Card ″DiaClon Anti-Jk^{b}″	Id-n°: **50280**	1 x 12. REF 007331

INTRODUCTION

The Kidd system (ISBT number 009) consists of 3 antigens Jk^{a} (JK1), Jk^{b} (JK2) and Jk3 (JK3).

According to Mollison [1], the phenotype frequency of the antithetical Kidd antigens Jk^{a} and Jk^{b} in the Caucasian population is as follows:


Phenotype	**Frequency**
Jk (a+b-)	26%
Jk (a+b+)	50%
Jk (a-b+)	24%
Jk (a-b-)	exceedingly rare

The antibodies of the Kidd system are often associated with severe hemolytic transfusion reactions (HTRs), notably delayed reactions. The transfusion of compatible blood in the Kidd system, where antibodies have previously been detected (or are suspected), is therefore important.

The ID-Cards ″DiaClon Anti-Jk^{a}″ and ″DiaClon Anti-Jk^{b}″ may be used to determine the Kidd phenotype of donor blood prior to crossmatching, for confirmation of Kidd phenotype after antibody identification, or for paternal/fetal antigen typing in alloimmunised pregnancies.

REAGENTS IVD

- ID-Card ″DiaClon Anti-Jk^{a}″ contains 6 microtubes with monoclonal anti-Jk^{a} (cell line MS-15) within the gel matrix. Preservative : < 0.1% NaN_3.
- ID-Card ″DiaClon Anti-Jk^{b}″ contains 6 microtubes with monoclonal anti-Jk^{b} (cell line MS-8) within the gel matrix. Preservative: < 0.1% NaN_3.

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1. Dispense 0.5 ml of ID-Diluent 1 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed red cells, mix gently.
3. Incubate the red cell suspension for 10 minutes at room temperature.

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1. Identify the appropriate microtube of the ID-Card ˝DiaClon Anti-Jk[a]˝ and/or ˝DiaClon Anti-Jk[b]˝ with the patient's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Dispense 10 or 12.5 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for Jk[a] and Jk[b]
- Positive reactions of + to ++++ indicate presence of the corresponding antigen. Strong ++++ reactions are very rare.
- Anti-Jk[a] may react more strongly with Jk (a+b-) than with Jk (a+b+) red blood cells (dosage effect).
- Negative reactions indicate absence of the corresponding antigen.

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DiaClon Anti-Le[a] and DiaClon Anti-Le[b]

English

B007221 v.01 2021-02-15

monoclonal
Determination of the Lewis antigens

Product-Identification:

ID-Card ″DiaClon Anti-Le[a]″	Id-n°: **50242**	1 x 12. REF 007221
ID-Card ″DiaClon Anti-Le[b]″	Id-n°: **50250**	1 x 12. REF 007231

INTRODUCTION

In the Lewis system (ISBT number 007) the antigens are adsorbed onto the red cell surface from the plasma, rather than being an inherent part of the membrane.

The Lewis gene (*LE* or *FUT3*) interacts with the secretor gene (*SE* or *FUT2*). When the *LE* gene is present and the *SE* gene is absent, the resulting phenotype is Le (a+b-). Where both genes are present, the Le(a) (LE1) antigen is converted to Le(b) (LE2) antigen, giving rise to the Le (a-b+) phenotype. Individuals lacking the LE gene group as *Le* (a-b-).

There is also a complex interaction with the ABO system, due to the commonality of precursor substance.
Lewis antigens are absent or only poorly expressed on red cells from newborns.

According to Mollison [1], the Lewis phenotype frequency in the Caucasian population is as follows:

Phenotype	**Frequency**
Le (a-b+)	75%
Le (a+b-)	20%
Le (a-b-)	5%
Le (a+b+)	very rare *


* more common in East Asia, South East Asia, the Pacific region and Australasia.

The ID-Cards ″DiaClon Anti-Le[a]″ and ″DiaClon Anti-Le[b]″ allow testing of 6 samples simultaneously.

REAGENTS IVD

- ID-Card ″DiaClon Anti-Le[a]″ with 6 microtubes containing monoclonal anti-Le(a) (cell lines P3N20V3 + LM112/161) within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ″DiaClon Anti-Le[b]″ with 6 microtubes containing monoclonal anti-Le(b) (cell line LM129/181 [LB1]) within the gel matrix. Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: modified bromelin solution for red cell suspension.

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1. Identify the appropriate microtube of the ID-Card ″DiaClon Anti-Le[a]″ and/or ″DiaClon Anti-Le[b]″ with the patient's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for the determination of the antigens Le[a] and Le[b]
- A positive reactions of + to ++++ indicate the presence of the corresponding antigen.
- A negative reaction indicates the absence of the corresponding antigen.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performance evaluation studies have been realized with ″DiaClon Anti-Le[a]″ and ″DiaClon Anti-Le[b]″. Results obtained showed that performance characteristics of the reagent are in conformity with the application field and comparable to other approved products.

Reproducibility
Intraassay reproducibility (repeatability) and inter-assay reproducibility of the ″DiaClon Anti-Le[a]″ and ″DiaClon Anti-Le[b]″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reaction strenghts in positive samples were less than one reaction strength.

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REMARKS

1. Prior to testing for the presence of antigens, it must be assured that the red cells are free of enzyme reactive autoantibodies.
2. Repeated centrifugation of the ID-Cards does not affect the performance of the product.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 1 may modify the reactions.
f) Too heavy or too weak red cell suspensions can cause aberrant results.
g) Lewis antigens may be lost as result of disease. Lewis antigens are often greatly reduced during pregnancy.

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INTRODUCTION

The antithetical antigens Kp[a] (KEL3), Kp[b] (KEL4) and Kp[c] (KEL21) are part of the Kell system (ISBT, number 006). Kp[a] is present in approximately 2% of the Caucasian population while Kp[b] is a high frequency (public) antigen. The Kp[c] antigen has an occurrence of less than 0.01%, but up to 0.32% in Japanese.

According to Mollison [1], the phenotype frequency in the Caucasian population of the two most important antigens is as follows:

Phenotype	**Frequency**
Kp a+b-	rare
Kp a+b+	± 2.3%
Kp a-b+	± 97.7%

Anti-Kp[a] and anti-Kp[b] are less common than anti-K but have similar serological characteristics and are considered to be clinically significant. They may occur after transfusion or by feto-maternal immunization.

The ID-Cards ″Anti-Kp[a]″ and ″Anti-Kp[b]″ allow testing of 6 samples simultaneously.

REAGENTS IVD

- ID-Card ″Anti-Kp[a]″ with 6 microtubes containing polyclonal anti-Kp[a] antibodies from human serum within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ″Anti-Kp[b]″ with 6 microtubes containing polyclonal anti-Kp[b] antibodies from human serum within the gel matrix. Preservative: < 0.1% NaN_3.

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1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to the appropriate microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for Kp[a] and Kp[b]
- Positive reactions of + to ++++ indicate presence of the corresponding antigen.
- Positive reactions may appear as a double population. However, a double population may also indicate the presence of cells positive and negative for the corresponding antigen (e.g. Post-transfusion, if the antigen configuration of the transfused blood was different from the patient's).
- Negative reactions indicate absence of the corresponding antigen.

REMARKS

1. Prior to testing for the presence of antigens, it must be assured that the red cells are free of enzyme reactive autoantibodies.
2. Repeated centrifugation of the ID-Card does not affect the performance of the product.

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LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either *in vivo* or *in vitro*, may react with all human antisera. Further investigation of such reactions is required.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions for red cells other than ID-Diluent 1 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant reactions.

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Anti-Lu^a and Anti-Lu^b

English B007351 v.01 2020-09-25

human
Determination of the Lu^a and Lu^b Antigens

Product-Identification:

ID-Card ″Anti-Lu^a″	Id-n°: **50310**	1 x 12. REF 007351
ID-Card ″Anti-Lu^b″	Id-n°: **50320**	1 x 12. REF 007361

INTRODUCTION

According to Mollison [1], the phenotype frequency of the Lutheran system in the Caucasian population is as follows:

Phenotype	Frequency
Lu (a+b-)	0.15%
Lu (a+b+)	7.50%
Lu (a-b+)	92.35%
Lu (a-b-)	very rare

Anti-Lu^b has been reported to be implicated in mild to moderate transfusion reactions and as the cause of mild HDN. The antigen-antibody reaction produces small and loose agglutinates, interspersed among unagglutinated red cells, causing characteristic mixed-field appearance [3].


Antibodies against antigens in the Lutheran system (ISBT terminology symbol LU, number 005) are not often encountered. Anti-Lu^a and anti-Lu^b are produced in response to pregnancy or transfusion. The antigens are poorly developed at birth, consequently, anti-Lu^a has not been found to be the cause of hemolytic disease of the newborn (HDN) and is not considered to be of clinical significance.

The ID-Cards ″Anti-Lu^a″ and ″Anti-Lu^b″ allow testing of 6 samples simultaneously.

REAGENTS IVD

- ID-Card ″Anti-Lu^a″ with 6 microtubes containing polyclonal anti-Lu^a antibodies from human serum, within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ″Anti-Lu^b″ with 6 microtubes containing polyclonal anti-Lu^b antibodies from human serum, within the gel matrix. Preservative: < 0.1% NaN_3.

Caution: The source materials from which these products were manufactured, were found non-reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

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1. Dispense 0.5 ml of ID-Diluent 1 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.
3. Incubate the red cell suspension 10 minutes at room temperature (18–25 °C).

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Allow ID-Cards to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to the appropriate microtubes.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the reactions.

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B) Reactions for Lu[a] and Lu[b]

- Positive reactions of + to ++++ indicate presence of the corresponding antigen.
- Positive reactions may appear as a double population. However, a double population may also indicate the presence of cells positive and negative for the corresponding antigen (e.g. Post-transfusion, if the antigen configuration of the transfused blood was different from the patient's).
- Strong ++++ reactions are very rare.
- Negative reactions indicate absence of the corresponding antigen.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the polyclonal antibodies present in the ID-Cards ″Anti-Lu[a]″ and ″Anti-Lu[b]″ are verified by tests with positive and negative samples to ensure conformity with product specifications. The results obtained with the ID-Cards ″Anti-Lu[a]″ and ″Anti-Lu[b]″ were compared directly with an established device already CE marked:

ID-Card	Total number of samples	Sensitivity	Specificity
Anti-Lu[a], polyclonal	342	100%	100%
Anti-Lu[b], polyclonal	342	100%	100%

Results obtained showed that performance characteristics of each test of the ID-Cards ″Anti-Lu[a]″ and ″Anti-Lu[b]″ are in conformity with the application field and comparable to other approved products.

REMARKS

Repeated centrifugation of the ID-Card does not affect the performance of the product.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Cells that have become polyagglutinable, due to cryptoantigen exposure, e.g. T antigen, either in vivo or *in vitro*, may react with all human antisera. Further investigation of such reactions is required.
c) Bacterial or other contamination of materials used can cause false positive or false negative results.
d) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions for red cells other than ID-Diluent 1 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant reactions.

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DiaClon Anti-M and DiaClon Anti-N

English B007011 v.01 2021-02-15

monoclonal
Determination of the M (MNS1) and N (MNS2) Blood groups

Product-Identification:

ID-Card ″DiaClon Anti-M″	Id-n°: **50212**	1 x 12 REF 007011
ID-Card ″DiaClon Anti-N″	Id-n°: **50221**	1 x 12 REF 007111

INTRODUCTION

M (MNS1) and N (MNS2) antigenic determinants are carried on glycophorin A (GPA), one of the major sialic acid containing structure of the red cell membrane.

According to Mollison [1], the phenotype frequencies of the MN system in the Caucasian population are as follows:

MM..... 28%
MN..... 50%
NN 22%

Interpretation of results with anti-M and anti-N sera require particular care. Discrepancies in the reactivity may be observed between results obtained with test sera of human origin, from rabbits, lectins and monoclonal reagents.

Glycophorin B (GPB) carries the S (MNS3), s (MNS4) antigenic determinants and also expresses the ″N″ antigen. Some anti-N sera may therefore cross react with red cells carrying double dose of the M antigen. Unexpected results may also occur when variants of the M and N antigens are present, for example, the rare Mg (MNS11) antigen which is encoded by a gene which produces virtually no M or N antigen. Red cells can therefore be erroneously classified NN or MM when they are genotypically NM_g or MM_g. M^c (MNS8) is another rare allele of M and N which reacts with the majority of anti-M sera and some anti-N sera. Me (MNS13) is a determinant shared by M and He antigens and some anti-M sera also contain an anti-Me component. These potential sources of discrepant results of MN typing can be of particular importance in parentage testing.

M and N typing is also undertaken in order to help to confirm allo-antibody specificity. Anti-M and anti-N are not normally considered to be of clinical significance unless reactive at 37 °C. In this case, MN typing of donor blood is important so that appropriate antigen negative units may be selected before crossmatching.

The ID-Cards ″DiaClon Anti-M″ and ″DiaClon Anti-N″ allow testing of 6 samples simultaneously.

REAGENTS IVD

- ID-Card ″DiaClon Anti-M″ contains 6 microtubes with monoclonal anti-M antibodies (cell line 9211 9D5 1G10) within the gel matrix (containing bovine albumin). Preservative: < 0.1% NaN_3.
- ID-Card ″DiaClon Anti-N″ contains 6 microtubes with monoclonal anti-N antibodies (cell lines 1422 C 7 + MN879) within the gel matrix (containing bovine albumine). Preservative: < 0.1% NaN_3.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

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B) Reactions for M and N antigens
- Positive reactions of + to ++++ indicate presence of the corresponding antigen. Strong (++++) reactions are very rare.
- Positive reactions may appear as a double population. However, a double population may also indicate the presence of cells positive and negative for the corresponding antigen (e.g. Post-transfusion, if the antigen configuration of the transfused blood was different from the patient's).
- Negative reactions indicate absence of the corresponding antigen.

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LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential.
e) The equipment should be checked regularly according to GLP procedures.
f) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
g) Too heavy or too weak red cell suspensions can cause aberrant results.

REMARKS

Repeated centrifugation of the ID-Cards does not affect the performance of the product.

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English

B007531 v.01 2021-02-15

Product-Identification: 09410
˝ID-DiaClon Anti-D˝

1 x 5 ml REF 007531

INTENDED USE

The ID-DiaClon Anti-D is intended to be used for the confirmation of RhD negative results and complementary characterization of D variants (weak D and partial D).

INTRODUCTION

The increasing use of monoclonal anti-D's from different cell lines can create problems in the interpretation of D antigen testing, in particular the differentiation of weak D and partial D (D variant).

In the ID-Card System, a +++ to ++++ reaction indicates the presence of the D antigen and weaker reactions indicate weak D. These weak positive reactions can vary depending on the monoclonal cell line used.


Cells giving negative reactions in the direct test and positive by IAT are considered as weak D or partial D (D variant).

A supplementary verification can be performed in the ID-System with the aid of a monoclonal anti-D formulated to characterize weak D's and DVI using the ID-Card ˝Coombs Anti-IgG˝ or the ID-Card ˝LISS/Coombs˝.

REAGENTS IVD

˝ID-DiaClon Anti-D˝ containing monoclonal anti-D IgG (cell line ESD1), supplied as a ready-to-use reagent in vials of 5 ml. Preservatives: trimethoprime and sulfamethoxazole.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Card ˝Coombs Anti-IgG˝ with 6 microtubes containing anti-IgG (rabbit), within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Card ˝LISS/Coombs˝ with 6 microtubes containing polyspecific AHG (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix. Preservative: < 0.1% NaN_3.
- ID-Diluent 2: modified LISS for red cell suspension.

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CONTROLS

Known weak D positive and negative samples should be included to validate the results.
As false positive results may occur in cases where the direct anti-globulin test (DAT) is positive, all positives should be confirmed negative in the DAT using the ID-Card ″Coombs Anti-IgG″ or the ID-Card ″LISS/Coombs″.

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1. Identify ID-Card with donor's/neonate's number/details as appropriate (see REMARKS).
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Add 50 µl of the red cell suspension to the appropriate microtube.
4. Add 50 µl of the ″ID-DiaClon Anti-D″ for confirmation of weak D by IAT to the microtube.
5. Incubate the ID-Card 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the result.

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* In order to distinguish between variants (weak D and partial-D), it is advised to further investigate the sample using the ID-Partial RhD Typing Set (REF 001451) or Extended Partial RhD Typing Set (REF 001461) (see related package inserts).

Differentiation between RhD and Rh weak D:	**Reactions**	**Reactions**
Direct test using an ID-Card with anti-D	+++ to ++++	++ to neg.
Test with ″ID-DiaClon Anti-D″ and the ID-Card ″Coombs Anti-IgG″ or the ID-Card ″LISS/Coombs″.	++++	+++ to +
Results	RhD positive	Rh weak D/DVI

REMARKS

The anti-D contained in this vial has been selected to react with the DVI phenotype.

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PERFORMANCE CHARACTERISTICS

Performances of ID-DiaClon Anti-D have been evaluated during internal studies against a panel of variant samples. This reagent was able to detect weak D's Type 1, 2, 3, 4.0, 4.1, 5, 6, 7, 8, 9 ,10 11, 15, 17.
The cell line included in the ID-DiaClon anti-D is also known to detect the following variants; DII/DNU, DIII, DV (only DVa evaluated), DCS, DVI, DVII, DOL, DFR, DAR and DAR-E, DKH/DAU-4 but not the following ones : DIV, DBT and DHAR (3).

BIBLIOGRAPHY

1. Mollison P. L., Engelfriet C. P. and Contreras M.: Blood Transfusion in Clinical Medicine. 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109–113.
3. Robb JS an Allan JC. (2005). Evaluation of a New Kit for the Identification of partial RhD Types by Heamagglutination. Transfusion Medicine, 15, Supplement 1: p51.

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1. Mollison P. L., Engelfriet C. P. and Contreras M.: Blood Transfusion in Clinical Medicine. 10th ed. 1997; Blackwell Scientific Publications, Oxford.
2. Lapierre Y., Rigal D., Adam J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Robb JS an Allan JC. ((2005). Evaluation of a New Kit for the Identification of partial RhD Types by Haemagglutination. Transfusion Medicine, 15, Supplement 1: p51.

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Set REF 001451
1 x 12 ID-Cards
6 x 0.4 ml ID-Anti-D 1–6

INTRODUCTION

The D antigen of the Rh system, first described in 1939, is now considered a mosaic of epitopes. Individuals whose cells lack part of the D mosaic may, when exposed to a complete D antigen, make antibody to the epitopes which are lacking on the individual's cells.


Investigations have shown that there are now more than 30 recognisable epitopes but for practical purposes the 9-epitope model of the D antigen is satisfactory. Selected panels of monoclonal anti-D reagents can be used to distinguish the red cell categories of partial-D phenotypes, often in conjunction with antibodies to associated low frequency antigens for confirmation.

The ˝ID-Partial RhD Typing Set˝ consists of 6 monoclonal anti-D which can aid further characterisation of the RhD antigen after routine RhD typing. The sera can differentiate between categories II, IV, V, VI, VII DFR, DBT and R_O^{Har}. Category I is obsolete and category III red cells react with all anti-D reagents but can be distinguished from normal RhD positive cells by the presence of allo-anti-D in the serum.

Partial DVI remains the most important category to define. Clinically, patients with partial DVI red cells should be treated in the same way as RhD negative individuals both for transfusion purposes and during/after pregnancy with an RhD positive fetus, to prevent the production of anti-D to the missing epitopes with possible subsequent clinical sequelae.

REAGENTS IVD


- ID-Card ˝ID-Partial RhD Typing Set˝ with 6 microtubes containing polyspecific antihuman globulin [rabbit anti-IgG and monoclonal anti-C3d (cell line C139-9)], within the gel matrix.
 Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

- 6 monoclonal anti-D's (cell lines see interpretation table B001452), with addition of bovine albumin, labelled 1–6, in vials of 0.4 ml.
 Preservative: < 0.1% NaN_3.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow the anti-D reagents and samples to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2 Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to each microtube.
4. Add 25 µl of each anti-D (1–6) to the appropriate microtube.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for partial RhD characterisation
See ″REMARKS″ and also the interpretation table B001452 and ″REMARKS″.

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REMARKS

1. Red cells must be free from *in vivo* or *in vitro* sensitisation with IgG: a negative DAT validates the reactions obtained in the ID-Partial RhD Typing Test.
2. Weak D phenotypes posessing lower numbers of antigen sites per cell may not give a positive reaction with the IgM anti-D number 6.
3. Some examples of DVI cells may fail to react with antibody no.1. This may be due to the DVI type I phenotype. This partial D type is associated with a cDE haplotype and is thought to carry fewer accessible RhD proteins per cell than DVI types II and III [3].
4. Some examples of R_O^{Har} may react with anti-D number 6.
5. There are many other known partial-D types. Therefore, if the red cells cannot be clearly characterised as one of the partial-D types shown on the table, then the sample should be submitted for further testing to a specialised laboratory.
6. Some examples of weak D red cells may react more weakly with antibodies numbers 3 and 6.

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Fibrin residues in the serum or the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
f) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
g) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Dr Daniels G. Human Blood Groups. 1995. Blackwell Scientific Publications. Oxford.
2. Lomas C, McColl K, Tippett P. Further complexities of the Rh antigen D disclosed by testing category DII cells with monoclonal anti-D. Transfusion Medicine. 1993. 3: 67-69.
3. Wagner FF; Gassner C; Muller TH, Schönitzer D; Schunter F; Flegel WA. Three molecular structures cause rhesus D category VI phenotypes with distinct immunohaematologic features. Blood, 91 (6): 2157-68.1998.

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1. Dr Daniels G. Human Blood Groups. 1995. Blackwell Scientific Publications. Oxford.
2. Lomas C, McColl K, Tippett P. Further complexities of the Rh antigen D disclosed by testing category DII cells with monoclonal anti-D. Transfusion Medicine. 1993. 3: 67-69.
3. Wagner FF; Gassner C; Muller TH, Schonitzer D; Schunter F; Flegel WA. Three molecular structures cause rhesus D category VI phenotypes with distinct immunohaematologic features. Blood, 91 (6): 2157-68.1998.

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English B002418 v.01 2020-10-07

A, B, DVI-/ 3 x AHG

Confirmation of blood group A, B and D combined with antibody screening

Product-Identification: 50682
ID-Card ˝DiaClon Type + Screen˝

4 x 12	REF 002437
24 x 12	REF 002431
60 x 12	REF 002439
112 x 12	REF 002438

INTRODUCTION

Under certain rigorous conditions, compatibility tests before a blood transfusion may be replaced by ˝Type and Screen˝ (T/S) of the recipient [1].

However, before elimination of the compatibility test, it should be taken into consideration that modern procedures such as the compatibility test with the gel test (ID-System) are as fast, easy and economical as additional typing and screening of the patient's blood. The compatibility test still offers the highest safety for any individual patient needing transfusion.

T/S might be considered for cases of low transfusion probability whereas in the presence of irregular antibodies and history of transfusion complications a T/S procedure must be firmly excluded in favour of a complete compatibility test [2].

The ID-Card ˝DiaClon Type + Screen˝ is designed for the confirmation of the recipient's ABO blood group and RhD type, and verification of the lack of irregular antibodies with the indirect antiglobulin test (IAT, LISS/Coombs) with 3 test cell reagents (ID-DiaCell I-II-III).


The ˝Type and Screen˝ test must not be used to replace full crossmatch where the crossmatch is designated as a mandatory procedure by the national regulatory body.

REAGENTS IVD

ID-Card ˝DiaClon Type + Screen˝, with 3 microtubes containing monoclonal anti-A [cell line A5], monoclonal anti-B [cell line G½] and monoclonal anti-D [cell lines LHM59/20 (LDM3) + 175-2], and 3 microtubes containing polyspecific anti-human globulin (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within the gel matrix. Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ˝DiaClon Type + Screen˝ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

Caution: All reagents should be treated as potentially infectious.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-DiaCell I-II-III for antibody screening
- ID-Diluent 2: Modified LISS for red cell suspensions.

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When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use to avoid fibrin residues, which may interfere with the reaction pattern.

PREPARATION OF BLOOD SAMPLE

a) Red cell suspension (for ABO/D typing)
Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed red cells, mix gently.

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b) Plasma or serum for antibody screening
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

CONTROLS

Known positive and negative samples should be included in accordance to the relevant guidelines of quality assurance.

TEST PROCEDURE

ABD-Confirmation and antibody screening
Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

Allow the test cell reagents and samples to reach room temperature before use and resuspend gently.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 10 or 12.5 µl of the patient's red cell suspension into the first 3 microtubes of the ID-Card (anti-A, anti-B and anti-D).
4. Pipette 50 µl of each ″ID-DiaCell I-II-III″ to the microtubes 4, 5 and 6 (3 x AHG).
5. Add 25 µl of the patient's plasma or serum to the microtubes 4, 5 and 6.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

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INTERPRETATION OF THE RESULTS

A) Principle [3]

Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for blood groups ABO

Anti-A	Anti-B	Blood group
+++ to ++++	negative	**A**
negative	+++ to ++++	**B**
+++ to ++++	+++ to ++++	**AB**
negative	negative	**O**

Weaker reactions than +++ may indicate A or B subgroups and further investigations should be performed. For correct interpretation, a complete grouping test should be performed. In the presence of weak or very weakly expressed antigens the reaction can be negative.

C) Reactions for RhD

RhD positive	RhD weak	RhD negative
+++ to ++++	± to ++*	negative

* ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial D types as appropriate for the category of sample being tested.

The anti-D test serum in the ID-Card ˝DiaClon Type + Screen˝ was selected so as not to react with DVI-variants.Weak D may give a negative reaction. If all weak D's are required to be detected, all D negative results must be retested. If used for testing of patients samples, note that most guidelines do not recommend further testing for weak or partial D phenotypes.

D) Reactions for antibody screening

- A negative reaction indicates the absence of detectable irregular antibodies in the patient's serum or plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions observed on the antigen table. Verify that the lot number of the ˝ID-DiaCell I-II-III˝ test cells correspond to the lot number indicated on the antigen table.
- Where there is a positive reaction with all ˝ID-DiaCell˝ test cells it is advisable to perform an autocontrol. Proceed as follows:
 1. Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
 Allow the diluent to reach room temperature before use.
 1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
 2. Add 10 µl of packed red cells or 20 µl of whole blood, mix gently.
 The cell suspension may be used immediately.
 2. Identify the appropriate microtube of the ID-Card ˝LISS/Coombs˝ with patient's name or number.
 3. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
 4. Pipette 50 µl of the patient's red cell suspension to the appropriate microtube.
 5. Add 25 µl of the patient's plasma or serum.
 6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
 7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
 8. Read and record the results.
- If the autocontrol is negative, the positive reactions with the test cells may be considered valid and the usual further tests should be undertaken to identify the antibody.
- If the autocontrol shows an identical positive reaction as the 3 test cells, it suggests that these positive reactions are probably due to a non-specific reaction.
- If one or several of the test cells show a stronger positive reaction than the autocontrol, it should be assumed that, as well as the non-specific reaction, an antibody may be present and the patient sample should be submitted to further testing.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon Type + Screen˝, as well as the antibody screening, have been evaluated according to the Common Technical Specification (CTS) [4] for *in vitro* diagnostic medical devices.

ABO/RH1 grouping

Evaluation was performed with samples coming from donors, patients and newborns for which ABO/RH1 phenotypes have been previously determined by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: A5	4385	100%	100%
Anti-B, Cell line: G1/2	4385	100%	100%
Anti-D(DVI-), Cell lines: LHM59/20(LDM3),175-2	4679*	99.72%**	100%

* These samples include 81 samples with weak (Weak/variant) expression of the RH1 antigen.
** 62 weak/variant samples have been detected directly. The DVI variants were not detected.

Remark: The Anti-DVI- is not conceived to detect every weak/variant expressions of the RH1 antigen.

Antibody screening

Evaluation was performed with samples coming from donors and patients already characterized by a reference method:

Antibodies	Total number of samples	Sensitivity	Specificity
Anti-IgG + C3d	201	97.08%	100%

In addition, the test procedure specific for the ˝DiaClon Type + Screen˝ has been evaluated with a minimum of 500 samples coming from donors and patients. The evaluation resulted in 100% of sensitivity and specificity on the concerned wells Anti-A, Anti-B and Anti-DVI-.

REMARKS

1. The ABD-confirmation does not replace the complete ABO/Rh-determination. It should only be used as a confirmation of previously determined blood groups where, according to general rules, ABO-typing is confirmed with reverse grouping and ABO/Rh-Typing is validated by the presence of a negative control test.
2. The anti-B of monoclonal origin does not react with the acquired B antigen.

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
f) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
g) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

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BIBLIOGRAPHY

1. Fopp M., et al, Safety of ″Type and Screen″ programs, First Int. Symposium, 17-18 March 1993.
2. Holland P.V. ed.: Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington, V.A: American Association of Bloodbanks, 1984: 15-39.
3. Lapierre, Y., Rigal, D., Adam, J. et al,: The gel test: A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

GLOSSARY OF SYMBOLS

The following symbols **may** be used for labelling purpose.

REF	Catalog reference
LOT	Batch number
IVD	*In vitro* diagnostic
	Consult instructions for use
	Expiry date (YYYY-MM-DD)
	Storage temperature
	Legal manufacturer
⚠	Consult **downloads.bio-rad.com** to download the latest version of these instructions for use

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

These products may only be used with the following analyzers authorized by Bio-Rad: IH-1000, IH-500, Saxo ID-Reader and Swing TwinSampler.*

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1. Fopp M., et al, Safety of "Type and Screen" programs, First Int. Symposium, 17-18 March 1993.
2. Holland P.V. ed.: Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington, V.A: American Association of Bloodbanks, 1984: 15-39.
3. Lapierre, Y., Rigal, D., Adam, J. et al,: The gel test: A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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INTRODUCTION

Test cell reagents are routinely used in blood group serology to detect the presence or absence of anti-A and anti-B isoagglutinins (in reverse grouping).

For reverse grouping, test cell reagents of group A_1 and B / A_1, A_2 and B / A_1, B and O or A_1, A_2, B and O are used, following various requirements and guidelines.

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
USE OF THE ID-TEST CELL REAGENTS

- All Test cell reagents are for use with the ID-Cards of the ID-System only.
- **Strictly follow the test procedures as described in the specific package inserts of the ID-Cards to be used.**
- Always gently resuspend the red cells, by inverting the vial several times before use and also before placing the vials into a pipetting automat.
- Make sure that the test cells are at room temperature (18–25 °C) when in use.
- During the working procedures, check that the test cell reagents remain in suspension. If there is settling of the cells, resuspend again.
- **For the ID-System, precise pipetting is of importance.** Use the ID-Pipetors for serial pipetting.
- Avoid contamination of the test cell reagents.
- After use, close the vials and replace them in the refrigerator.

INTERPRETATION OF THE RESULTS

A) Principle
Agglutination indicate the presence of the agglutinins.
No agglutination indicate the absence of agglutinins.

B) Reaction with ID-DiaCell ABO

A_1	A_2	B	O	Isoagglutinins	Blood groups
-	-	+	-	Anti-B	**A**
+	+	-	-	Anti-A	**B**
-	-	-	-	None	**AB**
+	+	+	-	Anti-A and Anti-B	**O**

Note: The results of the reverse grouping must agree with the blood group observed from forward (antigen) typing. Discrepant results should be subjected to further investigation, preferably with a new sample.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performances of the test cells ″ID-DiaCell ABO″ have been evaluated according to the Common Technical Specification (CTS) [3] for *in vitro* diagnostics medical devices. Evaluation was performed with samples coming from donors and patients for which ABO phenotype has been previously determined by a reference method:

ID-DiaCell	Total number of samples	Sensitivity	Specificity
A_1	N=5262	100%	100%
A_2	N=2788	97.62%*	100%
B	N=5262	100%	100%
O	N=2788	N/A	100%

* Sensitivity of 98.59% with 10 minutes of incubation at room temperature and 93.21% without incubation.

REMARKS

1. The optimal reaction temperature for isoagglutinins is 4 °C. If weak or doubtful reactions are observed, repeat the test with an incubation at 2–8 °C for 15 minutes.
2. Atypical reactions require further studies.
3. In general, anti-A and anti-B iso-antibodies are not detectable in newborns and infants. Older people and patients with agammaglobulinaemia may also lack these iso-antibodies.

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BIBLIOGRAPHY

1. Technical Manual of the American Association of Blood Blanks, 13th edition, 1999.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990; 30:109-113.
3. Official Journal of the European Union L 318/25 : Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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NaCl, Enzyme Test and Cold Agglutinins

English

B005014 v.01 2020-12-18

Detection of saline and enzyme reactive antibodies (for reverse grouping, antibody screening and identification, compatibility tests)

Product-Identification: 50520

ID-Card ″NaCl, Enzyme Test and Cold Agglutinins″

4 x 12. REF 005014
24 x 12. REF 005017
60 x 12. REF 005016
112 x 12. REF 005015

INTRODUCTION

Saline techniques are used to detect antibodies that react predominantly at 4 °C or at 18 – 25 °C (room temperature), such as anti-M, -N, -P_1, -Le^a, -Le^b, -I. They are a lso used to define hemolytic anemia associated with cold antibodies [1].


Enzyme techniques are useful when increased sensitivity in antibody screening is desired. They enhance the reactions of certain antibodies, notably in the Rh, Kell and Kidd systems. Since proteolytic enzymes usually destroy some antigens such as M (MNS1), N (MNS2), S (MNS3), Fy(a) (FY1), Fy(b) (FY2), an enzyme technique should never be the only method employed.

It is not considered necessary to include the room temperature saline test in routine antibody screening procedures, but the saline test at 4 °C can be used to detect cold agglutinins [2].

The ID-Card ″NaCl, Enzyme Test and Cold Agglutinins″ is suitable for antibody screening and identification procedures, compatibility testing and reverse grouping.

REAGENTS IVD

ID-Card ″NaCl, Enzyme Test and Cold Agglutinins″ with 6 microtubes containing neutral gel suspension.
Preservative: < 0.1% NaN_3.

Do not store near any heat, air conditioning sources or ventilation outlets.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 1: modified bromelin solution for red cell suspension.
- ID-Diluent 2: modified LISS for red cell suspension.
- ID-Papain: Papain solution for enzyme tests.
- ID-DiaCell I-II-III for antiglobulin test and cold agglutinins.
- ID-DiaCell I-II-III P (papainized) for enzyme technique.
- ID-DiaCell I-II for AHG test and cold agglutinins.
- ID-DiaPanel for AHG test and cold agglutinins.
- ID-DiaPanel P (papainized) for enzyme technique.
- ID-DiaCell ABO for reverse grouping.
- ID-I-negative cell.

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SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anti-coagulant) may also be used.

PREPARATION OF BLOOD SAMPLE

a) Red cells for compatibility test and autocontrol in the saline and enzyme techniques
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

The cell suspension may be used immediately.

b) Plasma or serum
Where samples are not for immediate testing they should be stored at 2 – 8 °C after separation (see also under ˝SAMPLE MATERIAL˝) for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow ID-test cell reagents, ID-Diluents, ID-Papain and plasma or serum samples to reach room temperature before use (except for cold agglutinins).

A) Antibody screening
I. Enzyme test with papainized cells (˝ID-DiaCell I-II-III P˝)
Resuspend test cells gently before use.
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ˝ID-DiaCell I-II-III P˝ test cells (papainized) to the appropriate microtubes (marked 1 – 3).
4. When an autocontrol is to be included, pipette 50 µl of the patient's red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. For the autocontrol, add 25 µl ID-Papain to the appropriate microtube.
7. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
8. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
9. Read and record the results.

II. Saline test at 4 °C (˝ID-DiaCell I-II-III˝)
Resuspend test cells gently before use.
Before use, place the ID-Cards in the refrigerator (2–8 °C) for 2 hours.
Use refrigerated (cold) test cell reagents and samples.
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ˝ID-DiaCell I-II-III˝ test cells to the appropriate microtube (marked 1–3 as appropriate).
4. When an autocontrol is to be included, pipette 50 µl of the patient's red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 30 minutes at 2–8 °C.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

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B) Antibody identification

I. Enzyme test with papainized cells (˝ID-DiaPanel P˝)

Resuspend test cells gently before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ˝ID-DiaPanel P˝ test cells (papainized) to the appropriate microtubes (marked 1–11).
4. When an autocontrol is to be included, pipette 50 µl of the patient's red cell suspension to microtube 12.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. For the autocontrol, add 25 µl ID-Papain to microtube 12.
7. Incubate the ID-Cards for 15 minutes at 37 °C in the ID-Incubator.
8. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
9. Read and record the results.

II. Saline test at 4 °C (˝ID-DiaPanel˝)

Resuspend test cells gently before use.
Before use, place the ID-Cards in the refrigerator (2–8 °C) for 2 hours.
Use refrigerated (cold) test cell reagents and samples.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ˝ID-DiaPanel˝ test cells to the appropriate microtubes (marked 1–11).
4. When an autocontrol is to be included, pipette 50 µl of the patient's red cell suspension to microtube 12.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. Incubate the ID-Cards for 30 minutes at 2–8 °C.
7. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

C) Compatibility test

I. Enzyme test (papain or bromelin)

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor's red cell suspension to the appropriate microtubes.
4. When an autocontrol is to be included, pipette 50 µl of the patient's red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Add 25 µl of ID-Papain or ID-Diluent 1 to each microtube.
7. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
8. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
9. Read and record the results.

II. Saline test

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor's red cell suspension to the appropriate microtubes.
4. When an autocontrol is to be included, add 50 µl of the patient's red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at room temperature (18–25 °C).
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

D) Reverse grouping (˝ID-DiaCell ABO˝)

Resuspend test cells gently before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ID-DiaCell ABO test cells (A_1, A_2, B and O) to the appropriate microtubes (marked with the corresponding test cell).
4. Add 50 µl of the patient's or donor's plasma or serum to each microtube.
5. An incubation for 10 minutes at room temperature is recommended (see remarks, point 2).
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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INTERPRETATION OF THE RESULTS

A) Principle [3]
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for:

I. Antibody Screening

- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with some test cells and a negative autocontrol suggest the presence of a specific antibody.
- A positive reaction with all test cells and a positive autocontrol may be due to an auto-antibody.
- Where there is a positive reaction with all test cells and a positive autocontrol but with one or several test cells showing a stronger positive reaction than the autocontrol, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

II. Antibody Identification

- One or more positive reactions indicate the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaPanel″ correspond to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present can, in most cases, be identified (autocontrol must be negative).
- A positive reaction with all ″ID-DiaPanel″ test cells and a negative autocontrol may be due to nonspecific reactions or may indicate the presence of an allo-antibody directed against a high frequency antigen.
- A positive reaction with all ″ID-DiaPanel″ test cells and a positive autocontrol are probably due to an auto-antibody.
- A positive reaction with all ″ID-DiaPanel″ test cells and a positive autocontrol but with one or more test cells showing stronger reactions than the autocontrol, may indicate an underlying allo-antibody and further investigation should be undertaken.
- It should be noted that enzyme auto-antibodies may react less strongly with the patient's cells in the one-stage technique than against ″ID-DiaCell I-II-III P″ (two-stage technique).
- As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Bio-Rad Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).
- Strong cold antibodies can be the cause of positive reactions, especially in the enzyme test procedure. Complete antibody screening procedures should be undertaken.
- Anti-I is a common specificity in cold agglutinin disease. Confirmation of specificity may be undertaken with the ″ID-I negative cell″.

III. Compatibility test

- A positive reaction indicates possible incompatibility due to the presence of irregular antibodies present in the recipient's serum or plasma, and further investigation is required.
- A negative reaction indicates that there are no detectable antibodies in the recipient's plasma/serum which are directed against antigens present on the donor's red cells, by the technique used.
- When results indicates an incompatible crossmatch, antibody screening and identification tests should be performed and ABO and Rh blood groups of both recipient and donor should be retested.

IV. Reverse grouping

Reaction with

Test cell reagents				Blood group
A_1	**A_2**	**B**	**O**	
negative	negative	+ to ++++	negative	**A**
+ to ++++	± to ++++	negative	negative	**B**
negative	negative	negative	negative	**AB**
+ to ++++	± to ++++	+ to ++++	negative	**O**

- An autocontrol should be included when there is a discrepancy between the ABO forward typing result with test sera and the ABO result of the reverse grouping and/or if unexpected positive reactions are obtained with the screening cells.
- When there are discrepant results between ABO typing and reverse grouping, consult the Bio-Rad blood group ABO identification chart for appropriate procedures.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performance evaluation studies have been realized with ID-Cards ˝NaCl, Enzyme Test and Cold Agglutinins˝. Results obtained showed that performance characteristics of each application of the card are in conformity with the application field and comparable to other approved products.

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ˝NaCl, Enzyme Test and Cold Agglutinins˝ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

REMARKS

1. If hemolysis is observed in the upper part of the microtube in the absence of a haemolysed sample it should be interpreted as positive.
2. For reverse grouping, an incubation of at least 10 minutes at 18–25 °C prior to centrifugation will enhance the reactions and minimise repeat testing due to weak isoagglutinins.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Bacterial or other contamination of materials used can cause false positive or false negative results.
c) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
d) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
e) Use of suspension solutions other than ID-Diluent 1 and 2 may modify the reactions.
f) Use of test cells other than the ID-Test cells may modify the reaction patterns.
g) Too heavy or too weak red cell suspensions can cause aberrant results.

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LISS/Coombs

English B004014 v.01 2020-11-05

Indirect and direct antiglobulin test

Product-Identification: 50531
ID-Card ″LISS/Coombs″

4 x 12	REF 004014
24 x 12	REF 004017
60 x 12	REF 004016
112 x 12	REF 004015

INTRODUCTION

Polyspecific anti-human globulin (AHG) reagents are used for routine alloantibody detection and identification, compatibility tests and the direct antiglobulin test (DAT).

The most important function of the polyspecific AHG reagent is to detect the presence of IgG. The importance of anticomplement in the AHG reagent is debatable since antibodies detectable only by their ability to bind complement are rather rare. However, anti-C3d activity is important for the DAT in the investigation of autoimmune hemolytic anemia (AIHA) [4].
A positive DAT generally indicates that the red cells are coated *in vivo* with immunoglobulin and/or complement.


The microtubes of the ID-Card ″LISS/Coombs″ contain polyspecific AHG, to be used for antibody screening, antibody identification, crossmatch and the DAT. For the indirect antiglobulin test (IAT), labour intensive washing procedures are eliminated, due to the fact that the red cell suspension is added to the microtube before the plasma/serum, creating a barrier over the gel suspension, thus avoiding neutralization of the AHG by serum IgG proteins.

The anti-human globulin IgG used in the ID-Card ″LISS/Coombs″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

The ID-Card ″LISS/Coombs″ is suitable for the DAT, for the compatibility test, for antibody screening and identification with ″ID-DiaCell″ and ″ID-DiaPanel″.

REAGENTS IVD

ID-Card ″LISS/Coombs″ with 6 microtubes containing polyspecific AHG (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
- ID-DiaCell, ID-DiaPanel: Test cell reagents.

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SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

PREPARATION OF BLOOD SAMPLE

a) Red cell suspension for DAT or autocontrol
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

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b) Red cell suspension for crossmatch procedures
Prepare a 0.8% red cell suspension in ID-Diluent 2 as above. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

c) Plasma or serum for indirect antiglobulin test (IAT) procedures
Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under ″SAMPLE MATERIAL″) for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURES

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

I. Direct antiglobulin test (DAT)
1. Identify the appropriate microtubes of the ID-Card ″LISS/Coombs″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

II. Antibody screening (IAT)
Use ready-to-use test cell reagents ″ID-DiaCell″.
Allow the test cell reagents and samples to reach room temperature before use.
1. Identify the appropriate microtubes of the ID-Card ″LISS/Coombs″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each test cell reagent to the appropriate microtubes (marked with the corresponding test cell).
4. When an autocontrol is to be included, pipette 50 µl of the sample's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

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III. Antibody Identification (IAT)
Use the ready-to-use test cell reagent ″ID-DiaPanel″.
Allow the test cell reagents and samples to reach room temperature before use.
1. Identify two ID-Cards ″LISS/Coombs″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ″ID-DiaPanel″ test cell to the appropriate microtube (marked 1 to 11).
4. Pipette 50 µl of the sample's own red cell suspension to the 12th microtube (autocontrol).
5. Add 25 µl of the patient's or donor's plasma or serum to all 12 microtubes.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

IV. Compatibility test
1. Identify the appropriate microtubes of the ID-Card ″LISS/Coombs″ with recipient's and donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor red cell suspensions to the appropriate microtubes.
4. For the autocontrol, pipette 50 µl of the patient's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

INTERPRETATION OF THE RESULTS

A) Principle
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for:
I. Direct antiglobulin test (DAT)
- A negative reaction indicates absence of detectable IgG antibodies or C3d complement component on the red cells.
- A positive reaction (± to ++++) indicates that the patient's red cells are sensitized (red cells coated with IgG antibodies and/or C3d).

II. Antibody screening
- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaCell I-II″ or ″ID-DiaCell I-II-III″ corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more test cells and a negative autocontrol suggest the presence of a specific antibody (see ″REMARKS″).
- A positive reaction with all test cells and a positive autocontrol may be due to non-specific reactions.
- A positive reaction with all test cells and a positive autocontrol but with one or more test cells showing a stronger positive reaction than the autocontrol, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

III. Antibody Identification
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaPanel″ corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present can, in most cases, be identified (autocontrol must be negative).
- A positive reaction with all ″ID-DiaPanel″ test cells and a negative autocontrol may be due to non-specific reactions or may indicate the presence of an alloantibody directed against a high frequency antigen.
- A positive reaction with all ″ID-DiaPanel″ test cells and a positive autocontrol may be due to non-specific reactions.
- A positive reaction with all ″ID-DiaPanel″ test cells and the autocontrol but with one or more test cells showing stronger reactions than the autocontrol, may indicate an underlying allo-antibody and further investigation should be undertaken.

IV. Compatibility test
- A negative reaction indicates compatibility of the donor blood with the recipient.
- A positive reaction indicates incompatibility of the donor blood with the recipient, due to presence of antibodies directed against antigens on the donor red cells. Further investigation to identify the antibody specificity should be performed.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performance evaluation studies have been realized with ID-Cards ″LISS/Coombs″. Results obtained showed that performance characteristics of each application of the card are in conformity with the application field and comparable to other approved products.

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ″LISS/Coombs″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

REMARKS

As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Bio-Rad Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
d) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either in vivo or *in vitro* may react with all human sera. Further investigation of such reactions is required.
e) Bacterial or other contamination of materials used can cause false positive or false negative results.
f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
h) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
i) The use of other test cell reagents than ″ID-DiaCell″ or ″ID-DiaPanel″ may modify the reaction patterns.
j) Too heavy or too weak red cell suspensions can cause aberrant results.

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English

B004024 v.01 2020-12-18

Direct and indirect antiglobulin test

Product-Identification: 50540
ID-Card ″Coombs Anti-IgG″

4 x 12. REF 004024
24 x 12. REF 004027
60 x 12. REF 004026
112 x 12. REF 004025

INTRODUCTION

The most important function of the anti-human globulin (AHG) reagent is to detect the presence of IgG antibodies. The importance of anticomplement in the AHG reagent is debatable since clinically significant antibodies detectable only by their ability to bind complement are rare. Some workers therefore prefer to use anti-IgG AHG serum because there is no interference from complement components which may be non-specifically bound to the red cell.


For the indirect antiglobulin test (IAT) with the ID-Card ″Coombs Anti-IgG″, labor-intensive washing procedures are eliminated, due to the fact that the red cell suspension is added first to the microtube, creating a barrier over the gel suspension, thus avoiding neutralisation of the AHG by immunoglobulins in the patient's plasma/serum.

The ID-Card ″Coombs Anti-IgG″ allows the testing of 6 samples, for direct and indirect antiglobulin test (DAT/IAT), antibody screening, antibody identification and crossmatch.

REAGENTS IVD

ID-Card ″Coombs Anti-IgG″ with 6 microtubes containing anti-human globulin anti-IgG (rabbit), within the gel matrix.
Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″Coombs Anti-IgG″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.
When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

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PREPARATION OF BLOOD SAMPLE

a) Red cell
For crossmatch, autocontrol and direct antiglobulin test prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow reagents and plasma or serum samples to reach room temperature before use.

I. Direct antiglobulin test (DAT)
1. Identify the appropriate microtube of the ID-Card ″Coombs Anti-IgG″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

II. Antibody Screening (IAT)
Use the ready-to-use test cell reagents ″ID-DiaCell I, II″, ″ID-DiaCell I, II, III″ or ″ID-DiaScreen I, II, III, IV″.

1. Identify the appropriate microtube of the ID-Card ″Coombs Anti-IgG″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each test cell reagent to the appropriate microtubes (marked with the corresponding test cell).
4. When an autocontrol is to be included, pipette 50 µl of the sample's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

III. Antibody Identification (IAT)
Use the ready-to-use test cell reagent ″ID-DiaPanel″.

1. Identify two ID-Cards ″Coombs Anti-IgG″ and identify with the patient's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each ″ID-DiaPanel″ test cell to the appropriate microtube (marked 1 to 11).
4. Pipette 50 µl of the sample's own red cell suspension to the 12th microtube (autocontrol).
5. Add 25 µl of the patient's or donor's plasma or serum to all 12 microtubes.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

IV. Crossmatch (IAT)
1. Identify the appropriate microtubes of the ID-Card ″Coombs Anti-IgG″ with the recipient's and donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor red cell suspensions to the appropriate microtubes.
4. For the autocontrol, pipette 50 µl of the patient's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

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B) Reactions for:
I. Direct antiglobulin test (DAT)
- A positive reaction (± to ++++) indicates that the patient's red cells are IgG sensitised.
- A negative reaction indicates absence of detectable IgG antibodies on the red cells.

II. Antibody screening
- A negative reaction indicates the absence of detectable irregular antibodies in the patient's serum/plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions observed on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaCell I, II″, ″ID-DiaCell I, II, III″ or ″ID-DiaScreen I, II, III, IV″ corresponds to the lot number indicated on the antigen table. Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with some screening cells and a negative autocontrol indicate presence of a specific antibody.
- A positive reaction with all screening cells and a positive autocontrol are probably due to autoantibody or non-specific reactions.
- Where there is a positive reaction with all screening cells and a positive autocontrol but with one or several screening cells showing a stronger positive reaction than the autocontrol, it should be assumed that, besides the non-specific reaction, a specific antibody may be present and the patients sample should be submitted to further testing.

III. Antibody Identification
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions observed on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaPanel″ corresponds to the lot number indicated on the antigen table. Following the reaction pattern and the antigen configuration, the specificity of the antibody present can, in most cases, be identified (autocontrol must be negative).
- Positive reactions with all ″ID-DiaPanel″ test cells and a negative autocontrol are probably due to nonspecific reactions or may indicate the presence of an alloantibody directed against a high frequency antigen.
- If the autocontrol is also positive, this may be due to unbound autoantibodies or non-specific reaction.

IV. Crossmatch
- A negative reaction indicates compatibility of the donor blood with the recipient.
- A positive reaction indicates incompatibility of the donor blood with the recipient. ABO and RhD blood groups of both recipient and donor must be verified and antibody screening/identification be performed. For complete ABO and RhD blood grouping procedure see related instruction sheets.

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LISS/Coombs + Enzyme Test

English B004514 v.01 2020-02-15

3 x AHG tests/3 x enzyme tests
Antibody screening

Product-Identification: 50581
ID-Card ″LISS/Coombs + Enzyme Test″

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It is generally considered most effective to perform screening tests by both anti-human globulin (AHG) and enzyme test procedures.

The most important function of the polyspecific anti-human globulin (AHG) reagent is to detect the presence of IgG antibodies. The importance of anticomplement (C3d) in the AHG reagent is debatable as antibodies detectable only by their ability to bind complement are rather rare.


Enzyme techniques are useful when increased sensitivity in antibody screening is desired or where more than one antibody may be present. They enhance the reactions of certain antibodies, notably in the Rh, Kell and Kidd systems, whereas antibodies to enzyme-sensitive antigens may not be detected, notably Duffy, M (MNS1), N (MNS2) and S (MNS3).

The ID-Card ″LISS/Coombs + Enzyme Test″ offers complete antibody screening for one donor or patient in an easy single step, when used in conjunction with ″ID-DiaCell I-II-III″ test cell reagents for the indirect antiglobulin test (IAT) procedure and ″ID-DiaCell IP-IIP-IIIP″ papainized test cell reagents for the two-stage enzyme technique.

REAGENTS IVD

ID-Card ″LISS/Coombs + Enzyme Test″: with 3 microtubes containing polyspecific anti-human globulin (AHG) serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9), within the gel matrix, for the AHG test, and 3 microtubes containing neutral gel for the enzyme test.
Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″LISS/Coombs + Enzyme Test″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-DiaCell I-II-III for IAT.
- ID-DiaCell IP-IIP-IIIP for enzyme test.
- ID-Diluent 2: Modified LISS for red cell suspensions (autocontrol).
- ID-Papain: papain solution ready to use for autocontrol.

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SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.
When the use of serum instead of plasma is required, the serum must be well cleared by centrifugation at 1500 g for 10 minutes before use to avoid fibrin residues, which may interfere with the reaction pattern.

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b) Plasma or serum
Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under ″Sample material″) for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Take off the aluminium foil from as many microtubes as needed by holding the ID-Card in the upright position.
3. Resuspend the test cells gently and pipette 50 µl of ″ID-DiaCell I-II-III″ and ″ID-DiaCell IP-IIP-IIIP″ to the appropriate microtubes (1 to 3 and 4 to 6 respectively and marked to indicate the corresponding test cell).
4. Add 25 µl of the patient's or donor's plasma or serum to each mictotube.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for antibody screening
- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction by either IAT or the enzyme test indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents ″ID-DiaCell I-II-III″ or ″ID-DiaCell IP-IIP-IIIP″ corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more antibody screening cells by either IAT or the enzyme test and a negative autocontrol suggest the presence of a specific antibody (see REMARKS).
- A positive reaction with all antibody screening cells by either IAT or the enzyme test and a positive autocontrol by the same technique may be due to an autoantibody.
- Where there is a positive reaction with all antibody screening cells and a positive autocontrol but with one or more antibody screening cells showing a stronger positive reaction than the autocontrol, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

REMARKS

I. Autocontrol
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Take off the aluminium foil from as many microtubes as needed by holding the ID-Card in the upright position.
3. Pipette 50 µl of the patient's or donor's red cell suspension to the appropriate microtube.
4. Add 25 µl of the patient's or donor's plasma or serum to the appropriate microtube.
5. Add 25 µl of ID-Papain to the microtube(s) for enzyme test.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

- If the autocontrol is negative, the positive reactions with the test cell reagents may be considered valid and further tests should be performed to identify the antibody.
- It should be noted that enzyme autoantibodies may react less strongly with the patient's cells in the one-stage technique than against ″ID-DiaCell IP-IIP-IIIP″ (two-stage technique). An autoantibody may mask the presence of other alloantibodies.

II.
As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Some antigens are destroyed or weakened by enzyme treatment. As a result, antibodies against these antigens may not be detected by enzyme procedures {e.g. M (MNS1), N (MNS2), S (MNS3), Fy^a (FY1) and Fy^b (FY2)}.
c) Enzyme test procedures are known to detect more autoantibodies. It is therefore important to include autocontrols when enzyme testing is performed.
d) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
e) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
f) Bacterial or other contamination of materials used can cause false positive or false negative results.
g) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
h) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
i) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
j) Use of test cells other than the ″ID-DiaCell″ may modify the reaction patterns.
k) Too heavy or too weak red cell suspensions can cause aberrant results.

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Product-Identification: 50870
ID-Card ″DAT IgG-Dilution″ 1 x 12. REF 004033

INTRODUCTION

Red cells sensitised with IgG may be subject to accelerated clearance by phagocytosis and/or antibody-dependent cell-mediated cytotoxicity (ADCC), resulting in extravascular hemolysis in the spleen or liver.

As an adjunct to other laboratory tests in the investigation of a positive Direct Antiglobulin Test (DAT), it can be useful to assess the relative amount of IgG sensitising the cells.


The numbers of IgG molecules per cell influence the accelerated red cell destruction seen in Autoimmune Hemolytic Anemia (AIHA), Hemolytic Disease of the newborn (HDN) and Transfusion Reactions.

The ID-Card ″DAT IgG-Dilution″ provides an indication of the clinical importance of the DAT positive result.

REAGENTS IVD

ID-Card ″DAT IgG-Dilution″ with 5 microtubes containing 5 different dilutions (1:10, 1:30, 1:100, 1:300 and 1:1000) of Anti-IgG (rabbit), within the gel matrix and the negative control.
Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
- ID-DiaCell as appropriate (see Test procedure B), ready to use.
- 0.9% isotonic saline solution.

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PREPARATION OF BLOOD SAMPLE

a) Red cells for direct antiglobulin test (DAT)
The red cells should be washed 3 times prior to the preparation of the suspension (see Remarks 2).
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

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b) Plasma or serum for indirect antiglobulin test (IAT)
Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under ˝SAMPLE MATERIAL˝) for a maximum of 48 hours, thereafter at -20 °C.

TEST PROCEDURE

A) Determination of IgG by a positive direct antiglobulin test
Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow reagents to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to each microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

B) Determination of IgG in the case of a positive indirect antiglobulin test
Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow reagents to reach room temperature before use.

I. Screening procedure with the appropriate ID-DiaCell reagent possessing the corresponding antigen.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the ID-DiaCell into each microtube.
4. Add 25 µl plasma or serum of the patient or of donor to each microtube.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

II. Reference procedure for anti-D based on sensitivity cut off by 1000 IgG1 molecules/125 IgG3 molecules.

1. Prepare a 5% red cell suspension of O ccD.Ee cells by adding 25 µl of packed cells to 0.5 ml isotonic saline solution.
2. Pipette 0.2 ml of this cell suspension into a clean test tube.
3. Add 0.4 ml of patient or of donor serum (containing anti-D).
4. Incubate for 30 minutes at 37 °C.
5. Wash 3 times with physiological saline. Decant completely after last wash.
6. Resuspend cell button in 1.0 ml ID-Diluent 2.
7. Identify the ID-Card with the unique patient or donor number/details as appropriate.
8. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
9. Pipette 50 µl of the red cell suspension to each microtube of the ID-Card ˝DAT IgG-Dilution˝.
10. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
11. Read and record the results.

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B) Reactions for IgG in the ID-Card ˝DAT IgG-Dilution˝

Dilutions	1:10	1:30	1:100	1:300	1:1000	ctl
	2+	±	-	-	-	-

Titre anti-IgG 1:30, not relevant, determination of IgG subclasses not necessary.

Dilutions	1:10	1:30	1:100	1:300	1:1000	ctl
	3+	3+	2+	1+	-	-

Titre anti-IgG 1:300, clinically relevant, determination of IgG subclasses necessary, risk of hemolysis.

REMARKS

1. It is important to note that the risk calculation for hemolysis or HDN must be combined with the patient's clinical condition, medical history and other laboratory data.
2. Where a positive reaction occurs in the negative control (ctl) proceed as follows: Wash the red cells with warm (37 °C) 0.9% isotonic saline solution (or ID-Diluent 2) then proceed as described under ˝PREPARATION OF BLOOD SAMPLE˝.
3. If the negative control (ctl) becomes positive when determining IgG subclasses from a positive indirect antiglobulin test due to a strongy agglutinating antibody (for example anti-D) one should first perform a titration in a neutral card to establish the dilution where spontaneous aggregation no longer occurs. This dilution should be used for the determination of the IgG titer.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
f) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
g) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

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Product-Identification: 50890
ID-Card ″DAT IgG1/IgG3″ 1 x 12. REF 004043

INTRODUCTION


The risk of hemolysis by AIHA or HDN depends on the amount of IgG and/or complement coated on the red cells as well as on the IgG subclasses involved. IgG1 and IgG3 may activate complement, while IgG2 and IgG4 rarely do so. To induce phagocytosis, the cells must be coated with 1000–4000 IgG1 molecules or 135–500 IgG3 molecules. Because of the functional differences between the subclasses, determination of the IgG subclasses is important.

In the ID-Card ″DAT IgG1/IgG3″, 2 dilutions of both anti-IgG1 and IgG3 are added to the gel, to allow the differentiation between low and high risk for hemolysis. A positive reaction with the first dilution has a sensitivity of approximately 1000 IgG1 molecules/cell and 125 anti-IgG3 molecules respectively. A positive reaction with the second dilution indicates a high concentration on anti-IgG1 and/or IgG3 antibodies.

To further estimate the risk of hemolysis, the titration of IgG should be performed by using the ID-Card ″DAT IgG-Dilution″ *(see related package insert)*.

REAGENTS IVD

ID-Card ″DAT IgG1/IgG3″ with 6 microtubes containing monoclonal Anti-IgG1 (cell line M345/795) in 2 different dilutions, Anti-IgG3 (cell line M346/805) in 2 different dilutions, Anti-IgG 1:10 (rabbit), within the gel matrix and the negative control. Preservative: < 0.1% NaN_3.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: Modified LISS for red cell suspensions.
- ID-DiaCell as appropriate (see Test procedure B), ready to use.
- 0.9% isotonic saline solution.
- O R_2r (ccD.Ee) cell.

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When the use of serum instead of plasma is required, the serum must be well cleared by centrifugation at 1500 g for 10 minutes before use to avoid fibrin residues, which may interfere with the reaction pattern.

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PREPARATION OF BLOOD SAMPLE

a) Red cell suspension for direct antiglobulin test (DAT)
The red cells should be washed 3 times prior to the preparation of the suspension (see REMARKS 2).
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of washed packed red cells, mix gently.

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b) Plasma or serum for indirect antiglobulin test (IAT)
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow reagents to reach room temperature before use.

A. Determination of IgG by a positive direct antiglobulin test
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to each microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

B. Determination of IgG in the case of a positive indirect antiglobulin test
I. Screening procedure with the appropriate ID-DiaCell reagent possessing the corresponding antigen.
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the ID-DiaCell into each microtube.
4. Add 25 µl plasma or serum of the patient to each microtube.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

II. Reference procedure for anti-D based on sensitivity cut off by 1,000 IgG1 molecules/125 IgG3 molecules.
1. Prepare a 5% red cell suspension of O ccD.Ee cells by adding 25 µl of packed cells to 0.5 ml isotonic saline solution.
2. Pipette 0.2 ml of this cell suspension into a clean test tube.
3. Add 0.4 ml of patient or donor serum (containing anti-D).
4. Incubate for 30 minutes at 37 °C.
5. Wash 3 times with physiological saline solution. Decant completely after last wash.
6. Resuspend cell button in 1.0 ml ID-Diluent 2.
7. Identify the ID-Card with the unique patient or donor number/details as appropriate.
8. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
9. Pipette 50 µl of the red cell suspension to each microtube of the ID-Card ˝DAT IgG1/IgG3˝.
10. Centrifuge the ID-Cards for 10 minutes in the ID-Centrifuge.
11. Read and record the results.

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B) Reactions * for Anti-IgG1/Anti-IgG3 subclasses

	IgG1		IgG3		ctl	IgG
	1:1	1:100	1:1	1:100		1:10
reactions:	2+	-	-	-	-	3+
interpretation:	moderate risk of hemolysis (subclass IgG1, concentration low)					
reactions:	3+	2+	-	-	-	3+
interpretation:	high risk of hemolysis (subclass IgG1, concentration high)					
reactions:	3+	2+	2+	-	-	3+
interpretation:	high risk of hemolysis (subclass IgG1 concentration high, IgG3 concentration low)					
reactions:	-	-	3+	2+	-	3+
interpretation:	high risk of hemolysis (subclass IgG3, concentration high)					
reactions:	-	-	2+	-	-	3+
interpretation:	moderate risk of hemolysis (subclass IgG3, concentration low)					
reactions:	2+	-	3+	2+	-	3+
interpretation:	high risk of hemolysis (subclass IgG1 concentration low, IgG3 concentration high)					
reactions:	2+	-	2+	-	-	3+
interpretation:	moderate risk of hemolysis (subclasses IgG1 and IgG3, concentration low)					
reactions:	3+	2+	3+	2+	-	3+
interpretation:	high risk of hemolysis (subclasses IgG1 and IgG3, concentration high)					

* Reactions indicated in this table are generally the maximum grade of reaction expected but may be weaker.

REMARKS

1. It is important to note that the risk calculation for hemolysis or HDN must be combined with the patient's clinical condition, medical history and other laboratory data. Reaction grades may differ according to the concentrations of each IgG subclass present.
2. Where a positive reaction occurs in the negative control (ctl) proceed as follows: wash the red cells with warm (37 °C) 0.9% isotonic saline solution (or ID-Diluent 2) then proceed as described under ″PREPARATION OF BLOOD SAMPLE″.
3. If the negative control (ctl) becomes positive when determining IgG subclasses from a positive indirect antiglobulin test due to a strongly agglutinating antibody (for example anti-D) one should first perform a titration in a neutral card to establish the dilution where spontaneous aggregation no longer occurs. This dilution should be used for the determination of the IgG titer.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Fibrin residues in the serum or the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
f) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
g) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
h) Too heavy or too weak red cell suspensions can cause aberrant results.

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BIBLIOGRAPHY

1. Zupanska B. et al; Phagozytosis of erythrocytes sensitized with known amounts of IgG1 and IgG3 anti-Rh antibodies. Vox Sang 53 (1987) 96-101.
2. Lynen R. et al; Flow cytometric analyses of the subclasses of red cell IgG antibodies. Vox Sang 69 (1995) 126-130.
3. Lynen R., B.Bernbeck et al: Flow cytometric determination of IgG subclasses and quantitation of red cell allo- and autoantibodies. p.49-58 in: Gutensohn K, Sonneborn H-H, Kühnl P (eds.): Aspects of the Flow-Cytometric Analysis of Red Blood Cells. Clin Lab Publications, Heidelberg 1997.
4. Lynen R, Bernbeck B, Legler TJ, Köhler M: Studies on the sensitivity of IgG red-cell antibody determination using the gel centrifugation test and flow cytometry.
5. New Dimensions for Gel Test Procedures; Satellite Symposium during the ISBT Meeting in Frankfurt, October 1st 1997.
6. Lynen R. et al; ISBT Meeting in Vienna July 9-14 2000.

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English B004704 v.01 2020-02-15

Detection of irregular antibodies with the ID-Card ˝DiaScreen˝

Product-Identification: 50571
ID-Card ˝DiaScreen˝

4 x 12. REF 004704
24 x 12. REF 004707
60 x 12. REF 004706
112 x 12. REF 004705

INTRODUCTION

The reliability of antibody detection is largely dependent on the availability of test cells with appropriate antigens and on the sensitivity of the test methods used.

The requirements for antigen configuration are stringent: it must allow the safe detection of all clinically significant antibodies. For the Rh system, MNSs, Duffy and Kidd, the antigens must be in homozygous form. The Lewis antigens must be present, as should the rare antigen Kp^a.

It is generally considered most effective to perform screening tests by both anti-human globulin (AHG) and enzyme test procedures. Due to higher sensitivity of the indirect antiglobulin test (IAT) with new procedures such as the Bio-Rad ID-System, some scientists in various countries have formed the opinion that the enzyme test has become somewhat less important.

However, enzyme techniques are useful when increased sensitivity in antibody screening is desired or where more than one antibody may be present. They enhance the reactions of certain antibodies, notably in the Rh, Kell and Kidd systems, whereas antibodies to enzyme-sensitive antigens may not be detected, notably in the Duffy and MNS Systems.


Considering these requirements, Bio-Rad has developed a new combination of test cells for use with a single ID-Card ˝DiaScreen˝ consisting of 4 microtubes containing polyspecific AHG for the IAT and 2 microtubes containing neutral gel for the two-stage enzyme technique.

The ID-Card ˝DiaScreen˝ therefore maintains the enzyme procedure, whilst maximising the detection of IgG antibodies by the IAT.

REAGENTS IVD

ID-Card ˝DiaScreen˝ with 4 microtubes containing gel with polyspecific anti-human globulin (AHG) serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) for the IAT and 2 microtubes containing neutral gel for the enzyme test.
Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ˝DiaScreen˝ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension (Id-n°: 05761).
- D-Papain: Liquid papain solution for papainizing red cells for autocontrol (Id-n°: 06311).
- ID-Cards ˝LISS/Coombs˝ (Id-n°: 50531) and ˝NaCl, Enzyme test and cold agglutinins˝ for autocontrol (Id-n°: 50520).
- ˝DiaScreen˝ test cells from single donations, consisting of 6 red cell suspensions at 0.8% (±0.1%) in 10 ml vials (Id-n°: 45070).

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PREPARATION OF BLOOD SAMPLE

a) Red cell suspension for autocontrol
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

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b) Plasma or serum
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

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Allow test cell reagents to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Resuspend the ″DiaScreen″ test cells gently and pipette from each of the vials I-VI 50 µl to the appropriate microtube of the ID-Card (marked with the corresponding test cell number).
4. Add 25 µl of the patient's or donor's plasma or serum to each microtube of the ID-Card.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for antibody screening
- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction by either IAT or the enzyme test indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents ″DiaScreen″ I-VI corresponds to the lot number indicated on the antigen table. Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more test cells by either IAT or the enzyme test and a negative autocontrol suggest the presence of a specific antibody (see REMARKS).
- A positive reaction with all test cells by either IAT or the enzyme test and a positive autocontrol may be due to an autoantibody.
- Where there is a positive reaction with all test cells by either IAT or the enzyme test and a positive autocontrol but with one or more test cells showing a stronger positive reaction than the auto-control, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

REMARKS

I. Autocontrol

1. Identify the appropriate microtubes of the ID-Card ″LISS/Coombs″ and/or "NaCl, Enzyme test and cold agglutinins″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the patient's or donor's red cell suspension to the appropriate microtube.
4. Add 25 µl of the patient's or donor's plasma or serum to the appropriate microtube.
5. Add 25 µl of ″ID-Papain″ to the microtube of the ID-Card ″NaCl, Enzyme test and cold agglutinins″ for the enzyme test.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

- If the autocontrol is negative, the positive reactions with the test cell reagents may be considered valid and further tests should be performed to identify the antibody.
- It should be noted that enzyme autoantibodies may react less strongly with the patient's cells in the one-stage technique than against ″DiaScreen cells V and VI″ (two-stage technique). An autoantibody may mask the presence of other alloantibodies.

II. As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Bio-Rad Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

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LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Some antigens are destroyed or weakened by enzyme treatment. As a result, antibodies against these antigens may not be detected by enzyme procedures [e.g. M (MNS1), N (MNS2), S (MNS3), Fy^a (FY1) and Fy^b (FY2)].
c) Enzyme test procedures are known to detect more autoagglutinins. It is therefore important to include autocontrols when enzyme testing is performed.
d) Certain drugs are known to cause positive anti-human globulin reactions.
e) Some pathological conditions are also reported as causing positive anti-human globulin tests.
f) Bacterial or other contamination of materials used can cause false positive or false negative results.
g) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifuation.
h) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
i) Use of test cells other than the ˝DiaScreen˝ I-VI may modify the reaction patterns.
k) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
l) Too heavy or too weak red cell suspensions can cause aberrant results.

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B004851 v.01 2020-12-18

IgG, IgA, IgM, C3c, C3d, ctl
Direct antiglobulin test (DAT)

Product-Identification: 50830
˝DC-Screening I˝

1 x 12. REF 004851
24 x 12. REF 004857
60 x 12. REF 004856


INTRODUCTION

A positive direct antiglobulin test (DAT) with polyspecific antihuman globulin (AHG) generally indicates that the red cells are coated *in vivo* with immunoglobulin and/or complement. To differentiate the reaction, monospecific AHG reagents are used, such as anti-IgG, -IgA, -IgM, -C3, -C3b, -C3c, -C3d, and -C4.

The configuration of the ID-Card ˝DC Screening I˝ is structured to provide clinically important answers. It consists of anti-IgG, - IgA, -IgM, -C3c and -C3d. It allows this important investigation to be performed in one easy step.

REAGENTS IVD

ID-Card ˝DC-Screening I˝ consisting of five different monospecific AHG reagents: anti-IgG, anti-IgA, anti-IgM, anti-C3c (all rabbit) and anti-C3d (monoclonal, cell line C139-9), suspended in gel; and the negative control.

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PREPARATION OF BLOOD SAMPLE

The red cells should not be washed prior to the preparation of the suspension (see ″REMARKS 1″).

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Hemolytic Disease of the newborn (HDN) results in the neonatal rbcs being coated *in vivo* by maternally-derived IgG antibody.

It is important to note that a positive DAT can be due to other phenomena, and that the patient's clinical condition, medical history and other laboratory data must be considered in parallel with the DAT results.

A positive DAT does not always indicate shortened survival of the rbc's, whereas a negative DAT does not rule out the possibility of a hemolytic process, as the causative antibody may be of low-affinity, or may be due to cell-bound immunoglobulin being below the sensitivity threshold of the test system.

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REMARKS

1. Where there is a positive reaction in the negative control (ctl), proceed as follows:
 Wash the red cells with warm (37 °C) 0.9% isotonic saline solution (or ID-Diluent 2), then proceed as described under ˝PREPARATION OF BLOOD SAMPLE˝ and ˝TEST PROCEDURE˝.
2. Certain drugs are known to cause a positive DAT.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Fibrin residues in the serum or red cell aggregates in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
c) Use of suspension solutions for red cells other than ID-Diluent 2 may modify the reactions.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regulary according to GLP procedures.
f) Too heavy red cell suspensions can cause aberrant results.

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B004831 v.01 2020-12-18

IgG, C3d, ctl/IgG, C3d, ctl
Direct antiglobulin test (DAT)

Product-Identification: 50560
ID-Card ″DC-Screening II″

1 x 12. REF 004831
24 x 12. REF 004837
60 x 12. REF 004836

INTRODUCTION

Polyspecific antihuman globulin (AHG) reagents are used for routine alloantibody detection, compatibility tests and the direct antiglobulin test (DAT). Polyspecific AHG contains antibody to human IgG and the C3d component of human complement, and may also react weakly with IgA and IgM molecules due to anti-light chain activity.

A positive DAT generally indicates that the red cells are coated *in vivo* with immunoglobulin and/or complement.


The most important function of the AHG reagent is to detect the presence of IgG-antibodies. The importance of anti-complement in the AHG reagent is debatable since antibodies detectable only by their ability to bind complement are rather rare. However, anti-C3d activity is important for the direct antiglobulin test in the investigation of autoimmune hemolytic anemia (AIHA).

The configuration of the ID-Card ″DC-Screening II″ allows differentiation between IgG and C3d coating of red cells during investigation of a positive DAT with polyspecific AHG, in one easy step.

REAGENTS IVD

ID-Card ″DC-Screening II″ consisting of monospecific AHG reagents:
- 2 microtubes with anti-IgG (rabbit)
- 2 microtubes with anti-C3d (monoclonal, cell line C139-9), suspended in gel
- 2 microtubes for negative control

The anti-human globulin IgG used in the ID-Card ″DC-Screening II″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: Modified LISS for red cell suspensions.

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PREPARATION OF BLOOD SAMPLE

The red cells should not be washed prior to the preparation of the suspension (see REMARKS 1).

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1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl packed cells, mix gently.

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1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the red cell suspension to the appropriate three microtubes.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

Remark: If you wish to test only one sample, remove the aluminium foil from the first 3 microtubes only. The card may then be used for another test at a subsequent time.

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B) Monospecific reactions
Important: The negative control (ctl) must show a negative reaction.

The results obtained with the monospecific reagents may help to define a clinical condition. Patients with warm antibody auto-immune hemolytic anemia (WAIHA) most often have IgG coating their red cells (rbcs), however C3 complement components may also be detected in about 24% of WAIHA patients. Less commonly, complement alone may be demonstrated on the rbcs (about 7%). IgM and/or IgA immunglobulins may also be present.

In Cold Agglutinin Disease (CAD), the most common result is the detection of complement alone coating the rbcs, but IgM and/or IgA may be detected simultaneously. Paroxysmal Cold Hemoglobinuria (PCH) is also usually demonstrated initially by complement coating the rbcs unless the rbcs are washed with cold (4 °C) saline, when IgG autoantibody may also be detectable. The Donath-Landsteiner test is used to assist in differentiating between CAD and PCH. In mixed-type AIHA, both IgG and C3 are usually demonstrated.

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Hemolytic Disease of the newborn (HDN) results in the neonatal rbcs being coated *in vivo* by maternally-derived IgG antibody. It is important to note that a positive DAT can be due to other phenomena, and that the patient's clinical condition, medical history and other laboratory data must be considered in parallel with the DAT results. A positive DAT does not always indicate shortened survival of the rbc's, whereas a negative DAT does not rule out the possibility of a hemolytic process, as the causative antibody may be of low-affinity, or may be due to cell-bound immunoglobulin being below the sensitivity threshold of the test system.

For identification of IgG autoantibodies, proceed to elution techniques and test the eluate against an appropriate panel of red cell reagents (ID-DiaPanel).

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity
Performance evaluation studies have been realized with ID-Cards ″DC-Screening II″. Results obtained showed that performance characteristics of each application of the card are in conformity with the application field and comparable to other approved products.

Reproducibility
Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards ″DC-Screening II″ have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

REMARKS

1. Where there is a positive reaction in the negative control (ctl), proceed as follows:
 Wash the red cells with warm (37 °C) 0.9% isotonic saline solution (or ID-Diluent 2), then proceed as described under ″PREPARATION OF BLOOD SAMPLE″ and ″TEST PROCEDURE″.
2. Certain drugs are known to cause a positive DAT.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
c) Use of suspension solutions for red cells other than ID-Diluent 2 may modify the reactions.
d) Bacterial or other contamination of materials used can cause false positive or false negative results.
e) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regulary according to GLP procedures.
f) Too heavy or too weak red cell suspensions can cause aberrant results.

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Reverse Grouping with Antibody Screening

English B003014 v.01 2021-02-15

Reverse grouping combined with antibody screening

Product-Identification: 50510

ID-Card ″Reverse Grouping with Antibody Screening″	4 x 12. REF 003014
	24 x 12. REF 003017
	60 x 12. REF 003016
	112 x 12. REF 003015

INTRODUCTION

ABO blood group determination requires that reverse grouping using plasma or serum should be carried out to confirm the accuracy of ABO red cell typing (forward grouping).

Reverse grouping uses red cell reagents of known ABO antigen specificity to indicate the presence or absence of anti-A and anti-B isoagglutinins, the results of which determine the reverse group. Discrepancies between forward and reverse grouping require further investigation.

Good laboratory practice dictates that ABO grouping is incomplete without a screening procedure to detect the presence of irregular antibodies.


The ID-Card ″Reverse Grouping with Antibody Screening″ is used with the test cell reagents ″ID-DiaCell ABO / I-II″ or ″ID-DiaCell ABO / I-II-III″, allowing reverse grouping and antibody screening (indirect antiglobulin test at room temperature) in an easy single test procedure. It provides following combinations of test profiles, according to individual needs:

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REAGENTS IVD

ID-Card ″Reverse Grouping with Antibody Screening″ consisting of 3 microtubes containing neutral gel and 3 microtubes containing anti-human globulin serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix. Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″Reverse Grouping with Antibody Screening″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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PREPARATION OF BLOOD SAMPLE

a) Red cells (for autocontrol or DAT)
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

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CONTROLS

An internal quality control should be included in accordance with the relevant guidelines of quality assurance.

TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow the test cell reagents and samples to reach room temperature before use and resuspend gently.

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each test cell reagent into the appropriate microtubes.
4. Add 50 µl of the patient's serum or plasma to each microtube.
5. An incubation for 10 minutes at room temperature is recommended (see REMARKS 3.).
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

If the test combination is to include an autocontrol, pipette 50 µl of the sample red cell suspension into the appropriate microtube and add 50 µl of the sample plasma or serum. If the test combination is to include a direct antiglobulin test (DAT), pipette 50 µl of the sample red cell suspension into the appropriate microtube (without adding plasma or serum).

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B) Reactions for reverse grouping

Test cell reagents				Blood group
A_1	A_2	B	O	
negative	negative	+ to ++++	negative	**A**
+ to ++++	± to ++++	negative	negative	**B**
negative	negative	negative	negative	**AB**
+ to ++++	+ to ++++	+ to ++++	negative	**O**

Note: Full ABO typing requires that the forward (cell) group determination and reverse (serum / plasma) group show concordance. Where there is a discrepancy, further investigation should be undertaken to establish the cause and ascertain the true ABO-type. Reaction patterns other than those indicated require further investigation.

C) Reactions for antibody screening

- Negative reactions with all antibody screening cells indicate the absence of detectable antibodies.
- Positive reactions with one or more antibody screening cells suggest the presence of irregular antibodies, and further investigation is required to determine the antibody specificity.
- Where there are positive reactions with each antibody screening cell, it is advisable to perform an autocontrol and to proceed to further tests (complete antibody screening at 37 °C), (see REMARKS).
- Where an AHG-microtube was used to perform a DAT, and there is a positive reaction, further investigation to elucidate the cause should be undertaken.

REMARKS

1. A positive autocontrol suggests the presence of an autoantibody; however, appropriate laboratory procedures should be performed to exclude the possible presence of an underlying alloantibody in the patient's plasma / serum.
2. A negative autocontrol with positive antibody screening results indicates the presence of an alloantibody and antibody identification procedures should be undertaken in order to determine the specificity.
3. For reverse grouping, an incubation of at least 10 minutes at 18–25 °C prior to centrifugation will enhance the reactions and minimise repeat testing due to weak isoagglutinins.
4. Antibody screening at room temperature is not recommended for detection of clinically important antibodies in antenatal patients or transfusion recipients.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
b) The optimal reaction conditions may vary considerably from one type of antibody to another. Consequently, one single test method may be insufficient to detect all antibodies and combined test procedures should be applied.
c) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
d) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
e) Bacterial or other contamination of materials used can cause false positive or false negative results.
f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
h) Use of test cells other than the ID-DiaCells may modify the reaction patterns.
i) Too heavy or too weak red cell suspensions can cause aberrant results.

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English B004350 v.01 2020-12-18

Test cell reagents for antibody screening: ID-DiaCell I-II, ID-DiaCell I-II-III, ID-DiaCell IP-IIP-IIIP (papainized), ID-DiaCell I-II-III Asia, ID-DiaCell Pool; for antibody identification: ID-Panel, ID-Panel-P; special antigens: ID-Di^a (Diego) Positive, ID-I Negative Cell

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The test cell reagents are specially designed for the ID-System.

REAGENTS IVD

All test cell reagents are of human origin, in a buffered suspension medium at 0.8% (± 0.1%).
Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

For antibody screening, single donors, blood group O:

ID-DiaCell I-II	$R_1^wR_1$+R_2R_2 for IAT and NaCl test
ID-DiaCell I-II-III	$R_1^wR_1$+R_2R_2+rr for IAT and NaCl test
ID-DiaCell IP-IIP-IIIP	papainized, for enzyme technique
ID-DiaCell Pool	R_1R_1+R_2R_2 (2 pooled cells for donor screening)
ID-DiaCell I-II-III Asia	R_1R_1+R_2R_2+cell of the GP.MUR phenotype, for IAT and NaCl test

For antibody identification, single donors, blood group O:

ID-DiaPanel	11 test cells for IAT and NaCl test
ID-DiaPanel-P	11 test cells papainized, for enzyme technique


Special antigens:
These cells are for use as an adjunct with other routine antibody screening cell sets.

ID-Di^a (Diego) positive
ID-I negative cell

Shipment on standing order, every 4 weeks.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Card ″LISS Coombs + Enzyme Test″ 3 microtubes with polyspecific anti-human globulin AHG) serum and 3 microtubes containing neutral gel (Id-n° 50581).
- ID-Card ″LISS/Coombs″ 6 microtubes with polyspecific AHG serum (Id-n° 50531).
- ID-Card ″Coombs Anti-IgG″ 6 microtubes containing rabbit anti-IgG (Id-n° 50540).
- ID-Card ″NaCl, Enzyme test and cold agglutinins″ 6 microtubes containing neutral gel (Id-n° 50520).
- ID-Card ″Reverse Grouping with Antibody Screening″ 3 microtubes containing neutral gel and 3 microtubes containing polyspecific anti-human globulin (AHG) serum (Id-n° 50510).
- ID-Diluent 2: modified LISS for red cell suspensions (Id-n° 05761).

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CONTROLS

Controls should be included in accordance with the relevant guidelines of quality assurance.

USE OF THE ID-TEST CELL REAGENTS

- All test cell reagents are for use with the ID-Cards of the ID-System only.
- **Strictly follow the test procedures as described in the specific package inserts of the ID-Cards to be used.**
- Always gently resuspend the red cells, by inverting the vial several times before use and also before placing the vials into a pipetting automate.
- Make sure that the test cells are at room temperature (18–25 °C) when in use.
- During the working procedures, check that the test cell reagents remain in suspension. If there is settling of the cells, resuspend again.
- **For the ID-System, precise pipetting is of importance.** Use the ID-Pipetors for serial pipetting.
- Avoid contamination of the test cell reagents.
- When recording the reactions, ensure that the lot number of the antigen table corresponds with the lot number of the reagent vials.
- After use, close the vials and replace them in the refrigerator.

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Note: Test cell reagents for antibody screening containing pooled cells can show a double cell population appearance depending on the antibody present. This is considered to be a positive result.

B) Reactions
See package inserts of the corresponding ID-Cards.

LIMITATIONS

a) Bacterial or other contamination of materials used can cause false positive or false negative results.
b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.

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During antibody identification one very often requires extra cells of certain phenotypes to confirm/exclude possible antibodies. The ″ID-DiaPanel Plus 6″ offers such additional cells.

The product ″ID-DiaPanel plus 6″ is for use in antibody identification, as an addition to the ″ID-DiaPanel″.

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REAGENTS IVD


All test cell reagents are of human origin, in a buffered suspension medium at 0.8% ± 0.1.
Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

For antibody screening and identification, single donors, blood group O:

ID-DiaPanel Plus 6	R_1R_1 ($R_1^wR_1$) / R_1R_1 / R_2R_2 / R_2R_2 / R_ZR_1 and a cell negative for a high frequency blood group antigen or positive for a low frequency blood group antigen.

Shipment on standing order, every 8 weeks.

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BIBLIOGRAPHY

1. Technical Manual of the American Association of Blood Blanks, 13th edition, 1999.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990; 30:109-113.

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ID-DiaScreen

English B004714 v.01 2021-02-15

Detection of irregular antibodies with ˝ID-DiaScreen˝ test cells

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It is generally considered most effective to perform screening tests by both anti-human globulin (AHG) and enzyme test procedures. Due to higher sensitivity of the indirect antiglobulin test (IAT) with procedures such as the ID-System, some scientists in various countries have formed the opinion that the enzyme test has become somewhat less important.

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Considering these requirements, Bio-Rad has developed a combination of test cells for use with a single ID-Card ˝DiaScreen˝ consisting of 4 microtubes containing polyspecific AHG for the IAT and 2 microtubes containing neutral gel for the two-stage enzyme technique.

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
REAGENTS IVD

˝ID-DiaScreen˝ test cells from single donations, consisting of 6 red cell suspensions at 0.8% (± 0.1%) in 10 ml vials.
Preservatives: the antibiotics trimethoprime and sulfamethoxazole.

I..... $R_1^wR_1$
II..... R_2R_2
III..... rr
IV..... not specified
V..... $R_1^wR_1$ papainised
VI..... R_2R_2 papainised

Shipment on standing order every 4 weeks.

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b) Plasma or serum
Where samples are not for immediate testing they should be stored at 2–8 °C after separation for a maximum of 48 hours, thereafter at -20 °C or in accordance with local / national policies / guidelines.

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TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.
Allow test cell reagents to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Resuspend the ″ID-DiaScreen″ test cells gently and pipette from each of the vials I-VI 50 µl to the appropriate microtube of the ID-Card (marked with the corresponding test cell number).
4. Add 25 µl of the patient's or donor's plasma or serum to each microtube of the ID-Card.
5. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
6. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
7. Read and record the results.

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B) Reactions for antibody screening

- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction by either IAT or the enzyme test indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table.
 Verify that the lot number of the test cell reagents ″ID-DiaScreen″ I-VI corresponds to the lot number indicated on the antigen table. Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more test cells by either IAT or the enzyme test and a negative autocontrol suggest the presence of a specific antibody (see REMARKS).
- A positive reaction with all test cells by either IAT or the enzyme test and a positive autocontrol may be due to an autoantibody.
- Where there is a positive reaction with all test cells by either IAT or the enzyme test and a positive autocontrol but with one or more test cells showing a stronger positive reaction than the auto-control, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

REMARKS

I. Autocontrol

1. Identify the appropriate microtubes of the ID-Card ″LISS/Coombs″ and / or ″NaCl, Enzyme Test and Cold Agglutinins″ with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Pipette 50 µl of the patient's or donor's red cell suspension to the appropriate microtube.
4. Add 25 µl of the patient's or donor's plasma or serum to the appropriate microtube.
5. Add 25 µl of ″ID-Papain″ to the microtube of the ID-Card ″NaCl, Enzyme Test and Cold Agglutinins″ for the enzyme test.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

- If the autocontrol is negative, the positive reactions with the test cell reagents may be considered valid and further tests should be performed to identify the antibody.
- It should be noted that enzyme autoantibodies may react less strongly with the patient's cells in the one-stage technique than against ″ID-DiaScreen cells V and VI″ (two-stage technique). An autoantibody may mask the presence of other alloantibodies.

II. As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Anti-D reference reagent (REF 105500) provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

LIMITATIONS

a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
b) Some antigens are destroyed or weakened by enzyme treatment. As a result, antibodies against these antigens may not be detected by enzyme procedures [e.g. M (MNS1), N (MNS2), S (MNS3), Fy^a (FY1) and Fy^b (FY2)].
c) Enzyme test procedures are known to detect more autoagglutinins. It is therefore important to include autocontrols when enzyme testing is performed.
d) Certain drugs are known to cause positive anti-human globulin reactions.
e) Some pathological conditions are also reported as causing positive anti-human globulin tests.
f) Bacterial or other contamination of materials used can cause false positive or false negative results.
g) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifuation.
h) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
i) Use of test cells other than the ″ID-DiaScreen″ I-VI may modify the reaction patterns.
j) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
k) Too heavy or too weak red cell suspensions can cause aberrant results.

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English

B004330 v.01 2021-02-15

Test cell reagents for antibody screening and identification

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The requirements for antigen configuration in antibody screening cells are stringent: they must allow the safe detection of all clinically significant antibodies. For the Rh system, MNSs, Duffy and Kidd, the antigens must be present in double dose. The Lewis antigens and the rare antigen Kp^a are present, when possible.

If an antibody screen is carried out using ″ID-DiaCell I-II-III″ on pregnant women having received anti-D prophylaxis or other patients allo-immunised with anti-D, then 2 out of 3 cells will give a positive result as they are RhD positive.

For such patients, the ″ID-DiaScreen Prophylax″ can be used as a adjunct to ″ID-DiaCell I-II-III″. The presence of other relevant antibodies can be detected using the RhD negative cells. These cells meet the criteria for double dose expression of the Duffy, Kidd and MNS antigens.


REAGENTS IVD

″ID-DiaScreen Prophylax″ test cells from single donations, from group O donors, consisting of 6 red cell suspensions at 0.8% ±0.1 in 10 ml vials. Preservatives: the antibiotics trimethoprime and sulfamethoxazole.

For antibody screening and identification:
ID-DiaScreen Prophylax R_1R_2 / r′r / r″r / rr / rr / rr

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INTERPRETATION OF THE RESULTS

Principle
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

For further comments see package inserts of the corresponding ID-Cards.

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BIBLIOGRAPHY

1. Technical Manual of the American Association of Blood Blanks, 13th edition, 1999.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990; 30:109-113.
3. Shulman, IA., Calderon, C., Nelson, IM., Nakayama, R.: The routin use of Rh-negative reagent red cell for the identification of anti-D red cell antibodies. Transfusion 1994; 34: 666-670.

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INTENDED USE

IH-QC tubes are intended to control the ID-System manually and/or with instruments for the determination of the ABO system, Rh (RH) system, Kell (KEL) system, and/or further immunohematology tests as listed in table 1 (INTENDED USE).

Table 1

	IH-QC							
Analysis, intended use	**IH-QC1**	**IH-QC2**	**IH-QC3**	**IH-QC4**	**IH-QC5**	**IH-QC6**	**IH-QC7**	**IH-QC8**
ABO forward and/or reverse grouping	x	x	x	x	x			
RhD typing	x	x	x	x	x			
Rh/K phenotyping	x	x	x	x	x			
Antibody Screening / Identification, Indirect Antiglobulin Test (IAT)	x	x	x	x	x			
Antibody Screening / Identification, 2 stage papain technique		x	x		x			
Crossmatch, Indirect Antiglobulin Test (IAT)	x	x	x	x	x			
Direct Antiglobulin Test (DAT)	x	x	x	x	x	x	x	x
RhD typing and anti-D for weak D confirmation						x		

SUMMARY

Transfusion guidelines recommend regular checking of test materials, test methods, personnel working procedures and automated equipment/instruments used. The control sample should always have the same characteristics as a patient sample and therefore be treated identically. These activities are intended to ensure accuracy and safety of the blood grouping serological results. The controls should be carried out at regular intervals. Please refer to local/national guidelines.

PRINCIPLE OF THE TEST

Concerning the principle of the test(s); please refer to the respective instruction(s) for use of the product(s) to be controlled according to the listed applications of the different IH-QC tubes in table 1 (INTENDED USE).

REAGENT CHARACTERISTICS

- The tubes of IH-QC1 to IH-QC8 contain 6 ml of a human red blood cell suspension.
- The hematocrit for IH-QC1/IH-QC2/IH-QC3/IH-QC4/IH-QC5 is 15 ±2%.
- The hematocrit for IH-QC6/IH-QC7/IH-QC8 is 10 ±2%.
- The reagents contain bovine albumin.
- Preservatives: Trimethoprim and Sulfamethoxazole.

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Table 2	ABO	RH	K	FY	Antibodies		DAT
IH-QC1	A_1 [ABO:1,-2,3,4]	ddccee (rr) C^w- [RH:-1,-2,-3,4,5,-8]	K+ [KEL:1]	N/A	Anti-B [Anti-ABO2]	Anti-D ≤ 0.05 IU/ml [Anti-RH1]	Neg
IH-QC2	B [ABO:-1,2,3]	DCcEe (R_1R_2) C^w- [RH:1,2,3,4,5,-8]	K- [KEL:-1]	Fy(a-) [FY:-1]	Anti-A [Anti-ABO1]	Anti-Fy^a [Anti-FY1]	Neg
IH-QC3	AB [ABO:1,2,3]	DCCee (R_1R_1) C^w- [RH:1,2,-3,-4,5,-8]	K- [KEL:-1]	N/A	Anti-c [Anti-RH4]		Neg
IH-QC4	O [ABO:-1,-2,-3]	DccEE (R_2R_2) C^w- [RH:1,-2,3,4,-5,-8]	K- [KEL:-1]	N/A	Anti-A and Anti-B [Anti-ABO1 and Anti-ABO2]	Anti-K [Anti-KEL1]	Neg
IH-QC5	A_2 [ABO:1,-2,3,-4]	DCcee (R_1r) C^w- [RH:1,2,-3,4,5,-8]	K- [KEL:-1]	N/A	Anti-B [Anti-ABO2]		Neg
IH-QC6	N/A	Weak D [RH:W1]	N/A	N/A	N/A		Neg
IH-QC7	N/A	N/A	N/A	N/A	N/A		Pos [IgG]
IH-QC8	N/A	N/A	N/A	N/A	N/A		Pos [C3b (c/d)]

FURTHER REAGENT CHARACTERISTICS

IH-QC1: The anti-D ≤ 0.05 IU/ml is standardised against the WHO (World Health Organization) 01/572 International Standard which is controlled by an external accredited reference laboratory by Continuous Flow Assay.
IH-QC2: Anti-Fy^a gives a positive reaction with a titre of ≤ 4 against Fy(a+b+) [FY:1,2] cells in IAT.
IH-QC3: Anti-c gives a positive reaction of ≤ ++ reaction strength against C+c+ [RH:2,4] cells in IAT.
IH-QC4: Anti-K gives a positive reaction with a titre of ≤ 4 against K+k+ [KEL:1,2] cells in IAT.
IH-QC5: No further reagent characteristics.
IH-QC6: The weak D cell is selected based on serologically obtained reactions. In direct testing, the weak D cell can show a reaction strength from + to +++.
By the use of Indirect Antiglobulin Technique (IAT), usually an enhancement is observed.
IH-QC7: Prepared by sensitizing red blood cells with an IgG antibody.
IH-QC8: Prepared by coating human red blood cells with C3b (c/d) using a modified Fruitstone method [1].

STORAGE AND HANDLING

- Store at 2–8 °C.
- Store in an upright position.
- Stability: see expiry date on the label.
- Once opened and if handled in accordance with good laboratory practice (GLP) principles and stored correctly, each tube IH-QC1/IH-QC2/IH-QC3/IH-QC4/IH-QC5 may be used for a maximum of 7 days.
- Once opened and if handled in accordance with good laboratory practice (GLP) principles and stored correctly, each tube IH-QC6/IH-QC7/IH-QC8* may be used for a maximum of 28 days.

* For the IH-QC8, refer to ˝LIMITATIONS˝ for further details.

MATERIAL PROVIDED

- IH-QC1/IH-QC2/IH-QC3/IH-QC4/IH-QC5/IH-QC6/IH-QC7/IH-QC8, see also ˝Product Identification˝
- Worksheet, IH-QC Modular System

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- Material and equipment for the routine test procedures with the ID-System.
- For the use with other reagents and/or instruments, laboratories must follow their approved validation procedures.

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TEST PROCEDURE

- Please refer to the respective test procedure of the ID-System to be controlled.
- Allow the ID-System reagents and IH-QC samples to reach room temperature (18–25 °C) prior to use.
- Treat the IH-QC samples in exactly the same way as regular samples; e.g. in terms of centrifugation.
- After use, place the IH-QC samples back in the refrigerator at 2–8 °C.

INTERPRETATION OF THE RESULTS

- When properly stored and used according to ID-System test procedures, the IH-QC1/IH-QC2/IH-QC3/IH-QC4/IH-QC5/IH-QC6/IH-QC7/IH-QC8 reagents will demonstrate the appropriate characteristics listed in the table under ″REAGENT CHARACTERISTICS″ and in the ″Worksheet, IH-QC Modular System″.
- If the results do not correspond to the chapter ″REAGENT CHARACTERISTICS″ and the ″Worksheet, IH-QC Modular System″, the test methods, working procedures, automated equipment/instruments and materials used should be checked immediately. Repeat the controls after making appropriate corrections.
- Different batches of antibody screening and identification red blood cells can give different reaction strengths with IH-QC1/IH-QC2/IH-QC3/IH-QC4.

PERFORMANCE

Performance evaluation studies have been realized with ″IH-QC Modular System″ for each application claimed for the different IH-QC tubes in Table 1 under ″INTENDED USE″:
ABO forward and/or reverse grouping, RhD typing, Rh/K phenotyping, Antibody Screening / Identification - Indirect Antiglobulin Test (IAT), Antibody Screening / Identification - 2 stage papain technique, Crossmatch - Indirect Antiglobulin Test (IAT), Direct Antiglobulin Test (DAT), RhD Typing and anti-D for weak D confirmation.
Interassays and interlots reproducibility of ″IH-QC Modular System″ have been evaluated. Obtained results have demonstrated that ″IH-QC Modular System″ characteristics (according to Table 1) allow a reproducible control of ID-Systems performances in both manual and automated methods.
As part of the lot release process, each lot of ″IH-QC Modular System″ is tested to ensure the described reagent characteristics.

The final release testing is performed according to the product specific Standard Operating Procedures.

LIMITATIONS

Erroneous and abnormal results may be caused by:
- Non respect of the intended use or shelf life or stability after opening.
- Bacterial or chemical contamination of the ″IH-QC Modular System″.
- In case the IH-QC8 would be tested more than once a day over 28 days, a decrease of reactivity may be noticed.

WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to the GLP procedures.
- Do not use reagents beyond the expiry date.
- Do not freeze or expose reagents to excessive heat.
- Do not use if markedly hemolyzed.
- All blood products should be treated as potentially infectious. Source material, from which this product was derived was found negative for HbsAg, HCV and HIV (1+2) tested with licenced reagents. No known test method can offer assurance that products derived from human blood will not transmit infectious agents.
- The bovine albumin used for the production of this reagent is purchased from BSE-free sources.
- Used test material must be discarded as hazardous materials. Manage waste according to local, state and/or national regulations.

BIBLIOGRAPHY

[1] Fruitstone MJ. C3b sensitized erythrocytes. Transfusion 1978;18:125.

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ID-Internal Quality Control

English B009925 v.01 2020-09-25

Product-Identification: 45341
˝ID-Internal Quality Control˝

Set REF 009925
5 x 4 ml red cell reagents
3 x 3 ml sera

INTRODUCTION

The safety of a patient needing a blood transfusion depends directly on the accuracy and the reliability of the blood group serology results. Consequently, the test material, the procedures and the performance of the laboratory staff should be monitored regularly by an internal quality control system.

˝ID-Internal Quality Control˝ is designed to allow control of the reagents in use and the working procedures with ID-System. The control should be used regularly and by each member of the laboratory staff.

REAGENTS IVD

˝ID-Internal Quality Control˝ contains 5 vials of 4 ml human red cell reagent from single donors as a 3–5% suspension in a buffered medium, and 3 vials with 3 ml serum of human origin (all reagents containing bovine albumin):


- Red cell reagent 1(Cell 1).....AB, RhD positive
- Red cell reagent 2(Cell 2).....O, RhD negative
- Red cell reagent 3(Cell 3).....O, Kell positive, R_1R_2
- Red cell reagent 4(Cell 4).....O, RhD weak
- Red cell reagent 5(Cell 5).....DAT positive

Preservatives: the antibiotics trimethoprim and sulfamethoxazole

- Serum 1 with irregular antibody (blood group AB)
- Serum 2 with irregular antibody (blood group O)
- Serum 3 without irregular antibody (blood group O)

Preservatives: the antibiotics trimethoprim and sulfamethoxazole

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Stability: see expiry date on label.

FURTHER MATERIALS REQUIRED

Exactly the same as used in routine procedures.

PREPARATION OF QUALITY CONTROL SAMPLES

Allow the quality control red cell reagents and sera to reach room temperature before use.

A) Red cell reagents

1. Take the required vial of test cell reagents and gently resuspend the cells.
2. Transfer at least 500 µl to a centrifuge tube.
3. Centrifuge the tube for 1 minute at 125 g or for 20 seconds at 1000 g.
4. Decant the supernatant.
5. From the packed cells, make a cell suspension in ID-Diluent 1 or 2 as appropriate, according to the ID-Card being used for the patient test procedure.

B) Quality control sera

Ready-to-use for reverse grouping, for antibody detection and identification procedures.

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TEST PROCEDURE

Use the quality control samples in the same manner as patient samples in routine work.

INTERPRETATION OF THE RESULTS

Compare the reactions and the results obtained with those on the enclosed table. If the reactions and results do not correspond to those of the table, the working procedures should be verified, along with routine test methods and the materials used. The necessary corrective action should be taken and the tests repeated.

Please note: each lot of ˝ID-Internal Quality Control˝ can give different reactions and results. Verify that the lot numbers indicated on the table correspond to the lot numbers on the vial labels.

PERFORMANCE CHARACTERISTICS

As part of the lot release process, each lot of ˝ID-Internal Quality Control˝ is tested to ensure the described reagent characteristics. The final release testing is performed according to the product specific Standard Operating Procedures.

LIMITATIONS

a) Bacterial or other contamination of materials used can cause false positive or false negative results.
b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
c) Too heavy or too weak red cell suspensions can cause aberrant results.

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A, B, DVI-/Enz./2 x AHG
Compatibility Testing

Product-Identification: 50601

ID-Card ″DiaClon Complete Crossmatch″

4 x 12	REF 004614
24 x 12	REF 004617
60 x 12	REF 004616
112 x 12	REF 004615

INTRODUCTION

The methods used in compatibility tests should include those that will demonstrate ABO and RhD incompatibility and detect clinically significant unexpected antibodies.
The ID-Card ″DiaClon Complete Crossmatch ″ is configured to meet the above requirements in one easy step.

- ABD-Confirmation of both donor and recipient.
- Major compatibility in anti-human globulin (AHG) and enzyme test.
- Autocontrol (AHG test).


Compatibility testing should also include checking of previous transfusion records for the recipient's serological history. Results of current testing must be compared with interpretation of previous testing [1, 3, 4, 5].

REAGENTS IVD

ID-Card ″DiaClon Complete Crossmatch″ with 3 microtubes containing monoclonal anti-A (cell line LM297/628 [LA-2]), monoclonal anti-B (cell line LM306/686 [LB-2]) and monoclonal anti-D [cell lines LHM59/20[LDM3], 175-2] within the gel matrix, for the ABO/D confirmation of both recipient and donor; 1 microtube containing neutral gel for the enzyme compatibility test and 2 microtubes with polyspecific AHG serum (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) for the compatibility test in AHG and the autocontrol (ac).
Preservative: < 0.1% NaN_3.

The anti-human globulin IgG used in the ID-Card ″DiaClon Complete Crossmatch″ is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

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Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
- ID-Papain: Papain solution for enzyme test.
- ID-Diluent 1: modified bromelin solution for enzyme test.

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SAMPLE MATERIAL

Patient
For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

Donor
The donor sample must be from a segment, originally attached to the blood bag, well mixed.

PREPARATION OF BLOOD SAMPLE

a) Patient's and donor's red cells
Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Pipette 10 µl of packed red cells; mix gently. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

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b) Patient's plasma or serum
Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under ˝SAMPLE MATERIAL˝) for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURE

A) with ABO/D test of both donor and patient
Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

Allow diluent and samples to reach room temperature before use.

1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the patient's red cell suspension to microtubes 1, 2, 3 and 6 (A, B, D, ac).
4. Pipette 50 µl of the donor's red cell suspension to microtubes 1, 2, 3, 4 and 5 (A, B, D, Enz., AHG).
5. Add 25 µl of the patient's plasma or serum to microtubes 4 and 5 (compatibility test) and microtube 6 (autocontrol).
6. Add 25 µl of ID-Diluent 1 or ID-Papain to microtube 4 (enzyme test).
7. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
8. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
9. Read and record the results.

B) with ABO/D test of patient only (confirmation of blood group)
1. Identify the ID-Card with the unique patient or donor number/details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the patient's red cell suspension to microtubes 1, 2, 3 and 6 (A, B, D, ac).
4. Pipette 50 µl of the donor's red cell suspension to microtubes 4 and 5 (compatibility test).
5. Add 25 µl of the patient's serum or plasma to microtubes 4 and 5 (compatibility test) and microtube 6 (autocontrol).
6. Add 25 µl of ID-Diluent 1 or ID-Papain to microtube 4 (enzyme test).
7. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
8. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
9. Read and record the results.

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INTERPRETATION OF THE RESULTS

A) Principle [2]
Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for complete compatibility test
I. ABD confirmation of both donor and patient
- ABO/D compatibility is proven when the reaction gives only one cell population: ++++ or negative.
- ABO incompatibility has to be considered when the reaction with anti-A and/or anti-B shows two distinct cell populations: positive and negative (sedimentation of non-agglutinated cells) causing a double cell population phenomenon.
- RhD incompatibility has to be considered when the reaction with anti-D shows two distinct cell populations (positive and negative).
- When there is a possible incompatibility, ABO and Rh blood groups of both recipient and donor need to be retested. For complete ABO and RhD blood grouping procedures see related instruction sheets.

II. ABD confirmation of patient only

Anti-A	Anti-B	Anti-D	Blood group
+++ to ++++	negative		**A**
negative	+++ to ++++		**B**
+++ to ++++	+++ to ++++		**AB**
negative	negative		**O**
		+++ to ++++	**D positive**
		± to ++*	**D weak**
		negative	**D negative**

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Weaker reactions than +++ with anti-A or anti-B can indicate the presence of A or B subgroups. In this case interpretation can only be made after a complete ABO bloodgrouping procedure. See related instruction sheets. In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-D test sera were selected so as **not to react** with DVI variants.

III. Compatibility test
- A positive reaction (± to ++++) in one or both microtubes 4 and 5 (enzyme and AHG test) and a negative reaction in microtube 6 (autocontrol) indicates an incompatibility due to the possible presence of antibodies.
- A negative reaction in the microtubes 4 and 5 (enzyme and AHG test) and a negative reaction in microtube 6 indicates the absence of detectable antibodies, indicates compatibility.
- A positive reaction in microtube 6 only (autocontrol) may indicate a positive direct antiglobulin test (DAT) of the patients cells. The patient's sample should be submitted to further testing.
- The patient's sample should be submitted to further testing if a positive reaction in one or both microtubes 4 and 5 and a positive autocontrol is observed.

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PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performances of the monoclonal antibodies present in the ID-Cards ˝DiaClon Complete Crossmatch˝, as well as the applications compatibility test with enzyme, compatibility test on AHG well and autocontrol, have been evaluated according to the Common Technical Specification (CTS) [6] for *in vitro* diagnostic medical devices.

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Antibodies	Total number of samples	Sensitivity	Specificity
Anti-A, Cell line: LM297/628 (LA-2)	4566	100%	100%
Anti-B, Cell line: LM306/686 (LB-2)	4566	100%	100%
Anti-D(DVI-), Cell lines: LHM59/20(LDM3), 175-2	4679*	99.72%**	100%

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Compatibility test and autocontrol

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Applications	Total number of samples	Sensitivity	Specificity
Compatibility test enzyme	523	98.91%	100%
Compatibility test AHG	649	99.59%	99.75%
Autocontrol AHG	523	N/A	100%

In addition, the test procedures specific for the ˝DiaClon Complete Crossmatch˝ (A and B) have been evaluated with a minimum of 500 samples coming from donors and patients. The evaluation resulted in 100% of sensitivity and specificity on the concerned wells Anti-A, Anti-B and Anti-DVI-.

REMARKS

1. When there is an incompatibility, antibody screening and identification tests should be performed.
2. Strong cold antibodies can be the cause of weak positive reactions, especially in the enzyme test microtube.
3. The anti-B of monoclonal origin does not react with the acquired B antigen.

LIMITATIONS

a) ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
b) Some antigens are destroyed or weakened by enzyme treatment. As a result, antibodies against these antigens can not be detected by enzyme procedures.
c) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
d) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
e) Bacterial or other contamination of materials used can cause false positive or false negative results.
f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
h) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
i) Too heavy or too weak red cell suspensions can cause aberrant results.

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BIBLIOGRAPHY

1. Holland P. V. ed.; Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington, V.A: American Association of Blood Banks, 1984: 15-39.
2. Lapierre, Y.,Rigal,D.,Adam, J.et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
3. Technical Manual; American Association of Blood Banks. 11th ed. 1993: chapter 8, 175-187.
4. Issit, P. D.: Applied blood group serology, 3rd ed. Miami; Montgomery Scientific Publications, 1985.
5. Mollison, P.L., Engelfriet, C.P. und Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed 1997; Blackwell Scientific Publications, Oxford.
6. Official Journal of the European Union L 318/25: Commission decision of 27 November 2009 amending decision 2002/364/EC on common technical specifications for *in vitro* diagnostic medical devices. (2009/886/EC)

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English

B005652 v.01 2020-12-18

Stabilisation solution for red cells

Product-Identification: 05740
˝ID-CellStab˝

1 x 500 ml REF 005650
2 x 100 ml REF 005660


INTRODUCTION

It is often necessary to preserve red cells, for example for further investigation of anomalous results or for use in antibody screening and/or identification procedures. Therefore it is indispensable that the antigens on the red cells remain stable.

˝ID-CellStab˝ is a specially formulated glycine buffered saline which stabilises red cells at 0.8% concentration for use in the ID-System for 4 weeks after preparation.

REAGENTS IVD

˝ID-CellStab˝; glycine buffered saline.
Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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TEST PROCEDURE

Red Cells for Antibody Screening/Identification

1. Wash the red cells 3 times with ˝ID-CellStab˝.
2. Decant the supernatant completely after the last wash.
3. Resuspend in ˝ID-CellStab˝ to a 0.8% concentration as follows:
 - dispense 5 ml of ˝ID-CellStab˝ into a clean glass tube.
 - add 50 µl of the washed packed cells.
4. Store at 2–8 °C.

Note: To prepare papainised cells see package insert for ˝ID-Papain˝. For freezing of cells, see package insert for ˝DiaMed CellFreeze˝ and ˝DiaMed CellThaw˝.

REMARKS

Use sterilised vials and droppers for storage.

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LIMITATIONS

1. Red cells will remain stable for a minimum of 4 weeks after preparation, but may be viable for longer periods if storage conditions are optimal. Appropriate controls should be employed to ensure red cell antigenicity.
2. Bacterial or other contamination may interfere with the red cell stability.
3. Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.

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DiaMed CellFreeze / DiaMed CellThaw

English B005702 v.01 2021-02-15

Freezing solution for long term storage of red cells / Thawing solution for frozen red cells

Product-Identification:

DiaMed CellFreeze, *Freezing solution for long term storage of red cells*	Id-n°: **40230**	2 x 100 ml	REF 005702
DiaMed CellThaw, *Thawing solution for frozen red cells*	Id-n°: **40240**	2 x 100 ml	REF 005712
ID-CellStab, *Stabilisation solution for red cells*	Id-n°: **05740**	1 x 500 ml	REF 005650

INTRODUCTION


The presence of antibodies against high and low incidence antigens often create problems in antibody identification procedures due to the lack of appropriate red cells.

It is therefore important to be able to preserve rare red cells over a long period of time.

Bio-Rad has developed two special buffers for easy freezing and thawing of red cells, allowing the preservation of interesting antigen configurations so that the laboratory has rare red cells available at any time.

REAGENTS IVD

- ″DiaMed CellFreeze″ freezing solution, containing: glycine phosphate saline buffer, with sugars, glycerine and bovine albumin, in 100 ml bottles.
 Preservatives: the antibiotics trimetoprime and sulfamethoxazole.
- ″DiaMed CellThaw″ thawing solution, containing: glycine phosphate saline buffer, with sugars and bovine albumin, in 100 ml bottles.
 Preservatives: the antibiotics trimetoprime and sulfamethoxazole.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ″ID-CellStab″

(see related package insert)

FURTHER MATERIALS REQUIRED

- Glass bottles (100 ml)
- Glass tubes
- 5 ml glass vials with screw caps (sterilized)
- Vial labels
- Freezer (-25 °C to -85 °C)
- Centrifuge for tubes (1000 g)

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Freezing procedure

For long term storage of red cells.

1. Wash the sample red cells 3 times with ″ID-CellStab″, decant the supernatant after the last centrifugation.
2. Add the freezing buffer ″DiaMed CellFreeze″ in the ratio of 100 ml to 1.2 ml washed packed red cells, mix well (for example 25 ml buffer + 0.3 ml packed red cells).
3. Take 5 ml sterilized glass vials and fill them with 4 ml of the red cell suspension (step 2. from above), close the vial and mix well.
4. Label the vials with appropriate identification.
5. Place the vials in the freezer at -25 °C (or down to -85 °C).

Thawing procedure

1. Take the required vial from the freezer and allow it to remain at room temperature (18–25 °C) for 15 minutes.
2. Gently mix the thawed red cell suspension and transfer 2 ml into a clean glass tube. Place the closed vial in the freezer for future use.
3. Add to the tube 2 ml of the thawing buffer ″DiaMed CellThaw″ and mix well.
4. Centrifuge the tube for 1 minute at 1000 g.
5. Decant the supernatant.
6. Add again 2 ml of the thawing buffer, mix, centrifuge and decant the supernatant.
7. Add 2 ml of ″ID-CellStab″, mix well and centrifuge 1 minute at 1000 g.
8. If the supernatant is clear, the red cells are ready for use, otherwise wash once more with ″ID-CellStab″.

The red cell suspension is stable for at least 4 weeks when stored at 2–8 °C.

REMARKS

1. If the red cells are to be used for the ID-System, decant the supernatant after the last centrifugation and with the packed red cells, prepare the cell suspension following the corresponding package insert with ″ID-Diluent 2″ or ″ID-CellStab″.
2. Samples to be frozen should be no more than 3 days old. Older samples may exhibit increased haemolysis on thawing, reducing the amount of red cells recovered for use.

LIMITATIONS

a) The freezing temperature must be at least -25 °C. At higher temperature, i.e. -20 °C, the suspensions do not freeze.
b) Bacterial or other contamination of materials used will shorten the lifespan of the red cells.
c) For freezing, use only glass vials of 5 ml size or smaller, and no other closures but screw caps, to avoid water adsorption during the freezing process (causing haemolysis).
d) Use of other freezing, thawing and stabilisation solutions than described may affect the stability of the red cells.

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English

B009190 v.01 2020-12-18

Modified Bromelin solution for red cell suspensions and enzyme tests

Produkt-Identification: 05751

″ID-Diluent 1″	2 x 100 ml	REF 009160
″ID-Diluent 1″	1 x 500 ml	REF 009180
″ID-Diluent 1″ Rack for **IH**-Analyzers	10 racks with 60 x 700 µl	REF 009190

INTRODUCTION

Proteolytic enzymes modify red blood cell antigens and thus enhance the reactivity of some antigen/antibody systems and suppress others such as M, N, S, Fy^a and Fy^b.

The most commonly used enzymes in blood group serology are papain and bromelin. Papain is generally used for enzyme treatment of red blood cells prior to their use, whereas bromelin is often used as an additive reagent.


″ID-Diluent 1″ is a modified bromelin solution in which the enzyme activity is stabilised for a long period, specially prepared for the ID-System. ″ID-Diluent 1″ is used for preparing suspensions of red cells for blood grouping and as an additive for enzyme tests with untreated red cells for antibody detection and crossmatching.

REAGENTS IVD

″ID-Diluent 1″: Modified Bromelin solution for red cell suspensions and enzyme tests, in 100 and 500 ml bottles or racks of 60 x 700 µl for the **IH**-Analyzers.

Preservatives: the antibiotics trimetoprim and sulfamethoxazole.

Once opened and if handled in accordance with good laboratory practice (GLP) principles and stored correctly, the 100 and 500 ml bottles can be used until the date of expiry or for up to 6 months afterwards, whichever is the sooner.

Stability: see expiry date on label.

FURTHER MATERIALS REQUIRED (FOR 100 AND 500 ML BOTTLES)

- ID-Dispenser
- Pipette 25 and 50 µl
- D-Tips (pipetor tips)
- Suspension Tubes
- ID-Working Table

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PREPARATION OF BLOOD SAMPLE (FOR 100 AND 500 ML BOTTLES)

A) For blood group determinations

Prepare a 5% red cell suspension of the patient's red cells in ″ID-Diluent 1″ as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 0.5 ml of ″ID-Diluent 1″ into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.
3. Incubate the red cell suspension for 10 minutes at room temperature.

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B) For use as enzyme reagent additive (antibody detection, crossmatch)
Allow the diluent to reach room temperature before use.
After pipetting of the 0.8% red cell suspension and the serum or plasma sample to the appropriate microtubes of an ID-Card corresponding to the required test procedure, add 25 µl of ˝ID-Diluent 1˝ to each microtube.

Important: consult the related instruction sheet of the ID-Cards for precise working procedures.

PREPARATION OF BLOOD SAMPLE (RACKS WITH 60 X 700 µL)

Please refer to the user manual of the **IH**-Analyzers.

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ID-Diluent 2

English B009290 v.01 2020-10-07

Modified LISS for red cell suspensions

Product-Identification: 05761

″ID-Diluent 2″	2 x 100 ml	REF 009260
	1 x 500 ml	REF 009280
″ID-Diluent 2″ Rack for **IH**-Analyzers	10 racks with 60 x 700 µl	REF 009290

INTRODUCTION


Low ionic strength solution (LISS) increases the rate of antibody association and thus enhances antigen/antibody reactions.

″ID-Diluent 2″ is a modified low ionic strength solution made for the ID-System, for preparing 5% red cell suspensions for blood grouping as well as 0.8% red cell suspensions for crossmatching, autocontrol, direct anti-human globulin test, blood grouping of newborns and test cells prepared in the laboratory.

REAGENTS IVD

″ID-Diluent 2″: modified LISS for Red Cell Suspension, in 100 and 500 ml bottles or racks of 60 x 700 µl for the **IH**-Analyzers. Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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- ID-Dispenser
- Pipette 25 and 50 µl
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working Table

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A) For blood group determinations
Prepare a 5% red cell suspension of the patient's red cells in ″ID-Diluent 2″ as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 0.5 ml of ″ID-Diluent 2″ into a clean tube
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

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B) For crossmatching, direct anti-human globulin test, autocontrol in antibody screening and other tests as specified in the relevant box inserts
Prepare a 0.8% red cell suspension in ″ID-Diluent 2″ as follows:
Allow the diluent to reach room temperature before use.
1. Dispense 1.0 ml of ″ID-Diluent 2″ into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

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C) For preparation of test cell reagents in the laboratory

1. Wash the red cells with isotonic saline solution or with ″ID-Diluent 2″ three times or more, until the supernatant is clear.
2. Decant the supernatant and resuspend the packed cells to a 0.8% suspension with ″ID-Diluent 2″.
3. Under aseptic conditions place the suspension into sterile glass vials or tubes, close with sterile caps.
 Red cell reagents thus prepared and stored at 2–8 °C are stable for one day.

Important: *consult the related instruction sheet of the ID-Cards for precise working procedures.*

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LIMITATIONS

a) Bacterial or other contamination of materials used can cause false positive or false negative results.
b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
c) ″ID-Diluent 2″ has been formulated in order to reduce the incidence of detection of non-significant red cell bound complement components.
d) Some antibodies in the Kell system may react more weakly in LISS techniques.
e) LISS solutions may affect autoantibody activity and therefore may cause problems with certain samples.

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INTRODUCTION

To prevent bacterial or other contamination of the pipetting system of instruments, and to eliminate the risk of sample cross-contamination, it is essential that the distribution needle and the tubing are constantly washed and rinsed.

Wash Solution A and B should be used according to the instructions described in the user manual of the corresponding Bio-Rad Immunohematology instrument.
Wash Solution A contains isotonic buffer solution and Wash Solution B is a cleaning solution with surfactant agent.

REAGENTS IVD

- **Wash Solution A Concentrate, isotonic buffer solution, 100 ml vial.** Preservative: ≤ 2% NaN_3.

WARNING

Acute Tox. 4 **H302+H312** Harmful if swallowed or in contact with skin.
Aquatic Chronic 3 **H412** Harmful to aquatic life with long lasting effects.

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.
P264 Wash thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P301+P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- **Wash Solution B Concentrate, cleaning solution with surfactant agent, 100 ml vial.** Preservative: ≤ 2% NaN_3.

DANGER

Eye Dam. 1 **H318** Causes serious eye damage.


Acute Tox. 4 **H302+H312** Harmful if swallowed or in contact with skin.
Aquatic Chronic 3 **H412** Harmful to aquatic life with long lasting effects.
Skin Irrit. 3 **H316** Causes mild skin irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.
P301+P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P363 Wash contaminated clothing before reuse.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

The product and the packaging are labeled in accordance with the rules laid down in the international GHS standard. In the case of the EU countries, the hazards listed in Section 2 of the safety data sheet can differ from the labeling of the product. The EU countries classify the hazards in Section 2 according to the CLP standard.

Chemicals must be disposed of in compliance with respective national regulations.

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
PREPARATION OF WORKING SOLUTIONS

- Fill the instrument containers A and/or B with 2 litres of distilled or demineralised water (Corresponding at least to Grad 2 - ISO 3696 or Type II - ASTM D1193-91).
- Add the entire contents (100 ml) of concentrate A to the container A, and concentrate B to the container B.
- Close the containers with the corresponding closure, mix well and connect them to the instrument.

Stability of the working solutions: 7 days

Disposal of the prepared solution
If preparation of the working solution is carried out in accordance with the instructions, the weight percentages of the hazardous material lie below the concentration value requiring declaration in accordance with the provisions for the relevant hazard symbols. The prepared solution can be disposed of in compliance with the respective national regulations.

REMARKS

Storage of Washing solution A at low temperature can cause the formation of crystals. Do not use if crystals are visible.
To dissolve, incubate the washing solution at 37 °C.

Wash and rinse the instrument containers A and B before preparing new solutions.

LIMITATIONS

Use of wash solutions other than the described products may result in malfunction of the instruments (blocking of tubing, sample cross-contamination, haemolysis of blood samples).
Use of distilled or demineralised water corresponding to a grade lower than Grad 2 - ISO 3696 or Type II - ASTM D1193-91 will lead to decrease performance of the solution which may in turn result in cross contamination, erroneous results or damage to the instrument.

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INTRODUCTION


To prevent bacterial or other contamination of the pipetting system of the instrument, and to eliminate the risk of sample cross-contamination, it is essential that the distribution needle and the tubing are constantly washed and rinsed.

Two bottles are connected to the pumping system of the ID-Automate: bottle A contains isotonic buffer solution, bottle B is a cleaning solution with surfactant agent.

REAGENTS

- Wash Solution A, isotonic buffer solution, 10 litres, ready-to-use.
 Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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
Stability: see expiry date on label.

- Wash Solution B, cleaning solution with surfactant agent, 10 litres, ready-to-use.
 Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

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Stability: see expiry date on label.

REMARK

Use of wash solutions other than the described products may cause malfunction of the instruments (blocking of tubing, sample cross-contamination, hemolysis of blood samples).

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DiaMed GmbH