Box size: 7 x 23
Overall size: 18 x 25
Folds to: 3 x 5


COVIDIEN

EEA™ OrVil™ Auto Suture™
Transoral Circular Stapler Anvil
Advancing Proximal Guide Suture

For use with EEA™ XL Circular Stapler with DST Series™ Technology

1306505

Page 1 of 2

EN

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical technique.

Do not use this product unless familiar with the surgical patterns and risks involved. Review or inservicing of this device may lead to its failure and subsequent injury. Reprocessing and/or re-sterilization of this device may lead to the risk of contamination and patient infection. Do not reuse, reprocess or re-sterilize this device.

DESCRIPTION

The EEA™ OrVil™ transoral circular stapler anvil comes with the anvil head secured in the fold position. The anvil assembly is mounted on a 45 cm long PDS™ delivery tube and is secured to the tube with a retention suture. An advancing proximal guide suture is attached to the anvil head to aid in the proximal control. An advancing proximal guide suture is supplied to aid and aid can be displaced as needed. A pledget is provided to aid in the tying of the purse-string suture. A pledget is provided to reinforce a planned, immediate after-staple transection. The stapler knife blade moves the entire linear, cinching circular anastomosis. The diameter of the staple line is 21 mm or 25 mm, depending on the device selected.

The EEA™ OrVil™ anvil is permanently affixed to the center rod and cannot be removed from the center rod.

INDICATIONS

The EEA™ OrVil™ trans-luminal circular stapler anvil when used with the DST Series™ circular stapling device is indicated for use in the creation of an end-to-end anastomosis between two sections of the GI tract. It is intended for applications throughout the abdominal cavity that for the treatment of end-to-end, end-to-side and side-to-side anastomoses in both open and laparoscopic surgeries, including bariatric surgery.

CONTRAINDICATIONS

1. Refer to the instructions for use for the EEA™ XL stapler for specific indications, contraindications, warnings, precautions, and operating instructions.

WARNING AND PRECAUTIONS

1. Do not use this product unless the instructions below fully apply. Refer to the stapler specification sheet for important warnings, precautions and cautions for use in the EEA™ XL stapler instructions for use.
2. EEA™ staplers should not be used on tissue that does not approximate easily.
3. The transoral circular stapler anvil is compatible with the corresponding stapler for EEA™ XL stapler length products (780). The EEA™ OrVil™ anvil 21 mm (aqua) is compatible with the EEA™ XL stapler 21 mm (aqua) and 25 mm (white) anvil. The EEA™ OrVil™ anvil 25 mm (white) is compatible with the EEA™ XL stapler 25mm. They are not compatible with the 21 cm long products or any other EEA™ staplers.
4. Use of an incompatible implant with this device may cause device failure and/or patient injury.
5. Do not use if sterile barrier has been damaged.
6. Do not use if expiration date has been exceeded.
7. Do not attempt to resterilize this device. Resterilization may compromise the integrity of the staple line, resulting in leakage or dissection.
8. Improper cutting and removal of the advancing proximal guide suture can result in excess suture left in the cavity. The Anvil head and Sutures not should be inspected after removal to ensure the advancing proximal guide suture has been removed with the Suture not.
9. Use of the anvil assembly without the locking cap or alternate means to insure the anvil retaining suture has been removed will be defective.
10. Do not use if retaining and advancing proximal guide sutures are not intended for use as a surgical suture in this application, they are only Sutures are to rotate and guide the anvil assembly.

ADVERSE EVENTS

Adverse events such as but not limited to dysphagia, stricture, stenosis, leak, perforation, fistula, bleeding, pain, infection, abscess, and/or death could occur as a result of the surgical procedure.

WARNING: Improper cutting and removal of anvil retaining suture could leave the suture in place, which may result in tissue damage or obstruction of the GI tract.

1. Remove delivery tube with retention suture attached. Removal of the tube will pull the suture from the anvil (pending anvil not long off the suture has been cut.)

WARNING: The delivery tube should be carefully inspected after removal to ensure the retaining suture has been removed with the delivery tube.

A. Insure the suture end of the anvil assembly to the EEA™ XL stapler. Since attached, the anvil head will be +0° to allow for smooth insertion.

WARNING: When mating the anvil with the stapler, ensure that the anvil center rod is aligned with the integrated trocar of the stapler.

1. Lubricate the anvil shaft with a water soluble lubricant and assembly. Driving a cut-end the longitudinal and make it difficult to attach beneath the anvil to the integrated trocar of the stapler. Allow advancing the anvil to the instrument, hold the black twist knob on the handle of the instrument firmly to prevent the integrated trocar from moving back slightly into the head of the device. This is especially important when using an EEA™ 29 mm.
2. Remove the anvil delivery tube with the anvil retention suture from the head of the device. This is especially important when using an EEA™ 29 mm.
3. Fit the anvil delivery tube over the anvil center rod and secure by pushing the suture well. Pull the delivery tube off until it is fully seated. Ensure the fitting is fully seated and secure in the tube prior to use.

WARNING: Improper cutting and removal of the advancing proximal guide suture can result in excess suture left in the cavity. The Anvil head and Sutures not should be inspected after removal to ensure the advancing proximal guide suture has been removed with the Suture not.

Use of the white accessory fitting (on a shortened delivery tube):

1. Cut the shortened length of the delivery tube to desired, the tube may be cut to the appropriate length.
2. Secure the white accessory fitting at the cut end of the delivery tube by pushing the barbed section of the fitting into the tube until it is fully seated. Ensure the fitting is fully seated and secure in the tube prior to use.

WARNING: A shortened delivery tube is not for use in transoral applications.

14. Follow instructions 2 through 4 above.

STORE AT ROOM TEMPERATURE. AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

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