IT

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EN

MULTIPOLAR EXTENSION CABLES - STERILE
INSTRUCTION FOR USE

PRODUCT DESCRIPTION

F5406/4REDS - F5406/6REDS - F5406/10REDS are connection cables for FIAB endocavitary SPIKE REDEL leads (or leads with REDEL connector, series P with female contacts) placed in the heart chambers and a pacing or recording device away from the surgical field.

The cables are constituted by a coaxial polyurethane shielded cable of 150 cm length ending in a pentagonal ABS box on one side.

From the box come out:

- 4 cables for **F5406/4REDS** model;
- 6 cables for **F5406/6REDS** model;
- 10 cables for **F5406/10REDS** model;

The cables are made of copper insulated with polyurethane ending with 2 mm touch-proof type gold plated plugs (complying with E-DIN 42802-2), numebered according to the model.

On the other end a male circular REDEL connector is attached that is made according to the model:

- 6 pin which is compatible with all SPIKE REDEL models up to 4 poles for **F5406/4REDS** model;
- 6 pin which is compatible with all SPIKE REDEL models up to 6 poles for **F5406/6REDS** model;
- 10 pin which is compatible with all decapolar SPIKE REDEL models for **F5406/10REDS** model;

SPIKE REDEL are the following models: SPIKE R, SPIKE ULTRA R, SPIKE LC ULTRA R.

PACKING CONTENTS

The models are sterilized and packaged in single paper-PE pouches.

STERILISATION

The package of the product for sale is STERILE.

NB: Sterility is guaranteed until the expiry date indicated on the package if it has been correctly stored and it is not damaged.

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
The product can be re-sterilized up to 20 times with ethylene oxide (EO) or low-temperature hydrogen peroxide gas plasma process ("STERRAD" / "STERIS"), by always taking care however to check for any alterations of the components and the overall functioning of the product.. Do not autoclave; do not use other heat sterilization methods.
The sterilization process must be compliant with the regulations in force.

- The validation of the process for the following sterilizations will be under user's liability.
- It's suggested to keep track of the number of sterilisation processes carried out on the device.
- For sterilization with ethylene oxide the cleaned product must be packaged in a suitable gas-permeable pouch. It is recommended to allow a proper aereation time in order to disperse the EO residues.
- To clean the product use detergent solutions (non ionic, surface-active) or disinfectant solutions (0.1 – 0.2 % quaternary ammonium base components). Follow the relative instructions for using these products. After cleaning remove the detergent by rubbing the lead several times with gauze soaked in distilled water. Dry well.

INDICATIONS
Connection to polygraphs or external pacemakers (with E-DIN 42802-2 compatible 2mm sockets) to SPIKE REDEL models leads placed into the heart.
CONTRAINDICATIONS
No contraindications exist as long as the product is used according to the instructions for use issued by the manufacturer.
PRECAUTIONS AND WARNINGS

1. **F5406/10REDS model can be used ONLY with SPIKE REDEL DECAPOLAR models and viceversa.**
2. **F5406/6REDS is ONLY compatible with SPIKE REDEL models with 6 poles or less.**
3. **F5406/4REDS are compatible with SPIKE REDEL models with 4 poles or less. IF USED WITH LEADS NOT MANUFACTURED BY FIAB, BEFORE USE, CAREFULLY CHECK THE COMPATIBILITY OF THE DEVICES. The REDEL connector of the lead to be used must be a P series one with 6 female contacts and connections to the electrodes as depicted in img. a (with number 1 connected to the distal electrode).**

img. a

4. Check packing before use. If wet, broken or damaged, sterility can't be guaranteed.
5. Sterility is guaranteed up to the expiry date printed on the package. If the latter has been correctly stored, it will not be damaged. In case of damage or visible defects, do not utilise the product but return it to **FIAB**.
6. The devices connected to the cables must comply with regulations in force.
7. N.B. use only with CF type devices.
8. Before connecting, review the instructions for use, supplied with the device which will be used
9. All electrical devices for patient support must be regularly grounded through equipotential node.
10. The devices must be connected and operated only by qualified personnel.
11. Before use check the integrity of the insulation of the cable, of the plastic parts and connection. Don't use the cable in presence of damages.
12. Insert the touch-proof plugs by pushing beyong the first resistance encountered to obtain the complet mating of the pin with the metallic insert of the socket.
13. Before the procedure check the effective insertion of every pin and connection by gently tugging them, to feel any movement.
14. Pay attention to the correct sequence of number on the plugs during insertion.
15. The connection between the lead and the cable must be the last in the connecting sequence.
16. Before the beginning of the procedure check the compatibility of the lead (with touch-proof plugs) with the devices it will be connected to. Shouldn't the sockets be compatible with touch proof type plugs, FIAB can provide suitable adapters FIAB PG922/2TP**X** (**X** = number of adapters in the package). ***N.B. at the disconnection of the lead from the device, make sure that the adapter (if used) remains fixed in the sockets in order to avoid the risk of accidental shocks transmitted to the cardiac tissue through the bare pin.***
17. Work in sterile environment.
18. Do not wash the leads with a solvent; if necessary, during the procedure, use a sterile physiologic solution. Dry well. Even during cleaning and drying it is necessary to stay within the temperature limits indicated on the label.
19. Make sure a defibrillator is ready for use

MODE OF USE
First, insert the plugs in the related sockets of the external device for recording or stimulation.
If the device sockets aren't touch-proof type, use adapters like FIAB PG922/2TP**X** (**X** = number of adapters in the package) following the enclosed intrucions for use. ***N.B. at the disconnection of the lead from the device, make sure that the adapters (if used) remains fixed in the sockets in order to avoid the risk of accidental shocks transmitted to the cardiac tissue through the bare pin.***
Insert the REDEL connector into the mating on the lead until they match with a click sound.
POTENTIAL COMPLICATIONS
None.
STORAGE
The product must be stored in the original package at the environmental conditions (temperature and relative humidity) specified on the pouch's label. Putting external heavy weights on the package, the product could be damaged.
WASTE DISPOSAL
Refuses deriving from health structures must be disposed in according to the regulation in force.
GENERAL NOTES
If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

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EN

F7818/4BIOLC
MULTIPOLAR EXTENSION CABLES
INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION
Model **F7818/4BIOLC** is a quadripolar cable. Conductive wire is copper and insulation material is polyurethane. At one end of the cable there are two couples of red (positive pole) and black (negative pole) screw sockets. Provided with a jaw mechanism which allows to modify the internal diameter from 0.6 to 2.0 mm by rotating the sockets same (clockwise or tightening and anti-clockwise for widening).
At the other end there is a round Redel® connector for the connection to external stimulators Biotronik Model Reocor D.
Cable is intended specially for the dual chamber pacing and has two grey colored leads for atrial chamber and two blue colored leads for ventricular chamber.
The total length of the cable is 200 cm.

CONTENTS OF PACKAGE
The cable is packaged separately in a double envelope, PE inside, thermo welded PE-paper outside.

STERILIZATION
The packet of the product for sale is STERILE and disposable, as shown on the labelling.

INDICATIONS
Cable F7818/4BIOLC is intended for the transmission of cardiac signal and pacing pulses for cardiac therapy and diagnosis
Cable F7818/4BIOLC is intended for the connection to external stimulators Biotronik Model Reocor D

CONTRAINDICATIONS
There are no contraindications as long as the product is used as indicated by the manufacturer.

PRECAUTIONS AND WARNINGS
- Check the package before using the product. If wet or torne, its sterility is not guaranteed.
- The devices that are connected to the cables must comply with the law in force.
- ***Only use devices that are suitable for direct cardiac application (like CF).***
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- Plug the round connector into the BIOTRONIK input device.
- Attention: use one color for each chamber (blue or grey).
- Before each procedure, carefully check that the plugs have been properly inserted and that all connections are solid, pulling them gently to check eventual movements.
- Work in a sterile environment.
- Make sure that a defibrillator is ready for use.
- The product life is reported on the label.

INSTRUCTIONS FOR USE
- Connect the lead to the stimulator as indicated in the previous paragraph.
- These cables are intended specially for the dual chamber pacing and has two grey colored leads for atrial chamber and two blue colored leads for ventricular chemaber.
- Connect the terminal corresponding to the CAUDAL epimyocardiac electrode to the BLACK socket. Firmly lock the contact by rotating the socket clockwise.
- Connect the terminal corresponding to the CRANIAL epimyocardiac electrode to the RED socket. Firmly lock the contact by rotating the socket clockwise.
- Start the pacing procedure.

POTENTIAL COMPLICATIONS
A partial insertion of the plugs into the sockets can cause an intermittent or ineffective stimulation and can lead to a lack of protection towards dispersion currents, with the following risk to prime dangerous ventricular arrhythmia

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SHELF LIFE
4 years if preserved in the original sterile packaging.
WASTE DISPOSAL
Refuses deriving from health structures must be disposed in according to the regulation in force.
GENERAL NOTES
If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

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EN INSTRUCTION FOR USE F7840/MEDT e F7840/4/MEDT

CONNECTION CABLES FOR MEDTRONIC PROGRAMMER 2090

PRODUCTS DESCRIPTION

- Interface cable that allows the Medtronic Programmer 2090 to be used with sterile threshold cables terminating with touch proof connectors.
- Constituted by a polyurethane shielded cable of 150cm (F7840/MEDT) and 400cm (F7840/4/MEDT), ending in a pentagonal ABS box on one side. On the box 4 touh-proof compatible sockets (2 red, 2 black) are placed (complying with EN 60601-1). On the other end there is a circular connector for the connection to MEDTRONIC PROGRAMMER 2090.
- The cables is not sterile.

PACKING CONTENTS

All the cables are packaged individually in pouch.

STERILISATION

The package of the product for sale is NOT STERILE. **WARNING: Sterilize before use**.

The product can be re-sterilized up to 20 times with ethylene oxide (EO) or low-temperature hydrogen peroxide gas plasma ("STERRAD" / "STERIS"), by always taking care however to check for any alterations of the components and the overall functioning of the product, in particular with regard to the rotation of the sockets. Do not autoclave; do not use other heat sterilization methods.

The sterilization process must be compliant with the regulations in force.

- The validation of the process for the following sterilizations will be under user's liability.
- It's suggested to keep track of the number of sterilisation processes carried out on the device.
- For sterilization with ethylene oxide the cleaned product must be packaged in a suitable gas-permeable pouch. It is recommended to allow a proper aereation time in order to disperse the EO residues. To clean the product use

INDICATIONS

To be used when performing a threshold test on perminant pacing and defibrilation leads at time of implant.

CONTRAINDICATIONS

No contraindications exist as long as the product is used according to the instructions for use issued by the manufacturer.

PRECAUTIONS AND WARNINGS

ATTENTION: F7840/4/MEDT model has cable with length higher than the standard. The cable may be subject to electromagnetic disturbances that, when particularly intense, could affect the operation of the device to which it is connected. If necessary remove noise sources from the cable.

1. **Used ONLY with MEDTRONIC PROGRAMMER 2090.**
2. Check packing before use.
3. External cardiac stimulators and intravenous potentials recording devices must comply with regulations in force.
4. ***N.B. use only with CF type devices.***
5. Before connecting, review the instructions for use, supplied with the device which will be used
6. All electrical devices for patient support must be regularly grounded through equipotential node.
7. The devices must be connected and operated only by qualified personnel.
8. Before use check the integrity of the insulation of the cable, of the plastic parts and connection. Don't use the cable in presence of damages.
9. Insert the touch-proof pins of the threshold cable by pushing beyond the first resistance encountered to obtain the complete mating of the pin with the metallic insert of the socket
10. Before the procedure check the effective insertion of every pin and connection by gently tugging them, to feel any movement.
11. The connection between the device and the cable must be the first in the connecting sequence.
12. Before the beginning of the procedure check the compatibility of the lead (with touch-proof pins) with the devices it will be connected to.
13. Work in sterile environment.
14. Do not wash the cable with a solvent; if necessary, during the procedure, use a sterile physiologic solution. Dry well. Even during cleaning and drying it is necessary to stay within the temperature limits indicated on the label.
15. Make sure a defibrillator is ready for use
16. The shelf-life of the product is indicated on the label.

MODE OF USE

Plug the circular connector into the threshold cable input of the Medtronic Programmer 2090 (It will only fit the correct socket). Unpack the sterile threshold cable and remove the touch-proof connector terminals from the sterile field. Plug the touch-proof connectors into the pentagonal ABS box respecting the polarity. Operate according to the instruction of the Medtronic programmer.

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EN

INSTRUCTIONS FOR USE

The instructions for use report to the following models:

- **Sterile Kit FIAB models PG922/2TPSx (x=; number of adapters);**
- **Non-sterile Kit FIAB models PG922/2TPxWG (x=; number of adapters);**
- **Sterile Kit FIAB models PG922/2TPVSx (x=; number of adapters);**
- **Non-sterile Kit FIAB models PG922/2TPVx (x=; colour of adapters);**

DESCRIPTION OF THE PRODUCT

The Kit were designed to make possible the connection of leads with "touch-proof" type connection pins (complying with EN 60601-1 standard) to 2 mm sockets of medical devices (cardiac potentials recorders, external pacemakers both intravenous or transesophageal ones) manufactured before the enforcement of the above mentioned standard. The adapters are constituted of an external shell made in Polypropylene in black, red and grey, with an inner gold plated brass contact shaped in order to get a socket on one end (complying with E-DIN 42802-2 standard) for the connection to the lead and a bare 2 mm male pin to be fixed on the device on the other end. PG922/2TPVSx and PG922/2TPVx models have a screw to secure the adapter to the plug of the lead.

PACKAGING

Sterile Kit FIAB, PE pouch + Paper-PE pouch. In 2 piece pouches one is black and the other one is red and they will be fixed to the corresponding black and red sockets on the equipment to be upgraded. In 4 piece pouches, except black and red, the other 2 adapters are grey an the will be fixed to the other sockets of the device.

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Non-Sterile Kit FIAB, PE ziplock pouch, contains 2, 3, 4, 6 or 10 adapters (depending on the model). In 2 piece pouches one is black and the other one is red and they will be fixed to the corresponding black and red sockets on the equipment to be upgraded. In other models red and black adapters are still present; the packaging is completed by grey adapters, which can be fixed to the other sockets of the device.
Non-Sterile Kit FIAB Mod. PG922/2TPVx: every package, PE ziplock pouch, contains adapters red and/or black, depending on the model.

STERILISATION
Sterile Kit FIAB are sterilised with EO gas. Do not reuse. Reusing the product may lead to: possible presence of biological residues that may cause cross-infections; alteration of materials; loss of initial functional features of product.
Other models are supplied **Non-Sterile**.

INDICATIONS
Assure the compatibility of "touch-proof" connection pins (E-DIN 42802-2) with not protected 2mm sockets in emergency case. At the end of emergency the adapters can be also fixed to the equipment by means of cyanacrylic adhesive supplied by FIAB separately.
Attention: the adapter must be connected at first in the corresponding sockets on the equipment to be upgraded. In this way the PLUG (pin conductive part of lead) will be always protected from accidental contacts.
N.B. Inter-axes between sockets must be at least 10 mm.

CONTRAINDICATIONS
No contraindications are indicated if the product is used according to the intended destination.

PRECAUTIONS AND WARNINGS
- Please, do open the package only after a careful check of the integrity of the pouches: for the Sterile Kit FIAB sterilization would be compromised. It is guaranteed until the expiry date indicated on the package if it has been correctly stored and it is not damaged.
- During the use, make sure that the adapter stays fixed to the equipment socket.
- If necessary (for example if a new positioning is needed), extract the plug so that the adapter remains fixed on the socket.

EMERGENCY MODE OF USE
- Open the package.
- Connect the adapters on the electronic device following the right polarity.
- For PG922/2TPVSx and PG922/2TPVx models, tighten the screw to secure the adapter to the plug of the lead.
- Insert the touch-proof pin of the catheter or extension cable only onto the adapters already mounted on the device, following the right polarity.
- If necessary (for example if a new positioning is needed) and however once the procedure is over, extract the PLUG (pin conductive part of lead) so that the adapter remains fixed on the socket. In this way the plug will be always protected from accidental contacts

PRECAUTIONS IN CASE OF GLUEING
- **FOR THE FIXING PROCEDURE OF THE ADAPTERS USE CYANACRYLATE ADHESIVE DISPENSER PROVIDED BY FIAB**
- **THE FIXING PROCEDURE OF THE ADAPTERS ON THE ELECTRONIC DEVICE MUST BE CARRIED OUT BY THE COMPETENT TECHNICAL DEPARTMENT OF THE HOSPITAL OR OTHER ADEQUATELY QUALIFIED PERSONNEL**
- Carry on the procedure of fixing the adapters at least 1 hour before the connection to the lead.
- Apply the adhesive carefully. Before its application check the dispenser to assure that the aluminium seal has been removed and firmly secure the plastic cap.
- During the use of the adhesive be careful because skin and eyes can be glued in few seconds.
- Store away from children's reach.
- Don breathe the vapours.
- Carry on the procedure in a well aired environment.
- Don't cover the pins with the adhesive because it could lead to the loss of electrical contact.

METHOD OF GLUEING
1. Insert completely the adapter pin in the socket which will be upgraded. **FIGURE 1** and **FIGURE 2.**
2. Extract for 2-3 mm the pin from the socket. **FIGURE 3.**
3. Put a drop of cyanoacrylate adhesive in the space between the socket and the adapters body. **FIGURE 4**
4. **N.B.** The extraction of the pin exceeding the indicated length and an excessive quantity of adhesive could lead to an alteration of the electrical conductivity because of the presence of adhesive inside the socket.
 In case of a spreading of the adhesive on a vast surface of the metallic pin, DISCARD the adapter and use another one. Check that the glue hasn't blocked the hole of the socket or dripped to deeply in the electrical contact **FIGURE 6.** If so, don't use the device and contact Servicing.
5. Reintroduce the adapter pin completely in the socket, holding the adapter pressed for 1-2 minutes. **FIGURE 5.**
6. Repeat the procedure on the other sockets.
7. Check if the adapters are firmly fixed to the device, otherwise repeat the procedure from point 1.

STORAGE
The product must be stored in the original package at the environmental conditions (temperature and relative humidity) specified on the pouch's label. Putting external heavy weights on the package, the product could be damaged.

SHELF LIFE
FIAB Sterile Kit, in their package, expire 4 years after manufacturing date.
FIAB Non-Sterile Kit in its package, expires 5 years after manufacturing date.

GENERAL NOTES
If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

WASTE DISPOSAL
Refuses deriving from health structures must be disposed in according to the regulation in force.

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EN

F7816xx SERIES
REUSABLE BIPOLAR EXTENSION CABLE

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

F7816xx models, supplied both in a sterile (*F7816S/xx*) and non-sterile (*F7816/xx*) packaging, are bipolar copper leads, sheated with soft plastic covering. At one end, the cable terminates with an box that supports two bushes where the polarity is indicated. The bushes incorporate a mechanism that, by means of a rotating element, varies the internal profile from 0 to 2 mm (clockwise rotation to tighten or anticlockwise to enlarge).

At the other end, depending on the model, there are:

- for "**30**" suffix : two shrouded 2 mm plugs Touch Proof, black and red, model complying with EN60601-1, E-DIN 42802-2;
- for "**33**" suffix : two 2 mm plugs, back and red, with protection cap;
- for "**MED**" suffix : a round male connector for connection to PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 and Mod. 5348;
- for "**RED**" suffix : a circular connector for connection to devices Redel® PACEMAKER PACE MEDICAL Models n. 4170 and 4580;
- for "**BIO**" suffix : a circular connector for connection to devices Redel® PACEMAKER BIOTRONIK Mod. Reocor S;

The total length of the cable is 200 mm.

CONTENTS OF PACKAGE

The cable is singly packed in closed PE pouch (zip-lock type).

STERILIZATION

F7816/xx series is supplied NON-STERILE. **WARNING: Sterilize before use**.

F7816S/xx series is supplied already sterilized with Ethylene Oxide (EO) for the first use.

The product may be sterilised up to 25 times:

- ethylene oxide (EO);
- steam autoclave;

NB: the sterilization process used must be compliant with the regulations in force.

The validation of the process for the following sterilizations will be the user's responsibility. You are advised to keep a tally of the number of sterilizations carried out.

After each sterilization check for any alterations of the components and the overall functioning of the product.

EO sterilization

For sterilization with ethylene oxide the clean product must be packaged in a suitable gas-proof packet. It is recommended that a sufficient aeration time be allowed in order for the EO residues to be completely dispersed.

Steam autoclave sterilization

Maximum sterilization times for steam autoclave at the temperature indicated are as follows:

Cycle	Exposure	Drying
134°C (273°F)	134°C (273°F); 2,1÷ 2,2 bar for 11 min	17 min
121°C (250°F)	121°C (250°F); 1,1÷1,2 bar for 21 min	17 min

ATTENTION: Should the products come in contact with the metal supports within the autoclave, their life may be shortened. It is therefore necessary the use of appropriate means spacing or wrapping as gauze.

CLEANING

- Before sterilization, clean the cable to remove all trace of blood or other body fluids. Caution: during the cleaning procedure, don't apply excessive stress that could affect the cable's integrity.
- For cleaning, use detergent solutions (non ionic, surface-active) or disinfectant solutions (0.1 – 0.2 % quaternary ammonium base components). Follow the relative instructions for using these products. After cleaning remove the detergent by rubbing the lead several times with gauze soaked in distilled water. Dry well.
- The use of mechanical washing systems may reduce the product's life, therefore such system shall be previously validated. The following method has been validated by the Manufacturer:
 1. Cold rinse with softened water – No detergent – Temperature less than 35°C.
 2. Main detergent wash with ph balanced detergent – Hot water typically 55°C.
 3. Reverse Osmosis rinse with softened water to remove any residual chemicals – typically 60°C.
 4. Main Reverse Osmosis disinfection – up to 93°C for approximately 1 minute (temperature rises up from 60 to 93°C).
 5. Drying –120°C for approximately 15 minutes prior to sterilization.
- Do not use ultrasound baths.
- Do not use abrasive mechanical systems to remove biological residue: the product could be damaged.

INDICATIONS

F7816xx cables are intended to connect an endocavitary bipolar temporary electrode with an external cardiac stimulator:

- "**30**" suffix in case of TOUCH-PROOF connections, complying with EN60601-1 and E-DIN 42802-2 regulations;
- "**33**" suffix in case of female 2mm connections;
- "**MED**" suffix in case of TOUCH-PROOF connections compatibility with PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 e Mod. 5348;
- "**RED**" suffix in case of connections compatible with Redel® PACEMAKER PACE MEDICAL Mod. 4170 and Mod. 4580;
- "**BIO**" suffix in case of connections compatible with Redel® PACEMAKER BIOTRONIK Mod. Reocor S;

CONTRAINDICATIONS

There are no contraindications as long as the product is used as indicated by the manufacturer.

WARNING: Always check the compatibility with the connections of the electrostimulator.

WARNING: "**MED**" suffix models can be used only with **MEDTRONIC PACEMAKER Mod. 5392 and Mod. 5388 (bicameral) - 5348 (momocameral)**.

WARNING: "**30**" suffix can be used only with devices fitted with sockets complying with **EN60601-1 and E-DIN 42802-2 regulations.**

WARNING: when using the product, follow all the precautions and warnings reported below, especially those concerning models with "**33**" suffix since the connections that are far from the patient are uncovered once the protection caps are removed.

WARNING: "**RED**" suffix can be used only with **PACEMAKER PACE MEDICAL Mod. 4580 (bicameral) and 4170 (unicameral)**.


WARNING: "**BIO**" suffix can be used only with **PACEMAKER BIOTRONIK Mod. Reocor S**.

PRECAUTIONS AND WARNINGS

- ***Sterilise before use.***
- The devices that are connected to the cables must comply with the law in force.
- ***Only use devices that are suitable for direct cardiac application (like CF).***
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- Make sure that the body of the lead is intact, that there are no incrustations on the jack-plugs.
- The cable must be connected to the equipment as the first connection of the system. The extension cable must be connected to the electrode as the last connection of the system.

For models with "**33**" suffix

After removing the protection caps, the accidental contact with the metal part of the plugs (free ends far from the patient) may cause dangerous arrhythmias due to induced currents transmitted to the cardiac tissue through the uncovered pin; therefore, the cable must be connected to the epimyocardial electrode after connecting the cable to the equipment; conversely, the cable must be disconnected from the epimyocardial electrode before disconnecting the cable from the equipment (see figure).

For models with "**30**" suffix

Insert the touch-proof plugs to the very end, going further than the first resistance that you encounter, so as to obtain the complete coupling between the plug and the metal insert of the bush.

For models with "**MED**" suffix

Insert the round connector into the input of the Medtronic device.

For models with "**RED**" suffix

Plug the round connector into the Pace Medical input device.

For models with "**BIO**" suffix

Plug the round connector into the Biotronik input device.

- Before each procedure, carefully check that the plugs have been properly inserted and that all connections are solid, pulling them gently to check eventual movements.
- Work in a sterile environment.
- Make sure that a defibrillator is ready for use.
- The product life is reported on the label.

INSTRUCTIONS FOR USE

- Connect the lead to the stimulator as indicated in the previous paragraph.
- Introduce the most CAUDAL electrode terminal in the bush marked "-" (minus) on the box. Firmly tighten the contact by turning the roller clockwise.
- Introduce the most CRANIAL electrode terminal in the bush marked "+" (plus) on the box. Firmly tighten the contact by turning the bush clockwise.
- Start the electrostimulation.

POTENTIAL COMPLICATIONS

An incomplete blockage of the electrodes to the cable can cause an intermittent or ineffective stimulation and can provide incomplete protection against dispersion currents, with the ensuing risk of dangerous ventricular arrhythmias.

STORAGE

The product must be stored in the original package at the environmental conditions (temperature and relative humidity) specified on the pouch's label. Putting external heavy weights on the package, the product could be damaged.

SHELF LIFE

In case of sterile products, check the expiry date and the integrity of the package before the first use.

The life of the product is a consequence of the number of uses and the state of maintenance and conservation of the product itself. Products are reusable and can be EO or steam re-sterilized up to 25 times.

WASTE DISPOSAL

Refuses deriving from health structures must be disposed in according to the regulation in force.

GENERAL NOTES

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

The brands and commercial names in this document are property of the corresponding manufacturers.

FIAB SpA
President of the Board
Alberto Calabrò

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EN

F7818xx SERIES
BIPOLAR EXTENSION CABLES
REUSABLE BIPOLAR EXTENSION CABLES
MULTIPOLAR EXTENSION CABLES

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

F7818xx models are bipolar copper leads, sheated with soft plastic covering. At one end the lead ends with two special sockets, black(negative pole) and red (positive pole), provided with a jaw mechanism which allows to modify the internal diameter from 0.6 to 2.0 mm by rotating the sockets same (clockwise or tightening and anti-clockwise for widening).

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- for "**30**" suffix : two male 2 mm plugs, black and red, TOUCH-PROOF model complying with EN60601-1, E-DIN 42802-2.
- for "**33**" suffix : two male 2 mm plugs, black and red, with protection cap (standard plugs).
- for "**MED**" suffix : a multipolar, round, male connector for connection to PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 and Mod. 5348.
- for "**R**" suffix : a circular connector for connection to devices Redel® PACEMAKER PACE MEDICAL Models n. 4170 and 4580;
- for "**BIO**" suffix : a circular connector for connection to devices Redel® PACEMAKER BIOTRONIK Mod. Reocor S;
- for "**OSY**" suffix : two 2mm plugs with bushing and protective cap for OSYPKA devices with screw sockets.

The total length of the cable is 200cm.

CONTENTS OF PACKAGE

The cable is packaged separately in a double envelope, PE inside, thermo welded PE-paper outside.

STERILIZATION

The packet of the product for sale is STERILE.

"LC" models are disposable and sterilised with EO gas (as shown on the labelling). Do not reuse. Reusing the product may lead to: possible presence of biological residues that may cause cross-infections; alteration of materials; loss of initial functional features of product.

Other models can be re-sterilized up to as many as 20 times in ethylene oxide (EO) or low-temperature hydrogen peroxide gas plasma ("STERRAD" / "STERIS"), by always taking care however to check for any alterations of the components and the overall functioning of the product, in particular with regard to the rotation of the sockets. Do not autoclave; do not use other heat sterilization methods.

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The validation of the process for the following sterilizations will be the user's responsibility. You are advised to keep a tally of the number of sterilizations carried out. For sterilization with ethylene oxide the clean product must be packaged in a suitable gas-proof packet. It is recommended that a sufficient aeration time be allowed in order for the EO residues to be completely dispersed.

INDICATIONS

F7818xx cables are intended to connect an endocavitary bipolar temporary electrode or an epimyocardiac electrode with an external cardiac stimulator. Each model should be used for a specific connection:

- "**30**" suffix in case of TOUCH-PROOF connections, complying with EN60601-1 and E-DIN 42802-2 regulations;
- "**33**" suffix in case of female 2mm connections;
- "**MED**" suffix in case of TOUCH-PROOF connections compatibili con PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 e Mod. 5348;
- "**R**" suffix in case of connections compatible with Redel® PACEMAKER PACE MEDICAL Mod. 4170 and Mod. 4580;
- "**BIO**" suffix in case of connections compatible with Redel® PACEMAKER BIOTRONIK Mod. Reocor S;
- "**OSY**" suffix in case of OSYPKA devices with screw sockets;

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WARNINGS


- Always check the compatibility with the bushes of the electrostimulator.
- "**MED**" suffix models can be used only with **MEDTRONIC PACEMAKER Mod. 5392 and Mod. 5388 (bicameral) or 5348 (momocameral)**.
- "**30**" suffix models can be used only with devices fitted with sockets complying with **EN60601-1 and E-DIN 42802-2 regulations**.
- When using the product, follow all the precautions and warnings reported below, especially those concerning models with "**33**" suffix since the connections that are far from the patient are uncovered once the protection caps are removed.
- "**R**" suffix models can be used only with **PACEMAKER PACE MEDICAL Mod. 4580 (bicameral) and 4170 (unicameral)**.
- "**BIO**" suffix models can be used only with **PACEMAKER BIOTRONIK Mod. Reocor S**.
- "**OSY**" suffix models can be used only with **OSYPKA devices with screw sockets**.

PRECAUTIONS AND WARNINGS

- Check the package before using the product. If wet or torne, its sterility is not guaranteed.
- The devices that are connected to the cables must comply with the law in force.
- *Only use devices that are suitabe for direct cardiac application (like CF).*
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- Before using the re-sterilized product (if possible) make sure that the body of the lead is intact, that there are no incrustations on the jack-plugs.
- For cleaning, use detergent solutions (non ionic, surface-active) or disinfectant solutions (0.1 – 0.2 % quaternary ammonium base components). Follow the relative instructions for using these products. After cleaning remove the detergent by rubbing the lead several times with gauze soaked in distilled water. Dry well. Even during cleaning and drying it is necessary to stay within the temperature limits indicated on the label.
- The cable must be connected to the equipment as the first connection of the system. The extension cable must be connected to the electrode as the last connection of the system.

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For models with "**R**" suffix

Insert the round connector into the Pace Medical input device.

For models with "**BIO**" suffix

Insert the round connector into the Biotronik input device.

For models with "**OSY**" suffix

Insert the plugs into the input of the Osypka device and screw to secure the connection.

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INSTRUCTIONS FOR USE

- Connect the lead to the stimulator as indicated in the previous paragraph.
- Connect the terminal corresponding to the CAUDAL epimyocardiac electrode to the BLACK socket. Firmly lock the contact by rotating the socket clockwise.
- Connect the terminal corresponding to the CRANIAL epimyocardiac electrode to the RED socket. Firmly lock the contact by rotating the socket clockwise.
- Start the pacing procedure.

POTENTIAL COMPLICATIONS

A partial insertion of the plugs into the sockets can cause an intermittent or ineffective stimulation and can lead to a lack of protection towards dispersion currents, with the following risk to prime dangerous ventricular arrhythmia

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SHELF LIFE

4 years if preserved in the original sterile packaging. "LC" models are disposable, other models are reusable and can be re-sterilized up to 20 times.

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EN

F7817xx SERIES
BIPOLAR EXTENSION CABLES
REUSABLE BIPOLAR EXTENSION CABLES

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

F7817xx models are bipolar copper leads, sheated with soft plastic covering. At one end the lead ends with two 2mm sockets, TOUCH-PROOF model, complying with EN60601-1 and E-DIN 42802-2 regulations, black (negative pole) and red (positive pole).

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- for "**30**" suffix : two male 2 mm plugs, black and red, TOUCH-PROOF model complying with EN60601-1, E-DIN 42802-2.
- for "**33**" suffix : two male 2 mm plugs, black and red, with protection cap (for suffix "**33**" is fitted with standard plugs, while for suffix "**33E**" is fitted with spring loaded plugs).
- for "**BIO**" suffix : a circular connector for connection to devices Redel® PACEMAKER BIOTRONIK Mod. Reocor S.
- for "**MED**" suffix : a multipolar, round, male connector for connection to devices PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 and Mod. 5348.
- for "**OSY**" suffix: two 2mm plugs with bushing and protective cap for OSYPKA devices with screw sockets.

The total length of the cable is 200cm.

CONTENTS OF PACKAGE

The cable is packaged separately in a double envelope, PE inside, thermowelded PE-paper outside.

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INDICATIONS

F7817xx cables are intended to connect an endocavitary bipolar temporary electrode (or a threshold cable) with an external cardiac stimulator. Each model should be used for a specific connection:

- "**30**" suffix in case of TOUCH-PROOF connections complying with EN60601-1 and E-DIN 42802-2 regulations;
- "**33**" suffix in case of female 2mm connections;
- "**BIO**" suffix in case of Redel® connections compatible with PACEMAKER BIOTRONIK Mod. Reocor S;
- "**MED**" suffix in case of TOUCH-PROOF connections compatible with PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 e Mod. 5348;
- "**OSY**" suffix in case of OSYPKA devices with screw sockets;

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WARNINGS


- Always check the compatibility with the bushes of the electrostimulator.
- "**MED**" suffix models can be used only with **MEDTRONIC PACEMAKER Mod. 5392 and Mod. 5388 (bicameral) or 5348 (momocameral)**.
- "**30**" suffix models can be used only with devices fitted with sockets complying with **EN60601-1 and E-DIN 42802-2 regulations**.
- When using the product, follow all the precautions and warnings reported below, especially those concerning models with "**33**" suffix since the connections that are far from the patient are uncovered once the protection caps are removed.
- "**BIO**" suffix models can be used only with **PACEMAKER BIOTRONIK Mod. Reocor S**.
- "**OSY**" suffix models can be used only with **OSYPKA devices with screw sockets**.

PRECAUTIONS AND WARNINGS

- Check the package before using the product. If wet or torne, its sterility is not guaranteed.
- The devices that are connected to the cables must comply with the law in force.
- ***Only use devices that are suitabe for direct cardiac application (like CF).***
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- Before using the re-sterilized product (if possible) make sure that the body of the lead is intact, that there are no incrustations on the jack-plugs.
- For cleaning, use detergent solutions (non ionic, surface-active) or disinfectant solutions (0.1 – 0.2 % quaternary ammonium base components). Follow the relative instructions for using these products. After cleaning remove the detergent by rubbing the lead several times with gauze soaked in distilled water. Dry well. Even during cleaning and drying it is necessary to stay within the temperature limits indicated on the label.
- The cable must be connected to the equipment as the first connection of the system. The extension cable must be connected to the electrode as the last connection of the system.

For models with "**33**" suffix

After removing the protection caps, the accidental contact with the metal part of the plugs (free ends far from the patient) may cause dangerous arrhythmias due to induced currents transmitted to the cardiac tissue through the uncovered pin; therefore, the cable must be connected to the electrocatheter after connecting the cable to the equipment; conversely, the cable must be disconnected from the electrocatheter before disconnecting the cable from the equipment (see figure).

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For models with "**BIO**" suffix

Insert the round connector into the input of the Biotronik.

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INSTRUCTIONS FOR USE

- Connect the lead to the stimulator as indicated in the previous paragraph.
- Place the electro-catheter.
- Introduce the plug relevant to the distal pole of electrocatheter into the black socket of the leads
- Introduce the plug relevant to the proximal pole of electrocatheter into the red socket of the leads
- Start electrostimulation procedure.

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EN

F7819xx SERIES
BIPOLAR EXTENSION CABLES
INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The models of the **F7819xx** series are copper bipolar cables covered in a flexible plastic material. At one end, the cable terminates with an ABS box that supports two bushes where the polarity is indicated. The bushes incorporate a mechanism that, by means of a rotating element, varies the internal profile from 0 to 2 mm (clockwise rotation to tighten or anticlockwise to enlarge). At the other end, depending on the model, there are:

- for "**30**" suffix : two shrouded 2 mm plugs Touch Proof, black and red, model complying with EN60601-1, E-DIN 42802-2;
- for "**33**" suffix : two 2 mm plugs, back and red, with protection cap;
- for "**MED**" suffix : a multipolar, round, male connector for connection to PACEMAKER MEDTRONIC Mod. 5392, Mod. 5388 and Mod. 5348;
- for "**BIO**" suffix : a circular connector for connection to devices Redel® PACEMAKER BIOTRONIK Mod. Reocor S;
- for "**OSY**" suffix: two 2mm plugs with bushing and protective cap for OSYPKA devices with screw sockets.

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STERILIZATION

The product package is STERILE, SINGLE-USE and sterilised with EO gas (as shown on the labelling). Do not reuse. Reusing the product may lead to: possible presence of biological residues that may cause cross-infections; alteration of materials; loss of initial functional features of product.

INDICATIONS

The **F7819xx** series cables are used for the connection between an epimyocardial or bipolar endocavity temporary electrode and an external cardiac stimulator with:

- "**30**" suffix in case of TOUCH-PROOF connections complying with EN60601-1 and E-DIN 42802-2 regulations;
- "**33**" suffix in case of female 2mm connections;
- "**MED**" suffix in case of touch-proof connections compatible with PACEMAKER; MEDTRONIC Mod. 5392, Mod. 5388 and Mod. 5348;
- "**BIO**" suffix in case of connections compatible with PACEMAKER BIOTRONIK Mod. Reocor S;
- "**OSY**" suffix in case of OSYPKA devices with screw sockets;

WARNING: Always check the compatibility with the bushes of the electrostimulator.

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WARNINGS


- Always check the compatibility with the bushes of the electrostimulator.
- "**BIO**" suffix models can be used only with **PACEMAKER BIOTRONIK Mod. Reocor S**.
- "**MED**" suffix models can be used only with **MEDTRONIC PACEMAKER Mod. 5392 and Mod. 5388 (bicameral) or 5348 (momocameral)**.
- "**30**" suffix models can be used only with devices fitted with sockets complying with **EN60601-1 and E-DIN 42802-2 regulations**.
- When using the product, follow all the precautions and warnings reported below, especially those concerning models with "**33**" suffix since the connections that are far from the patient are uncovered once the protection caps are removed.
- "**OSY**" suffix models can be used only with **OSYPKA devices with screw sockets**.

PRECAUTIONS AND WARNINGS

- Check the package before using the product. If wet or torne, its sterility is not guaranteed.
- The devices that are connected to the cables must comply with the law in force.
- ***Only use devices that are suitabe for direct cardiac application (like CF).***
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- The cable must be connected to the equipment as the first connection of the system. The extension cable must be connected to the electrode as the last connection of the system.

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HOW TO USE

- Connect the cable to the external cardiac stimulator, according of the previous paragraph.
- Introduce the most CAUDAL electrode terminal in the bush marked "-" (minus) on the box. Firmly tighten the contact by turning the roller clockwise.
- Introduce the most CRANIAL electrode terminal in the bush marked "+" (plus) on the box. Firmly tighten the contact by turning the bush clockwise.
- Start the electrostimulation.

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PRODUCT LIFE

4 years from the production date shown on the label, which reports the traceability data.

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