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
BARD Biopsy

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ENGLISH

DS4001 Vacuum and Rinse Tubing Cassette
DS4002 Vacuum Tubing Cassette

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EnCor Enspire™ Vacuum and Rinse Tubing Cassette provides vacuum and the option for a saline rinse to be delivered to the sample trap of an EnCor® Probe during a biopsy procedure.

The EnCor Enspire™ Vacuum Tubing Cassette provides vacuum to an EnCor® Probe during a biopsy procedure.

INDICATIONS FOR USE

The EnCor Enspire™ Breast Biopsy System is indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities.

- It is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.
- It is intended to provide breast tissue for histologic examination with partial removal of a palpable abnormality.

The extent of a histologic abnormality cannot always be readily determined from palpation or imaged appearance. Therefore, the extent of removal of the palpated or imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality, e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

In instances when a patient presents with a palpable abnormality that has been classified as benign through clinical and/or radiological criteria (e.g. fibroadenoma, fibrocystic lesion), the EnCor Enspire™ Breast Biopsy System may also be used to partially remove such palpable lesions. Whenever breast tissue is removed, histological evaluation of the tissue is the standard of care. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

CONTRAINDICATIONS

1. This device is not intended for use except as indicated.
2. The EnCor Enspire™ Breast Biopsy System is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

WARNINGS

1. **Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**
2. **Do not reuse. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with microorganisms which may lead to infectious complications.**
3. **Use of accessories not compatible with the EnCor Enspire™ Breast Biopsy System may create potentially hazardous conditions.**

PRECAUTIONS

1. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

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2. This accessory should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
3. Inspect tubing connections to the vacuum canister and the vacuum or vacuum and rinse tubing cassette to ensure proper vacuum levels are achieved and maintained during use.
4. Inspect the vacuum canister to ensure the lid is secure and that no damage has occurred during shipping or installation. A heavily scratched canister can break during use.

POTENTIAL COMPLICATIONS

1. Potential complications are those associated with any percutaneous removal/biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, a non-healing wound, pain and tissue adherence to the biopsy probe while removing it from the breast.

EQUIPEMENT REQUIRED

EnCor Enspire™ Breast Biopsy System
EnCor®, EnCor® MRI, EnCor® 360 Probe
EnCor®, EnCor® MRI, EnCor® 360 Driver
Vacuum Canister
Saline (when using optional saline rinse feature)

DIRECTIONS FOR USE

1. Open the doors on the EnCor Enspire™ Breast Biopsy System console.
2. Install vacuum canister in the console. Ensure correct orientation for proper fit. Verify the lid is securely snapped onto the vacuum canister. Securely attach the center port to the vacuum tubing. Verify all unused lid port caps are pressed firmly over the unused ports. (Figure 1)

Figure 1: Install the Vacuum Canister

3. Install the cassette into the console, and ensure the assembly is locked in place. The release lever on the cassette will click when installed properly. (Figure 2)

Figure 2: Install the Cassette

4. Connect the cassette tubing to the vacuum canister. Do not twist tubing during installation. Verify the connection is secure to prevent a decrease in vacuum. (Figure 3)


Figure 3: Connect Cassette Tubing to Vacuum Canister

5. For vacuum and rinse tubing cassette only, connect the saline spike to the saline. Place the saline in the saline cavity of the console. Verify the saline tubing is not kinked. (Figure 4A & 4B)

Figure 4A: Connect the Saline Spike

Figure 4B: Install the Saline in the Saline Cavity

6. When ready to begin the biopsy procedure, connect the probe vacuum tubing to the cassette vacuum tubing connector. If using rinse, connect the probe rinse tubing to the cassette rinse tubing connector. Verify the connections are secure to prevent leakage. (Figure 5)


Figure 5: Connect Probe Tubing to Cassette Connectors

7. At the conclusion of the procedure, disconnect all tubing, remove the cassette and dispose of properly.
 NOTE: After use, these products may be potential biohazards. Handle and dispose of in accordance with acceptable medical practice and applicable local, state and federal laws and regulations.
8. Refer to the EnCor Enspire™ Breast Biopsy System Instructions for Use, for information on cleaning and maintenance of the EnCor Enspire™ Breast Biopsy System.

HOW SUPPLIED

The vacuum and rinse tubing cassette and vacuum tubing cassette are supplied sterile and are intended for single use only.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

Extended warranty service agreements are available. Consult a Bard representative for details regarding terms and conditions.

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An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Manufacturer:
SenoRx, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

TEL: 1-480-894-9515
1-800-321-4254
FAX: 1-480-966-7062
1-800-440-5376
www.bardbiopsy.com

CE 0086

EC REP
Authorised Representative in the European Community
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

PK1237100 Rev. 2 03/12

Always confirm use of the current revision of this document. BAW1468600, 1

ENGLISH

EnCor® Directional Vacuum-Assisted Breast Biopsy Device and Driver

REF	
	EPC0112G EPC0116GV
	ECP0110G ECP017GV
	EPC0117G DRENCOR
	ECP0112GV

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The EnCor® Breast Biopsy Device is indicated for use in acquiring tissue for diagnosing breast abnormalities. The EnCor® Breast Biopsy Probe and EnCor® Breast Biopsy Driver are sold and packaged separately. The EnCor® Breast Biopsy Probe and EnCor® Breast Biopsy Driver are used together with one of the EnCor® Breast Biopsy Systems which include EnCor®, EnCor Essentra®, and EnCor Ultra™. These systems provide vacuum and rotation (EnCor®, EnCor Ultra™). Refer to the Indications for Use for the specific EnCor® Breast Biopsy System being used.

The EnCor® Breast Biopsy Probe contains a sharp trocar tip at the distal end for insertion and oscillating tubular cutter for tissue aspiration. A detachable tissue collection chamber is located on the proximal end of the device. The tissue collection chamber incorporates both a syringe luer for aspiration and transport of tissue samples and a saline rinse line for optional rinse of acquired samples in the chamber. A syringe adapter is provided to facilitate the connection of a syringe to the luer. The EnCor® Breast Tissue Markers which are sold and packaged separately may be inserted through the proximal end of the device.

The EnCor® Breast Biopsy Driver is designed as a handheld unit for ultrasound-guided breast biopsies and for use in a stereotactic platform using Bard® adapters. These adapters are sold and packaged separately.

Indications for Use

The EnCor® Breast Biopsy Probe and EnCor® Breast Biopsy Driver are indicated to acquire tissue for diagnostic sampling of breast abnormalities. The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined by the gross pathologic appearance. Therefore the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality, e.g. malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

Contraindications

This device is not intended for use except as indicated. In addition, it is not intended for use under conditions where tissue tension is contraindicated.

Warnings

- **Care must be taken when positioning the device trocar tip near the chest wall or skin margin.**
- **The EnCor® Breast Biopsy Probe has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination. Medical devices — particularly those with long and small lumina, joints, and/or crevices between components — are difficult or impossible to clean fully. Risks of disease transfer with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can prevent the complete sterilization of the device with pyrogens or microorganisms which may lead to infectious complications.**
- **Do not resterilize. After resterilization, the material's biocompatibility is impaired because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and resterilization of the single-use medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by mechanical and thermal changes.**
- **Di(2-ethylhexyl)phthalate (DEHP) is a plasticizer used in some polyvinyl chloride medical products. DEHP has been known to produce a range of adverse effects in experimental animals, notably liver toxicity and testicular atrophy. Although the toxic and carcinogenic effects of DEHP have been well established in experimental animals, the risk to humans (if any) of the compound to produce adverse effects in humans is controversial. There is no evidence that neonates, infants, pregnant and breast feeding women exposed to DEHP experience any related adverse effects. However, a lack of evidence of causation between DEHP-PVC and any disease or adverse effect does not mean that no such link exists.**
- **This device is not recommended for use in patients with breast implants.**
- **Do not use in the presence of infection.**
- **After use, this product is a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.**

Precautions

- A scalpel should be used to incise the skin prior to inserting the trocar tip.
- This device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use the EnCor® Breast Biopsy Probe if the product sterile barrier system or its packaging is compromised.
- Carefully inspect the device prior to use to verify that device has not been damaged. Do not use if product damage is evident and/or needle is blunt.
- Do not autoclave or immerse the EnCor® Breast Biopsy Driver in liquid.

Potential Complications

Potential complications may include, but are not limited to hematoma, hemorrhage, infection, adjacent tissue injury, pain, allergic reaction, and tissue adherence to the biopsy probe during removal from the breast (as per routine biopsy procedures, it may be necessary to cut tissue adhering to the stylet or coaxial while removing it from the breast).

Equipment Required

The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- EnCor® Breast Biopsy Systems which include EnCor®, EnCor Essentra®, or EnCor Ultra®
- EnCor® Breast Biopsy Driver
- EnCor® Breast Biopsy Probe
- Bard® Breast Biopsy Markers (if applicable)
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary

How Supplied

The EnCor® Breast Biopsy Probe is supplied sterile and is intended for single use only.
The EnCor® Breast Biopsy Driver is reusable and is supplied non-sterile.

Cleaning and Maintenance

Care should be exercised prior to use to ensure that the method of cleaning or handling has not caused deterioration to components.

Clean the EnCor® Breast Biopsy Driver by wiping all external surfaces with a cloth soaked with Dispatch Hospital cleaner or other hospital approved cleaner/disinfectant.

Directions For Use

Refer to Figures 1 and 2.

Figure 1. EnCor® Breast Biopsy Probe

1. Prepare the EnCor® Breast Biopsy System (EnCor®, EnCor Essentra®, EnCor Ultra®) according to the respective Instructions for Use.
2. Turn on the EnCor® Breast Biopsy System (EnCor®, EnCor Essentra®, EnCor Ultra®) and allow it to initialize. The EnCor® Breast Biopsy System will indicate when the EnCor® Breast Biopsy Driver initialization was successful and when the system is ready for the EnCor® Breast Biopsy Probe installation.

 Note: The EnCor® Breast Biopsy Driver will not initialize with the EnCor® Breast Biopsy Probe installed.
3. For stereotactic table use, refer to the Instructions for Use provided with the appropriate table adapter for proper installation of the EnCor® Breast Biopsy Driver.
4. Inspect the EnCor® Breast Biopsy Probe package to ensure the package integrity has not been compromised. The product is sterile unless the package has been previously opened or damaged. Do not use the product if the sterility has been compromised.
5. Using standard aseptic technique, remove the EnCor® Breast Biopsy Probe from the package and check for damage. Do not use the probe if damage is evident. Align the flat on the probe hub to the EnCor® Breast Biopsy Driver by sliding the distal end into place and then pressing down on the proximal end to lock (see Figure 2). Confirm that the probe is secure and that a solid audible tone connection to the EnCor® Breast Biopsy System (EnCor®, EnCor Essentra®, EnCor Ultra®) according to the appropriate Instructions for Use. Remove the probe cover from the tip.

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Copyright © 2012 C.R. Bard, Inc. All rights reserved.
Copyright ©2012 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2012 C. R. Bard, Inc. Todos los derechos reservados.
Copyright ©2012 C.R. Bard, Inc. Alle rechten voorbehouden.