smith&nephew

DYONICS Disposable

Arthroscopic Blades

1060595-Y

EC REP

Smith & Nephew Orthopaedics GmbH
AlemannenstraBe 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

Endoscopy

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

CE 2797

08/2020 1060595 Rev. Y

DYONICS™ Disposable
Arthroscopic Blades

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Device Description

Smith & Nephew DYONICS Disposable Arthroscopic Blades are rotary cutting devices for use in arthroscopic and endoscopic surgical applications for the resection of soft tissue and bone. Resected tissue is removed from the surgical site by suction through the inner lumen of the blade. DYONICS Disposable Arthroscopic Blades are made of surgical stainless steel and are available with a variety of cutting window and burr configurations, as well as straight, curved, retractable sheath, and curved rotatable styles.

Intended Use

Smith & Nephew DYONICS Disposable Arthroscopic Blades are intended for resection of soft and osseous tissue.

Indications for Use

Smith & Nephew DYONICS Disposable Arthroscopic Blades are indicated for resection of soft and osseous tissues in both large and small articular cavities.

Contraindications

- The Smith & Nephew DYONICS POWER Control Systems and blades should not be used with patients exhibiting ankylosis and also lacking adequate joint space or distention for arthroscopic inspection. Abrasion arthroplasty may not be effective in treating heavy patients or those with ankylosis, instability, or expectations beyond the relief of pain.
- Varus or valgus deformity is not in itself a contraindication in patients with good range of motion and without gross instability or extreme malalignment (15° varus, 30° valgus).
- Intra cortical abrasion arthroplasty may be contraindicated in patients not qualifying for high tibial osteotomy or total knee replacement.
- Synovectomy is contraindicated when the disease has progressed beyond the phase of synovial proliferation and advanced rheumatoid arthritis when erosion of the articular cartilage is present.

Warnings

Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.

- Contents are sterile unless package is opened or damaged. DO NOT RESTERILIZE. For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- DYONICS Disposable Arthroscopic Blades contain magnets that may be attracted to other instrumentation when placed in direct contact.
- When using a footswitch, ensure that the footswitch is not inadvertently activated while inserting the blade or plugging the connector cable into the control unit.
- Periodic irrigation of the blade is recommended to provide adequate cooling of the blade and to prevent accumulation of excised materials in the surgical site. Ensure that suction of 128 mmHg minimum is flowing while the instrument is running.

Precautions

℞ ONLY U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to:
- Biological Contamination (bacteria, virus, and/or protein) which may lead to patient infection;
- Device Malfunction (device resects too slowly, or does not resect at all).
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- DYONICS Disposable Arthroscopic Blades are packaged as a set. They must be used as supplied. Do not interchange blade components.
- Do not allow the rotating portion of any blade or burr to touch any metallic object such as a cannula or arthroscope. Damage to both instruments is likely. Damage to the blade can range from a slight distortion or dulling of the blade edge to actual fracture of the tip in vivo. If such contact does occur, inspect the tip. If there are cracks, fractures or dulling, or if there is any other reason to suspect a blade is damaged, replace it immediately.
- Excessive "side-loading" on the blade during use does not improve cutting performance and in extreme cases may result in wear and degradation of the inner assembly.
- Do not tum the adjustment collet beyond the limits of travel, doing so may result in damage to the blade.
- Irreversible damage to blades or burrs will result if they are run without the flow of irrigation (dry).
- All Smith & Nephew DYONICS Disposable Arthroscopic Blades are preassembled, packaged sterile, and ready for use. Irreversible damage will result from any attempt to disassemble curved blades. Damaged blades will not be replaced by Smith & Nephew.
- Make sure the handpiece is off while changing blade tip position.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

Instructions for Use

To remove a disposable Arthroscopic blade from its sterile package, peel TYVEK™ seal off blade tray and proceed with instructions below.

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Blade Insertion and Removal

For use with hand-controlled and non-hand-controlled handpieces, POWERMINI® Handpieces, and UltraLight handpieces.

1. To insert a blade, orient the handpiece so the key slot (A) and the key (B) can be seen on the blade (Figure 1). Insert the blade into the handpiece so that the key goes into the slot. The key should be pushed in proximally as far as possible for a positive lock.


Figure 1.

2. To remove a blade, depress the key on the blade using a thumb or fingertip. This will disengage the blade hub from the drive shaft. Withdraw the blade from the handpiece.

Note: Refer to the Smith & Nephew DYONICS Handpieces Instructions for Use (REF 10600266) for detailed instructions for use.

Positioning Smith & Nephew DYONICS ORBIT' Curved Blades

CAUTION: *Make sure handpiece is off while changing blade tip position.*

1. Grasp blade hub or handpiece (if blade is already inserted) with one hand and black collet with the other (Figure 2). Rotate collet to desired blade curve position (Figure 3).


Figure 2.


Figure 3.

2. Lock position by applying pressure and pushing hub (or handpiece and hub) and collet together until a click is felt.
3. To change blade curve position, pull hub (or handpiece and hub) and collet apart until another click is felt. Repeat steps 1 and 2.

Note: Blade tip opening will rotate relative to handpiece rotation.

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Positioning Smith & Nephew DYONICS Retractable Sheath Blades

CAUTION: *Do not turn the adjustment collet beyond the limits of travel doing so may result in damage to the blade.*

1. Grasp the blade hub or the handpiece (if the blade is already inserted) and rotate the collet to the desired sheath position. The sheath position is secure when a click is felt.

2. To retract the sheath to expose the burr, rotate the collet in the direction indicated in Figures 4 and 5.


Figure 4


Figure 5

3. To advance the sheath to cover the burr, rotate the collet in the direction indicated in Figures 6 and 7.


Figure 6.


Figure 7

Notes:

- Smith & Nephew DYONICS Disposable Arthroscopic Blades may only be used with Smith & Nephew handpieces.
- Ultra series curved blades and the Ultra series TURBOWHISKER® Blade may be used in all Smith & Nephew shaver handpieces except the mini handpieces. Standard curved blades and the standard TURBOWHISKER Blade will not fit Smith & Nephew DYONICS UltraLight Handpieces.
- The Smith & Nephew DYONICS ELITE® Burrs will fit only the DYONICS ELITE Handpiece.
- DYONICS Mini Blades will fit only the DYONICS Mini Motor Drive.
- DYONICS POWERMINI Blades will fit only the DYONICS POWERMINI Handpieces.

Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

www.smith-nephew.com
+1 978 7491000
+1 978 7491108 Fax
+1 800 343 5717 U.S. Customer Service

Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office.
All trademarks acknowledged.
Covered by one or more of the following U.S. patent numbers:
5,270, 622; 5,510,070; 5,563,481; 5,602,449;
5,620,447; 5,630,826; 5,672,945; 5,712,543; 5,745,647;
5,749,885; 5,804,936; 5,833,692; 5,871,493; 5,913,867;
Des. 390,956.
©2001,2020 Smith & Nephew, Inc. All rights reserved.

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DYONICS BONECUTTER 1회용 관절경 칼날

www.smith-nephew.com
+1 978 749 1000
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

DYONICS® BONECUTTER® Disposable Arthroscopy Blades

STERILE R

Device Description

Smith & Nephew DYONICS® BONECUTTER® Disposable Arthroscopy Blades are stainless steel rotary cutting devices for use in arthroscopic and endoscopic surgical applications for the resection of soft tissue and bone. Resected tissue is removed from the surgical site by suction through the inner lumen of the blade.

Smith & Nephew DYONICS BONECUTTER blades that run up to 10,000 rpm forward speed fit only the DYONICS POWERMAX® ELITE Handpieces.

This device is provided sterile, for single use only.

Intended Use

Smith & Nephew DYONICS Disposable Arthroscopy Blades are intended for resection of soft and osseous tissues.

Indications for Use

Smith & Nephew DYONICS BONECUTTER Disposable Arthroscopy Blades are indicated for resection of soft and osseous tissues in large and small articular cavities. Specific applications may include knee and shoulder reconstructive procedures such as notchplasty and acromioplasty.

Contraindications

- The DYONICS POWER System and blades should not be used with patients exhibiting ankylosis, and also lacking adequate joint space or distention. Abrasion arthroplasty may not be effective in treating heavy patients or those with ankylosis, instability, or expectations beyond the relief of pain.
- Varus or valgus deformity is not in itself a contraindication in patients with good range of motion and without gross instability or extreme malalignment (15° varus, 30° valgus).
- Intracortical abrasion arthroplasty may be contraindicated in patients not qualifying for high tibial osteotomy or total knee replacement.
- Synovectomy is contraindicated when the disease has progressed beyond the phase of synovial proliferation, and advanced rheumatoid arthritis when erosion of the articular cartilage is present.

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- **Contents are sterile unless package is opened or damaged. DO NOT RESTERILIZE. For single use only. Discard any open, unused product. Do not use after the expiration date.**
- **It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.**
- **Read these instructions completely prior to use.**
- **DYONICS BONECUTTER Disposable Arthroscopy Blades are packaged as a set. They must be used as supplied. Do not interchange blade components.**
- **DYONICS BONECUTTER Disposable Arthroscopy Blades contain magnets that may be attracted to other instrumentation when placed in direct contact.**
- **When using a footswitch, ensure that the footswitch is not inadvertently activated while inserting the blade or plugging the connector cable into the control unit.**
- **Periodic irrigation of the tip of the burr or cutter is recommended to provide adequate cooling of the blade and to prevent accumulation of excised materials in the surgical site. Ensure that suction of 128 mmHg minimum is flowing while the instrument is running.**

Precautions

℞ U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Do not allow the rotating portion of any blade or burr to touch any metallic object such as a cannula or arthroscope. Damage to both instruments is likely. Damage to the blade can range from a slight distortion or dulling of the blade edge to actual fracture of the tip *in vivo*. If such contact does occur, inspect the tip. If there are cracks, fractures or dulling, or if there is any other reason to suspect a blade is damaged, replace it immediately.
- Excessive "side-loading" on the blade during use does not improve cutting performance, and in extreme cases may result in wear and degradation of the inner blade.
- Irreversible damage to blades or burrs will result if they are run without the flow of irrigation (dry).
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.


Instructions for Use

To remove an arthroscopy blade from its sterile package, peel Tyvek® seal off blade tray and proceed with instructions below.

Blade Insertion and Removal

For use with hand-controlled and non-hand-controlled motor drive units.

1. To insert a blade, orient the motor drive so that the key slot (A) and the key (B) on the blade are visible. Insert the blade into the motor drive so the key goes into the slot. The key should be pushed in proximally as far as possible for a positive lock.

2. To remove a blade, depress the key on the blade using a thumb or fingertip. This will disengage the blade hub from the driveshaft. Withdraw the blade from the motor drive.

For Use with Mini Motor Drives

1. To insert a blade, slide the release ring on the distal end of the motor drive toward the proximal end as shown. This will reveal the key slot (A) on the motor drive (Figure 1). Orient the motor drive so that you can see the key slot (A) and the key (B) on the blade. Insert the blade into the motor drive so that the key goes into the slot. Let the release ring slide back to its original position.
2. To remove a blade, slide the release ring toward the proximal end of the motor drive unit and simultaneously withdraw the blade.

Note: If suction is being applied you will perceive a slight force holding the components together.

Note: Smith & Nephew DYONICS BONECUTTER blades that run up to 10,000 rpm forward speed fit only the DYONICS POWERMAX ELITE Handpieces.

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®Trademark of Smith & Nephew, certain marks registered U.S. Patent & Trademark Office. All trademarks acknowledged.

Covered by one or more of the following U.S. Patent numbers: 5,270,622; 5,322,505; 5,510,070; 5,563,481; 5,602,449; 5,620,447; 5,630,826; 5,672,945; 5,707,350; 5,712,543; 5,745,647; 5,749,885; 5,804,936; 5,833,692; 5,871,493; 5,913,867; 6,090,122; 6,328,752; Des. 381,425; Des. 390,955; Des. 390,956; Re. 34,556.

©2001, 2020 Smith & Nephew, Inc. All rights reserved.

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