Information
Instructions For Use

ZEISS Prescription Spectacle Lens

Carl Zeiss Vision GmbH
ZEISS
Regulatory Affairs

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Table of Content

Intended use and description of the principles of operation of the device and its mode of action **3**

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Intended use and description of the principles of operation of the device and its mode of action

A prescription spectacle lens is a glass or plastic device that is a lens intended to be worn by a patient in a spectacle frame to provide refractive corrections in accordance with a prescription for the patient. Corrective lenses are typically prescribed to meet the visual needs of a patient. Objective and subjective refraction methods allow to determine the patient's ametropia. Typically, wearers of corrective lenses can improve their visual performance and compensate deficits from near-sightedness (Myopia), far-sightedness (Hyperopia) or degradation of accommodation mostly as age-related loss of flexibility of the eye's crystalline lens (Presbyopia) together with effects of astigmatism, squinting and similar mechanisms.
The prescription consists of all the specifications necessary to make the lens according to an individual patient's visual needs. Prescriptions typically include the power specifications of each lens (typically 2 lenses for both eyes of an individual patient).
Corrective spectacle lenses come in distinct design configurations that can be grouped into 3 major groups. Ophthalmic lenses may be single vision, multi-focal or varifocal (or progressive addition) lenses. Single Vision lenses have one visual zone, either for far or near distance. Multi-focal lenses have two or three visual zones, separated by visible sharp delineation. Variable power lenses consist of a continuum of focal points, with transition from distance correction to intermediate correction to near correction (no focal segment lines).

General Instructions

The product is to be used as per the instruction of the Eye Care Professional based on the patient's individual prescription.
The intended use of ZEISS ophthalmic spectacle lenses is vision correction. Unless specifically ordered as safety glasses, **they do not claim to protect the eye from physical impact** because they are neither shatterproof nor unbreakable and are not intended for use as personal protective equipment.

Do not clean spectacle lenses with clothing, towels, paper or household cleaners. Use ZEISS cleaning products or just water, a pH-neutral cleaning agent and a microfiber or soft cotton cloth. Dry cleaning with an ordinary cloth or household paper can cause damage to the lens surface. Only use a dedicated lens cleaning cloth.

ZEISS Lens Wipes

1.  Using the folded lens wipe gently clean away dust and dirt
2.  Unfold the lens wipe and clean the surface using circular movements
3.  Wait a short time period until the lens surface is dry

ZEISS Cleaning Spray

1.  Unfold a ZEISS Microfiber cloth and gently clean away dust and dirt
 → If lenses are heavily soiled, rinse under water before application of cleaning spray.

Carl Zeiss Vision GmbH Regulatory Affairs

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2.  Apply ZEISS Lens Cleaner to both surfaces of each lens

3.  Wipe dry with the microfiber cloth

Without ZEISS Cleaning Solutions

1.  Rinse the lens under cold running water

2.  Use a pH-neutral cleaning agent

3.  Wipe dry with a microfiber or cotton cloth

Plastic lenses should not be exposed to short-time high temperature fluctuations; this can cause stress within the material and coating damage may result.
Short exposure time to temperatures above 80 °C (e.g. sauna visit), or 60 °C of long exposure time (e.g. on the dashboard of a car) can cause damage to the lens / coating.
Lenses with less than 8% transmission are not suitable for driving.

 Not suitable for driving


Lenses with less than 25% transmission are not suitable for driving at night.

 Not suitable for driving at night

If treated carefully by avoiding damage to the lens, the service life (life-time, durability) of lenses is only limited by the potential need to replace the lenses due to a change of the vision correction prescription needed by the patient, resulting in the need for new prescription lenses.

Storage conditions

- do not store on the dashboard of a care
- do not use in a Sauna at temperatures above 80°C
 avoid extreme temperature differences

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Symbols on the label

caution, read the Instructions for Use

date of manufacture

Unique Device Identifier

Lot number

CE Mark

Product Owner

Medical Device

L, R

Left eye, right eye

The spectacle lenses are provided with the CE marking.

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