9.75 x 9.75 inches (247.65 mm x 247.65 mm)
2020/Sep/24 9:43 AM

EverFlex™

Self-expanding Peripheral Stent with Entrust™ Delivery System
Selvekspanderende, perifert stent med Entrust™ fremføringssystem
Selbstexpandierender peripherer Stent mit dem Einführsystem Entrust™
Σύστημα αυτοδιαστελλόμενου περιφερικού stent με σύστημα τοποθέτησης Entrust™
Sistema de stent periférico autoexpansible con sistema de liberación Entrust™
Itsestään laajeneva perifeerinen stentti ja Entrust™-sisäänviejäjärjestelmä
Stent périphérique auto-expansible avec système de pose Entrust™
Entrust™ bevezetőrendszerrel szerelt öntáguló perifériás sztent
Stent periferico autoespandibile con sistema di rilascio Entrust™
Zelfontplooiende perifere stent met Entrust™-plaatsingssysteem
Selvekspanderende perifer stent med Entrust™-innføringssystem
Obwodowy stent samorozprężalny z systemem wprowadzającym Entrust™
Stent periférico autoexpansível com sistema de entrega Entrust™
Саморасширяющийся периферический стент с системой доставки Entrust™
Självexpanderande perifer stent med införingssystemet Entrust™
Kendinden Genişleyen Periferik Stent ile Entrust™ Taşıma Sistemi

CE
2797

510531-001 Rev. B

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Instructions for use

Device description

The EverFlex™ self-expanding peripheral stent with Entrust™ delivery system is a self-expanding nitinol stent system intended for permanent implantation. The EverFlex stent is made of a nickel titanium alloy (nitinol) and comes pre-mounted on a 5 Fr, 0.035 in over-the-wire (OTW) delivery system. The stent is cut from a nitinol tube into an open lattice design. The stent has tantalum radiopaque markers at the proximal and distal ends. After deployment, the stent opens to its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The Entrust delivery system, shown in **Figure 1**, has a triaxial shaft design comprised of an inner shaft (1), silver retractable outer sheath (2), and a gold isolation sheath (3) for ease of use during deployment. The ergonomic deployment handle (8) has a removable safety locking pin (9), deployment wheel (10), and luer hub (6).

The distal portion of the delivery system (enlarged section) contains two radiopaque markers; one marker distal (4) and one marker proximal (5) to the constrained stent. The distal marker is mounted on the retractable outer sheath and the proximal marker is mounted on the inner shaft. The EverFlex stent is constrained within the outer sheath. The flexible catheter tip (7) is at the distal end of the outer sheath.

Before the procedure, the operator flushes saline through the guidewire lumen until it exits at the end of the catheter. Radiopaque markers on the EverFlex stent give guidance for positioning the stent before deployment. After the stent is correctly positioned, the operator removes the safety locking pin to unlock the deployment wheel and enable the outer sheath to retract. The operator retracts the outer sheath by slowly rotating the deployment wheel in the direction indicated by the arrow on the deployment handle. When the radiopaque marker on the outer sheath passes the proximal radiopaque marker on the inner shaft, stent deployment is complete.


Figure 1. Entrust delivery system

Indications for use

The EverFlex self-expanding peripheral stent with Entrust delivery system is indicated for use in occlusions, lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA); or lesions believed to be at high risk for restenosis following PTA in the common iliac, external iliac, superficial femoral, proximal popliteal, or subclavian arteries. Stenting is intended to improve and maintain artery luminal diameter.

Contraindications

The EverFlex self-expanding peripheral stent with Entrust delivery system is contraindicated under the following conditions:

- Patients for whom anticoagulant or antiplatelet therapy is contraindicated
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Aneurysm of the artery to be treated
- Patients with known hypersensitivity to nickel titanium
- Patients with uncorrected bleeding disorders
- Use in the carotid artery
- All of the customary contraindications for PTA

Warnings

- This device is provided **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- Do not force stent deployment if high resistance is met when initially rotating the thumbwheel. Carefully withdraw the stent system without deploying the stent.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system, the vessel lumen, or cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Verify that the red locking pin is still engaged with the handle. Do not use the stent or the delivery system if damage is detected.
- Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
- Do not use the stent if it was partially deployed in the package, or if it is partially deployed before starting the deployment procedure.
- Always use an introducer sheath during the implant procedure to protect both the vessel and the puncture site. Support from an introducer sheath is necessary to minimize stent lengthening or shortening during stent deployment.
- Using a 0.014 in (0.36 mm) or 0.018 in (0.46 mm) guidewire in a highly tortuous vessel, highly stenosed lesion, or calcified lesion could result in crossing, tracking, or deployment complications.
- In lesions with a vessel diameter less than 2 mm, crossing or deployment complications may occur.
- Do not torque (turn in multiple rotations) the delivery system. Torquing the delivery system may lead to difficulty or make it impossible to deploy or move the stent system over the guidewire.
- Do not allow slack in the delivery system during deployment. Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.
- Completely remove the safety locking pin and tube assembly from the handle, as depicted in **Figure 2**. If the safety locking pin and tube assembly are not removed, the EverFlex stent may partially deploy, fail to deploy, or deploy with increased deployment force.
- The stent system is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition.
- Failure to hold the handle in a fixed position may cause foreshortening or lengthening of the stent.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the EverFlex stent beyond its nominal diameter.

Potential adverse events

The potential adverse effects (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

Abrupt or subacute closure	Dissection or intimal flap	Pseudoaneurysm
Allergic reaction to device materials or procedure medications	Edema	Renal failure requiring dialysis
Allergic reaction to nitinol	Embolism	Renal insufficiency (new or worsening)
Amputation	Failure to deploy stent	Restenosis
Aneurysm	Fever	Sepsis
Angina	Gastrointestinal bleeding due to anticoagulation	Shock
Arrhythmia	Hematoma	Stent collapse or fracture
Arteriovenous fistula	Hypertension or hypotension	Stent migration
Artery perforation or rupture	Infection	Stent misplacement
Bleeding requiring transfusion	Inflammation	Stroke
Bruising	Intraluminal thrombus	Surgical or endovascular intervention
Contrast medium reaction / renal failure	Myocardial infarction	Thrombosis or occlusion of the stent
Death	Pain	Transient ischemic attack
Device breakage	Partial stent deployment	Venous thromboembolism
		Vessel spasm

Directions for use

Preparation

Warning: The device is provided sterile for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

Recommended items

- 5 to 10 ml syringe filled with heparinized saline
- 0.035 in (0.89 mm) guidewire
- Hemostatic sheath
- PTA balloon

Select stent size

Measure the diameter of the reference vessel (proximal and distal to lesion). Measure the length of the target lesion. Refer to **Table 1** to select a stent size according to the vessel diameter and a stent length sufficient to extend proximal and distal to the target lesion.

Table 1. Stent diameter and length sizing

Stent diameter (mm)	Reference vessel diameter (mm)	Stent length (mm)
5	3.5 – 4.5	20, 40, 60, 80, 100, 120, 150
6	4.5 – 5.5	20, 40, 60, 80, 100, 120, 150
7	5.5 – 6.5	20, 40, 60, 80, 100, 120, 150
8	6.5 – 7.5	20, 40, 60, 80, 100, 120, 150

Stent delivery system preparation

1. Carefully inspect the pouch to ensure that there is no damage to the sterile barrier. Open the pouch and remove the tray and contents. Place the tray on a flat surface and carefully remove the stent delivery system.
 Caution: Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Verify that the red locking pin is still engaged with the handle. Do not use the stent or the delivery system if damage is detected.
2. Attach a 5 to 10 ml syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until fluid passes from the catheter tip.
 Caution: Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
3. Examine the distal end of the catheter to ensure that the stent is fully contained within the tip of the catheter.
 Caution: Do not use the stent if it was partially deployed in the package, or if it is partially deployed before starting the deployment procedure.

Stent deployment

Sheath and guidewire insertion

1. Gain access using an introducer sheath with a hemostatic valve that is compatible with a 5 Fr (1.7 mm) delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial artery or axillary artery.
 Caution: Always use an introducer sheath during the implant procedure to protect both the vessel and the puncture site. Support from an introducer sheath is necessary to minimize stent lengthening or shortening during stent deployment.
2. Insert a 0.035 in (0.89 mm) guidewire of appropriate length across the target lesion.
 Caution: Using a 0.014 in (0.36 mm) or 0.018 in (0.46 mm) guidewire in a highly tortuous vessel, highly stenosed lesion, or calcified lesion could result in crossing, tracking, or deployment complications.

Lesion dilation

Predilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire.

Caution: In lesions where vessel lumen diameter is less than 2 mm, crossing or deployment complications may occur.

Stent delivery system introduction

Advance the stent delivery system over the guidewire through the hemostatic valve and introducer sheath.

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Warning: If resistance is encountered at any time during the insertion procedure, do not force passage of the stent system. Forced passage may damage the stent system, the vessel lumen, or cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Caution: Do not torque (turn in multiple rotations) the delivery system. Torquing the delivery system may lead to difficulty or make it impossible to deploy or move the stent system over the guidewire.

Stent deployment

1. Advance the delivery system until the distal (leading) radiopaque marker is distal to the target lesion.
2. Pull back on the delivery system until there is no slack in the delivery system and until the radiopaque markers extend distal and proximal to the target lesion.
 Caution: Do not allow slack in the delivery system during deployment. Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.
3. Completely remove the safety locking pin and tube assembly by gently pulling in the direction indicated by the arrow on the safety locking pin. Discard the safety locking pin and the attached tube assembly.
 Caution: Completely remove the safety locking pin and tube assembly (2) from the handle (1), as depicted in **Figure 2**. If the safety locking pin and tube assembly are not removed, the EverFlex stent may partially deploy, fail to deploy, or deploy with increased deployment force.


Figure 2. Entrust delivery system with safety locking pin and tube assembly

1. Handle
2. Locking pin and tube assembly
3. Pull down to remove the safety locking pin and tube assembly

4. Visually confirm that the radiopaque markers still extend distal and proximal to the target lesion. Reposition if necessary.
5. Initiate stent deployment by holding the handle in a fixed position and slowly rotating the thumbwheel in the direction indicated by the arrow on the handle, see **Figure 2**.
 Note: If a second hand is used to hold the delivery system, gently support the delivery system at the gold isolation sheath (1). Do not constrict the silver outer sheath (2) during deployment. See **Figure 3**.


Figure 3. Silver retractable outer sheath and gold isolation sheath

1. Gold
2. Silver

Caution: The stent system is not designed for recapturing at any time, or designed for repositioning after establishing vessel apposition.

6. When the distal radiopaque markers of the stent are approximately 1 to 3 mm distal to marker on the retractable outer sheath, use fluoroscopic guidance to reposition the stent as needed before achieving vessel apposition.
7. Keep the entire length of the flexible deployment system as straight as possible while the stent is being released. Hold the handle in a fixed position to ensure that no slack is introduced to the delivery system. Deployment is complete when the marker on the outer retractable sheath reaches the inner shaft marker and the stent is released.
 Warning: Do not force stent deployment if high resistance is met when initially rotating the thumbwheel. Carefully withdraw the stent system without deploying the stent.
 Caution: Failure to hold the handle in a fixed position may cause foreshortening or lengthening of the stent.
 Caution: The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
 Note: If more than one stent is needed, place the distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.

After stent deployment

1. After stent deployment, use fluoroscopic guidance to withdraw the entire delivery system as one unit, over the guidewire, into the introducer sheath and out of the body. Remove the delivery system from the guidewire.
2. Using fluoroscopy, verify that the EverFlex stent is fully deployed.
3. If the EverFlex stent is not completely expanded at any point along the lesion, postdeployment balloon dilation can be performed at the discretion of the physician.
 Caution: Use caution when crossing a deployed stent with any adjunct device.
 Caution: Do not expand the EverFlex stent beyond its nominal diameter.
4. To dilate the stent, select an appropriately sized PTA balloon catheter and dilate the EverFlex stent using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel.
5. Confirm that the stent is fully expanded, then remove the PTA balloon from the patient.
6. Remove the guidewire and introducer sheath from the body.
7. Close entry wound.
8. Discard the delivery system, guidewire, and sheath.

MRI information

MR Conditional

Nonclinical testing demonstrated that the EverFlex stent is MR Conditional. A patient with this device can be scanned safely, immediately after stent placement under the following conditions:

- Static magnetic field of 3.0 Tesla or 1.5 Tesla
- Maximum spatial gradient magnetic field of 8000-Gauss/cm (extrapolated) or less (80 T/m)
- Normal operating mode (maximum whole body averaged SAR of 2.0 W/kg) for 15 minutes of scanning

MRI related temperature rise

In nonclinical testing, the EverFlex self-expanding peripheral stent produced the following temperature rises during MRI performed for 15 minutes of scanning (that is, per pulse sequence) in 1.5 Tesla/64 MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3.0 Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Single, 8 mm x 80 mm	**1.5-T**	**3.0-T**
MR system reported, whole body averaged SAR (W/kg)	2.9	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.0	2.4
Single, 8 mm x 120 mm	**1.5-T**	**3.0-T**
MR system reported, whole body averaged SAR (W/kg)	2.9	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.9	3.7

These temperature changes do not pose a hazard to a patient under the identified test conditions. Patients are recommended to register conditions under which the implant may be scanned safely with the MedicAlert Foundation (www.medicalert.org) or equivalent organization.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent. Therefore, it may be necessary to optimize MR imaging parameters to compensate for the presence of this device.

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Instructions for use

Device description

The Protégé™ EverFlex™ self-expanding stent system is a self-expanding nitinol stent system intended for permanent implantation. The stent is made of a nickel titanium alloy (nitinol) and comes premounted on a 6 Fr (0.035 in /0.89 mm) over-the-wire (OTW) delivery system. The stent is cut from a nitinol tube in an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent. After deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The delivery system, as shown in **Figure 1**, is comprised of an inner shaft (1) and outer sheath (2), which are secured with a safety lock (3). The nylon inner shaft terminates distally in a flexible catheter tip (4) and originates proximally at the hub (5).

The distal portion of the delivery system for the 20 to 150 mm stents, as shown in **Figure 2**, is comprised of two radiopaque markers, one marker distal (6) and one retainer marker proximal (7) to the constrained stent, on the inner shaft.

The distal portion of the delivery system for the 200 mm stents, as shown in **Figure 3**, is comprised of the same components as those in **Figure 2** except for the following radiopaque markers: one retainer marker distal (13) and one holder marker proximal (14) to the constrained stent, on the inner shaft.

As shown in **Figure 1**, the outer sheath connects proximally to the manifold (8). The self-expanding stent is constrained within the space between the inner shaft and the outer sheath. This space is flushed before the procedure through the stopcock (9). The outer sheath has a radiopaque marker at its distal end (10), as shown in **Figures 2** and **3**.

Stent positioning at the target lesion (peripheral indication) or the stricture (biliary indication) is achieved before deployment by using the two radiopaque markers on the inner shaft, which mark the location of the constrained stent. Before stent deployment, turn the safety lock counterclockwise to unlock the outer sheath.

Retract the outer sheath by pulling the distal grip (11) toward the proximal grip (12). The stent is completely deployed when the radiopaque marker on the outer sheath passes the proximal radiopaque marker on the inner shaft.


Figure 1. Stent on delivery system


Figure 2. Distal portion of 20 to 150 mm delivery system


Figure 3. Distal portion of 200 mm delivery system

Indications for use

Peripheral

The stent is indicated for use in the following types of occlusions or lesions:

- Occlusions or lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA)
- Lesions that appear to be at high risk for restenosis following PTA in the common iliac, external iliac, superficial femoral, proximal popliteal, or subclavian arteries

Stenting is intended to improve and maintain the artery luminal diameter.

Biliary

The stent is indicated as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications

The Protégé EverFlex self-expanding peripheral stent system is contraindicated under the following conditions:

Peripheral

- Patients in whom anticoagulant and/or antiplatelet therapy is contraindicated
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Aneurysm of the artery to be treated
- Patients with known hypersensitivity to nickel-titanium
- Patients with uncorrected bleeding disorders
- All the customary contraindications for PTA

Biliary

- None known

Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.
- If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.
- If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.

Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
- Do not use the stent if it is partially deployed.
- Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
- Always use a sheath during the implant procedure to protect the vessel (peripheral indication) or the biliary tract (biliary indication), and the puncture site.
- If the lesion is not predilated, it may be difficult to properly position or remove the stent system.
- The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).
- Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- The effects of overlapping stents have not been evaluated.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the stent beyond its nominal diameter.

Potential adverse events

Potential adverse events from the stent implantation procedure include the usual complications associated with any percutaneous procedure, PTA, or stent placement.

Table 1. Potential adverse events

Risks associated with any percutaneous procedure	Risks associated with PTA and stent placement
• Arteriovenous fistula	• Abrupt or subacute closure
• Allergic reaction to device materials or procedure medications	• Amputation
• Aneurysm	• Persons with allergic reaction to nickel titanium (nitinol) may suffer an allergic response to this implant
• Artery perforation or rupture	• Device breakage
• Bleeding requiring transfusion	• Dissection or intimal flap
• Contrast medium reaction or renal failure	• Failure to deploy stent
• Death	• Gastrointestinal bleeding due to anticoagulation
• Duct rupture resulting from overstretching the duct (biliary only)	• Partial stent deployment
• Embolism	• Renal insufficiency (new or worsening)
• Fever	• Renal failure requiring dialysis
• Hematoma	• Restenosis
• Hypertension or hypotension	• Stent collapse or fracture
• Infection (for biliary placement, secondary to contamination of the stent may lead to cholangitis, hemobilia, peritonitis, or abscess)	• Stent migration
• Inflammation	• Stent misplacement
• Intraluminal thrombus	• Surgical or endovascular intervention
• Myocardial infarction	• Thrombosis or occlusion of the stent
• Pseudoaneurysm	• Vessel spasm
• Sepsis	
• Shock	
• Stroke	
• Transient ischemic attack	
• Venous thromboembolism	

Storage

Store the stent system away from sunlight at room temperature in a dry place.

Directions for use

Warning: This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

Selection and preparation

Required items

- 5 to 10 cc syringe filled with heparinized saline
- 0.035 in (0.89 mm) exchange guidewire
- Hemostatic sheath
- PTA balloon (peripheral only)

Select stent size for peripheral interventions

Measure the diameter of the reference vessel (proximal and distal to lesion). Refer to **Table 2** for stent diameter sizing. Measure the length of the target lesion. Choose a stent length that extends proximal and distal to the target lesion.

Preparation for biliary interventions

1. Inject contrast medium.
2. Perform a percutaneous cholangiogram using standard technique.
3. Evaluate and mark the stricture by using fluoroscopy, and observe the most distal level of the biliary stricture.

Select stent size for biliary interventions

Measure the diameter of the reference bile duct (proximal and distal to stricture). Refer to **Table 2** for stent diameter sizing. Measure the length of the target stricture. Choose a stent length that extends proximal and distal to the tumor to protect against obstruction from further tumor growth.

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Table 2. Stent diameter sizing

Stent diameter (mm)	Reference vessel diameter (mm)	Bile duct diameter (mm)
5	3.5–4.5	3.5–4.5
6	4.5–5.5	4.5–5.5
7	5.5–6.5	5.5–6.5
8	6.5–7.5	6.5–7.5

Preparation of stent delivery system

1. Open the shelf box to reveal the pouch containing the stent and delivery catheter.
2. After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the pouch and extract the tray with contents.
3. Set the tray on a flat surface. Carefully pull off the tray lid and remove the stent and delivery system. If sterility has been compromised or the device is damaged, do not use the device.
 Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
4. Verify that the device is locked by tightening the safety lock clockwise.
 Caution: Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
5. Attach a 5 to 10 cc syringe filled with heparinized saline to the stopcock on the manifold. Open the stopcock and inject saline into the annular space between the inner shaft and the outer sheath until it comes out the outer sheath.
6. Attach a 5 to 10 cc syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until it comes out the catheter tip.
7. Ensure that the stent is contained within the outer sheath by examining the distal end of the catheter. Do not use the stent system if the stent is partially deployed. If a gap exists between the catheter tip and the outer sheath, open the safety lock and gently pull the inner shaft in a proximal direction until the gap is closed. Secure the safety lock after making the adjustment by turning the lock clockwise.
 Caution: Do not use the stent if it is partially deployed.

Stent deployment procedure

Insertion of sheath and guidewire

1. To gain access to the appropriate site, use a sheath with a hemostatic valve that is compatible with a 6 Fr (2.0 mm) delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery.
 Caution: Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
2. Insert an 0.035 in (0.89 mm) exchange guidewire of appropriate length across the target lesion (peripheral indication) or the stricture (biliary indication) via the sheath.
 Caution: Always use a sheath during the implant procedure to protect the vessel (peripheral indication) or the biliary tract (biliary indication), and the puncture site.

Dilation of lesion

Peripheral

Predilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire.
Caution: If the lesion is not predilated, it may be difficult to properly position or remove the stent system.

Biliary

Predilation is not typically performed on malignant strictures. However, if predilation is necessary, use standard balloon dilation techniques. Remove the balloon catheter from the patient while maintaining stricture access with the guidewire.

Introduction of stent delivery system

Advance the device over the guidewire through the hemostatic valve and sheath.
Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may cause damage to the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Stent deployment

1. Advance the delivery system until the distal (leading) radiopaque inner shaft marker is distal to the target lesion (peripheral indication) or the stricture (biliary indication).
 Note: When deploying the 200 mm stent delivery system, the stent moves back approximately 5 mm from the distal retainer after the initial release.
2. Pull back on the delivery system until there is no slack in the delivery system and the radiopaque inner shaft markers extend distal and proximal to the target lesion (peripheral indication) or the stricture (biliary indication).
3. Open the safety lock by turning it counterclockwise.
4. Initiate stent deployment by pinning down the inner shaft (proximal grip) in a fixed position and pulling the outer sheath (distal grip) toward the proximal grip as shown in **Figure 4**.
5. Once initial deployment is visible and before achieving vessel apposition (peripheral indication) or duct apposition (biliary indication), reposition the stent as needed using the radiopaque markers.
 Note: To ensure that there is no relative movement between the grips during repositioning, secure the safety lock.
 Caution: The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).
6. During release of the stent, keep the entire length of the flexible deployment system as straight as possible. To ensure that no slack is introduced into the delivery system, hold the proximal grip stationary and fixed. Deployment is complete when the outer sheath marker passes the proximal inner shaft stent marker and the stent is released.
 Warning: If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.


Figure 4. Stent deployment

Caution: Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.
Caution: The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
Note: If a second stent is needed, place the more distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.
Caution: The effects of overlapping stents have not been evaluated.

After stent deployment

1. While using fluoroscopy, withdraw the entire delivery system as one unit, over the guidewire, into the access sheath and out of the body. Remove the delivery system from the guidewire.
 Warning: If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.
2. Using fluoroscopy, visualize the stent to verify full deployment.
3. If the stent is not completely expanded at any point along the lesion (peripheral indication) or stricture (biliary indication), balloon dilation (standard PTA technique) can be performed after deploying the stent. To dilate the stent, select an appropriately sized PTA balloon catheter and dilate the stent by using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel. Remove the PTA balloon from the patient.
 Caution: Use caution when crossing a deployed stent with any adjunct device.
 Caution: Do not expand the stent beyond its nominal diameter.
4. Remove the guidewire and the sheath from the body.
5. Close the entry wound as appropriate.
6. Discard the delivery system, guidewire, and sheath.
 Note: Physician experience and discretion determines the appropriate drug regimen for each patient.

MRI Conditional

Stent migration and image artifact

Comparison to published test results indicate that the Protégé EverFlex stent will not exhibit substantial magnetic field interactions or migrate in an MR environment of 3.0 Tesla or less. MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent.

RF heating

Radio frequency (RF) heating in an MR environment was evaluated for the Protégé EverFlex stent. Testing was conducted on two stent/MR scanner scenarios: a 200 mm length stent with a 1.5 Tesla scanner and a 120 mm stent with a 3.0 Tesla scanner. These combinations were chosen since they represent the theoretical worst-case conditions for RF heating.
Protégé EverFlex Stent, single configuration, 8 x 200 mm length, SFA position (or overlapped stents totaling approximately 200 mm in length) was tested for RF heating at 1.5 Tesla. 3.0 Tesla RF heating testing was not performed on the 8 x 200 mm.

Recommendations:

If performing an MR procedure with a patient with a 200 mm Protégé EverFlex stent (or overlapping stents approximately 200 mm in length), exclude the stent area from RF exposure.
The local body SAR should be limited to ≤ 0.9 W/kg if an MR body coil must be used (the recommended SAR level is based on theoretical extrapolation to maintain localized temperature increase ≤ 2.0 °C based on in vitro test results, valid for test configuration only). Note that the whole body averaged SAR displayed by the MR system is inappropriate to scale exact local temperature increases.
Gradient magnetic fields: a stimulation level of 63.80% was used (1.5 T Siemens) in respect to heating issues. No tests have been performed on possible nerve or other tissue stimulation.
Protégé EverFlex stent, single configuration, 8 x 120 mm length, SFA position, was tested for RF heating at 3.0 Tesla. 1.5 Tesla RF heating testing was not performed on the 8 x 120 mm.

Recommendations:

Local-body SAR ≤ 3.1 W/kg with MR body coil (the recommended SAR level is based on theoretical extrapolation to maintain localized temperature increase ≤ 2.0 °C based on in vitro test results, valid for test configuration only). Note that the whole-body-averaged SAR displayed by the MR system is inappropriate to scale exact local temperature increases.
Gradient magnetic fields: stimulation level of 74.19% was used (3T Siemens) in respect to heating issues. No tests have been performed on possible nerve or other tissue stimulation.

Warranty

This product is subject to standard Medtronic warranty terms.

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Instructions for use

Device description

The Protégé™ GPS™ self-expanding nitinol stent and delivery technology is a self-expanding nitinol stent system intended for permanent implantation. The stent is made of a nickel titanium alloy (nitinol) and comes premounted on a 6 Fr (2.0 mm) over-the-wire (OTW) delivery system. The stent is cut from a nitinol tube in an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent. After deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The delivery system, as shown in Figure 1, is comprised of an inner shaft (1) and outer sheath (2), which are secured with a safety lock (3). The nylon inner shaft terminates distally in a flexible catheter tip (4) and originates proximally at the hub (5). Two radiopaque markers, one marker distal (6) and one retainer marker proximal (7) to the constrained stent, are mounted on the inner shaft.

The outer sheath connects proximally to the manifold (8). The self-expanding stent is constrained within the space between the inner shaft and outer sheath. This space is flushed before the procedure through the stopcock (9). The outer sheath has a radiopaque marker at its distal end (10).

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- Occlusions or lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA)
- Lesions that appear to be at high risk for restenosis following PTA in the common iliac, external iliac, or subclavian arteries

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Biliary

The stent is intended as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications

The Protégé GPS self-expanding stent system is contraindicated under the following conditions:

Peripheral

- Patients in whom anticoagulant or antiplatelet therapy is contraindicated
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Aneurysm of the artery to be treated
- Patients with known hypersensitivity to nickel titanium
- Patients with uncorrected bleeding disorders
- All the customary contraindications for PTA

Biliary

- None known

Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing or resterilizing could increase the risks of patient infection and compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.
- If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.
- If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.

Precautions

- Carefully inspect the sterile package and device before using the system to verify that no damage occurred during shipment. Do not use the device if the packaging or the device is damaged.
- Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
- Do not use the stent if it is partially deployed.
- Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
- Always use a sheath during the implant procedure to protect both the vessel and puncture site.
- If the lesion is not predilated, it may be difficult to properly position or remove the stent system.
- The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).
- Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- The effects of overlapping stents have not been evaluated.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the stent beyond its nominal diameter.

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Directions for use

Warning: This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing or resterilizing could increase the risks of patient infection and compromised device performance.

Selection and preparation

Required items

- 5 to 10 cc syringe filled with heparinized saline
- 0.035 in (0.89 mm) exchange guidewire
- Hemostatic sheath

Recommended items

- PTA balloon (peripheral only)

Stent size selection for peripheral interventions

1. Measure the diameter of the reference vessel (proximal and distal to lesion). Refer to **Table 2** for stent diameter sizing.
2. Measure the length of the target lesion.
3. Choose a stent length that extends proximal and distal to the target lesion.

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Stent size selection for biliary interventions

1. Measure the diameter of the reference bile duct (proximal and distal to the stricture). Refer to **Table 2** for stent diameter sizing.
2. Measure the length of the target stricture.
3. Choose a stent length that will extend proximal and distal to the tumor to protect against obstruction from further tumor growth.

Table 2. Stent diameter sizing

Stent diameter (mm)	Reference vessel diameter (mm)	Bile duct diameter (mm)
9	7.5–8.5	7.5–8.5
10	8.5–9.5	8.5–9.5
12	9.5–11.0	9.5–11.0
14	11.5–13.0	11.5–13.0

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1. Open the shelf box to reveal the pouch containing the stent and delivery catheter.
2. After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the pouch and extract the tray with contents.
3. Set the tray on a flat surface. Carefully pull off the tray lid and remove the stent and delivery system. If sterility has been compromised or the device is damaged, do not use the device.
 Caution: Carefully inspect the sterile package and device before using the system to verify that no damage occurred during shipment. Do not use the device if the packaging or the device is damaged.
4. Verify that the device is locked by tightening the safety lock clockwise.
 Caution: Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
5. Attach a 5 to 10 cc syringe filled with heparinized saline to the stopcock on the manifold. Open the stopcock and inject saline into the annular space between the inner shaft and the outer sheath until the saline comes out the outer sheath.
6. Attach a 5 to 10 cc syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until the saline comes out the catheter tip.
7. Ensure that the stent is contained within the outer sheath by examining the distal end of the catheter. Do not use the stent system if the stent is partially deployed.
 Note: If a gap exists between the catheter tip and the outer sheath, open the safety lock and gently pull the inner shaft in a proximal direction until the gap is closed. Secure the safety lock after making the adjustment by turning the lock clockwise.
 Caution: Do not use the stent if it is partially deployed.

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Stent deployment procedure

Insertion of sheath and guidewire

1. To gain access to the appropriate site, use a sheath with a hemostatic valve that is compatible with a 6 Fr (2.0 mm) delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery.
 Caution: Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
2. Insert an exchange guidewire of appropriate length across the target lesion via the sheath.
 Caution: Always use a sheath during the implant procedure to protect both the vessel and the puncture site.

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Caution: If the lesion is not predilated, it may be difficult to properly position or remove the stent system.

Biliary

Predilation is not performed on malignant strictures. However, if predilation is necessary, use standard balloon dilation techniques. Remove the balloon catheter from the patient while maintaining stricture access with the guidewire.

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Advance the device over the guidewire through the hemostatic valve and sheath.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Stent deployment

1. Advance the delivery system until the distal (leading) radiopaque inner shaft marker is distal to the target lesion (peripheral indication) or the stricture (biliary indication).
2. Pull back on the delivery system until there is no slack in the delivery system and the radiopaque inner shaft markers extend distal and proximal to the target lesion (peripheral indication) or the stricture (biliary indication).
3. Open the safety lock by turning it counterclockwise.
4. Initiate stent deployment by pinning down the inner shaft (proximal grip) in a fixed position and pulling the outer sheath (distal grip) toward the proximal grip as shown in Figure 2.
5. Once initial deployment is visible and before achieving vessel apposition (peripheral indication) or duct apposition (biliary indication), reposition the stent as needed using radiopaque markers.

Note: To ensure that there is no relative movement between the grips during repositioning, secure the safety lock.

Caution: The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).

6. During release of the stent, keep the entire length of the flexible deployment system as straight as possible. To ensure that no slack is introduced into the delivery system, hold the proximal grip stationary and fixed. Deployment is complete when the outer sheath marker passes the proximal inner shaft stent marker and the stent is released.

Warning: If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.


Figure 2: Stent deployment

Caution: Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.

Caution: The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.

Note: If a second stent is needed, place the more distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.

Caution: The effects of overlapping stents have not been evaluated.

After stent deployment

1. While using fluoroscopy, withdraw the entire delivery system as one unit, over the guidewire, into the access sheath out of the body. Remove the delivery system from the guidewire.
 Warning: If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.
2. Using fluoroscopy, visualize the stent to verify full deployment.
3. If the stent is not completely expanded at any point along the lesion (peripheral indication) or stricture (biliary indication), balloon dilation (standard PTA technique) can be performed after deploying the stent. To dilate the stent, select an appropriately sized PTA balloon catheter and dilate the stent by using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel. Remove the PTA balloon from the patient.
 Caution: Use caution when crossing a deployed stent with any adjunct device.
 Caution: Do not expand the stent beyond its nominal diameter.
4. Remove the guidewire and the sheath from the body.
5. Close the entry wound as appropriate.
6. Discard the delivery system, guidewire, and sheath.
 Note: Physician experience and discretion will determine the appropriate drug regimen for each patient.

MRI compatibility

The stent is MRI Safe and does not interfere with, nor is affected by, the operation of an MRI device.[1]

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[1] MR Imaging Artifacts, Ferromagnetism, and Magnetic Torque of Intravascular Filters, Stents, and Coils, Radiology 1988, 166: 657-664.

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Instructions for Use

Device description

The Visi-Pro™ balloon-expandable peripheral stent is intended as a permanently implanted device. The stent is made from a 316L stainless steel tube that is cut into an open lattice design and mounted onto a non-compliant balloon catheter. The Visi-Pro stent is deployed and expanded by inflating the balloon.

Indications for use

Peripheral

The Visi-Pro balloon-expandable peripheral stent is indicated for use in occlusions, lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA) or lesions believed to be at high risk of stenosis following PTA in the following arteries:

- common and external iliac, subclavian
- renal

Stenting is intended to improve and maintain artery luminal diameter.

Biliary

The Visi-Pro balloon-expandable peripheral stent is indicated for use as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications

Peripheral

The Visi-Pro stent is contraindicated for patients who share the following symptoms or factors:

- Patients for whom antiplatelet therapy, anticoagulants, or thrombolytic drugs are contraindicated.
- Patients who exhibit persistent acute intraluminal thrombus of the targeted lesion site, post thrombolytic therapy.
- Perforation at the angioplasty site evidenced by extravasation of contrast media.
- Aneurysm of the artery to be treated.
- All of the customary contraindications for PTA.

Biliary

- There are no known contraindications for the Visi-Pro stent in biliary interventions.

Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- Select a balloon inflation diameter that approximates the diameter of the lumen at the intended implant site.
- The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
- Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.
- If resistance is encountered during movement through the introducer sheath, carefully withdraw the stent system.
- Do not retract the stent into the sheath. Attempting to retract the stent into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.
- Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.
- Implanting the stent across a side branch may result in occlusion of the artery and prevent further access for future transluminal interventions.
- The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the 316L stainless steel stent should not be greater than artifacts caused by metal surgical clips. An MRI scan should not be performed until the implanted stent has had time to heal (estimated to be 8 weeks), to minimize the risk of stent migration under a strong magnetic field.

Precautions

- Carefully inspect the sterile package and the stent before use. Do not use the device if the packaging or the device is damaged.
- Pre-dilation may be required in cases of high grade stenosis to ensure that there is adequate lumen diameter for the stent to pass through the lumen.
- Always use an introducer sheath during the implant procedure to protect the vessel (peripheral use), the biliary tract (biliary use), and the puncture site. Select the introducer sheath using sizing recommendations on the label.

Potential adverse events

Potential adverse events from the stent implantation procedure include the usual complications with any percutaneous procedure, PTA, or stent placement.

The following potential adverse events are associated with peripheral interventions:

- Artery perforation or rupture
- AV fistula
- Bleeding requiring transfusion
- Cardiovascular events
- Cerebrovascular events
- Contrast media reaction/renal failure
- Death
- Embolism
- Hematoma
- Infection at the puncture site
- Intraluminal thrombus
- Pseudoaneurysm

The following complications are associated specifically with PTA or stent placement for occlusive or stenotic disease:

- Allergic response in persons with known allergic reactions to 316L stainless
- Amputation
- Dissection or intimal flap
- Failure to deploy the stent
- Restenosis
- Stent misplacement
- Stent migration
- Stent thrombosis or occlusion

The following complications are specifically associated with renal artery stent placement:

- Nephrectomy
- Renal infarction

The following potential adverse events are associated specifically with biliary interventions:

- Infection secondary to stent contamination, leading to cholangitis, hemobilia, peritonitis, or abscess
- Rupture resulting from overstretching the duct

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Recommended items

Peripheral indication

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile saline
- Appropriately sized exchange guidewire (refer to product labeling for sizing information)
- Appropriately sized hemostatic introducer sheath (refer to product labeling for sizing information)
- Appropriately sized guide catheter (refer to product labeling for sizing information)
- Renal use: appropriately sized long curved sheath or renal guide catheter (refer to product labeling)
- Inflation system

Biliary indication

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile saline
- Appropriately sized exchange guidewire (refer to product labeling for sizing information)
- Appropriately sized introducer sheath (refer to product labeling for sizing information)
- Chiba needle
- Inflation device
- Drainage catheter
- Cannula

Preparation

1. Select a stent system with a nominal balloon diameter matching the reference lumen diameter of the artery (peripheral use) or the obstructed duct and intended implant site (biliary use).

 Warning: Select a balloon inflation diameter that approximates the diameter of the lumen at the intended implant site.
2. Remove the stent system from the packaging. Visually inspect the stent to ensure that it is loaded between the proximal and distal balloon markers. Gently tug on the stent to ensure that it is crimped securely in place. If it is not securely attached, do not use the device.

 Caution: Carefully inspect the sterile package and the stent before use. Do not use the device if the packaging or the device is damaged.
3. Flush the balloon catheter lumen with saline, then load the stent delivery system onto the guidewire.
4. Connect a stopcock to the catheter inflation port.
5. Connect a syringe to the stopcock, open the stopcock, and pull the syringe plunger back to induce a negative pressure to remove air from the balloon and lumen. Close the stopcock and remove the syringe.
6. Connect an angioplasty inflation system to the stopcock. Open the stopcock and slowly fill the inflation lumen with diluted contrast medium.

Directions for peripheral stent implantation

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.

Caution: Pre-dilation may be required in cases of high grade stenosis to ensure that there is adequate lumen diameter for the stent to pass through the lumen.

1. If the iliac or renal artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery. Cross the lesion with the guidewire and advance the guide catheter or introducer sheath atraumatically over the guidewire across the lesion.

 Caution: Always use an introducer sheath during the implant procedure to protect the vessel and the puncture site. Select the introducer sheath using sizing recommendations on the label.
2. Advance the stent system through the introducer sheath or the guide catheter to the lesion site.

 Warning: If resistance is encountered during movement through the introducer sheath, carefully withdraw the stent system.

 Warning: Do not retract the stent system into the sheath. Attempting to retract the stent system into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.
3. Position the stent at the lesion site, using contrast medium to enable fluoroscopic visualization. For ostial lesions, ensure that the ostium of the lesion is covered by positioning the stent to protrude approximately 1 to 2 mm into the aorta. For iliac stenting, the entire stent should be placed within the iliac artery.
4. Under fluoroscopy, expand the stent by inflating the balloon until the stent is fully expanded to the desired diameter.

 Warning: The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
5. Deflate the balloon.
6. Gently rotate the deflated balloon delivery catheter to ensure that the stent is free and properly deployed.
7. Remove the delivery catheter.

 Warning: Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
8. Using angiography, determine whether the stent is appropriately expanded and positioned. If the stent does not cover the lesion, implant a second stent. Place the second stent overlapping the initially implanted stent with approximately 1 to 2 mm of overlap.

 Warning: Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.

Stent withdrawal

If resistance is encountered after the stent has exited the introducer sheath or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent and delivery assembly. Use the following steps for withdrawal.

1. Using fluoroscopic guidance, retract the stent to the end of the introducer sheath.
2. Inflate the delivery balloon 1 atm to reduce the likelihood of stent slippage or embolization.
3. Retract the introducer sheath and the stent system while preserving the guidewire position.
4. Withdraw the introducer sheath and the stent system assembly together as a unit.

Directions for biliary stent implantation

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.

A guidewire is first introduced through the stricture, and an internal or external biliary drain is positioned through the stricture. Next, the transhepatic tract is accessed through the routine transhepatic approach (a lateral incision is placed between the tenth and eleventh intercostal space, then punctured upward under fluoroscopy aiming at the ninth or tenth vertebral body). No predilation is performed with malignant strictures, but a large lumen catheter may be inserted. Subsequent dilation generally is performed for only a few seconds.

Biliary stent implantation

1. Select an introducer sheath using sizing recommendations on the label. Use the introducer sheath to gain access, and advance the stent system through the introducer sheath and over the guidewire to the correct position.

 Caution: Always use an introducer sheath during the implant procedure to protect both the liver tract and puncture site. Select the introducer sheath using sizing recommendations on the label.

 Warning: If resistance is encountered during movement through the sheath, carefully withdraw the stent system.

 Warning: Do not retract the stent into the sheath. Attempting to retract the stent into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.
2. Position the stent across the target stricture. The distal portion of the stent should be positioned so that it completely covers the distal end of the stricture.

 Note: Take care to position the stent so that it covers areas that are proximal and distal to an inoperable malignant tumor. Fully covering proximal and distal margins protects against obstruction from additional tumor growth.
3. Pull back on the sheath to expose the stent completely. Expand the stent by inflating the balloon catheter to the nominal inflation pressure. Use fluoroscopy to determine whether the balloon and stent are expanded appropriately. Do not exceed the recommended maximum balloon inflation pressure.

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Warning: The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.

4. Expand the stent to the lumen diameter of the bile duct.

 Warning: The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
5. Deflate the balloon.
6. Ensure that the stent is free and deployed properly by rotating the deflated balloon delivery catheter.

 Warning: Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
7. Remove the delivery catheter.

 Warning: Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.

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© 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
©2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd