39.5 x 27.25 inches (1003.3 mm x 692.15 mm)
2019/Sept/18 2:26 PM

Protégé™ RX

Self-expanding Peripheral and Carotid Stent System

CE 2797

English

Self-expanding Peripheral and Carotid Stent System

Instructions for use

1. Device description

The Protégé™ Rx self-expanding peripheral and carotid stent system is a self-expanding nitinol stent system intended for permanent implantation. The self-expanding stent is made of a nickel titanium alloy (nitinol) and comes pre-mounted on a 0.014 inch (0.36 mm) over-the-wire delivery system. It is from a nitinol tube in an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent. After deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to maintain vessel patency.

The delivery system, as shown in **Figure 1**, is comprised of an inner shaft and outer sheath, which are locked together with a tuohy-borst valve (2). The inner shaft terminates distally in a flexible pathway tip (3) and originates proximally in the tuohy-borst valve (1). The stent is positioned on the marker distal (3) and over the retainer marker proximal (4) to the constrained stent (5), are mounted on the inner shaft between the pathway tip and the retainer marker. A proximal radiopaque marker that identifies where the stent diameter transition begins.

Tapered stents are mounted onto the delivery catheter with the smallest diameter toward the distal end of the catheter. The outer sheath is retracted proximally by the tuohy-borst (6). The self-expanding stent is constrained within the space between the inner shaft and outer sheath. This space is flushed before the procedure through the side port of the Y-connector (7).

The delivery is back-loaded on a 0.36 mm (0.014 in) wire through the distal tip and exits through the Rx port (8).

Stent positioning at the target lesion is achieved before deployment by using the two radiopaque markers (3 and 10), which mark the location of the constrained stent. For stent deployment, the tuohy-borst valve (counterclockwise) to unlock the outer sheath.

The stent is deployed by retracting the outer sheath by pulling the Y-connector (6) toward the proximal grip (11).

2. Indications for use

2.1. Peripheral

The stent is indicated for use in the following types of occlusions or lesions:

- Occlusions or lesions that are at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA)
- Lesions that appear to be at high risk for restenosis following PTA, such as common iliac, femoral or subclavian arteries

Stenting is intended to improve and maintain artery luminal diameter.

2.2. Carotid

The stent is indicated for treatment of stenoses of the common carotid artery (CCA), internal carotid artery (ICA), and carotid bifurcation.

3. Contraindications

The Protégé RX self-expanding stent system is contraindicated under the following conditions:

- Patients with known hypersensitivity to nickel-titanium
- Peripheral vascular disease or anatomy, which preclude safe sheath, guide catheter or stent system insertion
- Total occlusion of the carotid artery
- Patients for whom antiplatelet therapy, anticoagulants, or thrombolytic drugs are contraindicated
- Patients with no evident persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Presence of intracranial tumor(s), arteriovenous malformations (AVMs), ipsilateral intracranial stenosis or occlusion, aneurysm, or aneurysm treatment
- Current radiation treatment for cerebral carcinoma or sarcoma presenting with occluded or subinimed vessels (carotid only)
- All the customary contraindications for PTA

4. Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

8.1. Selection and preparation

8.1.1. Recommended items

- 5 cc syringe filled with heparinized saline
- 0.014 in (0.36 mm) exchange guidewire
- Embolic protection device (carotid indication only)
- Hemostatic introducer sheath

8.1.2. Select stent size

Measure the diameter of the reference vessel (proximal and distal to lesion). Refer to **Table 2** for stent diameter sizing. Measure the length of the target lesion. Choose a stent that adequately covers the proximal and distal to the target lesion.

There is no foreshortening of the stent length during deployment.

Note: Properly size the tapered stent diameters to the lumen diameters.

Table 2. Stent diameter sizing

Stent diameter (mm)	Reference vessel diameter (mm)
6	4.5–5.5
7	5.5–6.5
8	6.5–7.5
9	7.5–8.5
10	8.5–9.5
8–6 (tapered)*	(6.5–7.5)-(4.5–5.5)
10–7 (tapered)*	(8.5–9.5)-(5.5–6.5)

*Tapered

9. Preparation of stent delivery system

a. Open the shelf box to reveal the pouch containing the stent and delivery catheter.

b. After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the outer pouch and extract the tray with contents.

c. Remove the lid from the tray. Sterile package and the device within are sterile. Do not use the device if the package or the device is damaged.

d. Set the tray on a flat surface. Carefully pull off the tray lid and remove the stent and delivery system. If sterility has been compromised or the device is damaged, do not use the device.

e. Examine the distal end of the catheter to ensure the stent is constrained within the outer sheath. Do not use the stent system if the stent is partially deployed. If a partially deployed stent is observed, the physician may open the tuohy-borst valve and gently pull the inner shaft in a proximal direction until the gap is closed. Secure the tuohy-borst valve after the procedure by turning it clockwise.

f. Flush the delivery system with heparinized saline through the hemostasis valve/port on the Y-connector. Ensure proper inject the saline. Apply constant, steady pressure to the syringe to ensure the delivery system is fully flushed.

g. The system is fully flushed when saline exits the distal end of the outer sheath and the Rx port while applying constant steady pressure.

h. Attach a 5 cc syringe filled with heparinized saline to the Y-connector side port.

i. Loosen the tuohy-borst valve (counterclockwise) and insert valve until it seeps from the proximal end of the tuohy-borst valve.

j. Close the tuohy-borst valve by turning it clockwise.

k. Continue to flush system by injecting saline until it comes out between the distal end of the outer sheath and the tip.

l. Cover the distal end of the outer sheath using fingers and continue to inject saline until it exits the Rx port.

m. Before loading the catheter onto the guidewire, submerse the distal tip in heparinized saline.

10. Stent Deployment Procedure

Warning: Careful procedures require the use of an embolic protection device. Refer to Medtronic product labeling for detailed instructions of the embolic protection device.

10.1. Insertion of introducer sheath and guidewire

a. To gain access at the appropriate site, use an introducer sheath with a hemostatic valve (2). Refer to the product labeling for the required sheath size. If the carotid or subclavian artery is used, a long sheath is recommended.

b. Insert a 0.014 in (0.36 mm) guidewire of appropriate length across the target lesion via the introducer sheath.

10.2. Dilation of lesion

Per physician discretion, pre-dilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire and the introducer sheath.

Caution: During dilation, do not expand the balloon to the point that bleeding or dissection could occur.

10.3. Introduction of stent delivery system

Advance the stent over the guidewire or the embolic protection device to the target lesion. Refer to Figure 1 for support. Use other hand to advance the delivery system through the hemostatic valve and guide catheter, or through the introducer sheath to the site of the lesion.

Warning: If resistance is encountered at any time during stent delivery, determine the cause before proceeding. Do not damage the stent or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Warning: If resistance is felt during movement through the tortuous anatomy, withdraw the stent system without deploying the stent.

Caution: Always use an introducer sheath during the implant procedure to protect both the vessel and the puncture site.

10.4. Stent deployment

a. Advance the delivery system until the radiopaque inner marker bands (distal and proximal ends of the stent) are proximal and distal to the target lesion. For tapered stents, ensure the middle radiopaque inner shaft marker is positioned at the location in which the artery begins to taper.

b. Loosen the tuohy-borst valve to unlock the system.

c. Deploy stent deployment by retracting the outer sheath (pulling the Y-connector) while holding the inner shaft (proximal grip) in a fixed position. To ensure the delivery system is straight as possible. To ensure the proximal shaft is not moved forward. Slightly more tension on the delivery catheter by using the proximal grip. Deployment is complete when the outer sheath has been retracted so that the distal stent marker and the stent is released from the retainer.

Caution: The stent is not designed for recapturing at any time and is not designed for repositioning after establishing vessel apposition.

Note: When more than one stent is required to open the lesion, place the more distal stent first. If overlap of stents is necessary, the minimal necessary, keep the amount of overlap to a minimum.

11. After stent deployment

1. While using fluoroscopy, withdraw the entire delivery system and the guidewire. If resistance is felt, withdraw the protection device, into the catheter across stents and out of the body. Remove the delivery device from the guidewire or the embolic protection device.

 Note: If resistance is met during delivery system withdrawal, advance the outer sheath so that the outer sheath marker contacts the catheter tip and withdraw the system as one unit.

2. Using fluoroscopy, visualize the stent to verify full deployment.

 Caution: Do not completely expanded at any point along the lesion, balloon dilation (standard rapid exchange PTA technique) can be performed.

 Caution: Use caution when crossing a deployed stent with any adjunct device.

 Note: Do not expand the stent beyond its nominal diameter.

 a. To dilate the stent, select an appropriately sized rapid exchange PTA balloon catheter and dilate the stent using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel.

4. Remove the guidewire and sheath from the patient.
5. Remove the guidewire and sheath from the body.
6. Close the entry wound as appropriate.
7. Discard the delivery system, guidewire, embolic protection device, and sheath.

Note: Physician use knowledge and discretion determines the appropriate drug regimen for each patient.

11.1. MRI compatibility

Nonclinical testing demonstrated that the stent is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- The maximum whole-body specific absorption rate (SAR) shall be limited to 3.0 W/kg for 20 minutes of scanning (per pulse sequence)

12. Warranty

This product is subject to standard Medtronic warranty terms. Please go to www.medtronic.com for standard terms and conditions of supply for Australia and New Zealand.

©2018 Medtronic. All rights reserved. Medtronic and the Medtronic logo are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

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1 Device description

The SpiderFX embolic protection device is designed to capture and remove dislodged debris during interventional procedures. This distal embolic protection device is delivered percutaneously. The SpiderFX capture wire is used as the 0.36 mm (0.014 in) wire for interventional device delivery. The SpiderFX device can be delivered over any 0.36 mm (0.014 in) or 0.46 mm (0.018 in) primary guidewire used to gain access to the lesion site.

The SpiderFX embolic protection device package contains the following components:

- One **capture wire** composed of a nitinol mesh filter with a distal floppy tip. The capture wire is mounted on a polytetrafluoroethylene (PTFE) coated 0.36 mm (0.014 in) stainless steel wire that is available in 190 cm or convertible 320/190 cm lengths. The filter has a heparin coating intended to maintain patency during filter deployment. This capture wire acts as the primary guidewire for other interventional devices compatible with a 0.36 mm (0.014 in) wire. The convertible 320/190 cm capture wire is scored, allowing the wire to be snapped to a 190 cm usable length for use with rapid exchange systems, if desired. The distal portion of the capture wire for rapid exchange (RX) use has a gold PTFE coating. The proximal portion for standard over-the-wire (OTW) use is black.
- One dual-ended **SpiderFX catheter** used to exchange the primary access guidewire with the capture wire, deploy the capture wire at the desired location, and recover the capture wire at the end of the procedure.
- One 23-gauge **blunt needle** for attachment to a syringe with luer connector for priming and flushing the catheter.

2 Indications for use

The SpiderFX embolic protection device provides distal embolization protection during general vascular use, including peripheral, coronary, and carotid interventions.

3 Contraindications

Use of the SpiderFX embolic protection device is contraindicated under these circumstances:

- Do not use in patients that have experienced heparin-induced thrombocytopenia.
- Do not use in patients with hypersensitivity to heparin.
- Do not deploy the filter within a previously implanted stent.
- Do not use in vessels with excessive tortuosity.

4 Warnings

- This device is sterilized using ethylene oxide and is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing increases the risk of patient infection and can compromise device performance.
- Do not grip the catheter at the distal white tip. Gripping the catheter's distal white tip can cause the embedded radiopaque marker band or inner lumen of the catheter to deform.
- Failure to adequately close off the primary guidewire exit port when flushing the delivery end of the SpiderFX catheter can result in air embolism.
- Do not undersize or oversize the filter in relation to the selected vessel diameter. Undersizing or oversizing can result in inadequate vessel wall apposition, incomplete deployment of the filter, slow-flow or no-flow, or debris embolization. (Refer to *Table 2* Sizing guide.)
- Minimize torquing of the capture wire after it has been introduced into the patient. Torquing the capture wire is likely to result in wire whipping and tangling.
- Do not apply excessive force while manipulating the capture wire. Excessive force can lead to distal embolization of debris, vessel damage, or device damage.
- Never withdraw or move an intravascular device against any resistance until the cause of the resistance is determined. Advancing against resistance can lead to debris embolization, vessel damage, or device damage.
- Do not place the SpiderFX capture wire tip in a vessel that does not allow for wire movement. Placing the capture wire in a vessel that does not allow for wire movement can lead to debris embolization, vessel damage, or device damage. Refer to *Table 2* for device sizing information.
- Minimize movement of the capture wire after initial placement. Excessive movement can cause debris embolization, vessel damage, or device damage.
- Always use the SpiderFX catheter when repositioning or removing the capture wire. Repositioning without using the SpiderFX catheter can cause debris embolization, vessel damage, or device damage.
- Observe the SpiderFX capture wire frequently under fluoroscopy, and monitor the patient to verify that the filter has not become occluded with debris, resulting in slow-flow or no-flow. Recover the filter if it becomes occluded or if the flow is compromised. See *Section 9.3* for recovery steps.
- Do not exceed 60 minutes of dwell time for the SpiderFX capture wire. Exceeding 60 minutes of dwell time can cause occlusion, resulting in slow-flow or no-flow.
- Exercise caution when using the partial recovery method to withdraw a full filter through a deployed stent. If resistance is met, advance the filter and recovery end of the catheter distally, and draw the filter more fully into the catheter before attempting withdrawal again.
- Do not use the capture wire with mechanical rotational thrombectomy, or with atherectomy or laser devices that require use of a specific guidewire.

5 Precautions

- Do not bend or reshape the distal floppy tip. Bending or reshaping can damage the device. The capture wire is not intended for use as a steerable guidewire.
- Do not use the SpiderFX embolic protection device if the sterile barriers of the package have been damaged or compromised.
- Ensure that the filter and distal floppy tip are visible within the clear segment of the catheter.
- Improperly bending the catheter at the primary guidewire exit port can damage the catheter.
- Do not allow the primary guidewire to advance proximally past the primary guidewire exit port. Contact with the distal floppy tip or the filter can damage the device.
- Clearly identify the capture wire and primary guidewire while advancing. Do not completely rotate the SpiderFX catheter. Rotating can cause the wires to wrap around the catheter.
- Position the filter sufficiently distal to the lesion so that the interventional device can be advanced without coming within 2 cm of the proximal end of the filter.

6 Potential adverse events

Potential adverse events associated with the use of the SpiderFX embolic protection device may include, but are not limited to:

- Access site adverse event (for example, arteriovenous fistula, hematoma, hemorrhage, pseudoaneurysm, or puncture site infection)
- Adjunct device entanglement
- Adverse reaction to antiplatelet/anticoagulation agents or contrast media
- Allergic reaction to device materials
- Amaurosis fugax (carotid indication only)
- Amputation (lower extremity indication only)
- Aneurysm
- Angina (coronary indication only)
- Arrhythmia
- Cardiac tamponade (coronary indication only)
- Death
- Device deformation, collapse, fracture, or rupture
- Device thrombosis (acute and subacute)
- Embolization of air, debris, plaque, or thrombus (evidenced by slow-flow or no-flow) from mechanical disruption
- Embolization or migration of the interventional device
- Emergency surgery
- Gastrointestinal bleeding due to anticoagulation
- Hemodynamic compromise (e.g., prolonged hypotension requiring treatment with intravenous medications)
- Heparin-induced thrombocytopenia
- Infection
- Intracerebral bleed (carotid indication only)
- Ischemia
- Myocardial infarction
- Renal failure requiring dialysis
- Repeat intervention to treatment site
- Seizure (carotid indication only)
- Significant cardiac arrhythmia requiring treatment with medications and/or transvenous pacing
- Thrombosis (acute and subacute)
- Transient ischemic attack (TIA) or stroke
- Vessel dissection, perforation, rupture, intimal flap, or vasospasm

7 Storage

Store the device away from sunlight at room temperature in a dry place.

8 Device components


The SpiderFX embolic protection device package contains one SpiderFX catheter, capture wire, and flushing needle. The SpiderFX capture wire and catheter are provided in a single hoop.

Note: See *Table 1* for product specifications.

Capture wire – The capture wire has a nitinol mesh filter with a gold tungsten proximal mouth indicator and a distal floppy tip, mounted on a PTFE-coated 0.36 mm (0.014 in) stainless steel wire that is available in 190 cm or convertible 320/190 cm lengths. The filter has a heparin coating intended to maintain filter patency during filter deployment. The capture wire is preloaded in the delivery end of the SpiderFX catheter with the filter extending from the catheter's distal tip.

Caution: Do not bend or reshape the distal floppy tip of the capture wire. Bending or reshaping can damage the device. The capture wire is not intended for use as a steerable guidewire.

Figure 1. Filter


1 Proximal mouth indicator
2 Filter
3 Proximal radiopaque marker band
4 Distal radiopaque marker band
5 Distal floppy tip

SpiderFX catheter – The SpiderFX catheter is dual-ended, with a delivery end (green) and a recovery end (blue) on opposite ends of the catheter.

- The **delivery end** (green) of the SpiderFX catheter has an embedded radiopaque marker at the catheter's distal tip, a primary guidewire exit port (indicated by a white stripe), a clear segment, and a capture wire exit port. The capture wire is preloaded into the delivery end of the catheter, with the filter portion extending from the catheter's distal tip.
- The center **white shaft** is used to manipulate the SpiderFX catheter while delivering and recovering the SpiderFX capture wire.
- The **recovery end** (blue) of the SpiderFX catheter has an embedded radiopaque marker at its distal tip, and a capture wire exit port.

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Figure 2. SpiderFX catheter

1 Recovery end (blue)
2 Delivery end (green)
3 Primary guidewire exit port
4 Radiopaque marker
5 Capture wire exit ports
6 Radiopaque marker

The SpiderFX embolic protection device does not contain latex.

Table 1. SpiderFX embolic protection device specifications[a]

SpiderFX capture wire					
Filter diameter (mm)	3.0	4.0	5.0	6.0	7.0
Capture wire length RX (cm)	190	190	190	190	190
Capture wire length OTW/RX (cm)	320/190	320/190	320/190	320/190	320/190
Wire diameter mm (in)	0.36 (0.014)	0.36 (0.014)	0.36 (0.014)	0.36 (0.014)	0.36 (0.014)
Distal wire tip length (cm)	1.2	1.2	1.2	1.2	1.2
SpiderFX catheter - delivery end (green)					
Nominal diameter mm (in)	1.02 (0.040)	1.02 (0.040)	1.02 (0.040)	1.02 (0.040)	1.02 (0.040)
Functional length (cm)	140	140	140	140	140
SpiderFX catheter - recovery end (blue)					
Diameter mm (Fr)	1.4 (4.2)	1.4 (4.2)	1.4 (4.2)	1.4 (4.2)	1.4 (4.2)
Functional length (cm)	140	140	140	140	140
Total catheter length (cm)	180	180	180	180	180

[a] Nominal

9 Directions for use

9.1 Preparation

1. Administer anticoagulation medications and monitor activated clotting time per standard institutional guidelines. Antiplatelet medications are recommended when using this device for intraluminal interventions.
2. Performing angiography, identify the location, approximately 3 cm distal to the lesion site within the vessel, where the SpiderFX filter will be deployed. Measure the vessel diameter at this location in millimeters, and select the appropriate filter diameter. Confirm that there is a vessel length of 4 to 5 cm distal to the lesion, to allow for adequate filter deployment.
 Warning: Do not undersize or oversize the filter in relation to the selected vessel diameter. Undersizing or oversizing can result in inadequate vessel wall apposition, incomplete deployment of the filter, slow-flow or no-flow, or debris embolization. (Refer to *Table 2* Sizing guide.)
3. Using the size specifications in *Table 2*, ensure that the guide catheter or sheath will accommodate the SpiderFX catheter. The specifications shown in *Table 2* are nominal values.

Table 2. Sizing guide


Target vessel size (mm)	Filter size[a] (mm)	Primary guidewire size mm (in)	Guide catheter or sheath min I.D.[b] mm (in)
2.0–3.0	3.0	0.36–0.46 (0.014–0.018)	1.68 (0.066)
3.1–4.0	4.0	0.36–0.46 (0.014–0.018)	1.68 (0.066)
4.1–5.0	5.0	0.36–0.46 (0.014–0.018)	1.68 (0.066)
4.5–6.0	6.0	0.36–0.46 (0.014–0.018)	1.68 (0.066)
5.5–7.0	7.0	0.36–0.46 (0.014–0.018)	1.68 (0.066)

[a] Choose the largest filter size based on the maximum vessel diameter at the deployment site.

[b] Select a guide catheter or sheath that provides coaxial alignment and stable backup support.

4. Using sterile technique, remove the hoop containing the SpiderFX device components from the packaging, and place in a sterile work area.
 Caution: Do not use the SpiderFX embolic protection device if the sterile barriers of the package have been damaged or compromised.
5. Detach the colored portion of the hoop to expose the filter and the delivery end (green) of the SpiderFX catheter.
6. Grasp the catheter near the primary guidewire exit port and gently remove the colored portion of the hoop to uncover the catheter. The remaining length of the catheter can remain in the hoop or can be removed from the hoop to expose the recovery end (blue) of the catheter.
 Note: If using the 320/190 cm convertible version of the SpiderFX device as a rapid exchange wire, firmly grasp the capture wire on each side of the score of the wire (see *Figure 3*) and snap the wire to break off the black proximal portion of the wire. Leave the black portion of the wire in the hoop or discard.


Figure 3. Capture wire

1 Score
2 Proximal (black)
3 Distal (gold)

7. Holding the catheter at its distal end, submerge the filter in heparinized saline to saturate and remove air from the filter.
 Warning: Do not grip the catheter directly at the distal white tip. Gripping the catheter's distal white tip can cause the embedded radiopaque marker bands or inner lumen of the catheter to deform.
8. Continue holding the catheter near the distal tip to provide support to the catheter shaft. The hydrophilic coating is slippery when it is wet. To avoid pinching the distal end, use gauze to hold the catheter. Pull the capture wire proximally until the filter portion is visible in the clear segment of the catheter.
 Caution: Ensure that both the filter (1) and distal floppy tip are visible within the clear segment of the catheter.

Figure 4. Filter in proper position in clear segment


9. Flush heparinized saline through the distal tip until all air is removed and fluid passes from the primary guidewire exit port (22.5 cm proximal to the distal tip of the catheter). The 23-gauge blunt needle included in the packaging can be used, if desired.
10. Gently apply pressure to block the primary guidewire exit port. Flush until all air is removed and fluid passes from the capture wire exit port (40 cm proximal to catheter distal tip).
 Warning: Failure to adequately close off the primary guidewire exit port when flushing the delivery end of the SpiderFX catheter can result in air embolism.

9.2 Delivery and deployment

1. Using standard technique to access the vessel, advance the primary guidewire at least 5 cm beyond the distal edge of the lesion being treated.
2. Ensure that the opening on the primary guidewire exit port is visible. The exit port is indicated by a white stripe on the catheter.
 Note: When loading the primary guidewire, gently bend the catheter at the primary guidewire exit port to create an approximately 2.5 cm radius, with the primary guidewire exit port at or near the apex of the curve.
 Caution: Improperly bending the catheter at the primary guidewire exit port can damage the catheter.
3. Load the distal tip of the delivery end (green) onto the proximal end of the primary guidewire. Advance the delivery end (green) until the primary guidewire exits at the exit port.
 Caution: Do not allow the primary guidewire to advance proximally past the primary guidewire exit port. Contact with the distal floppy tip or filter can damage the device.
4. Using imaging, gently advance the delivery end (green), carrying the capture wire over the primary guidewire and into the guide catheter or sheath. Continue advancement until the radiopaque marker at the distal tip of the delivery end is at least 4 to 5 cm beyond the distal edge of the lesion.
 Caution: Clearly identify the capture wire and primary guidewire while advancing. Do not completely rotate the SpiderFX catheter. Rotating can cause the wires to wrap around the catheter.
 Warning: Minimize torquing of the capture wire after it has been introduced into the patient. Torquing the capture wire is likely to result in wire whipping and tangling.
 Warning: Do not apply excessive force while manipulating the capture wire. Excessive force can lead to distal embolization of debris, vessel damage, or device damage.
 Warning: Never withdraw or move an intravascular device against any resistance until the cause of the resistance is determined. Advancing against resistance can lead to debris embolization, vessel damage, or device damage.
5. Holding the catheter stationary, withdraw the primary guidewire, leaving the delivery catheter and capture wire in place.
6. Holding the catheter stationary, gently advance the capture wire until the distal radiopaque marker band on the filter aligns with the radiopaque marker on the catheter distal tip.
 Warning: Do not place the SpiderFX capture wire in a vessel that does not allow wire movement. Placing the capture wire in a vessel that does not allow for wire movement can lead to debris embolization, vessel damage, or device damage. Refer to *Table 2* for device sizing information.
7. Using imaging, observe the filter position to ensure that the proximal radiopaque marker band is positioned at least 2 cm distal to the anticipated distal end of the interventional device that will be used to treat the lesion. Gently advance the catheter and capture wire as a unit if necessary so the proximal radiopaque marker band is at least 2 cm beyond the lesion distal edge.

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Figure 5. Positioning the SpiderFX catheter relative to interventional device.


1 Lesion
2 Anticipated distal tip of interventional device
3 SpiderFX capture wire
4 Proximal radiopaque marker band
5 SpiderFX catheter
6 2 cm

Caution: Position the filter sufficiently distal to the lesion so that the interventional device can be advanced without coming within 2 cm of the proximal end of the filter.

8. Holding the capture wire stationary, gently withdraw the catheter to expose and deploy the filter in the vessel. To reposition the filter (if required), advance the catheter back over the filter, adjust the position of the catheter and capture wire together, then redeploy the filter. Never attempt to move the filter when it is outside of the catheter.

Figure 6. Positioning the filter relative to the interventional device

1 Lesion
2 Anticipated distal tip of interventional device
3 Proximal radiopaque marker band
4 Filter
5 2 cm
6 Distal radiopaque marker band

Warning: Always use the SpiderFX catheter when repositioning or removing the capture wire. Repositioning without using the SpiderFX catheter can cause debris embolization, vessel damage, or device damage.

9. The gold radiopaque proximal mouth indicator will be visible using imaging when the filter is properly opened. If the proximal mouth indicator is not visible, attempt to redeploy the filter or to remove the filter, capture wire, primary guidewire, and the delivery catheter together. Do not reintroduce or reuse the capture wire after it has been removed from the patient.
10. Use imaging to confirm that the filter and proximal radiopaque marker band are in apposition to the vessel wall by including side view imaging as shown in *Figure 6*.
11. After the filter is deployed in the desired location, carefully remove the SpiderFX catheter from the patient. Keep the catheter in the sterile field. **Do not discard the catheter.** The recovery end (blue) of the catheter is required for capture wire recovery.
12. Place the SpiderFX catheter in sterile heparinized saline or in a sterile field, and save the catheter for recovering the SpiderFX capture wire. Do not bend or kink the catheter.
13. The 0.36 mm (0.014 in) capture wire can be used as the primary guidewire for other interventional devices. Always keep the capture wire stationary during use.

 Warning: Do not exceed 60 minutes of dwell time for the SpiderFX capture wire. Exceeding 60 minutes of dwell time can cause occlusion, resulting in slow-flow or no-flow.

 Warning: Minimize movement of the capture wire after initial placement. Excessive movement can cause debris embolization, vessel damage, or device damage.

 Warning: Do not use the capture wire with mechanical rotational thrombectomy, or with atherectomy or laser devices that require use of a specific guidewire.


 Warning: Observe the SpiderFX capture wire frequently under fluoroscopy, and monitor the patient to verify that the filter has not become occluded with debris, resulting in slow-flow or no-flow. Recover the filter if it becomes occluded or if the flow is compromised. See *Section 9.3* for recovery steps.

9.3 Recovery and removal

1. After the interventional procedure is completed, wet the distal 20 cm of the recovery end (blue) of the catheter with heparinized saline. Remove all air from the recovery end by flushing the distal tip until fluid passes from the capture wire exit port (25 cm from the tip).

 Warning: Do not grip the catheter directly at the distal white tip. Gripping the catheter's distal white tip can cause the embedded radiopaque marker bands or inner lumen of the catheter to deform.
2. Load the recovery end (blue) of the catheter onto the capture wire and advance the catheter through the guide catheter or sheath.
3. Using fluoroscopy, advance the catheter until the distal tip radiopaque marker aligns with the proximal radiopaque marker band on the filter.
4. Use one of the following methods to recover the filter:
 a. **Partial enclosure recovery:** Gently advance the recovery end of the catheter over the filter, until the proximal portion of the filter is inside the catheter (1), as shown in *Figure 7*. The recovery end of the filter is inside the catheter when the gold radiopaque proximal mouth indicator of the filter is fully compressed against the capture wire and is positioned proximally to the radiopaque marker on the recovery catheter.


Figure 7. Partial enclosure recovery

Warning: Exercise caution when using the partial recovery method to withdraw a full filter through a deployed stent. If resistance is met, advance the filter and recovery end of the catheter distally and draw the filter more fully into the catheter before attempting withdrawal again.

 b. **Full enclosure recovery:** Gently advance the recovery end over the filter until all radiopaque markers on the filter are within the catheter (1), proximal to the catheter distal marker, as shown in *Figure 8*. The distal floppy tip can remain outside the catheter.

Figure 8. Full enclosure recovery

Warning: Never withdraw or move an intravascular device against any resistance until the cause of the resistance is determined. Advancing against resistance can lead to debris embolization, vessel damage, or device damage.

5. Carefully remove the catheter and capture wire together. To allow the SpiderFX catheter to be removed without resistance, open the hemostatis valve on the guide catheter or sheath. Use care to avoid interaction with the intervention site. Discard the SpiderFX embolic protection device after removal from the patient.

 Warning: This device is sterilized using ethylene oxide and is supplied sterile for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing increases the risk of patient infection and can compromise device performance.

10 Warranty

This product is subject to standard Medtronic warranty terms. For standard terms and conditions of supply for Australia and New Zealand, go to www.medtronic.com.

M016462C001 A Medtronic Confidential
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
Composed: 2021-04-08 05:00:45
PSJ - Size-selectable package manual 14-SEP-2018

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (“TM*”) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
ev3, Inc.