Instruction for Use

GUtta Percha Points

(English Version)

Hand-Rolled Endodontic Points
Used to fill a prepared root canal prior to tooth restoration.

Composition
Gutta Percha, Zinc Oxide, Barium Sulfate, Coloring Agent

Application
1. Check the product if it is damaged, presence of impurities, and shape is suitable for use.
2. Check if the product is within the expiration date.
3. After treating the root canal of the damaged tooth, use the proper size gutta percha to fill.
4. This product is disposable and discarded after use.

Intended Purpose
1. Intended Use: Used to fill a prepared root canal prior to tooth restoration.
2. Indication: Used to fill a prepared root canal prior to tooth restoration.
3. Contraindication: Do not use for patient who is allergic to this product.
4. Patient Population: All age group
5. Intended user: Dentist
6. Side Effects/Adverse Effects: There are no side effects or adverse effects reported.
7. Potential risk if re-used: Re-use of used Gutta percha point may cause inflammation.
8. Can be used with sealer or cement and it shall be used under the control of the dentist according to the clinical indication.
9. Single Use Device

Precaution
1. Do not use other than dental specialist.
2. Do not use if the product is expired.
3. Before use, check for infection inside the root canal.
4. Do not use for patient who is allergic to this product.


Storage
Store at room temperature away from direct sunlight

Disposal : Dispose of contents/container according to the related regulations of the country.

Shelf Life : 4 years from the date of manufacture

Radiopaque : Equivalent to 6mm aluminium stepwedge

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OBELIS S.A (Authorized Representative)
Bd. Général Wahis, 53, 1030 Brussels, Belgium

■ In case that any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

■ Size Designation and Length

(Dimensions in millimeters)

Size			Tip Diameter	Length (minimum)
Taper .02	Taper .04	Taper .06		
08.02	-	-	0.08	28
10.02	-	-	0.10	
15.02	15.04	15.06	0.15	
20.02	20.04	20.06	0.20	
25.02	25.04	25.06	0.25	
30.02	30.04	30.06	0.30	
35.02	35.04	35.06	0.35	
40.02	40.04	40.06	0.40	
45.02	45.04	45.06	0.45	
50.02	50.04	50.06	0.50	
55.02	55.04	55.06	0.55	
60.02	60.04	60.06	0.60	
70.02	70.04	70.06	0.70	
80.02	80.04	80.06	0.80	
90.02	-	-	0.90	
100.02	-	-	1.00	
110.02	-	-	1.10	
120.02	-	-	1.20	
130.02	-	-	1.30	
140.02	-	-	1.40	

Size	Tip Diameter	Length (minimum)
X-Fine	0.10	28
F-Fine	0.12	
M-Fine	0.14	
Fine	0.13	
F-Medium	0.16	
Medium	0.17	
M-Large	0.15	
Large	0.23	
X-Large	0.18	
F-Large	0.20	25
XX-Fine	0.21	21
L-Medium	0.11	
Small	0.14	15
F-Small	0.17	17

F1	0.20	28
F2	0.25	
F3	0.30	
F4	0.40	
F5	0.50	

※ Please contact us for other special sizes.

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Instruction for Use

Gutta Percha Points(Obturator)

Used to fill a prepared root canal prior to tooth restoration.

Composition
Gutta Percha, Zinc Oxide, Zirconium Oxide, Barium Sulfate, Coloring Agent

Application
1. Check the product if it is damaged, presence of impurities, and shape is suitable for use.
2. Check if the product is within the expiration date.
3. After treating the root canal of the damaged tooth, fill the root canal with the product by using obturation gun.
4. This product is disposable and discarded after use.

Intended Purpose
1. Intended Use: Used to fill a prepared root canal prior to tooth restoration.
2. Indication: Used to fill a prepared root canal prior to tooth restoration.
3. Contraindication:
 Do not use for patient who is allergic to this product.
4. Patient Population: All age group
5. Intended user: Dentist
6. Side Effects/Adverse Effects: There are no side effects or adverse effects reported.
7. Potential risk if re-used: Re-use of used Gutta percha obturator may cause inflammation.
8. Can be used with sealer or cement and it shall be used under the control of the dentist according to the clinical indication.
9. Single Use Device

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
■ **In case that any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

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Instruction For Use

Absorbent Paper Points

Sterilized endodontic points
Used to remove moisture from the root canal or to introduce medicaments.

Composition
Paper 100%

Application
1. Check if there is any abnormality in the appearance of the product such as foreign substances or cracks.
2. Check if the product is within the expiration date.
3. After forming the root canal, select the size that fits the file and dry the root canal, and use it with medication after.
4. This product is a gamma ray sterilized disposable medical device.

Intended Purpose
1. Intended Use: Used to remove moisture from the root canal or to introduce medicaments.
2. Indication: Used to remove moisture from the root canal or to introduce medicaments.
3. Contraindication: No contraindication has been reported that can cause health problems.
4. Patient Population: All age group
5. Intended user: Dentist
6. Side Effects/Adverse Effects: There are no side effects or adverse effects reported.
7. Potential risk if re-used: Re-use of used Paper point may cause inflammation.
8. Can be used with medicaments and tools and it shall be used under the control of the dentist according to the clinical indication.
9. Single Use Device

Precaution
1. Do not use other than dental specialist.
2. Use disinfected tools.
3. This product is a gamma-ray sterilized disposable medical device.

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