Box Size: 15.5 x 5.25
Overall Size: 27 x 26
Folds to: 3.38 x 6.5

COVIDIEN™

Endo GIA™

Reinforced Reload with Tri-Staple™ Technology

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IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or sterilization of this device may create the risk of contamination and patient infections. Do not reuse, reprocess or resterilize this device.

DESCRIPTION

The Endo GIA™ reinforced reload with Tri-Staple™ technology is pre-attached with a polyglycolic acid (PGA) staple line reinforcement. The reload includes a non-woven absorbable mesh prepared of polyglycolic acid and dyed with D&C Green No. 6 on both the anvil and cartridge. The dissection fibres of tissue that can be compressed by 0.5 mm. The PGA reinforcement material design was a combination of hydrolytic and enzymatic pathways. The PGA staple line reinforcement material has been found to be biocompatible and non-antigenic, with a minimal tissue reaction or sex-dependant variation and other implantation devices. The strength profile for the reinforcement material shows mechanical strength through the critical wound healing period of 10 to 14 days post-surgery. The absorbable strength for the staple line reinforcement material is 52% at 7 days and 14% at 14 days and the reinforcement material is fully absorbed in 35 weeks.

The PGA staple line reinforcement material is secured to the anvil and cartridge of the Endo GIA™ reinforced reload with Tri-Staple™ technology and a Maxon™ synthetic absorbable anchoring suture. During firing, the Maxon™ anchoring sutures release the PGA reinforcement material from the anvil to allow its placement of the PGA reinforcement material on the tissue. The Maxon™ suture remains with the reload and is discarded after use.

The Endo GIA™ reinforced reload with Tri-Staple™ technology stays in a triple-staggered row of titanium staples and two layers of absorbable reinforcement material on either side of the cut line. As the staples are deployed, the tissue and staple line reinforcement material are simultaneously divided. The tissue clamp on a stapler is determined by the selection of the reload. The Endo GIA™ reinforced reload with Tri-Staple™ Technology is available in articulating 45 mm and 60 mm lengths in purple and black cartridges. The purple reload has height progressive over 4.0 mm, 4.5 mm and 4.8 mm staples in staples on either side of the cut line. The black reload has height progressive over 4.0 mm, 4.5 mm and 5.0 mm titanium staples on either side of the cut line.

The Endo GIA™ reinforced reload with Tri-Staple™ Technology can be used with the GIA™ universal staplers, Endo GIA™ universal staplers, Endo GIA™ Ultra universal staplers and the Endo GIA™ Ultra or a purpose-directype Drive™ Ultra with adaptative Drive Gold™ adapter.

Refer to the instructions for use for the Endo GIA™ universal stapler, Endo GIA™ universal stapler, Endo GIA™ Ultra universal stapler and the Drive™ Ultra powered handle with Endo GIA™ adapter for specific indications, contraindications, warnings, precautions, and operating instructions.

INDICATIONS

Endo GIA™ reinforced reload with Tri-Staple™ Technology preloaded with polyglycolic acid staple line reinforcement has applications in abdominal, gynecologic, pediatric, and thoracic surgery for resection, transection of tissue and creation of anastomoses. They may be used for transection and resection of liver substance, hepatic vasculature and biliary structures, and for transection and resection of pancreas.

CONTRAINDICATIONS

Do not use for the patch reconstruction of cardiovascular defects such as cardiac, great vessel and peripheral vascular anomies or vein.

Do allow for the reinforcement material thickness of 0.300 mm or unload, the tissue thickness ranges are as follows:

For the purple reload: Do not use on tissue that does not apply compress to 1.25 mm–2.25 mm.

For the black reload: Do not use on tissue that does not apply compress to 1.95 mm–2.7 mm.

WARNINGS

1. Use of this product or applications other than those indicated has the potential for serious complications, such as inadequate reinforcement strength, staple pullout, infection, abscesses, migration and tension.
2. When positioning the staples on the appropriate tissue, do not clamp tissue such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
3. Do not use the Endo GIA™ reinforced reload with Tri-Staple™ Technology in conjunction with other reinforcement materials. This may cause staple malfunction and poor hemostasis resulting in leaks and tissue ischemia.
4. In the event there is partial or no tissue compression, the material can elect a slight acute angle and edge, remove this material with scissors.
5. Do not attempt to remove the yellow shipping wedge until the Endo GIA™ reinforced reload with Tri-Staple™ Technology is loaded into the instrument. Failure to do so may result in loading and jaws not opening after clamping and firing on tissue. The shipping wedge should be retained until the completion of the surgical procedure. Please refer to the applicable Endo GIA™ universal stapler and Endo GIA™ Ultra universal stapler instructions for use for detailed loading instructions.
6. Do not clamp instrument prior to removing shipping wedge.

ADVERSE REACTIONS

Possible adverse reactions may include, but are not limited to: infection, inflammation, adhesions, and hemostasis.

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
PT00004274

EN

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques. This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or resterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize this device.

DESCRIPTION

The Tri-Staple™ 2.0 reinforced reload is pre-attached with a polyglycolic acid (PGA) staple line reinforcement. The reload includes a non-woven absorbable mesh prepared of polyglycolic acid and dyed with D&C Green No. 6 on both the anvil and cartridge. This decreases the thickness of tissue that can be compressed by 0.3 mm. The PGA reinforcement material degrades via a combination of hydrolytic and enzymatic pathways. The PGA staple line reinforcement material has been found to be biocompatible and non-antigenic, with a history of use as bio-absorbable sutures and other implantable devices. The strength profile for the reinforcement material shows mechanical strength through the critical wound healing period of 10 to 14 days post-surgery. The residual measurable strength for the reinforcement material is 32% at 7 days and 1.4% at 14 days and the reinforcement material is bio-absorbed in 15 weeks.

The PGA staple line reinforcement material is secured to the anvil and cartridge of the Tri-Staple™ 2.0 reinforced reload and a Maxon™ synthetic absorbable anchoring system. During firing, the Maxon™ anchoring sutures release the PGA reinforcement material from the reload to allow its placement of the PGA reinforcement material on the tissue. The Maxon™ suture remains with the reload and is discarded after use.

The Tri-Staple™ 2.0 reinforced reload places a triple-staggered row of titanium staples and two layers of absorbable reinforcement material on either side of the cut line. As the staples are deployed, the tissue and staple line reinforcement material are simultaneously divided. The size of the staples is determined by the selection of the reload. The Tri-Staple™ 2.0 reinforced reload is available in articulating 45 mm and 60 mm lengths in purple and black configurations.

The purple reload has height progressive rows of 3.0 mm, 3.5 mm and 4.0 mm titanium staples on either side of the cut line. The black reload has height progressive rows of 4.0 mm, 4.5 mm and 5.0 mm titanium staples on either side of the cut line.

The Tri-Staple™ 2.0 reinforced reload is provided STERILE for single use only. The Tri-Staple™ 2.0 reinforced reload is sterilized by ethylene oxide gas. The product should be stored at room temperature.

The Tri-Staple™ 2.0 reinforced reload will contain an intelligence chip. The intelligence chip will have the ability to communicate with Covidien™ powered stapling handles that have a compatible communications interface.

STAPLER HANDLE COMPATIBILITY

The Tri-Staple™ 2.0 reinforced reload can be used with the GIA™ universal stapler, Endo GIA™ universal staplers, Endo GIA™ Ultra universal staplers and Covidien™ powered stapler handles used with Endo GIA™ adapters.

Refer to the instructions for use for the GIA™ universal stapler, Endo GIA™ universal stapler, Endo GIA™ Ultra universal stapler and Covidien™ powered stapler handles used with Endo GIA™ adapters for specific indications, contraindications, warnings, precautions, and operating instructions.

INDICATIONS

Tri-Staple™ 2.0 reinforced reload preloaded with polyglycolic acid staple line reinforcement has applications in abdominal, gynecologic, pediatric and thoracic surgery for resection, transection of tissue and creation of anastomosis. They may be used for transection and resection of liver substance, hepatic vasculature and biliary structures, and for transection and resection of pancreas.

CONTRAINDICATIONS

1. Do not use for the patch reconstruction of cardiovascular defects such as cardiac, great vessel and peripheral vascular arteries or vein.

2. To allow for the reinforcement material total thickness of 0.300 mm on a reload, the tissue thickness ranges are as follows:

For the purple reload: Do not use on tissue that does not easily compress to 1.2 mm-1.95 mm.

For the black reload: Do not use on tissue that does not easily compress to 1.95 mm-2.7 mm.

WARNINGS

1. Use of this product in applications other than those indicated has the potential for serious complications, such as inadequate reinforcement strength, staple pullout, infection, abrasion, migration and erosion.

2. When positioning the stapler on the application site, ensure that no obstructions, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

3. Do not use the Tri-Staple™ 2.0 reinforced reload in conjunction with other reinforcement materials. This may cause staple malfunction and poor hemostasis resulting in leaks and tissue ischemia.

4. In the event there is extra material with or without staples protruding beyond the tissue and cut edge, remove this material with scissors.

5. Do not attempt to remove the yellow shipping wedge until the Tri-Staple™ 2.0 reinforced reload is loaded into the instrument. Failure to do so may result in loading unit jaws not opening after clamping and firing on tissue. The shipping wedge should be retained until the completion of the case. Please refer to the applicable Endo GIA™ universal stapler, Endo GIA™ Ultra universal stapler and Covidien™ powered stapler handles used with Endo GIA™ adapters instructions for use for detailed loading instructions.

6. Do not clamp instrument prior to removing shipping wedge.

ADVERSE REACTIONS

Possible adverse reactions and potential complications may include but are not limited to seroma/hematoma, bleeding/hemorrhage/leak, ischemia, chronic pain, infection, allergic reactions, inflammatory reaction, visceral adhesions, fistula, nerve entrapment, tissue erosion and small bowel perforation, need for future reoperation.

PRECAUTIONS

1. If the Tri-Staple™ 2.0 reinforced reload is inserted precutaneously, ensure the incision is wide enough to accommodate the instrument so that excessive force is not placed on the jaws of the instrument.

2. When using the Endo GIA™ Ultra universal stapler to grasp or manipulate tissue multiple times, the staple line reinforcement material secured on the anvil and cartridge may move and/or dislodge.

3. In the event that the reinforcement material shifts upon clamping, unclamp the stapler and reposition on tissue.

4. When using the Tri-Staple™ 2.0 reinforced reload with a Covidien™ powered stapler handle with Endo GIA™ adapter in challenging applications, the firing sequence may stop prematurely. The propensity to stop prematurely may increase if the reload is in the articulated position. If the firing does stop prematurely, release the blue CLOSE/FIRE button, assess the tissue for obstructions, and wait approximately 15 – 30 seconds to allow the tissue to compress. After waiting 15-30 seconds, press the CLOSE/FIRE button to continue the firing stroke. The knife blade will stop automatically at the distal end of the reload when the firing stroke is complete. A complete firing stroke can be confirmed by visual indicators on the powered stapler handle and the reload. In the event that the firing cannot be completed, the knife blade can be retracted and the jaws opened by pressing the silver OPEN button. Refer to Precaution #5 for instructions to release the distal Maxon™ anchoring sutures in the event of an incomplete firing.

5. In the event of an incomplete firing of the Tri-Staple™ 2.0 reinforced reload, the distal Maxon™ anchoring sutures used to secure the staple line reinforcement material to the anvil and cartridge may not be released, the inability to release the distal anchor suture will prevent the jaws from opening fully when the device is unclamped. The remaining intact staple line reinforcement material on the reload will need to be removed from the device by gentle traction and/or by cutting the anchor suture with scissors.

6. After complete firing, inspect the staple line for hemostasis; minor bleeding may be controlled by electrocautery or manual sutures. Additionally, small pieces of the Maxon™ sutures that secured the Tri-Staple™ 2.0 reinforced reload with the PGA staple line reinforcement material may be dislodged and observed in the surgical field.

7. Securing the target tissue prior to and during firing with an atraumatic grasper may reduce tissue slipping during the firing stroke. When closing the Tri-Staple™ 2.0 reinforced reload, the tissue may extrude distally and laterally as you clamp and/or fire when in thick tissue and/ or crossing intersecting staple lines. To minimize this, the use of an atraumatic grasper is recommended to place and secure the tissue within the jaws. As the instrument is fired, use gentle proximal traction with the grasper placed next to the reload until the deployed staples provide acceptable support.

8. As with any temporary implanted prosthetic, the material can elicit a slight acute inflammatory response in tissues, and its implantation may result in inflammation or partial extrusion of material with late absorption or mechanical irritation and may result in a pH decrease by accumulation of dissolved unmetabolized degraders and inflammation caused from the pH decrease.

9. As with any foreign body, prolonged contact of this or any other reinforcement material with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

10. To avoid damage or contamination, always use clean gloves and grease care is taken when handling the Tri-Staple™ 2.0 reinforced reload, so as not to damage the staple line reinforcement material.

11. Dispose of used instruments and used reloads in accordance with the end-user's medical and biological waste disposal requirements.

12. Store at room temperature. Avoid prolonged exposure to elevated temperatures.

13. If the desiccant packaged in the foil pouch has been compromised, discard the Tri-Staple™ 2.0 reinforced reload.

14. Before attempting to load a previously loaded but unused reload, reattach the shipping wedge ensuring that the raised post on the shipping wedge is properly seated in the anvil of the reload. Then follow the loading procedure as described in the instructions for use for the GIA™ universal stapler, Endo GIA™ universal stapler, Endo GIA™ Ultra universal stapler, and iDrive™ Ultra powered handle with Endo GIA™ adapter.

15. The titanium staples in this device are nonferromagnetic and may undergo MRI and/or NMR procedures at forces of 3.0 Tesla or less after implantation.

1 SCHEMATIC VIEW

A) SINGLE USE Tri-Staple™ 2.0 reinforced reload (45 or 60mm lengths)

B) ANVIL

C) STAPLE CARTRIDGE

D) SHIPPING WEDGE

E) INTELLIGENCE CHIP

F) STAPLE LINE REINFORCEMENT

NOTE: For complete staple instructions, always refer to the package insert included with the Endo GIA™ Ultra universal, Endo GIA™ universal, GIA™ universal single use stapler or Covidien™ powered stapler handles used with Covidien™ Endo GIA™ adapters for specific instructions, INDICATIONS, CONTRAINDICATIONS, WARNINGS, and PRECAUTIONS.

INSTRUCTIONS FOR USE

LOADING

1. The reload is packaged in the open position. Do not attempt to close the reload.

WARNING: Always select a reload with the appropriate staple size for the tissue thickness. Overly thick or thin tissue may result in unacceptable staple formation.

CAUTION: Do not attempt to remove the shipping wedge until the reload is loaded into the instrument. The shipping wedge should be retained until the completion of the case.

CAUTION: Do not clamp instrument prior to removing shipping wedge.

2. Depending upon the stapler handle used, refer to the applicable instructions for use for the GIA™ universal stapler, Endo GIA™ universal staplers, Endo GIA™ Ultra universal staplers and Covidien™ powered stapler handles used with Endo GIA™ adapters for specific indications, contraindications, warnings, precautions, and operating instructions.

NOTE: Dipping the reload in sterile saline prior to positioning the reload at the target tissue may help flatten the material on the cartridge.

UNLOADING

1. Dependent upon the stapler handle used, refer to the applicable instructions for use for the GIA™ universal stapler, Endo GIA™ universal staplers, Endo GIA™ Ultra universal staplers and Covidien™ powered stapler handles used with Endo GIA™ adapters for specific indications, contraindications, warnings, precautions, and operating instructions.

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