Instructions

On the use of cast glass bottles produced by SGD S.A., France.

Vials are made of clear and amber glass 1, 2, 3 hydrolytic types, with a capacity of 5.0-3000.0 ml. Vials made of molded glass according to international ISO standards. Vials are used for filling liquid, bulk and lyophilic medicines.

Each batch (series) of bottles is accompanied by a quality certificate of the manufacturer and corresponds to the quality specified in the certificate.

Package:

The vials are packed in a shrink film in the form of rectangular "seals", the sizes of which are standard for each type and capacity of vials according to material packaging descriptions. Euro-pallets with a certain amount of each type and capacity of bottles are formed from the "seals", which are marked with the standard label of the manufacturer. The label indicates the name of the plant, product name, batch number, quantity, date of manufacture.

Preparation of vials for use:

Before use, the vials undergo quality control in the QCD laboratory, after which the QCD laboratory gives permission for use. Before loading the technological line, the vials are visually inspected, the vials with identified quality deviations are rejected. Then the vials are placed in special metal cassettes, which are loaded into the technological line.

Note: After loading the cassettes into the technological line, the preparation of bottles, filling, capping is carried out automatically, in accordance with the preset mode.

Brief description of the technological process:

The vials, placed in cassettes, are washed with circulating water heated to 40-60 ° C, and blown through with compressed purified air. Then they are rinsed with purified water, followed by blowing with compressed purified air, then rinsed with water for injection, blown and dried with compressed purified air.

Sterilization of vials begins with warming, high-temperature frying in a sterilization tunnel and ends with cooling.

The sterilized vials are automatically fed to the dosing units for filling and then to the capping unit, where the bottle is closed with a stopper and rolled up with an aluminum (alumoplastic) cap. Vials with lyophilic preparations are sent for pre-closure with special stoppers for lyophilic preparations, then for freeze drying, then for sealing and seaming with an aluminum (alumoplastic) cap.

Correctly sealed vials with a medical product are sent to quarantine for quality control. After confirming the quality by the laboratory, the vials with the medical product are sent to the labeling and packaging line in a pack.

Please note that SGD Pharma is only responsible for the items listed in the introduction and "Package" paragraphs of this letter. It is under the responsibility of our customers to prepare the vials for use and to choose the technological process according to their final needs.

Laurent ZUBER – CCO
On the 31st of March, 2022 – Puteaux, France

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