PRODUCT: STERRAD NX (ASP GmbH Version)
PRODUCT DESCRIPTION: Cassette Outsert
PRODUCT CODE: 10133
LANGUAGES: English, French, German, Italian, Spanish, Dutch, Portuguese, Danish, Finnish, Norwegian, Swedish, Greek, Japanese, Polish, Hungarian, Slovak, Czech, Slovenian, Croatian, Turkish, Russian, Romanian, Estonian, Latvian, Lithuanian, Bulgarian

SIZE: 6 x 7 panels
COLORS: Black, PMS 185, PMS 109, PMS 226, PMS 2695 C
DATE: December 22, 2020
PART NUMBER: 101335-50
BARCODE: 2D datamatrix

Pantone Black | Pantone 185 | Pantone 109 | Pantone 226 | Pantone 2695 C

ASP
REF 10133
STERRAD NX™ Cassette
① STERRAD NX® Cassette
② STERRAD NX® Wirkstoff-Kassette
③ STERRAD NX® Cassette
④ STERRAD NX® Cassette
⑤ STERRAD NX® Cassette
⑥ STERRAD NX® Cassette
⑦ STERRAD NX® Cassette
⑧ STERRAD NX® Cassette
⑨ STERRAD NX® Cassette
⑩ STERRAD NX® Cassette
⑪ STERRAD NX® Cassette
⑫ STERRAD NX® Cassette
⑬ STERRAD NX® Cassette
⑭ STERRAD NX® Cassette
⑮ STERRAD NX® Cassette
⑯ STERRAD NX® Cassette
⑰ STERRAD NX® Cassette
⑱ STERRAD NX® Cassette
ADVANCED STERILIZATION PRODUCTS, INC.
33 Technology • Irvine, CA 92618, USA
Produced by:
ASP Global Manufacturing GmbH,
Switzerland under license from:
Advanced Sterilization Products, Inc., USA
949-581-5799
CE
0123
ASP, The Netherlands BV
Antonio Vivaldistraat 150
Eindhoven, The Netherlands
© ASP 2021
2020-12
101335-50

1

en

INSTRUCTIONS FOR USE
STERRAD NX® Cassette

Fits (S) Cycles/Cassette

Your Personal Protective Equipment (e.g. Chemical Resistant Gloves) at all times when handling cassettes.

PROCEDURE:
- Remove plastic wrap.
- Remove cassette from box.
- Insert cassette into STERRAD NX Sterilizer.
- The cassette will be accepted automatically.

WARNING: Each cell contains less than 1750 µl (1,800 µl, 450 µl of 58% hydrogen peroxide which is a strong oxidizer, irritant, and corrosive. If the yellow indicator has turned red, the cassette might have been damaged and do not remove the plastic wrap. Contact ASP Customer Service for instructions.

Refer to STERRAD NX Sterilizer Operator's Manual for additional information and warnings.

For Use in STERRAD NX Sterilizer Only.

Store at room temperature between 15°-30°C (59°- 86°F).

This product is intended for use in the processing of critical and semi-critical medical devices that are to be sterilized.

For a glossary of symbols used, go to: qr.asp.com.

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PRODUCT: STERRAD 100NX (ASP GmbH Version)
PRODUCT DESCRIPTION: Cassette Outsert
PRODUCT CODE: 10144
LANGUAGES: English, French, German, Italian, Spanish, Dutch, Portuguese, Danish, Finnish, Swedish, Estonian, Greek, Japanese, Polish, Hungarian, Slovakian, Czech, Slovenian, Croatian, Turkish, Russian, Romanian, Latvian, Lithuanian, Taiwan, Korean, Vietnamese, Bulgarian

SIZE: 6 x 7 panels
COLORS: Black, Pantone 185, Pantone 176, Pantone 209, Pantone 2239C
DATE: December 22, 2020
PART NUMBER: 10812-50
BARCODE: 2D datamatrix

Pantone Black | Pantone 185 | Pantone 176 | Pantone 209 | Pantone 2239 C

ASP
REF 10144
STERRAD® 100NX Cassette

- STERRAD® 100NX Cassette
- STERRAD® 100NX Kassette
- STERRAD® 100NX Cassette
- Cassette STERRAD® 100NX
- STERRAD® 100NX Cassetta
- STERRAD® 100NX Kassett
- STERRAD® 100NX Kasetti
- Kontio STERRAD® 100NX
- Κασέτα STERRAD® 100NX
- Kaseta STERRAD® 100NX
- Kazeta STERRAD® 100NX
- Kaseta STERRAD® 100NX
- STERRAD® 100NX kasete
- STERRAD® 100NX
- Kasszetta STERRAD® 100NX
- Cassett STERRAD® 100NX
- Kazeta STERRAD® 100NX
- STERRAD® 100NX 盒式磁带
- STERRAD® 100NX 盒式磁帶
- STERRAD® 100NX Bång cassette
- Kaseta STERRAD® 100NX

ADVANCED STERILIZATION PRODUCTS, INC
33 Technology Drive, Irvine, CA 92618, USA
or
ASP Global Manufacturing GmbH,
Switzerland and Ireland Branch
Advanced Sterilization Products, Inc., USA
CE 0086
ASP The Netherlands BV
EC REP BIO 1, 5657 DB,
Eindhoven, The Netherlands

© ASP 2021 | 2020-12 | 10812-50

EN

Caution: Rx only

Intended Use
Sterilization of heat and moisture sensitive medical devices.

Warnings
Do not use if package is damaged.
Do not use if sterilization indicator has changed color.
Do not use if sterilization indicator is not present.
Do not use if expiration date has passed.
Do not use if cassette is damaged.
Do not use if cassette is wet.
Do not use if cassette is not fully inserted.
Do not use if cassette is not locked.
Do not use if cassette is not unlocked.
Do not use if cassette is not removed.
Do not use if cassette is not removed after the cycle.
Do not use if cassette is not removed after the cycle is aborted.
Do not use if cassette is not removed after the cycle is completed.
Do not use if cassette is not removed after the cycle is interrupted.
Do not use if cassette is not removed after the cycle is paused.
Do not use if cassette is not removed after the cycle is resumed.
Do not use if cassette is not removed after the cycle is restarted.
Do not use if cassette is not removed after the cycle is stopped.
Do not use if cassette is not removed after the cycle is terminated.
Do not use if cassette is not removed after the cycle is ended.
Do not use if cassette is not removed after the cycle is finished

PRODUCT: STERRAD 100S [ASP GmbH Version]
PRODUCT DESCRIPTION: Cassette Outsert
PRODUCT CODE: 10113
LANGUAGES: English, French, German, Italian, Spanish, Dutch, Portuguese, Danish, Finnish, Norwegian, Swedish, Greek, Japanese, Polish, Hungarian, Slovak, Czech, Slovenian, Croatian, Turkish, Russian, Romanian, Estonian, Latvian, Lithuanian, Bulgarian

SIZE: 6 x 7 panels
COLORS: Black, PMS 185, PMS 109, PMS 226
DATE: April 20, 2021
PART NUMBER: 108433-50
BARCODE: 2D datamatrix

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