Amplatzer™ Talisman™ Delivery Sheath

Instructions for Use

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Device Description and Purpose

The Amplatzer™ Talisman™ Delivery Sheath is designed to deliver an Amplatzer™ Talisman™ PFO Occluder. The Amplatzer™ Talisman™ Delivery Sheath includes the sheath and dilator.

The body of the sheath is radiopaque for visibility under fluoroscopy.

Refer to the figure and table in Appendix A: Supplemental Information for details about the delivery sheath. The following components are identified in Figure 1:

A. Dilator – Eases penetration through tissue and minimizes vessel trauma

B. Sheath – Provides a pathway through which the Amplatzer™ Talisman™ PFO Occluder is delivered

Intended Use

The Amplatzer™ Talisman™ Delivery Sheath is intended to facilitate the delivery and deployment of an Amplatzer™ PFO Occluder.

Indication for Use

The Amplatzer™ Talisman™ Delivery Sheath is indicated for patients recommended for patent foramen ovale closure with an Amplatzer™ PFO Occluder.

Contraindications

None known.

Warnings

- This device was sterilized with ethylene oxide and is for single use only. Do not reuse or resterilize this device. Attempts to resterilize this device can cause a malfunction, insufficient sterilization, or harm to the patient.
- Do not use this device if the sterile package is open or damaged. Inspect all components before use. Do not use if the package or items appear to be damaged or defective.
- DO NOT use the Amplatzer™ Talisman™ Delivery Sheath after the Use-by date stated on the package label.
- This device should be used only by physicians who are trained in standard transcatheter techniques. The physician should determine which patients are candidates for procedures that use this device.
- Use a hemostasis valve to impede blood backflow during the implant procedure.

Precautions

- The physician should exercise clinical judgment in situations that involve the use of antithrombotic drugs before, during, and/or after the use of the delivery sheath.
- Use caution and rely on imaging guidance when advancing the sheath and dilator to minimize the risk of cardiovascular injury or interference with previously implanted medical devices. Whenever possible, advance the sheath and dilator over a guidewire.
- Do not attempt to use a guidewire larger than the maximum diameter specified in the Directions for Use.
- Do not use a power injection system to put contrast solution through the sheath.
- Remove the dilator and sheath from the patient slowly to prevent an ingress of air.

Limitations

Certain conditions may require special consideration when using this product. These may be, but are not limited to, the presence of a membrane in the right atrium such as Chiari network (CN), Eustachian Valve (EV), or cases with previously implanted medical devices (e.g., transvenous pacing leads).

Potential Adverse Events

Potential adverse events that may occur during or after a procedure using this delivery sheath may include, but are not limited to:

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- Allergic reaction/toxic effects due to hypersensitivity to contrast agent, anesthesia, device materials, or drugs used to minimize blood clot formation
- Arrhythmia
- Arteriovenous fistulae
- Bleeding
- Cardiac perforation
- Cardiac tamponade
- Cardiovascular injury
- Death
- Dissection
- Embolism (air, foreign body, and peripheral)
- Hematoma
- Infection
- Myocardial infarction
- Pericardial effusion
- Stroke
- Thromboembolism
- Thrombosis
- Vascular access site injury

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Performance Characteristics

The Amplatzer™ Talisman™ Delivery Sheath consists of the following:

- Sheath
- Dilator

With a working length of 80 cm and an inner diameter of 8 or 9 French, the sheath, via femoral vein access, provides flexible positioning within the cardiac anatomy for placement of the Amplatzer™ Talisman™ PFO Occluder.

A dilator facilitates sheath placement. When positioned within the sheath, the dilator provides an inner lumen compatible with .035 guidewires and a low-profile tip, which allows the sheath to be advanced over a guidewire until it reaches the desired cardiac structure.

Device Compatibility

The Amplatzer™ Talisman™ Delivery Sheath is compatible with the Amplatzer™ Talisman™ PFO Occluders shown in Table 1 in Appendix A: Supplemental Information.

Storage and Handling Information

- There are no special storage requirements beyond what is typical for medical devices.
- The instructions for use are recyclable. Dispose of all packaging materials appropriately.

Directions for Use

Materials recommended for use with the Amplatzer™ Talisman™ Delivery Sheath:

0.035-inch Amplatzer™ Guidewire (9-GW-002)

Procedure

CAUTION: When placing a device using an Amplatzer™ Talisman™ Delivery Sheath, refer to the instructions for use provided with the device.

General instructions for the Amplatzer™ Talisman™ Delivery Sheath are provided below.

1. Select the appropriate delivery sheath for use with the device being delivered. Refer to Appendix A: Supplemental Information for sizing information.
2. Prepare the delivery sheath for use:
 - Inspect the sterile pouch and verify that it is unopened and undamaged. Do not use the components if the sterile barrier has been compromised.
 - Gently open the sterile pouch and inspect the components for damage. Do not use damaged or kinked components.
 - Flush all components with sterile saline.
 - Wipe the dilator and sheath with sterile gauze moistened with sterile saline to remove any foreign material.
3. Access the desired vessel.

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4. Place the guidewire according to the Amplatzer™ Talisman™ PFO Occluder's instructions for use.
5. Insert the dilator into the sheath. You may encounter resistance as the dilator reaches the distal end of the sheath because the tip of the sheath is tapered.
6. Turn the rotating luer on the dilator clockwise to lock the components together.
7. Advance the dilator and sheath over the guidewire until the sheath is positioned according to the Amplatzer™ Talisman™ PFO Occluder's instructions for use.
8. Turn the rotating luer on the dilator counterclockwise to unlock the dilator. Slowly remove the dilator and guidewire from the sheath.

 WARNING: Remove the dilator and guidewire slowly to prevent an ingress of air.
9. Allow blood backflow to purge any air from the sheath.
10. Connect the sheath to the loader according to the Amplatzer™ Talisman™ PFO Occluder's instructions for use.
11. Position, deploy, and detach the device according to the Amplatzer™ Talisman™ PFO Occluder's instructions for use.
12. When the procedure is complete, slowly remove the sheath.

 WARNING: Remove the sheath slowly to prevent an ingress of air.

Disposal

After use, this device, its accessories and packaging should be appropriately classified for disposal, e.g., biohazard, sharps, non-hazardous waste etcetera, and carefully disposed of in compliance with facility procedures and applicable laws and regulations.

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Amplatzer™ Talisman™ PFO Occluder


EN: Instructions for Use

CAUTION: Read all instructions carefully. Failure to follow these instructions, warnings, and precautions may lead to patient injury or device damage.

Device Description and Purpose

The Amplatzer™ Talisman™ PFO Occluder is a self-expandable, double-disc device made from a Nitinol wire mesh. The two discs are linked together by a short connecting waist that allows free motion of each disc. The occluder has radiopaque marker bands on its distal and proximal ends, and also contains an end screw on the proximal end to facilitate delivery and deployment. The device is sterilized with ethylene oxide.

Figure 1. Amplatzer™ Talisman™ PFO Occluder

1	Right atrial disc
2	Left atrial disc
3	Connecting waist

Appendix A: Supplemental Information contains additional information about the occluder, delivery sheath, and dilator:

- Table 1 shows the device and delivery sheath dimensions.
- Figure 1 shows the Amplatzer™ Talisman™ PFO Occluder pre-attached to the Amplatzer™ Trevisio™ delivery cable. The illustration has call-outs for the following components:

1	Delivery cable
2	Hemostasis valve
3	Extension tube with three-way stopcock
4	Plastic vise
5	PFO occluder
6	Distal cable attached to the occluder
7	Loader

- Figure 2 shows intracardiac echocardiography (ICE) images.
- Figure 3 shows transesophageal echocardiography (TEE) images.

Intended Use

The Amplatzer™ Talisman™ PFO Occluder is a percutaneous, transcatheter occlusion device intended to close a patent foramen ovale.

Indication for Use

The Amplatzer™ Talisman™ PFO Occluder is indicated for percutaneous transcatheter closure of a patent foramen ovale (PFO) to reduce the risk of recurrent ischemic stroke in patients who have had an ischemic stroke due to a presumed paradoxical embolism.

Performance Characteristics

The Amplatzer™ Talisman™ PFO Occluder consists of the following:

- PFO occluder
- Amplatzer™ Trevisio™ delivery cable
- Loader
- Hemostasis valve with extension tube and three-way stopcock
- Plastic vise

It is intended to be delivered with the Amplatzer™ Talisman™ Delivery Sheath.

The Amplatzer™ Talisman™ PFO Occluder is a self-expandable, double-disc device made from a nitinol wire mesh available in four sizes to meet a variety of PFO anatomies (as described in Table 1 in this chapter). It is pre-attached to a 120-cm long delivery cable that facilitates placement of the device in a PFO via femoral vein access. The loader

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reduces the size of the Amplatzer™ Talisman™ PFO Occluder, allowing introduction of the device into the sheath. The hemostasis valve with extension tube and stopcock allows flushing of the delivery system and controls blood backflow during the procedure.

Contraindications

- Presence of thrombus at the intended site of implant, or documented evidence of venous thrombus in the vessels through which access to the defect is gained
- Patients with intra-cardiac thrombus, mass, vegetation, or tumor
- Patients whose vasculature, through which access to the PFO is gained, is inadequate to accommodate the appropriate sheath size
- Patients with anatomy in which the required Amplatzer™ Talisman™ PFO device size would interfere with other intracardiac or intravascular structures, such as valves or pulmonary veins
- Patients with another source of right-to-left shunts, including an atrial septal defect or a fenestrated septum
- Patients with active endocarditis or other untreated infections
- Patients who are unable to tolerate intra-procedural anticoagulation or post-procedural anti-platelet therapy

Warnings

- Do not use an open or damaged pouch; do not use a damaged device.
- Patients who are at increased risk for venous thromboembolic events should be managed with a thromboembolic risk reduction regimen after the PFO closure following standard of care.
- The safety and effectiveness of the Amplatzer™ Talisman™ PFO Occluder has not been established in patients with a hypercoagulable state.
- Prepare for situations that require percutaneous or surgical removal of this device. This includes availability of a surgeon and access to an operating room.
- Embolized devices must be removed as they may disrupt critical cardiac functions. Do not remove an embolized occluder through intracardiac structures unless the occluder is fully recaptured inside a catheter or sheath.
- The Amplatzer™ Talisman™ PFO Occluder device consists of a nickel-titanium alloy, which is generally considered safe. However, in vitro testing has demonstrated that nickel is released from this device for a minimum of 60 days. Patients who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may also develop an allergy to nickel if this device is implanted.
- Transient hemodynamic compromise may be encountered during device placement, which may require fluid replacement or other medications as determined by the physician.
- Prior to device detachment, evaluate the position of the device relative to the free atrial wall and the aortic root using echocardiography (see Figure 2 and Figure 3 in Appendix A: Supplemental Information).
- Use echocardiography to ensure that the device does not impinge on the free atrial wall or aortic root.
- Do not release the device from the delivery cable if the device does not conform to its original configuration, or if the device position is unstable or if the device interferes with any adjacent cardiac structure (such as Superior Vena Cava (SVC), Pulmonary Vein (PV), Mitral Valve (MV), Coronary Sinus (CS), aorta (AO)). If the device interferes with an adjacent cardiac structure, recapture the device and redeploy. If still unsatisfactory, recapture the device and either replace with a new device or consider alternative treatments. See Table 1.
- DO NOT use the Amplatzer™ Talisman™ PFO Occluder after the Use-by date stated on the package label.
- This device was sterilized with ethylene oxide and is for single use only; never reuse or re-sterilize the system. Use of expired, reused, or re-sterilized devices may result in infection.
- This device should be used only by physicians who are trained in standard transcatheter techniques.

Precautions for Special Populations

- Pregnancy: The safety and effectiveness of this occluder has not been established during pregnancy. Fluoroscopic X-ray guidance may be used during placement of the device. The risk of increased X-ray exposure for patients who are pregnant must be weighed against the potential benefits of this technique.
- Nursing Mother: The safety and effectiveness of this occluder has not been established in lactating mothers. There has been no quantitative assessment for the presence of leachables in breast milk.
- Pediatric Population: The safety and effectiveness of this occluder has not been established in a pediatric population.

Precautions

- Aspirin (325 mg/day) (or alternative antiplatelet/anticoagulant, if patient has aspirin intolerance) is recommended to be started at least 24 hours prior to the procedure.
- Antibiotics should be administered peri-procedurally.
- Patients should be fully heparinized throughout the procedure using adequate dosing so as to keep the activated clotting time (ACT) greater than 200 seconds.

Patient Counseling Information

Physicians should review the following information when counseling patients about the Amplatzer™ Talisman™ PFO Occluder and the implant procedure:

- PFO closure with the Amplatzer™ Talisman™ PFO Occluder can only reduce the risk for a recurrent stroke due to a paradoxical embolism through a PFO.
 - With aging there is an increased likelihood that non-PFO related risks for stroke may develop and cause a recurrent ischemic stroke independent of PFO closure.
- For procedural risks associated with the Amplatzer™ Talisman™ PFO Occluder, refer to Potential Adverse Events.
- The need for adherence to a defined adjunctive antithrombotic therapy following implantation of the Amplatzer™ Talisman™ PFO Occluder.
- Patients with a history of DVT or PE may benefit from continuation or resumption of anticoagulation therapy following implantation of the Amplatzer™ Talisman™ PFO Occluder to reduce the risk of recurrent DVT or PE.
- It is recommended that the medical team (neurologist and cardiologist) and the patient engage in a shared decision-making process and discuss the risks and benefits of PFO closure in comparison to using antithrombotic therapy alone, while taking into account the patient's values and preferences.

MR Conditional

Non-clinical testing has demonstrated that the Amplatzer™ Talisman™ PFO Occluder is MR Conditional. A patient with the Amplatzer™ Talisman™ PFO Occluder can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system-reported whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

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Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 3ºC after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 16 mm from the Amplatzer™ Talisman™ PFO Occluder when imaged with a gradient echo pulse sequence in a 3.0T MR system.

Potential Adverse Events

Potential adverse events that may occur during or after a procedure using this device may include, but are not limited to:

- Air embolus
- Allergic reaction/toxic effect due to: anesthesia, contrast media, medication, or metal
- Arrhythmia
- Arteriovenous fistulae
- Bleeding
- Cardiac perforation
- Cardiac tamponade
- Chest pain
- Death
- Deep vein thrombosis
- Device embolization
- Device erosion
- Endocarditis
- Esophagus injury
- Fever
- Headache/migraine
- Hematoma
- Hypertension/hypotension
- Infection
- Myocardial infarction
- Pacemaker placement secondary to PFO device closure
- Pain
- Pericardial effusion
- Pericarditis
- Peripheral embolism
- Pseudoaneurysm
- Pulmonary embolism
- Reintervention for residual shunt/device removal
- Stroke
- Transient ischemic attack
- Thrombus formation
- Valvular regurgitation
- Vascular access site injury
- Vessel perforation

Directions for Use

Materials for Use with this Device

- 0.035-inch Amplatzer™ Guidewire (9-GW-002)
- Amplatzer™ Talisman™ Delivery Sheath

For more information about the delivery sheath, see the Amplatzer™ Talisman™ Delivery Sheath instructions for use. See Table 1 in Appendix A: Supplemental Information for information about device sizes and compatible delivery sheathes; see Table 1 below for sizing guidelines.

The following materials are pre-packaged with the Amplatzer™ Talisman™ PFO Occluder:

- Amplatzer™ Trevisio™ delivery cable attached to occluder
- Loader
- Hemostasis valve with extension tube and three-way stopcock
- Plastic vise

Table 1. Device Sizing Guidelines

PFO Morphology	Example Anatomical Characteristics	Suggested Amplatzer™ Talisman™ PFO Occluder Size (mm)
Simple PFO or PFO with a non-prominent ASA PFO where a secure device position and effective PFO closure can be achieved when using the 25 mm device size	1. Absence of ASA, long tunnel, and thickened septum secundum 2. Non-prominent ASA (<20 mm total excursion) without a long tunnel (≥10 mm length) and without a thickened septum secundum (≥10 mm thickness)	25
Complex PFO PFO with one or more anatomical characteristic that may complicate the ability to achieve a secure device position and effective PFO closure when using the 25 mm device size	1. ASA (≥10 mm excursion) with long tunnel (≥10 mm length) 2. ASA (≥10 mm excursion) with thickened septum secundum (≥10 mm thickness) 3. Prominent ASA with excessive mobility (≥20 mm total excursion) 4. Lipomatous hypertrophy of septum secundum (≥15 mm thickness)	30 or 35
PFO with small anatomy Anatomy not suitable for 25 mm device size secondary to interference with adjacent cardiac structures	1. Septal primum length <20 mm	18

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Table 1. Device Sizing Guidelines (Continued)

PFO Morphology	Example Anatomical Characteristics	Suggested Amplatzer™ Talisman™ PFO Occluder Size (mm)
NOTE: Evaluate the position of the device after deployment, but before detachment. Use echocardiography to ensure that the device does not impinge on the free atrial wall or aortic root. If the device interferes with an adjacent cardiac structure (such as free atrial wall or aortic root), recapture the device and redeploy. If device position remains unsatisfactory, recapture the device and either replace with a smaller device (18 mm or 25 mm) or consider alternative treatments. Use of the 35 mm device should only be considered in cases where a secure device position cannot be achieved with the 30 mm device. The presence of an ASA alone does not necessarily prevent successful PFO closure with a 25 mm device size. In the RESPECT trial[a], 180 patients (36%) in the device closure group had an ASA. The 25 mm device size was used in the majority of patients with an ASA (77%) to close the PFO, and at six months post-implant, effective closure was achieved in 95% of these patients. There were no cases of device embolization in any patient in the study.		

a. Saver, Jeffrey L., John D. Carroll, David E. Thaler, Richard W. Smalling, Lee A. MacDonald, David S. Marks, and David L. Tirschwell. "Long-term outcomes of patent foramen ovale closure or medical therapy after stroke," New England Journal of Medicine 377, no. 11 (2017): 1022–1032.

Pre-Procedure Care

- Aspirin (325 mg/day) (or alternative antiplatelet/anticoagulant, if patient has aspirin intolerance) is recommended to be started at least 24 hours prior to the procedure.
- Antibiotics can be administered periprocedurally at operator's discretion.
- Patients should be fully heparinized throughout the procedure using adequate dosing so as to keep the activated clotting time (ACT) greater than 200 seconds.

Procedure

CAUTION: Intracardiac echocardiography (ICE) or transesophageal echocardiography (TEE) is recommended as an aid in evaluating the PFO and placing the Amplatzer™ Talisman™ PFO Occluder. If TEE is used, the patient's esophageal anatomy must be adequate for placement and manipulation of the probe.

CAUTION: Be cautious when using fluoroscopic X-ray guidance, which may be used during placement of the device.

CAUTION: Do not use a power injection system to put contrast solution through the sheath.

1. Puncture the femoral vein and perform a standard right-heart catheterization.
2. Perform an angiogram to demonstrate the PFO:
 a. Catheterize the left atrium using a 45° LAO position and cranial angulation of 35°-45°.
 b. Inject contrast medium into the right upper lobe pulmonary vein.

 NOTE: Occluder size and placement are based on the morphology and location of the PFO.
3. Use the J-tip guidewire to gain access through the PFO.
4. Use ICE or TEE to assess PFO morphology.

 NOTE: PFO morphology is best assessed using a short axis view, where the PFO is visualized adjacent to the aortic valve, and the septum primum is oriented in a near horizontal position (ICE septal short-axis view; TEE angles 30°, 45°, and 60°). Using this view, the excursion of the septum primum, length of the PFO tunnel, and the thickness of the septum secundum are assessed. These features should also be assessed using a long axis bicaval view (ICE septal bicaval view; TEE angles 90°, 105°, and 120° — see Figure 2 and Figure 3 in Appendix A: Supplemental Information).

 Atrial Septal Aneurysm (ASA): An atrial septal aneurysm is a redundancy or saccular deformity of the septum primum with increased mobility and bulging into the right atrium or the left atrium. ASA is defined as an excursion of the septum primum of ≥ 10 mm from the plane of the atrial septum into the right atrium or left atrium or a combined total excursion into the right atrium and left atrium of ≥ 15 mm. A prominent ASA is defined as having a combined excursion of ≥ 20 mm.

 PFO Tunnel Length: The PFO tunnel length is defined as the maximum overlap between the septum primum and septum secundum. A long tunnel is defined as a tunnel length ≥ 10 mm.

 Septum Secundum Thickness: The septum secundum thickness is defined as the maximum thickness of the septum secundum within 10 mm from the PFO. A thickened septum secundum is defined as thickness ≥ 10 mm. Lipomatous hypertrophy of the septum secundum is defined as having a thickness of ≥ 15 mm.

 NOTE: If TEE imaging is used, Steps 4 – 6 can be performed prior to femoral access, if desired.
5. Use ICE or TEE to measure the distance from the PFO to the aortic root and the distance from the PFO to superior vena cava orifice (see Figure 3B and Figure 3C in Appendix A: Supplemental Information).

ICE Imaging

2A	PFO as assessed by ICE in short axis view with measurement of septum primum excursion (a), PFO tunnel length (b), and septum secundum thickness (c).
2B	PFO as assessed by ICE in long axis bicaval/SVC view with measurement of septum primum length (a).
2C	PFO as assessed by ICE in septal view with measurement of septum primum excursion into left atrium (a) and right atrium (b).

TEE Imaging

3A	Measurement of septum primum excursion (a), PFO tunnel length (b), and septum secundum thickness (c) in the short axis 50-degree view (range 30 to 60 degrees to visualize the aortic root en face).
3B	Distance from the PFO to aortic root (a) in the 30-degree view (range 0 to 45 degrees to visualize the aortic root en face)
3C	Distance from the PFO to the orifice of the superior vena cava (SVC) (a) in the bicaval 90-degree view (range 80 to 125 degrees to minimize foreshortening of the SVC).

NOTE: If imaging cannot clearly localize the PFO, place a wire through the PFO to help with identification.

6. Size the device based on the PFO morphology while taking into consideration the distances from the PFO to the aortic root and SVC.

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Most PFOs can be adequately covered with the 25 mm device. Consider the use of a larger device (30-35 mm) if the PFO has an atrial septal aneurysm (ASA) with other features (e.g., long tunnel or thickened septum secundum) that may result in a substantial residual shunt or an unsecured position across the septum. Use of the 35 mm device should only be considered in cases where a secure device position cannot be achieved with the 30 mm device.


Use echocardiography and a push-pull test (see Step 20) to evaluate the position of the device after device deployment, but before detachment. If the device position after deployment is not secure, consider the use of a larger device (30 mm or 35 mm).

If the device interferes with an adjacent cardiac structure (such as free atrial wall or aortic root), recapture the device and redeploy. If device position remains unsatisfactory, recapture the device and either replace with a smaller device or consider alternative treatments. Refer to Table 1 in Appendix A: Supplemental Information and Table 1 (above) for device selection and sizing guidelines. Check SVC, inferior vena cava (IVC) and coronary sinus (CS) flows after device deployment, but before detachment.


CAUTION: Care should be taken not to entrap right atrial Chiari networks or large Eustachian valves under the right atrial side of the device.

WARNING: Use echocardiography to ensure that the device does not impinge on the free atrial wall or aortic root.

7. Prepare the Amplatzer™ Talisman™ PFO Occluder for use:
 a. Inspect the sterile package. Do not use the device if the sterile package is open or damaged.
 b. Open the sterile package and inspect the PFO occluder, Amplatzer™ Trevisio™ delivery cable, loader, and hemostasis valve. Do not use if damage is observed or if the occluder is detached from the cable.
 c. Remove the device with attached delivery accessories from the package (delivery cable, loader, and hemostasis valve) and ensure the connection is secure between the loader and the hemostasis valve.
 d. Ensure the PFO occluder is appropriately attached to the delivery cable:
 - Rotate the occluder one full turn counter-clockwise:

 - Rotate occluder clockwise until resistance is felt

 - Rotate occluder counter-clockwise 1/8 of one turn to facilitate device release upon deployment.

 NOTE: If the occluder detaches from the delivery cable, re-attach it to the distal end of the delivery cable by rotating the device clockwise approximately 4 to 5 turns until it is secure. Then, rotate the device counter-clockwise 1/8 of one turn to facilitate occluder release during implantation.
8. Load the occluder in the loader.
 a. Flush the loader and hemostasis valve with sterile saline.
 b. With the occluder and distal end of the loader submerged in sterile saline, pull back on the delivery cable to retract the device inside the loader until the distal marker band of the occluder reaches the loader tip.
9. Flush sterile saline through the hemostasis valve until no air bubbles are visible on the loaded occluder.

 NOTE: The attachment of the occluder to the cable may be checked after loading by rotating the cable clockwise until resistance is felt, then rotating the cable counter-clockwise 1/8 turn.
10. Prepare the Amplatzer™ Talisman™ Delivery Sheath according to the Amplatzer™ Talisman™ Delivery Sheath instructions for use.
11. Insert the dilator into the delivery sheath and tighten the rotating luer.
12. Advance the dilator and delivery sheath over the guidewire, through the PFO, and into the left atrium, confirming correct movement via echo and/or fluoroscopy.

 WARNING: Do not advance the dilator and delivery sheath if resistance is felt.
13. Slowly remove the guidewire and dilator to prevent ingress of air; allow blood back-flow to purge all air from the system.
14. While maintaining a wet-to-wet condition, attach the loader to the delivery sheath and tighten the rotating luer to lock the components together.
15. Advance the delivery cable and occluder through the delivery sheath until the occluder reaches the tip of the sheath. A mark on the cable indicates the occluder is nearing the end of the sheath.

 CAUTION: Do not rotate the cable while advancing the occluder; do not advance the delivery cable and occluder if resistance is felt.
16. Using echocardiography and fluoroscopy for guidance, hold the delivery cable in place while retracting the delivery sheath to deploy the left atrial disc and part of the connecting waist.
17. Pull the left atrial disc gently against the atrial septum. This can be felt or observed by echocardiography or fluoroscopy.
18. Maintain a slight tension on the delivery cable while retracting the delivery sheath approximately 5–10 cm to deploy the right atrial disc.
19. Use echocardiography or angiography to confirm that the device is in place and evaluate for residual shunt or valve insufficiency.
20. Prior to detaching the occluder, perform a gentle to and fro motion (push-pull test) with the delivery cable to assure a secure occluder position across the septum.
21. If the occluder position is unstable or interferes with any adjacent cardiac structure (such as SVC, PV, MV, CS, AO), do not release it from the delivery cable and perform the following steps:
 a. Recapture the occluder by stabilizing the delivery cable and advance the delivery sheath until the occluder is completely within the sheath.
 b. Reposition the occluder and deploy it again, or remove the occluder from the patient.

 NOTE: The occluder can be deployed a maximum of three times. If the position is still unsatisfactory, then remove and replace both the device and the sheath.
22. If the occluder does not conform to its original shape, recapture the device and replace it with a new device.
23. Place the plastic vise on the proximal end of the delivery cable.

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24. Detach the occluder from the delivery cable by turning the delivery cable counterclockwise (indicated by the arrow on the plastic vise).

 In the unlikely event that this should not be possible, advance the delivery sheath against the right atrial disc to secure the occluder and to facilitate detachment.

 WARNING: When the procedure is complete, slowly remove the delivery cable and delivery sheath from the patient. Remove the sheath slowly to prevent an ingress of air.

Post-Procedure Care

- Post-procedure follow-up should be per institution standard of care. Prior to discharge, patient should have a TTE.
- Patients with any observed pericardial effusion following the device implantation should be closely monitored with serial echocardiograms performed to ensure the pericardial effusion is not enlarging or until it resolves.
- Clinical follow-up with a cardiologist, including an echocardiography evaluation, is recommended per the institution's standard of practice and established guidelines.
- Patients who experience a stroke after device implantation should undergo TEE for evaluation of thrombus formation.
- If a left-sided thrombus is identified following device implant, the patient should be evaluated for a hypercoagulable state and therapeutic anticoagulation (warfarin) should be initiated for a minimum of three months. A TEE should be performed following anticoagulation treatment to confirm resolution of the device-related thrombus. Thrombolysis or surgical removal of the device should be considered if the patient does not respond to anticoagulant therapy.
- Patients should be educated to seek immediate medical attention that includes an echocardiogram, if they develop signs or symptoms of hemodynamic instability such as chest pain, arrhythmia, fainting, or shortness of breath.
- Patients should take appropriate endocarditis prophylaxis for six months following device implantation. The decision to continue endocarditis prophylaxis beyond six months is at the discretion of the physician.
- Patients should be treated with antithrombotic therapy, such as aspirin for six months post implant. The decision to continue antithrombotic therapy beyond six months is at the discretion of the physician.
- Patients should be instructed to avoid strenuous activity for a minimum of one month post-device implant or as directed by physician. Strenuous activities may lead to the increased risk of adverse events including erosion. Patients should be reminded that if they experience any symptoms of shortness of breath or chest pain at any time, and especially after strenuous activity, they should seek medical care immediately.
- For patients with a history of PE or DVT, chronic anticoagulation rather than antiplatelet therapy should be considered post-device implantation.

Reporting Device Incidents

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer.

Storage and Handling Information

There are no special storage requirements beyond what is typical for medical devices.

This instructions for use is recyclable.

Disposal

Amplatzer™ Talisman™ PFO Occluder accessories should be appropriately classified as biohazards and disposed of in compliance with applicable facility procedures and local and country laws and regulations. Return the Amplatzer™ Talisman™ PFO Occluder accessories to Abbott Medical at the end of their operating life.

Warranty

Abbott Medical warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Abbott Medical's obligation under this warranty is limited to replacing or repairing at its option, at its factory, this product if returned within the warranty period to Abbott Medical and after confirmed to be defective by the manufacturer.

EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

See the Terms and Conditions of Sale for further information.

Australian Warranty

This section is only provided in English.

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical at the address provided with the product and after Abbott Medical confirms that the goods are defective.

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED. See the Terms and Conditions of Sale for further information.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.
© 2022 Abbott. All Rights Reserved.