COVIDIEN™

REF CNN/SNN

INVOS™

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INVOS™

Infant Regional Saturation Sensor

For use with INVOS™ system technology.

INTENDED USE

The OxyAlert NIRSensor disposable sensor Model CNN/SNN has a small adhesive pad and a gentle hydrocolloid adhesive that has been designed for use with the INVOS™ cerebral/somatic monitoring systems for monitoring of site-specific regional oxygen saturation (rSO_2) in pediatric, infant, and neonatal patients weighing < 5 kilograms.

INDICATIONS FOR USE

The INVOS™ OxyAlert NIRSensor disposable sensor Model CNN/SNN is indicated for single patient use when cerebral/somatic monitoring of site-specific regional oxygen saturation (rSO_2) is required in patients weighing < 5 kilograms. This sensor is only intended to be used with the INVOS™ Near Infrared Spectroscopy (NIRS) technology including monitoring systems and devices integrated with INVOS™ NIRS technology. For additional information regarding setup and use of the INVOS™ System including indications for use, contraindications, warnings and cautions, consult the Monitoring System Operator's Manual.

CONTRAINDICATIONS

None known.

WARNINGS

- **The sensor is designed for single patient use only and may not be used on multiple patients due to an increased risk of cross-contamination among patients. Do not autoclave or gas sterilize the sensors.**
- **The sensor is designed for external use only as described in the instructions. Do not use the sensors internally for any reason.**
- **Do not immerse the INVOS™ System or sensor in any liquids as they may cause electric shock hazard or damage the device.**
- **Do not use the INVOS™ System in the presence of flammable anesthetics or in other flammable environments.**
- **Improper sensor placement can lead to inaccurate readings. Follow site selection and placement instructions in this document.**
- **The prospective clinical value of data from the INVOS™ System has not been demonstrated in disease states. The INVOS™ System should not be used as the sole basis for diagnosis or therapy.**
- **Carefully route cabling to reduce the possibility of strangulation.**
- **Choking hazard — The sensor cables include a strain-relief clip that, if detached, may pose a choking hazard.**
- **Environments with excessive ambient light such as bright sunlight or strong operating room lighting may require loosely covering the area of the sensor with an opaque drape.**
- **Do not use the sensor on patients who exhibit allergic reactions to the adhesive tape.**
- **Disconnect the sensor from the patient during magnetic resonance imaging (MRI) scanning. Induced current may cause burns.**

CAUTIONS

- INVOS™ readings represent a small volume of tissue beneath the sensor and may not reflect oxygenation disturbances that occur elsewhere.
- Make sure all connectors are fully engaged and free from moisture. Moisture intrusion may cause inaccurate readings, erratic readings or no readings at all.
- Different INVOS™ System sensors (adult, pediatric and infant/neonatal) cannot be used simultaneously on the same monitor.
- Use only Covidien recommended or provided accessories. Use of the INVOS™ System with any other sensor will compromise accuracy.
- Use care when placing or removing sensor. Do not place on broken or undeveloped skin.
- Do not place the sensor on regions with severe tissue edema to reduce the possibility of skin lesions.
- Use of an electrosurgical/electrocautery instrument in the vicinity of the INVOS™ System may interfere with the signal and result in no readings.
- To avoid pressure sores do not apply pressure (e.g. headbands, wraps, tape) to sensor.
- Sensors do not need to be removed for X-ray or CT Scan. However, sensors will appear on image.
- Sensors must be removed for MRI.
- When two sensors are placed, they must be connected into the same preamplifier.
- Do not use a damaged sensor or cable. Do not use a sensor with exposed optical components.

INSTRUCTIONS FOR USE

Pre-use Check: Remove the sensor from the package and examine for visual signs of damage. If any signs of damage are observed, select another sensor.

To achieve optimum adhesion, the patient's skin must be clean and dry. Dry skin with a gauze pad. Warm the sensor in your hands or an incubator to ease placement.

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Note: Refer to illustrations on page 2 throughout this section.

Cerebral Site Selection: Select sensor site on the right or left side of forehead. Placement of the sensor in other cerebral locations, or over hair, may cause inaccurate readings, erratic readings, or no readings at all. Do not place the sensor over nevi, sinus cavities, the superior sagittal sinus, subdural or epidural hematomas or other anomalies such as arteriovenous malformations, as this may cause readings that are not reflective of brain tissue or no readings at all.

EEG and level of consciousness sensors can be used simultaneously with the INVOS™ infant rSO_2 sensors, overlapping sensor pads if necessary, as long as the three optical windows on the INVOS™ sensor are not blocked.

Somatic Site Selection: Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over fatty deposits, hair or bony protuberances. Do not place the sensor over nevi, hematomas or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted on this Instruction For Use. Placements may include:

- Posterior flank (T10-L2, right or left of midline) 3
- Abdomen 4
- Forearm
- Calf
- Upper arm
- Chest
- Upper leg

Patient Prep: To achieve optimum adhesion, the patient's skin must be clean and dry. Dry skin with a gauze pad. Warm the sensor in your hands or an incubator to ease placement 1.

Sensor Placement: Holding the cable, peel back 2 the two halves of the white liner from the clear adhesive layer beginning at the center of the sensor pad. Apply to the skin 5. Continue applying the sensor by smoothing it to the skin from the center outward. Ensure edges of the sensor are sealed.

Monitoring: Plug the sensor cable into the preamplifier connector 6. When two sensors are placed, they must be connected into the same preamplifier. Secure the sensor cable to a fixed object to avoid strain on the sensor-to-skin interface using strain-relief clips. Ensure the cable is properly inserted into the preamplifier (via the reusable infant sensor adapter cables (RIC) on the 7100 system) 7. Calibration is automatic and monitoring will begin in seconds.

Status messages on the INVOS™ System display will appear if monitoring conditions are compromised. Check the sensor periodically to assess skin integrity (irritation or injury) and the security of sensor placement. Always ensure sensor is properly sealed to skin to reduce light interference. For extended monitoring, it is recommend to use a new sensor every 24 hours or if adhesive is inadequate to seal the sensor to the skin.

If present, the following may cause poor performance:

- "Cardiogreen,""Indigo Carmine,""Methylene Blue" or other intravascular dyes
- Carboxyhemoglobin or other dyshemoglobins
- Hemoglobinopathies
- Conjugated hyperbilirubinemia (direct)
- Myoglobin (Mb), hemoglobin from muscular tissue, in the blood.
- Dark skin pigment
- Externally applied coloring agents (dye, pigmented cream)

Cleaning/Disinfection/Sterilization

The sensor is intended for single use only and cannot be adequately cleaned, disinfected, or sterilized for safe reuse. Attempts to clean, disinfect, or sterilize the sensor for reuse may result in product failure and/or safety risk to the patient.

Disposal

Dispose of the sensor after a single use. Follow local governing ordinances regarding disposal of biologically hazardous waste.

Environmental Specification

- Operational temperature range 10°C to 35°C
- Storage temperature range 10°C to 30°C
- Transport temperature range -10°C to 50°C up to 72 hours
- Operational relative humidity range 15% to 95%
- Storage relative humidity range 10% to 95%
- Atmospheric pressure range 61.6 kPa to 107.5 kPa
- Storage relative pressure range 46.5 kPa to 107.5 kPa

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INVOS™

Infant Regional Saturation Sensor

For use with INVOS™ system technology.

INTENDED USE

The OxyAlert NIRSensor disposable sensor Model CNN/SNN has a small adhesive pad and a gentle hydrocolloid adhesive that has been designed for use with the INVOS™ cerebral/somatic monitoring systems for monitoring of site-specific regional oxygen saturation (rSO_2) in paediatric, infant and neonatal patients weighing < 5 kilograms.

INDICATIONS FOR USE

The INVOS™ OxyAlert NIRSensor disposable sensor Model CNN/SNN is indicated for single patient use when cerebral/somatic monitoring of site-specific regional oxygen saturation (rSO_2) is required in patients weighing < 5 kilograms. This sensor is only intended to be used with the INVOS™ Near Infrared Spectroscopy (NIRS) technology including monitoring systems and devices integrated with INVOS™ NIRS technology. For additional information regarding setup and use of the INVOS™ System including indications for use, contraindications, warnings and cautions, consult the Monitoring System Operator's Manual.

CONTRAINDICATIONS

None known.

WARNINGS

- **The sensor is designed for single patient use only and may not be used on multiple patients due to an increased risk of cross-contamination among patients. Do not autoclave or gas sterilise the sensors.**
- **The sensor is designed for external use only as described in the instructions. Do not use the sensors internally for any reason.**
- **Do not immerse the INVOS™ System or sensor in any liquids as they may cause electric shock hazard or damage the device.**
- **Do not use the INVOS™ System in the presence of flammable anaesthetics or in other flammable environments.**
- **Improper sensor placement can lead to inaccurate readings. Follow site selection and placement instructions in this document.**
- **The prospective clinical value of data from the INVOS™ System has not been demonstrated in disease states. The INVOS™ System should not be used as the sole basis for diagnosis or therapy.**
- **Carefully route cabling to reduce the possibility of strangulation.**
- **Choking hazard – The sensor cables include a strain-relief clip that, if detached, may pose a choking hazard.**
- **Environments with excessive ambient light such as bright sunlight or strong operating room lighting may require loosely covering the area of the sensor with an opaque drape.**
- **Do not use the sensor on patients who exhibit allergic reactions to the adhesive tape.**
- **Disconnect the sensor from the patient during magnetic resonance imaging (MRI) scanning. Induced current may cause burns.**

CAUTIONS

- INVOS™ readings represent a small volume of tissue beneath the sensor and may not reflect oxygenation disturbances that occur elsewhere.
- Make sure all connectors are fully engaged and free from moisture. Moisture intrusion may cause inaccurate readings, erratic readings or no readings at all.
- Different INVOS™ System sensors (adult, paediatric and infant/neonatal) cannot be used simultaneously on the same monitor.
- Use only Covidien recommended or provided accessories. Use of the INVOS™ System with any other sensor will compromise accuracy.
- Use care when placing or removing sensor. Do not place on broken or undeveloped skin.
- Do not place the sensor on regions with severe tissue oedema to reduce the possibility of skin lesions.
- Use of an electrosurgical/electrocautery instrument in the vicinity of the INVOS™ System may interfere with the signal and result in no readings.
- To avoid pressure sores do not apply pressure (e.g. headbands, wraps, tape) to sensor.
- Sensors do not need to be removed for X-ray or CT Scan. However, sensors will appear on image.
- Sensors must be removed for MRI.
- When two sensors are placed, they must be connected into the same preamplifier.
- Do not use a damaged sensor or cable. Do not use a sensor with exposed optical components.

INSTRUCTIONS FOR USE

Pre-use Check: Remove the sensor from the package and examine for visual signs of damage. If any signs of damage are observed, select another sensor.

To achieve optimum adhesion, the patient's skin must be clean and dry. Dry skin with a gauze pad. Warm the sensor in your hands or an incubator to ease placement.

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Note: Refer to illustrations on page 2 throughout this section.

Cerebral Site Selection: Select sensor site on the right or left side of forehead. Placement of the sensor in other cerebral locations, or over hair, may cause inaccurate readings, erratic readings, or no readings at all. Do not place the sensor over nevi, sinus cavities, the superior sagittal sinus, subdural or epidural haematomas or other anomalies such as arteriovenous malformations, as this may cause readings that are not reflective of brain tissue or no readings at all.

EEG and level of consciousness sensors can be used simultaneously with the INVOS™ infant rSO_2 sensors, overlapping sensor pads if necessary, as long as the three optical windows on the INVOS™ sensor are not blocked.

Somatic Site Selection: Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over fatty deposits, hair or bony protuberances. Do not place the sensor over nevi, haematomas or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted on this Instruction for Use. Placements may include:

- Posterior flank (T10-L2, right or left of midline) 3
- Abdomen 4
- Forearm
- Calf
- Upper arm
- Chest
- Upper leg

Patient Prep: To achieve optimum adhesion, the patient's skin must be clean and dry. Dry skin with a gauze pad. Warm the sensor in your hands or an incubator to ease placement 1.

Sensor Placement: Holding the cable, peel back 2 the two halves of the white liner from the clear adhesive layer beginning at the centre of the sensor pad. Apply to the skin 5. Continue applying the sensor by smoothing it to the skin from the centre outwards. Ensure edges of the sensor are sealed.

Monitoring: Plug the sensor cable into the preamplifier connector 6. When two sensors are placed, they must be connected into the same preamplifier. Secure the sensor cable to a fixed object to avoid strain on the sensor-to-skin interface using strain-relief clips. Ensure the cable is properly inserted into the preamplifier (via the reusable infant sensor adapter cables (RIC) on the 7100 system) 7. Calibration is automatic and monitoring will begin in seconds.

Status messages on the INVOS™ System display will appear if monitoring conditions are compromised. Check the sensor periodically to assess skin integrity (irritation or injury) and the security of sensor placement. Always ensure sensor is properly sealed to skin to reduce light interference. For extended monitoring, it is recommended to use a new sensor every 24 hours or if adhesive is inadequate to seal the sensor to the skin.

If present, the following may cause poor performance:

- 'Cardiogreen', 'Indigo Carmine', 'Methylene Blue' or other intravascular dyes
- Carboxyhaemoglobin or other dyshaemoglobins
- Haemoglobinopathies
- Conjugated hyperbilirubinaemia (direct)
- Myoglobin (Mb), haemoglobin from muscular tissue, in the blood.
- Dark skin pigment
- Externally applied colouring agents (dye, pigmented cream)

Cleaning/Disinfection/Sterilisation

The sensor is intended for single use only and cannot be adequately cleaned, disinfected or sterilised for safe reuse. Attempts to clean, disinfect or sterilise the sensor for reuse may result in product failure and/or safety risk to the patient.

Disposal

Dispose of the sensor after a single use. Follow local governing guidelines regarding disposal of biologically hazardous waste.

Environmental Specification

- Operational temperature range 10°C to 35°C
- Storage temperature range 10°C to 30°C
- Transport temperature range -10°C to 50°C up to 72 hours
- Operational relative humidity range 15% to 95%
- Storage relative humidity range 10% to 95%
- Atmospheric pressure range 61.6 kPa to 107.5 kPa
- Storage relative pressure range 46.5 kPa to 107.5 kPa

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	Non-Sterile
	Single use
	Do not use if package is damaged
	Consult instructions for use

EC REP	Authorized representative in the European Community
REF	Catalogue number
	Manufacturer
LOT	Lot number
	Date of manufacture

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Rx ONLY

DEHP
Not made with DEHP

Caution, consult accompanying documents

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Part No. PT00133881 Rev A 2021-01

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG.
™* brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

© 2020 Covidien.

Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

EC REP
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands.

www.covidien.com

COVIDIEN™

REF SPFB

INVOS™

Pediatric Regional Saturation Sensor

en-GB	(English)	Instructions For Use
cs	(Čeština)	Návod k použití
da	(Dansk)	Brugervejledning
de	(Deutsch)	Gebrauchsanweisung
es	(Español)	Modo de empleo
fr	(Français)	Instructions d'emploi
it	(Italiano)	Istruzioni per l'uso
hu	(Magyar)	Használati útmutató
nl	(Nederlands)	Gebruiksaanwijzing
no	(Norsk)	Bruksanvisning
pl	(Polski)	Instrukcja użytkowania
pt	(Português)	Instruções de Utilização
ro	(Română)	Instrucţiuni de utilizare
sk	(Slovenčina)	Návod na použitie
fi	(Suomi)	Käyttöohje
sv	(Svenska)	Bruksanvisning
tr	(Türkçe)	Kullanım Talimatları
el	(Ελληνικά)	Οδηγίες χρήσης
ru	(Русский)	Инструкция по применению

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore, Ireland.