CHEST DRAIN P.A.S.

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Cardiopulmonary Post-operative Autotransfusion System VENICE P.A.S. and ROME P.A.S.

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Chest Drain P.A.S.
Cardiopulmonary Post-operative Autotransfusion System
VENICE P.A.S. and ROME P.A.S.

The instructions for use given below are not exhaustive for the *VENICE P.A.S and ROME P.A.S* devices.
Also refer to the specific instructions for use of the *VENICE* and *ROME* devices for thoracic drainage management.

WARNING
- **Carefully read the instructions before use.**
- **The product must only be used by adequately trained and qualified personnel.**
- **The product must be used according to these instructions and multi-chamber thoracic drainage *VENICE* and/or *ROME* devices instructions. EUROSETS cannot be held responsible for damage deriving from incorrect use of the product.**
- **The product must not be subjected to any further processing.**
- **The product has been sterilized by ethylene oxide. Sterility is guaranteed only if the packaging is undamaged.**
- **The device must be handled aseptically**
- **It is recommended to carefully inspect the product before use.**
- **The product is non-pyrogenic.**
- **The product is disposable and to be used only once on one patient.**
- **Do not resterilize or clean for reuse.**
- **Do not subject to further treatments.**
- **Improper reuse may cause cross-infections and alteration of the device performance.**
- **If disinfectants type Betadine® and/or antibiotics not for parenteral use are utilised during surgery, the wound must thoroughly be washed with sterile physiological solution before suture.**
- **Avoid conditions of use which expose the blood to a temperature of higher than 37°C.**
- **Reinfuse only the blood collected within the first 6 (six) hours from suture.**
- **Do not store the blood collected, reinfuse immediately.**
- **Always reinfuse the blood with a filter for micro-aggregates.**
- **Do not reinfuse with devices under pressure to reduce the risks of gaseous embolism.**
- **The measurements made by the device are approximate.**
- **Remove any protective caps on the connectors before making any connection.**
- **Any serious incident occurring in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

CAUTION
- Store in a dry place between 0°C/32°F and 40°C/104°F.
- Do not use solvents such as alcohol, ether, acetone, etc., as contact may cause damage to the device.
- After use dispose of the device in accordance with the regulations applicable for biological hazardous waste in force in the country of use.

PRODUCT DESCRIPTION
The ***Chest Drain P.A.S.*** *Post-operative Autotransfusion System* is composed of:
1. ***VENICE P.A.S.*** *or* ***ROME P.A.S. multi-chamber thoracic drainage device.***
2. Accessories for ***P.A.S.*** reinfusion, such as:
 - ***P.A.S.*** reinfusion bag
 - Microaggregate filter
 - ***P.A.S.*** connectors and connection lines.

In addition, the following accessories are required for use of the system:
- Drainage catheter/catheters.
- Male 3/8 connector for drainage catheter connection(s)
- A reusable connection tube to the vacuum source with an inside diameter between 4.8 mm (3/16") and 8.7 mm (11/32").

INTENDED PURPOSE
The ***Chest Drain P.A.S.*** *Post-operative Autotransfusion System* has been optimised to respond to the specific thoracic drainage requirements with simultaneous blood recovery from the cardio surgical wound drainages.
The ***Chest Drain P.A.S.*** *Post-operative Autotransfusion System* has been designed to allow direct reinfusion of the recovered blood by means of the closed-circuit technique.

Max use period:
Reinfuse only the blood collected within the first 6

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(six) hours from suture.

PATIENT POPULATION
Adult patients.

INDICATIONS FOR USE
Treatment of patients undergoing cardiac surgery.

Contraindication
*** Reinfusion of the recovered blood is contraindicated:**

a) when the recovery is performed in the presence of bacterial and/or malignant contamination;

b) in the presence of massive haemolysis, in cases of significant coagulopathy or liver or renal failure. In these cases the blood may be reinfused subject to treatment with cell separators;

c) treatment of blood with cell separators is advisable if there is a large blood mass to be reinfused (>800-1000 ml).

Warning: The cell separator may be used in connection with ***VENICE P.A.S - ROME P.A.S*** exclusively with a suitable connector supplied on request. This connector, specific for each machine, is to be connected to the "lock connection" of the "emptying tube". Then follow the procedure in the chapter "Blood reinfusion".

Warning: The blood recoverable for purposes of autotransfusion with the ***VENICE P.A.S or ROME P.A.S device is only that collected in the 1st subchamber with a capacity of 1350 ml.*** The blood collected in the 2nd subchamber (more than the 1350 ml capacity of the collection chambers) is not recoverable for autotransfusion purposes.


INTENDED USER
VENICE P.A..S and ROME P.A.S. devices are intended to be used by thoracic surgery, cardiac surgery, trauma surgery and Intensive Care Unit (ICU) physician or by nurse under the physician direction.

TECHNICAL FEATURES


	Venice P.A.S.	Rome P.A.S.
Technical features		
Capacity	2500 ml	
Transfer blood capacity	1350 ml	
Blood filters	Macro filter inside chamber + 40 µ microaggregate filter (where present or supplied as accessory)	
Graduated scale in Collection chamber	Single collection chamber divided into 2 sub-chambers: The first sub-chamber is graduated in 10ml increments up to the maximum capacity of 1350ml. The second sub-chamber is graduated in 10ml increments up to the maximum capacity of 2500ml.	
Suction Control	From -5 to -25 cmH_20	From -10 to -40 cmH_20
Body size	Height (excluding tubing): 360mm Depth: 76mm Width: 270mm	
	Venice P.A.S.	**Rome P.A.S.**
Lines and connections		
Drainage connection line	Single	
Connection to the catheter	Depending on the version: • Connector for thoracic drainage catheters with an inside diameter (ID) between 6,36 mm (1/4") and 12,7 mm (1/2") or Connector for thoracic drainage catheters with an inside diameter (ID) of 6,36 mm (1/4") or • ale connector 3/8" (supplied on request).	
Connection to vacuum source	Connector for a tube of connection to the vacuum source with an inside diameter (ID) between 4.8 mm (3/16") and 8.7 mm (11/32") is required.	

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Instructions for use of the *P.A.S.* accessories in connection with the *VENICE P.A.S. or ROME P.A.S.* device for blood reinfusion

1. Carefully inspect the protective wrapping of the accessories. If it is open or damaged, do not use them.
2. Quantify the blood spilled from the drains directly in the blood collection chamber of the drainage device.
3. Set the drainage/suction device into standby (drainage and vacuum tubes closed/vent open) by pulling the "tube clamp plate" to the left.

4. Once the blood has passed through:
a) Close the clamp.
b) Pull the "tube-clamp plate" to the right to restore the normal drainage/suction functions.
5. Proceed with blood reinfusion using the microaggregate filter (microfiltration ≤ 40µm).

Warning
Do not disconnect the reinfusion set from the reservoir as long as there is a possibility of further reinfusion.
6. For multiple reinfusion, do not disconnect the microaggregate filter from the reinfusion bag or the patient; only close the upstream clamp (on the bag) and then operate as described in steps 2 to 5.

Warning

If there are difficulties in the passage of blood from the reservoir to the reinfusion bag, act as follows:
a) Check that you have correctly carried out the procedure described above.
b) Increase the difference in height between the unit and the reinfusion bag.
c) Lightly shake the blood in the collection chambers in order to remove any aggregations which may have formed on the bottom.

WARRANTY CONDITIONS

This limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.
EUROSETS guarantees that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.
EUROSETS guarantees that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.
However, EUROSETS cannot guarantee that the user will use the device correctly, nor that an incorrect diagnosis or therapy and/or that particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.
EUROSETS, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences

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arising directly or indirectly from improper use of this device.
EUROSETS will replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by EUROSETS up to the time of delivery to the end user unless such defect has been caused by mishandling by the purchaser.
The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose.
No person, including any representative, agent, dealer, distributor or intermediary of EUROSETS or any other industrial or commercial organisation is authorized to make representation or warranty concerning this medical device except as expressly stated herein.
EUROSETS disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein.
The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with EUROSETS, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.
The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and Jurisdiction. The court chosen is the Court of Modena (Italy).

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MICROAGGREGATES REMOVAL FILTER

WARNINGS

- Store in a dry place between 32°F/0°C and 104°F/40°C.
- Carefully read the instructions before using the product.
- The product must only be used by adequately trained and qualified personnel.
- The product must be used according to these instructions. EUROSETS is not responsible for damage deriving from incorrect use of the product.
- Before using the product it is advisable to carefully inspect it. Shipping and handling could cause structural and functional damage to the device
- Atoxic - STERILE non-pyrogenic
- Sterile only if protective bag is not opened, damaged or broken
- Ethylene oxide sterilized
- Do not expose to the heat
- The product is disposable and to be used only once on one patient.
- Do not resterilize or clean for reuse.
- Do not subject to further treatments.
- Improper reuse may cause cross-infections and alteration of the device performance.
- After use dispose of the device in accordance with the regulations in force in the country of use.
- The circuit must be used immediately after removing the protective caps

INTENDED USE

The "GOCCIA MICRO" device is intended to be used for transfusion of whole blood and blood products in order to remove the microaggregates.


PRECAUTION

- If transfusion is performed with pressurised devices:
 a) Do not exceed 400mmHg
 b)Be extremely careful not to infuse "air" as there is a risk of "gaseous embolism".
- The device has been tested to allow transfusion of several units of whole blood and/or blood products.

INSTRUCTIONS FOR USE FOR VERSIONS WITH PRE-ASSEMBLED DRIP CHAMBER


Warning: Do not use if the packaging has been damaged or the caps has fallen off.

1


Align the roller clamp with the drip chamber and close it while turning it upside down. Insert the spike of the device in the bag outlet connector.

2

Open the roller clamp; with regular pressure on the bag slowly fill the filter chamber and 2/3 of the drip chamber with blood. Close the roller clamp

3

Hang up the bag and bleed the air from the drip chamber line before connecting the device to the patient. Use the roller clamp to adjust to the desired drip rate.

4

- To change blood bag, close the roller clamp, so that the Filter remains filled.
- WARNING: If the filter body empties out, the device is no longer usable.

INSTRUCTIONS FOR USE FOR VERSIONS WITHOUT PRE-ASSEMBLED DRIP CHAMBER

1) Remove the lower cap from the filter and connect it to the drip chamber.
2) Continue as for the versions with pre-assembled drip chamber.

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TECHNICAL FEATURES

	Venice & Venice V.A.L.M. Single chamber	Venice & Venice V.A.L.M. Double chamber	Venice V.A.L.M. Pediatric
Technical features			
Capacity	2500 ml		200 ml
Graduated scale in Collection chamber	Single collection chamber divided into 3 sub-chambers: The first sub-chamber is graduated in 1ml increments up to the 100 ml capacity and in 2 ml increments from the 100 ml capacity up to the maximum capacity of 200 ml. The second sub-chamber is graduated in 10 ml increments up to the maximum capacity of 1350 ml. The third sub-chamber is graduated in 10ml increments up to the maximum capacity of 2500 ml.	Double collection chambers. The first chamber is divided into two sub-chamber: The first sub chamber is graduated in 1ml increments up to the 100 ml capacity and in 2 ml increments from the 100 m capacity up to the maximum capacity of 200 ml. The second sub-chamber is graduated in 10 ml increments up to the maximum capacity of 1350 ml. The second chamber is graduated in 10ml increments up to the maximum capacity of 1150 ml.	Single collection chamber. The collection chamber is graduated in 1ml increments up to the 100 ml capacity and in 2 ml increments from the 100 ml capacity up to the maximum capacity of 200 ml.
Suction Control	From -5 to -25 cm H_20		
Maximum negative pressure	-60 cm H_20		
Body size	Height (excluding tubing): 360 mm Depth: 76 mm Width: 270 mm		
	Venice & Venice V.A.L.M. Single chamber	**Venice & Venice V.A.L.M. Double chamber**	**Venice V.A.L.M. Pediatric**
Lines and connections			
Drainage connection line	Single	Double	Single
Connection to the catheter	Connector for thoracic drainage catheters with an inside diameter (ID) between 6,36 mm (1/4") and 12,7 mm (1/2").		Connector for thoracic drainage catheters with an inside diameter (ID) between 3,2 mm (1/8") and 11,1 mm (7/16").
Connection to vacuum source	Connector for a tube of connection to the vacuum source with an inside diameter (ID) between 4.8 mm (3/16") and 8.7 mm (11/32") is required.		

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INSTALLATION

Warning

- Carefully inspect the protective wrapping of the device. If it is opened or damaged, do not use the device.

1. Open the packaging using aseptic techniques and hand the contents to the sterile operator.
2. In the sterile field, remove the inner wrapping, open the floor stand and start with:

a) "Water seal" chamber filling

Using either the funnel provided or a common syringe, inject 45 ml of sterile saline solution through the "vacuum source connection line". The water will turn light blue.

Warnings

- Do not fill above a level of 2 cm.
- Remove any excess water using a syringe fitted with a needle through the pierceable port located at the base of the "Water Seal" chamber.

b) "Suction Control" chamber filling

Remove the "holed cap" and pour in water until the desired suction pressure level is reached. Close the "holed cap".

Warning

- Do not obstruct the holes in the cap when using suction pressures $\leq$ - 25 cm H_2O.

c) Patient connection

Connect the drainage catheter(s) with the drainage collection lines using the connector, if provided, or a suitable connector.

Warning

- In VENICE and VENICE Volumetric Air leak Meter (V.A.L.M.) "Dual collection" version there are two "drainage connection lines", one ending with red cap and connected with the first collection chamber and the other ending with blue cap and connected with the second chamber.

d) Drainage activation with suction or gravity

Place the device outside the sterile field; open the floor stand and place the thoracic drainage on the floor or hang it on the patient bed using the appropriate support hooks.

For suction drainage, connect a suction source to the "Vacuum Source Connection" line using a reusable tube.

Use a suction source controlled by a vacuum regulator. Once the drainage unit has been connected, slowly increase the suction pressure until gentle and constant bubbling in the "Suction Control" chamber is obtained. Small adjustments may be made to reduce bubbling by turning the "macro vacuum regulation" knob on the drainage unit.


Warning

- The drainage device must always be located below the level of the patient's wound in an upright position.

e) Blood drainage / suction operation verification

- "Suction Control": to obtain the desired suction pressure a vacuum source must be connected. The "Suction Control" chamber must be filled with water and maintained at the level prescribed by the physician over time; water must bubble gently and continuously.
- "Water Seal": to guarantee the unidirectional sealing functions the "Water Seal" chamber must be filled with water and maintained at a level of 2 cm over time. It is an indication of correct operation if, at least at the moment of activation of the suction, there is bubbling in the "Water seal" chamber.

f1) Detection of air leaks from the circuit and/ or patient via the instantaneous air leakage detector (monitoring can be performed using the VENICE model or the VENICE Volumetric Air Leak Meter (V.A.L.M) model with the stopcock in "BUBBLING MODE")


If bubbling at the base of the "Water Seal" chamber persists after suction is turned on, it is an indication of air suction from the thoracic cavity (the frequency and intensity of bubbling are proportional to the sucked air quantity).

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Warning

- Bubbling may be caused by accidental air suction through the drainage catheter connections with the "drainage connection line". To check this, clamp the drains upstream of the connection. If bubbling continues, tighten the connections; otherwise the air is sucked by the pleural-mediastinal cavity.

f2) Detection of air leaks from the circuit and/ or patient via the cumulative/instantaneous air leakage detector (monitoring can only be performed using the VENICE Volumetric Air Leak Meter (V.A.L.M) model with the stopcock turned to the "VOLUMETRIC MODE")

When the stopcock is turned into the "VOLUMETRIC MODE" position (lever horizontal) the air drained from the thoracic cavity is progressively collected in the central part of the "Water Seal" chamber, causing the water column to move from the center to the left part of the seal chamber.
The "cumulative/volumetric detection of air leakage" can be performed more than once by returning the dedicated stopcock to the "BUBBLING MODE" position (accumulation is zeroed).
The "Water Seal" maintains its functionality during "VOLUMETRIC MODE" operation. Anyway, once the "cumulative/volumetric detection" is complete, return the stopcock to the "BUBBLING MODE" position.
Air leaks can be measured with VENICE Volumetric Air Leak Meter (V.A.L.M.) both under controlled suction and by gravity.

Warning

- In order to obtain reliable results, when Volumetric Air Leak Meter is reading do not vary the suction pressure of the drainage device and do not switch from suction to gravity drainage and viceversa.

VENICE Volumetric Air Leak Meter can also be used to quantify the air flow drained from the thoracic cavity in a time unit using a stopwatch and the conversion scale on the label applied to the device floor stand.

1. When the stopcock is turned in "BUBBLING MODE", check that the water in the "Water Seal" chamber is at the "START/ 2 cm"; if not, add or remove water as needed using the pierceable port located on the back of the device.
2. Start the stopwatch while turning the tap from the "BUBBLING MODE" position to the "VOLUMETRIC MODE" position.
3. Stop the stopwatch when the water reaches the "END" level on the left part and/or if the central part of the "Water Seal" chamber empties completely.
4. Based on the time measured, the conversion scale will provide and indicative value of the drained air per time unit.

Warning

- The value of drained air per time unit, is purely indicative and cannot be used as the only element for drawing diagnostic conclusions.

g) Detection of the patient's intrathoracic pressure
The patient's intrathoracic pressure is detected by the level of the water column in the right part of the "Water Seal" chamber. If the device is operating under suction, the intrathoracic pressure is given by the algebraic sum of the value set on the "Suction Control" chamber added to the value read on the right part of the "Water Seal" chamber [for example, this addition would give the following result: (-20 cm H_2O) + (+2 cm H_2O) = -18 cm H_2O]. If, on the other hand, the device is gravity-fed (drainage by gravity), the intrathoracic pressure is equivalent to that of the water level in the right-hand portion of the "Water Seal" chamber.

h) Automatic floating discharge valve (high negative pressure valve)
The "Water Seal" features a maximum negative pressure valve; this valve comes into operation when in the "Water Seal" chamber the water column reaches the full scale (-20 cm H_2O) and prevents depression increases beyond this value. The cause of high negative pressure can be maneuvers such as squeezing (milking) the tubes (peristaltic pump effect).

i) Restoring the physiological intrathoracic pressure
When pulling the "sliding clamp" to the left, air

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inlets into the collection chamber, bringing the pressure inside this chamber to atmospheric zero ("Water Seal" column down to zero).
This operation does not in any way place the thoracic cavity in communication with the atmosphere since the "tube-clamp plate" system in one single operation allows:
1) opening the vent, letting air in the collection chamber through the septum;
2) clamping the drainage connection line(s);
3) closing the vacuum source, if active.
When pulling the "sliding clamp" to the right, the previous drainage mode is restored, by either suction or gravity.

Warning
- In particular, under gravity drainage, it may be necessary to move the "tube-clamp plate" repeatedly in order to bring the intrathoracic pressure close to atmospheric zero.

j) Automatic positive pressure valve
The positive pressure relief valve, located in the upper part of the drainage in connection with the "Water seal" chamber, opens automatically to release any positive pressure that may develop under gravity drainage.

k) Water levelling in the "Suction Control" and "Water Seal" chambers
Position the drainage/suction system in "Standby" by pulling the "sliding clamp" to the left.
- It is possible to add water in the "Suction Control" chamber through the dedicated hole plugged with the holed cap. To remove water, insert a needle through the pierceable port on the back of the "Suction Control" chamber and aspirate with a syringe.
- It is possible to add or remove water in the "Water Seal" chamber: insert a needle into the pierceable port of the chamber and with a syringe, aspirate or insert water.

Once the leveling operations have been carried out, restore the normal drainage functions by pulling the "sliding clamp" to the right.

l) Sampling of fluid drained from the thoracic cavity
In order to take a sample of fluid, with a suitable aseptic technique, insert a needle in the pierceable port on the back of the first collection chamber and drawn a sample with a syringe.

m) Use of suction pressure above -25 cm H_2O
If it is necessary to apply to the thoracic cavity a suction pressure greater than -25 cm H_2O, occlude the holes on the "holed cap" using plasticine or non-porous adhesive tape and set the desired vacuum through a regulator and a wall-mounted pressure gauge.

Caution
- The device has been designed to support vacuums not greater than -60 cm H_2O.

Warning in case of difficulty in evacuating air and fluids from thoracic cavity:
a) Check and return the "sliding clamp" to the drainage/aspiration position (pulled to the right);
b) Check and, if necessary, eliminate the loops containing the drained fluid on the "drainage connection line" (e.g. a loop of 20 cm high counteracts a suction pressure of -20 cm H_2O);
c) Check if drains and/or drainage connection line(s) are clogged by clots and, if necessary, free them by squeezing (milking) the connection drainage line or the drains.
d) Check if the quantity of water in the "Water Seal" chamber is at the prescribed level; if necessary, restore the level. Repeated overturning of the device may cause water to flow from the "Suction Control" chamber to the "Water Seal" chamber. Each cm over the prescribed level counteracts suction by one cm of water column, or creates a resistance of one cm of water column for gravity drainage.
e) Check if collection chamber is completely full of drained fluid; if so, replace the device.
f) If using suction drainage, check the water level in the "Suction Control" chamber and, if necessary, restore the desired level. Also, check if there is constant bubbling in the chamber. If not, increase the suction source intensity.
Bubbling in the "Suction Control" chamber may fail for different reasons: a) a drop in the suction source flow rate as a result of clamping a tube or problems with the centralized system, b) an increase in air drained from the thoracic cavity of the patient.
g) If using under gravity drainage, check that the maximum positive pressure valve is functioning (the ball must be free to move inside its housing) and disconnect the device from the vacuum source tube, completely unscrewing the "macro regulation" knob; this will allow any positive pressure from the valve to be released as well as from the tube left open.

n) Disconnection of the system
Clamp the "Drainage Connection Line" and all thoracic catheters before disconnecting the

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thoracic drainage from the patient.

DEVICE REPLACEMENT

During drainage procedure, a spare thoracic drainage device must always be available. Replace the device when the collection chamber is completely full with drained fluid or if the thoracic drainage is damaged.

WARRANTY CONDITIONS

This limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.
EUROSETS guarantees that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.
EUROSETS guarantees that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.
However, EUROSETS cannot guarantee that the user will use the device correctly, nor that an incorrect diagnosis or therapy and/or that particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.
EUROSETS, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from improper use of this device.
EUROSETS will replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by EUROSETS up to the time of delivery to the end user unless such defect has been caused by mishandling by the purchaser.
The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose.
No person, including any representative, agent, dealer, distributor or intermediary of EUROSETS or any other industrial or commercial organisation is authorized to make representation or warranty concerning this medical device except as expressly stated herein.
EUROSETS disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein.
The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with EUROSETS, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.
The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and Jurisdiction. The court chosen is the Court of Modena (Italy).

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