KAN TRACHEOSTOMY TUBES

INSTRUCTIONS FOR USE

Pictograms


REF	catalogue number		manufacturer
LOT	serial No.	STERILE EO	Sterilized by ethylene oxide
	date of manufacture	LATEX	latex-free
	use by		consult instructions for use
CE 1011	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment.		Caution: read all warnings and precautions in instructions for use
PHT DEHP	free from phthalates (incl. DEHP)		keep dry
	do not use if packaging is damaged		keep away from sunlight
2 STERILIZE	do not resterilize	40°C 5°C	store at 5-40 °C

Introduction
The instructions are intended for doctors, nursing staff and the patient (hereinafter referred to as the users), in order to ensure safe and correct use of KAN tracheostomy tubes both in healthcare facilities and in domestic conditions.

Read the instruction for use before first use of the device!
The packaging of the tracheostomy tube must be kept for the entire period the tube is used, due to the identification data and the device data on the label.
Users have to be trained by professional medical staff on the use of the device.

KAN tracheostomy tubes and KAN tracheostomy short tubes from SUMI are intended for use only by one patient.
Multiple uses of the device for one and the same patient are allowed.
The maximum usage time of one tube will not exceed 29 days.
It is recommended that the tubes be replaced and cleaned frequently, especially internal cannulas. Preferably by medical staff.
If damages, cracks or deformation occur during handling, replace the tube or damaged component.
Damaged device components or tubes must not be used.
In order to avoid confusion with tube selection and easier identification, each tube size is marked with a suitable color print on the flange. The device is latex free.
The device in its original, undamaged packaging is sterilized by ethylene oxide. The device must not be resterilized.

This instruction for use applies only to the following sets of tracheostomy tubes:
KAN tracheostomy tube without cuff, REF 47-xx00, 47-xx01
KAN tracheostomy tube without cuff, fenestrated, REF 47-xx10, 47-xx11
KAN tracheostomy tube short, REF 47-xx20, 47-xx21
KAN tracheostomy tube short, fenestrated, REF 47-xx30, REF 47-xx31
KAN tracheostomy tube with cuff, REF 48-xx00, REF 48-xx01
KAN tracheostomy tube with cuff, fenestrated, REF 48-xx10, REF 48-xx11
KAN tracheostomy tube with suction lumen, REF 49-xx10,

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Composition of sets

REF	KAN tracheostomy tube	KAN tracheostomy tube, fenestrated	KAN tracheostomy tube with cuff	KAN tracheostomy tube with cuff, fenestrated	KAN tracheostomy tube with suction lumen	KAN tracheostomy tube short	KAN tracheostomy tube short, fenestrated	Inner cannula with collar	Inner cannula with collar, fenestrated	Inner cannula with 15 mm connector	Inner cannula with 15 mm connector, fenestrated	Obturator	Tube holder	Connecting plug	15 mm connector	Protective plug	Decannulation cap	Speaking valve
47-XX00	1							2		1		1	1	1	1	1		
47-XX01	1							2				1	1	1				
47-XX10		1							1	1	1	1	1	1	1	1	1	1
47-XX11		1							2			1	1	1				1
48-XX00			1					1		2		1	1	1	1	1		
48-XX01			1							2		1	1					
48-XX10				1					1	1	1	1	1	1	1	1	1	1
47-XX20						1		2		1		1	1	1	1	1		
47-XX21						1		2					1	1				
47-XX30							1		1	1	1	1	1	1	1	1	1	1
47-XX31							1		2			1	1					1
49-XX10					1			1		2		1	1	1	1	1		

Indications
The KAN tracheostomy tubes are intended for patients after tracheotomy.
Contraindications
Do not use if the patient is allergic to the materials used.
Attention!!!
For mechanical ventilation, it is recommended to use tubes with cuff.
Fenestrated tubes should be used during ventilation only after consultation with a doctor.
Tubes with a speaking valve must not be used in patients after laryngectomy (without a larynx).
Product description
The set of sterile KAN tracheostomy tubes and short tracheostomy tubes contains the following elements:

Drawing	Name	Application and properties
	KAN tracheostomy tube without cuff	Outer tracheostomy tube made of medical quality PVC
	KAN tracheostomy tube short	Tracheostomy tube short is a version of the KAN tracheostomy tube without cuff with a smaller total length.
	KAN tracheostomy tube with cuff	KAN tracheostomy tube with a low-pressure cuff
	Inner cannula with collar	To be used when there is no need to use a 15 mm connector
	Inner cannula with 15 mm connector	To be used when the use of a 15 mm connector is needed (e.g. for respirator or HME filter connection)
	Obturator	The obturator closes the distal end of the tube, ensuring non-traumatic positioning of the tube in the stoma.
	15 mm connector	Enables the tube to be connected to the respirator without an inner cannula
	Protective plug	Protects the tube against aspiration of foreign bodies.
	Connecting plug.	For fixation of the inner cannula
	Speaking valve	Allows the patient to talk
	Tube holder	It provides comfortable attachment of tube flange to the patient's neck.
	Red decannulation cap	It is used to weaning the patient from using the tube.

AKTUALNY
data 230822
Director Representative
for the Quality System
Zbigniew Milewski

Instructions for use

Tube insertion

The tracheotomy procedure should be performed by a surgeon or a laryngologist.
The use of sterile dust-free disposable gloves is recommended during the entire procedure.

1. Remove appropriately sterile tracheostomy tube from intact sterile packaging.
2. Ensure that the set is complete and functions properly.
3. Check that no damage, crack or separation of individual components occurred during transport and storage. Damaged tubes must not be used.
4. Check if the connector matches the respiratory tube, using a cannula with a 15 mm connector or just the 15 mm connector.
5. Insert the obturator into the tube and check if it is easily removed from the tube.
6. Insert the cannula into the tube and check if it is easily removed from the tube.
7. In the case of tubes with cuff, check the cuff and the pilot balloon and the valves of each tube by filling the cuff, using a suitable instrument, with a test pressure of 50 cm H_2O and observing it for 1 minute. After establishing the cuff is airtight, empty it completely.
8. Before inserting the tube into stoma, check the position of the obturator inside the tube, so that its end closes the distal end of the tube. Keep the obturator in the correct position using a thumb during the entire tube insertion process. The obturator is also locked in the tube with a ring.
9. After inserting the tube into the trachea and verifying its correct positioning, remove the obturator.
10. Insert the cannula with the 15 mm connector in the tube and attach it with the connecting plug. The cannula may also be omitted, in order not to reduce the inner diameter of the tube, in which case use only the 15 mm connector on the pipe collar, securing it with a locking ring.
11. If required by the procedure, connect the tube to the respirator.
12. Ensure that the tube rests in the stoma without undue excessive tension.
13. Inflate the cuff to a minimum pressure that ensures a seal. Under typical conditions it is in the range of 20 cm H_2O to 30 cm H_2O.
14. After inflating the cuff disconnect the inflation instrument from the valve.
15. Place a sterile gauze dressing or a pad under the tube flange.
16. Use the neck strap to attach the tube to the patient's neck. The neck straps must be tied tight enough so that the tube does not fall, and loose enough that the neck is not compressed.

Attention!!!

The cuff must not be overinflated, as this may damage the delicate walls of the trachea.
The cuff pressure typically changes over a longer period, due to the diffusion of gases through the cuff walls. It is therefore recommended to monitor the cuff pressure on a regular basis.

Removing the secretions

Should viscous secretions collect in the internal cannula, and it is not possible to suction them or remove them by coughing, then you should remove the inner cannula, insert the spare cannula, and then clean the removed cannula and prepare it for the next use. There is no need to replace the entire tube in this case.

Maintaining open fenestrated holes

If a fenestrated tube with a replaceable cannula remains in the patient for a longer period, it is necessary to regularly check that the fenestrated holes are not obstructed by mucus. If necessary, replace the entire tube.

Speaking

In order to allow the patient to talk, a fenestrated tube with fenestrated cannula and a speaking valve should be used.

Attention!!!

Do not use fenestrated tubes with a speaking valve in patients with laryngeal stenosis, paralysis of vocal cords, tracheal stenosis, airway obstruction and respiratory tract infections.
Only conscious patients who can breathe spontaneously can use tubes with speaking valve.
Always remove the speaking valve before sleep.
In case of a strong cough or a shallow breath, remove the speaking valve to allow the patient to breathe freely.

Closing the fenestrated tube in order to weaning the patient from the tube

Attention!!! The closing of the tube with a decannulation cap or decannulation screw cap will only be carried out by trained medical personnel.
Before attempting to dishabituate the patient from using the tube, make sure that the upper respiratory tract is clear and that the patient can breathe spontaneously.
When the tube is being closed, the patient's condition must be monitored closely. The decannulation cap is placed on the cannula or on the plug with 15 mm connector.
The decannulation screw cap is screwed onto a tube equipped with a cannula with a cone.
If the patient exhibits symptoms of shallow, short breathing, it is absolutely required to open the tube.
Attention!!! The use of decannulation caps with non-fenestrated tubes is not allowed.

Removing the tube

Before removing the tube, remove the additional equipment such as the speaking valve or the protective plug, heat-moisture exchanger.
In the case of an unstable stoma a second spare tube should be prepared before removing the tube, and inserted immediately to avoid the collapse of the stoma.
It might be necessary to suction the secretions before removing the tube.
The suction is performed with a catheter introduced through the tube, while the tube is being pulled out of the trachea.
The removal of the tube should be performed with the patient's head slightly tilted back.
In the case of tubes with cuff, it is recommended to insert a suction catheter through the tracheostomy tube and then to completely remove the air from the cuff using a syringe. In this way, the accumulated secretions may be suctioned without any effort and without causing a cough. When the air is completely removed from the cuff, the tube can be safely removed.

Cleaning

1. Thoroughly wash the tube or cannula with an aqueous solution of sodium bicarbonate or an aqueous solution of Schulke Perform, in accordance with manufacturer's instructions, removing thoroughly any visible internal contamination with a single-use disposable cleaning swab, consisting of a shaft and a lint-free polyurethane sponge.
2. Contaminants which are difficult to remove should be cleaned using a soft bristle brush.

Attention!!!

The SUMI products may not be washed with brushes with sharp bristles, due to the danger of "scratching" of the tube, which may lead to increased bacterial growth on the tube walls.
The SUMI products must not be washed with commercially available household detergents, with denture cleaning products nor with alcohol.
It is not allowed to subject it to boiling or steam sterilizing (wet heat).

3. Afterwards the device should be rinsed with potable water.
4. When washing tubes with cuff, special caution must be taken when cleaning the delicate cuff. Do not allow water, cleaning agent or disinfectant penetrate into the filling valve or the cuff itself.
5. Thoroughly dry the tube and the cannula. Do not wipe with cloths and any other materials leaving lint.
6. Check if the device is not damaged; if damaged - dispose.

Disinfecting

The manufacturer does not recommend the use of disinfectants, because the upper respiratory tract, also in healthy humans, is not free of bacteria.
If it turns out that disinfection is necessary, then it should be carried out with, e.g. Schulke Perform or Opticit, in accordance with the manufacturer's recommendations. Afterwards the device should be rinsed with sterile saline solution.

Attention!!!

Under any circumstances do not use disinfectants releasing chlorine, containing strong alkalis or phenol derivatives.
It is forbidden to heat the tubes, cannulas and accessories above 65°C.
It is recommended that spare tubes and sterile internal cannulas are available at all times in their original packaging or washed, stored in suitable packaging to prevent contamination.
Spare tubes and internal cannulas must always be kept clean and dry.
Residues of cleaning and disinfecting agents on the tracheostomy tube may result in irritation of the mucous membranes or other negative impact on the patient's health.
The use of unsuitable cleaning and disinfecting agents may result in the damage to the device.

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Storage

Sterile, spare tubes and currently unused washed tubes will be stored in a dry environment, protected against dust, sunlight, at temperature 5°C–40°C.

Precautions

1. The tracheostomy patient may use the tracheostomy tube only when it has been prescribed by a doctor.
2. The patient's doctor selects the correct tube shape and size for the patient and decides on what additional accessories may be used.
3. The insertion of the tracheostomy tube through the tracheostoma to the trachea always increases air resistance (because of the position of the tracheostomy tube in the trachea and its small inner diameter).
4. If the dyspnea persists, the patient should immediately call a nurse or a doctor. In most cases strong coughing of the secretion or suctioning of the respiratory tract is effective.
5. Do not, in any case, plug the end of the tube with a speaking valve or a sealing plug, or a similar object without ensuring that the patient can breathe spontaneously.
6. The tube closure with the plug to weaning the patient from using the tube can only be carried out under medical supervision and at the doctor's responsibility.
7. In any case the use of set components is not allowed if the body shaft is damaged (e.g. surface cracks) and when the connection between the collar and the body loosened or separated.
8. Avoid touching the tubes with electrosurgery electrodes or laser beam, as PVC will release toxic fumes or cause ignition.
9. When the tracheostomy tube is to be used outside the hospital, the doctor or nurse should instruct the user on the safe use of the device.
10. If an allergic reaction occurs, contact a doctor immediately.
11. Never wet the pads with alcohol-based disinfectants.
12. If the patient inserts the tube by himself or herself, the procedure may be facilitated by performing it in front of a mirror.

Warning:

Due to the danger of respiratory tract infections, current approved medical procedures and techniques must be observed.
Observing current approved medical techniques and having the procedure performed by qualified medical personnel reduces the risk of complications such as: tracheal and laryngeal stenosis, wound bleeding or infection, pneumothorax, mediastinal emphysema, tracheoesophageal fistula, blocking of the tracheostomy tube, inflammation and decubitus ulcer in the area of trachea, air leakage around the tracheostomy tube, respiratory tract infection, breathing difficulties and breathing effort, muscle lesions, tracheal hemorrhage.

Remarks

1. Do not use if the package is damaged.
2. Do not use after expiration date.
3. The device is not intended for resterilization.
4. If the tube is used correctly and in accordance with the instructions, the injury of the mucous membrane or cartilage may be significantly limited.
5. Any modification of the tubes, especially shortening, perforation and repairs of the tubes by the patient or user are unacceptable.

Responsibility – legal guidelines

The KAN tracheostomy tubes manufactured by SUMI are subjected to physical and chemical testing and bio-compatibility testing. This allows their safe use and ensures highest quality.
Carefully read the instructions for use before using the tracheostomy tubes.
Tracheostomy tubes made by SUMI are intended for use by qualified medical staff and trained patients (users).
The SUMI company will not be liable for the consequences caused by improper selection of tube size or by the use of tubes contrary to the intended purpose and instructions for use.

In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Photographic documentation is also recommended.

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Ed. 02 01.04.2021.

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Director Representative
for the Quality System
Zbigniew Milewski

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EN

SUMINI TRACHEOSTOMY TUBES – INSTRUCTIONS FOR USE

These instructions for use refer to the following products:
3N-XXXX, 3P-XXXX, 2N-XXXX, 2P-XXXX

Sumini tracheostomy tubes are intended only for one patient.
The maximum usage time of one tube should not exceed 29 days.
It is recommended that the tubes be replaced frequently, as recommended by your doctor.

Description:
Sumini tracheostomy tubes are intended for infants and children. The tubes are made of medical grade polymer, which ensures adapting of tube and flange shape in patient's body temperature and provides visibility under X-rays. All Sumini tracheostomy tubes are equipped with 15 mm connector complying with ISO 5356 standard. Sumini tracheostomy tubes with cuff are equipped with low-pressure sealing cuff made of soft plastic that minimizes pressure on the trachea wall. Set contains a cotton neck tape which allows for patient-friendly fastening to the patient neck.

Kit components:
• SUMINI tracheostomy tubes • obturator • cotton neck tape

Indications:
Sumini tracheostomy tubes are intended for infants and children after tracheotomy surgery.

Instructions for use:
1. Prior to use make sure that the set is complete.
2. The procedure of tracheostomy should be performed by surgeon or otolaryngologist.
3. Remove a sterile tracheostomy tube from packaging.
4. In the case of tubes with cuff, check the cuff and the pilot balloon and the check valve of tube by filling the cuff before use. Connect a syringe with luer connector to check valve and fill the cuff with air.
5. After filling the cuff and establishing the cuff is airtight, deflate it completely.
6. Before inserting the tube into stoma, check the position of the obturator inside the tube, so that its end closes the distal end of the tube. Keep the obturator in the correct position using a thumb during the entire tube insertion procedure. The procedure may be facilitated by putting a small dose of lubricating water-soluble agent on the end of the tube and obturator.
7. After inserting the tube into the trachea and verifying its correct positioning, remove the obturator.
8. Ensure that the tube rests in the stoma without undue excessive tension.
9. If tube has a cuff, inflate the cuff with a minimal amount of air that will cause an effective sealing.
10. Place a sterile gauze dressing or a pad under the tube's flange.
11. Use the neck strap tied to a tube's flange to attach the tube to the patient's neck. The neck straps must be tied tight enough so that the tube does not fall, and loose enough so that the neck is not compressed.

Precautions:
1. The tracheostomy patient may use the tracheostomy tube only when it has been prescribed by a doctor.
2. The patient's doctor selects the correct tube shape and size for the patient and decides on what additional accessories may be used.
3. The insertion of the tracheostomy tube through the tracheostoma to the trachea always increases air resistance (because of the position of the tracheostomy tube in the trachea and its small inner diameter).
4. If the dyspnea persists, the patient should immediately call a nurse or a doctor. In most cases strong coughing to clear secretions out or suctioning of the respiratory tract is effective.
5. Under no circumstances the end of the tracheostomy tube should be occluded, for example with speaking valve, decannulation cap or a similar object.
6. Protect the cuff from damage due to contact with sharp edges.
7. Avoid touching the tubes with electrosurgery electrodes or laser beam, as PVC will release toxic fumes or cause ignition.
8. Do not use during the Magnetic Resonance Imaging (MRI).
9. Under no circumstances you should continue to use kit components if the tracheostomy tube or the cuff has been damaged (e.g. cracks on the surface, a leaky cuff) or when there is loosening of connections between the flange and the tube, or if the flange separates from the tube.
10. If an allergic reaction occurs, contact a doctor immediately.
11. When the Sumini tracheostomy tube is to be used outside the hospital, the doctor or nurse should instruct the patient's parent on the safe use of the device.

Warnings:
1. Due to the danger of respiratory tract infections, current approved medical procedures and techniques must be observed.
2. The device is not intended for resterilization
3. Observing current approved medical techniques and having the procedure performed by qualified medical personnel reduces the risk of complications such as: tracheal and laryngeal stenosis, wound bleeding or infection, pneumothorax, mediastinal emphysema, tracheoesophageal fistula, blocking of the tracheostomy tube, inflammation and decubitus ulcer in the area of trachea, air leakage around the tracheostomy tube, respiratory tract infection, breathing difficulties and breathing effort, muscle lesions, tracheal hemorrhage.
4. Do not use after expiration date or if the package is damaged.

Storage:
Sterile tubes should be stored in a dry environment, protected against dust, sunlight, at temperature 5–40°C.

Responsibility – legal guidelines:
The KAN tracheostomy tubes manufactured by SUMI are subjected to physical testing, chemical testing and biocompatibility testing. This allows for their safe use and ensures highest quality. Read the instructions for use carefully before using the tracheostomy tubes. The SUMI company will not be liable for the consequences caused by improper selection of tube size or by the use of tubes contrary to the intended purpose and instructions for use.
In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Photographic documentation is also recommended.

REF	catalogue number		manufacturer
LOT	batch no.	STERILE EO	sterilized using ethylene oxide
	date of manufacture	LATEX	latex-free
	use-by date		consult instructions for use
CE 1011	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment		caution consult precautions and warnings in instructions for use
PHT DEHP	phthalates free (DEHP free)		keep dry
	do not use if package is damaged		keep away from sunlight
2 STERILIZE	do not resterilize	40°C 5°C	store at 5-40 °C

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Director of Operations
for the Quality System
Zbigniew Milewski

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EN

CRICOTHYROTOMY SET – INSTRUCTIONS FOR USE

These instructions for use refer to product REF 30-0001

Description:
Cricothyrotomy is an invasive technique used in the emergency conditions as a quick and immediate manner of restoring patency of the airways, which have been subjected to the airway obstruction at the height or above the rima glottidis. The laryngological procedure consists in an incision of the cricothyroid ligament, located between an inferior margin of lamina of the larynx thyroid cartilage and a superior margin of the arch of the larynx cricoid cartilage.
Cricothyrotomy is a procedure of choice if there are the emergency indications in case of a lack of the possibilities of a ventilation or a conventional intubation of the patient, used when all other methods of the airway management failed to be effective.
Cricothyrotomy set is a single-use device.

Kit components:
• Tracheostomy tube size 6.0 • Bougie introducer • 10 ml syringe
• Scalpel N°10 • Strap for tracheostomy tube

Indications:
Inability to ventilate and oxygenate the patient in life-threatening conditions.

Instructions for use:
Before use, ensure whether a set is complete. Cricothyrotomy should be performed by an anaesthesiologist familiar with the procedure. It is recommended to use the pictorial instruction included in the set during the procedure.

Warning!
During the whole procedure you should remain calm and try to oxygenate the patient.

1A. Impalpable cricothyroid membrane

- Stand on the left side of the patient
- Place a pillow under the patient's shoulders and position his head to extend neck.
- Using left hand to stabilize larynx, tense the place of incision between your fingers
- Make a 8-10cm vertical incision, caudad to cephalad
- Use blunt dissection with fingers of both hands to separate tissues until you identify the cricothyroid membrane

1B. Impalpable cricothyroid membrane

- Stand on the left side of the patient
- Place a pillow under the patient's shoulders and position his head to extend neck.
- Identify the cricothyroid membrane with left hand

2. Using left hand to stabilize larynx, make a transverse stab incision through cricothyroid membrane
3. Turn blade through 90° (sharp edge caudally)
4. Hold the scalpel with your left hand. Insert the bougie introducer to 10-15cm, caudad .
5. Slide the tracheostomy tube along the guide, confirm the correct placement of the tube and secure its position with the strap.

Precautions:
1. Under no circumstances the end of the tracheostomy tube should be occluded, for example with speaking valve, decannulation cap or a similar object.
2. Protect the cuff from damage due to contact with sharp edges
3. Avoid a contact with the electrosurgical electrodes or laser radiation, because PVC will emit the toxic fumes or will cause an ignition of the air enriched with oxygen.
4. The medical device is not proper for the patients during the radiotherapy.
5. Do not use during the Magnetic Resonance Imaging (MRI).
6. Under no circumstances you should continue to use kit components if the tracheostomy tube or the cuff has been damaged (e.g. cracks on the surface, a leaky cuff) or when there is loosening of connections between the flange and the tube, or if the flange separates from the tube.

Warnings:
1. Because of a risk of infection of airways with some hospital strains that are treatment resistant, the current approved medical procedures and techniques should be observed.
2. The tube is exclusively for single use.
3. Observance of the current approved medical techniques and performing of a procedure by qualified medical personnel reduces down to a minimum possibility of complications e.g.: stenosis of larynx, trachea, bleeding or infection of a wound mediastinal emphysema, tracheal and esophageal fistula, occlusion of the tracheostomy tube, infection and decubitus ulcer within a windpipe , leak of air around the tracheostomy tube, infection of respiratory system, dyspnoea and respiratory effort, inflam[illegible] of trachea and bronchial cells, haemorrhages from larynx.

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Notes:
1. The set should be stored in 5-40°C, in a dark, dry place.
2. Do not use in if package is damaged.
3. Do not use after the expiry date.

Responsibility:
Products manufactured by SUMI undergo physical, biological and chemical tests. This fact allows to use them safely and guarantees top quality. It is recommended to thoroughly acquaint the INSTRUCTIONS FOR USE prior to use of the cricothyrotomy set. SUMI shall not bear responsibility for consequences caused by using the cricothyrotomy set contrary to their intended use or instructions for use. In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Photographic documentation is also recommended.

Pictograms:

REF	Catalogue number	LATEX	Latex free
LOT	Batch number		Prior to use consult instructions for use
	Date of manufacture		Do not re-use
	Use-by date		Do not use if package is damaged
CE 1011	Product compliant with European Union requirements / number of a notified body participating in the evaluation of device conformity	PHT DEHP	Phthalates free (DEHP free)
	Manufacturer		Keep dry
STERILE EO	Sterilized using ethylene oxide		Keep away from sunlight

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(refers to: ref 45-0093, ref 45-0094)

Pictograms

REF	catalogue number		manufacturer	PHT DEHP	free from phthalates (incl. DEHP)
LOT	serial No.	STERILE EO	Sterilized by ethylene oxide		do not use if packaging is damaged
	date of manufacture	LATEX	latex-free	2	Do not re-use
	use by		consult instructions for use		keep dry
CE 1011	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment.		Caution: read all warnings and precautions in instructions for use		keep away from sunlight
				40°C 5°C	store at 5-40 °C

Assortment
Tracheostomy pads are produced in following sizes:
ref 45-0093 8,0cm x 9,0cm ± 1,0cm
ref 45-0094 6,5cm x 6,5cm ± 1,0cm

Physical characteristics:
Density (Net weight) 77,9 kg/m^3
Thickness 5 mm
Dimension of cells 208µm
Retention 62%
Absorption time 7 sec.
Length increase after absorption 20,5%
Width increase after absorption 18,7%
Thickness increase after absorption 21,1%
Retention capacity under pressure 12.5 gr/gr
Retention capacity under pressure of 4,5 mm 4.392 gr/m^2
Absorption capacity 1600%± 200
MVTR: 1300 ± 200 g/m² in 24 h

Indications
1. In case of decubitus ulcer and to prevent it.
2. Prevention of leaking secretions spreading.
3. Protection of skin around stoma opening from irritation by a tracheostomy tube.

Instructions for use
1. Rinse wound in conformance with local procedure.
2. Prepare and clean skin around tracheostomy tube by removal of all impurities.
3. Prior to use of product, remove separating film.
4. Adjust product to the stoma opening.
5. Place product in place with the semi-absorbent layer to the patient neck.
6. In case of need, product may be cut or formed with use of sterile scissors.
7. During the initial period of curing, the stoma opening should be nurtured often.
8. Pad should be changed according to condition of the stoma opening, once per week at minimum.

Precautions
1. Pad should be used by medical staff and for treatment purposes, only.
2. The product is sterile. Open package right before use.
3. Do not use the product when packing is damaged or opened.
4. Change pad once product is soaked in 90% towards the edges.
5. In case of a haemorrhage use the product after it is stopped
6. Do not use the product in case of symptoms of a wound infection or the surrounding area

Contraindications
Cease usage immediately after notice symptoms of allergy or infection.
Use proper treatment in accordance with medical staff indications.

Warnings
Due to the danger of respiratory tract infections, current approved medical procedures and techniques must be observed.

Notices
1. Do not use in case of damaged package.
2. Do not use after the expiry date.
3. Store in room temperature, protect from high temperature and humidity.
4. Protect product from direct sunlight.
5. Do not re-use (single use product). Dispose after single use.

Responsibility
Tracheostomy pads produced by SUMI are subjected to physical, biological and chemical tests. It allows for safe use of the product and guarantees highest quality. It is recommended to read instruction for use carefully, before use of the pads.
Tracheostomy pads manufactured by SUMI may be used by qualified medical staff, only. In case of a complaint pads and its immediate packaging must be kept.

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HEAT AND MOISTURE EXCHANGER – INSTRUCTIONS FOR USE
(refers to ref 45-0050, 45-0055, 45-0051, 45-0053, 45-0054, 45-0056, 45-0057)

Description
Exchanger consists of plastic housing with two or one PU foam filters, oxygen port, and a safety valve.

Application
The exchangers are designated for spontaneously breathing patients.
It heats and moistures the inhaled air.
The heat and moisture exchanger capture heat and moisture from exhaled air and return it to the patient with the inhaled air. It prevents formation of secretions inside of trachea and reduce the risk of infection. Exchanger filters breathed air, thus preventing particles from getting inside the airway.

Technical features (in accordance with DIN EN ISO 9360-2)
Heat and moisture exchanger 45-0050, 45-0055
Tidal volume (VT) min 500ml
Internal volume 12 ml
Humidity loss 13.9 mg/l

Pressure drop
0,5 l/s 0.40 cm H_2O
1,0 l/s 0, 90 cm H_2O
1,5 l/s 1,50 cm H_2O

Weight approx. 7.0 g
Length 33 mm
Connection 15mm connector in accordance with ISO, for connecting tubes, oxygen port used for connecting oxygen

Heat and moisture exchanger 45-0051, 45-0053, 45-0054, 45-0056, 45-0057
length 22,0 mm (38.0 mm z with connecting plug)
diameter 24,0 mm
weight 2,8 g
inner volume. 3,8 ml
humidity loss 9,7 mg/l
Pressure drop
Inspiration 1,0 l/s 0,30 cm H2O
2,0 l/s 0,69 cm H2O
3,0 l/s 1,12 cm H2O
Expiration 1,0 l/s 0,34 cm H2O
2,0 l/s 0,78 cm H2O
3,0 l/s 1,28 cm H2O
Connection 15mm connector in accordance with ISO, for connecting tubes, oxygen port used for connecting oxygen

Instructions for use
Remove sterile exchanger out of its package and connect it to tracheostomy or tracheal tube, through its 15mm connector

Attention!
In case of use of a MINI exchanger designated for KAN tubes, connect the exchanger with a connecting plug, and then mount it on the tracheostomy tube and secure it in place by twisting the locking system.
Place special care on correct connection of exchanger with a tube.

Precautions
Do not use heat and moisture exchanger together with hot water humidifier, sprayer and do not pour water into the exchanger due to risk of resistance increase.

Warnings
1. Due to the risk of cross infection of airways medical procedures should be followed, Re-sterilization and reuse may pose risk for the user.
2. Product is single use, only.
3. Heat and moisture exchanger may be used for up to 24 hours.
4. In case of secretions overflow, exchanger should be replaced more often.

Notes
1. Store product in room conditions.
2. Do not use if the package is damaged.
3. Do not use after the expiry date.

Responsibility
The exchangers produced by SUMI Company undergo physical, biological and chemical tests. This allows them to be used safely and guarantees top quality.
It is recommended to thoroughly understand the instructions for use prior to use of the exchangers.
SUMI shall not bear responsibility for consequences caused by inappropriate selection and use of exchangers.
In the event of a claim one should keep the product together with direct packaging.

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EN **FEEDING TUBE - INSTRUCTIONS FOR USE**

This instructions for use refers to the following medical devices:
50-XXXX, 51-XXXX, 52-XXXX,

Pictograms

REF	catalogue number		manufacturer	PHT DEHP	free from phthalates (incl. DEHP)
LOT	serial No.	STERILE EO	sterilized by ethylene oxide		do not use if packaging is damaged
	date of manufacture	LATEX	latex-free	2	do not re-use
	use by		consult instructions for use		keep dry
CE 1011	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment.	⚠	caution: read all warnings and precautions in instructions for use		keep away from sunlight
				40°C 5°C	store at 5-40 °C

Description
The Catheters for nasal feeding is of the form of a transparent tube of circular section, On the one side of it there is stuck in a cap, the colour of which specifies a size of the catheter. The other end of the catheter is smoothly melted. On the side surface thereof, there are made two alternate holes. Catheters can be equipped with a green line providing contract in X-rays. There may be visible on the drain a code of length of the form of single marks or of each cm.

Indication
Administering of drugs, nutrition of unconscious patients, feeding of infants.

Precautions
Due to a danger of the upper respiratory tract infection mainly with the drug resistant strains in hospitals you should adhere to the actual accepted medical procedures and techniques. Re-sterilisation or reuse may constitute a thread for a patient. Observing the actual accepted medical techniques and performing a medical intervention by the qualified medical personnel reduce to a minimum a possibility of an occurrence of the complications such as: Inserting the tube to the respiratory tract, rhinitis and sinusitis, pneumonia, sore throat and earache, dry and ruptured oral mucous membrane, necrosis of nose, oesophagus, stomach soft tissues, haemorrhage.

Warnings
The medical device is intended only for its single use.

Comments
Store the packages with the medical devices in the room conditions, in a temperature between 5-40˚C
Do not use when a package is damaged.
Do not use after the expiry date.

Responsibility:
Products manufactured by SUMI undergo physical, biological and chemical tests. This fact allows to use them safely and guarantees top quality. It is recommended to thoroughly acquaint the INSTRUCTIONS FOR USE prior to use of the feeding tube. SUMI shall not bear responsibility for consequences caused by using the feeding tube contrary to their intended use or instructions for use. In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Photographic documentation is also recommended.

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