IFU10407/30012021

Evermine 50™ — Instructions for use - English

Everolimus Eluting Coronary Stent System

1. Device Description:
The Evermine 50™ Everolimus Eluting Coronary Stent System comprises of following components-
- A balloon-expandable L605 Cobalt Chromium Coronary Stent
- A stent coating that consists of a blend of anti-proliferative drug and polymers
 1. Immunosuppressant drug Everolimus
 2. Bio-compatible, bio-degradable co-polymer coating which acts as drug reservoir and drug release platform
- A rapid-exchange stent delivery PTCA balloon catheter
- The stent is pre mounted on balloon catheter & placed between two platinum-iridium radio opaque markers bands.
- Evermine 50™ is permanent implantation device.

1.1. Device Components Description:

1.1.1 Available stent lengths & diameters:
Available stent lengths & diameters (87 configurations) are shown in table-1 below

Table-1 Evermine 50™ Stent System Product Description

Available Stent Diameters (mm)	Available Stent Lengths (mm)										
	8	13	16	19	24	29	32	37	40	44	48
2.00	-	EVF20013	EVF20016	EVF20019	EVF20024	EVF20029	EVF20032	EVF20037	EVF20040	EVF20044	EVF20048
2.25	EVF22508	EVF22513	EVF22516	EVF22519	EVF22524	EVF22529	EVF22532	EVF22537	EVF22540	EVF22544	EVF22548
2.50	EVF25008	EVF25013	EVF25016	EVF25019	EVF25024	EVF25029	EVF25032	EVF25037	EVF25040	EVF25044	EVF25048
2.75	EVF27508	EVF27513	EVF27516	EVF27519	EVF27524	EVF27529	EVF27532	EVF27537	EVF27540	EVF27544	EVF27548
3.00	EVF30008	EVF30013	EVF30016	EVF30019	EVF30024	EVF30029	EVF30032	EVF30037	EVF30040	EVF30044	EVF30048
3.50	EVF35008	EVF35013	EVF35016	EVF35019	EVF35024	EVF35029	EVF35032	EVF35037	EVF35040	EVF35044	EVF35048
4.00	EVF40008	EVF40013	EVF40016	EVF40019	EVF40024	EVF40029	EVF40032	EVF40037	EVF40040	EVF40044	EVF40048
4.50	EVF45008	EVF45013	EVF45016	EVF45019	EVF45024	EVF45029	EVF45032	EVF45037	EVF45040	EVF45044	EVF45048


Table-2

1.1.2	Stent Material	Electropolished L605 Cobalt Chromium alloy, laser-cut from seamless tubing in a hybrid design pattern.
1.1.3	Stent delivery balloon catheter system	Name of Delivery System: Xpedient™ Rx PTCA Balloon Dilatation Catheter - UNS Lineage Semi-compliant Polyamide balloon, nominally 0.5 mm longer both side than the stent length. Mounted stent length & location is defined by two platinum-iridium swaged radioopaque markers under the balloon catheter. Two proximal delivery system shaft markers (90 cm and 100 cm proximal to distal tip) indicate the relative position of the delivery system to the end of brachial or femoral guiding catheter.
1.1.4	Delivery system usable length	142 cm
1.1.5	Guide wire lumen	Starts at the distal tip of the balloon catheter & ends approximately 25 cm from distal tip of the balloon catheter
1.1.6	Guide-wire rapid exchange (Rx) port	Starts at the distal tip of the balloon catheter emerges approximately 25cm from distal tip of the balloon catheter. A disposable stylet protects the distal catheter from an inadvertent kinking
1.1.7	Shaft outer profile	Proximal 2.13F Distal 2.7F
1.1.8	Stent dilatation / Balloon inflation pressures	Nominal Pressure: 8 atm, Rated Burst Pressure: 14/16 atm* RBP* 14 atm for 4.50 mm diameter, all lengths. RBP* 14 atm for all diameter with length > 40mm.
1.1.9	Guide catheter compatibility	5F (Min I.D. 0.056" / 1.42 mm)
1.1.10	Guide wire compatibility	0.014" (0.36 mm)
1.1.11	Product code format EVFxxxy	EVF = Evermine 50 xxx = nominal stent diameter (mm) yy = nominal stent length (mm) For e.g. EVF30016 xxx = 300 = Diameter 3.00 mm yy = 16 = Length 16 mm

1.2. Drug Component Description:
The drug component is coated on the stent. This coating consists of a blend of Everolimus drug (the active ingredient), and biodegradable polymers (the inactive ingredient).

1.2.1 Everolimus:
Everolimus is an active pharmaceutical ingredient. It is a novel semisynthetic macrolide immunosuppressant, synthesized by chemical modification of rapamycin (Sirolimus). The everolimus chemical name is 40-O-(2-hydroxyethyl)-rapamycin and the chemical structure is shown in Figure below:

Fig.1 Everolimus Drug Chemical Structure

The drug content on Evermine50™-Everolimus Eluting Coronary Stent System ranges 38.6 µg to 363.9 µg.

1.2.2 Polymer:
The inactive ingredient of the coating consists of a blend of Lactide and Glycolide based biodegradable polymers. These polymers control the drug release kinetics and they degrade as the drug is released from the stent.

2. How Supplied:
Sterile: This device is sterilized with Ethylene oxide (ETO) gas and is non-pyrogenic. It is intended for single use only. Do not resterilize. Do not use the device if the package is opened or damaged.
Contents: One (1) Evermine 50™ Everolimus Eluting Coronary Stent System housed in a protective circular hoop tray, one (1) Instructions for Use, two (2) Stent Implant Card.
Storage: Store between 15-25°C (59-77°F) temperature in a dry, dark, cool place. Protect from light.
Shelf Life: 3 years.

3. Indications:
The Evermine 50™ Everolimus Eluting Coronary Stent System is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease due to de novo & in-stent restenotic lesions (lengths ≤ 44 mm) in native coronary artery with a reference vessel diameter of 2.00 mm to 4.5 mm in patients eligible for Percutaneous Transluminal Coronary Angioplasty (PTCA) and Stenting procedures.

4. Contraindications:
Evermine 50™ Everolimus Eluting Coronary Stent System is contraindicated in the following patient types:
- Patients with hypersensitivity or allergies to aspirin, heparin, clopidogrel, ticlopidine, drug such as Everolimus or similar drugs or any analogue or derivative, cobalt, chromium, nickel, molybdenum, tungsten or any contrast media.
- Patients in whom anti-platelet and/or anti coagulant therapy are contraindicated.
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon.
- Transplant patients

5. Warnings:
- Judicious patient selection is necessary during use of this device since it carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events.
- Safety and effectiveness of stenting saphenous vein grafts has not been established.
- Never try to straighten a kinked hypotube. Straightening of a kinked metal may result in breakage of the shaft.
- Effect of multiple Stenting is not evaluated independently. Decision of multiple Stenting is at the physician's discretion.

6. Precautions:

6.1. General Precautions:
- Only physicians who have received adequate training should perform implantation of the stent.
- Stent placements should only be performed at hospitals where emergency coronary artery bypass graft surgery (CABG) is readily available.
- Subsequent blockage may require repeat dilatation of the arterial segment containing the stent. The long term outcome following repeat dilatation of the endothelialized stents is not well characterized.

6.2. Stent Handling Precautions:
- Do not use if the package has been opened or damaged.
- Use the device before the "Use By" date as specified on the product label.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.
- Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Remove the protective stylet from the guide wire lumen and discard.
- Do not remove the stent from the delivery system as removal may damage the stent and / or lead to stent embolization. The Evermine 50™ Everolimus Eluting Coronary Stent System is intended to perform as a system.
- The stent should not be removed for use in conjunction with other dilatation catheters.
- Special care must be taken not to handle or in any way disrupt the stent position on the delivery device. This is especially important during catheter removal from the packaging, placement of the guide wire, advancement through the rotating haemostatic valve adaptor and guiding catheter hub.
- Do not manipulate, touch or handle the stent with fingers or contact with liquids prior to the preparation and delivery as this may result in coating damage, contamination or dislodgement of the stent from the delivery balloon catheter.
- Do not expose or wipe the device with organic solvents such as alcohol or detergents.
- Use only the appropriate balloon inflation media. Do not use any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment, provide adequate support to shaft segments.

6.3. Stent Placement Precautions:
- Do not prepare or pre-inflate the balloon prior to stent deployment, other than as directed.
- Do not induce vacuum (negative pressure) on the delivery balloon catheter before reaching the target lesion.
- Implantation of a stent may lead to dissection of the vessel distal and / or proximal to the stented portion and may cause acute closure of the vessel requiring additional intervention (e.g. CABG, further dilatation or placement of additional stents).
- Do not expand the stent if it is not properly positioned in the vessel.
- Long term outcome following repeat dilatation of endothelialized coronary stents is unknown at present.
- Placement of stents has the potential to compromise side branch patency.
- Do not exceed Rated Burst Pressure as indicated on labelling. Use of pressures higher than those specified on the product label may result in a ruptured balloon and potential trauma to the vessel wall.
- Guiding catheters used must have lumen sizes that are suitable to accommodate the introduction of Evermine 50™ Stent (see Table 2).
- Stent retrieval methods (use of additional wires, snares or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma or pseudoaneurysm.
- To avoid the possibility of disimilar metal corrosion, do not implant the stents of different material in tandem or overlap or contact the stents.
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the risk of dislodging the proximal stent.
- The safety and effectiveness of Evermine 50™ coronary stent in patient with prior brachytherapy of the target lesion have not been established.
- The safety and effectiveness of using mechanical artherectomy devices or laser angioplasty catheters in conjunction with Evermine 50™ Everolimus Eluting Coronary Stent implantation have not been established.
- The Tyvek pouch is the sterile barrier. Therefore only the contents of the sealed Tyvek pouch should be considered sterile. Do not remove the contents from the Tyvek pouch until immediately prior to use.
- During withdrawal of the delivery system, hold saline-soaked gauze around the exposed catheter shaft and pull the catheter through the gauze to remove any excess contrast media.
- If reinserting the catheter, flush the guide wire lumen using the flushing needle before insertion.
- Additional expansion of the deployed stent may cause a flow limiting dissection. This may be treated by implantation of another stent. When multiple stents are implanted, the ends should overlap slightly.

6.4. Stent/System Removal Precautions:
- Should any unusual resistance be felt at any time during either lesion access or removal of the stent delivery system pre-stent implantation, the entire system must be removed as a single unit.
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer

11.3. Preparation:

11.3.1Guidewire Lumen Flush:
- Remove the protective stylet from the guide wire lumen and discard.
- Flush the guide wire with HepNS until the fluid exits the guide wire exit port approximately 25 cm distal to catheter distal tip.

Caution: Avoid manipulation of stent during flushing of guide wire lumen, as this may disrupt the placement of the stent on the balloon.

11.3.2Delivery system Preparation:
- Prepare an inflation device with diluted contrast medium
- Attach inflation device to stopcock; attach to hub (balloon inflation port)

Caution: Do not apply negative or positive pressure to balloon at this time.
- With the tip down orient the delivery system vertically
- Open stopcock to stent delivery system
- Leave inflation device on neutral
- Purge the inflation device of all air

11.3.3Delivery Procedure:
- Prepare vascular access site according to standard practice.
- Prepare lesion site according to standard practice. Predilate the lesion with a PTCA balloon.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See Stent/System Removal Precautions, section 6.4 for specific stent delivery system removal instructions.
- Tighten the rotating haemostatic valve. The stent is now ready to be deployed.

11.3.4Deployment Procedure:

Caution: Refer to product label for in-vitro stent inner diameter and RBP.
- Before deployment, reconfirm the correct position of the stent relative to target lesion via the radiopaque balloon markers.
- Attach the inflation device (only partially filled with contrast media) to a stopcock (three-way minimum) and apply negative pressure to purge the balloon of air.
- Turn the stopcock to the off position to the catheter and purge the inflation device of air. Close the side port of the stopcock.
- Under fluoroscopic visualization, inflate the balloon to deploy the stent but do not exceed the labelled stent burst pressure. Maintain inflation pressure for up to 30 seconds to fully unfold the stent against the vessel wall, with the stent internal diameter matching the size of the reference vessel diameter. Stent wall contact should be verified through routine angiography or intravascular ultrasound.
- Deflate the balloon by pulling a vacuum with the inflation device. Make sure the balloon is fully deflated before any attempted movement of the catheter.
- Confirm adequate stent expansion by angiographic injection through to the guiding catheter.

11.3.5Further Dilatation of the Stented Segments:
All efforts should be taken to assure that the stent is not under dilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. If the initial angiographic results are suboptimal, the stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent.

Note: Post dilatation is recommended for stent length > 40 mm

Caution: Do not dilate the stent beyond the following limits

Nominal Stent Diameter	Dilatation Limits
2.00 mm - 2.50 mm	3.00 mm
2.75 mm - 3.50 mm	4.00 mm
4.00 mm - 4.50 mm	5.00 mm

11.3.6Removal Procedure:
- Ensure that the balloon is fully deflated
- While maintaining slight tension and negative pressure on the inflation device, withdraw the stent delivery system.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See Stent/System Removal Precautions, section 6.4 for specific stent delivery system removal instructions.
- Repeat angiography to assess the stented area. If an adequate expansion as not been obtained, exchange back to the original delivery catheter or another proper sized balloon to the vessel wall.
- Final internal stent diameter should match reference vessel. **ASSURE THAT THE STENT IS NOT UNDERDILATED.**

11.3.7Disposal of used devices:
Dispose devices used during this procedure as per local regulatory requirements for medical device waste disposal.

12. Antiplatelet Regimen:
Physician should use the information from the current Drug Eluting Stent literature, guideline and expert opinion to select the antiplatelet/anticoagulation regime to be used for their patients in general practice.
Current guidelines for the DAPT discontinuation should be followed and are recommended. The decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition. In case an unanticipated interruption or discontinuation of DAPT is required any time after one month following DES coronary stent implantation, data from published literature show low stent thrombosis rates and no observed increased risk for stent thrombosis.

It is very important that the patient is compliant with the post procedure antiplatelet recommendation. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infraction or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventionalist and patient should carefully consider whether a Drug Eluting Stent and its associated recommended antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risk and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant bleeding should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

13. Disclaimer of Warranty and Limitation of Remedy:
There is no express or implied warranty, including without limitation any implied warranty of merchantability of fitness for a particular purpose, on the Meril Life Sciences Pvt. Ltd. product(s) described in this publication. Under no circumstances shall Meril Life Sciences Pvt. Ltd. be liable for any direct, indirect, incidental or consequential damages. Some jurisdictions do not allow the exclusion of or limitations on an implied warranty. For jurisdictions which do not allow the exclusion or limitations of incidental or consequential damages, some of the above exclusions may not apply. The patient may also have other rights, which vary from jurisdiction to jurisdiction.
Descriptions or specifications in Meril Life Sciences Pvt. Ltd. printed matter, including this publications, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.
Meril Life Sciences Pvt. Ltd. will not be responsible for any direct, indirect, incidental, or consequential damages resulting from use of the product.

Symbols used in labeling

Stent inner diameter		Stent length	
Temperature Limitation		Contains one unit	
MGCID		Keep dry	
Min guide catheter I.D.			
Do not use if package is damaged		Serial number	SN
Reference number	REF	Lot number	LOT
Consult instructions for use		Max. guide wire diameter	
Do not resterilize		Non-pyrogenic	
Sterilized using ethylene oxide	STERILE EO	Use-by date	
Do not re-use		Caution	
Manufacturer		Date of manufacture	
Keep away from sunlight		MR Conditional	

7. Adverse Effects:
Undesirable effects/adverse events (in alphabetical order) that may be associated with the implantation of a coronary stent in native coronary arteries include but are not limited to:
- Abrupt stent closure
- Acute myocardial infarction
- Allergic reactions
- Aneurysm
- Angina
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT)
- Arteriovenous fistula
- Cardio tamponade
- Coronary Artery Occlusion
- Cardiogenic shock
- Death
- Dissection
- Drug reactions to antiplatelet agents / anticoagulation agents / contrast media
- Emboli, distal (air, tissue or thrombotic emboli)
- Embolization, stent
- Emergency Coronary Artery Bypass Graft Surgery (CABG)
- Failure to deliver the stent at the intended site
- Fever
- Fistulization
- Heart Failure
- Hematoma
- Hemorrhage
- Hypotension / Hypertension
- Incomplete Stent Apposition
- Infection, including infection and/or pain at the access site
- Myocardial Infarction
- Myocardial Ischemia
- Perforation or rupture
- Pericardial effusion
- Prolong Angina
- Pseudoaneurysm
- Renal failure
- Respiratory failure
- Restenosis of stented segment
- Rupture of native and bypass graft
- Shock / Pulmonary edema
- Spasm
- Stent compression
- Stent migration
- Stroke / cerebrovascular accident / TIA
- Stent thrombosis (acute, subacute, or late)/occlusion
- Ventricular fibrillation
- Vessel perforation
- Vessel spasm
- Vessel trauma requiring surgical repair or reintervention

Potential adverse events, not captured above, that may be related to Everolimus following oral administration:
- Abdominal pain (including upper abdominal pain)
- Anemia
- Angioedema
- Anorexia
- Arthralgias
- Diarrhea
- Constipation
- Cough
- Delayed wound healing/fluid accumulation
- Dyspnea
- Dysgeusia
- Dysuria
- Dry skin
- Edema (peripheral)
- Epistaxis
- Fatigue
- Headache
- Hematuria
- Hypercholesterolemia
- Hypersensitivity, including anaphylactic/anaphylactoid type reactions
- Hypertriglyceridemia
- Hyperlipidemia
- Hyperkalemia
- Hypertension
- Hypokalemia
- Hypomagnesemia
- Hypophosphatemia
- Increased serum creatinine
- Infections
- Interstitial lung disease
- Insomnia
- Interaction with strong inhibitors and inducers of CYP3A4 or PgP
- Leukopenia
- Lymphoma and other malignancies
- Male infertility (azospermia and / or oligospermia)
- Mucosal inflammation (including oral ulceration and oral mucositis)
- Myalgia
- Neutropenia
- Non-infectious pneumonitis
- Pain: extremity, incision site and procedural, and back
- Proteinuria
- Pruritus
- Pyrexia
- Rash
- Stomatitis
- Thrombotic microangiopathy (TMA) / Thrombotic thrombocytopenic purpura (TTP) /Hemolytic uremic syndrome (HUS)
- Tremor
- Thrombocytopenia
- Urinary tract infection
- Vomiting

There may be other potential adverse events that are unforeseen at this time.

8. Recommended Drug Regimen:
Antiplatelet or anticoagulant therapy is recommended as per institional practices for coronary stenting.

9. Individualization of Treatment:
The risk and benefits should be considered for each patient before use of Evermine 50™ Everolimus Eluting Coronary Stent System. Patient selection factors should include a judgement regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active or gastritis or peptic ulcers disease.
Pre-morbid conditions that increase the risk of a poor initial result or the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure and severe obesity), should be reviewed.
A review of the vessel location, reference vessel size, lesion length, qualitative target lesion characteristics, and the amount of myocardium in jeopardy from acute or subacute thrombosis must also be considered.
Thrombosis following stent implantation is affected by several baseline angiographic and procedural factors. These include vessel diameter less than 3 mm, intra-procedural thrombus and dissection following stent implantation. In patients who have undergone coronary stenting, the persistence of a thrombus or dissection should be considered a marker for subsequent thrombotic occlusion. These patients should be monitored very carefully during the first month after stent implantation.

10. Use in Special Populations:
The safety and effectiveness of the Evermine 50™ Everolimus Eluting Coronary Stent System has not been established in the following patient populations:
- Patient with unresolved vessel thrombus at the lesion site
- Patient with coronary artery reference vessel diameter < 2.00 mm.
- Patients with unprotected left main coronary stenosis in the left main coronary artery.
- Patients with brachytherapy treatment of the target lesion

Pregnant Patients: There are no adequate and well controlled studies in pregnant women or men intending to father children. Effective contraception should be initiated before implanting Evermine 50™ Everolimus Eluting Coronary Stent System and for 12 weeks after implantation. The Evermine 50™ Everolimus Eluting Coronary Stent System should be used during pregnancy only if the potential benefit outweighs the potential risk to the embryo or fetus.

Lactation: It is not known whether Everolimus is distributed in human breast milk. Because similar drugs are known to be excreted in human milk, and because of the risk of adverse reactions in nursing infants, a decision should be made whether to discontinue nursing infants, a decision should be made whether to discontinue nursing or implant the stent, taking into account the importance of the stent to the mother.
The safety and effectiveness of using brachytherapy treatment, mechanical artherectomy devices (directional artherectomy catheters, rotational artherectomy catheters) or laser angioplasty catheters to treat in stent stenosis of a Evermine 50™ Everolimus Eluting Coronary Stent + Everolimus Eluting Coronary Stent has not yet been established.
The safety and effectiveness of Evermine 50™ has not been evaluated in pediatric subjects below 18 years.

11. Clinical Use Information:

11.1. Inspection prior to use:
- Carefully inspect the sterile package before opening
- Do not use if the package has been damaged or opened
- The product should not be used after the "Use By" Date
- If the sterile package appears intact, carefully remove the system from the package and inspect for bends, kinks and other damage.
- Tear open the sterile pouch to carefully remove the product and pass on or drop the contents into the sterile field using aseptic technique.
- Verify that the stent is located between the radiopaque markers
- Do not use if any defects are noted

11.2. Materials Required:
- Appropriate guiding catheter(s)
- 2-3 syringes (10-20 cc)
- 1000 u/500 cc Normal heparinised saline (HepNS)
- 0.014" (0.36 mm) diameter guide wire, 175 cm minimum length
- Rotating haemostatic valve with an appropriate internal diameter (if applicable).
- Contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer

11.3. Preparation:

11.3.1Guidewire Lumen Flush:
- Remove the protective stylet from the guide wire lumen and discard.
- Flush the guide wire with HepNS until the fluid exits the guide wire exit port approximately 25 cm distal to catheter distal tip.

Caution: Avoid manipulation of stent during flushing of guide wire lumen, as this may disrupt the placement of the stent on the balloon.

11.3.2Delivery system Preparation:
- Prepare an inflation device with diluted contrast medium
- Attach inflation device to stopcock; attach to hub (balloon inflation port)

Caution: Do not apply negative or positive pressure to balloon at this time.
- With the tip down orient the delivery system vertically
- Open stopcock to stent delivery system
- Leave inflation device on neutral
- Purge the inflation device of all air

11.3.3Delivery Procedure:
- Prepare vascular access site according to standard practice.
- Prepare lesion site according to standard practice. Predilate the lesion with a PTCA balloon.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See Stent/System Removal Precautions, section 6.4 for specific stent delivery system removal instructions.
- Tighten the rotating haemostatic valve. The stent is now ready to be deployed.

11.3.4Deployment Procedure:

Caution: Refer to product label for in-vitro stent inner diameter and RBP.
- Before deployment, reconfirm the correct position of the stent relative to target lesion via the radiopaque balloon markers.
- Attach the inflation device (only partially filled with contrast media) to a stopcock (three-way minimum) and apply negative pressure to purge the balloon of air.
- Turn the stopcock to the off position to the catheter and purge the inflation device of air. Close the side port of the stopcock.
- Under fluoroscopic visualization, inflate the balloon to deploy the stent but do not exceed the labelled stent burst pressure. Maintain inflation pressure for up to 30 seconds to fully unfold the stent against the vessel wall, with the stent internal diameter matching the size of the reference vessel diameter. Stent wall contact should be verified through routine angiography or intravascular ultrasound.
- Deflate the balloon by pulling a vacuum with the inflation device. Make sure the balloon is fully deflated before any attempted movement of the catheter.
- Confirm adequate stent expansion by angiographic injection through to the guiding catheter.

11.3.5Further Dilatation of the Stented Segments:
All efforts should be taken to assure that the stent is not under dilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. If the initial angiographic results are suboptimal, the stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent.

Note: Post dilatation is recommended for stent length > 40 mm

Caution: Do not dilate the stent beyond the following limits

Nominal Stent Diameter	Dilatation Limits
2.00 mm - 2.50 mm	3.00 mm
2.75 mm - 3.50 mm	4.00 mm
4.00 mm - 4.50 mm	5.00 mm

11.3.6Removal Procedure:
- Ensure that the balloon is fully deflated
- While maintaining slight tension and negative pressure on the inflation device, withdraw the stent delivery system.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See Stent/System Removal Precautions, section 6.4 for specific stent delivery system removal instructions.
- Repeat angiography to assess the stented area. If an adequate expansion as not been obtained, exchange back to the original delivery catheter or another proper sized balloon to the vessel wall.
- Final internal stent diameter should match reference vessel. **ASSURE THAT THE STENT IS NOT UNDERDILATED.**

11.3.7Disposal of used devices:
Dispose devices used during this procedure as per local regulatory requirements for medical device waste disposal.

12. Antiplatelet Regimen:
Physician should use the information from the current Drug Eluting Stent literature, guideline and expert opinion to select the antiplatelet/anticoagulation regime to be used for their patients in general practice.
Current guidelines for the DAPT discontinuation should be followed and are recommended. The decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition. In case an unanticipated interruption or discontinuation of DAPT is required any time after one month following DES coronary stent implantation, data from published literature show low stent thrombosis rates and no observed increased risk for stent thrombosis.

It is very important that the patient is compliant with the post procedure antiplatelet recommendation. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infraction or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventionalist and patient should carefully consider whether a Drug Eluting Stent and its associated recommended antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risk and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant bleeding should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

13. Disclaimer of Warranty and Limitation of Remedy:
There is no express or implied warranty, including without limitation any implied warranty of merchantability of fitness for a particular purpose, on the Meril Life Sciences Pvt. Ltd. product(s) described in this publication. Under no circumstances shall Meril Life Sciences Pvt. Ltd. be liable for any direct, indirect, incidental or consequential damages. Some jurisdictions do not allow the exclusion of or limitations on an implied warranty. For jurisdictions which do not allow the exclusion or limitations of incidental or consequential damages, some of the above exclusions may not apply. The patient may also have other rights, which vary from jurisdiction to jurisdiction.
Descriptions or specifications in Meril Life Sciences Pvt. Ltd. printed matter, including this publications, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.
Meril Life Sciences Pvt. Ltd. will not be responsible for any direct, indirect, incidental, or consequential damages resulting from use of the product.

Symbols used in labeling

Stent inner diameter		Stent length	
Temperature Limitation		Contains one unit	
MGCID		Keep dry	
Min guide catheter I.D.			
Do not use if package is damaged		Serial number	SN
Reference number	REF	Lot number	LOT
Consult instructions for use		Max. guide wire diameter	
Do not resterilize		Non-pyrogenic	
Sterilized using ethylene oxide	STERILE EO	Use-by date	
Do not re-use		Caution	
Manufacturer		Date of manufacture	
Keep away from sunlight		MR Conditional	

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