Roadsaver™

INSTRUCTIONS FOR USE

TERUMO

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
English
ROADSAVER™ Carotid Artery Stent
Instructions for Use

DEVICE DESCRIPTION

The Terumo ROADSAVER Carotid Artery Stent system consists of two basic components: a self expanding stent (see Figure 2) and a 5 French Rapid Exchange (RX) delivery catheter (See Figure 1). The Stent component is a dual-layer braided closed cell design (See Figure 2). The stent outer layer is a woven closed cell structure with flared ends. The stent inner layer is a braided closed cell structure with micromesh sized pores. Both the outer and inner stents are composed of a nickel titanium (Nitinol) material. The stent is designed to expand to a pre-determined diameter when deployed by the delivery catheter. Upon exiting the delivery catheter at the target lesion, the implant expands to the vessel lumen diameter. The ROADSAVER system can be used with a compatible .014" wire based EPD (Embolic Protection Device) as well as .014" guidewires. The delivery system is 143 cm length.

Figure 1. ROADSAVER System Set-up

1. Knob
2. RHV
3. Guidewire Lumen Mandrel
4. Proximal RX Port
5. Distal Tip

Figure 2. ROADSAVER Stent (Unconstrained)



Figure 3. ROADSAVER Stent (Implanted/Constrained)

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INDICATIONS

The ROADSAVER Carotid Stent System is indicated for use in patients with carotid arterial atherosclerotic disease.

CONTRAINDICATIONS

The ROADSAVER Carotid Stent System is contraindicated for use in:

- Patients in whom anticoagulant, antiplatelet therapy or thrombolytic drugs is contraindicated
- Patients with known hypersensitivity to nickel-titanium
- Patients with severe vascular tortuosity or anatomy that would preclude the safe introduction of a guide catheter, sheath, embolic protection system or stent system
- Patients with uncorrected bleeding disorders
- Lesions in the ostium of the common carotid artery

WARNINGS

The ROADSAVER stent should be implanted in vessels with the referenced diameters according to the product label.

Should unusual resistance be felt at any time during access or removal, the sheath introducer/guide catheter and ROADSAVER system should be removed as a single unit. Applying excessive force during delivery or retrieval of the ROADSAVER system can potentially result in loss or damage to the device and delivery components.

If an EPD device is used, allow for and maintain adequate distance between the filter, the stent delivery system or deployed stent to avoid potential entanglement.

Do not reposition the ROADSAVER system without fully retrieving the device. The ROADSAVER stent MUST be retrieved into the delivery system and redeployed at the desired target location or removed completely from the patient.

Do not attempt to reposition the ROADSAVER implant after detachment.

The use of a guiding sheath or guiding catheter with a fixed hemostasis valve may cause the distal tip of the delivery catheter to detach at the hemostasis valve upon removal if the valve is not adequately opened.

When multiple stents are required, stent materials should be of similar composition.

Persons allergic to nickel titanium (Nitinol) may suffer an allergic reaction to this implant.

It is not recommended that the stent be used in patients with the following additional characteristics:

- Patients with poor renal function and in the opinion of physician may be at risk for a reaction to contrast medium
- Pregnant patients
- Patients with perforated vessels evidence by extravasation of contrast medium
- The presence of extensive atherosclerotic disease involving the aortic arch and proximal common carotid artery that would preclude the safe introduction of a guiding catheter, sheath or EPD

PRECAUTIONS

Only physicians trained and familiar with carotid stent placement and the associated complications, side effects and hazards should use this device.

The stent should not be sized to the vessel beyond the reference vessel diameter range listed on the product label.

Carefully inspect the sterile package and the ROADSAVER system device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components, or if the package is opened or damaged.

See the product label for shelf life. Do not use the ROADSAVER system beyond the labeled use by date.

Carefully remove the guide wire lumen mandrel (see Figure 1 for Proximal Rx port location) and fully flush the guide wire lumen. Do not use the delivery system if flush is not observed exiting at the stent Rx delivery system port locations (proximal and distal).

Exercise caution when crossing the deployed/detached ROADSAVER system with adjunctive devices such as guidewires, catheters or balloon catheters to avoid disrupting the device geometry and device placement.

The ROADSAVER system is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose in accordance with hospital, administrative and/or local government policy.

Store in a cool and dry place.

POTENTIAL COMPLICATIONS

Possible complications include but are not limited to the following:

- Arrhythmia
- Aneurysm and pseudoaneurysm formation
- Abrupt vessel closure
- Allergic reactions (including to antiplatelet agents, contrast medium or stent materials)
- Arteriovenous fistula
- Bleeding from anticoagulation/antiplatelet medication
- Bradycardia and hypotension
- Carotid artery spasm
- Coronary ischemia
- Death
- Disseminated intravascular coagulation
- Emboli (air, tissue, plaque, thrombus, device or other)
- Emergency artery bypass graft surgery
- Hematoma
- Hemorrhagic or embolic stroke/TIA
- Infection and/or pain at insertion site/Sepsis
- Intimal tear/dissection
- Myocardial Infarction (MI)
- New or worse encephalopathy
- Renal failure/insufficiency
- Respiratory arrest
- Stent misplacement
- Tissue necrosis
- Vessel injury/dissection/perforation/rupture/trauma
- Vessel occlusion or thrombosis
- Vessel spasm or recoil

MR Information

The ROADSAVER Carotid Stent has been determined to be **MR conditional** according to the terminology specified by the American Society for Testing Materials (ASTM).

Non-clinical testing demonstrated that the device is MR conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla or less.
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg for 15 minutes of scanning.
- First Level Controlled Operating Mode for the MR system.

MRI-Related Heating

In non-clinical testing, the device produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	1.5-Tesla	3-Tesla
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.8°C	+2.0°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Image Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5-mm relative to the size and shape of this implant.

Optimum MR Imaging Parameters				
Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	564 mm^2	125 mm^2	955 mm^2	155 mm^2
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Terumo Corporation recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation or equivalent organization. A ROADSAVER patient implant card is included in the package, which should be completed and provided to the patient.

CLINICIAN USE INFORMATION

Materials

The following parts are required to use the ROADSAVER system:

Accessories for performing a procedure but NOT supplied with the system should be selected based on the physician's experience and preferences:

- Appropriate Guiding sheath or catheter with a minimum internal diameter (I.D.) of .074 in. (1.9 mm)

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- .014" guide wire or embolic protection device compatible with the ROADSAVER System
- Saline solution/heparin-saline solution continuous flush set
- 5 cc syringe for flushing the ROADSAVER System
- Contrast solution
- Rotating Hemostatic Valve (RHV)
- Pressurized sterile Infusion solutions – IV stand
- Femoral artery access device, sterile needle, guide wire
- The ROADSAVER system does not contain latex or PVC materials

PACKAGING AND STORAGE

The ROADSAVER System is placed inside a protective packaged which includes a tray, sterile sealed pouch and unit carton. The ROADSAVER System will remain sterile unless the package is opened, damaged, or the expiration date has passed. Store at a controlled room temperature in a dry place.

SHELF LIFE

See the product label for the device shelf life. Do not use the device beyond the labeled use by date.

PREPARATION FOR USE

Device and Delivery System Selection

Appropriate selection of the ROADSAVER Carotid Stent System is important for patient safety. In order to choose the optimal ROADSAVER stent size for any given lesion, examine pre-treatment angiograms for correct and accurate vessel measurements. Reference the Product Specification chart for the ROADSAVER Carotid Stent System with the in-service guide Instructions (not included within IFU) or see the device box label. The stent diameter should be slightly (1-2 mm) larger than the diameter of the reference vessel to obtain secure placement. Choose the unconstrained diameter and dual layer stent length as indicated on the outer box label to correspond with the implanted dual layer length, at a minimum, as long as the target lesion.

Directions for Use

Delivery System Preparation

Note: You can prepare and flush the delivery system within the packaging tray or you can remove it and flush the delivery system outside the packaging.

1. Remove the pouched tray enclosing the ROADSAVER System from the carton.
2. Inspect the pouch for any signs of damage to the sterile barrier.

Warning: If it is suspected that the sterile barrier seal has been opened or compromised, do not use the carotid delivery system and return it to the manufacturer.

The delivery system should be flushed using a 5-ml syringe.

1. Peel open the pouch, remove the tray with delivery system from the pouch and inspect the delivery system for any signs of damage.
2. Remove the guidewire lumen mandrel from the proximal Rx port and flush the delivery system (See Figure 1 for Proximal Rx Port location).
3. Attach a 5-ml syringe filled with sterile heparinized saline solution to the RHV luer and flush the delivery system (See Figure 1 for RHV). Apply positive pressure until the saline weeps or comes out in drops from the valve seal of the RHV. Close the seal of the RHV, continue to apply positive pressure until the saline weeps or comes out in drops from the proximal RX port (See Figure 1). Pinch the proximal RX port and continue to apply positive pressure until the saline weeps or comes out in drops from the distal RX port or distal tip of the delivery catheter.
4. Carefully remove the delivery system from the tray and remove the catheter from the tubing within the tray without kinking the catheter.

Pre-procedure

The use of an Embolic Protection Device (EPD) should be considered. Ensure that the EPD is compatible with the ROADSAVER Carotid Delivery System and moves freely within the wire delivery section of the catheter. Gain vascular access according to standard angiographic practice and perform a diagnostic angiogram to determine the appropriate stent diameter and length for the target lesion. Refer to the Product Specification chart for the ROADSAVER Carotid Stent System with the in-service guide Instructions (not included within IFU) or see the device box label for choosing the stent size.

Introduction of the Stent Delivery System

Refer to the outer box label for choosing the proper guide sheath or guiding catheter for use with the ROADSAVER Carotid Stent System.

Access the treatment site using the appropriate accessory equipment and insert an EPD usable as a guidewire. If an EPD is not preferred, use a sufficiently long .014" guidewire and place it across the lesion.

Warning: Always use a guidewire or EPD when advancing or withdrawing the delivery system. If necessary perform a pre-dilatation at the lesion target site using a standard PTA technique using the recommended guide sheath or guiding catheter on the outer box label. The guide wire or EPD device should be advanced past the target lesion site. Then the delivery system should be advanced over the guide wire or EPD to the lesion site using fluoroscopy.

Warning: If necessary perform a pre-dilatation using a standard PTA procedure and then introduce the system through the sheath or guiding catheter.

Caution: Do not advance the delivery system against significant resistance. The cause of the resistance should be determined by fluoroscopy and remedial action taken. Withdraw the system and use a new one.

Slack Removal

Ensure that the delivery system catheter outside the patient remains flat and straight.

Warning: Slack in the delivery system catheter either outside or within the patient may result in misplacement of the stent beyond the lesion.

Stent Deployment

Note: The ROADSAVER stent should be sized to the vessel according to the reference vessel diameter range listed on the product label. The ROADSAVER stent foreshortens as it expands to the diameter of the native vessel. The stent foreshortening should be taken into account when sizing and deploying the ROADSAVER Carotid Stent System. If the ROADSAVER stent positioning is not satisfactory across the target lesion, the stent may be recaptured and repositioned if it is not fully deployed. The implant may be recaptured up to 50% of its deployment from the catheter delivery system (**See Figure 5 and 6**). Recapturing of stent can be accomplished by removing the slack in the delivery catheter shaft and then maintaining the inner catheter (handle) position and advancing or pushing forward the outer catheter. After recapturing the stent into the delivery system then you can redeploy the stent across the lesion, taking into consideration the foreshortening of the stent.

Ensure the RHV is open. The stent is deployed by outer catheter retraction. Deployment is complete by maintaining inner catheter position (holding handle) while retracting the outer catheter and allowing the stent to expand (**See Figure 4**).


Figure 4. Stent Deployment

Figure 5. Stent Retraction

Figure 6. Stent Recapture: Approximately 50% Deployment

Withdrawing Stent Delivery System

After the stent has been completely deployed, using fluoroscopy, withdraw the delivery system carefully, leaving the guidewire or EPD in place. Perform standard post-procedural angiography.

Post-deployment Stent Dilatation

If the stent is not completely expanded within the length of the lesion, post-dilation with the balloon catheter inside the stent may be done with standard PTA technique. The inflated nominal diameter of the PTA balloon used for post-dilation should be approximately the same as the diameter of the referenced native vessel.

Post Stent Placement

Upon final angiography confirming a satisfactory result, the guidewire or EPD should be removed along with the guide sheath or guiding catheter. Hemostasis of the puncture site should be established.

HOW SUPPLIED

Sterile: This device is sterilized using Ethylene Oxide. Non-pyrogenic

Contents: One (1) ROADSAVER system

Storage: Store product in a dry, cool place.

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WARRANTY DISCLAIMER

Terumo warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage, cleaning, and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure, and other matters beyond Terumo's control directly affect the device and the results obtained from its use. Terumo's obligation under this warranty is limited to the repair or replacement of this device through its expiration date. Terumo shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. Terumo neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Terumo assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications, and model availability are subject to change without notice.

© Copyright 2016 Terumo Corporation. All rights reserved.

Terumo® is a registered trademark of Terumo Corporation.

MicroVention® is a registered trademark of MicroVention, Inc.

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© Copyright 2016 Terumo Corporation. Tous droits réservés.
© Copyright 2016 Terumo Corporation. Alle Rechte vorbehalten.
© Copyright 2016 Terumo Corporation. Reservados todos los derechos.
© Copyright 2016 Terumo Corporation. Tutti i diritti riservati.
© Copyright 2016 Terumo Corporation. Todos direitos reservados.
© Copyright 2016 Terumo Corporation. Alle rettigheder forbeholdes.