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Intraoral Digital Impression Instrument

Product Specification

（manual）

Model：PANDA P3

CE

Ziyang Freqty Medical Equipment Co., Ltd.

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Contents

1 General Information....3
2 Product Information....3
3 Product Composition....3
4 Main Dimension....4
5 Intended Use and Contraindication....4
5.1 Intended Use....4
5.2 Users....5
5.3 Contraindication....5
6 Environmental Requirements....5
7 Working Power Requirements....5
8 Safety Information....6
8.1 Prerequisites....6
8.2 Mechanical Hazards....6
8.3 Explosion Hazards....7
8.4 Electrical Safety....7
8.5 Laser Safety....7
8.6 Cautions....8
9 Product Hardware Installation Instructions....8
10 Product Software Description....9
10.1 Software Operation Configuration Requirements....9
10.2 Software Basic Information....9
10.3 Software Installation Method....10
10.4 Main Software Interface....10
11 Application Method....10
11.1 Operating Steps....10
11.2 Scan technique....10
11.3 Calibration....11
12 Care and Maintenance Methods....11

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17.1 Liability of the Manufacturer..... 15

17.2 Warranty Description.....15

18 About EMC Descriptions and Risk Warning.....15

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CAUTION

Be sure to observe all warnings!

Please observe all safety information and warnings to prevent personal injury material damage or damage to your instrument.Safety information and warnings are highlighted in this guide using the words WRNING, CAUTION.

The symbols used in this document imply the following:

WARNING

Warnings regarding situations where a risk of injury to person exists if the information is not observed.

CAUTION

Safety information where hazards such as:loss of data, invalidation of warranty or service contract, risk of property damage, damage to the instrument exist if the information is not observed.

2 Product Information

Product Name: Intraoral Digital Impression Instrument

Model: PANDA P3

Manufacturer Name: Ziyang Freqty Medical Equipment Co., Ltd.

Manufacturer Address: Floor 2-3, unit 7, building 3, No. 222, West Section 3, outer ring road, Yanjiang District, Ziyang City, 641300,Sichuan Province,P.R. China

Manufacturer contact: Tel: +86-028-26577388 E-mail: sales@freqty.com

Classification: The product is energized from an external electrical power source is classified as CLASS I. The APPLIED PARTS is classified as TYPE B.


3 Product Composition

The product is composed of the probe, the power adapter, the calibrator, electric switch box, the supporting software (release version:P3V1) and the probe bracket. The probe includes the probe body and P2 probe head assembly (include normal probe, D probe and M probe).

The P2 probe head assembly is applied part.

The structure of PANDA P3 is showed as below Fig 1.

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5 Intended Use and Contraindication

5.1 Intended Use

This product uses the optical scanning method to obtain the three-dimensional shape feature data of the surfaces of teeth, gums and other tissue. It output the three-dimensional digital impression data which can be used in the CAD / CAM denture design and processing.

⚠	**WARNING** Unintended use of instrument can results in physical injury to patients, operators and damage to the product.

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5.2 Users

The product may be operated only by trained dental professionals and qualified personnel.

5.3 Contraindication

Patients with the following contraindications are not suitable for intraoral digital impressions. These contraindications include but are not limited to:

1. Patients have oral mucosal disease; patients have mental illness; patients have Parkinson's disease; patients have ADHD (Attention Deficit and Hyperactivity Disorder); patients have epilepsy.

2. Spray optical shading powder on smoked areas when patients have very severe black smoke stains that are not conducive to optical scanning. If the special shading powder need to be sprayed, these diseases need to avoid dust are contraindications, mainly including Not limited to: allergic or multi-drug allergic; severe respiratory diseases, asthma patients, etc.

3. It should not be used on patients who have or have had photobiological reactions (including those with excessive sun exposure or porphyria) or who have been treated with photosensitive drugs (including methoxsalen or chlortetracycline).

6 Environmental Requirements

Operating temperature: 5 ℃ ~ 30 ℃

Storage temperature: -10 ℃ ~ 55 ℃

Operating humidity: ≤80%

Storage humidity: ≤93%

Atmospheric pressure: 700hPa ~ 1060hPa

Operating Environment: Home healthcare environment and professional healthcare facility environment. Indoor operation, prevent direct sunlight and strong lights, and keep away from electromagnetic sources, cold and heat sources, and vibration sources

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8 Safety Information

8.1 Prerequisites

	**CAUTION** Read all instructions carefully including all warnings and cautions. You must comply with the warnings in the manual to prevent injury to persons and damage to equipment. Proper functionality and safety can only be guaranteed if the safety precautions in this safety guide and on the instrument are observed.

Preventive inspection before use of the system.

	**CAUTION** Please examine the instrument for any mechanical damage on: ● All enclosures; ● All cables Safety can only be guaranteed if NO DAMAGE on the instrument is observed.

Modification

	**WARNING** No modification of this equipment is allowed.

Approved software only

	**CAUTION** Install only approved software to prevent interference with the runtime reliability of the system and programs within it.

In case of equipment failure

	**WARNING** If at any time the instrument malfunctions, or if you suspect in any way that the instrument is not working correctly: ● Remove the scanner from contact with the patient. ● Unplug the probe and make sure it cannot be used before it is checked. ● Contact your reseller. ● DO NOT attempt to open any covers on the instrument.

8.2 Mechanical Hazards

Dropped or damaged equipment

	**CAUTION** If you drop a P2 probe head assembly on the floor, you MUST dispose of it immediately and NOT use the same assembly again for scanning. There is high risk that the mirror in the P2 probe head assembly has become dislodged and can fall out.
	**CAUTION** If the probe body is dropped or bumped it should immediately be calibrated before further use. If calibration fails, please contact your technical service provider. See this guide instructions on Calibrating the probe.

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8.3 Explosion Hazards

Environment

	**WARNING** The product is not designed to be used in environments that are potentially explosive such as in close proximity to flammable liquids or gases or in oxygen-enriched atmospheres.

Distance to the patient

	**WARNING** There is a potential explosion hazard if the product is operated in the presence of flammable anesthetic.

8.4 Electrical Safety

Ensure a grounded/earth connection

	**WARNING** To avoid the risk of electric shock, the product must only be connected to a supply mains with protective earth.

Electrical shock

	**WARNING** There is a risk of electrical if you attempt to access the inside of any part of the product. Only authorized and qualified service personnel may access the inside of any part of the system.

Stress on cables

	**CAUTION** All externally connected cables such as power cable, data cable must never be subjected to pulling stress.

Spilled Liquids

	**WARNING** Do not bring liquids such as beverages near the product. Do not spill liquids on the product.

Disconnected from mains

	**WARNING** The power ON/OFF switch is on the electric switch box.Do not position the product so that it is difficult to operate the electric switch box.

8.5 Laser Safety

This product uses a visible laser light.

	**WARNING** Do not look into the visible laser beam in the process of use, and prohibit the beam of the scan window (laser window) from directly hitting the operator and the patient's eyes.

The laser wavelengths used by the product are 450nm and 520nm.Both beam divergence angles(Parallel) are

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9°,beam divergence angle(Perpendicular) are 44° and 49°. The pulse width is 60μs. The repetition frequency is 15Hz, and the maximum power is less than 0.5mW. The warning label related laser product has been affixed to the external surface of the product which can be clearly seen. After use, please place the probe on the probe bracket with the scanning window facing down.

WARNING
The patient should wear goggles before starting scanning.

The requirements for goggles are as follows: protection wavelength 200-540nm, OD4+ (transmittance of 0.01%), visible light transmittance of 60%. Goggles should be kept properly after use, for example, put them in a glasses case.

8.6 Cautions

This product is a precision optical measuring instrument and must not be impacted during use.

This product meets the requirements for electromagnetic compatibility of medical devices in use, but it is not recommended to use it in environments with strong magnetic fields, strong switches and strong light sources, otherwise it may affect the performance of the product.

This product can only be connected to the USB interface of UL/CSA 60950-1 (or GB4943.1) certified computer equipment.

After the product is used at the end of its life, the product should be disposed of in accordance with the requirements of local laws and regulations, or contacted by the manufacturer for recycling and centralized disposal in accordance with local laws and regulations.

The product label has two parts, including the label on the body of the product and the label on the electric switch box, the serial number on the product label (e.g. IS027E50820ZY) and the number on the electric switch box label (e.g. IS027E50820ZY-1) is actually the same, but in order to distinguish the two labels, -1 is added after the number of the electric control switch box as the distinction mark, -1 is only used for the distinction between the two and has no other meaning.

9 Product Hardware Installation Instructions

This product is a precision optical equipment. Manufacturers and distributors shall not be liable for the loss of product safety and reliability and performance if the operator do not operate in accordance with the instructions, or if they do not use the product in a collision and fall due to improper use. After falling, please check the product function and calibrate the product with a calibrator. If the calibration fails, please contact the manufacturer for repair.

When the probe is not scanned, it should be placed on the probe bracket and on the horizontal operating table as a whole to avoid falling damage caused by improper placement.
Installation step:

- Connect the power supply cord to the adapter.
- Insert the adapter into the power adapter socket on the electric switch box.
- Insert the probe power supply line into the probe power socket on the electric switch box.

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- Connect the USB 3.0 of the probe data cable to the USB 3.0 port of the user's computer.
- Plug the power cord into the SUPPLY MAINS.
- Turn on the power switch on the electric switch box, and the green indicator light will turn on.
- Run the scanning software and scan according to the requirements of the scanning operation. During normal scanning, light is projected from the scanning window.
- After scanning, power should be turned off.

10 Product Software Description

10.1 Software Operation Configuration Requirements

This product can only be used by installing software on the computer. The requirement for the recommended configuration of computer hardware is no less than the following configuration:

CPU:	Intel i7-8700（3.2GHz）and above
GPU:	NVIDIA GTX1060 (6G) and above
Memory:	DDR2400 16G and above
Hard disk:	Solid state drive SSD 120G and above
Display:	Resolution 1920×1080 and above


10.2 Software Basic Information

Software Name: Intraoral Digital Impression Instrument Scan Software

Release version: P3V1

Software security level: A

Network security: User access control can choose to use user name and password for identity authentication. The user type is ordinary user. Ordinary users can use the instrument normally and view data results. The login interface is shown in the figure below:

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
Data saving format: standard STL, PLY format and PTY format defined by our company

10.3 Software Installation Method

See Intraoral Digital Impression Instrument Software Operation Manual.

10.4 Main Software Interface

The specific method of use of the software is provided by the company with training materials and operation manuals.

11 Application Method

11.1 Operating Steps

Follow the instructions for Product Hardware Installation in Chapter 9.


Open the software and scan after the startup is complete. During the scanning process, the function button can be used to control the start, pause and end of the scan.

11.2 Scan technique

Hold the probe body in the same way as a pen while scanning, Due to the limitation of the actual space in the mouth, it is necessary to ensure that the head window of the probe is as close to the tooth surface as possible (it is recommended to keep it within 2mm) for scanning, and the operation mode of suddenly far and suddenly near should be avoided.

In general, a normal probe is used for scanning, which can be used for partial or full scanning operations. For some special cases (such as insufficient mouth opening, or severe occlusion between teeth), D probe can be used to supplement the missing area of the distal surface, or the M probe can be used to supplement the missing area of the mesial surface. Please refer to the following figure to select the probe.

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Axial drag of the probe was the main scanning method, and radial drag of the probe was used in the scanning of the front teeth and occlusion points. Start scanning from the end teeth, first scan the oeclusal data, then scan part of the buccal and lingual data, and drag from oeclusal to mesial to scan the next tooth, and follow the same operation to complete the frame scan of the posterior tooth area. When entering the lingual surface of the anterior tooth area, drag the probe radial direction left and right to scan the lingual surface and incisal data, and then scan part of the labial surface data after the lingual surface is completed.

11.3 Calibration

According to the usage, it is recommended to use the calibrator to calibrate the product once a week. The product has not been used for three months, it is recommended to calibrate before use. When the device is impacted, or the product is moved or vibrated greatly, or in order to maintain the accuracy of the scanning accuracy, the scanner needs to be calibrated. Refer to "Operation Manual" for the calibration method.

⚠	**CAUTION** The calibrator of the product should be properly kept. Once the calibrator is defaced, the performance of the product will be degraded.

12 Care and Maintenance Methods

The product is not expected to have long and frequent oral contact with patients. It is recommended to sterilize the P2 probe head assembly by means of moist heat steam sterilization (121 °C, 15 min or 134 °C, 6 min).

⚠	**WARNING** In order to ensure the normal performance of the product, it is recommended that the times of repeated sterilization of the P2 probe head assembly shall not more than 20.

The P2 probe head should be replaced when the appearance is damaged or the times of sterilizations is 20. The P2 probe head assembly can be purchased separately from the seller or manufacturer.

Recommended sterilization method:

- Clean the P2 probe head assembly with soapy water and a soft brush, then place it under running water for rinsing.

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- Wipe the water stain on the surface of the P2 probe head assembly with medical gauze and wipe it thoroughly with absolute alcohol. Pay special attention to whether there are stains or water stains on the head mirror. If there is, use another medical gauze to draw the absolute alcohol and carefully wipe the head mirror. The sample was allowed to stand for two minutes after wiping.
- Place the P2 probe head assembly which had been cleaned into 90* 260mm Self-sealing sterilization pouch (materials: Medical high-temperature dialysis paper and medical CPP/PET complex film) and seal the sterilization pouch. Then place the packaged P2 probe head assembly into sterilizing instrument tray.
- Place the sterilizing instrument tray into a small pressure steam sterilizer and set the sterilization parameters according to the instructions of the small steam sterilizer: temperature 121 ° C, 15 min, or temperature 134 ° C, 6 min.

Keep the outside of the product clean.

If the probe head reflective glass smudging, can dipping a small amount of anhydrous alcohol with skimmed cotton, from the center to gently wipe the rotation. If the glass is scratched, it needs to be replaced.

In the course of equipment, software errors and warnings can be self-healed by the software. Serious problems may require shutting down the software and restarting. General hardware errors can be restored by turning off the power and then turning the power back on. If something cannot be recovered, contact the manufacturer or the seller.

The maintenance personnel must take laser protective measures during the inspection process, such as wearing goggles.

	**WARNING** During the inspection, ensure that there is no person in the direction of laser irradiation..

Replacement equipment parts must be obtained from the manufacturer or manufacturer approved dealer.

	**CAUTION** The parts that not supplied by the manufacturer may reduce the accuracy and safety of the equipment.

Disclaimer: We can provide the necessary information for maintenance equipment to the users with corresponding maintenance qualifications.

13 Transport and Storage Life

Transport conditions: temperature -10 ℃ ~ 55 ℃, relative humidity ≤93% and an atmospheric pressure 700hPa ~ 1060hPa, note rain, drop.

Storage: should be stored at ambient temperature of -10 ℃ ~ 55 ℃, relative humidity not exceed 93% and an atmospheric pressure 700hPa ~ 1060hPa, well-ventilated, non-corrosive gas chamber. Prevent moisture, corrosion, avoid direct sunlight.

Expected service life: five years. Over period of use, the degradation of the product's main electronic and optical components may reduce product performance.

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14 Parts List

Parts Name	Quantity
Probe body	1
Probe bracket	1
Calibrator	1
Power adapter	1
Electric switch box	1
Power supply cord	1
Calibrator cable	1
P2 Probe head assembly	5
Specification	1
Qualified label	1
Warranty card	1

15 Revision History

Revision	Revision content	Date
A.00	First release	2022.03.01

16 Legend of Labels and Symbols

Caution


General Warning

Type B application part

Refer to instruction manual/ booklet

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Serial number


The device should be sent to the special agencies according to local regulations for separate collection after its useful life.

CE marking in conformity with Regulation (EU) 2017/745

Indicates a medical device that needs to be protected from moisture.

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17 Other Content

17.1 Liability of the Manufacturer

The installation, adjustment, modification, and repair of this product are performed by persons or organizations approved by the manufacturer or distributor. And the manufacturer must be able to ensure the safety of the product in accordance with the electrical, environmental, storage, maintenance, and operation requirements of the manual. Responsibility for reliability and performance.

17.2 Warranty Description

The warranty period for this product is one year, calculated from the date of sale (according to the date of sales invoice).

One of the following cases the company provides free maintenance:

(1) Non-user subjective destruction within the warranty period, product failure caused only by product quality;

(2) Product failure caused by force majeure (such as earthquake, flood, typhoon, etc.) during the warranty period.

One of the following cases the company provides paid maintenance:

(1) Failure of the product due to non-subjective damage by the user during the warranty period;

(2) Failure of the product after the warranty period but within the service period.

We no longer provides maintenance in one of the following situations:

(1) The product has been used for more than six years;

(2) Product failure caused by subjective destruction of users

18 About EMC Descriptions and Risk Warning

This product has passed the electromagnetic compatibility test and meets the requirements of EN 60601-1-2 Medical electrical equipment–Part 1-2: General requirements for basic safety and essential performance–Collateral Standard: Electromagnetic disturbances – Requirements and tests.

The following application requirements shall be strictly observed during use, otherwise it may cause electromagnetic interference to other devices or reduce the anti-electromagnetic interference capability of the therapeutic device, or even lose the basic performance.

This product belongs to the Group 1 Class B equipment specified in IEC/CISPR 11, non-permanent installation equipment, non-living support equipment, and belongs to equipment that is expected to be directly connected to the public power grid.

WARNING

Portable RF communications equipment(including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm(12inches) to any part of the product including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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The cable information of this product is as shown in the following table. If there is a fault in the connection cable, please contact our company for maintenance or replacement. Other-wise it may cause excessive electromagnetic interference. If there is something wrong with this product, please contact our company promptly. Do not repair or replace the components yourself, or it may cause excessive electromagnetic interference.

No.	Name	Cable length (m)	If the shield	Remarks
1	Power supply cord	1.6	No	\
2	DC power supply line	1.5	No	With a magnetic loop
3	Data cable	2.0	Yes	\
4	Probe power supply line	0.3	No	\

	**WARNING** The use of accessories or cables outside of the regulations together with equipment and systems may result in increased emissions or reduced immunity of the equipment or system..
	**WARNING** This product should not be used near or stacked with other devices. If it must be used close to or stacked, it should be observed and verified to work properly under its configuration.

Pass and fail criteria: Free from distortion in an image or error of a displayed. In the test, the communication was normal during the continuous scanning process, and the image of the tooth model could be obtained normally.

Test method: The device is powered on, connected to the test software, set to continuous scan mode, a dental plaster model is placed on the front end of the probe for continuous scanning.

Work mode: continuous scan mode. After the device is connected to the power test software, for continuous scanning.

Trouble Shooting:

Issue	Solution
No image display in 2D image area	-Make sure the device's USB interface is properly connected to the computer's USB 3.0 interface. -Restart the software and scanning device to check if the image can be displayed normally.
2D image flicker	-Check if the modulator is connected properly. -Replace the USB port of the device with the computer. -Connect your computer to the Internet.
Scans are easily interrupted and not smooth	-Inappropriate scan brightness. For plaster model scanning, choose 1/2, for resin model scanning, choose 3, for the intraoral scanning, choose 4, 5 is suitable for patients with darker teeth in the mouth. -During scanning, confirm that A above the image area is blue. If it is black, use the keyboard A key to switch. -Standardize scanning methods. Ensure coverage of scanned data with existing data.
Out-sync of data between 2D and	-Confirm whether the computer configuration meets the requirements (higher

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3D	than or equal to our recommended configuration). -Delays caused by too many scans (single jaw scans should be completed within 3 minutes). -Uninstall antivirus software or add scanning software to the whitelist of antivirus software. -Check the status of windows update. If the update is in progress or has failed, please restart the computer after the update is completed before using the scanning software.
Difficulty for scan relocation	-Ensure that the scanning direction is consistent with the previous scanning when repositioning -Avoid long scans.
No 3D data when scanning	Recalibration
Abnormal interrupt during scanning	-Check the status of windows update. If the update is in progress or has failed, please restart the computer after the update is completed before using the scanning software. -Check whether the remaining storage space of drive C is sufficient. -Turn off or uninstall anti-virus software.

This product declarations to meet Table 1, Table 2, Table 3of Contents.

Table1

Manufacturer's Declaration - Electromagnetic Emissions	
The product intended for use in the electromagnetic environment specified below. The customers or users should ensure that it is used in such an environment.	
Emission measurement	Conformity
RF emission CISPR 11	Group 1
RF emission CISPR 11	Class B
Harmonic emission IEC 61000-3-2	Complies
voltage fluctuations / flicker emission IEC 61000-3-3	Complies

Table 2

Manufacturer's Declaration - Electromagnetic Immunity		
The product intended for use in the electromagnetic environment specified below. The customers or users should ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level guidelines	Compliance level
electrostatic discharge IEC 61000-4-2	Contact: ± 8 kV. Air: ± 2kV, ±4kV, ±8kV, ± 15 kV	Contact: ± 8 kV. Air: ± 2kV, ±4kV, ±8kV, ± 15 kV.
Conducted disturbances included by RF fields IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

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Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines.	± 2 kV for power supply lines
Surge IEC 61000-4-5	± 1 kV line(s) to line(s). ± 2 kV line(s) to earth.	± 1 kV line(s) to line(s). ± 2 kV line(s) to earth.
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Power frequency magnetic field (50Hz) IEC 61000-4-8	30A/m.	30A/m.
Power input line voltage dips, short interruptions and voltage variations IEC 61000-4-11	0% U_T, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T, 1 cycle and 70%U_T,25/30 cycles Single phase:at 0° 0% U_T, 250/300 cycles	0% U_T, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T, 1 cycle and 70%U_T,25/30 cycles Single phase:at 0° 0% U_T, 250/300 cycles

Table 3

Manufacturer's declaration-electromagnetic immunity

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

The PANDA P3 is intended for use in the electromagnetic environment (for professional healthcare) specified below. The customer or the user of the PANDA P3 should assure that it is used in such an environment.

Test frequency (MHz)	**Band** [a] (MHz)	**Service** [a]	**Modulation** [b]	**Maximum power** (W)	**Distance** (m)	**IMMUNITY TEST LEVEL** (V/m)	**Compliance LEVEL** (V/m) (For professional healthcare)
385	380-390	TETRA 400	Pulse modulation [b] 18 Hz	1,8	0,3	27	27
450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation	2	0,3	28	28
710	704-787	LTE Band 13, 17	Pulse modulation [b] 217 Hz	0,2	0,3	9	9
745							
780							
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850,	Pulse modulation [b] 18 Hz	2	0,3	28	28
870							
930							
1720	1700-	GSM 1800;	Pulse	2	0,3	28	28

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NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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Manufacturer Information

Manufacture:Ziyang Freqty Medical Equipment Co., Ltd.


Address:Floor 2-3, unit 7, building 3, No. 222, West Section 3, outer ring road, Yanjiang District, Ziyang City, 641300,Sichuan Province,P.R. China

Tel: +86-028-26577388

Email: sales@freqty.com

Website: http://www.freqty.com/

European Authorized Representative

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Ziyang Freqty Medical Equipment Co., Ltd.

Intraoral Digital Impression Instrument Software Operation Manual

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Contents

1 Equipment Introduction..... 6
1.1 Product use..... 6
1.2 Product composition..... 6
2 Hardware product..... 7
2.1 Configuration requirements..... 7
2.1.1 Computer configuration table..... 7
2.1.2 Power Supply..... 7
2.1.3 Environmental Requirements..... 7
2.2 Equipment Installation..... 7
3 Software Instructions..... 9
3.1 Installation Instructions..... 9
3.2 Software Interface..... 10
3.2.1 Management Page..... 11
3.2.2 Treatment Page..... 15
3.2.3 Acquisition..... 16
3.2.4 Processor..... 24
3.2.5 Storage..... 36
3.2.6 Setting..... 37
3.3 Shortcut key..... 42
4 Equipment sterilization and routine maintenance..... 43
4.1 Probe head sterilization..... 43
4.2 Routine maintenance..... 44
5 Software and hardware common problems and solutions..... 44
5.1 Software startup issues..... 44
5.1.1 The software cannot be opened normally..... 44
5.1.2 Software starts too slowly..... 44
5.1.3 Software icon on the desktop turns white..... 45
5.2 Problems connecting devices..... 45
5.2.1 Camera connection failed..... 45
5.3 Image display problems..... 45

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5.3.1 No image display in 2D image area....45
5.3.2 2D image flicker....45
5.4 Scanning issues....46
5.4.1 Scans are easily interrupted and not smooth....46
5.4.2 Out-sync of data between 2D and 3D....46
5.4.3 Difficulty for scan relocation....46
5.4.4 No 3D data when scanning....46
5.5 Abnormal interrupt during scanning....46
5.6 Problems with calibration....46
5.6.1 Failed to add position or calibration failed....46
5.7 Other problems....47
5.7.1 Computer restarts repeatedly....47
5.7.2 The device cannot be powered on normally....47
6 Maintenance and technical support....48

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Copyright

Opyright ©Ziyang Freqty Medical Equipment Co., Ltd. 2015-2025. All rights reserved.

Any part of this manual, including words, pictures, graphics, etc., belongs to Ziyang Freqty Medical Equipment Co., Ltd. (hereinafter referred to as "the Company" or "Freqty Technology"). No unit or individual may extract, reproduce, translate or modify all or part of this manual in any way without written permission. Unless otherwise agreed, the Company does not provide any express or implied representations or warranties to this manual.

Trademark

频泰科技 and FREQTY are registered trademarks of Ziyang Freqty Medical Equipment Co., Ltd. Other trademarks referred to in this manual are owned by their respective owners.

Statement

This manual is used as a guide. The photos, graphics and illustrations provided in the manual are for explanatory and illustrative purposes only and may differ from specific products. Due to product version upgrade or other needs, the company may update this manual, if you need the latest version of the manual, please log on the company's official website. Freqty Technology recommends that you use this manual under the guidance of professionals.

When using this product, please strictly follow the applicable law. If this product is used to infringe upon the rights of a third party or other improper use, the company will not bear any responsibility.

If the contents of this manual conflict with the applicable law, the provisions of the law shall prevail.

Safety Warning

First of all, thank you for choosing the intraoral digital impression instrument of Freqty Technology.

Please take a few minutes to read the safety warning before using it to avoid various hazards and damage to the equipment.

Do not touch equipment and equipment plugs with hands or other objects with water.

Please connect the equipment to a power supply with protective grounding.

Please avoid pulling, knotting, trampling and other operations on the cable, so as not to damage the power cord.

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When there are problems in the process of use, do not disassemble the equipment privately, please contact us in time.

After each use on the patient, the probe head must be disinfected and sterilized.

When using the equipment, please handle it gently, do not drop or smash the equipment, so as not to cause damage to the equipment.

Do not block the heat sink under the equipment.

Do not use this equipment in an overheated, wet or cold environment.

Do not touch the reflector of the probe head with other objects, so as not to dirty and damage the reflector.

Do not touch the lens with your hands or other objects, so as to avoid getting dirty or damaging the lens.

Do not shine a laser directly on eyes of any people.

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1 Equipment Introduction


1.1 Product use

This product uses the optical scanning method to obtain the three-dimensional geometric data of dentition directly, and provides the digital three-dimensional model of CAD / CAM denture design and processing, which can be used in the fields of dental repair, orthodontics and implant.

Compared with the traditional method of crown production, the Freqty impression scanner has great advantages.

- No powder: improve the comfort of patients.
- Timely feedback: the corresponding 3D data can appear in the computer immediately, so that doctors and patients can understand the oral situation at the first time.
- High precision: this product has ultra-high accuracy, the data can be directly sent to the processing factory for designing and processing.

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2.1.1 Computer configuration table

Table 1 Necessary computer configuration

Computer Parts	Configuration
CPU	i7-7700k (3.6GHz; quad-core; 8-thread) or above
RAM	DDR 2400 16G or above
Hard Disk	SSD 240G or above
Graphics Card	GTX 1060 Memory 6G or above
Operating System	Windows10 (2004 and above)
Graphics Driver Version	NVIDIA* Driver 461.92 - WHQL version and above

2.1.2 Power Supply

Table 2 External power supply requirements

Power Parameters	Requirements
Rated Voltage	AC100V-240V
Rated Current	1.1A
Rated Frequency	50Hz-60Hz

2.1.3 Environmental Requirements

- Working environment temperature: 5℃～30℃
- Working environment relative humidity: ≤80%
- Working environment: indoor, avoid any other strong light, away from strong electromagnetic interference source.
- Indoor storage conditions: -10℃～50℃, keep dry.

2.2 Equipment Installation

After removing the packaging, put the probe body❸ and probe bracket❺ on a safe and stable horizontal table. If the table is unstable or uneven, there will be a large sloshing in

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the operation of the equipment, which may cause deterioration of image data.

Insert the probe body❸ into P2 probe head assembly❶, the probe head can be replaced.

The probe body❸ and data cable❹ are integral and cannot be disassembled. If you pull out the cable or check the plug and unplug connection by yourself, the cable will be damaged.

Plug the power cord into the probe power socket ❽and turn on the power.

Plug the USB 3.0❻ into the computer's USB3.0 port (blue).


Clean the scanner regularly to keep it clean.

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3 Software Instructions


3.1 Installation Instructions

Open the installation package (>> Select "No" to overwrite the local configuration file) >> carefully read the license agreement and select agree >> click "next" according to the prompts >> finish ——the installation is complete.


(This interface will not appear if the user has never installed the company's software.)

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3.2 Software Interface

>>> Main Page

① Main menu: consist of Orders, Planning, Acquisition, Processor, Storage and Parameter Setting, etc, each of which opens the corresponding page.

② Display area: show the contents of each page.

③ Status bar: show the feedback information of each move and hints in the process of operation.

Settings button: click Settings to open the settings interface, you can set the language, format, save path, etc.

Create button: click to enter the creation interface, you can create or modify doctor-patient-related information.

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Administration button: click to enter the administration interface, you can create / modify / view doctor or patient medical record information.

Import button: click to load the medical record from the computer folder.

Practice button: you can directly enter the scan interface for exercise without creating doctor and patient information.

Recent button: click to open the recent case interface to view or open the selected case.

Minimize button: click to minimize the current software interface to a button on the taskbar.

Close button: click to close the software.


3.2.1 Management Page

>>> Case

- Method 1 for creating patients and their cases


Input doctor and patient's name and ID >> click  **to complete data >> click OK >> next step——create successfully**

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- Create case of existing patient

Click button **>> select a patient >> next step —— create new case of the patient successfully**

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
Drop-down button: display the name and ID of the patient of the last five cases.

fter selecting any existing patient in the drop-down list, if any information is modified in this interface, it is considered to create a new patient. (If it is modified in the right extension window, it is considered to modify the patient.)

>>> Manage

- Method 2 for creating patients and their cases

Click New Dentist >> fill in doctor information and apply >> fill in patient information and apply —— create successfully


The "Tab" key is to switch the edit column.

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- Open the existing medical record information

Select a doctor >> click to view the patient >> click Open —— Open successfully

① Search: search the doctor / patient information through the keyword (Doctor name / number).

② Patient basic information: including patient name, ID and creation time.

③ Patient medical record: including patient scanning scheme, doctor and creation time of this diagnosis and treatment information.

④ Detailed function area: load the two-dimensional data in the patient medical record for viewing.

Delete button: click to delete the selected doctor / patient / medical record information.

Expand button: click to view the details of the changed patient or doctor.

Open **Open button:** click to enter the processing interface to view the existing 3D data.

CAD **Design button:** view data and design in EXO CAD

Package **Package button:** package the data and open the location of the folder where is located.


Upload **Upload button:** upload the current data to the network. (This function is to be open)

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Click Import >> Select the folder of medical records to view >> open the corresponding IFTSCAN file —— successfully imported


>>> Recent

Used to view cases within the last 7 days/month/year.

Uncompleted cases after quick scanning can be automatically processed by clicking "Processing" in the "Recent" interface.

The "red dot" indicates the number of cases that need to be batched after a quick scan.

The Handle button is hidden when there are no cases to handle.

3.2.2 Treatment Page

>>> Planning selection

Select tooth number >> select scanning type >> fill in information >> click next

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——**Successful selection**

Empty **Empty button:** click to empty all the selected teeth (click the selected tooth to cancel the selection).

Orthodontics **Orthodontic mode:** just need to fill in the note.

Restoration **Restoration mode:** fill material, color, remarks. You can choose the pre-preparation scan function.

Implant **Implant mode:** fill the implanting brand, material, color and remarks.

Pre-preparation scan: after selection, the dentition before preparation and the abutment after preparation can be scanned. This procedure is only used in cases where the shape after repair is required to be consistent with that before repair.

Note: After the backup tooth scanning is selected, the occlusal displayed on the scanning image acquisition interface is the occlusal displayed before the backup tooth, and the occlusal displayed on the model processing interface is the occlusal displayed after the backup tooth


3.2.3 Acquisition

>>> Acquisition Page

- Scan the upper and lower jaw

 Click open (S) to open light >> click start (S) to scan >> click end (Enter)-scan ——completed

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1. Detailed explanation of areas marked with pictures

① 2D real-time display area: display real-time scanning in the mouth, prompting the number of scanning frames and scanning time.

There is a hint of whether the “A” intelligent scanning function is enabled：

- Blue when scanning: indicates that “A” intelligent scanning is running. Only enabled in mode 4 or 5.
- Black when scanning: indicates that “A” smart scanning is not running. You can pause the scan, then click the keyboard "A" key to enable the A function, and then start the scan.
- Gray when scanning: indicates that it cannot be used. When it exists in mode 1 to 3, the mode has no “A” intelligent function.

② Scan control area: start, pause, and end of the scan

③ Scan display area: display the data of the scan in the mouth. The operations are as follows.

- Data control: After clicking the data, the left mouse button controls the data rotation, the mouse wheel controls the data zoom, and the right mouse button controls the data translation.
- Data Supplement Scan: when the data is missing, return to the scanning interface, click Start, and then scan the missing data location.

④ Auxiliary scanning tool area: provides auxiliary tools for data inspection and modification.

⑤ Dental jaw selection area: select the area to scan the position of the corresponding jaw.


⑥ Prompt area: users can follow the prompts. If you want to close the interface, need to go to the settings interface to close it.

⑦ Brightness selection area: 1-5;1, 2 for plaster, 3 for resin, 4, 5 for real teeth.

⑧ Depth of field selection area: 1-3; Respectively for shallow, medium and deep gear position; If you need to scan more deeply, use the deep gear. Image retrieval based on scanning needs self-adjustment.

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
2、 Upper limit of scanning frames: 6000 frames are set as the upper limit for scanning single jaw. After 6000 frames, supplementary scanning is not allowed after pause or stop.

3. Out of memory:

(1)Out of memory was indicated during image scanning, as shown in Figure below.


The model needs to be processed immediately and then come back to continue the scan.

(2) Lack of memory is indicated during supplementary scanning, as shown in Figure below.

After selecting "OK", the jaw will automatically process the model and clear the cached data, and the cleared jaw will not be re-swept.

After selecting "Cancel", the current jaw scan cannot be continued

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4. Scanning music: the music plays normally during scanning, the music becomes smaller when the scanning is interrupted, and the music stops when the scanning is paused or finished

Note:

-H black smoke tooth enhanced scan mode is only used to supplement data.
After normal contour scanning, supplementary scanning should be carried out for the black teeth/smoke teeth where data cannot be swept out.
- Smoke teeth without black teeth do not require H mode.

Delete button: delete the whole data, which cannot be recovered again.

Brush button: after clicking, the mouse becomes a brush. And delete excess data as a brush.

Polygon button: after clicking, the mouse becomes a multilateral shape, and the closed area that is circled is deleted.

Lock button: after clicking, you can circle a closed area, and the circled area will not be affected during the scan.

Undo button: click once to restore the previous state of this operation, and click again to restore all.

Photo button: saves the image of the current two-dimensional imaging frame.

Color button: click to make the current 3D data switch between true color and gypsum color.

Upper jaw button: click to enter the Upper jaw scanning interface for related **operations.**

Mandibular button: click to enter the mandibular scanning interface for related operations.

Bite button: click to enter the bite scanning interface for related operations.

Swap button: click to swap the upper and lower jaw.


Next button: complete the scan operation and enter the model processing.

- Lock function

It is used to protect the area that you do not want to change during the scanning process. It's only used when scanning is paused.

Start scan >> scan pause >> click "Lock" >> circle Target area >> continue scan >> end scan

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Radio button: you can circle one area at a time, circle multiple times to save only the last selected area.

Multiple selection button: select the area with circle where all selections are saved multiple times.


Remove button: check the locked area to unlock the corresponding position.

Note:
- Function can only be used when the scanning is paused.
- Select the locked area and start scanning.
- The unlock button and the end button can release the locked area.

- Scan byte

 Click start (S) to scan >> scan bite point 1 (>> scan bite point 2) >> bite success-scan — scan finished

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Delete button: click to delete the current data.

Bite point 1: click to switch to scan bite point 1.

Bite point 2: click to switch to scan bite point 2.

Quick scan button ： It can be used when a quick scan is needed. Click it to enter the interface of creating a new case. After that, the user can conduct batch model processing in the "Recently" interface.

Note: - The quick fusion function is only supported when occlusal data is available.

- When there is no occlusal data, click the Quick Fusion button to realize the next function.

- If you need to process fast-scanning cases in batches, users can click the "Recent" button in the left navigation bar of the home page of the software. After entering the interface of recent cases, users can click "Process" to automatically process the cases.

>>> Scanning method

- Principles of operation in the mouth
- Scanning distance is constant, avoid far and near.
- In the operation of the mouth, the fulcrum is used to ensure the uniform and stable scanning speed.
- Dental scan first scans the occlusal surface to determine the dental arch.
- The current scan and the scanned data are multiplexed to confirm that the data is complete

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- Upper and lower jaw scanning:

- Recommended grip: pen-style
- Scanning distance: 2mm-5mm
- Initial scanning position: scan from the occlusal plane, parallel to the occlusal plane, and scan along the arc of the dental arch;

 Note: when you start scanning, the probe head is facing the occlusal surface. Do not scan the probe at an angle!
- Scan path: occlusal surface → tongue tip / lingual surface → cheek tip / lip surface
- **Scanning steps:**


① Place the scanning probe head directly on the posterior occlusal surface and start scanning.

② When reaching the position of tooth No. 3, the probe head is tilted toward the buccal side and scans the anterior lip data and incisor data at the same time.


③ After crossing the midline, the probe head slowly turns to the palatal side, scans the lingual / palatal side, turns from the palatal side to the opposite side, and scans the contralateral buccal data to the posterior teeth.

④ Turn to the lingual /palatal side to scan data.

⑤ Turn to the buccal side and scan. After reaching the midline, remove the scanning probe head from the mouth and place it on the opposite contralateral posterior occlusal surface. After repositioning, slowly turn to the buccal side to scan the buccal data.

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For scanning methods, please refer to the matching scanning training video.

Note: - The saliva and blood stains should be cleaned before scanning.

- The lips should be removed when scanning the anterior teeth of the upper and lower jaw.
- When scanning the buccal and lingual side, the data of the incisal teeth should be covered to avoid repeated scanning in the same position.

- Bite scanning: when scanning the bite, it is necessary to scan 1 to 2 buccal points to determine the bite relationship. If scanning the full mouth, scan the two points on the left and right buccal sides; if it is half mouth, scan one point on the left or right buccal side.

① Clamp the patient's upper and lower jaw (swallow), keep still.

② Place the probe head on the patient's buccal side, and slowly move forward with the “S” path between the upper and lower jaw. Take 4 teeth and match the upper and lower jaw. Then click the “S” key to pause and switch to the buccal point 2, move slowly between the upper and lower jaw in an “S” shaped path, take 4 teeth and match the upper and lower jaw, click Enter.

- Point switching: after scanning one bite relation, pause the scan and switch the buccal point 2 to scan. Switching the buccal point can be achieved through the shortcut key "P".

- Implant scanning: when the treatment planning is implant, the jaw selection area will appear in two types as is shown in the following figures: jaw scan and body scan.

- Implant working jaw: first remove the healing cap and scan the working jaw where the cuff is located.
- Implant scan body: put on the scan body, and then only need to scan the data of the scan body

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completely, without scanning other data.

After the scan is completed, a small green check mark will appear in the lower right corner of the icon; when not scanning, the icon remains unchanged. When data are missed, a red exclamation mark will appear in the lower right corner of the icon.

Note: - The bite relationship of the patient must remain the same during the buccal scan, and the swallowing action is good for the bite relationship.
- Before processing model, please use editing tool to remove extra data.


3.2.4 Processor

After the scan is completed, click button  or the “Enter” to enter the model processing interface, and the scanned model will be processed.

>>> Processor Page

- Upper/lower jaw interface

Click the maxillofacial selection button below the interface to switch the maxillofacial surface.

Edit button: click to enter the editing function area, including circle selection, restore, flip, return and other functions.


Polygon button: after clicking, the mouse will become polygonal, and the circled closed area will be deleted.

Reback button: click once to restore the previous state of this operation, and click again to restore all.

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Flip button: convert impression data into normal 3D data.

Clean up button: automatically clean up the redundant data, used in the end of scanning state.

Return button: exit the current editing ribbon.

Capture button: capture the three-dimensional data of the current interface and save.

Margin button: click to mark the cervical margin of the selected tooth position.

Hollow button: click to check the concave situation from various directions.

Base button: add the base for the data model; click the button again, cancel the base.

Measure button: select any two points on the teeth and jaw to measure the distance between them.(Only in orthodontic mode)

Color button: Click to switch the current 3D data between true color and plaster color.

Compare button:You can select target data to compare with current model data and calculate r.If the data to be compared is the case being scanned, click the orthodontic (repair/implantation) button in the upper right corner to save the case and then compare the data and calculate the error


/

/


Save button: this save function only exists in the model processing interface;After loading and batch processing

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Note: - Hold down the left mouse button to rotate the current 3D model.
- Hold the right mouse button to move the current 3D model.
- Scroll the mouse wheel to zoom in and out the current 3D model.

- Bite interface

Align **Registration button:** manually register the bite relationship.

Distance **Distance button:** check the occlusal distance of the current model data.

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The cervical margin confirms the position of the shoulder edge by drawing the margin line, which provides a guidance basis for subsequent data processing, as is shown in the figure:

Adjustment button: click to adjust, the cervical margin line is displayed in the form of small dots, drag the small dot to change the cervical margin line.

Clear button: click to clear, remove the cervical margin mark of the selected tooth position.


Step 1

Hold down the left mouse button to draw the neck edge line in the abutment shoulder. position.

Step 2

Click the "Trim" button, the mark line appears a lot of points, the left mouse button to select the need to adjust the position of the point, you can fine-tune the position of the mark line.

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- Click the "clear" button to clear the currently marked neck line.
- When multiple teeth are restored, directly click the tooth number ID to switch the teeth position.

The data is saved and the corresponding neck line data file is generated in the default save path. When the data is designed, the marked neck line can be displayed after the file and the scan data are imported into the dental CAD software at the same time.

Note ：
- Tip: double-click the model data to enlarge the position clicked.
- Only by selecting teeth in restoration mode can you use the neck margin function.
- Click the shortcut key C to switch to plaster color to help the neck margin mark.

- Hollow

Click the “Hollow” to check the undercut at different perspectives. The undercut area will be marked with gradient color.

Step 1 ：

Adjust the data model to the point of view that needs to be observed.

Step 2 ：


Click the “Hollow” button to view the concave condition. Hold down the left mouse button to rotate the model data. Release the mouse and click the left button again to turn the data model to adjust the observation direction. Double-click to zoom in on the tooth. The direction of the arrow represents the current observation direction.

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- **Base**

Preview function to add a base to the data model, as is shown in the figure. Clicking again to cancel the base. When the data center module is saved, the model data saved after adding the base is selected, which is consistent with the base style of the current preview.


- Compare


There are two ways for selecting the models to be compared.

– The first way is selecting the current historical case of the patients to get the data to compare. Double click the treatment history in Medical Case， and click the "OK" button to make a comparison.

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
- The second way is directly selecting the data from the storage path, and click "Confirm" after selecting.

Step 1

Click compare, select the object to be compared (patient-related medical records or upper and lower jaw storage path), and click "Confirm". As is shown below.

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
Step 2

Enter the model alignment interface, using direct automatic alignment or three-point alignment, and click “OK”. As is shown below.

Direct alignment: the alignment interface is determined directly, and the software automatically carries out.

Three points alignment： select 3 distinctive points in corresponding places in the current and the standard data. And click “OK” to start.


Unselect points： for unselecting some selected point, double click while pressing "Ctrl". For unselecting all the selected point, click “Reset” to clear all the points.


Step 3

After selecting, users can check the comparison results, as is shown below. The middle is the result of the alignment between the current model and the selected standard model. The left side is the maximum distance setting area. The different distances within the maximum distance are displayed in different colors. The lower left corner of the page is the calculated coverage and the accuracy.

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
- Align

Recalculate the current occlusion relationship.

As is shown below, after the selection is completed, the middle zone is the bite registration data display area. The left side is the selection of upper / lower jaw, and the currently selected upper/lower jaw model data in the lower left corner.

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Step 1

Click the “Align” button to enter the alignment interface.

Step 2

After selecting "Upper jaw ", use the three-point alignment method or click “Align” directly.

Step 3

After selecting "lower jaw", use the three-point alignment method or click “Align” directly.

Step 4

Click the "check" button, and the display area will display the bite data after alignment.

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Note:


- Direct alignment: select the “Align” button on the right side of the upper and lower jaw for automatic alignment.
- Three-point alignment: select three points with obvious features in the corresponding place between the current data and the standard data, click the “Align” button for alignment. The “Reset” button can clear the selected points.

- Distance

After the alignment, click “Distance” to calculate the occlusal distance and show the result. As is shown below, the left side is the rainbow diagram for distance reference, and the middle zone is the calculated results of the two jaws.


The darker the occlusal area is, the smaller the occlusal distance is. The lighter the occlusal area is, the larger the occlusal distance is. When the occlusal area is dark blue, it indicates that the distance is negative. When there is no rainbow picture distribution in the occlusal face, it indicates that the occlusal distance between the upper and lower jaw in this area is greater than that of 2mm.

When it is found that the occlusion relationship of the scanned data is significantly different from the actual occlusion state, you can use the "Adjust" button to automatically optimize the occlusal jaw.


Double-click to separate the upper and lower jaw to see the distance, as is shown below, double-click again to restore the occlusal state.

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● Measurement


Select any two points on the jaw and measure the distance between the two points

Double-click to cancel all measurement distances. The measurement function can only be used in "Orthodontic" mode.

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3.2.5 Storage

① Information display area: displays the relevant information of the current case, including the patient's name, number, treatment doctor, treatment plan description and tips for the current data preservation format.

② Design / save area: including choosing whether to add a base to the saved data, editing and displaying the remarks bar, and CAD software to open the data, save the data, package the data, upload the data on the network (this function is to be opened) and other ways to save the data.

③ Select save data type: according to their own needs to set the data type to save.

④ Copy: Store the current case as backup data. After operating on the current case, the original data will also exist. It is suitable for taking multiple bite scenes.

Add the base of the model: save after selection, and the stored data has a base, which is consistent with the preview of the base in the processor page.

Design button: open the related dental CAD software directly, then create the project to process the scanned 3D data and design.

Save as: save the current case to the specified path.

Package: compress and package the scanned data and open the folder.


Upload button： The scanned data will be uploaded to the cloud. During the upload process, the operation of data uploading can be stopped by clicking the Cancel button

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
3.2.6 Setting

Software settings include system, network, registration and calibration. Click on the upper left corner  to open the settings interface, as is shown in the following figure:

>>> ommon

Set some default options when scanning, including language, Frame rate prompt, data format, etc.


- Language: Select different languages as you needed，only English and Chinese are supported now
- Tip: select whether to display the prompt box. “Yes”: a prompt box is displayed in the lower left corner during scanning. “No”: there is no prompt box during scanning.

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
- Music type: select the type of music when scanning. Turn off the music.
- Data format supported formats of the output data are STL, PLY and PTY (*).
- Frame rate: You can choose different frame rates according to your needs. There are two options by default.
- Save path name format: including [time-number], [time-name], [scan type-time-number], [scan type-time-name], [doctor//scan type_time_number], [doctor //Scan Type_Time_Name] Six formats.
- Department:
 - Orthodontics: no need of treatment setting.
 - Restoration: needs Treatment setting.
 - Implant: needs Treatment setting.
- Mode:
 - Grade 1 or 2 is for scanning white plaster models.
 - Grade 3 is suitable for resin molds, and grade 2 or 3 is suitable for scanning colored plaster models.
 - Grade 4 is for scanning resin artificial teeth.
 - Grade 4 or 5 are for scanning real teeth.
- Package path: the path to save the current compressed packaged project for the convenient transmission of the data. The default path: “D: /data/compress”.
- Project path: the path to save the project data. The default path: “D: /data”.

If users need to restore the system default settings, click "Default" to pop up and the confirmation box will pop up as is shown in the following figure.

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
- Color


Color function is used to calibrate the color of the data. As is shown in the following picture click “White Balance” to turn it on. Set the correct grade, put the scanner close to the white paper, and adjust the brightness to gear 1. Then click Start. After the progress, the information feedback area will display the information of the success or failure. **In general, do not change it.**

- 3D

When the scanner has undergone strong impact, severe temperature change or long transport, it should be calibrated immediately. Calibration should be taken regularly every 5-7 days for maintaining the accuracy of the scanner.

After turning on the scanner, connect the calibrator, click Calibration Mode icon and start calibration mode. And the user can follow the prompts in the information feedback area. As is shown below.

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- **Calibration operation flow**


- **Automatic calibrator calibration**

(1) Remove the probe head assembly of this product, and align the silver protrusion of the probe body with the groove of the calibrator and insert it to make it tightly connected.

(2) Click the position button to start calibration.

>>> Account


This module displays the currently logged on network disk space account information. Includes account name, account space address, system port, and system backup and system space automatic optimization functions. As is shown below.

>>> Help

As is shown below. If you need to feedback the experience of using the software or upload Bug, you can send E-mail to support@freqty.com, or leave us a message on Wechat official account. We will contact you as soon as possible.

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
>>> About

This software must be licensed before it can be used. As is shown below. If the software does not have an authorization number and expiration date, it can be considered as expired or pirated software, please stop using it and contact us.

This software will be updated and maintained from time to time. Please pay attention to the software version number for details. You can get the latest version of the software by following Update the version number.

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Note. The software version in the figure is an example of the software version, which is not the latest and current version

3.3 Shortcut key

Open an external case	Ctrl + O	Save the current project	Ctrl + S
Minimize	Ctrl + M	**Close**	Ctrl + Q
Switch between the bite points	P	**Save images in real time**	T
Color check	C	**Start/Pause**	S
Finish scan	Enter	**Switch between the jaws**	E/R

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4 Equipment sterilization and routine maintenance

4.1 Probe head sterilization

As the contamination of the probe head mirror will seriously affect the scanning effect, in order to ensure the normal use of the scanner, please strictly follow the high temperature and high-pressure disinfection steps below.

Step 1 ：

Rinse the probe head with water and wipe it dry with absorbent gauze in time. Do not wash the equipment body with water. Use absolute alcohol after drying, clean the probe head mirror carefully after soaking the absorbent gauze, keep the lens surface free of dirt.


Step 2 ：

Then fold the clean absorbent gauze to the long strips which is similar to lens size, please wear gloves or finger cots to prevent contamination of gauze.

Step 3 ：

Then cover the surface of the probe head mirror with absorbent gauze. Please take good care of contact and coverage.

Put the head of the probe covered with gauze into the sterilization bag. Do not shift or drop the gauze during placing processes.

Step 4 ：

For final disinfection, place the probe head window face down. Suggest to sterilize with high-temperature steam for a few minutes, such as134 °C for at least 6 minutes or 121 °C for at least 15 minutes.

The bag with probe head must be maintained vacuum during sterilization.

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- The absorbent gauze used in each step of the operation cannot be reused.
- Do not use absorbent cotton instead of absorbent gauze. The former is prone to residual cotton wool and difficult to clean.
- Use absolute alcohol when wiping the mirror. Do not use medical alcohol.
- The head of the probe, especially the lens, must be wiped dry after cleaning in the first step. Be careful not to scratch the mirror surface.
- When the disinfection is repeated for many times, which causes the mirror surface to be stained and affect the scanning effect, please replace the probe with a new one in time.

4.2 Routine maintenance

- **Lens use:**

It is prohibited to blow the inside of the probe with a three-purpose air gun to prevent dust from being blown into the photoelectric system inside the device.

- **Equipment placement:**

The base is placed on a stable table, and the handle is placed in the groove of the base after use. Pay attention to protection to avoid unnecessary vibration or bumps.

- **Accuracy calibration:**

Under normal circumstances, the calibration is performed once a week. If the scan is easily interrupted, the calibration should be performed immediately.

5 Software and hardware common problems and solutions

5.1 Software startup issues

5.1.1 The software cannot be opened normally

- Check if the USB-KEY is connected normally and within the validity period.
- If prompted “Workflow initial error” or “Unable to open file directory”, please reinstall the software.

5.1.2 Software starts too slowly

- Check if the computer power is connected normally.
- Confirm that the scanning software is running in “Administrator Mode” of Windows operating system.

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- Check if the Windows operating system is being updated. If it is being updated, please complete the update before using the software。
- Check if other software of the computer can start normally.

5.1.3 Software icon on the desktop turns white

- Check if there's any anti-virus software, if so, either put the software into the white list or uninstall the anti-virus software.
- Right-click the icon and select "Open the file location" to see if the software icon is in the folder turns white. If the software can be opened normally and the icon in the folder is normal, delete the icon on the desktop and send the icon in the folder to the desktop. If the icon in the folder turns white, contact the customer service for help
- The icon turns white. If you cannot run the program by double-clicking, please reinstall the scanning software.

5.2 Problems connecting devices

5.2.1 Camera connection failed

- Check whether the device power is on normally, whether the computer power is tightly connected, and whether the USB interface is correctly connected to the computer's USB3.0 interface.
- Replace interface or connect Hub, etc.
- If possible, replace the computer or device and try again.

5.3 Image display problems

5.3.1 No image display in 2D image area

- Make sure the device's USB interface is properly connected to the computer's USB 3.0 interface.
- Restart the software and scanning device to check if the image can be displayed normally.

5.3.2 2D image flicker

- Check if the modulator is connected properly.
- Replace the USB port of the device with the computer.
- Connect your computer to the Internet.

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5.4 Scanning issues

5.4.1 Scans are easily interrupted and not smooth

- Inappropriate scan brightness. For plaster model scanning, choose 1/2, for resin model scanning, choose 3, for the intraoral scanning, choose 4, 5 is suitable for patients with darker teeth in the mouth.
- During scanning, confirm that A above the image area is blue. If it is black, use the keyboard A key to switch.
- Standardize scanning methods. Ensure coverage of scanned data with existing data.

5.4.2 Out-sync of data between 2D and 3D

- Confirm whether the computer configuration meets the requirements (higher than or equal to our recommended configuration).
- Delays caused by too many scans (single jaw scans should be completed within 3 minutes).
- Uninstall antivirus software or add scanning software to the whitelist of antivirus software.
- Check the status of windows update. If the update is in progress or has failed, please restart the computer after the update is completed before using the scanning software.

5.4.3 Difficulty for scan relocation

- Ensure that the scanning direction is consistent with the previous scanning when repositioning
- Avoid long scans.

5.4.4 No 3D data when scanning

- Recalibration.

5.5 Abnormal interrupt during scanning

- Check the status of windows update. If the update is in progress or has failed, please restart the computer after the update is completed before using the scanning software.
- Check whether the remaining storage space of drive C is sufficient.
- Turn off or uninstall anti-virus software.

5.6 Problems with calibration

5.6.1 Failed to add position or calibration failed.

- If there are obvious stains or debris on the calibration plate in the 2D image display area of the software, please invert the calibrator first to make the debris out.

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- Do not blow the inside of the calibrator directly with your mouth.

5.7 Other problems

5.7.1 Computer restarts repeatedly

- Reinstall NVIDIA graphics driver.
- If you still have problems, please replace your computer.

5.7.2 The device cannot be powered on normally

- Check if the adapter indicator is on and the adapter is powered on normally.
- Check if the current socket is powered.
- Replace the adapter and power cord.

Kind reminder: If the above operations do not solve your problem, please contact technical support in time.

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6 Maintenance and technical support

Do not open the device without permission, if the device has any problems that you cannot solve. If you have any questions or requests, please contact us in the following ways.

198 183 628（Overseas）

813 576 869（China）

(+86)181 6123 1898

(+86)028 8531 2603

market@freqty.com

Ziyang Freqty Medical Equipment Co., Ltd.

Company address: Floor 2-3, unit 7, building 3, No. 222, West Section 3, outer ring road, Yanjiang District, Ziyang City, Sichuan Province, China

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