BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: N50 621 rev. 003	Date: 07-SEP-2021 Proofread No.: 1
Designer: CDS	Page: 1 of 36
Colour Reference: BLACK	Material Specification: HEC1129061

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en

Symbols are referenced from ISO 15223-1 *Medical Devices: Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General Requirements*

DEFINITION OF SYMBOLS

Symbol and Symbol Number	Symbol Title	Description
5.1.1	Manufacturer	Indicates the medical device manufacturer, i.e. the natural or legal person with responsibility for the design, manufacture, packaging and labelling of a device before it is placed on the market under his own name.
5.1.3	Date of Manufacture	Indicates the date when the medical device was manufactured. Format should be YYYY-MM-DD.
5.1.4	Use-by-date	Indicates the date after which the medical device is not to be used. Format should be YYYY-MM-DD.
LOT 5.1.5	Batch code	Indicates the manufacturer's batch code. Synonyms are "lot number" and "batch number".
REF 5.1.6	Catalogue number	Indicates the manufacturer's catalog number. Synonyms are "reference number" and "reorder number".
5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
5.3.6	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed.
5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use. Synonym for "Consult instructions for use" is "Consult operating instructions".
MD	Medical Device	Indicates that the device is a medical device.

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en

GENERAL INFORMATION - REST OF THE WORLD

INDICATIONS FOR USE

The **U9000 Plus Ultrafilter** is indicated for purification of incoming water for dialysis fluid as well as purification of dialysis fluid to obtain *microbiologically high quality dialysis fluid**. The risk of exposure to bacteria and endotoxins can thereby be minimized.

The **U9000 Plus Ultrafilter** is intended to be used in conjunction with a water treatment system.

WARNING! The **U9000 Plus Ultrafilter** can only be used with Gambro/Baxter/Hospal dialysis machines equipped with holder for an ultrafilter.

CONTRAINDICATIONS

There are no known Contraindications for the use of **U9000 Plus Ultrafilter.**

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

WARNING! Is used to alert the operator not to take a specific action that can cause potential hazard which, if ignored, can result in an adverse reaction, injury or death.

CAUTION! Is used to alert the operator to take a specific action to protect against possible hazard which, if ignored, could have an adverse effect on the patient or equipment.

GENERAL WARNINGS AND PRECAUTIONS

WARNING! If a leakage is observed, exchange the filter.

WARNING! Water used for dialysis has to obtain microbiological quality according to international standards (e.g. ISO 23500).

WARNING! Do not use if package and/or protection caps are damaged or not in place.

CAUTION! Federal law (USA) restricts this device to sale by or on the order of a physician.

ADVERSE REACTIONS

If the dialysis fluid or the incoming water is heavily contaminated or the handling of **U9000 Plus Ultrafilter** is incorrect, there is a risk of bacterial and/or pyrogenic overload in the Ultrafilter with consequent risk of breakthrough. The patient may then be exposed to bacteria or pyrogens sufficient to induce a pyrogenic reaction with symptoms of shivering and raised body temperature. Pyrogenic reaction may be accompanied by severe reduction in blood pressure.
Patient exposure to high concentrations of disinfectants due to inadequate rinsing of **U9000 Plus Ultrafilter** after chemical disinfection may lead to symptoms of rash, chest pain and hemolysis. It is therefore imperative to validate the rinsing procedure after each disinfection procedure (residual test).

WARRANTY AND LIMITATION OF LIABILITY

The manufacturer warrants that the **U9000 Plus** has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **U9000 Plus** becoming apparent before the expiry date.

Expiry date: Refer to information on the unit container label.

The warranty above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or fitness or other warranties, which extend beyond those described.

The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **U9000 Plus** and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **U9000 Plus**, whether as a result of any defect therein or otherwise.


The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings, directions and instructions in the labelling, damage arising from events after the manufacturer's release of the **U9000 Plus** before use in order to insure that the **U9000 Plus** is in proper condition, or any warranty given by independent distributors or dealers.

PERFORMANCE

** microbiologically high quality dialysis fluid*
U9000 Plus Ultrafilter is installed to the dialysis fluid circuit of the dialysis machine as described in the respective operator's manual; Concentrates and dialysis water as specified in ISO 23500 have to be used for the preparation of dialysis fluid.

Depending on the fluid circuit setup and the sequence of Ultrafilters, **standard dialysis fluid, ultrapure dialysis fluid and online-prepared substitution fluid used for haemodiafiltration or haemofiltration as described in ISO 23500 can be achieved.** For information on which dialysis fluid quality can be achieved in a given setup, refer to the operator's manual of the respective dialysis machine.

The following reduction factors for bacteria and endotoxins are typically achieved by filtration with **U9000 Plus Ultrafilter**:
- Concentration of bacteria is reduced by a factor of > 10^7 (Brevundimonas diminuta ATCC19146)
- Concentration of endotoxins is reduced by a factor of > $10^{3.5}$ (E. coli 055:B5)

Details of the test methodologies are available to the user upon request.


SPECIFICATIONS

REF	U9000 Plus
Maximum Q_{UF} (l/min)	1.2
Surface area (m^2)	2.4
Storage conditions:	Below + 30°C (+ 86°F).
Materials Membrane Potting Housing Gasket	 PAES/PVP Polyurethane Polyphenylsulfone Silicone rubber

DIRECTIONS FOR USE - REST OF THE WORLD

INSTALLATION / CHANGE OF ULTRAFILTER

U9000 Plus Ultrafilter is connected to the dialysis fluid circuit as illustrated in the figure below:

For detailed information refer to operator's manual of the dialysis machine.

OPERATION

After installation, the **U9000 Plus Ultrafilter** is operated completely by the dialysis machine.

WARNING! The **U9000 Plus Ultrafilter** is delivered non sterile. The **U9000 Plus Ultrafilter** has to be disinfected prior to first use.
For further information see Operator's manual of the dialysis machine.

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FILTER LIFE TIME / FREQUENCY OF CHANGE

The **U9000 Plus Ultrafilter** has to be exchanged regularly.

The frequency for changing the Ultrafilter can vary between once a month and once every three months according to the information given in the operator's manual of the dialysis machine.

After exchange, handle and dispose the product according facility procedures and applicable local, state and federal laws and regulations.

WARNING! U9000 Plus Ultrafilter has to be exchanged according table below.

Maximum time of installation	Maximum number of disinfection (in total)	Maximum number of sodium hypochlorite based disinfections
90 days	150	12

DISINFECTION / CLEANING

WARNING! Disinfection according to this IFU is essential to ensure safety and performance of the product.

U9000 Plus Ultrafilter is included in all dialysis machine cleaning and disinfection procedures. The filter remains mounted in the dialysis machine for storage between treatments.

The following types of disinfection can be performed with **U9000 Plus Ultrafilter**:

- Heat disinfection
- Chemical disinfection with following agents:
 - peracetic acid (≤ 0.1 %)
 - sodium carbonate (≤ 0.5 %)
 - sodium hypochlorite (≤ 0.5 %)
 - citric acid (≤ 2.0 %)

For detailed information refer to operator's manual of the dialysis machine.

en

GENERAL INFORMATION - USA

INDICATIONS FOR USE

The **U9000 Plus Ultrafilter** is indicated for purification of incoming water for dialysis fluid as well as purification of dialysis fluid to obtain standard dialysis fluid in accordance with international standards requirements and local regulations. The risk of exposure to bacteria and endotoxins can thereby be minimized.

The **U9000 Plus Ultrafilter** is intended to be used in conjunction with a water treatment system.

WARNING! The **U9000 Plus Ultrafilter** can only be used with AK98 dialysis machines equipped with holder for an ultrafilter.

CONTRAINDICATIONS

There are no known Contraindications for the use of **U9000 Plus Ultrafilter.**

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

WARNING! Is used to alert the operator not to take a specific action that can cause potential hazard which, if ignored, can result in an adverse reaction, injury or death.

CAUTION! Is used to alert the operator to take a specific action to protect against possible hazard which, if ignored, could have an adverse effect on the patient or equipment.

GENERAL WARNINGS AND PRECAUTIONS

WARNING! If a leakage is observed, exchange the filter.

WARNING! Water used for dialysis has to obtain microbiological quality according to international standards (e.g. ISO 23500).

WARNING! Do not use if package and/or protection caps are damaged or not in place.

CAUTION! Federal law (USA) restricts this device to sale by or on the order of a physician.

ADVERSE REACTIONS

If the dialysis fluid or the incoming water is heavily contaminated or the handling of **U9000 Plus Ultrafilter** is incorrect, there is a risk of bacterial and/or pyrogenic overload in the Ultrafilter with consequent risk of breakthrough. The patient may then be exposed to bacteria or pyrogens sufficient to induce a pyrogenic reaction with symptoms of shivering and raised body temperature. Pyrogenic reaction may be accompanied by severe reduction in blood pressure.

Patient exposure to high concentrations of disinfectants due to inadequate rinsing of **U9000 Plus Ultrafilter** after chemical disinfection may lead to symptoms of rash, chest pain and hemolysis. It is therefore imperative to validate the rinsing procedure after each disinfection procedure (residual test).

WARRANTY AND LIMITATION OF LIABILITY

The manufacturer warrants that the **U9000 Plus** has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **U9000 Plus** becoming apparent before the expiry date.

Expiry date: Refer to information on the unit container label.

The warranty above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or fitness or other warranties, which extend beyond those described.

The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **U9000 Plus** and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **U9000 Plus**, whether as a result of any defect therein or otherwise.


The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings, directions and instructions in the labelling, damage arising from events after the manufacturer's release of the **U9000 Plus** before use in order to insure that the **U9000 Plus** is in proper condition, or any warranty given by independent distributors or dealers.

PERFORMANCE

U9000 Plus Ultrafilter is installed to the dialysis fluid circuit of the dialysis machine as described in the respective operator's manual; Concentrates and dialysis water as specified in ISO 23500 have to be used for the preparation of dialysis fluid.

The following reduction factors for bacteria and endotoxins are typically achieved by filtration with **U9000 Plus Ultrafilter**:

- Concentration of bacteria is reduced by a factor of > 10^7 (Brevundimonas diminuta ATCC19146)
- Concentration of endotoxins is reduced by a factor of > $10^{3.5}$ (E. coli 055:B5)

Details of the test methodologies are available to the user upon request.


SPECIFICATIONS

REF	U9000 Plus
Maximum Q_{UF} (l/min)	1.2
Surface area (m^2)	2.4
Storage conditions:	Below + 30°C (+ 86°F).
Materials Membrane Potting Housing Gasket	 PAES/PVP Polyurethane Polyphenylsulfone Silicone rubber

DIRECTIONS FOR USE - USA

INSTALLATION / CHANGE OF ULTRAFILTER

U9000 Plus Ultrafilter is connected to the dialysis fluid circuit as illustrated in the figure below:

For detailed information refer to operator's manual of the dialysis machine.

OPERATION

After installation, the **U9000 Plus Ultrafilter** is operated completely by the dialysis machine.

WARNING! The **U9000 Plus Ultrafilter** is delivered non sterile. The **U9000 Plus Ultrafilter** has to be disinfected prior to first use.

For further information see Operator's manual of the dialysis machine.

FILTER LIFE TIME / FREQUENCY OF CHANGE

The **U9000 Plus Ultrafilter** has to be exchanged regularly.

The frequency for changing the Ultrafilter can vary between once a month and once every three months according to the information given in the operator's manual of the dialysis machine.

After exchange, handle and dispose the product according facility procedures and applicable local, state and federal laws and regulations.

WARNING! U9000 Plus Ultrafilter has to be exchanged according table below.

Maximum time of installation	Maximum number of disinfection (in total)	Maximum number of sodium hypochlorite based disinfections
90 days	150	12

DISINFECTION / CLEANING

WARNING! Disinfection according to this IFU is essential to ensure safety and performance of the product.

U9000 Plus Ultrafilter is included in all dialysis machine cleaning and disinfection procedures. The filter remains mounted in the dialysis machine for storage between treatments.

The following types of disinfection can be performed with **U9000 Plus Ultrafilter**:

- Heat disinfection
- Chemical disinfection with following agents:
 - peracetic acid (≤ 0.1 %)
 - sodium hypochlorite (≤ 0.5 %)
 - citric acid (≤ 2.0 %)

For detailed information refer to operator's manual of the dialysis machine.

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Gambro, Hospal, Baxter and U9000 Plus Ultrafilter are trademarks of Baxter International Inc. or its subsidiaries.