Cranberry®

# Instruction for Use

1. **Product Name:** Dental Dam (Latex and Non Latex)

2. **Product Specifications :**
This product contains Latex Dental Dam or Non Latex Dental Dam, refer to table below for Dental Dam Product Specifications.
Table 1 - Product Measurement and Colour.

| Dental Dam | Measurement | Colour |
|---|---|---|
| Latex | 127mm×127mm (5″×5″) | Green |
| | 127mm×127mm (5″×5″) | Blue |
| | 152mm×152mm (6″×6″) | Green |
| | 152mm×152mm (6″×6″) | Blue |
| | 152mm×152mm (6″×6″) | Black |
| Non Latex | 152mm×152mm (6″×6″) | Purple |

3. **Composition:**
Latex Dental Dam is made from Natural Rubber Latex. Non Latex Dental Dam is made from Synthetic Polyisoprene.

4. **Intended Use:**
Dental Dams are thin sheets of latex or non latex used to isolate specific teeth during dental procedures.
The domain of usage of dental dam is to ensure the best use of dental treatment, their benefits includes:
- Infection control barrier
- Improve the operation way/field
- Increase operational efficiency
- Isolate a tooth to obtain a dry working field
- Patient management

5. **Storage Conditions:**
Latex dental dam and non latex dental dam is an easy to aging, and the condition of preservation of the product is very sensitive. Many variable factors (especially temperature) can affect this material. Therefore the storage temperature should be lower than 104°F or 40°C.

6. **Warranty Period:**
In order to achieve the best usage effect, dental dam should be used before the expiry date as stated on colour box. Expiration of latex dental dams is 5 years of the manufacturing date while for non latex dental dams is 3 years of the manufacturing date.

7. **Duration of use:**
It is intended for continuous use for less than 60 minutes during dental procedures.

8. **Basic principles which medical device works on and their features:**
(1) Place the Dental Dam in a flat surface.
(2) Determine the center of dam and punch a hole approximately 2mm diameter.
Choose the right size of puncher according to the contacted teeth.
(3) Use a forceps to enlarge the hole and place a metal clamp to hold the enlarged hole.
(4) Place a dam and the clamp into the affected tooth.
(5) Place the dental dam frame outside the mouth to hold excess dam materials.

9. **Caution:** This product contains natural rubber latex which may cause allergic reactions.

10. **Precautions:** (1) Do not use damaged dental dam. (2) For single use only.

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