Cranberry®

Instruction for Use

Product Name
Powder-Free Nitrile Examination Gloves

Description
Powder-Free Nitrile Examination Gloves, Non-Sterile, Single Use Only

Intended Use
A powder-free nitrile examination glove is a disposable glove product that is only applicable for examination, medical treatments, and handling contaminated medical devices to prevent cross contamination between users and patients. It is strictly not suitable for surgical use.

Caution
This product contains nitrile rubber which may cause allergic reactions in some individuals. Do not expose these gloves to any person with a known or suspected sensitivity to nitrile rubber. If allergic reactions occur, discontinue use immediately and consult a physician. This product is not made with natural rubber latex.

Sterilization
Non-sterile

Harsh Chemicals
These gloves should not be used for mechanical, thermal or harsh chemical protection. Although good quality nitrile gloves provide an excellent barrier, they are not intended for application involving prolonged, direct exposure to harsh chemicals, where heavy duty industrial gloves are required. Variability in material thickness, glove integrity, chemical concentration, temperature, and length of exposure to chemicals may affect performance. User is responsible for determining glove suitability for individual application. These gloves are for single and transient use only.

Storage
Store in a dry place at room temperature. Exposed product should be shielded from direct sunlight, intense artificial light, x-ray machines, and other sources of ozone.

Disposal
Gloves are biodegradable products where no special decommissioning or disposal is required. If the glove is contaminated with a toxic compound or biological material that is covered by any disposal regulations, the gloves must be handled in the same way as the toxic material itself. If gloves are not contaminated or have been properly decontaminated, either landfill or incineration is a satisfactory means of disposal.

Shelf Life
Three years from manufacturing date

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