
ChM®

EN	REUSABLE ORTHOPAEDIC AND SURGICAL INSTRUMENTS	3
PL	NARZĘDZIA CHIRURGICZNE I ORTOPEDYCZNE WIELOKROTNEGO UŻYTKU	9
RU	ХИРУРГИЧЕСКИЕ И ОРТОПЕДИЧЕСКИЕ ИНСТРУМЕНТЫ МНОГОРАЗОВОГО ИСПОЛЬЗОВАНИЯ	15
ES	INSTRUMENTOS QUIRÚRGICOS Y ORTOPEDICOS REUTILIZABLES	23
DE	WIEDERVERWENDBARE CHIRURGISCHE UND ORTHOPÄDISCHE INSTRUMENTE	29
CS	OPAKOVANĚ POUŽÍVANÉ ORTOPEDICKÉ A CHIRURGICKÉ NÁSTROJE	35
IT	STRUMENTI ORTOPEDICI E CHIRURGICI RIUTILIZZABILI	41
UK	ХІРУРГІЧНІ ІНСТРУМЕНТИ БАГАТОРАЗОВІ	47
LV	ATKĀRTOTI LIETOJAMIE ORTOPĒDISKIE UN ĶIRURĢISKIE INSTRUMENTI	53
FR	OUTILS CHIRURGICAUX ET ORTHOPÉDIQUES RÉUTILISABLE	59
BG	ОРТОПЕДИЧНИ И ХИРУРГИЧЕСКИ ИНСТРУМЕНТИ ЗА МНОГОКРАТНА УПОТРЕБА	65

IFU-I-001/2023-08-10

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EN

INSTRUCTIONS FOR USE

REUSABLE ORTHOPAEDIC AND SURGICAL INSTRUMENTS

1 INDICATIONS *4*

2 DESCRIPTION *4*

3 MATERIALS *4*

4 WARNINGS AND PRECAUTIONS *4*

5 CLEANING, DISINFECTION, STERILIZATION *5*

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1 INDICATIONS

1. Surgical and orthopaedic instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.

2 DESCRIPTION

1. The unit package contains one piece of the product in non-sterile condition. Clear plastic bags are a typical packaging material. The products may also be supplied as a complete set (*arranged on palettes and placed into specially designed sterilization containers*).
2. The package is equipped with the product label. The label (*as a primary label*) contains, among others:
 1) Logo **ChM** and the address of the manufacturer.
 2) Catalogue number (*REF*), e.g.: 40.XXXX.XXX, and device name and size.
 3) Production batch number (*LOT*), e.g. XXXXXXX.
 4) NON-STERILE sign - indicates non-sterile product.
 5) Information symbols (*described in the footer of this Instructions For Use*).
 6) CE conformity mark.
3. Depending on the size or type of the product, the following information may be marked on its surface: manufacturer's logo, production batch no. (*LOT*), catalogue no. (*REF*), type of material and device size.

3 MATERIALS

1. For the production of instruments, **ChM** uses mainly: steel, aluminum alloys and plastics, approved for use in surgical instruments and in accordance with applicable procedures.
2. Instruments are produced of corrosion-resistant steel. The protective layer (*passive layer*) against corrosion is formed on the surface of the device due to high content of chromium.
3. Devices produced of aluminium are mainly stands, palettes, cuvettes and some parts of instruments such as e.g. handles. The protective oxide layer which may be dyed or stays in natural colour (*silvery-grey*) is formed on the aluminium as an effect of electrochemical treatment of its surface.
4. Devices made of aluminium with processed layer have good corrosion resistance. However, the contact with strong alkaline cleaning and disinfecting agents, solutions containing iodine or some metal salts, due to chemical interference with the processed aluminium surface, shall be avoided.
5. Devices produced of plastics are mainly stands, palettes, cuvettes and some parts of instruments such as e.g. handles. Plastics used in the manufacture of instruments are mainly: PPSU (*Polyphenylsulfone*), PEEK (*Polyetheretherketone*), teflon (*PTFE - Polytetrafluoroethylene*) and silicone. The above-mentioned materials can be processed (*washed, cleaned, sterilized*) at temperatures not higher than 140°C. They are stable in aqueous solution of washing-disinfecting agents with a pH value from 4 to 10.8.
6. Steel surgical instruments with a hardened insert are more durable than steel products. The advantage of the product is the sintered carbide insert placed in the working part of the instrument. This insert is characterized by great hardness and abrasion resistance.
7. If the material of the device cannot be specified, please contact **ChM** representative.

4 WARNINGS AND PRECAUTIONS

1. Instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.
2. Improper, careless and inconsistent with the recommendations provided below handling of the instruments can lead to their chemical, electrochemical or mechanical damage which can adversely affect corrosion resistance and shorten the service life of the devices.
3. Instruments are intended only for specific procedures and must be used strictly according to their intended purpose. Use of instruments not in accordance with their intended purpose may lead to malfunction, accelerated wear and, in consequences, damage to the instrument.
4. The surgeon should be familiar with all components of the device before use and should personally verify if all components and instruments are present before the surgery begins.
5. Before the procedure begins, all instruments should be carefully inspected for their condition and proper functioning. They should be undamaged and without any signs of corrosion. Blades and cutting edges should be sharp and undamaged. Damaged or corroded instruments should be immediately replaced. The use of bent, damaged or corroded instruments is not allowed.
6. Tissue structures close to the operative site must be protected.

7. Collision of the instrument with metal operating equipment, retractor or other device may cause damage that necessitates intraoperative replacement of that instrument.
8. Do not apply excessive force when using the instrument – it may lead to its permanent damage and, in consequences, to malfunction of the device.
9. Instruments are subject to constant wear processes. While rare, intraoperative fracture or breakage of the instrument can occur. Instruments which have been subjected to prolonged use or excessive forces are more susceptible to fractures, depending on care taken during surgery and the number of procedures performed. Should breakage occur, the instrument parts must be removed and disposed of immediately in accordance with valid medical facility procedures.
10. In order to confirm the removal of all undesired metal fragments from the surgical field, intraoperative X-Ray examination is recommended.
11. In the case of suspected or documented allergy or intolerance to metallic materials, surgeon shall find out if the patient develops allergic reaction to the instrument material by ordering appropriate tests.
12. It is extremely important to follow the calibration deadline which is permanently marked on the torque instruments (*see CALIBRATION*). Use of a torque instrument with an overstepped calibration date may lead to potential injury, implant or device damage, or loss of correction. If there appear any irregularities in device operation, e.g. due to heavy usage, prior to next calibration date, the instrument should be immediately sent to the manufacturer for its re-calibration.
13. Instrument which had contact with tissues or body fluids of another patient cannot be re-used prior to its reprocessing due to a potential risk of cross-infection caused by viruses, bacteria and prions.
14. Middle and working part of the surgical devices with hardened insert shall be used during the surgical procedure. Improper or inconsistent with the intended purpose usage of the product may lead to damage of the working part e.g. damage to the inserts.

5 CLEANING, DISINFECTION, STERILIZATION

1. Prior to use of a non-sterile device, the following rules apply:
 1) The device must undergo cleaning, disinfection and sterilization procedures.
 2) Effective cleaning is a complicated procedure depending on the following factors: the quality of water, the type and the quantity of used detergent, the technique of cleaning (*manual, automated*), the proper rinsing and drying, the proper preparation of the device, the time, the temperature and carefulness of the person conducting this process, etc.
 3) The hospital facility remains responsible for the effectiveness of the conducted cleaning, packaging and sterilization processes with the use of existing equipment, materials and properly trained personnel.
2. Preparation at the place of use.
 1) Immediately after use, remove from instrument blood and other contaminants with disposable cloth or paper towels. Additionally, it is recommended to rinse the instrument under running water or to place it in the aqueous disinfectant solution. Do not let blood, tissues, body fluids or other biological impurities dry out on the surface of the device.
 2) In order to prevent blood and debris from drying out on the instrument surface, transport the product to the processing area in a closed container or covered with a damp cloth.
 3) In order to avoid contamination during transportation, the dirty instruments should be separated from the clean ones.
3. Preparation for washing and disinfection (*for all methods*).
 1) The used instruments should be reprocessed as soon as possible.
 2) If the instrument can be disassembled, it must be done before cleaning processes.
 3) Rinse under running water and remove surface debris using a disposable cloth, paper towel or plastic brushes (*nylon brushes are recommended*). Particular attention should be paid to openings and places difficult to be cleaned. Very dirty devices should be soaked in an aqueous solution of a detergent or a washing-disinfecting agent, e.g. neodisher® MediClean forte, at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 4) CAUTION: It is forbidden to use brushes made of metal, bristles or materials which could damage the product.
4. Cleaning and disinfection process.
 1) This Instructions for Use describes two **ChM**-approved cleaning and disinfection methods: manual with ultrasound cleaning and automated method. It is recommended to use automated cleaning and disinfection procedures (*in a washer-disinfector*).
 2) The chosen washing and disinfecting agents must be suitable and approved for use with medical devices. It is important to follow the instructions and restrictions specified by the producer of those cleaning agents. It is recommended to use aqueous solutions of washing-disinfecting agents with a pH value between 10.4 and 10.8. **ChM** used the following materials during the validation process of the described recommendations for cleaning and disinfection. It is allowed to use other materials than those listed below which may also give a comparable effect:

a) detergent - Dr.Weigert (*producer*) neodisher® MediClean forte (*name of the detergent*);
b) disinfectant - Dr.Weigert (*producer*) neodisher® Septo Active (*name of disinfectant*).

3) To prevent product damage (*pitting, rust, discoloration*), do not use aggressive cleaning agents (*NaOH, NaOCl*), saline solutions and unsuitable cleaning agents.
4) Where possible, it is recommended to use demineralized water to avoid the formation of spots and stains caused by chlorides and other compounds present in ordinary water.
5) Manual with ultrasound cleaning.
 a) Equipment and materials: a device for ultrasound cleaning, soft, lint-free cloths, plastic brushes, syringes, aqueous solutions of cleaning agent.
 b) Manual cleaning: Initial manual cleaning must be performed prior to ultrasound cleaning.
 c) Rinse under running water until the product is visually clean. Use plastic brushes to remove heavy or large debris.
 d) Soak the product for at least 10 minutes in an aqueous solution of a detergent at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 e) Rinse the product under cold water for at least 2 minutes, paying particular attention to the holes and places difficult to be cleaned.
 f) Prepare fresh washing solution. Clean the surfaces and gaps of the product, carefully. Use suitable brushes to clean the holes. Clean the product immersed in the solution.
 g) Rinse the product thoroughly under warm running water for at least 2 minutes, paying special attention to the gaps, blind holes, hinges and joints. When cleaning, use brushes and perform multiple reciprocating movements on the surface of the product.
 h) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-h until the product is visually clean.
 i) Ultrasound cleaning: prepare an aqueous cleaning solution at a temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the cleaning agent, in respect of temperature, concentration, exposure time and water quality*). Immerse fully the product in the aqueous cleaning solution and have it washed in ultrasounds for 15 minutes.
 j) Rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 k) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-k until the product is visually clean.
 l) Use demineralized water for final rinsing of the device.
 m) Dry the device thoroughly using disposable, soft, lint-free cloth or compressed air.
 n) Prepare an aqueous solution of disinfecting agent at a temperature of 20+/-2°C using 20g of the agent per 1 liter of water. Immerse the product in the solution, exposure time – 15 min (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 o) After the exposure time, rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 p) The cannulated instruments should be treated using a compressed air or air supplied from the syringe.
 q) Dry the device thoroughly. It is recommended to dry the product in a dryer at a temperature ranging from 90°C to 110°C, for 60 min.
 r) Visually inspect the entire surface of the device.
 s) CAUTION: If the obstruction in the cannula cannot be removed as indicated in the Instructions for Use, the device should be considered at the end of its useful life and should be discarded in accordance with facility procedures and guidelines.
6) The automated method using a washer - disinfector.
 a) Equipment and materials: a washer - disinfector, aqueous solutions of cleaning agent.
 b) Cleaning in the washer-disinfector must be preceded by a manual and ultrasound cleaning, following the procedure described in subsections c-h of paragraph 5.
 c) CAUTION: The equipment used for washing/disinfection should meet the requirements of ISO 15883. Procedure of washing in the washer-disinfector shall be performed according to internal hospital procedures, recommendations of the washer-disinfector manufacturer, and instructions for use prepared by the washing-disinfecting agent manufacturer.
 d) The device should undergo the process of machine washing in the washer-disinfector using the following cycle parameters: (*1*) - pre-washing in cold tap water, duration – 2 min; (*2*) - washing in an aqueous solution of cleaning agent at 55+/-2°C and pH of 10.4 - 10.8, duration – 10 min; (*3*) - rinsing under demineralized water, duration – 2 min; (*4*) - thermal

disinfection in demineralised water at 90°C, minimal duration – 5 min; (*5*) - drying at the temperature ranging from 90°C to 110°C, duration - 40 min.

5. Inspection
 1) Each time before re-use and re-sterilization, all medical devices should be inspected.
 2) All parts of the product should be checked for visible dirt and corrosion. Particular attention should be paid to:
 a) Places where dirt can be found, such as joints, latches, etc.
 b) Holes, grooves and gaps the debris could have been pressed into during use.
 3) Generally unmagnifed visual inspection under good light conditions is sufficient.
 4) Each time before re-use and re-sterilization, the functional check of the product should be performed, consisting of:
 a) Verifying the connections in the mating instruments, such as tips, shafts and quick coupling devices.
 b) Verifying the correct functioning of mechanisms, e.g. screw, ratchet, snap mechanism, etc.
 c) Verifying all rotating devices for straightness (*this can be simply achieved by rolling the device on a flat surface*).
 d) Verifying cutting edges for sharpness.
 e) Verifying instruments for damage to material structure (*cracks, dents, peels, etc.*).
 5) Damaged or defective product cannot be approved for further use.
 6) Prior to storage, the instrument must be checked for dryness.
 7) CAUTION:
 a) The **ChM** does not define the maximum number of uses appropriate for re-usable medical instruments. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful and proper use reduces the risk of damage to the product and extends its serviceable life.
 b) The manufacturer does not recommend using any preservatives on medical devices.
6. Packaging
 1) Washed and dried devices shall be stored (*if possible*) in suitable stands placed in special sterilization containers. Separate items should be packed in a packaging intended for the recommended steam sterilization. Sterilization containers, item packaging and packaging process itself have to meet the requirements of ISO 11607 standards. The packaging procedure must be performed in controlled purity conditions. The device must be packed so that during its removal from the packaging, when used, there is no risk for its re-contamination
7. Sterilization
 1) The instrument is washed, disinfected and dried to undergo a sterilization process. The recommended sterilization procedure for non-sterile medical devices validated by **ChM** is autoclave vacuum steam sterilization (*with water vapor under overpressure*) with the following parameters:
 a) temperature: 134°C,
 b) minimum exposure time: 5 min,
 c) minimum drying time: 20 min.
 2) In order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease (*especially for surgical instruments that could come into contact with the central nervous system*), national health authorities in some countries recommend sterilization according to the following parameters:
 a) temperature: 134°C,
 b) minimum exposure time: 18 min,
 c) minimum drying time: 20 min.
 3) CAUTION:
 a) The sterilization process used in the medical facility must be validated and routinely controlled in accordance with the requirements of EN ISO 17665-1.
 b) Sterilization must be effective and in accordance with requirements of the EN 556-1 standard to ensure the required level of guaranteed sterility SAL 10^{-6} (*where SAL stands for Sterility Assurance Level*).
 c) Device must not be sterilized in the packaging in which it was delivered, except specially designed sterilization containers.
 d) The method of sterilization using ethylene oxide, gas plasma and dry heat should not be used, unless the Instructions for Use for the product contains sterilization recommendations using these methods.
 e) The sterilization temperature for plastic products (*PPSU, PEEK, PTFE, silicone*) cannot be higher than 140°C.

6 STORAGE

1. The devices should be properly stored. When storing surgical instruments, it is recommended that they never be stacked together. It may lead to damage of cutting edges (*nick or dull*) and/or initiation of corrosion centers. Instruments should be

stored in a clean and dry room, at room temperature and off the direct sunlight. If possible, instruments should be stored in suitable palettes placed into specially designed sterilization containers.

7 CALIBRATION

1. Regular calibration is required in case of torque wrenches, handles and connectors. Torque instruments are factory-calibrated, the nominal torque of a calibrated instrument is marked on the device (*e.g. 4 Nm*). To maintain a high level of safety and accuracy of operation of a torque instrument, it is necessary to follow the calibration deadline which is marked on the device.
2. Instrument calibration is performed by the manufacturer. Any attempt of unauthorized modifications to the construction or factory setting of the torque devices can lead to a potential injury or damage to the product and is prohibited.

8 COMPATIBILITY

1. **ChM** specialist instrument sets are designed for insertion of **ChM** implants. A specific, illustrated operating technique that describes the proper use of instruments included in the instrument set that is designed for particular implant system, is provided together with such instrument set. It is not allowed to combine **ChM** instruments with products from other manufacturers. The physician bears all responsibility for the use of the **ChM** instruments together with implants and instruments from other manufacturers.

If these instructions appear unclear, please contact the manufacturer, who shall provide all required explanations.
Updated INSTRUCTIONS FOR USE are available at the following website: **ifu.chm.eu**

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1 INDICATIONS

1. Cannula for bone cement is intended for use as an intermediary element during the injection of bone cement into a patient's bone. One end of the cannula is connected to a compatible screw through which cement is delivered into the bone. The other end is equipped with a standardized connector, suitable for a bone cement delivery set.
2. Surgical and orthopaedic instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.

2 DESCRIPTION

1. The package is equipped with the product label. The label (*as a primary label*) contains, among others:
 1) Logo **ChM** and the address of the manufacturer.
 2) Catalogue number (*REF*), e.g.: 40.XXXX.XXX, and device name and size.
 3) Production batch number (*LOT*), e.g. XXXXXXX.
 4) NON-STERILE sign - indicates non-sterile product.
 5) CE conformity mark.
 6) Information symbols (*described in the footer of this Instructions For Use*).
2. The unit package contains one piece of the product in non-sterile condition. Clear plastic bags are a typical packaging material. The products may also be supplied as a complete set (*arranged on palettes and placed into specially designed sterilization containers*).
3. Depending on the size or type of the product, the following information may be marked on its surface: manufacturer's logo, production batch no. (*LOT*), catalogue no. (*REF*), type of material and device size.

3 MATERIALS

1. Instruments are produced of corrosion-resistant steel. The protective layer (*passive layer*) against corrosion is formed on the surface of the device due to high content of chromium.
2. If the material of the device cannot be specified, please contact **ChM** representative.

4 WARNINGS AND PRECAUTIONS

1. Instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.
2. The surgeon should be familiar with all components of the device before use and should personally verify if all components and instruments are present before the surgery begins.
3. Before the procedure begins, all instruments should be carefully inspected for their condition and proper functioning. They should be undamaged and without corrosion. Damaged or corroded instruments should be immediately replaced. The use of bent, damaged or corroded instruments is not allowed.
4. Tissue structures close to the operative site must be protected.
5. Collision of the instrument with metal operating equipment, retractor or other device may cause damage that necessitates intraoperative replacement of that instrument.
6. Do not apply excessive force when using the instrument – it may lead to its permanent damage and, in consequences, to malfunction of the device.
7. Should breakage occur, the instrument parts must be removed and disposed of immediately in accordance with valid medical facility procedures.
8. In order to confirm the removal of all undesired metal fragments from the surgical field, intraoperative X-Ray examination is recommended.
9. In the case of suspected or documented allergy or intolerance to metallic materials, surgeon shall find out if the patient develops allergic reaction to the instrument material by ordering appropriate tests.
10. Improper, careless and inconsistent with the recommendations provided below handling of the instruments can lead to their chemical, electrochemical or mechanical damage which can adversely affect corrosion resistance and shorten the service life of the devices.
11. Instruments are intended only for specific procedures and must be used strictly according to their intended purpose. Use of instruments not in accordance with their intended purpose may lead to malfunction, accelerated wear and, in consequences, damage to the instrument.
12. Only use bone cement offered or recommended by **ChM**.
13. Consult the instructions for use delivered with the bone cement, cannulas and cement mixing/delivery system for warnings associated with the use this equipment.

14. Only use cannulas for cement delivery provided by **ChM** and intended for use with **ChM** implants. The use of cannulas from other manufacturers can result in cement leakage and associated complications.
15. Only use cement mixing/delivery devices that are provided or recommended by **ChM**. Assure that the components of the mixing/delivery device are compatible with cannula for bone cement delivery and are firmly connected prior to cement injection. Improperly secured connections could result in the unintended disconnection of components and cement leakage.
16. Do not attempt to force the injection of cement if excessive resistance is felt. Always determine the cause of the resistance and take appropriate action.
17. When using cement augmentation, attention must be paid not to exceed the working time of the cement prior to completion of cement delivery. When the cement working time is close to completion, a new cement package should be used.
18. After cement introduction is complete, the cannula for bone cement must be removed immediately from the screw. The cement setting while the cannula is still connected to the screw may lead to difficulty in removal of the device.
19. CAUTION: The cannula for bone cement delivery as well as the mixing/delivery device are designed for use with one package of bone cement only. If a second package of cement is needed, use a new cannula and mixing/delivery device.
20. CAUTION: has the cannula been used for cement delivery and working time for cement been exceeded, the hardened material can no longer be removed from the cannula, and therefore, the product should be considered worn out and disposed of in accordance with procedures and guidelines of a given facility.

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1. Prior to use of a non-sterile device, the following rules apply:
 1) The device must undergo cleaning, disinfection and sterilization procedures.
 2) Effective cleaning is a complicated procedure depending on the following factors: the quality of water, the type and the quantity of used detergent, the technique of cleaning (*manual, automated*), the proper rinsing and drying, the proper preparation of the device, the time, the temperature and carefulness of the person conducting this process, etc.
 3) The hospital facility remains responsible for the effectiveness of the conducted cleaning, packaging and sterilization processes with the use of existing equipment, materials and properly trained personnel.
2. Preparation for washing and disinfection (*for all methods*).
 1) If the instrument can be disassembled, it must be done before cleaning processes.
 2) Rinse under running water and remove surface debris using a disposable cloth, paper towel or plastic brushes (*nylon brushes are recommended*). Particular attention should be paid to openings and places difficult to be cleaned. Very dirty devices should be soaked in an aqueous solution of a detergent or a washing-disinfecting agent, e.g. neodisher® MediClean forte, at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 3) CAUTION: It is forbidden to use brushes made of metal, bristles or materials which could damage the product.
3. Cleaning and disinfection process.
 1) This Instructions for Use describes two **ChM**-approved cleaning and disinfection methods: manual with ultrasound cleaning and automated method. It is recommended to use automated cleaning and disinfection procedures (*in a washer-disinfector*).
 2) The chosen washing and disinfecting agents must be suitable and approved for use with medical devices. It is important to follow the instructions and restrictions specified by the producer of those cleaning agents. It is recommended to use aqueous solutions of washing-disinfecting agents with a pH value between 10.4 and 10.8. **ChM** used the following materials during the validation process of the described recommendations for cleaning and disinfection. It is allowed to use other materials than those listed below which may also give a comparable effect:
 a) detergent - Dr.Weigert (*producer*) neodisher® MediClean forte (*name of the detergent*);
 b) disinfectant - Dr.Weigert (*producer*) neodisher® Septo Active (*name of disinfectant*).
 3) To prevent product damage (*pitting, rust, discoloration*), do not use aggressive cleaning agents (*NaOH, NaOCl*), saline solutions and unsuitable cleaning agents.
 4) Where possible, it is recommended to use demineralized water to avoid the formation of spots and stains caused by chlorides and other compounds present in ordinary water.
 5) Manual with ultrasound cleaning.
 a) Equipment and materials: a device for ultrasound cleaning, soft, lint-free cloths, plastic brushes, syringes, aqueous solutions of cleaning agent.
 b) Manual cleaning: Initial manual cleaning must be performed prior to ultrasound cleaning.
 c) Rinse under running water until the product is visually clean. Use plastic brushes to remove heavy or large debris.
 d) Soak the product for at least 10 minutes in an aqueous solution of a detergent at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature,*

concentration, exposure time and water quality).

e) Rinse the product under cold water for at least 2 minutes, paying particular attention to the holes and places difficult to be cleaned.

f) Prepare fresh washing solution. Clean the surfaces and gaps of the product, carefully. Use suitable brushes to clean the holes. Clean the product immersed in the solution.

g) Rinse the product thoroughly under warm running water for at least 2 minutes, paying special attention to the gaps, blind holes, hinges and joints. When cleaning, use brushes and perform multiple reciprocating movements on the surface of the product.

h) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-h until the product is visually clean.

i) Ultrasound cleaning: prepare an aqueous cleaning solution at a temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the cleaning agent, in respect of temperature, concentration, exposure time and water quality*). Immerse fully the product in the aqueous cleaning solution and have it washed in ultrasounds for 15 minutes.

j) Rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.

k) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-k until the product is visually clean.

l) Use demineralized water for final rinsing of the device.

m) Dry the device thoroughly using disposable, soft, lint-free cloth or compressed air.

n) Prepare an aqueous solution of disinfecting agent at a temperature of 20+/-2°C using 20g of the agent per 1 liter of water. Immerse the product in the solution, exposure time – 15 min (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).

o) After the exposure time, rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.

p) The cannulated instruments should be treated using a compressed air or air supplied from the syringe.

q) Dry the device thoroughly. It is recommended to dry the product in a dryer at a temperature ranging from 90°C to 110°C, for 60 min.

r) Visually inspect the entire surface of the device.

s) CAUTION: If the obstruction in the cannula cannot be removed as indicated in the Instructions for Use, the device should be considered at the end of its useful life and should be discarded in accordance with facility procedures and guidelines.

6) The automated method using a washer - disinfector.

a) Equipment and materials: a washer - disinfector, aqueous solutions of cleaning agent.

b) Cleaning in the washer-disinfector must be preceded by a manual and ultrasound cleaning, following the procedure described in subsections c-h of paragraph 5.

c) CAUTION: The equipment used for washing/disinfection should meet the requirements of ISO 15883. Procedure of washing in the washer-disinfector shall be performed according to internal hospital procedures, recommendations of the washer-disinfector manufacturer, and instructions for use prepared by the washing-disinfecting agent manufacturer.

d) The device should undergo the process of machine washing in the washer-disinfector using the following cycle parameters: (*1*) - pre-washing in cold tap water, duration – 2 min; (*2*) - washing in an aqueous solution of cleaning agent at 55+/-2°C and pH of 10.4 - 10.8, duration – 10 min; (*3*) - rinsing under demineralized water, duration – 2 min; (*4*) - thermal disinfection in demineralised water at 90°C, minimal duration – 5 min; (*5*) - drying at the temperature ranging from 90°C to 110°C, duration - 40 min.

4. Inspection

1) Generally unmagnifed visual inspection under good light conditions is sufficient.

2) All parts of the product should be checked for visible dirt and corrosion. Particular attention should be paid to:

a) Places where dirt can be found, such as joints, latches, etc.

b) Holes, grooves and gaps the debris could have been pressed into during use.

3) Damaged or defective product cannot be approved for further use.

4) Prior to storage, the instrument must be checked for dryness.

5) CAUTION:

a) The manufacturer does not recommend using any preservatives on medical devices.

5. Packaging

1) Washed and dried devices shall be stored (*if possible*) in suitable stands placed in special sterilization containers. Separate items

should be packed in a packaging intended for the recommended steam sterilization. Sterilization containers, item packaging and packaging process itself have to meet the requirements of ISO 11607 standards. The packaging procedure must be performed in controlled purity conditions. The device must be packed so that during its removal from the packaging, when used, there is no risk for its re-contamination

6. Sterilization
 1) Washed, disinfected, and dried device shall undergo the sterilization process. The recommended method of sterilization is vacuum-type steam sterilization (*with water vapor under overpressure*):
 a) temperature: 134°C,
 b) minimum exposure time: 5 min,
 c) minimum drying time: 20 min.
 2) CAUTION:
 2) In order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease (*especially for surgical instruments that could come into contact with the central nervous system*), national health authorities in some countries recommend sterilization according to the following parameters:
 a) temperature: 134°C,
 a) The sterilization process must be validated and routinely monitored in accordance with the requirements of EN ISO 17665-1.
 b) Sterilization must be effective and in accordance with requirements of the EN 556-1 standard to ensure the required level of guaranteed sterility SAL 10^{-6} (*where SAL stands for Sterility Assurance Level*).
 b) minimum exposure time: 18 min,
 c) Device must not be sterilized in the packaging in which it was delivered, except specially designed sterilization containers.
 c) minimum drying time: 20 min.
 d) The method of sterilization using ethylene oxide, gas plasma and dry heat should not be used, unless the Instructions for Use for the product contains sterilization recommendations using these methods.

6 STORAGE

1. The devices should be properly stored. When storing surgical instruments, it is recommended that they never be stacked together. It may lead to damage of cutting edges (*nick or dull*) and/or initiation of corrosion centers. Instruments should be stored in a clean and dry room, at room temperature and off the direct sunlight. If possible, instruments should be stored in suitable palettes placed into specially designed sterilization containers.

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8 RE-STERILIZATION

1. It is permitted to re-sterilize a device, when its sterile packaging has been damaged or opened. In this case, the product should be washed and sterilized in the manner described in the chapter CLEANING, DISINFECTION, STERILIZATION.
2. It is strictly forbidden to re-sterilize and re-use cement constituents, the cannula for cement delivery and device for cement mixing/delivery.

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1. INTENDED USE 11
2. CLINICAL BENEFITS 11
3. PATIENT TARGET GROUP 11
4. INTENDED USERS 11
5. SPECIFICATIONS AND EQUIPMENT 11
6. ACCESSORIES 12
7. CONSTRUCTION 13
8. USAGE 13
9. CLEANING AND DISINFECTION 15
10. STORAGE 15
11. WARRANTY AND SERVICE 15
12. THE LIFETIME AND DISPOSAL 15
13. DISPOSAL 16

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Read the Instructions for Use prior to the first use.

1. INTENDED USE

The pump is used to supply compressed air to tourniquets used to control the blood flow in the extremities in order to produce bloodless surgical field. The pump is a stand-alone device that does not need to be powered by electricity or by compressed air from central systems of medical gases distribution.

2. CLINICAL BENEFITS

Clinical benefits of using a pump with a tourniquet include blood loss reduction during surgery and making the procedure easier to perform.

3. PATIENT TARGET GROUP

The device can be used in surgical procedures performed on infants, children, adolescents and adult patients.

4. INTENDED USERS

The devices are designed and intended for use by professional users - specialists in orthopedics and traumatology who are familiar with relevant surgical principles and techniques and who have practical skills in using **ChM**'s compatible instrumentation.

5. SPECIFICATIONS AND EQUIPMENT

Output pressure	0÷750mmHg
Dimensions *(H/W/D)*	39cm x 6.5cm x 3.5cm
Weight	0.5kg

Equipment supplied with the device.

	name	Pcs.
	Pump air hose cat. no. 30.0077.200	1
	Pneumatic tube male cat. no. 30.0022.100	1

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
Tourniquets are additional equipment of the pump and are offered as separate devices.

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7. CONSTRUCTION

1. Manometer
2. Knob of pressure reducing valve
3. Tourniquet connection port
4. Pump handle
5. Pump air hose
6. Pneumatic tube male

8. USAGE

NOTE: Before each use, the pump should be carefully checked. Should any defects or doubts about the technical condition appear, the pump must not be used. Contact your representative or **ChM** directly.

NOTE: It is unacceptable, during pumping, to:

- plug *(e.g. with a finger)* the hole of the tourniquet connection port *(3)*
- bend or clamp the spiral air hose *(5)*, or plug its end when connected to the pump.

Such actions may cause excessive pressure increase in the pneumatic system of the pump *(that exceeds the permissible range of 750mmHg)* and damage to the manometer.

Place the chosen **ChM** tourniquet on an operated limb as intended. Install the tourniquet air hose to the pump by screwing the threaded connection up into the tourniquet connection port *(3)* in the pump.

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Should a non-**ChM** tourniquet be used, the tourniquet hose must be provided with the pneumatic tube male **[30.0022.100]** *(6)* supplied with the device.


NOTE: **ChM** shall not be liable for the correctness of the connection of the pneumatic tube male 30.0022.100 to a non-**ChM** tourniquet.

INFORMATION: If the length of the tourniquet air hose is insufficient, use the pump air hose **[30.0077.200]** supplied with the device.

Turn the knob *(2)* up clockwise to close the pressure reducing valve. Do not use excessive force.

Perform gentle reciprocating movements with pump handle *(4)* and inflate the tourniquet to the desired pressure, the value of which will be indicated on the monometer *(1) (maximum 750mmHg)*.

To reduce or completely release the pressure from the tourniquet, turn the knob *(2)* counter-clockwise.

INFORMATION: The pressure drop depends on the angle by which the knob will be turned.

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9. CLEANING AND DISINFECTION

The outer surfaces of the device should be cleaned with a soft, moistened *(not dripping)* cloth with a mild detergent solution then wiped dry. In the case of contaminating the device with infectious agents *(blood, body fluids, etc.)*, disinfection using disinfectants with a neutral pH is recommended. In such cases, soak a soft cloth in disinfectant, wipe the device and allow it to dry.

NOTE: Do not allow the detergents or disinfectants to enter inside the device. Should the detergents or disinfectants enter inside the device, send it to the manufacturer for cleaning processes.

NOTE: The pump should NOT be sterilized.

10. STORAGE

The device should be stored in a clean and dry room, off the direct sun, at the temperature of 15 ÷ 30°C and humidity of 20 ÷ 70%.

11. WARRANTY AND SERVICE

When servicing, warranty and post-warranty repairing is concerned, please contact the manufacturer's representative or contact **ChM** directly:

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, Poland
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu

NOTE: Do not replace parts of the pump by yourself! The changes made will render the warranty void.

The warranty does not include product damage due to: negligence, inappropriate usage, usage not in accordance with the intended purpose, mechanical damage.
Accessories and spare parts are presented in **ChM**'s catalogues.
Warranty and non-warranty services are to be handled solely by the manufacturer - **ChM**.

NOTE: It is recommended to perform technical inspection and calibration of the device at the manufacturer's site at least once every 12 months.

12. THE LIFETIME AND DISPOSAL

The **ChM** sp. z o.o. does not define the maximum number of uses for the pump. The useful life of the device depends on many factors including the method and duration of each use, applied cleaning and disinfection and the handling between uses.

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1. INDICATIONS FOR USE

Tourniquet, in combination with an external source of compressed air or nitrogen, is used to control blood flow in the limb and produce a bloodless surgical field. Tourniquet control unit or hand pump (*compressed air*) produced by **ChM** may be examples of external source of compressed air or nitrogen.

2. CONTRAINDICATIONS

- Open fracture of the lower limb.
- Post-traumatic hand reconstruction.
- Severe crush injuries.
- Elbow surgery (*with excessive swelling*).
- Severe hypertension.
- Skin grafts in which all bleeding points must be distinguished.
- Compromised circulation e.g. peripheral arterial disease.
- Diabetes.
- Sickle cell disease is a relative contraindication against use of a tourniquet (*See PRECAUTIONS*).
- Avoid using a tourniquet in patients undergoing secondary or delayed treatment after fixation.

NOTE:
The abovementioned list of contraindications is not exhaustive.

3. CLINICAL BENEFITS

- Clinical benefits of using a tourniquet include blood loss reduction during surgery and making the procedure easier to perform.

4. PATIENT TARGET GROUP

- The devices can be used for treatment of infants, children, adolescents and adult patients.

5. INTENDED USERS

- The devices are designed and intended for use by professional users who are familiar with relevant surgical principles and techniques and who have practical skills in using **ChM**'s instrumentation.

6. USE

Depending on the area of application and the type of surgical procedure, tourniquets of different structure may be used so that they are more effective. **ChM** offers several types of tourniquets differing in shape, length and width, allowing the surgeon to choose a tourniquet appropriate to the location and type of procedure to be performed.

6.1. PLACE OF USE

The tourniquet shall be applied on a limb at the place of a large muscle mass which helps to protect vessels and nerves from damage.

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The typical location of the tourniquet:

- arms and thighs - where the limb has the the largest circumference closest to the operation field
- forearm - the central part of the forearm
- shank - the proximal edge of the tourniquet should be placed on the largest circumference of the calf
- foot - 1/3 of the lower part of the shank with the distal edge of the tourniquet passing above the ankle

6.2. SINGLE AND DUAL TOURNIQUETS

Depending on the surgical procedure, single or dual tourniquets may be used. Single tourniquets are mostly used with general anesthesia, brachial plexus blocks and central blocks.

Dual tourniquets are mostly used with local intravenous anesthesia (*also called Bier's block*). Thanks to the independent inflation of the proximal and distal section of the tourniquet, the surgical procedure is more patient-friendly, especially during longer operations.

6.3. LENGTH

The length of a tourniquet is chosen on the basis of the length of the limb circumference measured at the point where the tourniquet is to be applied. Knowing the circumference of the limb and the place of application, select the tourniquet from the table (*Tab. 1*).

NOTE:
Use of excessively short or long tourniquets can lead to its unstable positioning on a limb, problems with full occlusion of blood vessels and damage to the skin.

6.4. WIDTH

As wide as possible tourniquets should be applied. A wider tourniquet requires lower occlusion pressure than a narrow one what reduces the risk of complications associated with its use.

6.5. TOURNIQUET PRESSURE

Tourniquet pressure used during surgery depends on many variables which include: the patient's age, skin condition, a place the tourniquet is applied, the shape and size of the tourniquet. According to the publicly available literature data, occlusion pressure:

- for the upper limb should be about 75÷100mmHg above systolic blood pressure,
- for the lower limb should double the value of systolic blood pressure.

6.6. PREPARING A PATIENT FOR SURGERY

Each patient that is scheduled for surgery with the use of a tourniquet, should undergo a preoperative evaluation including, inter alia: determination of the physical condition of the patient, a medical history and medication review. On the basis of the preoperative evaluation, possible complications that may occur during the surgery are determined and whether the patient is allowed to be treated with the use of the tourniquet. The measurement of the blood pressure and the length of the limb at the point where the tourniquet will be applied should also be performed preoperatively.

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6.7. TOURNIQUET USE

Carefully check the condition and tightness of the tourniquet prior to each use. Should a defect or doubts as to the condition of the tourniquet appear, it is unacceptable to continue the use of the tourniquet. Service and warranty repairs may be performed only by the manufacturer - **ChM**.

Apply disposable protective cuff on the selected preoperatively tourniquet that would protect the tourniquet from blood, body fluids or other infectious agents contamination. The skin under the tourniquet should be protected from direct exposure to that tourniquet by applying padding underneath, keeping in mind though that excessively thick padding may increase occlusion pressure. Moreover, do not allow any fluids to get under the tourniquet because they can cause skin injuries. When applying the tourniquet, make sure that the surfaces that come into direct contact with the padding (*skin*) were not wrinkled because it can lead to skin damage after inflating the tourniquet. Applied tourniquet should adhere to the limb with its entire perimeter. The tourniquet feeder should be placed outside of the limb and operating field. Avoid accidental pressure of the feeder on nerves or its bend. Connect the tourniquet to a tourniquet control unit and start inflating until the established occlusion pressure is reached. Inflating the tourniquet should proceed quickly and smoothly so that arteries and veins are closed at the same time.

As a general accepted rule, the maximum time of the tourniquet being inflated (*without its deflating*), in the case of a healthy adult patient, should not exceed 120 minutes. During the surgery, to reduce the risk of complications, intraoperative control of blood pressure, the pressure in the tourniquet and the time the tourniquet is applied should be conducted.

Before deflating the tourniquet, it is recommended to close the wound. Deflating the tourniquet before wound closure increases blood loss. Deflating the tourniquet should proceed quickly and smoothly. After deflating, immediately remove the tourniquet and the padding from the limb so that they will not block the return of venous blood and resultant blood accumulation in the operating field.

7. SIDE EFFECTS

After tourniquet use, dull ache in the limb may appear.
Prolonged ischemia caused by use of the tourniquet can lead to temporary or permanent damage to tissues, blood vessels and nerves and can cause changes in blood clotting.

Excessively high pressure in the tourniquet can cause paralysis of limbs.
Excessively low pressure in the tourniquet can cause passive congestion of limbs with possible irreversible loss of function.

Intraoperative bleeding may be caused by:

- blood remaining in the limb due to insufficient emptying the limb of blood
- the use of improper occlusion pressure
- blood penetrating from the nourishing vessels of the long bones.

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8. CONSTRUCTION

The tourniquet is made of polyamide fabric.

The latex bladder connected with feeder is placed in the tourniquet pocket.

The quick-fit port connecting the tourniquet to the tourniquet control unit is placed at the end of the feeder. The tourniquet is reinforced by silicone insert to obtain greater stiffness.

There are special Velcro and protective tapes on the inner side of the cuff to protect the tourniquet against loosening when applied on the limb and inflated.

9. DISINFECTION AND CLEANING

The tourniquet should be cleaned by hand with a soft brush damped in warm water and a small amount of soap causing no foaming, unless the tourniquet got contaminated with infectious agents, e.g. blood. Only exterior surfaces of the rubber feeder and cuff should be cleaned, avoiding the excessive soaking of the tourniquet material. It is not allowed to: immerse the tourniquet in the water or aqueous solutions of washing and disinfecting agents, machine-wash the tourniquet or to pull out the inner latex bladder from the cuff. Wipe the washed tourniquet with a damp lint-free cloth and air-dry (*do not roll*). When dried, the tourniquet may be wiped with cloths soaked in isopropyl alcohol.

Should the tourniquet have contact with infectious agents, e.g. blood, for cleaning and disinfection use mild washing and disinfecting agents approved for use with medical devices. Follow the instructions provided by the manufacturer of these agents. It is recommended to use aqueous solutions of washing or washing-disinfecting agents with a pH value between 7 and 10.8.

1. Remove blood and/or dirt from the product using disposable cloths or paper towels.
2. Prepare a fresh detergent solution using e.g., Neodisher MediClean Forte from Dr.Weigert at a temperature of 40°C ± 2, using 15ml of the agent per 1l of water. Thoroughly clean the surface of the cuff with a soft, lint-free cloth.
3. Dry the device thoroughly using disposable, soft, lint-free cloth.
4. Prepare an aqueous solution of disinfectant, e.g. Neodisher Septo Active from Dr.Weigert of 2% concentration (*20g of agent per 1 liter of water*) and temperature 20°C ± 2. Thoroughly clean the surface of the cuff with a soft, lint-free cloth.
5. After disinfecting, clean the tourniquet with a damp cloth to remove residual solution that may have an allergic and irritating effect and shorten the life of the product. Dry the tourniquet in the open air - do not roll.
6. Clean the tourniquet with cloths soaked in isopropyl alcohol.

NOTE:
The tourniquet, when used in surgery, should always be placed in a protective cuff which protects it from contamination by infectious agents.

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10. STORAGE

The tourniquet shall be stored unwound or slightly twisted in a roll, in a dry, not in direct sunlight place. The tourniquet storage should be preceded by its washing and thorough drying.

11. WARRANTY AND SERVICE

The manufacturer provides annual warranty for the usage of the tourniquet provided that the above principles are observed.

Unauthorised changes introduced by the user might render the warranty void.

ChM sp. z o.o. does not specify the maximum number of cycles for use of tourniquets. The tourniquet life is dependent on many factors, including the duration of each use, the pressures used, re-processing.

Each time before re-use, inspect the functionality of the product. In case of damage, the tourniquet must not be used for surgery.

12. DISPOSAL

1. If the device is damaged, it should be removed immediately and safely disposed of in accordance with the procedures of the medical facility.
2. The medical devices are classified as medical waste. The safe disposal of medical waste should be carried out in accordance with applicable national and local legislation.

Contact details:

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu
www.chm.eu

CAUTION: Any serious incident involving the device shall be reported to the manufacturer and to the competent authority of the member state in which the user or patient resides.

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10. INTENDED USE

The stand is used as a mobile platform for the **ChM** family of tourniquet control units, as well as the accessories needed to be used with the device, such as tourniquets, skin protection padding, power cords and other accessories.

11. CLINICAL BENEFITS

Mobile stand facilitates surgical procedures that require the use of tourniquet systems.

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13. INTENDED USERS

The device is designed and intended for use by professional users who are familiar with relevant surgical principles and techniques and who have practical skills in using **ChM**'s instrumentation.

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The stand is equipped with the basket and hooks for storing additional equipment used with the tourniquet control unit.

The stand wheels with rubber tread make the stand move smoothly and quietly. Wheel locks provided allow for stand immobilization to prevent its uncontrolled movement. The wheel base and stand body are made of light, yet very durable aluminum alloys, so that the stand is both light and stable.

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15. TOURNIQUET CONTROL UNIT INSTALLATION AND STAND HANDLING

There are two screws **1** (provided with the stand) that are used to connect the tourniquet control unit with the stand. Insert the screws into the holes of the stand top and then in the bottom threaded holes of the tourniquet control unit.

NOTE:
When installing the tourniquet control unit onto the stand, hold the device firmly until at least one of the screws is tightened.

Move the stand by holding the handle at the top of the tourniquet control unit.

In order to immobilize the stand, push down the locks of the wheels (the locks are provided for 2 out of 5 wheels). Push the locks up to unlock the stand.

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16. CLEANING AND DISINFECTION

The stand should be cleaned with a soft cloth soaked (not dripping) with a mild detergent solution then wiped dry. In the case of contaminating the stand with infectious agents (blood, body fluids, etc.), wipe the stand with a soft cloth soaked with neutral pH disinfectant, then with a cloth soaked with mild detergent and finally wipe dry.

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18. SERVICE

When servicing, warranty and post-warranty repairing is concerned, please contact the manufacturer's representative or **ChM** directly:

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu

NOTE:
Do not repair or replace parts of the device by yourself!
The changes made will render the warranty void.
Any serious incident involving the device shall be reported to the manufacturer and to the competent authority of the member state in which the user or patient resides..

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EN STERILIZATION CONTAINERS SYSTEM

1. INTENDED USE OF **CHM** STERILIZATION CONTAINERS SYSTEM 4

2. COMPONENTS OF **CHM** STERILIZATION CONTAINERS SYSTEM 4

2.1. Box 5

2.2. Perforated lid 5

2.3. Single-use paper filter 5

3. MAINTENANCE AND CONTROL 6

4. CLEANING AND DISINFECTION 9

4.1. Disassembly of components prior to cleaning 10

4.2. Cleaning and disinfection process 11

4.2.1 Manual with ultrasound cleaning. 11

4.2.2 Automated method using a washer - disinfector. 12

5. STERILIZATION 13

6. TRANSPORTATION AND STORAGE 14

7. STORAGE CONDITIONS 14

8. TROUBLESHOOTING 14

Symbols for the user

Before first use, read the Instructions for Use.

This symbol stands for a warning against a threat to life or health.

This symbol stands for a warning against the risk of injury or product damage.

Before using the product, read carefully the Instructions for Use to avoid improper use, and product damage which may lead to loss of product warranty.

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1. INTENDED USE OF CHM STERILIZATION CONTAINERS SYSTEM

Prior to the first use, read carefully this manual and keep it in a place accessible to all users of ChM sterilization containers system.

ChM sterilization containers system is designed for transportation, packaging, sterilization and storage *(before and after sterilization)* of medical devices manufactured by **ChM** are intended for use in aseptic conditions.

The sterilization container system shall be used for vacuum steam sterilization *(overpressed steam)* of implants and instruments.

ChM sterilization containers system shall only be used as recommended in this manual. The provided principles of safe use and maintenance of the product must be strictly observed. Damaged, worn-out and expired constituent elements should be replaced immediately with new original ones. **ChM** sterilization containers system with all accessories can only be used by experienced and qualified users. Keep these operating instructions in a place accessible to related personnel.


Do not use **ChM** sterilization containers system in sterilization units with gravity venting.

2. COMPONENTS OF CHM STERILIZATION CONTAINERS SYSTEM

Description

1) Box
2) Perforated lid
(for systems with paper filter)
3) Handle
4) Window / place for label
5) Lid closing lock
6) Lid closing flexible mechanism
7) Silicone lid gasket
8) Paper filter
9) Filter holder
10) Filter holder locking mechanism

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ChM sterilization containers systems *(box and lid)* are made of aluminum alloy. Aluminum, by electrochemical surface treatment *(anodizing)*, produces a resistant oxide coating that can be colored or stay natural *(a silver-gray colour)*.
Products made of anodized aluminum prove good corrosion resistance properties; however, contact with strongly alkaline cleaning and disinfecting agents as well as solutions containing iodine or certain metal salts should be avoided since they can cause chemical damage and/or destruction of treated aluminum surface.

2.1. BOX

The non-perforated boxes, for sterilization process, must be used with perforated lids. The non-perforated box with non-perforated lid can only be used for transportation.

2.2. PERFORATED LID

2.3. SINGLE-USE PAPER FILTER

Single-use elements must be replaced before each sterilization cycle.

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
3. MAINTENANCE AND CONTROL

Prior to each use *(sterilization process)* **of sterilization containers system, check the following:**

a) Condition of the gasket in the lid,
b) Container locking system,
c) Filter holder gaskets,
d) Single-use paper filter

3.a) Lid gasket control


Should the lid gasket be mechanically damaged *(e.g. deformed, cracked, torn)* or should there be no gasket at all; or should the box or lid be damaged within their connection *(cracked or bent sheet)*, the effectiveness of the sterilization process and the maintenance the sterility in storage until use may be impaired. If any damage to any of the above-mentioned elements occurs, secure immediately the container against re-use and contact a **ChM** representative for further proceedings.

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3.b) Control of sterilization container locking system

Before each sterilization, check the function of the sterilization container locking system. Place the lid on the box. Close the locks.

Close the locks on both sides of the lid, simultaneously. Uneven application of force on both sides of the lid may prevent the locks from being closed properly and the sealing being not secure, what results in incorrect sterilization process.

If the closing system of the sterilization container is not working properly, such a container cannot be used in the sterilization process.

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Make sure the lid is tightly closed so that it does not move over the box.

If the lid moves on the box, secure immediately the container against re-use and contact a **ChM** representative for further proceedings.

3.c) Control of the gasket of filter holder


Check that the filter holders of the sterilization containers and filter holder locking mechanisms are functioning correctly before each use.

As show above, remove the filter holder from the lid by pressing the latch on it.

Check the gaskets on the inside edges of the filter holder. Do not use the filter holders with torn or deformed gaskets.

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Press your finger on the button until a *"click"* sound is heard. Try to rotate the filter holder by hand. If the filter holder does not rotate, the mechanism is seated.

3.d) Control of single-use paper filter

Disposable filters *(paper filters)* must be replaced before each sterilization cycle. They must be clean, dry and free of foreign bodies, they must not be worn, expired or damaged.

4. CLEANING AND DISINFECTION

Prior to cleaning the sterilization containers, remove the stand*(s)* with medical devices from the inside. Products made of anodized aluminum prove good corrosion resistance properties; however, contact with strongly alkaline cleaning and disinfecting agents as well as solutions containing iodine or certain metal salts should be avoided since they can cause chemical damage and/or destruction of treated aluminum surface. Effective cleaning is a complicated procedure depending on the following factors: the quality of water, the type and the quantity of used detergent, the technique of cleaning *(manual, automated)*, the proper rinsing and drying, the proper preparation of the device, the time, the temperature and carefulness of the person conducting this process, etc.

The hospital facility remains responsible for the effectiveness of the conducted cleaning, packaging and sterilization processes with the use of existing equipment, materials and properly trained personnel.

Immediately after use, remove from sterilization containers blood and other contaminants with disposable cloth or paper towels. Additionally, it is recommended to rinse the sterilization containers under running water or to place it in the aqueous disinfectant solution. Do not let blood, tissues, body fluids or other biological impurities dry out on the surface of the device.

Preparation for washing and disinfection *(for all methods)***.**

The used sterilization containers should be reprocessed as soon as possible.
Before cleaning, disassemble the components of sterilization containers as indicated in section 4.1.
Rinse under running water and remove surface debris using a disposable cloth, paper towel or plastic brushes *(nylon brushes are recommended)*. Particular attention should be paid to places difficult to be cleaned. Very dirty devices should be soaked in an aqueous solution of a detergent or a washing-disinfecting agent,

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
e.g. neodisher® MediClean forte, at temperature of 40+/- 2°C and pH of 10.4-10.8 *(follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality)*.

It is forbidden to use brushes made of metal, bristles or materials which could damage the product.

4.1. DISASSEMBLY OF COMPONENTS PRIOR TO CLEANING

Remove all detachable parts of the container prior to cleaning.

- Remove the lid from the box,

- Remove paper labels and metal badges from the box,

- Remove the filter holders and filters from the lid.

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4.2. CLEANING AND DISINFECTION PROCESS

These Instructions for Use describes two **ChM**-approved cleaning and disinfection methods: manual with ultrasound cleaning and automated method. It is recommended to use automated cleaning and disinfection procedures *(in a washer-disinfector)*.
The chosen washing and disinfecting agents must be suitable and approved for use with medical devices. It is important to follow the instructions and restrictions specified by the producer of those cleaning agents. It is recommended to use aqueous solutions of washing-disinfecting agents with a pH value between 10.4 and 10.8. **ChM** used the following materials during the validation process of the described recommendations for cleaning and disinfection. It is allowed to use other materials than those listed below which may also give a comparable effect:

a) detergent - Dr.Weigert *(producer)* neodisher® MediClean forte *(name of the detergent)*;
b) disinfectant - Dr.Weigert *(producer)* neodisher® Septo Active *(name of disinfectant)*.

To prevent product damage *(pitting, rust, discoloration)*, do not use aggressive cleaning agents *(NaOH, NaOCl)*, saline solutions and unsuitable cleaning agents.
Where possible, it is recommended to use demineralized water for rinsing to avoid the formation of spots and stains caused by chlorides and other compounds present in ordinary water.

4.2.1 MANUAL WITH ULTRASOUND CLEANING.

First, remove all dismountable elements of the sterilization container *(see section 4.1)*.
Equipment and materials: a device for ultrasound cleaning, soft, lint-free cloths, plastic brushes, syringes, aqueous solutions of cleaning agent.
Manual cleaning: Initial manual cleaning must be performed prior to ultrasound cleaning.

a) Rinse under running water until the product is visually clean. Use plastic brushes to remove heavy or large debris.
b) Soak the product for at least 10 minutes in an aqueous solution of a detergent at temperature of 40+/- 2°C and pH of 10.4-10.8 *(follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality)*.
c) Rinse the product under cold water for at least 2 minutes, paying particular attention to the holes and places difficult to be cleaned.
d) Prepare fresh washing solution. Clean the surfaces and gaps of the product, carefully. Use suitable brushes to clean the holes. Clean the product immersed in the solution.
e) Rinse the product thoroughly under warm running water for at least 2 minutes, paying special attention to the gaps, hinges and joints. When cleaning, use brushes and perform multiple reciprocating movements on the surface of the product.
f) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections a-f until the product is visually clean.
g) Ultrasound cleaning: prepare an aqueous cleaning solution at a temperature of 40 +/- 2°C and pH of 10.4 - 10.8 *(follow the information contained in the instructions prepared by the manufacturer of the cleaning agent, in respect of temperature, concentration, exposure time and water quality)*.

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h) Immerse fully the product in the aqueous cleaning solution and have it washed in ultrasounds for 15 minutes.

i) Rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleVisually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-j until the product is visually clean. Use demineralized water for final rinsing of the device.

j) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-j until the product is visually clean. Use demineralized water for final rinsing of the device.

Dry the device thoroughly using disposable, soft, lint-free cloth or compressed air.
Prepare an aqueous solution of disinfecting agent at a temperature of 20+/-2 °C using 20g of the agent per 1 liter of water. Immerse the product in the solution, exposure time – 15min *(follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality)*. After the exposure time, rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned. Dry the device thoroughly. It is recommended to dry the product in a dryer at a temperature ranging from 90°C to 110°C, duration - 60min. Visually inspect the entire surface of the device.

To remove heavy or large debris during initial washing, do not use cleaning tools that can scratch the surface. Only use brushes with soft bristles. Scratches on the surface of the sterilization container may damage the anodized surface of the product.

4.2.2 AUTOMATED METHOD USING A WASHER - DISINFECTOR.

The automated cleaning of the sterilization container is only recommended if the washer-disinfector has a special program for cleaning aluminum containers.

Cleaning in the washer-disinfector must be preceded by a manual and ultrasound cleaning, following the procedure described in subsections a-j of paragraph 4.2.1. Before automated cleaning, disassemble the components of sterilization containers as indicated in section 4.1. Equipment and materials: a washer - disinfector, aqueous solutions of cleaning agent.

CAUTION: The equipment used for washing/disinfection should meet the requirements of ISO 15883. Procedure of washing in the washer-disinfector shall be performed according to internal hospital procedures, recommendations of the washer-disinfector manufacturer, and instructions for use prepared by the washing-disinfecting agent manufacturer.

The device should undergo the process of machine washing in the washer-disinfector using the following cycle parameters: *(1)* - pre-washing in cold tap water, duration – 2min; *(2)* - washing in an aqueous solution of cleaning agent at 55+/-2°C and pH of 10.4 - 10.8, duration – 10min; *(3)* - rinsing under demineralized water, duration – 2min; *(4)* - thermal disinfection in demineralised water at 90°C, minimal duration – 5min;

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(5) - drying at the temperature ranging from 90°C to 110°C, duration - 40min.

Sterilization containers should not be washed together with other medical devices. In addition, the box should be placed upside down in the washer-disinfector to avoid water accumulation inside the box.

5. STERILIZATION

Before sterilization process, check all the elements of the sterilization container system according to section 3 of this manual. Make sure the expiry date of the disposable paper filter is not exceeded and that the openings in the lid are neither obstructed nor clogged as this may damage or even destroy the sterilization container.

The recommended sterilization method - vacuum steam sterilization *(overpressed steam)*:

a) temperature: 134°C;

b) minimum exposure time: 5 min;

c) minimum drying time: 20 min.


The sterilization process should comply with EN 285 and EN ISO 17665-1.

After the sterilization process, allow the sterilization container to cool down to the room temperature. Transportation of the hot container, due to the large temperature difference, may result in condensate accumulation inside the sterilization container.

Caution should be taken due to the risk of burns. Immediately after sterilization, the containers are hot.

To ensure an effective sterilization process, unobstructed flow of sterilizing agent inside the container should be guaranteed. Products sterilized in the sterilization container must not come into contact with the lid.

After effective sterilization, the color of the chemical indicator changes from blue to black.

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6. TRANSPORTATION AND STORAGE

Always hold sterilization containers by handles, use protective gloves, if possible. Do not hold the sterilization container by the closing locks of the closing mechanism.

Holding the sterilization container by the closing locks of the closing mechanism may result in loss of tightness between the lid and the box, and consequently loss of sterility of the products inside.

7. STORAGE CONDITIONS

Sterilized medical devices can be stored in the sterilization container until they are used, but no longer than for a period of 6 months, as long as the storage conditions of these containers *(after sterilization)* are met:

a) storage temperature 15-26°C;

b) relative humidity: 30-50%.

8. TROUBLESHOOTING

Problem	Possible causes	Solution
Excessive moisture inside the sterilisation container after sterilization	The temperature of the products before sterilization is too low.	Sterilize products at room temperature, not lower than 20°C.
	The drapes are very damp.	The drapes to be placed in the sterilization container should be dry.
	Sterilization containers have not been properly placed in the sterilization unit.	The total height of the charge should be 45 cm, with the heavier sterilization containers at the bottom.
	Sterilization unit does not meet the EN 285 requirements.	When using the sterilization unit, consider the manufacturer's recommendations. Check the drying time and the steam quality. Calibrate if necessary.
	The idle loop and vacuum test are not performed daily before sterilization.	Perform idle cycle and vacuum tests before sterilization. Make sure the sterilization unit is working properly
Color of indicator has not changed after sterilization.	Incorrect sterilization process.	Have the manufacturer repair the sterilization unit
Deformation of the sterilization container	Filter holes have been obstructed during sterilization	Do not block the filter holes in the lid

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Lid lock failure	Sterilization container, instead of handles, has been held by the closing locks of the lid.	Contact company representative. A sterilization container with a damaged lid lock cannot be used for sterilization
	Closing locks of the lid has been mechanically damaged	

Any modifications introduced to the product may render the warranty void.

For warranty, please contact the company representative or directly the manufacturer of the device - **ChM**.

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, Poland
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu

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