**Pacjent**

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everall7

# VILLACRYL REPAIR FREE

Cold-curing acrylic resin for denture repairs

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**Cold-curing acrylic resin for denture repairs**

**Product characteristics**

Acrylic resin for cold-curing polymerization in powder and liquid form [illegible] indirect relining of removable dentures. Without MMA (methyl methacrylate).

**Classification according to ISO 20795-1: Type 2, Class 1**

**Composition**

Powder: polymethyl methacrylate (PMMA), benzoyl peroxide (BPO), [illegible]henyl-barbituric acid (BPBA), titanium dioxide, pigments, silica, viscose fibers*.

Liquid: ethyl methacrylate (EMA), ethylene glycol dimethacrylate (E[illegible]imethyl-p-toluidine (DMPT), hydroquinone (HQ).

*Depending on the color variety, the composition of the product may v[illegible]terms of quantity and quality.

**Operating parameters:** The technical data necessary for proper pro[illegible]rials are given in Table 1.

**Intended use of the products**

**Indications:** repairs of complete, partial dentures, acrylic parts of f[illegible]ures and acrylic orthodontic appliances.

**Contraindications**

Do not use the product to persons allergic to methyl methacrylate [illegible]e for traditional MMA-based denture repair materials.

Acrylic monomers and benzoyl peroxide are allergenic substances. The[illegible]out of the denture in the oral environment and affect the patient's body. Do not use the products to [illegible]diagnosed immunodeficiency concerning allergens.

If you notice any allergic symptoms, stop using the device and consult

Do not perform repairs, reline or other work directly in the patient's or

**Target group(s) of patients and intended users of the devices**

**Users:** Dental technicians making dentures.

**Patients:** Denture users. Users of acrylic dental appliances, children, youth, adults.

**Side effects**

Irritation, redness or signs of hypersensitivity, including the occurrence of anaphylactic shock, may occur in the case of allergy to any of the components of the products.

Oral irritation and inflammation caused by fungi and bacteria adjacent to the repaired denture are possible in the absence of proper hygiene.

**Clinical benefits**

Restoring the biting and chewing function of dentures damaged during use.

Improving the comfort and well-being of denture users who have suffered mechanical damage during use.

Restoring the correct functioning of dental appliances damaged during use, in order to correct the dentition or prevent secondary movement of the teeth being corrected.

**Range***

**Color shades**

V2 – milk pink veined T4 – pink
V3 – dark pink veined 0 - transparent
V4 – pink veined

**Packaging**

| | Powder + liquid set | Powder | Liquid |
|---|---|---|---|
| **Villacryl Repair Free** | 100 g + 50 ml | 100 g<br>1 kg | 50 ml<br>200 ml<br>500 ml<br>1 l |

*Please contact your distributor to obtain information about the availability of the range.

**Method of Use (Step by step instructions)**

**Initial operations**

Read the instructions carefully before using the product.

Before repairing the denture, make sure that the patient for whom the denture is repaired does not have a diagnosed immunodeficiency concerning allergens.

It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering face mask suitable for liquids for respiratory protection, protective glasses and work clothing.

**Cleaning and disinfecting the denture before repair**

Remove possible deposits of tartar and bacterial plaque from the surface of the repaired denture. Wear protective gloves and masks. It is recommended to use **Denture Cleaner** by **Everall7** for cleaning the denture. Wash and disinfect the cleaned denture.

**Preparing**

**Preparing the denture for repair**

Use acrylic or composite teeth when attaching teeth to partial dentures. The teeth should have the mucous surface part ground during setting.

Prepare the repaired denture base thoroughly using the standard method.

Clean the damaged denture base and moisten it with a monomer.

**Material preparation**

Shake the powder and liquid containers 2 or 3 times as shown on the packaging and product labels to facilitate powder dispersion in the monomer.

**Measuring and mixing of ingredients**

**Recommended mixing ratio: (Table 1).**

The stated mixing ratio must be strictly observed. Failure to comply with these values may lead to the formation of a porous structure at the denture repair site, a reduction in its mechanical properties and an increase in the residual monomer content.

**Component mixing**

Pour an appropriate amount of liquid into a chemically resistant dish (ceramic, porcelain, glass or stainless steel) and carefully measure the amount of powder to be used. Mix the powder and liquid thoroughly with a spatula and leave in a closed container for 1-2 minutes at 20-23°C.

**Application**

The waiting time in Table1. refers to an ambient temperature of 20-23°C. A higher temperature shortens and a lower temperature extends the setting time. The suitability of the material for application at the site of repair should be checked from time to time, especially when the ambient temperature is above 20-23°C. If the dough dries too quickly (high ambient temperature), moisten the surface of the dough applied to the denture with liquid.

**Polymerization**

The product belongs to the group of materials chemically curable in a water environment under pressure. The reaction is initiated as a result of the reaction of catalysts contained in the powder and the liquid.

The correct polymerization process parameters for the material are given in Table 1. They guarantee obtaining a repaired denture with appropriate parameters and ensure high patient safety related to the use of the repaired denture.

**Polymerization in an aqueous environment**: the use of pressure equipment with a minimum working pressure of 2 bar is required. The material should be polymerized in water according to the parameters given in Table 1.

**Residual monomer content**: the product does not contain MMA Table 2. Technical data. Nevertheless, traces of applied monomers may remain in the material after polymerization, which can be rinsed out of the denture and have harmful effects on the patient.

Observe the recommended mixing ratio, temperature and setting time of the acrylic dough and the parameters of the polymerization process provided for the specific material.

**Final treatment of the repaired denture**

Use the standard method of finishing. Polishing **with Pumeks** and **Pasta Polerska** by **Everall7.**

**Packaging and transportation of the repaired denture to the dentist's office**: the repaired denture should be protected against mechanical damage and high temperatures and handed over to the dentist.

**Hygiene, sterilization, disinfection**

Product

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

The prosthesis made of materials

There is a risk of transmission of infection. The repaired prosthesis should be cleaned and disinfected by the dentist before being replaced in the patient's mouth.

For the disinfection of dentures made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water.

Provide the appropriate information to the dentist.

**Additional Notes**

Keep the powder mixed with the liquid in a closed container.

The preparation time of the acrylic dough depends on the ambient temperature.

Moisten the repair site with liquid to ensure proper connection.

The denture repair site must be thoroughly polished. Scratches left in the denture are a place of dirt accumulation.

Do not overheat the material during the processing and polishing process – permanent changes in the material structure may occur, leading to an increase in the presence of a residual monomer in the repaired denture.

In case of any allergic reaction, stop working immediately and consult your dermatologist.

**Application parameters (1)**

| | Application parameters | Villacryl Repair Free |
|---|---|---|
| I | Powder / liquid mixing ratio | 10 g / 5 g (5.3 ml) liquid |
| II | Waiting time (23±2°C) | 1-2' |
| III | Polymerization Pressure method | 15-20' - 50-60°C 2 bar |

**Technical data (2)**

| | Technical Data: ISO 20795-1 | Villacryl Repaire Free |
|---|---|---|
| I | Classification | Type 2, Class 1 |
| II | Porosity | None |
| III | Bending strength | >60 MPa* |
| IV | Elastic modulus | >1500 MPa* |
| V | Sorption capacity | < 32 µg/mm3* |
| VI | Solubility | < 8.0 µg/mm3* |
| VII | Residual MMA Monomer | Without MMA |

* Requirements of PN-EN ISO 20795-1:2013-07E "Dentistry - Base plate polymers - Part 1: Polymers for denture plates".

**Warning**

**Information for dental technicians**

**Physical properties of products:** Highly flammable liquid and vapor. Keep away from heat/sparks/open flames/hot surfaces. - Smoking is not allowed.

**Skin contact:** Liquid, polymethyl methacrylate powder and dusts resulting from the processing of the product may be irritating to skin. Wear protective gloves/protective clothing/eye protection/face protection.

**Allergic reactions:** The liquid, benzoyl peroxide contained in the polymer, and the dust produced during processing and grinding of the repaired denture may cause an allergic skin reaction. IN CASE OF CONTACT WITH THE SKIN: Wash with plenty of soap and water.

**Airway contact:** Liquid, powder and the resulting dust can cause respiratory irritation. Avoid breathing in dust/smoke/gas/mist/vapor/sprayed liquid.

**Eye contact:** Accidental direct contact of the product components and the dust generated by the denture repair with the eyes may lead to their irritation. In the event of contact with the eyes, rinse with water. Contact with an ophthalmologist is recommended in case of irritation.

**Ingestion:** Do not eat. In case of accidental consumption of an ingredient, contact a doctor.

**Possibility of infection transmission:** the procedure being used enables the transmission of the infection. Apply the principles of good laboratory practice in technical and dental laboratories. Denture submitted for repair shall be washed and disinfected. Remove possible deposits of tartar and bacterial plaque from the surface of the repaired denture. The denture should be disinfected after repair before it is handed over to the patient.

**Safety measures**

**Information for dental technicians**

Do not use the product if you are diagnosed with immunodeficiency concerning allergens.

Flammable liquid (contains ethyl methacrylate). Keep away from sources of fire and radiation.

Do not eat or smoke when working with the material.

Standard protective measures used in technical and dental laboratories: protective masks, gloves and aprons are sufficient to safely handle the product.

In addition, it is advisable to use forced ventilation of the rooms and exhaust systems in order to reduce the risk of harmful effects of chemical substances contained in the products on the technician's body.

In case of a contact sensitization, stop using the device and consult an allergist. In case of anaphylactic shock, take immediate rescue action.

Do not use the products in patients with diagnosed immunodeficiency concerning allergens, including methyl methacrylate. In such situations, use other material that is appropriate for the patient.

Avoid inhalation of liquid vapor, direct contact of product ingredients and moist acrylic dough with skin and eyes.

In case of accidental consumption of an ingredient in the product, contact a specialist or poison center.

Do not pour liquids into the sewers.

**Patient**

Patients with diagnosed immunodeficiency concerning allergens: using the product may cause various types of sensitization reactions including anaphylactic shock. In such situations, use other material that is appropriate for the patient.

**Information to be given to patients by the dentist**

Odorless soap, various preparations for cleaning and care of dentures and ultrasonic cleaners can be used to maintain the hygiene of the denture repaired with the material.

The material is not resistant to the long-term effects of alcoholic beverages with a concentration above 40%.

The denture repair site may change color or lose its original physical properties in the event of contact with chemicals or misuse.

In case of any allergic reaction to the material of the repaired denture, immediately discontinue use and consult a dentist.

Report any repeated damage to the denture during the period of use to the dentist. Repairs may only be carried out by a dental technician.

If a piece of the damaged denture is swallowed, see a doctor immediately. The material is not visible on X-rays.

**Storing and handling the product**

The expiry date refers to unopened properly stored packaging.

The products are intended for repeated use.

The liquid is sensitive to temperature and UV radiation, under the influence of which an uncontrolled polymerization reaction may occur. Each time before using the product again, check that the liquid retains its properties. If there are visible signs of a change in density, opacification or change in color, stop using the product.

Store product components (powder and liquid) in original packaging, in a dark and well-ventilated room, at a temperature of 5-25°C in a place inaccessible to third parties.

Do not use the product after the expiry date indicated on the marketing packaging and component labels.

**Waste treatment**

Dispose of in accordance with current waste regulations. Cardboard packaging, labels, instructions – recyclable, disposable as non-hazardous municipal waste. Powder components and powder packaging are not hazardous waste. Dispose of as non-hazardous municipal waste. In case of packaging, it is possible to be recycled. Liquid ingredients and liquid packaging – hazardous waste. Hand it over to an authorized waste collection entity.

**Important points**

**Material to be used only in technical and dental laboratories**

Information provided in any way, even during a demonstration, does not invalidate the instructions for use. Users are required to check that the product is suitable for the intended use. The manufacturer shall not be liable for damage, including damage to third parties, resulting from failure to observe the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied.

Report all serious incidents involving the medical device to the manufacturer and the relevant authorities.

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everall7

# VILLACRYL BASE FREE

Heat-curing acrylic resin
for denture bases

EN | PL | RU |

Last Update: 2023-02

**Everall7 Sp. z o.o.**
ul. Augustówka 14 | 02-981 Warszawa; Polska
T +48 22 858 82 72 | info@everall7.pl
www.everall7.pl

CE

**Heat-curing acrylic resin for denture bases**
**Product Characteristics**
**Villacryl Base Free:** Heat-curing acrylic resin for denture bases. Material in powder and liquid form for making removable dentures. Without MMA (methyl methacrylate).
**Classification according to ISO 20795-1: Type 1, Class 1**
**Composition**
**Villacryl Base Free**
Powder: polymethyl methacrylate (PMMA), titanium dioxide, pigments, silica, viscose fi
Liquid: ethyl methacrylate (EMA), ethylene glycol dimethacrylate (EGDMA), hydroquino
*Depending on the color variety, the composition of the product may vary slightly in ter y and quality.
**Operating parameters:** The technical data necessary for proper processing of materia Table 1.
**Intended use of the product**
**Indications**
**Villacryl Base Free:** complete dentures (edentulous); partial dentures; relining of remo in dental lab.
**Contraindications**
Do not use the product to persons allergic to ethyl methacrylate as a substitute for trad based denture base materials. Acrylic monomers and benzoyl peroxide contained in the polymer are stances. They can be rinsed out of the denture in the oral environment and affect the patient's body. e products to patients with diagnosed immunodeficiency concerning allergens. If you notice any all ns, stop using the device and consult your doctor.
**Target group(s) of patients and intended users of the devices**
**Users**: Dental technicians making dentures.
**Patients**: Persons of different ages with partial or complete lack of teeth, eligible fo eatment with removable dentures.
**Side effects**
Irritation, redness or signs of hypersensitivity, including the occurrence of anaphylactic shock, may occur in the case of allergy to any of the components of the products. Oral irritation and inflammation caused by fungi and bacteria adhering to dentures made of materials are possible in the absence of proper hygiene.
**Clinical benefits**
Restoration of lost chewing and biting functions as a result of partial or total loss of teeth, resulting from irreversible conditions of the oral cavity, natural teeth, or from unfortunate random events. Improving the comfort and well-being of people who have lost their natural teeth. Improving the comfort and well-being of denture users who have lost the adhesion to the palate.
**Range***
**Color shades**

V2 – milk pink veined
V3 – dark pink veined
V4 – pink veined
T4 – pink
0 - transparent

**Packaging**

| | Powder + liquid set | Powder | Liquid |
|---|---|---|---|
| **Villacryl Base Free** | 750 g + 400 ml<br>300 g + 150 ml | 750 g<br>300 g<br>4 kg<br>2 kg | 400 ml<br>150 ml<br>1 l |

*Please contact your distributor to obtain information about the availability of the range.
**Method of use (Step by step instructions)**
**Initial operations**
Read the instructions carefully before using the product. Before making the denture, make sure that the patient for whom the denture is made does not have a diagnosed immunodeficiency concerning allergenic factors. It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering face mask suitable for liquids for respiratory protection, protective glasses and work clothing.
**Preparing**
**Tooth preparation**
Use acrylic or composite teeth. The teeth should have the mucous surface ground during setting. You can use porcelain teeth.
**Preparation of the plaster molds**
In order to obtain dental stone model in polymerization flask, use dental stone of II class (**Stodent II**), dental stone III class to be used for working models (especially **Stodent III**). After scalding out, the wet dental stone model should be covered with **Izo-Sol** isolating liquid **twice**. Wait a minute until the isolating liquid forms a shiny film on the model.
Clean the mucosal surface of acrylic teeth very thoroughly from the insulator and moisten it with a monomer. In special cases, follow the recommendations of the tooth manufacturer.
**Material preparation**
Shake the powder and liquid containers 2 or 3 times as shown on the packaging and product labels to facilitate powder dispersion in the monomer.
**Measuring and mixing of ingredients**
**Recommended mixing ratio**: (Table 1)
The stated mixing ratio must be strictly observed, especially for very thick lower dentures. Failure to comply with these values may lead to the formation of a porous denture structure, a decrease in its mechanical properties and an increase in the residual monomer content of ethyl methacrylate in the finished denture.
**Component mixing**
Pour an appropriate amount of liquid into a chemically resistant dish (ceramic, porcelain, glass or stainless steel) and carefully measure the amount of powder to be used. Mix the powder and liquid thoroughly with a spatula and leave in a closed container.
**Dough setting and kneading**
Check from time to time whether the material is ready for kneading (it does not stick to gloves and the walls of the container). The dough setting time is given in Table 1. The setting time for the material refers to an ambient temperature of 20-23 °C. A higher temperature shortens and a lower temperature extends the setting time.
If the material does not stick, it should be thoroughly kneaded for 1-2 minutes. Special attention should be paid to the very thorough kneading of the dough in its last phase, especially in cases of thick lower dentures. When kneading the acrylic dough, the hands (gloves) must not be moist and must not be lubricated with foreign substances such as cosmetic creams, talc, etc.
**Pressing**
Place the prepared material in a polymerization flask and press slowly until a pressure of 5000 kg is achieved. Leave under load for at least 15 minutes. In the case of thick lower dentures, extend the pressing time under load to 30-40 minutes.

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**Polymerization**
The product belongs to the group of materials thermally cured in an aqueous environment. For polymerization, use standard equipment for heat-curing polymerization. The correct polymerization process parameters for the material are given in Table 1. They guarantee obtaining a denture with appropriate parameters and ensure high patient safety related to the use of the denture. In the case of thick dentures use extended polymerization time.
**Cooling and removal from the flask**
The way the material is cooled after polymerization influences the dimensional stability of the denture and its adhesion to the palate. The correct post-polymerization cooling process parameters for the material are given in Table 1.
**Standard method polymerization**: Remove the flask in frame from the polymerizer and cool to room temperature. Open a completely cooled polymerization flask and remove the denture from the plaster.
**Polymerization extended over time**: Leave the flask in frame after polymerization in the polymerizer to cool down spontaneously to ambient temperature. Remove from water, carefully open and free the denture from plaster.
**Residual monomer content**: The product does not contain MMA. (Table 2. Technical data). Nevertheless, traces of applied monomers may remain in the denture after polymerization, which can be rinsed out of the denture and have harmful effects on the patient. Observe the recommended mixing ratio, temperature and setting time of the acrylic dough and the parameters of the polymerization process. To further reduce the residual monomer content, the finished denture should be soaked in water at room temperature for 48 hours before being given to the dentist.
**Final treatment of the denture**
Use the standard method of treatment. Polishing **with Pumeks** and **Pasta Polerska by Everall7.**
**Corrections and repairs of damaged dentures**
Corrections and repairs of damaged dentures can be performed with the **Villacryl Repair Free** acrylic material after the corrected or repaired element has been moistened with liquid.
**Packaging and transportation of the denture to the dentist's office**: Protect the finished denture from mechanical damage and heat and hand it over to the dentist.
**Hygiene, sterilization, disinfection**
Product
The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.
The prosthesis made of materials.
There is a risk of transmission of infection. The prosthesis should be cleaned and disinfected by the dentist before being placed in the patient's mouth for the first time.
For the disinfection of dentures made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.
It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water. Provide the appropriate information to the dentist.
**Additional Notes**
The dentures must be thoroughly polished. Scratches left in the denture are a place of dirt accumulation. Do not overheat the material during the processing and polishing process – permanent changes in the material structure may occur, leading to an increase in the presence of a residual monomer in the finished denture. In case of any allergic reaction, stop working immediately and consult your dermatologist.

**Application parameters (1)**

| | Application parameters | Villacryl Base Free |
|---|---|---|
| I | Powder / liquid mixing ratio | 23 g / 10 g (10, 6 ml) |
| II | Dough preparation time (23±2°C) | 20 - 25' |
| III | Shelf life of the dough | 30' |
| IV | Polymerization<br>Standard method<br>Extended method | 50'- 100°C<br>Cold water → 100°C |
| V | Cooling<br>Standard method<br>Extended method | 30' → 25°C<br>Leave in the polymerizer to cool down spontaneously to ambient temperature. |

**Technical data (2)**

| | Technical Data: ISO 20795-1 | Villacryl Base Free |
|---|---|---|
| I | Classification | Type 1, Class 1 |
| II | Porosity | None |
| III | Bending strength | > 65 MPa * |
| IV | Elastic modulus | >2000 MPa* |
| V | Sorption capacity | < 32 µg/mm3* |
| VI | Solubility | < 1.6 µg/mm3* |
| VII | Residual MMA Monomer | 0% |

* Requirements of PN-EN ISO 20795-1:2013-07E "Dentistry - Base plate polymers - Part 1: Polymers for denture plates."

**Warning:**
**Information for dental technicians:**
**Physical properties of products**: Highly flammable liquid and vapor. Keep away from heat/sparks/open flames/hot surfaces. - Smoking is not allowed.
**Skin contact:** Liquid, polymethyl methacrylate powder and dusts resulting from the processing of the product may be irritating to skin. Wear protective gloves/protective clothing/eye protection/face protection.
**Allergic reactions:** The liquid, benzoyl peroxide contained in the polymer, and the dust produced during processing and grinding of the denture may cause an allergic skin reaction. IN CASE OF CONTACT WITH THE SKIN: Wash with plenty of soap and water.
**Airway contact:** Liquid, powder and the resulting dust can cause respiratory irritation. Avoid breathing in dust/smoke/gas/mist/vapor/sprayed liquid.
**Eye contact:** Accidental direct contact of the product components and the dust generated by the denture with the eyes may lead to their irritation. In the event of contact with the eyes, rinse with water. Contact with an ophthalmologist is recommended in case of irritation.
**Ingestion:** Do not eat. In case of accidental consumption of an ingredient, contact a doctor.
**Possibility of infection transmission:** The procedure being used enables the transmission of the infection. Apply the principles of good laboratory practice in technical and dental laboratories. The denture should be washed and disinfected before being placed in the patient's mouth.
**Safety measures**
**Dental technicians**
Do not use the product if you are diagnosed with immunodeficiency concerning allergens. Flammable liquid (contains ethyl methacrylate). Keep away from sources of fire and radiation. Do not eat or smoke when working with the material. Standard protective measures used in technical and dental laboratories: protective masks, gloves and aprons are sufficient to safely handle the product. In addition, it is advisable to use forced ventilation of the rooms and exhaust systems in order to reduce the risk of harmful effects of chemical substances contained in the products on the technician's body. In case of a contact sensitization, stop using the device and consult an allergist. In case of anaphylactic shock, take immediate rescue action. Do not use the products in patients with diagnosed immunodeficiency concerning allergens, including methyl methacrylate. In such situations, use other material that is appropriate for the patient. Avoid inhalation of liquid vapor, direct contact of product ingredients and moist acrylic dough with skin and eyes. In case of accidental consumption of an ingredient in the product, contact a specialist or poison center. Do not pour liquids into the sewers.
**Patient**
Patients with diagnosed immunodeficiency concerning allergens: using the product may cause various types of sensitization reactions including anaphylactic shock. In such situations, use other material that is appropriate for the patient.
**Information to be given to patients by the dentist**
Wash the denture after each meal with a soft brush under cool, running water. Odorless soap, various preparations for cleaning and care of dentures and ultrasonic cleaners can be used for cleaning. It is recommended to wash the denture from time to time in a vibrating washer filled with a special denture cleaner. Wash unused dentures and store them in a sealed container filled with clean, cool water. Do not leave the denture in the sun or near heat sources. It is forbidden to place the denture in an open fire – flammable material. The denture is not resistant to the long-term effects of alcoholic beverages with a concentration above 40%. The denture may change color or lose its original physical properties in the event of contact with chemicals or misuse. In case of any allergic reaction to the denture material, consult your dentist immediately and stop using it. During use, the patient should have the denture periodically checked by a dentist. Report any damage to the denture during the period of use to the dentist. Repairs may only be carried out by a dental technician. If a piece of the damaged denture is swallowed, see a doctor immediately. The material is not visible on X-rays. It is recommended that the denture plate be lined or rebased in case of alveolar atrophy.
**Storing and handling the product**
The expiry date refers to unopened properly stored packaging. The product is designed for repeated use. The liquid is sensitive to temperature and UV radiation, under the influence of which an uncontrolled polymerization reaction may occur. Each time before using the product again, check that the liquid retains its properties. If there are visible signs of a change in density, opacification or change in color, stop using the product. Store product components (powder and liquid) in original packaging, in a dark and well-ventilated room, at a temperature of 5-25 °C in a place inaccessible to third parties. Do not use the product after the expiry date indicated on the marketing packaging and component labels.
**Waste treatment**
Dispose of in accordance with current waste regulations. Cardboard packaging, labels, instructions – recyclable, disposable as non-hazardous municipal waste. Powder components and powder packaging are not hazardous waste. Dispose of as non-hazardous municipal waste. In case of packaging, it is possible to be recycled. Liquid ingredients and liquid packaging – hazardous waste. Hand it over to an authorized waste collection entity.
**Important points**
**Material to be used only in technical and dental laboratories**
Information provided in any way, even during a demonstration, does not invalidate the instructions for use. Users are required to check that the product is suitable for the intended use. The manufacturer shall not be liable for damage, including damage to third parties, resulting from failure to observe the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied. Report all serious incidents involving the medical device to the manufacturer and the relevant authorities.

Please send any comments concerning our products to our address.

MD This symbol means a medical device as defined in the EU Regulation 2017/745.

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