Box size: 2.25 x 5
Overall size: 18 x 26
Folds to: 2.25 x 4.33

COVIDIEN™

Maxon™ and Maxon™ CV

Monofilament Absorbable Suture

PTD0105321 ISOP0105321

en

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

DESCRIPTION

Maxon™ and Maxon™ CV monofilament polyglycolate synthetic absorbable sutures are prepared from a copolymer of glycolide and trimethylene carbonate. Maxon™ and Maxon™ CV synthetic absorbable sutures are dyed with D&C Green No. 6. Maxon™ sutures are also available in clear (undyed). Maxon™ and Maxon™ CV synthetic absorbable sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP) for synthetic absorbable surgical sutures, except for minor variations in suture diameter.

Such variations are:

Maximum Suture Oversize in Diameter (mm) from USP and EP			
USP/EP Size	Designation in Tadlip (EP mm)	Dejacement maximal (mm)	
7-0/0.5	0.050-0.069	0.025	
6-0/0.7	0.070-0.099	0.025	
5-0/1	0.100-0.149	0.035	
4-0/1.5	0.15-0.199	0.040	
3-0/2	0.20-0.249	0.055	
2-0/3	0.30-0.339	0.055	
0/3.5	0.35-0.399	0.070	
1/4	0.40-0.499	0.070	

INDICATIONS

Maxon™ and Maxon™ CV sutures are indicated for use as absorbable sutures in general soft tissue approximation and/or ligation, in pediatric cardiovascular tissue, where growth is expected to occur, and in peripheral vascular tissue.

Maxon™ and Maxon™ CV sutures are not indicated for use as adult cardiovascular tissue, ophthalmic surgery, neurosurgery, and microsurgery.

ACTIONS

Maxon™ and Maxon™ CV synthetic absorbable sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the sutures by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Maxon™ and Maxon™ CV synthetic absorbable sutures occur by nonenzymatic hydrolysis. Absorption is complete when subsequently absorbed and excreted by the body.

Studies indicate that tensile strength for Maxon™ and Maxon™ CV sutures is approximately 80% of the original suture strength at one week post-implant. Approximately 75% of the original suture strength remains at two weeks, 65% at three weeks, 50% at four weeks and 25% at six weeks post-implant. Absorption of Maxon™ and Maxon™ CV synthetic absorbable sutures is minimal until about the 90th day after implantation. Absorption is substantially complete within six months.

CONTRAINDICATIONS

1. The use of this suture is contraindicated in patients with known sensitivities or allergies to its components.
2. This suture, being absorbable, should not be used where extended approximation of tissue is required or for permanent cardiovascular prostheses or synthetic grafts.

WARNINGS

1. Rinse, reprocessing and/or reintroduction of this device may create the risk of contamination, patient infection, permanent impairment or life-threatening injury. Do not rinse, reprocess or reintroduce this device. Use a new, sterile device for each procedure. Do not use if package is opened or damaged. Store suture in cool, dry place at room temperature. Avoid prolonged exposure to elevated temperature.
2. In surgery of the urinary and biliary tracts, care should be taken to avoid prolonged contact of this, or any other, suture with salt solutions, as catalysis formation may result.

PRECAUTIONS

1. Under sterile technique with a surgical procedure and techniques involving absorbable sutures before employing Maxon™ and Maxon™ CV absorbable sutures for wound closure, or the risk of wound dehiscence may vary with the site of application and the human reaction.
2. Acceptable surgical practice must be followed with respect to drainage and closure of contaminated or infected wounds.
3. The use of this suture may be inappropriate for patients with any condition which, in the opinion of the surgeon, may cause or contribute to delayed wound healing.
4. The use of supplemental absorbable sutures should be considered by the surgeon in the closure of sites subject to expansion, stretching or distention, or requiring additional support.

PRECAUTIONS

1. Under some circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon.
2. In handling this or any other suture material care should be taken to avoid damage from handling. Avoid crushing or overtightening due to application of surgical instruments such as forceps or needle holders.
3. Adequate tissue security requires the accepted surgical technique of flat, square knot, with additional throws as warranted for surgical circumstances and the experience of the surgeon.
4. The use of additional throws may be particularly appropriate when using monofilament sutures.
5. Consideration should be given to selecting sutures used in continuous mass closure that will provide adequate wound support for the anticipated healing period of the patient.
6. Sterile packaging is not guaranteed if the package is damaged.

ADVERSE REACTIONS

1. Tissue drag is a common complication encountered in the use of absorbable sutures.
2. Any of the following complications may be encountered with the use of absorbable sutures: local tissue reaction, infection, extrusion, tissue necrosis, and/or sinus formation.

HOW SUPPLIED

Maxon™ sutures are available in 0.9 and 0.7 cm (1.9 Metric) through 7-0 (0.1 Metric). Maxon™ CV sutures are available in 0.9 cm (0.5 Metric) through 6-0 (0.7 Metric). Maxon™ and Maxon™ CV sutures are dyed with D&C Green No. 6. Maxon™ sutures are also available in clear (Undyed). The sutures are supplied sterile, in a single-unit foil package, packaged in an offset to various needle types, using both permanent and removable needle attachment techniques. The sutures are available in box quantities of one, two and three dozen.

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Surgipro™ and Surgipro™ II*

Monofilament Polypropylene

*Referenced as Surgipro™ polypropylene sutures throughout the Instructions for Use.

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DESCRIPTION

Surgipro™ monofilament sutures (clear or pigmented) are inert, nonabsorbable, sterile sutures composed of an isotactic crystalline stereoisomer of polypropylene (a synthetic linear polyolefin) and polyethylene. The suture is pigmented blue to enhance visibility.

Surgipro™ polypropylene sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopeia (EP) for nonabsorbable surgical sutures, except for size 12 of 2 and 8-0, which differ from USP and EP maximum diameter requirements by up to 0.007 mm.

INDICATIONS

Surgipro™ polypropylene sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, and ophthalmic tissues.

ACTIONS

Surgipro™ polypropylene sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Surgipro™ polypropylene sutures are not absorbed, nor is any significant change in strength retention known to occur in vivo.

CONTRAINDICATIONS

None known.

WARNINGS

1. As with any foreign body, prolonged contact of any suture with salt solutions may result in calculus formation.
2. Reuse, reprocessing and/or resterilization of this device may create the risk of contamination, patient infection, permanent impairment or life-threatening injury. Do not reuse, reprocess or resterilize this device. Sterile, unless packaging has been opened or damaged. Discard open, unused sutures. Dispose of used instruments in accordance with the end-user's medical and hospital waste disposal requirements. Store at room temperature. Avoid prolonged exposure to elevated temperatures.
3. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Surgipro™ polypropylene sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
4. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
2. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS

Adverse effects associated with the use of this product, include: wound dehiscence, calculi formation when prolonged contact with salt solutions occurs, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED

Surgipro™ polypropylene sutures are available either undyed (clear) or dyed with copper phthalocyanine blue. Surgipro™ sutures are available in sizes 2 (5 Metric) through 8-0 (0.4 metric) and 10-0 (0.2 metric).

The sutures are supplied sterile; affixed to various types of needles, using both permanent and removable needle attachment techniques. The sutures are available in box quantities of one, two and three dozen.

Surgipro™ polypropylene sutures are also available with PTFE (polytetrafluoroethylene) pledgets for use as a pad between the suture area and the tissue surface to increase the load bearing area. Surgipro™ polypropylene sutures are available in presentations containing bead and collar components to anchor the end of the suture for subcuticular closure or for use as tendon sutures.

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COVIDIEN™

PTFE Polymer Pledgets

Nonabsorbable Pre-punched Pledgets

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IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reuse, reprocessing and/or re sterilization of this device may create the risk of contamination, patient infection, permanent impairment or life threatening injury. Do not reuse, reprocess or re-sterilize this device.

DESCRIPTION

PTFE Pledgets are nonabsorbable, undyed fabrics composed of poly (tetrafluoroethylene).

INDICATIONS

PTFE Pledgets are indicated for use as nonabsorbable suture buttresses when there is a possibility of sutures tearing through tissue.

ACTIONS

PTFE Pledgets elicit a minimal acute inflammatory reaction in tissue, followed by gradual encapsulation of the fabric by fibrous connective tissue. PTFE Pledgets are not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.

CONTRAINDICATIONS

The use of this pledget is contraindicated in patients with known sensitivities or allergies to its components.

WARNINGS

1. Sterile unless packaging has been opened or damaged. Discard open unused pledgets.

2. Radiation should not be used on these materials.

3. Users should be trained in surgical procedures and techniques involving suture pledgets before employing PTFE Pledgets.

4. Acceptable surgical practice must be followed with respect to the management of contaminated or infected wounds.

PRECAUTIONS

1. In handling this or any other pledget material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

2. Thread suture through the prepunched holes in order to ensure adequate margins and intervening bridge.

3. Dispose of contaminated device and packaging materials utilizing standard hospital procedures and universal precautions for biohazardous waste.

4. For use with non-absorbable suture size 2 USP (metric size 5) or smaller.

ADVERSE REACTIONS

Adverse reactions and potential complications associated with bulk pledgets include leak, separation from suture material, obstruction, extrusion, infection/abscess formation, erosion, migration, fistula formation, encapsulation/granuloma, allergic reaction, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED

PTFE Pledgets are available sterile, prepunched, square, oval or rectangular in shape and in box quantities of 3 dozen.

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DESCRIPTION

Ti-Cron™ braided polyester sutures are nonabsorbable, sterile, surgical sutures composed of Polyethylene Terephthalate. They are prepared from fibers of high molecular weight, long chain, linear polyesters having recurrent aromatic rings as an integral component. Ti-Cron™ braided sutures are available coated uniformly with silicone to increase surface lubricity, thereby enhancing handling characteristics, ease of passage through tissue, and knot run-down properties. Ti-Cron™ sutures are inert and elicit only minimal local tissue reaction. Ti-Cron™ sutures are available dyed with D&C Blue No. 6 colorant or undyed (white).

Ti-Cron™ sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopeia (EP) for nonabsorbable surgical sutures.

INDICATIONS

Ti-Cron™ polyester sutures are indicated for use in general soft tissue approximation and/or ligation.

ACTIONS

Ti-Cron™ polyester sutures elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Ti-Cron™ polyester sutures are not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.

CONTRAINDICATIONS

None known.

WARNINGS

1. As with any foreign body, prolonged contact of any suture with salt solutions may result in calculus formation.
2. Reuse, reprocessing and/or resterilization of this device may create the risk of contamination, patient infection, permanent impairment or life threatening injury. Do not reuse, reprocess or resterilize this device. Sterile unless packaging has been opened or damaged. Discard open, unused sutures. Store at room temperature. Avoid prolonged exposure to elevated temperatures.
3. Users should be familiar with surgical procedures and techniques involving sutures before employing Ti-Cron™ polyester nonabsorbable sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
4. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
2. Adequate knot security requires the acceptable surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Dispose of used instruments in accordance with the end-user's medical and biological waste disposal requirements.

ADVERSE REACTIONS

Adverse effects associated with the use of this product are: wound dehiscence, calculi formation when prolonged contact with salt solutions occurs, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED

Ti-Cron™ polyester sutures are available undyed (white) or dyed with D&C Blue No. 6 and braided, coated, in sizes 5 (7 Metric) though 7-0 (0.5 Metric). The sutures are supplied sterile, in pre-cut lengths and ligating reels, non-needled or affixed to various needle types using both permanent and removable needle attachment techniques. The sutures are available in box quantities of one, two and three dozen. Ti-Cron™ polyester sutures are also available with PTFE (polytetrafluoroethylene) pledgets for use as a pad between the suture area and the tissue surface to increase the load-bearing area.

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STERILE EO

Rx ONLY

Caution, consult accompanying documents

CE 0123

© 2011 Covidien.

Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore, Ireland.

www.covidien.com

+1 800 633 8766

+1 763 514 4000

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, ™* brands are trademarks of their respective owner.

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