Box size: 2.25 x 5
Overall size: 18 x 12.5
Folds to: 2.25 x 4.1667

COVIDIEN™

Steel

Monofilament Stainless Steel

1300380

EN

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.
DESCRIPTION
Stainless monofilament steel sutures are composed of 316L stainless steel conforming to ASTM Standard F138 Grade 2 ("stainless steel" or "stainless steel for surgical implants").
Steel sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP) for nonabsorbable surgical sutures.
INDICATIONS
Steel sutures are indicated for use in abdominal wound closure, sternotomy assurances, hernia repair and sternal closure.
ACTIONS
Steel sutures elicit a minimal acute inflammatory reaction in tissue, and are not absorbed.
CONTRAINDICATIONS
The use of these sutures are contraindicated in patients with known sensitivities or allergies to steel and prior to prolonged exposure to chromium and nickel. Additionally, the presence of steel may interfere with certain radiodiagnostics and its use is contraindicated where radiotherapeutics of suture material is required.
WARNINGS
Do not resterilize. Sterile unless packaging has been opened or damaged. Discard open, unused sutures. Store at room temperature. Avoid prolonged exposure to elevated temperatures.
Care should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing steel sutures for wound closure, as the ends of steel sutures may be sharp. Use care when tying steel sutures and the suture material used.
Acceptable surgical practice must be followed with respect to drainage and closure of infected wounds or infected wounds.
PRECAUTIONS
Acceptable surgical practice should be followed to place steel sutures sufficiently distal to wound edges to reduce the incidence of the suture cutting through tissue when the wound is placed under extreme stress. Additionally, care should be taken to cut the "tail" of steel surgical knots to allow for the bending of that end away from surrounding tissue.
ADVERSE REACTIONS
Adverse effects, which may be associated with the use of this product, include: wound dehiscence; allergic response in patients with known sensitivities to steel or to prolonged exposure to chromium and nickel; prolonged contact with salt solutions occurs, enhanced bacterial infection), minimal acute inflammatory reaction, and transitory local irritation.
HOW SUPPLIED
Steel sutures are available in monofilament strands in sizes 3-0 (3/0metric) through 7 (4 Metric). The sutures are supplied sterile in pre-cut lengths, straight or attached to needles of various types of needles in various needle attachment techniques. The sutures are available in two quantities of one and three dozen.

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